

Technical Publications

Direction DOC2148033

Revision 1

OEC One Workstation Software Release WKS-2.x DICOM CONFORMANCE STATEMENT

Copyright© 2018 by General Electric Co.

Do not duplicate

GE Healthcare

REVISION HISTORY

| REV | VER | DATE | REASON FOR CHANGE |
|------------|------------|-------------|---------------------------------------------------|
| 1 | 1 | May 14 2018 | Draft for review |
| 1 | 2 | May 21 2018 | Update according to the review record from Roland |
| 1 | 3 | May 28 2018 | Update document by Support Center Review Comments |

THIS PAGE LEFT INTENTIONALLY BLANK

CONFORMANCE STATEMENT OVERVIEW

This DICOM Conformance Statement captures the DICOM capabilities of the OEC One Workstation identified below. The content is consistent with the 2015c version of the DICOM Standard.

Table 0.1 provides an overview of the network services supported by OEC One Workstation system.

Table 0.1 – NETWORK SERVICES

| SOP Classes | User of Service (SCU) | Provider of Service (SCP) |
|------------------------------------------------------|-----------------------|---------------------------|
| Transfer | | |
| Computed Radiography Image Storage | No | Yes |
| Secondary Capture Image Storage | No | Yes |
| X-Ray Angiographic Image Storage | Yes | Yes |
| X-Ray Radiofluoroscopic Image Storage | No | Yes |
| X-Ray Radiation Dose SR | Yes | No |
| Query/Retrieve | | |
| Patient Root Query/Retrieve Information Model – FIND | Yes | No |
| Patient Root Query/Retrieve Information Model – MOVE | Yes | No |
| Study Root Query/Retrieve Information Model – FIND | Yes | No |
| Study Root Query/Retrieve Information Model – MOVE | Yes | No |
| Print Management | | |
| Basic Film Session SOP Class | Yes | No |
| Basic Film Box SOP Class | Yes | No |
| Basic Grayscale Image Box SOP Class | Yes | No |
| Basic Grayscale Print Management Meta SOP Class | Yes | No |
| Printer SOP Class | Yes | No |
| Workflow Management | | |
| Storage Commitment Push Model SOP Class | Yes | No |
| Modality Performed Procedure Step SOP Class | Yes | No |
| Modality Worklist Information Model – FIND SOP Class | Yes | No |

Table 0.2 provides an overview of the Media Storage Application Profiles supported by Workstation system.

Table 0.2 - MEDIA SERVICES

| Media Storage Application Profile | Write Files (FSC or FSU) | Read Files (FSR) |
|-----------------------------------|--------------------------|------------------|
| USB | | |
| X-Ray Angiographic Image Storage | Yes | Yes |
| X-Ray Radiation Dose SR | Yes | No |

TABLE OF CONTENTS

| | |
|------------------------------------------------------------------------------|-----------|
| 1. INTRODUCTION | 11 |
| 1.1 Overview | 11 |
| 1.2 Overall DICOM Conformance Statement Document Structure..... | 12 |
| 1.3 Intended Audience | 13 |
| 1.4 Scope and Field of Application | 13 |
| 1.5 Important Remarks | 13 |
| 1.6 References..... | 14 |
| 1.7 Definitions..... | 14 |
| 1.8 Symbols and Abbreviations..... | 16 |
| 2. NETWORK CONFORMANCE STATEMENT | 19 |
| 2.1 Introduction..... | 19 |
| 2.2 Implementation Model..... | 19 |
| 2.2.1 Application Data Flow Diagram | 19 |
| 2.2.2 Functional Definition of AE's..... | 21 |
| 2.2.3 Sequencing of Real-World Activities..... | 22 |
| 2.3 AE Specifications | 22 |
| 2.3.1 Workstation AE Specification..... | 22 |
| 2.3.1.1 Association Establishment Policies..... | 23 |
| 2.3.1.1.1 General..... | 23 |
| 2.3.1.1.2 Number of Associations..... | 23 |
| 2.3.1.1.3 Asynchronous Nature..... | 23 |
| 2.3.1.1.4 Implementation Identifying Information..... | 23 |
| 2.3.1.2 Association Initiation Policy | 23 |
| 2.3.1.2.1 Real-World Activity Verify | 24 |
| 2.3.1.2.1.1 Associated Real-World Activity | 24 |
| 2.3.1.2.1.2 Proposed Presentation Context Table | 24 |
| 2.3.1.2.2 Real-World Activity Schedule Update..... | 24 |
| 2.3.1.2.2.1 Associated Real-World Activity | 24 |
| 2.3.1.2.2.2 Proposed Presentation Context Table | 24 |
| 2.3.1.2.3 Real-World Activity Start/End Exam | 26 |
| 2.3.1.2.3.1 Associated Real-World Activity | 26 |
| 2.3.1.2.3.2 Proposed Presentation Context Table | 26 |
| 2.3.1.2.4 Real-World Activity Image/RDSR Send and Storage Commit..... | 26 |
| 2.3.1.2.4.1 Associated Real-World Activity | 26 |
| 2.3.1.2.4.2 Proposed Presentation Context Table | 26 |
| 2.3.1.2.5 Real-World Activity Image Query/Retrieve and Local Store SCP | 29 |
| 2.3.1.2.5.1 Associated Real-World Activity | 29 |
| 2.3.1.2.5.2 Proposed Presentation Context Table | 29 |

| | | |
|-------------|--------------------------------------------------------------------|-----------|
| 2.3.1.2.6 | Real-World Activity “DICOM Print” | 31 |
| 2.3.1.2.6.1 | Associated Real-World Activity | 31 |
| 2.3.1.2.6.2 | Proposed Presentation Context Table | 31 |
| 2.3.1.2.7 | Real-World Activity Report MPPS..... | 39 |
| 2.3.1.2.7.1 | Associated Real-World Activity | 39 |
| 2.3.1.2.7.2 | Proposed Presentation Context Table | 39 |
| 2.3.1.3 | Association Acceptance Policy | 40 |
| 2.3.1.3.1 | Real-World Activity Verification SCP | 40 |
| 2.3.1.3.1.1 | Associated Real-World Activity | 40 |
| 2.3.1.3.1.2 | Accepted Presentation Context Table | 41 |
| 2.3.1.3.1.3 | Presentation Context Acceptance Criterion | 41 |
| 2.3.1.3.1.4 | Transfer Syntax Selection Policies | 41 |
| 2.3.1.3.2 | Real-World Activity Local Store SCP | 41 |
| 2.3.1.3.2.1 | Associated Real-World Activity | 41 |
| 2.3.1.3.2.2 | Accepted Presentation Context Table | 41 |
| 2.3.1.3.2.3 | Presentation Context Acceptance Criterion | 42 |
| 2.3.1.3.2.4 | Transfer Syntax Selection Policies | 42 |
| 2.3.1.3.3 | Real-World Activity- Storage Commit – SCU role | 42 |
| 2.3.1.3.3.1 | Associated Real-World Activity | 42 |
| 2.3.1.3.3.2 | Accepted Presentation Context Table | 42 |
| 2.3.1.3.3.3 | Presentation Context Acceptance Criterion | 42 |
| 2.3.1.3.3.4 | Transfer Syntax Selection Policies | 43 |
| 2.4 | Communication Profiles | 43 |
| 2.4.1 | Supported Communication Stacks | 43 |
| 2.4.2 | Physical Media Support | 43 |
| 2.4.3 | Additional Protocol | 43 |
| 2.4.4 | IPv4 and IPv6 Support | 43 |
| 2.5 | Extensions / Specializations/ Privatizations | 43 |
| 2.6 | Configuration | 43 |
| 2.6.1 | AE Title/Presentation Address Mapping..... | 43 |
| 2.6.2 | Configurable Parameters..... | 43 |
| 2.7 | Support of Extended Character Sets..... | 45 |
| 2.8 | Codes and Controlled Terminology..... | 45 |
| 2.8.1 | Fixed Coded Terminology | 45 |
| 2.8.2 | Mapped Coded Terminology | 46 |
| 2.9 | Security Profiles | 46 |
| 3. | MEDIA STORAGE CONFORMANCE STATEMENT | 47 |
| 3.1 | Introduction..... | 47 |
| 3.1.1 | Application Data Flow Diagram | 47 |
| 3.1.2 | Functional Definition of AE’s..... | 47 |
| 3.1.3 | Sequencing of Real-World Activities..... | 47 |
| 3.1.4 | File Meta Information Options (See PS3.10)..... | 48 |
| 3.2 | AE Specifications | 48 |
| 3.2.1 | Workstation AE Specification..... | 48 |
| 3.2.1.1 | File Meta Information for the Workstation Application Entity | 48 |
| 3.2.1.2 | Real-World Activities for the Workstation Application Entity | 48 |

| | | |
|-------------|-----------------------------------------------------------------------|-----------|
| 3.2.1.2.1 | Real-World Activity Storage | 48 |
| 3.2.1.2.1.1 | Media Storage Application Profile for the RWA Update File Set | 48 |
| 3.2.2 | Media Viewer AE Specification | 49 |
| 3.2.2.1 | File Meta Information for the Media Viewer Application Entity | 49 |
| 3.2.2.2 | Real-World Activities for the Media Viewer Application Entity | 49 |
| 3.2.2.2.1 | Real-World Activity Directory Listing | 49 |
| 3.2.2.2.1.1 | Media Storage Application Profile for the RWA Directory Listing | 50 |
| 3.2.3 | Media Viewer AE Specification | 50 |
| 3.3 | Augmented and Private Application Profiles | 50 |
| 3.4 | Extensions / Specializations / Privatizations | 50 |
| 3.5 | Configuration | 50 |
| 3.6 | Support of Extended Character Sets..... | 51 |
| 4. | X-RAY ANGIOGRAPHIC INFORMATION OBJECT IMPLEMENTATION | 52 |
| 4.1 | Introduction..... | 52 |
| 4.2 | Workstation Mapping of DICOM Entities..... | 52 |
| 4.3 | IOD Module Table | 52 |
| 4.4 | Information Module Definitions | 53 |
| 4.4.1 | Patient Entity Modules | 53 |
| 4.4.1.1 | Patient Module | 53 |
| 4.4.2 | Study Entity Modules..... | 55 |
| 4.4.2.1 | General Study Module | 55 |
| 4.4.2.2 | Patient Study Module | 56 |
| 4.4.3 | Series Entity Modules | 56 |
| 4.4.3.1 | General Series Module | 56 |
| 4.4.4 | Equipment Entity Modules..... | 58 |
| 4.4.4.1 | General Equipment Module | 58 |
| 4.4.5 | Image Entity Modules | 59 |
| 4.4.5.1 | General Image Module..... | 59 |
| 4.4.5.2 | Image Pixel Module | 60 |
| 4.4.5.3 | Contrast/Bolus Module | 61 |
| 4.4.5.4 | Cine | 61 |
| 4.4.5.5 | Multi-Frame | 62 |
| 4.4.5.6 | Mask..... | 62 |
| 4.4.5.7 | X-Ray Image Module..... | 62 |
| 4.4.5.7.1 | Image Type | 63 |
| 4.4.5.8 | X-Ray Acquisition Module | 64 |
| 4.4.5.9 | XA Positioner Module | 65 |
| 4.4.5.10 | Overlay Plane Module..... | 65 |
| 4.4.5.11 | VOI LUT module..... | 65 |
| 4.4.5.12 | SOP Common Module | 66 |
| 4.5 | Standard Extended and Private Data Attributes..... | 69 |
| 4.5.1 | Standard Attributes | 69 |
| 4.6 | Standard Extended and Private Context Groups | 69 |

5. X-RAY RADIATION DOSE STRUCTURED REPORT INFORMATION OBJECT IMPLEMENTATION70

- 5.1 Introduction..... 70**
- 5.2 Workstation Mapping of DICOM Entities 70**
- 5.3 IOD Module Table 70**
- 5.4 Information Module Definitions 71**
 - 5.4.1 Patient Entity Modules 71
 - 5.4.1.1 Patient Module 71
 - 5.4.2 Study Entity Modules..... 73
 - 5.4.2.1 General Study Module 73
 - 5.4.2.2 Patient Study Module 73
 - 5.4.3 Series Entity Modules 74
 - 5.4.3.1 SR Document Series Module 74
 - 5.4.4 Equipment Entity Modules..... 75
 - 5.4.4.1 General Equipment Module 75
 - 5.4.4.2 Enhanced General Equipment Module..... 75
 - 5.4.5 Document Entity Modules 75
 - 5.4.5.1 SR Document General Module 75
 - 5.4.5.2 SR Document Content Module 77
 - 5.4.5.2.1 SR Document Content Descriptions 78
 - 5.4.5.2.1.1 Content Template..... 78
 - 5.4.5.3 SOP Common Module 78
- 5.5 Standard Extended and Private Data Attributes 80**
- 5.6 Standard Extended and Private Context Groups 80**
 - 5.6.1 Standard Extended Context Groups 80
 - 5.6.2 Configurable Context Groups 80
- 5.7 Standard, Standard Extended and Private Templates 81**
 - 5.7.1 Standard Templates..... 81
 - 5.7.1.1 Template ID 10001 X-Ray Radiation Dose 81
 - 5.7.1.2 Template ID 10002 Accumulated X-Ray Dose..... 82
 - 5.7.1.3 Template ID 10003 Irradiation Event X-Ray Data 84
 - 5.7.1.4 TID 10003A Irradiation Event X-Ray Detector Data 87
 - 5.7.1.5 TID 10003B Irradiation Event X-Ray Source Data 87
 - 5.7.1.6 Template ID 10004 Accumulated Projection X-Ray Dose 90
 - 5.7.1.7 TID 10007 Accumulated Total Projection Radiography Dose 91
 - 5.7.1.8 Template ID 1002 Person Observer Context 92
 - 5.7.1.9 Template ID 1002 Device Observer Context 93
 - 5.7.1.10 TID 1003 Person Observer Identifying Attributes 93
 - 5.7.1.11 Template ID 1004 Device Observer Identifying Attributes 93
 - 5.7.2 Standard Extended Templates..... 94

6. MODALITY WORKLIST QUERY IMPLEMENTATION95

- 6.1 Introduction..... 95**
- 6.2 Workstation Mapping of DICOM Entities 95**
- 6.3 Worklist Query Module Table 95**

| | | |
|------------|----------------------------------------------------------------------|------------|
| 6.4 | Worklist Query Module Definitions | 96 |
| 6.4.1 | Common Scheduled Procedure Step Entity Modules..... | 96 |
| 6.4.1.1 | SOP Common Module | 96 |
| 6.4.1.1.1 | Specific Character Set..... | 96 |
| 6.4.1.2 | Scheduled Procedure Step Module | 97 |
| 6.4.1.2.1 | Scheduled Station AE Title..... | 99 |
| 6.4.2 | Common Requested Procedure Entity Modules..... | 99 |
| 6.4.2.1 | Requested Procedure Module | 99 |
| 6.4.2.1.1 | Study Instance UID..... | 101 |
| 6.4.3 | Common Imaging Service Request Entity Modules..... | 101 |
| 6.4.3.1 | Imaging Service Request Module | 101 |
| 6.4.4 | Common visit Entity Modules..... | 102 |
| 6.4.4.1 | Visit Identification | 102 |
| 6.4.4.2 | Visit Status | 102 |
| 6.4.4.3 | Visit Relationship..... | 102 |
| 6.4.4.4 | Visit Admission..... | 102 |
| 6.4.5 | Common Patient Entity Modules | 102 |
| 6.4.5.1 | Patient Relationship | 102 |
| 6.4.5.2 | Patient Identification | 102 |
| 6.4.5.3 | Patient Demographic | 103 |
| 6.4.5.4 | Patient Medical | 104 |
| | | |
| 7. | MODALITY PERFORMED PROCEDURE STEP IMPLEMENTATION | 105 |
| 7.1 | Introduction..... | 105 |
| 7.2 | Relationship between Scheduled and Performed Procedure Steps | 105 |
| 7.3 | Modality Performed Procedure Step Module Table..... | 105 |
| 7.4 | Modality Performed Procedure Step Module Definitions..... | 105 |
| 7.4.1 | SOP Common Module | 106 |
| 7.4.2 | Performed Procedure Step Relationship Module | 106 |
| 7.4.3 | Performed Procedure Step Information Module | 108 |
| 7.4.4 | Image Acquisition Results Module | 109 |
| 7.4.5 | Radiation Dose Module | 110 |
| 7.4.6 | Billing and Material Management Codes Module | 110 |
| 7.5 | Standard Extended and Private Data Attributes | 111 |
| 7.5.1 | Standard Attributes | 111 |
| 7.6 | Standard Extended and Private Context Groups | 111 |
| 7.6.1 | Standard Extended Context Groups | 111 |
| 7.6.1.1 | Context ID 9301 Modality PPS Discontinuation Reasons | 111 |
| | | |
| 8. | STORAGE COMMITMENT PUSH MODEL IMPLEMENTATION..... | 113 |
| 8.1 | Storage Commitment Push Model Information Object Definition..... | 113 |
| 8.1.1 | STORAGE COMMITMENT MODULE FOR N-ACTION | 113 |
| 8.1.2 | STORAGE COMMITMENT MODULE FOR N-EVENT-REPORT | 113 |
| 8.1.2.1 | Processing of Failure Reason when received in a N-Event-Report | 114 |
| | | |
| 9. | BASIC DIRECTORY INFORMATION OBJECT IMPLEMENTATION | 115 |

9.1 IOD Module Table 115

9.2 Information Module Definitions 115

9.2.1 File Set identification Module..... 116

9.2.2 Directory Information Module 116

9.2.3 Definition of Specific Directory Records..... 117

9.2.3.1 Patient Directory Record Definition..... 117

9.2.3.2 Study Directory Record Definition 117

9.2.3.3 Series Directory Record Definition..... 118

9.2.3.3.1 Icon Image Key Definition 118

9.2.3.4 Image Directory Record Definition..... 119

9.2.3.4.1 Icon Image Key Definition 119

9.2.3.5 SR Document Directory Record Definition 120

9.3 Private Data Dictionary..... 120

10. PRINT MANAGEMENT IMPLEMENTATION 121

10.1 Introduction..... 121

10.2 Basic Film Session SOP Class 121

10.2.1 Basic Film Session N-Create Attributes 121

10.2.2 Basic Film Session N-Delete..... 121

10.3 Basic Film Box SOP Class 121

10.3.1 Basic Film Box N-Create Attributes 121

10.3.2 Basic Film Box N-Action Attributes 123

10.3.3 Basic Film Box N-Delete 123

10.4 Basic Grayscale Image Box SOP Class 123

10.4.1 Basic Grayscale Image Box Pixel N-Set Attributes 123

10.5 Printer SOP Class 124

10.5.1 Printer N-Event-Report Attributes 124

10.5.2 Printer N-Get Attributes 125

11. QUERY IMPLEMENTATION 126

11.1 Workstation Mapping of DICOM Entities 126

11.2 Information Model Keys 126

11.2.1 Common Query Keys..... 126

11.2.1.1 Q/R Common Attribute Descriptions 126

11.2.1.1.1 Specific Character Set..... 126

11.2.2 Patient Level 127

11.2.3 Study Level – Patient Root 128

11.2.4 Study Level – Study Root 129

11.2.5 Series Level..... 131

11.2.6 Image Level..... 132

1. INTRODUCTION

1.1 OVERVIEW

This DICOM Conformance Statement is divided into Sections as described below:

Section 1 (Introduction), which describes the overall structure, intent, and references for this Conformance Statement

Section 2 (Network Conformance Statement), which specifies the GEHC equipment compliance to the DICOM requirements for the implementation of Networking features.

Section 3 (Media Storage Conformance Statement), specifies the Workstation equipment compliance to the DICOM requirements for the implementation of Media Storage features.

Section 4 (X-Ray Angiographic Information Object Implementations), specifies the Workstation Medical Systems equipment compliance to DICOM requirements for the implementation of an X-Ray Information Objects.

Section 5 (X-Ray Radiation Dose SR Information Object Implementations), specifies the Workstation Medical Systems equipment compliance to DICOM requirements for the implementation of X-Ray Radiation Dose SR Information Object Implementation feature.

Section 6 (Modality Worklist Query Information Model), specifies the Workstation Medical Systems equipment compliance to DICOM requirements for the implementation of the Modality Worklist service.

Section 7 (Modality Performed Procedure Step Implementation), specifies the Workstation Medical Systems equipment compliance to DICOM requirements for the implementation of the Modality Performed Procedure Step SOP Classes.

Section 8 (Storage Commitment Push Model Implementation), specifies the Workstation Medical Systems equipment compliance to DICOM requirements for the implementation of the Storage Commitment Push feature.


Section 9 (Basic Directory Information Object Implementation), specifies the Workstation Medical Systems equipment compliance to DICOM requirements for the implementation of a Basic Directory Information Object.

Section 10 (Print Management Implementation), specifies the Workstation Medical Systems equipment compliance to DICOM requirements for the implementation of the Basic Print Meta SOP Classes (Grayscale).

Section 11 (Query/Retrieve Information Model), specifies the Workstation Medical Systems equipment compliance to DICOM requirements for the implementation of the Patient, Study Root Query/Retrieve service.

1.2 OVERALL DICOM CONFORMANCE STATEMENT DOCUMENT STRUCTURE

The Documentation Structure of the GEHC DICOM Conformance Statements is shown in the Illustration below.



This document specifies the DICOM implementation. It is entitled:

OEC One Workstation Version WKS-2.x
Conformance Statement for DICOM
1941526

This DICOM Conformance Statement documents the DICOM Conformance Statement and Technical Specification required interoperating with the GEHC network interface.

The GEHC Conformance Statement, contained in this document, also specifies the Lower Layer communications which it supports (e.g., TCP/IP). However, the Technical Specifications are defined in the DICOM Part 8 standard.

For more information regarding DICOM, copies of the Standard may be obtained on the Internet at <http://medical.nema.org>. Comments on the Standard may be addressed to:

DICOM Secretariat
NEMA
1300 N. 17th Street, Suite 1752
Rosslyn, VA 22209
USA
Phone: +1.703.841.3200

1.3 INTENDED AUDIENCE

The reader of this document is concerned with software design and/or system integration issues. It is assumed that the reader of this document is familiar with the DICOM Standard and with the terminology and concepts which are used in that Standard.

1.4 SCOPE AND FIELD OF APPLICATION

It is the intent of this document to provide an unambiguous specification for GEHC implementations. This specification, called a Conformance Statement, includes a DICOM Conformance Statement and is necessary to ensure proper processing and interpretation of GEHC medical data exchanged using DICOM. The GEHC Conformance Statements are available to the public.

The reader of this DICOM Conformance Statement should be aware that different GEHC devices are capable of using different Information Object Definitions. For example, a GEHC CT Scanner may send images using the CT Information Object, MR Information Object, Secondary Capture Object, etc.

Included in this DICOM Conformance Statement are the Module Definitions which define all data elements used by this GEHC implementation. If the user encounters unspecified private data elements while parsing a GEHC Data Set, the user is well advised to ignore those data elements (per the DICOM standard). Unspecified private data element information is subject to change without notice. If, however, the device is acting as a "full fidelity storage device", it should retain and re-transmit all of the private data elements which are sent by GEHC devices.

1.5 IMPORTANT REMARKS

The use of these DICOM Conformance Statements, in conjunction with the DICOM Standards, is intended to facilitate communication with GE imaging equipment. However,

by itself, it is not sufficient to ensure that inter-operation will be successful. The user (or user's agent) needs to proceed with caution and address at least four issues:

- **Integration** - The integration of any device into an overall system of interconnected devices goes beyond the scope of standards (DICOM v3.0), and of this introduction and associated DICOM Conformance Statements when interoperability with non-GE equipment is desired. The responsibility to analyze the applications requirements and to design a solution that integrates GE imaging equipment with non-GE systems is the **user's** responsibility and should not be underestimated. The **user** is strongly advised to ensure that such an integration analysis is correctly performed.
- **Validation** - Testing the complete range of possible interactions between any GE device and non-GE devices, before the connection is declared operational, should not be overlooked. Therefore, the **user** should ensure that any non-GE provider accepts full responsibility for all validation required for their connection with GE devices. This includes the accuracy of the image data once it has crossed the interface between the GE imaging equipment and the non-GE device and the stability of the image data for the intended applications.

Such a validation is required before any clinical use (diagnosis and/or treatment) is performed. It applies when images acquired on GE imaging equipment are processed/displayed on a non-GE device, as well as when images acquired on non-GE equipment is processed/displayed on a GE console or workstation.

- **Future Evolution** - GE understands that the DICOM Standard will evolve to meet the user's growing requirements. GE is actively involved in the development of the DICOM Standard. DICOM will incorporate new features and technologies and GE may follow the evolution of the Standard. The GEHC protocol is based on DICOM as specified in each DICOM Conformance Statement. Evolution of the Standard may require changes to devices which have implemented DICOM. **In addition, GE reserves the right to discontinue or make changes to the support of communications features (on its products) described by these DICOM Conformance Statements.** The **user** should ensure that any non-GE provider, which connects with GE devices, also plans for the future evolution of the DICOM Standard. Failure to do so will likely result in the loss of function and/or connectivity as the DICOM Standard changes and GE Products are enhanced to support these changes.
- **Interaction** - It is the sole responsibility of the **non-GE provider** to ensure that communication with the interfaced equipment does not cause degradation of GE imaging equipment performance and/or function.

1.6 REFERENCES

NEMA PS3 Digital Imaging and Communications in Medicine (DICOM) Standard, available free at <http://medical.nema.org/>

1.7 DEFINITIONS

Informal definitions are provided for the following terms used in this Conformance Statement. The DICOM Standard is the authoritative source for formal definitions of these terms.

Abstract Syntax – the information agreed to be exchanged between applications, generally equivalent to a Service/Object Pair (SOP) Class. Examples: Verification SOP Class, Modality Worklist Information Model Find SOP Class, Computed Radiography Image Storage SOP Class.

Application Entity (AE) – an end point of a DICOM information exchange, including the DICOM network or media interface software; i.e., the software that sends or receives DICOM information objects or messages. A single device may have multiple Application Entities.

Application Entity Title – the externally known name of an *Application Entity*, used to identify a DICOM application to other DICOM applications on the network.

Application Context – the specification of the type of communication used between *Application Entities*. Example: DICOM network protocol.

Association – a network communication channel set up between *Application Entities*.

Attribute – a unit of information in an object definition; a data element identified by a *tag*. The information may be a complex data structure (Sequence), itself composed of lower level data elements. Examples: Patient ID (0010,0020), Accession Number (0008,0050), Photometric Interpretation (0028,0004), Procedure Code Sequence (0008,1032).

Information Object Definition (IOD) – the specified set of *Attributes* that comprise a type of data object; does not represent a specific instance of the data object, but rather a class of similar data objects that have the same properties. The *Attributes* may be specified as Mandatory (Type 1), Required but possibly unknown (Type 2), or Optional (Type 3), and there may be conditions associated with the use of an Attribute (Types 1C and 2C). Examples: MR Image IOD, CT Image IOD, Print Job IOD.

Joint Photographic Experts Group (JPEG) – a set of standardized image compression techniques, available for use by DICOM applications.

Media Application Profile – the specification of DICOM information objects and encoding exchanged on removable media (e.g., CDs)

Module – a set of *Attributes* within an *Information Object Definition* that are logically related to each other. Example: Patient Module includes Patient Name, Patient ID, Patient Birth Date, and Patient Sex.

Negotiation – first phase of *Association* establishment that allows *Application Entities* to agree on the types of data to be exchanged and how that data will be encoded.

Presentation Context – the set of DICOM network services used over an *Association*, as negotiated between *Application Entities*; includes *Abstract Syntaxes* and *Transfer Syntaxes*.

Protocol Data Unit (PDU) – a packet (piece) of a DICOM message sent across the network. Devices must specify the maximum size packet they can receive for DICOM messages.

Security Profile – a set of mechanisms, such as encryption, user authentication, or digital signatures, used by an *Application Entity* to ensure confidentiality, integrity, and/or availability of exchanged DICOM data

Service Class Provider (SCP) – role of an *Application Entity* that provides a DICOM network service; typically, a server that performs operations requested by another *Application Entity* (*Service Class User*). Examples: Picture Archiving and Communication System (image storage SCP, and image query/retrieve SCP), Radiology Information System (modality worklist SCP).

Service Class User (SCU) – role of an *Application Entity* that uses a DICOM network service; typically, a client. Examples: imaging modality (image storage SCU, and modality worklist SCU), imaging workstation (image query/retrieve SCU)

Service/Object Pair (SOP) Class – the specification of the network or media transfer (service) of a particular type of data (object); the fundamental unit of DICOM interoperability specification. Examples: Ultrasound Image Storage Service, Basic Grayscale Print Management.

Service/Object Pair (SOP) Instance – an information object; a specific occurrence of information exchanged in a *SOP Class*. Examples: a specific x-ray image.

Tag – a 32-bit identifier for a data element, represented as a pair of four digit hexadecimal numbers, the “group” and the “element”. If the “group” number is odd, the tag is for a private (manufacturer-specific) data element. Examples: (0010,0020) [Patient ID], (07FE,0010) [Pixel Data], (0019,0210) [private data element]

Transfer Syntax – the encoding used for exchange of DICOM information objects and messages. Examples: *JPEG* compressed (images), little endian explicit value representation.

Unique Identifier (UID) – a globally unique “dotted decimal” string that identifies a specific object or a class of objects; an ISO-8824 Object Identifier. Examples: Study Instance UID, SOP Class UID, SOP Instance UID.

Value Representation (VR) – the format type of an individual DICOM data element, such as text, an integer, a person’s name, or a code. DICOM information objects can be transmitted with either explicit identification of the type of each data element (Explicit VR), or without explicit identification (Implicit VR); with Implicit VR, the receiving application must use a DICOM data dictionary to look up the format of each data element.

1.8 SYMBOLS AND ABBREVIATIONS

| | |
|-------|------------------------------------------------|
| AE | Application Entity |
| AET | Application Entity Title |
| CAD | Computer Aided Detection |
| CD-R | Compact Disk Recordable |
| CR | Computed Radiography |
| CT | Computed Tomography |
| DB | Database |
| DHCP | Dynamic Host Configuration Protocol |
| DICOM | Digital Imaging and Communications in Medicine |
| DX | Digital X-ray |
| FSC | File-Set Creator |

| | |
|--------|-------------------------------------------------|
| FSU | File-Set Updater |
| FSR | File-Set Reader |
| HL7 | Health Level 7 Standard |
| HLF | High Level Fluoro |
| IOD | Information Object Definition |
| IPv4 | Internet Protocol version 4 |
| IPv6 | Internet Protocol version 6 |
| ISO | International Organization for Standards |
| JPEG | Joint Photographic Experts Group |
| LUT | Look-up Table |
| MPPS | Modality Performed Procedure Step |
| MR | Magnetic Resonance Imaging |
| MWL | Modality Worklist |
| O | Optional (Key Attribute) |
| PACS | Picture Archiving and Communication System |
| PDU | Protocol Data Unit |
| R | Required (Key Attribute) |
| RDSR | Radiation Dose Structured Reports |
| RF | Radiofluoroscopy |
| SC | Secondary Capture |
| SCP | Service Class Provider |
| SCU | Service Class User |
| SOP | Service-Object Pair |
| SPS | Scheduled Procedure Step |
| SR | Structured Reporting |
| TCP/IP | Transmission Control Protocol/Internet Protocol |
| U | Unique (Key Attribute) |
| VR | Value Representation |

OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT

GE Healthcare
2148033 REV 1

XA X-ray Angiography

2. NETWORK CONFORMANCE STATEMENT

2.1 INTRODUCTION

This section of the DICOM Conformance Statement specifies the Workstation compliance to DICOM requirements for **Networking** features.

The Workstation runs on industrial computer. It provides the following DICOM functionalities:


- Send Echo message to DICOM Verification SCP.
- Export DICOM images and RDSR to a DICOM Storage SCP.
- Query DICOM Modality Worklist from a DICOM Worklist SCP.
- Querying for saved examinations and retrieving saved examinations from a DICOM Query/Retrieve SCP
- Sending exam start and end messages to a DICOM Modality Performed Procedure Step SCP.
- Sending storage commitment requests to a DICOM Storage Commitment SCP.
- Print acquired images to a DICOM Printer.

2.2 IMPLEMENTATION MODEL

2.2.1 Application Data Flow Diagram

The network application model for the Workstation is shown in the following Illustration :

ILLUSTRATION 2-1
WORKSTATION NETWORK APPLICATION MODEL AND DATA FLOW DIAGRAM



The Workstation supports activity “Storage” which performs images/RDSR storage and storage commitment functions; The results of that functionalities are communicated to Remote Store Server for DICOM image/RDSR storage and storage commitment.

The Workstation supports activity “Print” which performs images print functions; The results of that functionalities are communicated to Remote Printer/Server for DICOM image print functions.

The Workstation supports activity “Verify” which interacts with Remote Verify SCP for DICOM Echo purpose.

The Workstation supports activity “Modality Worklist Query” which performs Modality worklist query, from remote Modality Worklist SCP.

The Workstation supports activity “Query”: Initiate an association with a Remote AE to query for images on the remote host. Once the association has been established, it sends a sequence of requests (Patient, then Study then Series Level requests) to the Remote AE. After all responses are received, Workstation AE will issue a Series-Level C-FIND-RQ request to get the series for a study in the list.

The Workstation supports activity “Retrieve”: Send a C-MOVE-RQ request to a Remote AE after successful association establishment. The Workstation AE’s Storage SCP will receive the images over a separate association.

The Workstation supports activity “Report MPPS”: When the user begins the image acquisition process and generates the first image, the Workstation AE sends N-CREATE message to the configured MPPS SCP to indicate that the image acquisition process has been started for the requested procedure. The operator can close the acquisition session either by completing the acquisition process or discontinuing the ongoing scan. On closing the acquisition session, the Workstation AE sends N-SET message to the configured MPPS SCP to indicate the acquisition state of the requested procedure, with appropriate MPPS status (COMPLETED/DISCONTINUED).

2.2.2 Functional Definition of AE's

Application Entity “OEC One Workstation AE ” initiates the following functionalities :

- DICOM Host Verification: Initiate an association and send a C-ECHO-RQ message to the remote DICOM AE; the remote DICOM server will send back a C-ECHO-RSP message with a status of “success” if the operation is successful.
- DICOM Image Storage: Initiate an association to a remote AE to send image(s). If the remote AE accept the presentation context applicable to the image(s), the WorkstationAE will send the image(s) by invoking C-STORE-RQ operation for each image on the same association. The remote DICOM server will send back C-STORE-RSP with status “success” if the operation is successful.
- Storage Commitment: Initiates a DICOM Storage Commitment request once a Store request has completed successfully. An N-ACTION request specifying all the successfully transferred images will be issued to the remote AE. The remote AE will reply accepting or refusing the request on the same association. If accepted the remote AE will issue an N-EVENT response specifying which images have been committed and which, if any, have not.
- Modality Worklist Query: Initiate an association with a remote AE to query for scheduled exams, a C-FIND-RQ request will be sent to the remote AE. C-FIND-RSP with query results will be received.
- Query: Initiate an association with a Remote AE to query for images on the remote host. Follows the Patient->Study->Series sequence C-FIND-RQ request (the request more than once) will be sent to the Remote AE once an association has been established. Once all the responses have been received, the operator needs to select an exam in the local database browser; on selection of the exam the DICOM Server AE will issue a Series-Level C-FIND-RQ request to get the series for a study in the

list. Similarly the Image-Level C-FIND-RQ will be issued for the series selected from the series list.

- Retrieve: Send a C-MOVE-RQ request to a Remote AE for retrieve of images after successful establishment. At this time the Local Store SCP listener will start, not until after receiving complete images the Local Store SCP listener will stopped.
- DICOM Print: Initiate an association with a remote AE to print grayscale images, please refer to Section 2.3.1.2.6.2.1 for details

2.2.3 Sequencing of Real-World Activities

This sequence is only applicable for Remote AE where Storage Commitment Option is allowed in Network Manager.

1. The user selects the images and sends them to a remote host.
2. If the remote DICOM AE is associated with a Storage Commitment Provider AE and if the images are successfully sent to the DICOM AE, then an N-ACTION-RQ request is sent automatically to the associated Storage Commitment Provider AE. The Storage Commitment Provider AE can be configured independently from the remote DICOM AE with network address, port and AE title.

2.3 AE SPECIFICATIONS

2.3.1 Workstation AE Specification

The Workstation Application Entity provides Standard Conformance to the following DICOM SOP Classes as an **SCU** and/or as an **SCP**:

| SOP Class Name | SOP Class UID | SCU | SCP |
|------------------------------------------------------|-------------------------------|-----|-----|
| Verification SOP Class | 1.2.840.10008.1.1 | Yes | Yes |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | No | Yes |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | No | Yes |
| X-Ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 | Yes | Yes |
| X-Ray Radiofluoroscopic Image Storage | 1.2.840.10008.5.1.4.1.1.12.2 | No | Yes |
| X-Ray Radiation Dose SR | 1.2.840.10008.5.1.4.1.1.88.67 | Yes | No |
| Patient Root Query/Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.1.1 | Yes | No |
| Patient Root Query/Retrieve Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.1.2 | Yes | No |
| Study Root Query/Retrieve Information Model - FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Yes | No |
| Study Root Query/Retrieve Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.2.2 | Yes | No |
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 | Yes | No |

| | | | |
|-------------------------------------------------|-------------------------|-----|----|
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Yes | No |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Yes | No |
| Basic Film Session SOP Class | 1.2.840.10008.5.1.1.1 | Yes | No |
| Basic Film Box SOP Class | 1.2.840.10008.5.1.1.2 | Yes | No |
| Basic Grayscale Image Box SOP Class | 1.2.840.10008.5.1.1.4 | Yes | No |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9 | Yes | No |
| Printer SOP Class | 1.2.840.10008.5.1.1.16 | Yes | No |

2.3.1.1 Association Establishment Policies

2.3.1.1.1 General

The DICOM Application Context Name (ACN), which is always proposed, is:

| | |
|---------------------------------|------------------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|---------------------------------|------------------------------|

The maximum length PDU receive size for the Workstation AE is:

| | |
|---------------------------|--------------|
| Maximum Length PDU | 16384 |
|---------------------------|--------------|

2.3.1.1.2 Number of Associations

The Workstation AE will initiate a maximum of 1 simultaneous association to remote nodes.

The Workstation AE will support a maximum of 1 simultaneous associations initiated by remote nodes.

2.3.1.1.3 Asynchronous Nature

Asynchronous mode is not supported. All operations will be performed synchronously.

2.3.1.1.4 Implementation Identifying Information

The Implementation UID for this DICOM Implementation is:

| | |
|--------------------------------------------|-----------------------------|
| OEC One Implementation UID | 1.2.840.113619.6.432 |
| OEC One Implementation Version Name | GEHC_OEC_200 |

2.3.1.2 Association Initiation Policy

When the Workstation Application Entity initiates an Association for any Real-World Activity, it will propose the Presentation Contexts for all Real-World Activities; i.e., there is only a single, comprehensive Presentation Context Negotiation proposed for the AE.

The Workstation AE proposes only a single Transfer Syntax in each Presentation Context; i.e., for each Abstract Syntax in the following Presentation Context Tables, the AE proposes one Presentation Context for each specified Transfer Syntax.

2.3.1.2.1 Real-World Activity Verify

2.3.1.2.1.1 Associated Real-World Activity

The user may initiate a DICOM Verify Request in the configuration screen for each respective remote SCP configuration (Print, Store/Store Commit, Worklist Query, MPPS, SC and Q/R).

A valid response from the SCP will result in a “The DICOM server has verified successfully” screen displayed on the screen. In the event that the SCP does not response for some reason, the operation will timeout and the Workstation AE will close the association and results in “The DICOM Verify Fails” screen displayed to user.

Note: The default timeout is 30 seconds for DICOM Echo and this value could be configured from system setup screen.

2.3.1.2.1.2 Proposed Presentation Context Table

| Presentation Context Table – Proposed by Workstation AE for Activity Verify | | | | | |
|------------------------------------------------------------------------------------|-------------------|---------------------------|---------------------|-------------|-----------------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Verification SOP Class | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

2.3.1.2.1.2.1 SOP Specific DICOM Conformance Statement for Verify SOP Classes

The message “The DICOM server has verified successfully” is displayed when the DICOM Verify is successfully performed, otherwise message “The DICOM Verify Fails” is displayed.

2.3.1.2.2 Real-World Activity Schedule Update

2.3.1.2.2.1 Associated Real-World Activity

Upon a request by the operator, the current scheduled exam list should be updated from the remote Modality Worklist Server.

2.3.1.2.2.2 Proposed Presentation Context Table

| Presentation Context Table – Proposed by Workstation AE for Activity Modality Worklist | | |
|-----------------------------------------------------------------------------------------------|------------------------|-------------|
| Abstract Syntax | Transfer Syntax | Role |

| Name | UID | Name List | UID List | |
|--------------------------------------------|------------------------|---------------------------|---------------------|-----|
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | |

2.3.1.2.2.2.1 SOP Specific DICOM Conformance Statement for the Modality Worklist Information Model - FIND SOP Class

The Workstation AE includes matching keys in the Modality Worklist queries as described in Section 6.

The results of the Worklist query are displayed to the user in the Scheduled Exam Screen. Currently, the maximal capability of received scheduled exams is 500 and not configurable.

The Workstation AE checks the Specific Character Set tag (0008,0005) contained in a C-FIND RSP, ISO_IR 100 or blank will be accepted, and all other value will be discarded. Responses without this tag will also be accepted.

A C-FIND CANCEL will be sent when the maximal capability of received scheduled exams(500 exams) is reached or user cancels operation by clicking the Cancel button on the progressing screen.

Following are the status codes that are more specifically processed when receiving messages from a **Modality Worklist SCP** equipment:

| Service Status | Status Code | Further Meaning | Application Behavior When Receiving Status Code |
|----------------|-------------|------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Failure | A700 | Refused: Out of resources | Terminates the association and operation. Logs failure message along with error comment returned from SCP and displays failure message. |
| | A900 | Error: Identifier does not match SOP Class | |
| | C000-CFFF | Error: Unable to process | |
| | 0122 | SOP Class Not Supported | |
| Cancel | FE00 | Matching terminated due to cancel | The progress screen is disappeared, the received scheduled exams will be displayed. |
| Success | 0000 | Matching is complete - No final identifier is supplied | The progress screen is disappeared, all the received scheduled exams will be displayed. |
| Pending | FF00 | Matches are continuing - Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys. | The progress screen is still displayed. |
| | FF01 | Matches are continuing - Warning that one or more | The progress screen is still displayed. |

| | | | |
|---|---|--------------------------------------------------------------------|-------------------|
| | | Optional Keys were not supported for existence for this Identifier | |
| * | * | Any other status code. | Deemed a failure. |

2.3.1.2.3 Real-World Activity Start/End Exam

2.3.1.2.3.1 Associated Real-World Activity

The system will start collecting X-ray dose information for the selected exam as soon as the X-ray activation switch is pressed. This starts the exam.

The user is provided a control to end/abandon an exam on the Image and Saved Exams screens. This will create the RDSR and optionally send the RDSR to an RDSR SCP.

2.3.1.2.3.2 Proposed Presentation Context Table

See Image RDSR Send and Storage Commit section [2.3.1.2.4](#) to see RDSR conformance.

2.3.1.2.4 Real-World Activity Image/RDSR Send and Storage Commit

2.3.1.2.4.1 Associated Real-World Activity

The user is able to store Images and RDSRs from the Images screen. Images can be stored at any time during or after the exam. RDSRs are only available after an exam is ended with a End Exam or Complete Exam.

2.3.1.2.4.2 Proposed Presentation Context Table

| Presentation Context Table – Proposed by Workstation AE for Activity Storage | | | | |
|-------------------------------------------------------------------------------------|-------------------------------|----------------------------------------------------------------------------------|-----------------------------------------------------------------|-------------|
| Abstract Syntax | | Transfer Syntax | | Role |
| Name | UID | Name List | UID List | |
| X-Ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 | Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 | SCU |
| X-Ray Radiation Dose SR | 1.2.840.10008.5.1.4.1.1.88.67 | Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 | SCU |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 | SCU |

2.3.1.2.4.2.1 SOP Specific DICOM Conformance Statement for All Storage SOP Classes

The Workstation AE includes optional data elements in the SOP Instances as described in Sections 4.

Following are the status codes that are more specifically processed when receiving messages from a **Storage** SCP equipment : If store command completed normally, with a status of success or some warning then the image was accepted.

| Service Status | Status Code | Further Meaning | Application Behavior When Receiving Status Code |
|----------------|-------------|------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Failure | A700-A7FF | Refused: Out of resources | Terminate the operation. Logs store failed message along with error comment returned from SCP. Displays store failure message. |
| | A900-A9FF | Error: Data Set does not match SOP Class | |
| | C000-CFFF | Error: Cannot Understand | |
| | 0122 | SOP Class Not Supported | |
| Warning | B000 | Coercion of Data Elements | Ignored by application and current operation continues |
| | B006 | Elements Discarded | |
| | B007 | Data Set does not match SOP Class | |
| Success | 0000 | DICOM storage operation is success | The DICOM operation progress screen disappearance indicates current operation is success. |
| * | * | Any other status code. | Upon receiving a C-STORE confirmation containing a Refused status, this implementation will terminate the association. |

2.3.1.2.4.2.2 SOP Specific DICOM Conformance Statement for the X-Ray Radiation Dose Storage SOP Classes

See Section 2.3.1.2.4.2.1 “SOP Specific DICOM Conformance Statement for All Storage SOP Classes” for details on general Storage Service SCU processing also applicable to the X-Ray Radiation Dose Storage SOP Classes.

The Workstation AE supports creation and transmission of X-Ray Radiation Dose SOP Instances referencing Instances of the following Storage SOP Classes:

| SOP Class Name | SOP Class UID |
|----------------------------------|-------------------------------|
| X-Ray Radiation Dose SR | 1.2.840.10008.5.1.4.1.1.88.67 |
| X-Ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 |

2.3.1.2.4.2.3 SOP Specific DICOM Conformance Statement for the Storage Commitment Push Model SOP Class SCU

The Workstation AE uses DICOM network storage services to transfer SOP Instances which are to be committed.

The Workstation AE may request Storage Commitment for Instances of any of the Composite SOP Classes it supports as an SCU (see Section 2.3.1.2.4.2.1).

The Storage Commitment Information Object is described in Section 8.

Upon receiving a Storage Commitment N-EVENT-REPORT (Storage Commitment Result on the same association), the Workstation AE will return a Success status.

Note: N-EVENT-REPORT is only expected on the same association.

Following are the status codes that are more specifically processed when receiving N-Action responses from Storage Commitment SCP equipment:

| Service Status | Status Code | Further Meaning | Application Behavior When Receiving Status Code |
|-----------------------|------------------------|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Failure | 0119 | Class-instance conflict | Terminates the association and operation. Logs failure message along with error comment returned from SCP and displays failure message. |
| | 0210 | Duplicate invocation | |
| | 0115 | Invalid argument value | |
| | 0117 | Invalid SOP Instance | |
| | 0212 | Mistyped argument | |
| | 0123 | No such action | |
| | 0114 | No such argument | |
| | 0118 | No such SOP Class | |
| | 0112 | No such SOP Instance | |
| | 0110 | Processing failure | |
| | 0213 | Resource limitation | |
| 0211 | Unrecognized operation | | |
| Success | 0000 | | The progress screen is disappeared. |
| * | * | Any other status code. | Deemed a failure. |

2.3.1.2.5 Real-World Activity Image Query/Retrieve and Local Store SCP

2.3.1.2.5.1 Associated Real-World Activity

Upon a request from the user, a C-FIND-RQ command will be issued to the Query/Retrieve SCP to query for exams that match the user-defined query criteria. The association will be released upon the receipt of the C-FIND-RSP confirmation. To retrieve an exam, the user can either use the information returned from a query or the information stored in the local database when the exam was archived. The C-MOVE-RQ command is used to retrieve the exam from the Query/Retrieve SCP. The C-STORE-RQ from the SCP is expected on another association.

After a DICOM file is retrieved, the attributes are extracted and saved onto the local system.

2.3.1.2.5.2 Proposed Presentation Context Table

| Presentation Context Table – Proposed by Workstation AE for Activity Modality Worklist | | | | |
|-----------------------------------------------------------------------------------------------|-----------------------------|----------------------------------------------------------------------------------|-----------------------------------------------------------------|-------------|
| Abstract Syntax | | Transfer Syntax | | Role |
| Name | UID | Name List | UID List | |
| Study Root Query/Retrieve Information Model - FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 | SCU |
| Study Root Query/Retrieve Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.2.2 | Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 | SCU |
| Patient Root Query/Retrieve Information Model - FIND | 1.2.840.10008.5.1.4.1.2.1.1 | Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 | SCU |
| Patient Root Query/Retrieve Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.1.2 | Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 | SCU |

2.3.1.2.5.2.1 SOP Specific DICOM Conformance Statement for the Patient Root Query/Retrieve Information Model - FIND, Study Root Query/Retrieve Information Model - FIND SOP Classes

The Workstation AE includes matching keys in the queries as described in Sections 11.

The Acquisition AE evaluates a C-FIND RSP containing the Specific Character Set tag (0008,0005). Values of ISO_IR 100 or blank will be accepted, and all others will be rejected.

Responses without this tag will also be accepted.

A C-FIND CANCEL will be sent when the display limit for results is reached.

The Workstation does not support Extended Negotiation for C-FIND.

Following are the status codes that are more specifically processed when receiving messages from a Query SCP equipment:

| Service Status | Status Code | Further Meaning | Application Behavior When Receiving Status Code |
|----------------|-------------|--------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Failure | A700 | Refused: Out of resources | Terminates the association and operation. Logs failure message along with error comment returned from SCP and displays failure message. |
| | A900 | Error: Identifier does not match SOP Class | |
| | C000-CFFF | Error: Unable to process | |
| | 0122 | SOP Class Not Supported | |
| Cancel | FE00 | Matching terminated due to cancel | The progress screen is disappeared, the received scheduled exams will be displayed |
| Success | 0000 | Matching is complete - No final identifier is supplied | The progress screen is disappeared, all the received scheduled exams will be displayed. |
| Pending | FF00 | Matches are continuing - Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys. | The progress screen is still displayed. |
| | FF01 | Matches are continuing - Warning that one or more Optional Keys were not supported for existence and/or matching for this Identifier | The progress screen is still displayed. |
| * | * | Any other status code. | Deemed a failure. |

2.3.1.2.5.2.2 SOP Specific DICOM Conformance Statement for the Patient Root Query/Retrieve Information Model - MOVE , Study Root Query/Retrieve Information Model - MOVE SOP Classes

The C-MOVE-RQ will use the AE Title of the Workstation Application Entity as the Move Destination AE Title.

A C-FIND CANCEL message is sent if the user clicks the Cancel button during the move.

Following are the status codes that are more specifically processed when receiving messages from a **Retrieve** SCP equipment:

| Service Status | Status Code | Further Meaning | Application Behavior When Receiving Status Code |
|----------------|-------------|-----------------------------------------------------------|--------------------------------------------------------------------------------------------------|
| Failure | A701 | Refused: Out of resources - Unable to calculate number of | Terminates the association and operation. Logs failure message along with error comment returned |

| | | | |
|---------|-----------|--------------------------------------------------------------|-----------------------------------------------------------------------------------------|
| | | matches | from SCP and displays failure message. |
| | A702 | Refused: Out of resources - Unable to perform sub-operations | |
| | A801 | Refused: Move Destination Unknown | |
| | A900 | Error: Identifier does not match SOP Class | |
| | C000-CFFF | Error: Unable to process | |
| | 0122 | SOP Class Not Supported | |
| Cancel | FE00 | Sub-operations terminated due to a Cancel indication | The progress screen is disappeared, the received scheduled exams will be displayed |
| Warning | B000 | Sub-operations Complete - One or more Failures. | Ignored by application and current operation continues. |
| Success | 0000 | Sub-operations Complete - No Failure. | The progress screen is disappeared, all the received scheduled exams will be displayed. |
| Pending | FF00 | Sub-operations are continuing - | The progress screen is still displayed. |
| * | * | Any other status code. | Deemed a failure. |

2.3.1.2.6 Real-World Activity “DICOM Print”

2.3.1.2.6.1 Associated Real-World Activity

Upon a request by the operator the selected images will be sent to the configured DICOM Printer (DICOM Print SCP).

2.3.1.2.6.2 Proposed Presentation Context Table

| Presentation Context Table – Proposed by AE Workstation for Activity DICOM Print | | | | | |
|-----------------------------------------------------------------------------------------|-----------------------|---------------------------|---------------------|-------------|-----------------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Basic Film Session SOP Class | 1.2.840.10008.5.1.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Basic Film Box SOP Class | 1.2.840.10008.5.1.1.2 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

**OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | | | |
|-------------------------------------------------|------------------------|---------------------------|---------------------|-----|------|
| Basic Grayscale Image Box SOP Class | 1.2.840.10008.5.1.1.4 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Printer SOP Class | 1.2.840.10008.5.1.1.16 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

2.3.1.2.6.2.1 SOP Specific DICOM Conformance Statement for Basic Grayscale Print Management SOP Classes

The Workstation AE uses the following DIMSE services of the supported SOP Classes:

| SOP Class | SOP Class UID | DIMSE Service Element | SCU Usage |
|----------------------------------|------------------------|-----------------------|------------------|
| Basic Film Session | 1.2.840.10008.5.1.1.1 | N-CREATE | Used (Mandatory) |
| | | N-SET | Not Used |
| | | N-DELETE | Used |
| | | N-ACTION | Not Used |
| Basic Film Box | 1.2.840.10008.5.1.1.2 | N-CREATE | Used (Mandatory) |
| | | N-ACTION | Used (Mandatory) |
| | | N-DELETE | Used |
| | | N-SET | Not Used |
| Printer | 1.2.840.10008.5.1.1.16 | N-EVENT-REPORT | Used (Mandatory) |
| | | N-GET | Used |
| Basic Grayscale Image Box | 1.2.840.10008.5.1.1.4 | N-SET | Used (Mandatory) |

When a manual print operation is initiated, the AE:

1. Initiates a DICOM association and negotiates Presentation Contexts
2. N-GETs printer status using the Printer SOP Class
3. N-CREATEs a Basic Film Session SOP Instance, and then followed by N-DELETE of the Basic Film Session SOP Instance.
4. N-CREATEs a Basic Film Box SOP Instance for each film
5. N-SETs the Image Box SOP Instance for each image on the film
6. Prints by N-ACTION on the Basic Film Box SOP Instances, then followed by N-DELETE of the Basic Film Box SOP Instance
7. Receives N-EVENT-REPORTs of the well know Printer SOP Instance indicating printer status if the SCP send a request.

8. Releases the DICOM association after printing is successful or failure has been signaled to the user

9. The DICOM operation progress screen disappearance indicates current operation is success.

The Print management IODs is defined in Section.

2.3.1.2.6.2.1.1 Basic Film Session SOP Class

Following are the status codes that are more specifically processed when receiving messages from a **Print** SCP equipment for the Basic Film Session SOP Class N-CREATE:

| Service Status | Status Codes | Further Meaning | Application Behavior When receiving Status Codes |
|-----------------------|------------------------|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Failure | 0119 | Class-instance conflict | Terminate the operation. Logs store failed message along with error comment returned from SCP. Displays print failure message. |
| | 0210 | Duplicate invocation | |
| | 0111 | Duplicate SOP Instance | |
| | 0106 | Invalid attribute value | |
| | 0117 | Invalid SOP instance | |
| | 0120 | Missing attribute | |
| | 0121 | Missing attribute value | |
| | 0212 | Mistyped argument | |
| | 0105 | No such attributes | |
| | 0118 | No such SOP Class | |
| | 0112 | No such SOP Instance | |
| | 0110 | Processing failure | |
| | 0213 | Resource limitation | |
| 0211 | Unrecognized operation | | |
| Warning | B600 | Memory allocation not supported | Ignored by application and current operation continues |
| Success | 0000 | Film session successfully created | The progress screen is still displayed. |
| * | * | Any other status code. | Deemed a failure. |

Following are the status codes that are more specifically processed when receiving messages from a **Print** SCP equipment for the Basic Film Session SOP Class N-DELETE:

| Service Status | Status Codes | Further Meaning | Application Behavior When receiving Status Codes |
|-----------------------|---------------------|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Failure | 0119 | Class-instance conflict | Terminate the operation. Logs store failed message along with error comment returned from SCP. Displays print failure message. |
| | 0210 | Duplicate invocation | |
| | 0117 | Invalid SOP instance | |
| | 0212 | Mistyped argument | |
| | 0118 | No such SOP Class | |
| | 0112 | No such SOP Instance | |
| | 0110 | Processing failure | |
| | 0213 | Resource limitation | |
| | 0211 | Unrecognized operation | |
| Success | 0000 | Film session successfully deleted | The progress screen is still displayed. |
| * | * | Any other status code. | Deemed a failure. |

2.3.1.2.6.2.1.2 Basic Film Box SOP Class

Following are the status codes that are more specifically processed when receiving messages from a **Print** SCP equipment for the Basic Film Box SOP Class N-CREATE:

| Service Status | Status Codes | Further Meaning | Application Behavior When receiving Status Codes |
|-----------------------|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Failure | C616 | There is an existing Film Box that has not been printed and N-ACTION at the Film Session level is not supported. A new Film Box will not be created when a previous Film Box has not been printed. | Terminate the operation. Logs store failed message along with error comment returned from SCP. Displays print failure message. |
| | 0119 | Class-instance conflict | |
| | 0210 | Duplicate invocation | |
| | 0111 | Duplicate SOP Instance | |
| | 0106 | Invalid attribute value | |

| | | | |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------|-----------------------------------------|
| Warning | 0117 | Invalid SOP instance | |
| | 0120 | Missing attribute | |
| | 0121 | Missing attribute value | |
| | 0212 | Mistyped argument | |
| | 0105 | No such attributes | |
| | 0118 | No such SOP Class | |
| | 0112 | No such SOP Instance | |
| | 0110 | Processing failure | |
| | 0213 | Resource limitation | |
| | 0211 | Unrecognized operation | |
| B605 | Requested Min Density or Max Density outside of printer's operating range. The printer will use its respective minimum or maximum density value instead. | Ignored by application and current operation continues. | |
| Success | 0000 | Film box successfully created | The progress screen is still displayed. |
| * | * | Any other status code. | Deemed a failure. |

Following are the status codes that are more specifically processed when receiving messages from a **Print** SCP equipment for the Basic Film Box SOP Class N-ACTION:

| Service Status | Status Codes | Further Meaning | Application Behavior When receiving Status Codes |
|-----------------------|---------------------|---------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Failure | C602 | Unable to create Print Job SOP Instance; print queue is full | Terminate the operation. Logs store failed message along with error comment returned from SCP. Displays print failure message. |
| | C603 | Image size is larger than image box size (by using the specified magnification value) | |
| | C604 | Image position collision : multiple images assigned to single image position | |
| | C613 | Combined Print Image size is larger than the Image Box size | |
| | 0119 | Class-instance conflict | |
| | 0210 | Duplicate invocation | |
| | 0115 | Invalid argument value | |
| | 0117 | Invalid SOP Instance | |
| | 0212 | Mistyped argument | |

| | | | |
|---------|------|-------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------|
| | 0123 | No such action | |
| | 0114 | No such argument | |
| | 0118 | No such SOP Class | |
| | 0112 | No such SOP Instance | |
| | 0110 | Processing failure | |
| | 0213 | Resource limitation | |
| | 0211 | Unrecognized operation | |
| Warning | B603 | Film Box SOP Instance hierarchy does not contain Image Box SOP Instances (empty page) | Ignored by application and current operation continues. |
| | B604 | Image size is larger than image box size, the image has been demagnified. | |
| | B609 | Image size is larger than the Image Box size. The Image has been cropped to fit. | |
| | B60A | Image size or Combined Print Image size is larger than the Image Box size. Image or Combined Print Image has been decimated to fit. | |
| Success | 0000 | Film accepted for printing; if supported, the Print Job SOP Instance is created | The progress screen is still displayed. |
| * | * | Any other status code. | Deemed a failure. |

Following are the status codes that are more specifically processed when receiving messages from a **Print** SCP equipment for the Basic Film Box SOP Class N-DELETE:

| Service Status | Status Codes | Further Meaning | Application Behavior When receiving Status Codes |
|-----------------------|---------------------|-------------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Failure | 0119 | Class-instance conflict | Terminate the operation. Logs store failed message along with error comment returned from SCP. Displays print failure message. |
| | 0210 | Duplicate invocation | |
| | 0117 | Invalid SOP instance | |
| | 0212 | Mistyped argument | |
| | 0118 | No such SOP Class | |
| | 0112 | No such SOP Instance | |
| | 0110 | Processing failure | |
| | 0213 | Resource limitation | |
| | 0211 | Unrecognized operation | |

| | | | |
|---------|------|-------------------------------|-----------------------------------------|
| Success | 0000 | Film box successfully deleted | The progress screen is still displayed. |
| * | * | Any other status code. | Deemed a failure. |

2.3.1.2.6.2.1.3 Basic Grayscale Image Box SOP Class

Following are the status codes that are more specifically processed when receiving messages from a **Print SCP** equipment for the Basic Grayscale Image Box SOP Class N-SET:

| Service Status | Status Codes | Further Meaning | Application Behavior When receiving Status Codes |
|-----------------------|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Failure | C603 | Image size is larger than image box size | Terminate the operation. Logs store failed message along with error comment returned from SCP. Displays print failure message. |
| | C605 | Insufficient memory in printer to store the image | |
| | C613 | Combined Print Image size is larger than the Image Box size | |
| | 0119 | Class-instance conflict | |
| | 0210 | Duplicate invocation | |
| | 0106 | Invalid attribute value | |
| | 0212 | Mistyped argument | |
| | 0117 | Invalid SOP instance | |
| | 0121 | Missing attribute value | |
| | 0105 | No such attributes | |
| | 0118 | No such SOP Class | |
| | 0112 | No such SOP Instance | |
| | 0110 | Processing failure | |
| | 0213 | Resource limitation | |
| 0211 | Unrecognized operation | | |
| Warning | B604 | Image size larger than image box size, the image has been demagnified. | Ignored by application and current operation continues. |
| | B605 | Requested Min Density or Max Density outside of printer's operating range. The printer will use its respective minimum or maximum density value instead. | |

| | | | |
|---------|------|-----------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------|
| | B609 | Image size is larger than the Image Box size. The Image has been cropped to fit. | |
| | B60A | Image size or Combined Print Image size is larger than the Image Box size. The Image or Combined Print Image has been decimated to fit. | |
| Success | 0000 | Image successfully stored in Image Box | The progress screen is still displayed. |
| * | * | Any other status code. | Deemed a failure. |

2.3.1.2.6.2.1.4 Printer SOP Class

The Workstation AE supports the Printer SOP Class to receive information on the status of the printer.

For the Workstation behavior description when receiving N-Event-Report request, refer to section 10.5.1 Printer N-Event-Report Attributes.

Also, following are the status codes the Application may send back in the **N-Event-Report** response command to the **Printer SOP Class** SCP Equipment that sent the N-Event-Report request:

| Service Status | Status Code | Further Meaning | Status Code Explanation | Related Fields Sent Back to the SCP |
|----------------|-------------|-----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------|
| Success | 0000 | | If an N-EVENT-REPORT DIMSE service is received when the association is active, Print SCU handles the relevant states but the data received is ignored and send back "Success". | None |

For the attribute list requested by this product and for the behavior of this product on each returned value, refer to Section 10.5.22 Printer N-Get Attributes

Following are the status codes that are more specifically processed when receiving messages from a **Print** SCP equipment for the Printer SOP Class N-GET:

| Service Status | Status Codes | Further Meaning | Application Behavior When receiving Status Codes |
|----------------|--------------|-------------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Failure | 0107 | Attribute list error | Terminate the operation. Logs store failed message along with error comment returned from SCP. Displays print failure message. |
| | 0119 | Class-instance conflict | |
| | 0210 | Duplicate invocation | |
| | 0117 | Invalid SOP instance | |
| | 0212 | Mistyped argument | |
| | 0118 | No such SOP Class | |

| | | | |
|---------|------|------------------------|-----------------------------------------|
| | 0112 | No such SOP Instance | |
| | 0110 | Processing failure | |
| | 0213 | Resource limitation | |
| | 0211 | Unrecognized operation | |
| Success | 0000 | Success | The progress screen is still displayed. |
| * | * | Any other status code. | Deemed a failure. |

2.3.1.2.7 Real-World Activity Report MPPS

2.3.1.2.7.1 Associated Real-World Activity

The Modality Performed Procedure Step N-CREATE is sent with status of “IN PROGRESS” when start exam is initiated by the operator. The Modality Performed Procedure Step N-SET is sent with status of “COMPLETED” when the exam is completed (all images acquired). If the exam is aborted, a status of “DISCONTINUED” is sent.

2.3.1.2.7.2 Proposed Presentation Context Table

| Presentation Context Table – Proposed by Workstation AE for Activity Modality Worklist | | | | | |
|-----------------------------------------------------------------------------------------------|-------------------------|----------------------------------------------------------------------------------|-----------------------------------------------------------------|-------------|-----------------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 | SCU | None |

2.3.1.2.7.2.1 SOP Specific DICOM Conformance Statement for Modality Performed Procedure Step SOP Class

MPPS will be created at the start of an image acquisition step; an MPPS will be set after the completion/discontinue of a post-processing step when the user selects “Complete”/“Discontinue”.

The Workstation AE includes attributes in the Modality Performed Procedure Step N-CREATE as described in Section 7.4

Following are the status codes that are more specifically processed when receiving an N-CREATE or N-SET response from an MPPS SCP equipment:

| Service Status | Status Code | Further Meaning | Application Behavior When Receiving Status Code |
|-----------------------|--------------------|-------------------------|--------------------------------------------------------|
| Failure | 0119 | Class-instance conflict | Terminates the association and operation. Logs |

| | | | |
|---------|------|-------------------------|------------------------------------------------------------------------------------------|
| | 0210 | Duplicate invocation | failure message along with error comment returned from SCP and displays failure message. |
| | 0106 | Invalid attribute value | |
| | 0117 | Invalid SOP instance | |
| | 0120 | Missing attribute | |
| | 0121 | Missing attribute value | |
| | 0212 | Mistyped argument | |
| | 0105 | No such attributes | |
| | 0118 | No such SOP Class | |
| | 0112 | No such SOP Instance | |
| | 0110 | Processing failure | |
| | 0213 | Resource limitation | |
| | 0211 | Unrecognized operation | |
| Success | 0000 | | The progress screen is disappeared. |
| | 0111 | Duplicate SOP Instance | |
| * | * | Any other status code. | Deemed a failure. |

2.3.1.3 Association Acceptance Policy

If the Remote AE needs to “Query/Retrieve Images”, to the local system then it has to be configured in the Local system to do the same.

When the DICOM Server AE accepts an association for image storage, it will receive any images transmitted on that association and store the images on disk.

It will also respond to queries from Remote AEs by sending matching entries. Any Remote AE can request and receive a list of images on the local database. The Remote AE must be configured in the local config file list of Remote AE for it to be able to retrieve images from DICOM Server AE.

2.3.1.3.1 Real-World Activity Verification SCP

2.3.1.3.1.1 Associated Real-World Activity

An incoming Verification request will cause the Workstation AE to accept the association and respond with a verification response after the Retrieve DICOM Server AE accepts an sub-association for image storage.

2.3.1.3.1.2 Accepted Presentation Context Table

| Presentation Context Table - Accepted by Workstation AE for Activity Receive Images | | | | | |
|--------------------------------------------------------------------------------------------|-------------------|---------------------------|---------------------|-------------|-----------------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Verification SOP Class | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

2.3.1.3.1.3 Presentation Context Acceptance Criterion

The Workstation evaluates each Presentation Context independently, and accepts any Presentation Context that matches an Abstract Syntax for any Real-World Activity.

2.3.1.3.1.4 Transfer Syntax Selection Policies

Within each Presentation Context, the Workstation will accept the first proposed transfer syntax that it also supports for that Abstract Syntax.

2.3.1.3.2 Real-World Activity Local Store SCP

2.3.1.3.2.1 Associated Real-World Activity

The Real-World Activity associated with the Receive Images operation is the storage of the image on the disk drive.

2.3.1.3.2.2 Accepted Presentation Context Table

| Presentation Context Table - Accepted by Workstation AE for Activity Receive Images | | | | | |
|--------------------------------------------------------------------------------------------|------------------------------|---------------------------|---------------------|-------------|-----------------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| X-Ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| X-Ray Radiofluoroscopic Image Storage | 1.2.840.10008.5.1.4.1.1.12.2 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Computed Radiography Image | 1.2.840.10008.5.1.4.1.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |

| | | | | | |
|---------|--|------------------------|---------------------|--|--|
| Storage | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
|---------|--|------------------------|---------------------|--|--|

2.3.1.3.2.2.1 SOP Specific DICOM Conformance Statement for all Storage SOP Classes

The Workstation AE provides Level 2 (FULL) Conformance, and stores all standard and private data elements of received SOP Instances.

Following are the status codes the Application may send back to the SCU Equipment after performing the requested **Storage**:

| Service Status | Status Code | Further Meaning | Status Code Explanation | Related Fields Sent Back to the SCU |
|----------------|-------------|---------------------------|-------------------------|-------------------------------------|
| Failure | A700 | Refused: Out of resources | Out of Disk space | (0000,0902) |
| Success | 0000 | | | None |

2.3.1.3.2.3 Presentation Context Acceptance Criterion

The Workstation evaluates each Presentation Context independently, and accepts any Presentation Context that matches an Abstract Syntax for any Real-World Activity.

2.3.1.3.2.4 Transfer Syntax Selection Policies

Within each Presentation Context, the Workstation will accept the first proposed transfer syntax that it also supports for that Abstract Syntax.

2.3.1.3.3 Real-World Activity- Storage Commit – SCU role

2.3.1.3.3.1 Associated Real-World Activity

For Image Send with Storage Commit, the Real-World Activity of Image Send to a DICOM Store/Commit server alias on the Workstation will cause an association to be initiated from Storage Commitment SCP for the N-EVENT-REPORT.

2.3.1.3.3.2 Accepted Presentation Context Table

| Presentation Context Table - Accepted by AE Elite Workstation for Activity – Storage Commit - SCP Role | | | | | |
|---------------------------------------------------------------------------------------------------------------|----------------------|----------------------------------------------------------------------------------|-----------------------------------------------------------------|-------------|-----------------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 | SCU | None |

2.3.1.3.3.3 Presentation Context Acceptance Criterion

The Workstation evaluates each Presentation Context independently, and accepts any Presentation Context that matches an Abstract Syntax for any Real-World Activity.

2.3.1.3.3.4 Transfer Syntax Selection Policies

Within each Presentation Context, the Workstation will accept the first proposed transfer syntax that it also supports for that Abstract Syntax.

2.4 COMMUNICATION PROFILES

2.4.1 Supported Communication Stacks

The DICOM Upper Layer Protocol is supported using TCP/IP, as specified in DICOM PS3.8.

The TCP/IP stack is inherited from the Linux Operating System.

2.4.2 Physical Media Support

The Workstation is provided with a 10/100 Mb/s or 1Gb/s auto-sensing Ethernet interface. Additional or alternate network interfaces may be available.

2.4.3 Additional Protocol

The Workstation supports DHCP, user could configure DHCP Enable/Disable from System Setup screen. When DHCP is Enabled, then, the network information will be gotten automatically.

Default is Enabled.

2.4.4 IPv4 and IPv6 Support

The Workstation only supports IPv4.

2.5 EXTENSIONS / SPECIALIZATIONS/ PRIVATIZATIONS

Not Applicable.

2.6 CONFIGURATION

2.6.1 AE Title/Presentation Address Mapping

2.6.2 Configurable Parameters

The following fields are configurable:

The default value for these fields is blank unless otherwise specified.

Network

- Use DHCP or static IP (Default –DHCP).
- Local IP Address
- Local Subnet mask
- Default Gateway

- Local Listening Port Number (Available with Retrieve option default is 1004)

Configuration:

- Local AE Title (Default – GEHC OEC)
- Station Name
- Hospital Name

Every Remote DICOM AE SCP: (Store, Storage Commitment, MPPS, Print, Modality Worklist, Query/Retrieve)

- Server Alias
- Remote AE Title
- Remote IP Address
- Remote Port Number

USB Options (Media Exchange)

- Dose Summary (RDSR, SC. Default - SC)
- Remove Patient Information (Default - No).
- Include Viewer (Default - Include)
- Include Overlay(Default - Include)

Additional fields for Remote DICOM Store

- RDSR on Exam Completion or Discontinue(Default – Checked)
- Dose Summary (RDSR, SC. Default - RDSR)
- Store Overlay options (Always if Available, Never. Default - Never)

Additional fields for Remote DICOM Print

- Min Density for Printer
- Max Density for Printer
- Configuration Information
- Border Density
- Empty Density
- Number of Copies
- Print Priority
- Film Destination

- Medium Type
- Film Size
- Format
- Film Orientation

Additional fields for Remote DICOM Query Worklist Filter

- Patient Name
- Patient ID
- Requested Procedure ID
- Accession Number
- Exam Date: Today, Next 2 Days, Next 3 Days, Next 7 Days, Date Range
- Performing Physician Name
- Modality Type (RF, XA, CR or SC. Default - RF)

Additional fields for Remote DICOM Query/Retrieve Filter

- Date: Today, Next 2 Days, Next 3 Days...
- Patient Name
- Patient ID
- Accession Number
- Modality Type (RF, XA, CR or SC. Default - RF)

Note: All configurations must be performed by a GE Field Engineer.

2.7 SUPPORT OF EXTENDED CHARACTER SETS

As a Query SCU, it will similarly accept response items with ISO_IR 100 values of Specific Character Set or empty value of Specific Character Set (absence of (0008,0005) which mean usage of ASCII Chars only).

The Workstation user interface will allow the user to enter characters from the console keyboard that is within ISO 8859-1. Responses with non-compatible values of (0008,0005) will be discarded.

2.8 CODES AND CONTROLLED TERMINOLOGY

2.8.1 Fixed Coded Terminology

The product uses the fixed (non-configurable, non-extensible) coded terminology in Image SOP Instance and X-Ray Radiation Dose Structured Report attributes, as described in Sections 6.4.3 where the VT is CODE. In these cases the standard or extended value for the code is defined .

2.8.2 Mapped Coded Terminology

The product maps, without change, coded terminology values supplied in Modality Worklist Scheduled Procedure Steps into Image SOP Instance and X-Ray Radiation Dose Structured Report attributes, as described in Sections 5.7 for RDSR and 6.4.1.2, 6.4.2, 6.4.3, and 6.4.5 for images.

2.9 SECURITY PROFILES

The Workstation does not conform to any defined DICOM Security Profiles.

It is assumed that the Workstation is used within a secured environment. It is assumed that a secured environment includes at a minimum:

1. Firewall or router protections to ensure that only approved external hosts have network access to the Workstation.
2. Firewall or router protections to ensure that the Workstation only has network access to approved external hosts and services.
3. Any communications with external hosts and services outside the locally secured environment use appropriate secure network channels (such as a Virtual Private Network (VPN))

3. MEDIA STORAGE CONFORMANCE STATEMENT


3.1 INTRODUCTION

This section of the DICOM conformance statement specifies the Workstation AE compliance to DICOM requirements for **Media Interchange**. It details the DICOM Media Storage Application Profiles and roles that are supported by this product.

The Workstation AE is able to export images and RDSR to DICOM media.

3.1.1 Application Data Flow Diagram

The media interchange application model for the Workstation AE is shown in the following Illustration:



The Workstation AE supports the General Purpose USB Interchange profiles.

3.1.2 Functional Definition of AE's

The Workstation AE can perform these functions:

- Create a new DICOM File-Set on Media
- Update a DICOM File-Set by adding new SOP instances to the media which supports multi-session.

3.1.3 Sequencing of Real-World Activities

Storage/Export – Images are acquired and archived on Workstation; From the Image Directory screen, images/RDSR can be selected and copied to media; Images/RDSR can be viewed by any DICOM application that supports Media Exchange FSR and X-Ray Radiation Dose Structured Reports.

Directory Listing – During the Image Copy, user could choice copying image/RDSR with a Media Viewer, then the images/RDSR will be copied with a Media Viewer instance;

Images/RDSR can be viewed by this media viewer only on Windows OS (Independent of the Workstation AE).

3.1.4 File Meta Information Options (See PS3.10)

The File Meta-Information for this implementation is:

| | |
|--------------------------------------------|-----------------------------|
| File Meta-Information Version | 1 |
| OEK One Implementation UID | 1.2.840.113619.6.432 |
| OEK One Implementation Version Name | GEHC_OEK_200 |

3.2 AE SPECIFICATIONS

3.2.1 Workstation AE Specification

The Workstation Application Entity provides standard conformance to DICOM Interchange Option of the Media Storage Service Class. The supported Application Profiles and roles are listed below.

| Supported Application Profile | Real World Activity | Role | Option |
|--------------------------------------|----------------------------|-------------|---------------|
| STD-GEN-USB | Storage/ Export | FSC,FSU | Interchange |

3.2.1.1 File Meta Information for the Workstation Application Entity

Following are the values set in the File Meta Information for this AE Title:

| | |
|--------------------------------------------|-----------------------------|
| Source Application Entity Title | Not used |
| OEK One Implementation UID | 1.2.840.113619.6.432 |
| OEK One Implementation Version Name | GEHC_OEK_200 |

3.2.1.2 Real-World Activities for the Workstation Application Entity

3.2.1.2.1 Real-World Activity Storage

The Workstation AE acts as a FSC/FSU when requested to copy SOP instance from local database to the interchange media.

When user selects “Storage” function, the selected images should be saved as DICOM SOP instances and copied to media. The DICOM File set on the media should be created or updated accordingly.

3.2.1.2.1.1 Media Storage Application Profile for the RWA Update File Set

For the list of Application Profiles that invoke this AE for the Real-World Activity Storage, see the Table in Section 3.2.1 “Workstation AE Specification” where the table describing the Application Profiles and Real-World Activity is defined.

3.2.1.2.1.1.1 Options for STD-GEN-USB Application Profile

Following are the optional SOP Classes supported by this AE. All SOP Instances use the Explicit VR Little Endian Uncompressed Transfer Syntax, UID 1.2.840.10008.1.2.1.

| SOP Class | SOP Class UID |
|----------------------------------|------------------------------|
| X-Ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 |

Common DICOMDIR Directory Records created by this AE will include key attributes as described in Section 9. Following are the Additional DICOMDIR Keys supported for this profile:

ADDITIONAL DICOMDIR KEYS

| Key Attribute | Tag | Directory Record Type | Type | Notes |
|---------------|-------------|-----------------------|------|-------------------------------------------------------------------------------------------------------------------|
| Image Type | (0008,0008) | IMAGE | 1C | FSC/FSU: Present when the DICOMDIR file is created and this value should be extracted from referenced image file. |

3.2.2 Media Viewer AE Specification

The Media Viewer Application Entity provides standard conformance to DICOM Interchange of the Media Storage Service Class. The supported Application Profiles and roles are listed below.

| Supported Application Profile | Real World Activity | Role |
|-------------------------------|---------------------|------|
| STD-GEN-USB | Directory Listing | FSR |

3.2.2.1 File Meta Information for the Media Viewer Application Entity

Following are the values set in the File Meta Information for this AE Title:

| | |
|----------------------------------------|-----------------|
| Source Application Entity Title | Not used |
| Media Viewer Implementation UID | Not used |
| Implementation Version Name | Not used |

Since Media Viewer just used to read file-set and display image, so, above information could not be set by Media Viewer.

3.2.2.2 Real-World Activities for the Media Viewer Application Entity

3.2.2.2.1 Real-World Activity Directory Listing

The Media Viewer AE acts as a FSR when listing the directory information from interchange media.

When user load the media on Window OS, the Media Viewer will read the DICOM File-Set and display/list the DICOM directory entries for those instances in the File-Set on the media.

3.2.2.2.1.1 Media Storage Application Profile for the RWA Directory Listing

For the list of Application Profiles that invoke this AE for the Real-World Activity Directory Listing, see the Table in Section 3.2.3 “Media Viewer AE Specification” where the table describing the Application Profiles and Real-World Activity is defined.

3.2.2.2.1.1.1 STD-GEN-USB Application Profile

Following are the SOP Classes supported by this AE and only the Explicit VR Little Endian Uncompressed Transfer Syntax (UID 1.2.840.10008.1.2.1) is supported.

| SOP Class | SOP Class UID |
|----------------------------------|------------------------------|
| X-Ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 |

Common DICOMDIR Directory Records read by this AE will include key attributes as described in Section 9. Following are the Additional DICOMDIR Keys supported for this profile:

| Key Attribute | Tag | Directory Record Type | Type | Notes |
|---------------|-------------|-----------------------|------|------------|
| Image Type | (0008,0008) | IMAGE | 1C | FSR: Used. |

3.2.3 Media Viewer AE Specification

The Media Viewer Application Entity provides standard conformance to DICOM Interchange of the Media Storage Service Class. The supported Application Profiles and roles are listed below.

| Supported Application Profile | Real World Activity | Role |
|-------------------------------|---------------------|------|
| STD-GEN-USB | Directory Listing | FSR |

3.3 AUGMENTED AND PRIVATE APPLICATION PROFILES

No augmented/private profile is implemented.

3.4 EXTENSIONS / SPECIALIZATIONS / PRIVATIZATIONS

No Extensions / Specializations / Privatizations are implemented.

3.5 CONFIGURATION

Media Options (Media Exchange)

- Dose Summary (RDSR, SC. Default - SC)
- Remove Patient Information (Default - No).
- Include Viewer (Default - Include)

- Include Overlay(Default - Include)

3.6 SUPPORT OF EXTENDED CHARACTER SETS

The Workstation AE will support fully only the ISO_IR 100 as extended character sets for media storage.

4. X-RAY ANGIOGRAPHIC INFORMATION OBJECT IMPLEMENTATION

4.1 INTRODUCTION

This section specifies the use of the DICOM X-Ray Angiographic Image IOD to represent the information included in X-Ray Angiographic Images by this implementation. Corresponding attributes are conveyed using the module construct.

4.2 WORKSTATION MAPPING OF DICOM ENTITIES

The Workstation maps DICOM Information Entities to local Information Entities in the product's database and user interface.

TABLE 4-1
MAPPING OF DICOM ENTITIES TO WORKSTATION ENTITIES

| DICOM IE | Workstation Entity |
|-----------|--------------------|
| Patient | Patient |
| Study | Exam |
| Series | Series |
| Equipment | Equipment |
| Image | Image |

4.3 IOD MODULE TABLE

The X-Ray Angiographic Image Information Object Definition comprises the modules of the following table, plus Standard Extended and Private attributes. Standard Extended and Private attributes are described in Section 4.5.

TABLE 4-2
X-RAY ANGIOGRAPHIC IMAGE IOD MODULE

| Entity Name | Module Name | Usage | Reference |
|--------------------|------------------------|----------|-----------|
| Patient | Patient | Used | 4.4.1.1 |
| | Clinical Trial Subject | Not Used | N/A |
| Study | General Study | Used | 4.4.2.1 |
| | Patient Study | Used | 4.4.2.2 |
| | Clinical Trial Study | Not Used | N/A |
| Series | General Series | Used | 4.4.3.1 |
| | Clinical Trial Series | Not Used | N/A |
| Frame of Reference | Synchronization | Not Used | N/A |
| Equipment | General Equipment | Used | 4.4.4.1 |
| Image | General Image | Used | 4.4.5.1 |

| | | |
|---------------------|----------|----------|
| Image Pixel | Used | 4.4.5.2 |
| Contrast/Bolus | Used | 4.4.5.3 |
| Cine | Used | 4.4.5.4 |
| Multi-Frame | Used | 4.4.5.5 |
| Frame Pointers | Not Used | N/A |
| Mask | Used | 6.4.2 |
| Display Shutter | Not Used | N/A |
| Device | Not Used | N/A |
| Intervention | Not Used | N/A |
| X-Ray Image | Used | 4.4.5.7 |
| X-Ray Acquisition | Used | 4.4.5.8 |
| X-Ray Collimator | Not Used | N/A |
| X-Ray Table | Not Used | N/A |
| XA Positioner | Used | 4.4.5.9 |
| DX Detector | Not Used | N/A |
| Overlay Plane | Used | 4.4.5.10 |
| Multi-Frame Overlay | Not Used | N/A |
| Modality LUT | Not Used | N/A |
| VOI LUT | Used | 4.4.5.11 |
| SOP Common | Used | 4.4.5.12 |
| Frame Extraction | Not Used | N/A |

4.4 INFORMATION MODULE DEFINITIONS

Please refer to DICOM Part 3 (Information Object Definitions) for a description of each of the entities, modules, and attributes contained within the X-Ray Angiographic Information Object.

The following modules are included to convey Enumerated Values, Defined Terms, and Optional Attributes supported and/or expected. Type 1 & Type 2 Attributes are also included for completeness and to define what values they may take and where these values are obtained from when generating the instance as well as what are the expected values when loading such instance. It should be noted that they are the same ones as defined in the DICOM Standard Part 3 (Information Object Definitions). Also note that Attributes not present in tables are not supported.

4.4.1 Patient Entity Modules

4.4.1.1 Patient Module

**TABLE 4-3
PATIENT MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------|-------------|-------------|----------------------------------------------------------------------------------------------------------|
| Patient's Name | (0010,0010) | 2 | May be entered from user interface or got from Worklist. Supports 5 different components delimited by |

**OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | |
|---------------------------------------------------|-------------|----|---------------------------------------------------------------------------------------------------------------------------------------------------|
| | | | “^”. Supports a maximum length of 64 bytes including the delimiter. |
| Patient ID | (0010,0020) | 2 | May be entered from user interface or got from Worklist. Truncated to 64 bytes, end truncated. |
| Issuer of Patient ID | (0010,0021) | 3 | May be entered from user interface or got from Worklist. Truncated to 64 bytes, end truncated. |
| Patient's Birth Date | (0010,0030) | 2 | May be entered from user interface or got from Worklist Truncated to 8 bytes, end truncated |
| Patient's Sex | (0010,0040) | 2 | May be entered from user interface or got from Worklist. Truncated to 2 bytes, end truncated |
| Referenced Patient Sequence | (0008,1120) | 3 | Not Used |
| <i>>Include 'SOP Instance Reference Macro'</i> | | | |
| Patient's Birth Time | (0010,0032) | 3 | Not Used |
| Other Patient IDs | (0010,1000) | 3 | Not Used |
| Other Patient IDs Sequence | (0010,1002) | 3 | Not Used |
| >Patient ID | (0010,0020) | 1 | Not Used |
| >Issuer of Patient ID | (0010,0021) | 1 | Not Used |
| >Type of Patient ID | (0010,0022) | 1 | Not Used |
| Other Patient Names | (0010,1001) | 3 | Not Used |
| Ethnic Group | (0010,2160) | 3 | Not Used |
| Patient Comments | (0010,4000) | 3 | May be entered from user interface or got from Worklist From user interface truncated to 512 bytes. From worklist truncated to 10240 bytes. |
| Patient Species Description | (0010,2201) | 1C | Not Used |
| Patient Species Code Sequence | (0010,2202) | 1C | Not Used |
| <i>>Include 'Code Sequence Macro'</i> | | | |
| Patient Breed Description | (0010,2292) | 2C | Not Used |
| Patient Breed Code Sequence | (0010,2293) | 2C | Not Used |
| <i>>Include 'Code Sequence Macro'</i> | | | |
| Breed Registration Sequence | (0010,2294) | 2C | Not Used |
| >Breed Registration Number | (0010,2295) | 1 | Not Used |
| >Breed Registry Code Sequence | (0010,2296) | 1 | Not Used |
| <i>>>Include 'Code Sequence Macro'</i> | | | |
| Responsible Person | (0010,2297) | 2C | Not Used |
| Responsible Person Role | (0010,2298) | 1C | Not Used |
| Responsible Organization | (0010,2299) | 2C | Not Used |
| Patient Identity Removed | (0012,0062) | 3 | Used. |

**OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | |
|----------------------------------------|-------------|----|------------------------------------------------------------------------------------|
| | | | If Patient Information is removed then set to "YES", otherwise set to "NO" |
| De-identification Method | (0012,0063) | 1C | Used. If Patient Information is removed. Set to "Limited Data Set" in this case |
| De-identification Method Code Sequence | (0012,0064) | 1C | Not Used |
| >Include 'Code Sequence Macro' | | | |

4.4.2 Study Entity Modules

4.4.2.1 General Study Module

**TABLE 4-4
GENERAL STUDY MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|----------------------------------------------------|-------------|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Study Instance UID | (0020,000D) | 1 | The product will adopt the Study Instance UID which is returned in the MWL response. If there is no Study Instance UID returned from MWL server, then new Study Instance UID will be created. |
| Study Date | (0008,0020) | 2 | This value is set to the exam created date. |
| Study Time | (0008,0030) | 2 | This value is set to the exam created time. |
| Referring Physician's Name | (0008,0090) | 2 | May be entered from the user interface or got from Worklist. |
| Referring Physician Identification Sequence | (0008,0096) | 3 | Not Used |
| >Include 'Person Identification Macro' | | | |
| Study ID | (0020,0010) | 2 | Got from Worklist or generated by equipment. |
| Accession Number | (0008,0050) | 2 | May be entered from the user interface or got from Worklist. |
| Study Description | (0008,1030) | 3 | May be entered from the user interface or got from Worklist. The default value will be set as Requested Procedure Description when the tag got from MWL is not empty. Else this value will be set as Procedure Description on user interface. |
| Physician(s) of Record | (0008,1048) | 3 | Not Used |
| Physician(s) of Record Identification Sequence | (0008,1049) | 3 | Not Used |
| >Include 'Person Identification Macro' | | | |
| Name of Physician(s) Reading Study | (0008,1060) | 3 | Not Used |
| Physician(s) Reading Study Identification Sequence | (0008,1062) | 3 | Not Used |
| >Include 'Person Identification Macro' | | | |
| Referenced Study Sequence | (0008,1110) | 3 | Got from MWL. |
| >Include 'SOP Instance Reference Macro' | | | |

| | | | |
|--------------------------------|-------------|---|----------|
| Procedure Code Sequence | (0008,1032) | 3 | Not Used |
| >Include 'Code Sequence Macro' | | | |

4.4.2.2 Patient Study Module

**TABLE 4-5
PATIENT STUDY MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------------------|-------------|------|--------------------------------------------------------------------------------------------------|
| Admitting Diagnoses Description | (0008,1080) | 3 | Not Used |
| Admitting Diagnoses Code Sequence | (0008,1084) | 3 | Not Used |
| >Include 'Code Sequence Macro' | | | |
| Patient's Age | (0010,1010) | 3 | Not Used |
| Patient's Size | (0010,1020) | 3 | May be entered from user interface or got from Worklist. Truncated to 16 bytes, end truncated |
| Patient's Weight | (0010,1030) | 3 | May be entered from user interface or got from Worklist. Truncated to 16 bytes, end truncated |
| Occupation | (0010,2180) | 3 | Not Used |
| Additional Patient's History | (0010,21B0) | 3 | Not Used |
| Admission ID | (0038,0010) | 3 | Not Used |
| Service Episode ID | (0038,0060) | 3 | Not Used |
| Service Episode Description | (0038,0062) | 3 | Not Used |
| Patient's Sex Neutered | (0010,2203) | 2C | Not Used |

4.4.3 Series Entity Modules

4.4.3.1 General Series Module

**TABLE 4-6
GENERAL SERIES MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------------|-------------|------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| Modality | (0008,0060) | 1 | Defined by the user by setting the modality on the DICOM Server or Media Image Save. |
| Series Instance UID | (0020,000E) | 1 | Created by equipment, using prefix: 1.2.840.113619.2.432.+MAC+seriesIndicator(20)+TimeStamp_Exam MAC: Number string of device MAC address. |
| Series Number | (0020,0011) | 2 | A number that identifies this Series. |
| Laterality | (0020,0060) | 2C | The value is empty. |
| Series Date | (0008,0021) | 3 | This value is set to the first shot date. |
| Series Time | (0008,0031) | 3 | This value is set to the first shot time. |
| Performing Physicians' Name | (0008,1050) | 3 | Get from the Worklist scheduled performing physician, or the user can enter from the user interface. |

OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X

GE Healthcare

DICOM CONFORMANCE STATEMENT

2148033 REV 1

| | | | |
|--------------------------------------------------------|-------------|----|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Performing Physician Identification Sequence | (0008,1052) | 3 | Not Used |
| <i>>Include 'Person Identification Macro'</i> | | | |
| Protocol Name | (0018,1030) | 3 | When using MWL, if Scheduled Protocol Sequence Code Meaning is not empty then this value is Scheduled Protocol Sequence Code Meaning, else this value is Scheduled Procedure Step Description. May be entered from user interface or set to "Unknown Protocol" for local created exam. |
| Series Description | (0008,103E) | 3 | The value equal protocol Name. |
| Operators' Name | (0008,1070) | 3 | Not Used |
| Operator Identification Sequence | (0008,1072) | 3 | Not Used |
| <i>>Include 'Person Identification Macro'</i> | | | |
| Referenced Performed Procedure Step Sequence | (0008,1111) | 3 | Used. |
| <i>>Include 'SOP Instance Reference Macro'</i> | | | |
| Related Series Sequence | (0008,1250) | 3 | Not Used |
| >Study Instance UID | (0020,000D) | 1 | Not Used |
| >Series Instance UID | (0020,000E) | 1 | Not Used |
| >Purpose of Reference Code Sequence | (0040,A170) | 2 | Not Used |
| <i>>>Include 'Code Sequence Macro'</i> | | | |
| Body Part Examined | (0018,0015) | 3 | The value is empty. |
| Patient Position | (0018,5100) | 2C | Not Used |
| Smallest Pixel Value in Series | (0028,0108) | 3 | Not Used |
| Largest Pixel Value in Series | (0028,0109) | 3 | Not Used |
| Request Attributes Sequence | (0040,0275) | 3 | Used |
| >Requested Procedure ID | (0040,1001) | 1C | This value could got from MWL response, or the user can enter from the user interface. |
| >Accession Number | (0008,0050) | 3 | Not Used |
| >Study Instance UID | (0020,000D) | 3 | Not Used |
| >Referenced Study Sequence | (0008,1110) | 3 | Not Used |
| <i>>> Include 'SOP Instance Reference Macro'</i> | | | |
| >Requested Procedure Description | (0032,1060) | 3 | This value could got from MWL response, or the user can enter from the user interface. |
| >Requested Procedure Code Sequence | (0032,1064) | 3 | Sequence set with values from the MWL. No user input is used. |
| <i>>>Include 'Code Sequence Macro'</i> | | | |
| >Reason for the Requested Procedure | (0040,1002) | 3 | Not Used |
| >Reason for Requested Procedure Code Sequence | (0040,100A) | 3 | Not Used |
| <i>>>Include 'Code Sequence Macro'</i> | | | |
| >Scheduled Procedure Step ID | (0040,0009) | 1C | This value could got from MWL response, or the |

**OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | |
|------------------------------------------|-------------|---|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| | | | user can enter from the user interface. |
| >Scheduled Procedure Step Description | (0040,0007) | 3 | This value could got from MWL response, or the user can enter from the user interface. |
| >Scheduled Protocol Code Sequence | (0040,0008) | 3 | Sequence set with values from the MWL. No user input is used. |
| >>Include 'Code Sequence Macro' | | | |
| >>Protocol Context Sequence | (0040,0440) | 3 | Not Used |
| >>>Include 'Content Item Macro' | | | |
| >>>Content Item Modifier Sequence | (0040,0441) | 3 | Not Used |
| >>>>Include 'Content Item Macro' | | | |
| Performed Procedure Step ID | (0040,0253) | 3 | Set to the value of 'Scheduled Procedure Step ID' value from MWL. Or, Set to value entered by the user. Or, Set to value generated by the equipment. |
| Performed Procedure Step Start Date | (0040,0244) | 3 | Used |
| Performed Procedure Step Start Time | (0040,0245) | 3 | Used |
| Performed Procedure Step Description | (0040,0254) | 3 | Used |
| Performed Protocol Code Sequence | (0040,0260) | 3 | Not Used |
| >Include 'Code Sequence Macro' | | | |
| >Protocol Context Sequence | (0040,0440) | 3 | Not Used |
| >>Include 'Content Item Macro' | | | |
| >>Content Item Modifier Sequence | (0040,0441) | 3 | Not Used |
| >>>Include 'Content Item Macro' | | | |
| Comments on the Performed Procedure Step | (0040,0280) | 3 | Not Used |

4.4.4 Equipment Entity Modules

4.4.4.1 General Equipment Module

**TABLE 4-7
GENERAL EQUIPMENT MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|-------------------------------|-------------|------|---------------------------------------------------------------------------|
| Manufacturer | (0008,0070) | 2 | This value is set to "GE Hualun Medical Systems, Co. Ltd". |
| Institution Name | (0008,0080) | 3 | Configurable by user on setup screen (The field name is "Hospital Name"). |
| Institution Address | (0008,0081) | 3 | Not Used |
| Station Name | (0008,1010) | 3 | Configurable by user on DICOM setup screen. |
| Institutional Department Name | (0008,1040) | 3 | Not Used |
| Manufacturer's Model Name | (0008,1090) | 3 | This value is set to type and model of the system: "OEC One". |
| Device Serial Number | (0018,1000) | 3 | Configurable during system install. |
| Software Versions | (0018,1020) | 3 | This value is set to workstation software |

| | | | |
|--------------------------|-------------|----|----------|
| | | | version. |
| Gantry ID | (0018,1008) | 3 | Not Used |
| Spatial Resolution | (0018,1050) | 3 | Not Used |
| Date of Last Calibration | (0018,1200) | 3 | Not Used |
| Time of Last Calibration | (0018,1201) | 3 | Not Used |
| Pixel Padding Value | (0028,0120) | 1C | Not Used |

4.4.5 Image Entity Modules

4.4.5.1 General Image Module

**TABLE 4-8
GENERAL IMAGE MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|---------------------------------------------------|--------------|-------------|------------------------------------------------|
| Instance Number | (0020,0013) | 2 | This value is set to image number. |
| Patient Orientation | (0020,0020) | 2C | The value is empty. |
| Content Date | (0008,0023) | 2C | This value is set to the date image was taken. |
| Content Time | (0008,0033) | 2C | This value is set to the time image was taken. |
| Image Type | (0008,0008) | 3 | See 4.4.5.7.1 |
| Acquisition Number | (0020,0012) | 3 | This value is set to shot number. |
| Acquisition Date | (0008,0022) | 3 | Not Used |
| Acquisition Time | (0008,0032) | 3 | Not Used |
| Acquisition DateTime | (0008,002A)) | 3 | This value is set to shot DateTime. |
| Referenced Image Sequence | (0008,1140) | 3 | Not Used |
| <i>>Include 'SOP Instance Reference Macro'</i> | | | |
| >Purpose of Reference Code Sequence | (0040,A170) | 3 | Not Used |
| <i>>>Include 'Code Sequence Macro'</i> | | | |
| Derivation Description | (0008,2111) | 3 | Not Used |
| Derivation Code Sequence | (0008,9215) | 3 | Not Used |
| <i>>Include 'Code Sequence Macro'</i> | | | |
| Source Image Sequence | (0008,2112) | 3 | Not Used |
| <i>>Include 'SOP Instance Reference Macro'</i> | | | |
| >Purpose of Reference Code Sequence | (0040,A170) | 3 | Not Used |
| <i>>>Include 'Code Sequence Macro'</i> | | | |
| >Spatial Locations Preserved | (0028,135A) | 3 | Not Used |
| >Patient Orientation | (0020,0020) | 1C | Not Used |
| Referenced Instance Sequence | (0008,114A) | 3 | Not Used |
| <i>>Include 'SOP Instance Reference Macro'</i> | | | |
| >Purpose of Reference Code Sequence | (0040,A170) | 1 | Not Used |
| <i>>>Include 'Code Sequence Macro'</i> | | | |
| Images in Acquisition | (0020,1002) | 3 | Not Used |
| Image Comments | (0020,4000) | 3 | Not Used |

| | | | |
|----------------------------------------|-------------|---|----------|
| Quality Control Image | (0028,0300) | 3 | Not Used |
| Burned In Annotation | (0028,0301) | 3 | Not Used |
| Lossy Image Compression | (0028,2110) | 3 | Not Used |
| Lossy Image Compression Ratio | (0028,2112) | 3 | Not Used |
| Lossy Image Compression Method | (0028,2114) | 3 | Not Used |
| Icon Image Sequence | (0088,0200) | 3 | Not Used |
| <i>>Include 'Image Pixel Macro'</i> | | | |
| Presentation LUT Shape | (2050,0020) | 3 | Not Used |

4.4.5.2 Image Pixel Module

**TABLE 4-9
IMAGE PIXEL MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|---------------------------------------------|-------------|-------------|----------------------------------------------------------------------------------------|
| Samples per Pixel | (0028,0002) | 1 | The value is set to 1. |
| Photometric Interpretation | (0028,0004) | 1 | The value is set to "MONOCHROME2". |
| Rows | (0028,0010) | 1 | This value is set to 1000. Note: For dose summary snapshot, this value will be 960. |
| Columns | (0028,0011) | 1 | This value is set to 1000. Note: For dose summary snapshot, this value will be 960. |
| Bits Allocated | (0028,0100) | 1 | This value is set to 16. Note: For dose summary snapshot, this value will be 8. |
| Bits Stored | (0028,0101) | 1 | This value is set to 12. Note: For dose summary snapshot, this value will be 8. |
| High Bit | (0028,0102) | 1 | This value is set to 11. Note: For dose summary snapshot, this value will be 7. |
| Pixel Representation | (0028,0103) | 1 | This value is set to 0 for unsigned integer. |
| Pixel Data | (7FE0,0010) | 1 | Pixel data of the image. |
| Planar Configuration | (0028,0006) | 1C | Not Used |
| Pixel Aspect Ratio | (0028,0034) | 1C | Not Used |
| Smallest Image Pixel Value | (0028,0106) | 3 | Not Used |
| Largest Image Pixel Value | (0028,0107) | 3 | Not Used |
| Red Palette Color Lookup Table Descriptor | (0028,1101) | 1C | Not Used |
| Green Palette Color Lookup Table Descriptor | (0028,1102) | 1C | Not Used |
| Blue Palette Color Lookup Table Descriptor | (0028,1103) | 1C | Not Used |
| Red Palette Color Lookup Table Data | (0028,1201) | 1C | Not Used |

| | | | |
|---------------------------------------|-------------|----|----------|
| Green Palette Color Lookup Table Data | (0028,1202) | 1C | Not Used |
| Blue Palette Color Lookup Table Data | (0028,1203) | 1C | Not Used |
| ICC Profile | (0028,2000) | 3 | Not Used |
| Pixel Data Provider URL | (0028,7FE0) | 1C | Not Used |
| Pixel Padding Range Limit | (0028,0121) | 1C | Not Used |

4.4.5.3 Contrast/Bolus Module

**TABLE 4-10
CONTRAST/BOLUS MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|----------------------------------------------|-------------|-------------|----------------------------------|
| Contrast/Bolus Agent | (0018,0010) | 2 | Got from MWL, else set as empty. |
| Contrast/Bolus Agent Sequence | (0018,0012) | 3 | Not Used |
| <i>>Include 'Code Sequence Macro'</i> | | | |
| Contrast/Bolus Route | (0018,1040) | 3 | Not Used |
| Contrast/Bolus Administration Route Sequence | (0018,0014) | 3 | Not Used |
| <i>>Include 'Code Sequence Macro'</i> | | | |
| >Additional Drug Sequence | (0018,002A) | 3 | Not Used |
| <i>>>Include 'Code Sequence Macro'</i> | | | |
| Contrast/Bolus Volume | (0018,1041) | 3 | Not Used |
| Contrast/Bolus Start Time | (0018,1042) | 3 | Not Used |
| Contrast/Bolus Stop Time | (0018,1043) | 3 | Not Used |
| Contrast/Bolus Total Dose | (0018,1044) | 3 | Not Used |
| Contrast Flow Rate(s) | (0018,1046) | 3 | Not Used |
| Contrast Flow Duration(s) | (0018,1047) | 3 | Not Used |
| Contrast/Bolus Ingredient | (0018,1048) | 3 | Not Used |
| Contrast/Bolus Ingredient Concentration | (0018,1049) | 3 | Not Used |

4.4.5.4 Cine

**TABLE 4-11
CINE MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|--------------------------------|-------------|-------------|---------------------------------------------|
| Preferred Playback Sequencing | (0018,1244) | 3 | Value is 0 |
| Frame Time | (0018,1063) | 1C | Nominal time (in msec) per individual frame |
| Frame Time Vector | (0018,1065) | 1C | Not Used |
| Start Trim | (0008,2142) | 3 | Value is 1 |
| Stop Trim | (0008,2143) | 3 | Value equals number of frames |
| Recommended Display Frame Rate | (0008,2144) | 3 | Value equals number of frames per second |
| Cine Rate | (0018,0040) | 3 | Number of frames per second |
| Frame Delay | (0018,1066) | 3 | Not Used |
| Image Trigger Delay | (0018,1067) | 3 | Not Used |

| | | | |
|--------------------------------------------------------|-------------|----|--------------------------------------------------------------------------------------|
| Effective Duration | (0018,0072) | 3 | Total time in seconds that data was actually taken for the entire Multi-frame image. |
| Actual Frame Duration | (0018,1242) | 3 | Elapsed time of data acquisition in msec per each frame. |
| Multiplexed Audio Channels Description Code Sequence | (003A,0300) | 2C | Not Used |
| >Channel Identification Code | (003A,0301) | 1 | Not Used |
| >Channel Mode | (003A,0302) | 1 | Not Used |
| >Channel Source Sequence | (003A,0300) | 1 | Not Used |
| >>Include Table 8.8-1 "Code Sequence Macro Attributes" | | | |

4.4.5.5 Multi-Frame

**TABLE 4-12
MULTI-FRAME MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|-------------------------|-------------|-------------|-----------------------------------------|
| Number of Frames | (0028,0008) | 1 | Number of frames in a Multi-frame Image |
| Frame Increment Pointer | (0028,0009) | 1 | Value is "(0x0018, 0x1063)" |
| Stereo Pairs Present | (0022,0028) | 3 | Not Used |

4.4.5.6 Mask

**TABLE 4-13
MASK MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------------|-------------|-------------|------------------------------|
| Mask Subtraction Sequence | (0028,6100) | 1 | Used |
| >Mask Operation | (0028,6101) | 1 | Value is NONE |
| >Subtraction Item ID | (0028,9416) | 1C | Not Used |
| >Applicable Frame Range | (0028,6102) | 1C | Not Used |
| >Mask Frame Numbers | (0028,6110) | 1C | Not Used |
| >Contrast Frame Averaging | (0028,6112) | 3 | Not Used |
| >Mask Sub-pixel Shift | (0028,6114) | 3 | Not Used |
| >TID Offset | (0028,6120) | 2C | Not Used |
| >Mask Operation Explanation | (0028,6190) | 3 | Not Used |
| >Mask Selection Mode | (0028,9454) | 3 | Not Used |
| Recommended Viewing Mode | (0028,1090) | 2 | Value is NAT |

4.4.5.7 X-Ray Image Module

**TABLE 4-14
X-RAY IMAGE MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|------------------------------|-------------|-------------|------------------------------|
| Lossy Image Compression | (0028,2110) | 1C | Not Used |
| Image Type | (0008,0008) | 1 | See 4.4.5.7.1 |
| Pixel Intensity Relationship | (0028,1040) | 1 | Value is DISP |
| Samples per Pixel | (0028,0002) | 1 | Value is 1 |

| | | | |
|-----------------------------------------------|--------------|----|-----------------------------------------------------------------------|
| Photometric Interpretation | (0028,0004) | 1 | Value is MONOCHROME2 |
| Bits Allocated | (0028,0100) | 1 | Value is 16 Note: For dose summary snapshot, this value will be 8. |
| Bits Stored | (0028,0101) | 1 | Value is 12 Note: For dose summary snapshot, this value will be 8. |
| High Bit | (0028,0102) | 1 | Value is 11 Note: For dose summary snapshot, this value will be 7. |
| Pixel Representation | (0028, 0103) | 1 | Value is 0000H |
| Scan Options | (0018,0022) | 3 | Not Used. |
| Anatomic Region Sequence | (0008,2218) | 3 | Not Used. |
| > <i>Include 'Code Sequence Macro'</i> | | | |
| >Anatomic Region Modifier Sequence | (0008,2220) | 3 | Not Used |
| >> <i>Include 'Code Sequence Macro'</i> | | | |
| Primary Anatomic Structure Sequence | (0008,2228) | 3 | Not Used |
| > <i>Include 'Code Sequence Macro'</i> | | | |
| >Primary Anatomic Structure Modifier Sequence | (0008,2230) | 3 | Not Used |
| >> <i>Include 'Code Sequence Macro'</i> | | | |
| R Wave Pointer | (0028,6040) | 3 | Not Used |
| Reference Image Sequence | (0008,1140) | 1C | Not Used |
| >Referenced SOP Class UID | (0008,1150) | 1 | Not Used |
| >Referenced SOP Class UID | (0008,1155) | 1 | Not Used |
| >Referenced Frame Number | (0008,1160) | 1C | Not Used |
| >Referenced Segment Number | (0062,000B) | 1C | Not Used |
| >Purpose of Reference Code Sequence | (0040,A170) | 3 | Not Used |
| >> <i>Include 'Code Sequence Macro'</i> | | | |
| Derivation Description | (0008,2111) | 3 | Not Used |
| Acquisition Device Processing Description | (0018,1400) | 3 | Not Used |
| Frame Label Vector | (0018,2002) | 3 | Not Used |
| Frame Dimension Pointer | (0028,000A) | 3 | Not Used |
| Calibration Image | (0050,0004) | 3 | Not Used |

4.4.5.7.1 Image Type

Define the values of Image Type (0008,0008) that may be sent and under what circumstances, or refer to the more general description in the General Image Module.

Value 1 shall have the following Enumerated Values:

- ORIGINAL identifies an Original Image

- DERIVED identifies a Derived Image.
- Value 2 shall have the following Enumerated Value:
- PRIMARY identifies a Primary Image
 - SECONDARY identifies a Secondary Image.

Specify which Enumerated Values of Value 3 are created/supported:

- SINGLE PLANE

Note: This value will always is “ORIGINAL\PRIMARY\SINGLE PLANE” in the system.

4.4.5.8 X-Ray Acquisition Module

**TABLE 4-15
X-RAY ACQUISITION MODULE**

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------------------------|-------------|-------------|------------------------------------------------------------------------------------------------------------|
| KV | (0018,0060) | 2 | KVP of the shot. |
| Radiation Setting | (0018,1155) | 1 | All X-Ray modes send 'SC' except for Digital Spot which sends 'GR'. |
| X-Ray Tube Current | (0018,1151) | 2C | Used |
| X-Ray Tube Current in μ A | (0018,8151) | 3 | Used |
| Exposure Time | (0018,1150) | 2C | mSec of the shot. |
| Exposure Time in microS | (0018,8150) | 3 | Not Used |
| Exposure | (0018,1152) | 2C | Not Used |
| Exposure in microAs | (0018,1153) | 3 | Not Used |
| Grid | (0018,1166) | 3 | Not Used |
| Average Pulse Width | (0018,1154) | 3 | Only used when 0018,115A is PULSED. Pulse Width = 49.0ms |
| Radiation Mode | (0018,115A) | 3 | The value reported will be 'CONTINUOUS' for all X-Ray modes except pulsed mode which will report 'PULSED'. |
| Type of Filters | (0018,1161) | 3 | Not Used |
| Intensifier Size | (0018,1162) | 3 | Not Used |
| Field of View Shape | (0018,1147) | 3 | Not Used |
| Field of View Dimension(s) | (0018,1149) | 3 | Not Used |
| Imager Pixel Spacing | (0018,1164) | 3 | Not Used |
| Pixel Spacing | (0028,0030) | 1C | Not Used |
| Pixel Spacing Calibration Type | (0028,0A02) | 3 | Not Used |
| Pixel Spacing Calibration Description | (0028,0A04) | 1C | Not Used |
| Focal Spot | (0018,1190) | 3 | Not Used |
| Image and Fluoroscopy Area Dose Product | (0018,115E) | 3 | DAP of the shot. |

4.4.5.9 XA Positioner Module

**TABLE 4-16
XA POSITIONER MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|---------------------------------------------|-------------|-------------|------------------------------|
| Distance Source to Patient | (0018,1111) | 3 | Not Used |
| Distance Source to Detector | (0018,1110) | 3 | 1000mm. |
| Estimated Radiographic Magnification Factor | (0018,1114) | 3 | Not Used |
| Positioner Motion | (0018,1500) | 2C | Value is "STATIC". |
| Positioner Primary Angle | (0018,1510) | 2 | Send with empty. |
| Positioner Secondary Angle | (0018,1511) | 2 | Send with empty. |
| Positioner Primary Angle Increment | (0018,1520) | 2C | Not Used |
| Positioner Secondary Angle Increment | (0018,1521) | 2C | Not Used |
| Detector Primary Angle | (0018,1530) | 3 | Not Used |
| Detector Secondary Angle | (0018,1531) | 3 | Not Used |

4.4.5.10 Overlay Plane Module

The Workstation supports up to 1 overlay in an image.

**TABLE 4-17
OVERLAY PLANE MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|------------------------|-------------|-------------|--------------------------------------------------|
| Overlay Rows | (60xx,0010) | 1 | Value is 1000. |
| Overlay Columns | (60xx,0011) | 1 | Value is 1000. |
| Overlay Type | (60xx,0040) | 1 | Value is G. |
| Overlay Origin | (60xx,0050) | 1 | Value is 0\0. |
| Overlay Bits Allocated | (60xx,0100) | 1 | 1 |
| Overlay Bit Position | (60xx,0102) | 1 | 0 |
| Overlay Data | (60xx,3000) | 1 | Overlay pixel data |
| Overlay Description | (60xx,0022) | 3 | Value is "Annotation/Measurement Graphic Layer". |
| Overlay Subtype | (60xx,0045) | 3 | Not Used |
| Overlay Label | (60xx,1500) | 3 | Value is "Annotation/Measurement Graphic Layer". |
| ROI Area | (60xx,1301) | 3 | Not Used |
| ROI Mean | (60xx,1302) | 3 | Not Used |
| ROI Standard Deviation | (60xx,1303) | 3 | Not Used |

4.4.5.11 VOI LUT module

**TABLE 4-18
VOI LUT MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------|-------------|-------------|------------------------------|
| VOI LUT Sequence | (0028,3010) | 1C | Not Used |

| | | | |
|-----------------------------------|-------------|----|---------------------------------------------------------------------------------------------|
| >LUT Descriptor | (0028,3002) | 1 | Not Used |
| >LUT Explanation | (0028,3003) | 3 | Not Used |
| >LUT Data | (0028,3006) | 1 | Not Used |
| Window Center | (0028,1050) | 1C | Only one value is present. Set to $2^{\wedge}bitsStored/2$ when export to DICOM images |
| Window Width | (0028,1051) | 1C | Only one value is present. Set to $2^{\wedge}bitsStored - 1$ when export to DICOM images |
| Window Center & Width Explanation | (0028,1055) | 3 | Not Used |

4.4.5.12 SOP Common Module

**TABLE 4-19
SOP COMMON MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------------------------|-------------|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| SOP Class UID | (0008,0016) | 1 | SOP Class UID for the Class that the dataset contains, just XA UID. |
| SOP Instance UID | (0008,0018) | 1 | Created by equipment, using prefix: 1.2.840.113619.2.432.+MAC+imageIndicator(40)+TimeStamp_Exam MAC: Number string of device MAC address. |
| Specific Character Set | (0008,0005) | 1C | Always set to ISO_IR 100(= Latin Alphabet No. 1) |
| Instance Creation Date | (0008,0012) | 3 | Not Used |
| Instance Creation Time | (0008,0013) | 3 | Not Used |
| Instance Creator UID | (0008,0014) | 3 | Not Used |
| Related General SOP Class UID | (0008,001A) | 3 | Not Used |
| Original Specialized SOP Class UID | (0008,001B) | 3 | Not Used |
| Coding Scheme Identification Sequence | (0008,0110) | 3 | Not Used |
| >Coding Scheme Designator | (0008,0102) | 1 | Not Used |
| >Coding Scheme Registry | (0008,0112) | 1C | Not Used |
| >Coding Scheme UID | (0008,010C) | 1C | Not Used |
| >Coding Scheme External ID | (0008,0114) | 2C | Not Used |
| >Coding Scheme Name | (0008,0115) | 3 | Not Used |
| >Coding Scheme Version | (0008,0103) | 3 | Not Used |
| >Coding Scheme Responsible Organization | (0008,0116) | 3 | Not Used |
| Timezone Offset From UTC | (0008,0201) | 3 | Not Used |
| Contributing Equipment Sequence | (0018,A001) | 3 | Not Used |
| >Purpose of Reference Code Sequence | (0040,A170) | 1 | Not Used |
| >>Include 'Code Sequence Macro' | | | |

OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X

GE Healthcare

DICOM CONFORMANCE STATEMENT

2148033 REV 1

| | | | |
|----------------------------------------------|-------------|----|---------------|
| >Manufacturer | (0008,0070) | 1 | Not Used |
| >Institution Name | (0008,0080) | 3 | Not Used |
| >Institution Address | (0008,0081) | 3 | Not Used |
| >Station Name | (0008,1010) | 3 | Not Used |
| >Institutional Department Name | (0008,1040) | 3 | Not Used |
| >Manufacturer's Model Name | (0008,1090) | 3 | Not Used |
| >Device Serial Number | (0018,1000) | 3 | Not Used |
| >Software Versions | (0018,1020) | 3 | Not Used |
| >Spatial Resolution | (0018,1050) | 3 | Not Used |
| >Date of Last Calibration | (0018,1200) | 3 | Not Used |
| >Time of Last Calibration | (0018,1201) | 3 | Not Used |
| >Contribution DateTime | (0018,A002) | 3 | Not Used |
| >Contribution Description | (0018,A003) | 3 | Not Used |
| Instance Number | (0020,0013) | 3 | Image number. |
| SOP Instance Status | (0100,0410) | 3 | Not Used |
| SOP Authorization Date and Time | (0100,0420) | 3 | Not Used |
| SOP Authorization Comment | (0100,0424) | 3 | Not Used |
| Authorization Equipment Certification Number | (0100,0426) | 3 | Not Used |
| MAC Parameters Sequence | (4FFE,0001) | 3 | Not Used |
| >MAC ID Number | (0400,0005) | 1 | Not Used |
| >MAC Calculation Transfer Syntax UID | (0400,0010) | 1 | Not Used |
| >MAC Algorithm | (0400,0015) | 1 | Not Used |
| >Data Elements Signed | (0400,0020) | 1 | Not Used |
| Digital Signatures Sequence | (FFFA,FFFA) | 3 | Not Used |
| >MAC ID Number | (0400,0005) | 1 | Not Used |
| >Digital Signature UID | (0400,0100) | 1 | Not Used |
| >Digital Signature DateTime | (0400,0105) | 1 | Not Used |
| >Certificate Type | (0400,0110) | 1 | Not Used |
| >Certificate of Signer | (0400,0115) | 1 | Not Used |
| >Signature | (0400,0120) | 1 | Not Used |
| >Certified Timestamp Type | (0400,0305) | 1C | Not Used |
| >Certified Timestamp | (0400,0310) | 3 | Not Used |
| >Digital Signature Purpose Code Sequence | (0400,0401) | 3 | Not Used |
| >>Include 'Code Sequence Macro' | | | |
| Encrypted Attributes Sequence | (0400,0500) | 1C | Not Used |
| >Encrypted Content Transfer Syntax UID | (0400,0510) | 1 | Not Used |
| >Encrypted Content | (0400,0520) | 1 | Not Used |

OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X**GE Healthcare****DICOM CONFORMANCE STATEMENT**

2148033 REV 1

| | | | |
|-----------------------------------------------------------------------------------------------------------------------|-------------|----|------------------------------------------------------------------------|
| Original Attributes Sequence | (0400,0561) | 3 | Not Used |
| >Source of Previous Values | (0400,0564) | 2 | Not Used |
| >Attribute Modification DateTime | (0400,0562) | 1 | Not Used |
| >Modifying System | (0400,0563) | 1 | Not Used |
| >Reason for the Attribute Modification | (0400,0565) | 1 | Not Used |
| >Modified Attributes Sequence | (0400,0550) | 1 | Not Used |
| >>Any Attribute from the main data set that was modified or removed; may include Sequence Attributes and their Items. | | | |
| HL7 Structured Document Reference Sequence | (0040,A390) | 1C | Not Used |
| >Referenced SOP Class UID | (0008,1150) | 1 | Not Used |
| >Referenced SOP Instance UID | (0008,1155) | 1 | Not Used |
| >HL7 Instance Identifier | (0040,E001) | 1 | Not Used |
| >Retrieve URI | (0040,E010) | 3 | Not Used |
| Longitudinal Temporal Information Modified | (0028,0303) | 3 | Used. If Include Remove Patient Information then fix to "MODIFIED". |

4.5 STANDARD EXTENDED AND PRIVATE DATA ATTRIBUTES

Not applicable

4.5.1 Standard Attributes

Not applicable

4.6 STANDARD EXTENDED AND PRIVATE CONTEXT GROUPS

Not applicable

5. X-RAY RADIATION DOSE STRUCTURED REPORT INFORMATION OBJECT IMPLEMENTATION

5.1 INTRODUCTION

This section specifies the use of the DICOM X-Ray Radiation Dose SR IOD to represent results produced and/or received by this implementation. Corresponding attributes are conveyed using the module construct.

5.2 WORKSTATION MAPPING OF DICOM ENTITIES

The Workstation maps DICOM Information Entities to local Information Entities in the product's database and user interface.

TABLE 5-1
MAPPING OF DICOM ENTITIES TO WORKSTATION ENTITIES

| DICOM IE | Workstation Entity |
|-----------|--------------------|
| Patient | Patient |
| Study | Exam |
| Series | Series |
| Equipment | Equipment |
| Document | Document |

5.3 IOD MODULE TABLE

The X-Ray Radiation Dose Structured Report Information Object Definitions comprise the modules of the following tables, plus Standard Extended and Private Attributes. Standard Extended and Private Attributes are described in Section 5.5

The contents of the SR Document Content are constrained by the supported template, as identified in Section 5.4.5.2.1.1. Standard, Standard Extended and Private Templates are further described in Section 5.7.

TABLE 5-2
STRUCTURE REPORT IOD MODULES

| Entity Name | Module Name | Usage | Reference |
|-------------|-------------------------|----------|-----------|
| Patient | Patient | Used | 5.4.1.1 |
| | Specimen Identification | Not Used | N/A |
| | Clinical Trial Subject | Not Used | N/A |
| Study | General Study | Used | 5.4.2.1 |
| | Patient Study | Used | 5.4.2.2 |
| | Clinical Trial Study | Not Used | N/A |
| Series | SR Document Series | Used | 5.4.3.1 |
| | Clinical Trial Series | Not Used | N/A |

| | | | |
|--------------------|----------------------------|----------|---------------------------------------------------|
| Frame Of Reference | Synchronization | Not Used | N/A |
| Equipment | General Equipment | Used | 5.4.4.1 |
| | Enhanced General Equipment | Used | 5.4.4.2 |
| Document | SR Document General | Used | 5.4.5.1 |
| | SR Document Content | Used | Error! Reference source not found. |
| | SOP Common | Used | 5.4.5.3 |

5.4 INFORMATION MODULE DEFINITIONS

Please refer to DICOM Part 3 (Information Object Definitions) for a description of each of the entities, modules, and attributes contained within the SR Information Objects.

The following modules are included to convey Enumerated Values, Defined Terms, and Optional Attributes supported and/or expected. Type 1 & Type 2 Attributes are also included for completeness and to define what values they may take and where these values are obtained from when generating the instance as well as what are the expected values when loading such instance. It should be noted that they are the same ones as defined in the DICOM Standard Part 3 (Information Object Definitions). Also note that Attributes not present in tables are not supported.

5.4.1 Patient Entity Modules

5.4.1.1 Patient Module

**TABLE 5-3
PATIENT MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------|-------------|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Patient's Name | (0010,0010) | 2 | May be entered from user interface or got from Worklist. Supports 5 different components delimited by “^”. Supports a maximum length of 64 bytes including the delimiter. |
| Patient ID | (0010,0020) | 2 | May be entered from user interface or got from Worklist. Truncated to 64 bytes, end truncated. |
| Issuer of Patient ID | (0010,0021) | 3 | Got from Worklist. Truncated to 64 bytes, end truncated. |
| Patient's Birth Date | (0010,0030) | 2 | May be entered from user interface or got from Worklist Truncated to 8 bytes, end truncated |
| Patient's Sex | (0010,0040) | 2 | May be entered from user interface or got from Worklist. Truncated to 2 bytes, end truncated |

OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X

GE Healthcare

DICOM CONFORMANCE STATEMENT

2148033 REV 1

| | | | |
|-----------------------------------------|-------------|----|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Referenced Patient Sequence | (0008,1120) | 3 | Not Used |
| >Include 'SOP Instance Reference Macro' | | | |
| Patient's Birth Time | (0010,0032) | 3 | Not Used |
| Other Patient IDs | (0010,1000) | 3 | Not Used |
| Other Patient IDs Sequence | (0010,1002) | 3 | Not Used |
| >Patient ID | (0010,0020) | 1 | Not Used |
| >Issuer of Patient ID | (0010,0021) | 1 | Not Used |
| >Type of Patient ID | (0010,0022) | 1 | Not Used |
| Other Patient Names | (0010,1001) | 3 | Not Used |
| Ethnic Group | (0010,2160) | 3 | Not Used |
| Patient Comments | (0010,4000) | 3 | May be entered from user interface or got from Worklist From user interface truncated to 512 bytes. From worklist truncated to 10240 bytes. |
| Patient Species Description | (0010,2201) | 1C | Not Used |
| Patient Species Code Sequence | (0010,2202) | 1C | Not Used |
| >Include 'Code Sequence Macro' | | | |
| Patient Breed Description | (0010,2292) | 2C | Not Used |
| Patient Breed Code Sequence | (0010,2293) | 2C | Not Used |
| >Include 'Code Sequence Macro' | | | |
| Breed Registration Sequence | (0010,2294) | 2C | Not Used |
| >Breed Registration Number | (0010,2295) | 1 | Not Used |
| >Breed Registry Code Sequence | (0010,2296) | 1 | Not Used |
| >>Include 'Code Sequence Macro' | | | |
| Responsible Person | (0010,2297) | 2C | Not Used |
| Responsible Person Role | (0010,2298) | 1C | Not Used |
| Responsible Organization | (0010,2299) | 2C | Not Used |
| Patient Identity Removed | (0012,0062) | 3 | Used. If Patient Information is removed then set to "YES", otherwise set to "NO" |
| De-identification Method | (0012,0063) | 1C | Used. If Patient Information is removed. Set to "Limited Data Set" in this case |
| De-identification Method Code Sequence | (0012,0064) | 1C | Not Used |
| >Include 'Code Sequence Macro' | | | |

5.4.2 Study Entity Modules

5.4.2.1 General Study Module

**TABLE 5-4
GENERAL STUDY MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|----------------------------------------------------|-------------|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Study Instance UID | (0020,000D) | 1 | The product will adopt the Study Instance UID which is returned in the MWL response. If there is no Study Instance UID returned from MWL server, then new Study Instance UID will be created. |
| Study Date | (0008,0020) | 2 | This value is set to the exam created date. |
| Study Time | (0008,0030) | 2 | This value is set to the exam created time. |
| Referring Physician's Name | (0008,0090) | 2 | Got from Worklist. |
| Referring Physician Identification Sequence | (0008,0096) | 3 | Not Used |
| <i>>Include 'Person Identification Macro'</i> | | | |
| Study ID | (0020,0010) | 2 | Got from Worklist or generated by equipment. |
| Accession Number | (0008,0050) | 2 | May be entered from the user interface or got from Worklist. |
| Study Description | (0008,1030) | 3 | May be entered from the user interface or got from Worklist. The default value will be set as Requested Procedure Description when the tag got from MWL is not empty. Else this value will be set as Procedure Description on user interface. |
| Physician(s) of Record | (0008,1048) | 3 | Not Used |
| Physician(s) of Record Identification Sequence | (0008,1049) | 3 | Not Used |
| <i>>Include 'Person Identification Macro'</i> | | | |
| Name of Physician(s) Reading Study | (0008,1060) | 3 | Not Used |
| Physician(s) Reading Study Identification Sequence | (0008,1062) | 3 | Not Used |
| <i>>Include 'Person Identification Macro'</i> | | | |
| Referenced Study Sequence | (0008,1110) | 3 | Got from Worklist. |
| <i>>Include 'SOP Instance Reference Macro'</i> | | | |
| Procedure Code Sequence | (0008,1032) | 3 | Not Used |
| <i>>Include 'Code Sequence Macro'</i> | | | |

5.4.2.2 Patient Study Module

**TABLE 5-5
PATIENT STUDY MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|---------------------------------|-------------|-------------|------------------------------|
| Admitting Diagnoses Description | (0008,1080) | 3 | Not Used |

| | | | |
|-----------------------------------|-------------|----|--------------------------------------------------------------|
| Admitting Diagnoses Code Sequence | (0008,1084) | 3 | Not Used |
| >Include 'Code Sequence Macro' | | | |
| Patient's Age | (0010,1010) | 3 | Not Used |
| Patient's Size | (0010,1020) | 3 | May be entered from the user interface or got from Worklist. |
| Patient's Weight | (0010,1030) | 3 | May be entered from the user interface or got from Worklist. |
| Occupation | (0010,2180) | 3 | Not Used |
| Additional Patient's History | (0010,21B0) | 3 | Not Used |
| Admission ID | (0038,0010) | 3 | Not Used |
| Service Episode ID | (0038,0060) | 3 | Not Used |
| Service Episode Description | (0038,0062) | 3 | Not Used |
| Patient's Sex Neutered | (0010,2203) | 2C | Not Used |

5.4.3 Series Entity Modules

5.4.3.1 SR Document Series Module

**TABLE 5-6
SR DOCUMENT SERIES MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Use |
|----------------------------------------------|-------------|-------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| Modality | (0008,0060) | 1 | Value = SR |
| Series Instance UID | (0020,000E) | 1 | Created by equipment, using prefix: 1.2.840.113619.2.432.+MAC+seriesIndicator(21)+TimeStamp_Exam MAC: Number string of device MAC address. |
| Series Number | (0020,0011) | 1 | A number that identifies this Series. This number will be set to 1 for image series and set to 2 for RDSR series |
| Series Date | (0008,0021) | 3 | First shot date. |
| Series Time | (0008,0031) | 3 | First shot time. |
| Series Description | (0008,103E) | 3 | The value equal protocol Name. |
| Referenced Performed Procedure Step Sequence | (0008,1111) | 2 | Always send. Referenced SOP Class UID = MPPS SOP Class UID. Referenced SOP Instance UID = MPPS SOP Instance UID. |
| > 'Referenced SOP Class / Instance UIDs' | | | |

5.4.4 Equipment Entity Modules

5.4.4.1 General Equipment Module

**TABLE 5-7
GENERAL EQUIPMENT MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|-------------------------------|-------------|------|---------------------------------------------------------------------------|
| Manufacturer | (0008,0070) | 2 | This value is set to "GE Hualun Medical Systems, Co. Ltd". |
| Institution Name | (0008,0080) | 3 | Configurable by user on setup screen (The field name is "Hospital Name"). |
| Institution Address | (0008,0081) | 3 | Not Used |
| Station Name | (0008,1010) | 3 | Configurable by user on DICOM setup screen. |
| Institutional Department Name | (0008,1040) | 3 | Not Used |
| Manufacturer's Model Name | (0008,1090) | 3 | This value is set to type and model of the system: "OEC One". |
| Device Serial Number | (0018,1000) | 3 | Configurable during system installation. |
| Software Versions | (0018,1020) | 3 | This value is set to workstation software version. |
| Gantry ID | (0018,1008) | 3 | Not Used |
| Spatial Resolution | (0018,1050) | 3 | Not Used |
| Date of Last Calibration | (0018,1200) | 3 | Not Used |
| Time of Last Calibration | (0018,1201) | 3 | Not Used |
| Pixel Padding Value | (0028,0120) | 1C | Not Used |

5.4.4.2 Enhanced General Equipment Module

**TABLE 5-8
ENHANCED GENERAL EQUIPMENT MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|---------------------------|-------------|------|---------------------------------------------------------------|
| Manufacturer | (0008,0070) | 1 | This value is set to "GE Hualun Medical Systems, Co. Ltd". |
| Manufacturer's Model Name | (0008,1090) | 1 | This value is set to type and model of the system: "OEC One". |
| Device Serial Number | (0018,1000) | 1 | Configurable during system installation. |
| Software Versions | (0018,1020) | 1 | This value is set to workstation software version. |

5.4.5 Document Entity Modules

5.4.5.1 SR Document General Module

**TABLE 5-9
SR DOCUMENT GENERAL MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Use |
|-----------------|-------------|------|-------------|
| Instance Number | (0020,0013) | 1 | Value is 1. |
| Completion Flag | (0040,A491) | 1 | COMPLETE |

OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X

GE Healthcare

DICOM CONFORMANCE STATEMENT

2148033 REV 1

| | | | |
|------------------------------------------------------|-------------|----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Completion Flag Description | (0040,A492) | 3 | Not Used |
| Verification Flag | (0040,A493) | 1 | UNVERIFIED |
| Content Date | (0008,0023) | 1 | First shot date. |
| Content Time | (0008,0033) | 1 | First shot time. |
| Verifying Observer Sequence | (0040,A073) | 1C | Not Used. |
| >Verifying Observer Name | (0040,A075) | 1 | Not Used |
| >Verifying Observer Identification Code Sequence | (0040,A088) | 2 | Not Used |
| >>Include 'Code Sequence Macro' | | | |
| >Verifying Organization | (0040,A027) | 1 | Not Used |
| >Verification DateTime | (0040,A030) | 1 | Not Used |
| Author Observer Sequence | (0040,A078) | 3 | Used |
| >Include 'Identified Person or Device Macro' | | | |
| Participant Sequence | (0040,A07A) | 3 | Not Used |
| >Participation Type | (0040,A080) | 1 | Not Used |
| >Participation DateTime | (0040,A082) | 2 | Not Used |
| >Include 'Identified Person or Device Macro' | | | |
| Custodial Organization Sequence | (0040,A07C) | 3 | Not Used |
| >Institution Name | (0008,0080) | 2 | Not Used |
| >Institution Code Sequence | (0008,0082) | 2 | Not Used |
| >>Include 'Code Sequence Macro' | | | |
| Predecessor Documents Sequence | (0040,A360) | 1C | Not Used |
| >Include 'Hierarchical SOP Instance Reference Macro' | | | |
| Identical Documents Sequence | (0040,A525) | 1C | Not Used |
| >Include 'Hierarchical SOP Instance Reference Macro' | | | |
| Referenced Request Sequence | (0040,A370) | 1C | Used |
| >Study Instance UID | (0020,000D) | 1 | The product will adopt the Study Instance UID which is returned in the MWL response. If there is no Study Instance UID returned from MWL server, then new Study Instance UID will be created by system during exam creation. |
| >Referenced Study Sequence | (0008,1110) | 2 | Set as empty |
| >>Include 'SOP Instance Reference Macro' | | | |
| >Accession Number | (0008,0050) | 2 | The product will adopt the Accession Number which is returned in the MWL response. If there is no Accession Number returned from MWL server, then the value will keep empty. For local created exam, user could input value from screen during exam creation. |
| >Placer Order Number/Imaging Service Request | (0040,2016) | 2 | Value is empty. |
| >Filler Order Number/Imaging Service Request | (0040,2017) | 2 | Value is empty. |

| | | | |
|-------------------------------------------------------|-------------|----|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| >Requested Procedure ID | (0040,1001) | 2 | This value could got from MWL response, or the user can enter from the user interface. |
| >Requested Procedure Description | (0032,1060) | 2 | This value could got from MWL response, or the user can enter from the user interface. |
| >Requested Procedure Code Sequence | (0032,1064) | 2 | Sequence set with values from the MWL. No user input is used. |
| >>Include 'Code Sequence Macro' | | | |
| >Reason for the Requested Procedure | (0040,1002) | 3 | Not Used |
| >Reason for Requested Procedure Code Sequence | (0040,100A) | 3 | Not Used |
| >>Include 'Code Sequence Macro' | | | |
| Performed Procedure Code Sequence | (0040,A372) | 2 | Empty. |
| >Include Code Sequence Macro | | | |
| Current Requested Procedure Evidence Sequence | (0040,A375) | 1C | With the value(s) from the system of all saved images for the exam. This is the list of all images whether they have been stored to PACS or not. |
| >Include 'Hierarchical SOP Instance Reference Macro' | | | |
| Pertinent Other Evidence Sequence | (0040,A385) | 1C | Not Used |
| > Include 'Hierarchical SOP Instance Reference Macro' | | | |
| Referenced Instance Sequence | (0008,114A) | 1C | Not Used |
| >Include 'SOP Instance Reference Macro' | | | |
| >Purpose of Reference Code Sequence | (0040,A170) | 1 | Not Used |
| >>Include 'Code Sequence Macro' | | | |

5.4.5.2 SR Document Content Module

**TABLE 5-10
SR DOCUMENT CONTENT MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Use |
|--------------------------------------|-------------|-------------------------------------------|-----------------|
| Observation DateTime | (0040,A032) | 1C | Not Used |
| Content Template Sequence | (0040,A504) | 1C | Not Used |
| >Mapping Resource | (0008,0105) | 1 | Not Used |
| >Template Identifier | (0040,DB00) | 1 | Not Used |
| Value Type | (0040,A040) | 1 | CONTAINER |
| Continuity of Content | (0040,A050) | 1C | SEPARATE |
| Concept Name Code Sequence | (0040,A043) | 1C | Used |
| >Include 'Code Sequence Macro' | | | |
| Insert Concept Value attribute(s) | | (113701,DCM, X-Ray Radiation Dose Report) | |
| Content Sequence | (0040,A730) | 1C | Used |
| > Relationship Type | (0040,A010) | 1 | HAS CONCEPT MOD |
| > Referenced Content Item Identifier | (0040,DB73) | 1C | Not Used |

| | |
|-------------------------------------------|-------------------------------------------------------------------------------------------------------------------|
| > <i>Insert SR DocumentContent Module</i> | Recursive inclusion to create document content tree. See section 5.4.5.2.1 for the list of supported templates |
|-------------------------------------------|-------------------------------------------------------------------------------------------------------------------|

5.4.5.2.1 SR Document Content Descriptions

5.4.5.2.1.1 Content Template

The Workstation supports the following root Templates for SR SOP Instances created, processed, or displayed by the product.

**TABLE 5-11
SR ROOT TEMPLATES**

| SOP Class | Template ID | Template Name | Use |
|-------------------------|-------------|----------------------|----------------|
| X-Ray Radiation Dose SR | 10001 | X-Ray Radiation Dose | Create/Display |

Refer to section 5.7 for a detailed description of the supported templates.

5.4.5.3 SOP Common Module

**TABLE 5-12
SOP COMMON MODULE ATTRIBUTES**

| Attribute Name | Tag | Type | Attribute Description |
|---------------------------------------|-------------|------|------------------------------------------------------------------------------------------------------------------------------------------------|
| SOP Class UID | (0008,0016) | 1 | SOP Class UID for the Class that the dataset contains, just RDSR UID. |
| SOP Instance UID | (0008,0018) | 1 | Created by equipment, using prefix: 1.2.840.113619.2.432.+MAC+ srIndicator (50)+TimeStamp_Exam MAC: Number string of device MAC address. |
| Specific Character Set | (0008,0005) | 1C | Always set to ISO_IR 100(= Latin Alphabet No. 1) |
| Instance Creation Date | (0008,0012) | 3 | Used |
| Instance Creation Time | (0008,0013) | 3 | Used |
| Instance Creator UID | (0008,0014) | 3 | Used |
| Related General SOP Class UID | (0008,001A) | 3 | Not Used |
| Original Specialized SOP Class UID | (0008,001B) | 3 | Not Used |
| Coding Scheme Identification Sequence | (0008,0110) | 3 | Not Used |
| >Coding Scheme Designator | (0008,0102) | 1 | Not Used |
| >Coding Scheme Registry | (0008,0112) | 1C | Not Used |
| >Coding Scheme UID | (0008,010C) | 1C | Not Used |
| >Coding Scheme External ID | (0008,0114) | 2C | Not Used |
| >Coding Scheme Name | (0008,0115) | 3 | Not Used |
| >Coding Scheme Version | (0008,0103) | 3 | Not Used |
| >Coding Scheme Responsible | (0008,0116) | 3 | Not Used |

OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X

GE Healthcare

DICOM CONFORMANCE STATEMENT

2148033 REV 1

| | | | |
|----------------------------------------------|-------------|----|------------------|
| Organization | | | |
| Timezone Offset From UTC | (0008,0201) | 3 | Not Used |
| Contributing Equipment Sequence | (0018,A001) | 3 | Not Used |
| >Purpose of Reference Code Sequence | (0040,A170) | 1 | Not Used |
| <i>>>Include 'Code Sequence Macro'</i> | | | |
| >Manufacturer | (0008,0070) | 1 | Not Used |
| >Institution Name | (0008,0080) | 3 | Not Used |
| >Institution Address | (0008,0081) | 3 | Not Used |
| >Station Name | (0008,1010) | 3 | Not Used |
| >Institutional Department Name | (0008,1040) | 3 | Not Used |
| >Manufacturer's Model Name | (0008,1090) | 3 | Not Used |
| >Device Serial Number | (0018,1000) | 3 | Not Used |
| >Software Versions | (0018,1020) | 3 | Not Used |
| >Spatial Resolution | (0018,1050) | 3 | Not Used |
| >Date of Last Calibration | (0018,1200) | 3 | Not Used |
| >Time of Last Calibration | (0018,1201) | 3 | Not Used |
| >Contribution DateTime | (0018,A002) | 3 | Not Used |
| >Contribution Description | (0018,A003) | 3 | Not Used |
| Instance Number | (0020,0013) | 3 | Instance number. |
| SOP Instance Status | (0100,0410) | 3 | Not Used |
| SOP Authorization Date and Time | (0100,0420) | 3 | Not Used |
| SOP Authorization Comment | (0100,0424) | 3 | Not Used |
| Authorization Equipment Certification Number | (0100,0426) | 3 | Not Used |
| MAC Parameters Sequence | (4FFE,0001) | 3 | Not Used |
| >MAC ID Number | (0400,0005) | 1 | Not Used |
| >MAC Calculation Transfer Syntax UID | (0400,0010) | 1 | Not Used |
| >MAC Algorithm | (0400,0015) | 1 | Not Used |
| >Data Elements Signed | (0400,0020) | 1 | Not Used |
| Digital Signatures Sequence | (FFFA,FFFA) | 3 | Not Used |
| >MAC ID Number | (0400,0005) | 1 | Not Used |
| >Digital Signature UID | (0400,0100) | 1 | Not Used |
| >Digital Signature DateTime | (0400,0105) | 1 | Not Used |
| >Certificate Type | (0400,0110) | 1 | Not Used |
| >Certificate of Signer | (0400,0115) | 1 | Not Used |
| >Signature | (0400,0120) | 1 | Not Used |
| >Certified Timestamp Type | (0400,0305) | 1C | Not Used |
| >Certified Timestamp | (0400,0310) | 3 | Not Used |
| >Digital Signature Purpose Code | (0400,0401) | 3 | Not Used |

| | | | |
|-----------------------------------------------------------------------------------------------------------------------|-------------|----|------------------------------------------------------------------------|
| Sequence | | | |
| >>Include 'Code Sequence Macro' | | | |
| Encrypted Attributes Sequence | (0400,0500) | 1C | Not Used |
| >Encrypted Content Transfer Syntax UID | (0400,0510) | 1 | Not Used |
| >Encrypted Content | (0400,0520) | 1 | Not Used |
| Original Attributes Sequence | (0400,0561) | 3 | Not Used |
| >Source of Previous Values | (0400,0564) | 2 | Not Used |
| >Attribute Modification DateTime | (0400,0562) | 1 | Not Used |
| >Modifying System | (0400,0563) | 1 | Not Used |
| >Reason for the Attribute Modification | (0400,0565) | 1 | Not Used |
| >Modified Attributes Sequence | (0400,0550) | 1 | Not Used |
| >>Any Attribute from the main data set that was modified or removed; may include Sequence Attributes and their Items. | | | |
| HL7 Structured Document Reference Sequence | (0040,A390) | 1C | Not Used |
| >Referenced SOP Class UID | (0008,1150) | 1 | Not Used |
| >Referenced SOP Instance UID | (0008,1155) | 1 | Not Used |
| >HL7 Instance Identifier | (0040,E001) | 1 | Not Used |
| >Retrieve URI | (0040,E010) | 3 | Not Used |
| Longitudinal Temporal Information Modified | (0028,0303) | 3 | Used. If Include Remove Patient Information then fix to "MODIFIED". |

5.5 STANDARD EXTENDED AND PRIVATE DATA ATTRIBUTES

Not applicable

5.6 STANDARD EXTENDED AND PRIVATE CONTEXT GROUPS

The Workstation supports coded terminology using Standard Extended, Private, and Configurable Context Groups defined in the following sections.

5.6.1 Standard Extended Context Groups

The Workstation supports the following extensions to standard Context Groups for SR SOP Instances created by this product. Extensions are indicated by **bold text**.

5.6.2 Configurable Context Groups

The Workstation supports the following Configurable Context Groups for SR SOP Instances created by this product.

| Context Group | Default Value Set | Use |
|----------------------|--------------------------|------------|
| N/A | N/A | N/A |

5.7 STANDARD, STANDARD EXTENDED AND PRIVATE TEMPLATES

The Workstation supports the Standard Extended and Private Templates defined in the following sections.

5.7.1 Standard Templates

The Workstation supports the following standard templates for SOP Instances created by this product.

5.7.1.1 Template ID 10001 X-Ray Radiation Dose

**TID 10001
PROJECTION X-RAY RADIATION DOSE
Type: Extensible Order: Significant**

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|----|-----|-----------------|-----------|-----------------------------------------------------|-----|----------|-------------------------|-----------------------------------------------------------------|
| 1 | | | CONTAINER | EV (113701, DCM, "X-Ray Radiation Dose Report") | 1 | M | | Root Node |
| 2 | > | HAS CONCEPT MOD | CODE | EV (121058, DCM, "Procedure reported") | 1 | M | | DT (113704, DCM, "Projection X-Ray") |
| 3 | >> | HAS CONCEPT MOD | CODE | EV (G-C0E8, SRT, "Has Intent") | 1 | M | | (R-408C3,SRT,Diagnostic Intent) |
| 4 | > | CONTAINS | CODE | EV (122142, DCM, "Acquisition Device Type") | 1 | U | | (DCM, 113957, Fluoroscopy-Guided Projection Radiography System) |
| 5 | > | | INCLUDE | DTID (1002) Observer Context | 1-n | M | | |
| 6 | > | HAS OBS CONTEXT | CODE | EV (113705, DCM, "Scope of Accumulation") | 1 | M | | (113016, DCM, "Performed Procedure Step"). |
| 7 | >>> | HAS PROPERTIES | UIDREF | DCID (10001) UID Types | 1 | M | | (121126, DCM, "Performed Procedure Step SOP Instance UID") |
| 8 | > | CONTAINS | CODE | EV (113945, DCM, "X-Ray Detector Data Available") | 1 | U | DCID 230 "Yes-No" | (SRT, R-0038D, Yes) |
| 9 | > | CONTAINS | CODE | EV (113943, DCM, "X-Ray Source Data Available") | 1 | U | DCID 230 "Yes-No" | (SRT, R-0038D, Yes) |
| 10 | > | CONTAINS | CODE | EV (113944, DCM, "X-Ray Mechanical Data Available") | 1 | U | DCID 230 "Yes-No" | (SRT, R-00339, No) |
| 11 | > | CONTAINS | INCLUDE | DTID (10002) Accumulated X-Ray Dose | 1 | MC | IFF Single Plane system | \$Plane = EV (113622, DCM, "Single Plane" |
| 12 | > | CONTAINS | INCLUDE | DTID (10002) | 1 | MC | IFF Biplane system | Not Used |

**OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | | | | | | |
|----|---|----------|---------|------------------------------------------------|-----|----|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------|
| | | | | Accumulated X-Ray Dose | | | | |
| 13 | > | CONTAINS | INCLUDE | DTID (10002) Accumulated X-Ray Dose | 1 | MC | IFF Biplane system | Not Used |
| 14 | > | CONTAINS | INCLUDE | DTID (10003) Irradiation Event X-Ray Data | 1-n | M | IF any of the values of TID(10001) Row 18 are not (113858, DCM, "MPPS Content"), (113866, DCM, "Copied From Image Attributes") or (113867, DCM, "Computed From Image Attributes") | Used |
| 15 | > | CONTAINS | TEXT | EV (121106, DCM, "Comment") | 1 | U | | Patient Comment |
| 16 | > | CONTAINS | IMAGE | EV (121342, DCM, "Dose Image") | 1-n | U | | Not Used |
| 17 | > | CONTAINS | INCLUDE | DTID (1020) Person Participant | 1 | U | | Not Used |
| 18 | > | CONTAINS | CODE | EV (113854, DCM, "Source of Dose Information") | 1-n | M | | (113856, DCM, "Automated Data Collection") |

5.7.1.2 Template ID 10002 Accumulated X-Ray Dose

**TID 10002
ACCUMULATED X-RAY DOSE
Type: Extensible Order: Significant**

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|---|----|-----------------|-----------|-------------------------------------------------|-----|----------|-----------------------------------|--------------------------------------------------------------------------------------------|
| 1 | | | CONTAINER | EV (113702, DCM, "Accumulated X-Ray Dose Data") | 1 | M | | |
| 2 | > | HAS CONCEPT MOD | CODE | EV (113764, DCM, "Acquisition Plane") | 1 | M | | (113622, DCM, "Single Plane") |
| 3 | > | CONTAINS | CONTAINER | EV (122505, DCM, "Calibration") | 1-n | MC | IFF Calibration Data is available | |
| 4 | >> | HAS CONCEPT MOD | CODE | EV (113794, DCM, "Dose Measurement Device") | 1 | M | | (A-2C090, SRT, "Dosimeter") |
| 5 | >> | CONTAINS | DATETIME | EV (113723, DCM, "Calibration Date") | 1 | M | | The date of DAP Calibration performed. |
| 6 | >> | CONTAINS | NUM | EV (122322, DCM, "Calibration Factor") | 1 | M | | Units = EV (1, UCUM, "no units") From QC user interface otherwise default value is 1.0. |

**OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare

2148033 REV 1

| | | | | | | | | |
|----|----|----------|---------|-------------------------------------------------------------------------------|---|----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|
| 7 | >> | CONTAINS | NUM | EV (113763, DCM, "Calibration Uncertainty") | 1 | M | | Units = EV (% , UCUM, "Percent") From QC user interface otherwise default value is 25. |
| 8 | >> | CONTAINS | TEXT | EV (113724, DCM, "Calibration Responsible Party") | 1 | M | | From QC user interface otherwise default value is GE. |
| 9 | >> | CONTAINS | TEXT | EV (113720, DCM, "Calibration Protocol") | 1 | U | | From QC user interface otherwise default value is empty. |
| 10 | > | CONTAINS | INCLUDE | DTID 10004 "Accumulated Fluoroscopy and Acquisition Projection X-Ray Dose" | 1 | MC | IFF TID (10001) Row 4 = (113957, DCM, "Fluoroscopy-Guided Projection Radiography System") or TID (10001) Row 2 = (113704, DCM, "Projection X-Ray") and TID (10001) Row 4 is absent) | Used |
| 11 | > | CONTAINS | INCLUDE | DTID 10005 "Accumulated Mammography X-Ray Dose" | 1 | MC | IFF TID (10001) Row 2 = (P5-40010, SRT, "Mammography") | Not Used |
| 12 | > | CONTAINS | INCLUDE | DTID 10007 "Accumulated Total Projection Radiography Dose" | 1 | MC | IFF TID (10001) Row 4 = (113958, DCM, "Integrated Projection Radiography System") or TID (10001) Row 4 = (113957, DCM, "Fluoroscopy-Guided Projection Radiography System") or TID (10001) Row 2 = (113704, DCM, "Projection X-Ray") and TID (10001) Row 4 is absent) | Used |
| 13 | > | CONTAINS | INCLUDE | DTID 10006 "Accumulated Cassette-based Projection Radiography Dose" | 1 | MC | IFF TID (10001) Row 4 = (113959, DCM, "Cassette-based Projection Radiography System") | Not Used |
| 14 | > | CONTAINS | INCLUDE | DTID (1021) Device Participant | 1 | MC | Required if the irradiating device is | Not Used |

**OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | | | | | | |
|--|--|--|--|--|--|--|---------------------------------------------------------------------------|--|
| | | | | | | | not the recording device and the dose was accumulated on a single device. | |
|--|--|--|--|--|--|--|---------------------------------------------------------------------------|--|

5.7.1.3 Template ID 10003 Irradiation Event X-Ray Data

**TID 10003
IRRADIATION EVENT X-RAY DATA
Type: Extensible Order: Significant**

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|----|----|-----------------|-----------|--------------------------------------------------|----|----------|-------------------------------------------------------------------|--------------------------------------------------------------------------|
| 1 | | | CONTAINER | EV (113706, DCM, "Irradiation Event X-Ray Data") | 1 | M | | |
| 2 | > | HAS CONCEPT MOD | CODE | EV (113764, DCM, "Acquisition Plane") | 1 | M | | (113622, DCM, "Single Plane") |
| 3 | > | CONTAINS | UIDREF | EV (113769, DCM, "Irradiation Event UID") | 1 | M | | The UID value of the irradiation event from the system. |
| 4 | > | CONTAINS | TEXT | EV (113605, DCM, "Irradiation Event Label") | 1 | U | | Value is shot number. |
| 5 | >> | HAS CONCEPT MOD | CODE | EV (113606, DCM, "Label Type") | 1 | MC | IF the value of Row 4 is the value of an Attribute in the images. | (DCM, 113608, Acquisition Number) |
| 6 | > | CONTAINS | DATETIME | DT (111526, DCM, "DateTime Started") | 1 | M | | Date and time of each shot |
| 7 | > | CONTAINS | CODE | EV (113721, DCM, "Irradiation Event Type") | 1 | M | | (P5-06000, SRT, "Fluoroscopy"). |
| 8 | > | CONTAINS | TEXT | EV (125203, DCM, "Acquisition Protocol") | 1 | U | | The text value of the IOD tag (0018,1030) Protocol Name from the system. |
| 9 | > | CONTAINS | CODE | EV (T-D0005, SRT, "Anatomical structure") | 1 | U | | Not Used |
| 10 | >> | HAS CONCEPT MOD | CODE | EV (G-C171, SRT, "Laterality") | 1 | UC | If anatomy is bilateral | Not Used |

**OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|----|----|-----------------|------|--------------------------------------------------|-----|----------|--------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 11 | > | CONTAINS | CODE | EV (111031, DCM, "Image View") | 1 | U | | Not Used |
| 12 | >> | HAS CONCEPT MOD | CODE | EV (111032, DCM, "Image View Modifier") | 1-n | U | | Not Used |
| 13 | >> | CONTAINS | CODE | EV (113946, DCM, "Projection Eponymous Name") | 1 | U | | Not Used |
| 14 | > | CONTAINS | CODE | EV (113745, DCM, "Patient Table Relationship") | 1 | U | | Not Used |
| 15 | > | CONTAINS | CODE | EV (113743, DCM, "Patient Orientation") | 1 | U | | Not Used |
| 16 | >> | HAS CONCEPT MOD | CODE | EV (113744, DCM, "Patient Orientation Modifier") | 1 | M | | Not Used |
| 17 | > | CONTAINS | CODE | EV (123014, DCM, "Target Region") | 1 | M | | 1. The code from the system using DICOM Part 16 Context ID 4031 values and supporting the code (T-D0001, SRT, "Topography unknown") for an 'Unknown' value which is SNOMED term. 2. User could set this value in EM->Additional Screen. Default is (T-D0001, SRT, "Topography unknown") |
| 18 | > | CONTAINS | NUM | EV (122130, DCM, "Dose Area Product") | 1 | MC | IFF TID (10001) Row 2 = (113704, DCM, "Projection X-Ray") | UNITS = EV (Gy.m2, UCUM, "Gy.m2") |
| 19 | > | CONTAINS | NUM | EV (111634, DCM, "Half Value Layer") | 1 | U | | Not Used |
| 20 | > | CONTAINS | NUM | EV (111638, DCM, "Patient Equivalent Thickness") | 1 | U | | Not Used |

OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X

GE Healthcare

DICOM CONFORMANCE STATEMENT

2148033 REV 1

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|----|----|-----------------|---------|-------------------------------------------------------|-----|----------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------|
| 21 | > | CONTAINS | NUM | EV (111636, DCM, "Entrance Exposure at RP") | 1 | MC | IF TID (10001) Row 2 = (P5-40010, SRT, "Mammography") and (TID (10001) Row 9 is absent or value is (R-0038D, SRT, "Yes")) and (TID (10001) Row 10 is absent or value is (R-0038D, SRT, "Yes")) | Not Used |
| 22 | > | CONTAINS | TEXT | EV (113780, DCM, "Reference Point Definition") | 1 | MC | IF Row 21 is present and Row 23 is not present | Not Used |
| 23 | > | CONTAINS | CODE | EV (113780, DCM, "Reference Point Definition") | 1 | MC | IF Row 21 is present and Row 22 is not present | The value of code (113861, DCM, "30cm in Front of Image Input Surface). |
| 24 | > | CONTAINS | INCLUDE | DTID 4007 "Mammography CAD Breast Composition" | 1 | U | | Not Used |
| 25 | > | CONTAINS | TEXT | EV (121106, DCM, "Comment") | 1 | U | | Patient Comment |
| 26 | > | CONTAINS | INCLUDE | DTID 1020 "Person Participant" | 1-n | U | | Not Used |
| 27 | > | CONTAINS | INCLUDE | DTID 10003A "Irradiation Event X-Ray Detector Data" | 1 | MC | IFF TID (10001) Row 8 is absent or has a value of (R-0038D, SRT, "Yes") | Used |
| 28 | > | CONTAINS | INCLUDE | DTID 10003B "Irradiation Event X-Ray Source Data" | 1 | MC | IFF TID (10001) Row 9 is absent or has a value of (R-0038D, SRT, "Yes") | Used |
| 29 | > | CONTAINS | INCLUDE | DTID 10003C "Irradiation Event X-Ray Mechanical Data" | 1 | MC | IFF TID (10001) Row 10 is absent or has a value of (R-0038D, SRT, "Yes") | Not Used |

**OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare

2148033 REV 1

5.7.1.4 TID 10003A Irradiation Event X-Ray Detector Data

**TID 10003A
Irradiation Event X-Ray Detector Data
Type: Extensible Order: Non-Significant**

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|---|----|-----------------|---------|-------------------------------------------|-----|----------|-----------------------------------------------------------|----------------------|
| 1 | | | NUM | EV (113845, DCM, "Exposure Index") | 1 | MC | IF the value is displayable to the X-Ray system operator. | Not Used |
| 2 | | | NUM | EV (113846, DCM, "Target Exposure Index") | 1 | MC | IF the value is displayable to the X-Ray system operator. | Not Used |
| 3 | | | NUM | EV (113847, DCM, "Deviation Index") | 1 | MC | IF the value is displayable to the X-Ray system operator. | Not Used |
| 4 | | | INCLUDE | DTID 1021 "Device Participant" | 1 | U | | Not Used |
| 5 | | | IMAGE | EV (113795, DCM, "Acquired Image") | 1-n | MC | IFF Image Object is created for this irradiation event | All image UIDs. |

5.7.1.5 TID 10003B Irradiation Event X-Ray Source Data

**TID 10003B
Irradiation Event X-Ray Source Data
Type: Extensible Order: Non-Significant**

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|---|----|-----------------|------|------------------------------------------------|----|----------|--------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------|
| 1 | | | NUM | EV (113738, DCM, "Dose (RP)") | 1 | MC | IF TID (10001) Row 2 = (113704, DCM, "Projection X-Ray") AND any of the values of TID (10001) Row 18 are not (113858, DCM, "MPPS Content") | UNITS = EV (Gy, UCUM, "Gy") Dose applied at the Reference Point (RP) |
| 2 | | | TEXT | EV (113780, DCM, "Reference Point Definition") | 1 | MC | IF Row 1 is present and Row 3 is not present | Not Used |
| 3 | | | CODE | EV (113780, DCM, "Reference | 1 | MC | IF Row 1 is present and Row 2 is not | The value of code (113861, DCM, "30cm in Front of |

**OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare

2148033 REV 1

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|----|----|-----------------|------|--------------------------------------------|-----|----------|---------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| | | | | Point Definition") | | | present | Image Input Surface). |
| 4 | | | NUM | EV (111631, DCM, "Average Glandular Dose") | 1 | MC | IFF TID (10001) Row 2 = (P5-40010, SRT, "Mammography") | Not Used |
| 5 | | | CODE | EV (113732, DCM, "Fluoro Mode") | 1 | UC | IFF TID (10003) Row 7 value = (P5-06000, SRT, "Fluoroscopy") | 1. (113630, DCM, "Continuous") for continuous mode irradiation events. 2. (113631, DCM, "Pulsed") for pulse mode irradiation events. |
| 6 | | | NUM | EV (113791, DCM, "Pulse Rate") | 1 | MC | IFF Row 5 value = (113631, DCM, "Pulsed") | Units = EV ({pulse}/s, UCUM, "pulse/s") For Pulsed mode irradiation events. |
| 7 | | | NUM | EV (113768, DCM, "Number of Pulses") | 1 | MC | IFF Row 5 is not present or Row 5 is present and equals (113631, DCM, "Pulsed") | Units = EV (1, UCUM, "no units") 1. For Pulsed mode irradiation events. 2. Estimated: a. Manual: exposure time/pulse width b. Auto: exposure time/pulse width (margin of error might be bigger than manual) c. Pulse Width = 49.0ms |
| 8 | > | HAS CONCEPT MOD | CODE | EV (121401, DCM, "Derivation") | 1 | MC | IFF count of pulses in Row 7 is estimated | EV (R-10260, SRT, "Estimated") |
| 9 | | | NUM | EV (113793, DCM, "Pulse Width") | 1-n | U | | Units = EV (ms, UCUM, "ms") Value is 49.0ms |
| 10 | | | NUM | EV (113742, DCM, "Irradiation Duration") | 1 | U | | Not Used |
| 11 | | | NUM | EV (113733, DCM, "KVP") | 1-n | M | | UNITS = EV (kV, UCUM, "kV") |

**OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|----|----|-----------------|-----------|----------------------------------------------------|-----|----------|--------------------------------|-----------------------------------------------------------------------------|
| | | | | | | | | For Continuous/Pulsed mode irradiation events. |
| 12 | | | NUM | EV (113734, DCM, "X-Ray Tube Current") | 1-n | MC | IF Row 15 is not present | UNITS = EV (mA, UCUM, "mA") |
| 13 | | | NUM | EV (113767, DCM, "Average X-Ray Tube Current") | 1 | U | | UNITS = EV (mA, UCUM, "mA") The same with above "X-Ray Tube Current" |
| 14 | | | NUM | EV (113824, DCM, "Exposure Time") | 1 | MC | IF Row 15 is not present | UNITS = EV (ms, UCUM, "ms") |
| 15 | | | NUM | EV (113736, DCM, "Exposure") | 1-n | MC | IF Row 12 or 14 is not present | Not Used |
| 16 | | | NUM | EV (113766, DCM, "Focal Spot Size") | 1 | U | | Not Used |
| 17 | | | CODE | EV (111632, DCM, "Anode Target Material") | 1 | U | | Not Used |
| 18 | | | CONTAINER | EV (113771, DCM, "X-Ray Filters") | 1-n | U | | Not Used |
| 19 | > | CONTAINS | CODE | EV (113772, DCM, "X-Ray Filter Type") | 1 | U | | Not Used |
| 20 | > | CONTAINS | CODE | EV (113757, DCM, "X-Ray Filter Material") | 1 | U | | Not Used |
| 21 | > | CONTAINS | NUM | EV (113758, DCM, "X-Ray Filter Thickness Minimum") | 1 | U | | Not Used |
| 22 | > | CONTAINS | NUM | EV (113773, DCM, "X-Ray Filter Thickness Maximum") | 1 | U | | Not Used |

**OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare

2148033 REV 1

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|----|----|-----------------|---------|---------------------------------------------|-----|----------|-----------------------------------------------------------------|----------------------|
| 23 | | | NUM | EV (113790, DCM, "Collimated Field Area") | 1 | U | | Not Used |
| 24 | | | NUM | EV (113788, DCM, "Collimated Field Height") | 1 | U | | Not Used |
| 25 | | | NUM | EV (113789, DCM, "Collimated Field Width") | 1 | U | | Not Used |
| 26 | | | CODE | EV (111635, DCM, "X-Ray Grid") | 1-n | U | | Not Used |
| 27 | | | INCLUDE | DTID 1021 "Device Participant" | 1 | MC | Required if the irradiating device is not the recording device. | Not Used |

5.7.1.6 Template ID 10004 Accumulated Projection X-Ray Dose

**TID 10004
ACCUMULATED PROJECTION X-RAY DOSE
Type: Extensible Order: Significant**

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|---|----|-----------------|-----|----------------------------------------------------|----|----------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|
| 1 | | | NUM | EV (113726, DCM, "Fluoro Dose Area Product Total") | 1 | MC | IFF TID (10003) Row 7 value = (P5-06000, SRT, "Fluoroscopy") for at least one irradiation event | UNITS = EV (Gy.m2, UCUM, "Gy.m2") The numeric value of Total Radioscopic Dose Area Product. |
| 2 | | | NUM | EV (113728, DCM, "Fluoro Dose (RP) Total") | 1 | MC | IFF TID (10003) Row 7 value = (P5-06000, SRT, "Fluoroscopy") for at least one irradiation event AND any of the values of TID (10001) Row 18 are not (113858, DCM, "MPPS Content"). | UNITS = EV (Gy, UCUM, "Gy") The numeric value of Cumulative Radioscopic Air Kerma. |

OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X

GE Healthcare

DICOM CONFORMANCE STATEMENT

2148033 REV 1

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|---|----|-----------------|-----|---------------------------------------------------------|----|----------|-------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|
| 3 | | | NUM | EV (113730, DCM, "Total Fluoro Time") | 1 | MC | IFF TID (10003) Row 7 value = (P5-06000, SRT, "Fluoroscopy") for at least one irradiation event. | UNITS = EV (s, UCUM, "s") The numeric value of Total Radioscopic time to the second. |
| 4 | | | NUM | EV (113727, DCM, "Acquisition Dose Area Product Total") | 1 | M | | UNITS = EV (Gy.m2, UCUM, "Gy.m2") The numeric value of Total Radiographic Dose Area Product. |
| 5 | | | NUM | EV (113729, DCM, "Acquisition Dose (RP) Total") | 1 | MC | IF any of the values of TID (10001) Row 18 are not (113858, DCM, "MPPS Content"). May be present otherwise. | UNITS = EV (Gy, UCUM, "Gy") The numeric value of Cumulative Radiographic Air kerma. |
| 6 | | | NUM | EV (113855, DCM, "Total Acquisition Time") | 1 | M | | UNITS = EV (s, UCUM, "s") The numeric value of Total Radiographic pedal time in milliseconds. |

5.7.1.7 TID 10007 Accumulated Total Projection Radiography Dose

TID 10007

Accumulated Total Projection Radiography Dose

Type: Non-Extensible Order: Significant

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|---|----|-----------------|-----|---------------------------------------------|----|----------|--------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|
| 1 | | | NUM | EV (113722, DCM, "Dose Area Product Total") | 1 | M | | UNITS = EV (Gy.m2, UCUM, "Gy.m2") The numeric value of Total Dose Area Product. (Radioscopic) |
| 2 | | | NUM | EV (113725, DCM, "Dose (RP) Total") | 1 | MC | IF TID (10001) Row 4 = (113958, DCM, "Integrated Projection Radiography System") or any of the values of TID | UNITS = EV (Gy, UCUM, "Gy") The numeric value of Cumulative Air Kerma.(Radioscopic) |

OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X

GE Healthcare

DICOM CONFORMANCE STATEMENT

2148033 REV 1

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|---|----|-----------------|------|---------------------------------------------------------|----|----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|
| | | | | | | | (10001) Row 18 are not (113858, DCM, "MPPS Content"). | |
| 3 | | | NUM | EV (113737, DCM, "Distance Source to Reference Point") | 1 | U | | UNITS = EV (mm, UCUM, "mm") Value is 700 |
| 4 | | | NUM | EV (113731, DCM, "Total Number of Radiographic Frames") | 1 | U | | UNITS = EV (1, UCUM, "no units") The numeric value of Total Film and Digital Spot Frames. |
| 5 | | | CODE | EV (113780, DCM, "Reference Point Definition") | 1 | MC | IF any of (113725, DCM, "Dose (RP) Total"), (113728, DCM, "Fluoro Dose (RP) Total") or (113729, DCM, "Acquisition Dose (RP) Total") are present, and Row 6 is not present. | (113861, DCM, "30cm in Front of Image Input Surface"). |
| 6 | | | TEXT | EV (113780, DCM, "Reference Point Definition") | 1 | MC | IF any of (113725, DCM, "Dose (RP) Total"), (113728, DCM, "Fluoro Dose (RP) Total") or (113729, DCM, "Acquisition Dose (RP) Total") are present, and Row 5 is not present. | Not Used |

5.7.1.8 Template ID 1002 Person Observer Context

TID 1002

OBSERVER CONTEXT

Type: Non-Extensible Order: Significant

| | NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|---|----|-----------------|---------|----------------------------------------------------|----|----------|-------------------------------------------------------------|-------------------------|
| 1 | | HAS OBS CONTEXT | CODE | EV (121005,DCM, "Observer Type") | 1 | MC | IF Observer type is device | (121006, DCM, "Person") |
| 2 | | HAS OBS CONTEXT | INCLUDE | DTID (1003) Person observer identifying attributes | 1 | MC | IFF Row 1 value = (121006,DCM, "Person") or Row 1 is absent | Used |

OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X

GE Healthcare

DICOM CONFORMANCE STATEMENT

2148033 REV 1

| | | | | | | | |
|---|------------------------|---------|----------------------------------------------------------|---|----|---------------------------------------------|----------|
| 3 | HAS OBS CONTEX T | INCLUDE | DTID (1004) Device observer identifying attributes | 1 | MC | IFF Row 1 value = (121007,DCM, "Device") | Not Used |
|---|------------------------|---------|----------------------------------------------------------|---|----|---------------------------------------------|----------|

5.7.1.9 Template ID 1002 Device Observer Context

TID 1002

OBSERVER CONTEXT

Type: Non-Extensible Order: Significant

| NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|----|------------------------|---------|----------------------------------------------------------|----|-------------|-------------------------------------------------------------------|-------------------------|
| 1 | HAS OBS CONTEX T | CODE | EV (121005,DCM, "Observer Type") | 1 | MC | IF Observer type is device | (121007, DCM, "Device") |
| 2 | HAS OBS CONTEX T | INCLUDE | DTID (1003) Person observer identifying attributes | 1 | MC | IFF Row 1 value = (121006,DCM, "Person") or Row 1 is absent | Not Used |
| 3 | HAS OBS CONTEX T | INCLUDE | DTID (1004) Device observer identifying attributes | 1 | MC | IFF Row 1 value = (121007,DCM, "Device") | Used |

5.7.1.10 TID 1003 Person Observer Identifying Attributes

TID 1003

DEVICE OBSERVER IDENTIFYING ATTRIBUTES

Type: Extensible Order: Significant

| NL | Rel with Parent | VT | Concept Name | V M | Req Type | Condition | Value Set Constraint |
|----|--------------------|-----------|----------------------------------------------------------------------|--------|----------|-----------|---------------------------------------------------------------------------------|
| 1 | | PNAM E | EV (121008, DCM, "Person Observer Name") | 1 | M | | This value is set to the performing physician. |
| 2 | | TEXT | EV (121009, DCM, "Person Observer's Organization Name") | 1 | U | | Defaults to Institution Name (0008,0080) of the General Equipment Module. |
| 3 | | CODE | EV (121010, DCM, "Person Observer's Role in the Organization") | 1 | U | | ("J-0016E", "SRT", "Medical Practitioner") |
| 4 | | CODE | EV (121011, DCM, "Person Observer's Role in this Procedure") | 1 | U | | ("121094", "DCM", "Performing") |

5.7.1.11 Template ID 1004 Device Observer Identifying Attributes

TID 1004

DEVICE OBSERVER IDENTIFYING ATTRIBUTES

Type: Extensible Order: Significant

| NL | Rel with Parent | VT | Concept Name | VM | Req Type | Condition | Value Set Constraint |
|----|--------------------|--------|-------------------------------------------|----|-------------|-----------|----------------------|
| 1 | | UIDREF | EV (121012,DCM, "Device Observer UID") | 1 | M | | System Generated |

**OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare

2148033 REV 1

| | | | | | | | |
|---|--|------|----------------------------------------------------------------------------|-----|---|--|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| 2 | | TEXT | EV (121013,DCM, "Device Observer Name") | 1 | U | | Defaults to value of Station Name (0008,1010) in General Equipment Module. If 'Station Name' is blank the value for Local AE Title shall be used. |
| 3 | | TEXT | EV (121014,DCM, "Device Observer Manufacturer") | 1 | U | | Defaults to value of Manufacturer (0008,0070) in General Equipment Module |
| 4 | | TEXT | EV (121015,DCM, "Device Observer Model Name") | 1 | U | | Defaults to value of Manufacturer's Model Name (0008,1090) in General Equipment Module |
| 5 | | TEXT | EV (121016,DCM, "Device Observer Serial Number") | 1 | U | | Defaults to value of Device Serial Number (0018,1000) in General Equipment Module |
| 6 | | TEXT | EV (121017,DCM, "Device Observer Physical Location during observation") | 1 | U | | Hospital Name. |
| 7 | | CODE | EV (113876, DCM, "Device Role in Procedure") | 1-n | U | | (113859, DCM, "Irradiating Device") |

5.7.2 Standard Extended Templates

The Workstation supports the following extensions or restrictions to standard templates for SOP Instances created by this product. Extensions or restrictions are indicated by **bold text**.

6. MODALITY WORKLIST QUERY IMPLEMENTATION

6.1 INTRODUCTION

This section specifies the use of the DICOM Modality Worklist Information Model used to organize data and against which a Modality Worklist Query will be performed.

6.2 WORKSTATION MAPPING OF DICOM ENTITIES

The Workstation maps DICOM Information Entities to local Information Entities in the product's database and user interface.

TABLE 6-1
MAPPING OF DICOM ENTITIES TO WORKSTATION ENTITIES

| DICOM | Workstation Entity |
|--------------------------|--------------------|
| Scheduled Procedure Step | Exam |
| Requested Procedure | Exam |
| Imaging Service Request | Exam |
| Visit | Not Applicable |
| Patient | Patient |

6.3 WORKLIST QUERY MODULE TABLE

See DICOM PS 3.3 and PS 3.4 for a complete definition of the entities, modules, and attributes.

**TABLE 6-2
MODALITY WORKLIST INFORMATION MODEL MODULES**

| Entity Name | Module Name | Reference |
|--------------------------|--------------------------|------------------|
| Scheduled Procedure Step | SOP Common | 6.4.1.1 |
| | Scheduled Procedure Step | 6.4.1.2 |
| Requested Procedure | Requested Procedure | 6.4.2.1 |
| Imaging Service Request | Imaging Service Request | 6.4.3.1 |
| Visit | Visit Identification | 6.4.4.1 |
| | Visit Status | 6.4.4.2 |
| | Visit Relationship | 6.4.4.3 |
| | Visit Admission | 6.4.4.4 |
| Patient | Patient Relationship | 6.4.5.1 |
| | Patient Identification | 6.4.5.2 |
| | Patient Demographic | 6.4.5.3 |
| | Patient Medical | 6.4.5.4 |

6.4 WORKLIST QUERY MODULE DEFINITIONS

Please refer to DICOM Standard PS 3.3. (Information Object Definitions) for a description of each of the query key attributes contained within the Modality Worklist Information Model.

6.4.1 Common Scheduled Procedure Step Entity Modules

6.4.1.1 SOP Common Module

**TABLE 6-3
SOP COMMON MODULE ATTRIBUTES**

| Attribute Name | Tag | Expected Matching Key Type | Expected Returned Key Type | Mapped into Instance / MPPS | Note |
|------------------------|-------------|-----------------------------------|-----------------------------------|------------------------------------|-----------------------------|
| Specific Character Set | (0008,0005) | O | 1C | Yes/Yes | See Section 6.4.1.1.1 below |

6.4.1.1.1 Specific Character Set

The attribute Specific Character Set (0008,0005) will always be sent.

As a Query SCU, it will similarly accept response items with ISO_IR 100 values of Specific Character Set or empty value of Specific Character Set (absence of (0008,0005) which mean usage of ASCII Chars only).

The Workstation user interface will allow the user to enter characters from the console keyboard that is within ISO 8859-1 characters.

Please also refer to Section 2.7.

**OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare

2148033 REV 1

6.4.1.2 Scheduled Procedure Step Module

**TABLE 6-4
SCHEDULED PROCEDURE STEP MODULE ATTRIBUTES**

| Attribute Name | Tag | Expected Matching Key Type | Expected Returned Key Type | Mapped into Instance / MPPS | Note |
|--------------------------------------|-------------|-----------------------------------|-----------------------------------|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Scheduled Procedure Step Sequence | (0040,0100) | R | 1 | No/No | |
| >Scheduled Station AE Title | (0040,0001) | R | 1 | No/No | No Matching is performed. |
| >Scheduled Procedure Step Start Date | (0040,0002) | R | 1 | No/No | Using the date configured from the DICOM query filter screen, the following matching can be done: YYYY- year, MM- month, DD- day. Range matching can also be done YYYYMMDD-YYYYMMDD - all matches included between those dates. The “dash” is important -YYYYMMDD - all matches prior to and including this date. The beginning “dash” is important. YYYYMMDD- all matches after and including this date. The ending “dash” is important. Single matching where the dates entered is equal. The data entered on the user interface is passed to the Server as entered; if it does not work see the Server DICOM conformance statement. |
| >Scheduled Procedure Step Start Time | (0040,0003) | R | 1 | No/No | No Matching is performed. |
| >Scheduled Procedure Step End Date | (0040,0004) | O | 3 | No/No | Not Requested |
| >Scheduled Procedure Step End Time | (0040,0005) | O | 3 | No/No | Not Requested |
| >Modality | (0008,0060) | R | 1 | No/No | Matching is supported and is configured by the user interface. This value could be set to RF/XA/CR/OT/ALL. “ALL” means any modality. Note 1: Send the field as empty to request universal matching when "ALL" is selected from the UI. Note 2:User could only select RF,XA,CR or OT as Modality Type during image archiving. |

**OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | | | |
|----------------------------------------|-------------|---|------|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| >Scheduled Performing Physician's Name | (0040,0006) | R | 2 * | Yes/Yes | Matching is performed. The data entered on the user interface is passed to the server as it is. Single matching is supported for this data element. This value is mapped into (0008, 1050) in the instance. This value is mapped into (0008, 1050) in the MPPS. Supports 5 different components delimited by “^”. Supports a maximum length of 64 bytes including the delimiter. |
| >Scheduled Procedure Step Description | (0040,0007) | O | 1C * | Yes/Yes | No Matching is performed. This value is mapped into (0040,0275)->(0040,0007) & Protocol Name(0018,1030) & Series Description(0008,103E) in the image instance. This value is mapped into Series Description(0008,103E) & EV(125203, DCM, “Acquisition Protocol”) in the RDSR instance. This value is mapped into (0040,0270)->(0040,0007) & Protocol Name(0018,1030) & Series Description(0008,103E) in the MPPS. |
| >Scheduled Station Name | (0040,0010) | O | 2 | No/No | No Matching is performed. |
| >Scheduled Procedure Step Location | (0040,0011) | O | 2 | No/No | No Matching is performed. |
| >Scheduled Protocol Code Sequence | (0040,0008) | O | 1C | Yes/Yes | No Matching is performed. This value is mapped into (0040,0275)-> (0040,0008) in the instance. This value is mapped into (0040,0270)-> (0040,0008) in the MPPS. |
| >>Code Value | (0008,0100) | O | 1 | Yes/Yes | No Matching is performed. |
| >>Coding Scheme Designator | (0008,0102) | O | 1 | Yes/Yes | No Matching is performed. |
| >>Coding Scheme Version | (0008,0103) | O | 3 | Yes/Yes | No Matching is performed. |
| >>Code Meaning | (0008,0104) | O | 3 * | Yes/Yes | No Matching is performed. |
| >Pre-Medication | (0040,0012) | O | 2C | No/No | No Matching is performed. |
| >Scheduled Procedure Step ID | (0040,0009) | O | 1 * | Yes/Yes | No Matching is performed. This value is mapped into (0040,0275)-> (0040,0009) in the instance. This value is mapped into (0040,0270)-> (0040,0009) in the MPPS. |
| >Requested Contrast Agent | (0032,1070) | O | 2C | Yes/No | No Matching is performed. This value is mapped into (0018, 0010) in the instance. |

OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT

GE Healthcare
 2148033 REV 1

| | | | | | |
|-------------------------------------------|-------------|---|---|-------|---------------|
| >Scheduled Procedure Step Status | (0040,0020) | O | 3 | No/No | Not Requested |
| >Comments on the Scheduled Procedure Step | (0040,0400) | O | 3 | No/No | Not Requested |

Note: * in the *Expected Return Key Type* column indicates that this information is displayed on screen, if available

6.4.1.2.1 Scheduled Station AE Title

Always empty in the request, and not used or set for matching.

6.4.2 Common Requested Procedure Entity Modules

6.4.2.1 Requested Procedure Module

TABLE 6-5
REQUESTED PROCEDURE MODULE ATTRIBUTES

| Attribute Name | Tag | Expected Matching Key Type | Expected Returned Key Type | Mapped into Instance / MPPS | Note |
|---------------------------------|-------------|----------------------------|----------------------------|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Requested Procedure ID | (0040,1001) | O | 1 * | Yes/Yes | Matching is performed. The data entered on the user interface is passed to the server as it is. Single matching is supported for this data element. This value is mapped into (0020,0010) in the instance and is displayed as the Study ID of the exam if a value is returned. This value is mapped into(0040,0275)->(0040,1001) in the image instance. This value is mapped into (0040,a370)->(0040,1001) in the RDSR. This value is mapped into (0040,0270)->(0040,1001) in the MPPS. |
| Requested Procedure Description | (0032,1060) | O | 1C * | Yes/Yes | This value is mapped into (0040,0275)->(0032,1060) & Study Description(0008,1030) in the image instance. This value is mapped into (0040,A370)->(0032,1060) in the RDSR. This value is mapped into (0040,0270)->(0032,1060) & Study Description(0008,1030) in the MPPS. |

OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT

GE Healthcare
 2148033 REV 1

| | | | | | |
|-----------------------------------------|-------------|---|----|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Requested Procedure Code Sequence | (0032,1064) | O | 1C | Yes/Yes | No Matching is performed. This value is mapped into (0032, 1064) in the image instance. This value is mapped into (0032, 1064) in the RDSR. This value is mapped into (0008, 1032) in the MPPS. |
| >Code Value | (0008,0100) | O | 1 | Yes/Yes | No Matching is performed. |
| >Coding Scheme Designator | (0008,0102) | O | 1 | Yes/Yes | No Matching is performed. |
| >Coding Scheme Version | (0008,0103) | O | 3 | Yes/Yes | No Matching is performed. |
| >Code Meaning | (0008,0104) | O | 3 | Yes/Yes | No Matching is performed. |
| Study Instance UID | (0020,000D) | O | 1 | Yes/Yes | The product will adopt the Study Instance UID which is returned in the MWL response. If there is no Study Instance UID returned from MWL server, then new Study Instance UID will be system created. |
| Study Date | (0008,0020) | O | 3 | No/No | Not Requested |
| Study Time | (0008,0030) | O | 3 | No/No | Not Requested |
| Referenced Study Sequence | (0008,1110) | O | 2 | Yes/Yes | This value is mapped into (0008,1110) in the image instance. This value is mapped into (0008,1110) in the RDSR. This value is mapped into (0008, 1110) in the MPPS. |
| >Referenced SOP Class UID | (0008,1150) | O | 1C | Yes/Yes | |
| >Referenced SOP Instance UID | (0008,1155) | O | 1C | Yes/Yes | |
| Requested Procedure Priority | (0040,1003) | O | 2 | No/No | No Matching is performed. |
| Patient Transport Arrangements | (0040,1004) | O | 2 | No/No | Not Requested |
| Requested Procedure Location | (0040,1005) | O | 3 | No/No | Not Requested |
| Confidentiality Code | (0040,1008) | O | 3 | No/No | Not Requested |
| Reporting Priority | (0040,1009) | O | 3 | No/No | Not Requested |
| Names of Intended Recipients of Results | (0040,1010) | O | 3 | No/No | No Matching is performed. |
| Reason for the Requested Procedure | (0040,1002) | O | 3 | No/No | No Matching is performed. |
| Requested Procedure Comments | (0040,1400) | O | 3 | No/No | Not Requested |

Note: * in the *Expected Return Key Type* column indicates that this information is displayed on screen, if available

6.4.2.1.1 Study Instance UID

The Workstation will adopt the Study Instance UID which is returned in the MWL response. If there is no Study Instance UID returned from MWL server, then new Study Instance UID will be created.

6.4.3 Common Imaging Service Request Entity Modules

6.4.3.1 Imaging Service Request Module

**TABLE 6-6
IMAGING SERVICE REQUEST MODULE ATTRIBUTES**

| Attribute Name | Tag | Expected Matching Key Type | Expected Returned Key Type | Mapped into Instance / MPPS | Note |
|-----------------------------------------------|-------------|-----------------------------------|-----------------------------------|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Accession Number | (0008,0050) | O | 2 * | Yes/Yes | Matching is performed. The data entered on the user interface is passed to the server as it is. Single matching is supported for this data element. Truncated to 16 bytes, end truncated This value is mapped into (0008,0050) in the instance or RDSR. This value is mapped into (0008, 0050) in the MPPS. |
| Requesting Physician | (0032,1032) | O | 2 | No/No | No Matching is performed. |
| Referring Physician's Name | (0008,0090) | O | 2 * | Yes/No | No Matching is performed. This value is mapped into (0008,0090) in the instance or RDSR. |
| Requesting Service | (0032,1033) | O | 3 | No/No | Not Requested |
| Imaging Service Request Comments | (0040,2400) | O | 3 | No/No | Not Requested |
| Issue Date of Imaging Service Request | (0040,2004) | O | 3 | No/No | Not Requested |
| Issue Time of Imaging Service Request | (0040,2005) | O | 3 | No/No | Not Requested |
| Placer Order Number / Imaging Service Request | (0040,2016) | O | 3 | No/No | Not Requested |
| Filler Order Number / Imaging Service Request | (0040,2017) | O | 3 | No/No | Not Requested |
| Order entered by ... | (0040,2008) | O | 3 | No/No | Not Requested |
| Order Enterer's Location | (0040,2009) | O | 3 | No/No | Not Requested |
| Order Callback Phone Number | (0040,2010) | O | 3 | No/No | Not Requested |

Note: * in the *Expected Return Key Type* column indicates that this information is displayed on screen, if available

6.4.4 Common visit Entity Modules

6.4.4.1 Visit Identification

Not applicable

6.4.4.2 Visit Status

Not applicable

6.4.4.3 Visit Relationship

Not applicable

6.4.4.4 Visit Admission

Not applicable

6.4.5 Common Patient Entity Modules

6.4.5.1 Patient Relationship

Not applicable

6.4.5.2 Patient Identification

TABLE 6-7
PATIENT IDENTIFICATION MODULE ATTRIBUTES

| Attribute Name | Tag | Expected Matching Key Type | Expected Returned Key Type | Mapped into Instance / MPPS | Note |
|----------------------|-------------|----------------------------|----------------------------|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Patient's Name | (0010,0010) | R | 1 * | Yes/Yes | Matching is performed. Single or Wildcard matching is supported for this data element. Supports 5 different components delimited by “^”. Supports a maximum length of 64 bytes including the delimiter. |
| Patient ID | (0010,0020) | R | 1 * | Yes/Yes | Matching is performed. Single matching is supported for this data element. Truncated to 64 bytes, end truncated |
| Issuer of Patient ID | (0010,0021) | O | 3 | Yes/Yes | No Matching is performed. |
| Other Patient IDs | (0010,1000) | O | 3 | No/No | Not Requested |
| Other Patient Names | (0010,1001) | O | 3 | No/No | Not Requested |
| Patient's Birth Name | (0010,1005) | O | 3 | No/No | Not Requested |

**OEK ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | | | |
|-------------------------------|-------------|---|---|-------|---------------|
| Patient's Mother's Birth Name | (0010,1060) | O | 3 | No/No | Not Requested |
|-------------------------------|-------------|---|---|-------|---------------|

Note: * in the *Expected Return Key Type* column indicates that this information is displayed on screen, if available

6.4.5.3 Patient Demographic

**TABLE 6-8
PATIENT DEMOGRAPHIC MODULE ATTRIBUTES**

| Attribute Name | Tag | Expected Matching Key Type | Expected Returned Key Type | Mapped into Instance / MPPS | Note |
|-----------------------------------------------------|-------------|----------------------------|----------------------------|-----------------------------|-----------------------------------------------------------------------------------------------------------|
| Patients Birth Date | (0010,0030) | O | 2 * | Yes/Yes | No Matching is performed.. Truncated to 64 bytes, end truncated. |
| Patient's Sex | (0010,0040) | O | 2 * | Yes/Yes | No Matching is performed. Truncated to 2 bytes, end truncated |
| Patient's Weight | (0010,1030) | O | 2* | Yes/No | No Matching is performed. Mapped into image and RDSR instance. Truncated to 16 bytes, end truncated |
| Confidentiality constraint on patient data | (0040,3001) | O | 2 | No/No | No Matching is performed. |
| Patient's Size | (0010,1020) | O | 3* | Yes/No | No Matching is performed. Mapped into image and RDSR instance. Truncated to 16 bytes, end truncated |
| Patient's Address | (0010,1040) | O | 3 | No/No | Not Requested |
| Patient's Telephone Numbers | (0010,2154) | O | 3 | No/No | Not Requested |
| Patient's Age | (0010,1010) | O | 3 | No/No | Not Requested |
| Occupation | (0010,2180) | O | 3 | No/No | Not Requested |
| Patient's Birth Time | (0010,0032) | O | 3 | No/No | Not Requested |
| Patient's Insurance Plan Code Sequence | (0010,0050) | O | 3 | No/No | Not Requested |
| > 'Code Sequence Macro' | | | | | |
| Patient's Primary Language Code Sequence | (0010,0101) | O | 3 | No/No | Not Requested |
| > 'Code Sequence Macro' | | | | | |
| > Patient's Primary Language Code Modifier Sequence | (0010,0102) | O | 3 | No/No | Not Requested |
| >> 'Code Sequence Macro' | | | | | |

**OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | | | |
|--------------------------------|-------------|---|----|--------|--------------------------------------------------------------------------------------------------------------|
| Military Rank | (0010,1080) | O | 3 | No/No | Not Requested |
| Branch of Service | (0010,1081) | O | 3 | No/No | Not Requested |
| Country of Residence | (0010,2150) | O | 3 | No/No | Not Requested |
| Region of Residence | (0010,2152) | O | 3 | No/No | Not Requested |
| Ethnic Group | (0010,2160) | O | 3 | No/No | Not Requested |
| Patient's Religious Preference | (0010,21F0) | O | 3 | No/No | Not Requested |
| Patient Comments | (0010,4000) | O | 3* | Yes/No | No Matching is performed. Mapped into image and RDSR instance. Truncated to 10240 bytes, end truncated |

Note: * in the *Expected Return Key Type* column indicates that this information is displayed on screen, if available

6.4.5.4 Patient Medical

Not applicable

7. MODALITY PERFORMED PROCEDURE STEP IMPLEMENTATION

7.1 INTRODUCTION

This section specifies the use of the DICOM Modality Performed Procedure Step information to be communicated to the Hospital/Radiology information system.

This feature works in conjunction with DICOM Modality Worklist feature, if installed. However the conformance of this feature is independent of Modality Worklist feature. For information on conformance of Modality Worklist feature to DICOM standard please refer to the appropriate section in this document.

7.2 RELATIONSHIP BETWEEN SCHEDULED AND PERFORMED PROCEDURE STEPS

The Workstation supports following relationship between Scheduled Procedure Step and PPS:

1. Supports a one-to-one relationship between Scheduled Procedure Step and PPS
2. Supports one-to-multiple relationship: Append Case with new PPS based on Original SPS
3. Supports Unscheduled Case or Acquisition without MWL Data

7.3 MODALITY PERFORMED PROCEDURE STEP MODULE TABLE

See DICOM PS 3.3 and PS 3.4 for a complete definition of the entities, modules, and attributes.

TABLE 7-1
MODALITY PERFORMED PROCEDURE STEP MODULES

| Module Name | Reference |
|---------------------------------------|------------------|
| SOP Common | 7.4.1 |
| Performed Procedure Step Relationship | 7.4.2 |
| Performed Procedure Step Information | 7.4.3 |
| Image Acquisition Results | 7.4.4 |
| Radiation Dose | 7.4.5 |
| Billing and Material Management Codes | 7.4.6 |

7.4 MODALITY PERFORMED PROCEDURE STEP MODULE DEFINITIONS

Please refer to DICOM Standard PS 3.3. (Information Object Definitions) for a description of each of the attributes contained within the Modality Performed Procedure Step Information Object Definition.

7.4.1 SOP Common Module

**TABLE 7-2
SOP COMMON MODULE ATTRIBUTES**

| Attribute Name | Tag | Type for SCU N-CREATE | Type for SCU N-SET | Use |
|------------------------|-------------|--------------------------|-----------------------|------------|
| Specific Character Set | (0008,0005) | 1C | 1C | ISO_IR 100 |

7.4.2 Performed Procedure Step Relationship Module

**TABLE 7-3
PERFORMED PROCEDURE STEP RELATIONSHIP MODULE ATTRIBUTES**

| Attribute Name | Tag | Type for SCU - N-CREATE | |
|----------------------------------------------|-------------|------------------------------------------------------------------|-------------------------------|
| | | Acquisition without MWL Entry | Acquisition with MWL Entry |
| Scheduled Step Attributes Sequence | (0040,0270) | 1, Only 1 item | 1, Only 1 item |
| >Study Instance UID | (0020,000D) | 1, System generated. | 1, Got from MWL. |
| >Referenced Study Sequence | (0008,1110) | 2, Empty. | 2, Got from MWL, Only 1 item. |
| >>Referenced SOP Class UID | (0008,1150) | 1, Not sent | 1, Got from MWL. |
| >>Referenced SOP Instance UID | (0008,1155) | 1, Not sent | 1, Got from MWL. |
| >Accession Number | (0008,0050) | 2, User input. If user does not provide input values, then empty | 2, Got from MWL. |
| >Placer Order Number/Imaging Service Request | (0040,2016) | 3, Not sent | 3, Not sent |
| >Filler Order Number/Imaging Service Request | (0040,2017) | 3, Not sent | 3, Not sent |
| >Requested Procedure ID | (0040,1001) | 2, User input. If user does not provide input values, then empty | 2, Got from MWL. |
| >Requested Procedure Code Sequence | (0032,1064) | 3, Empty | 3, Not sent |
| >>Code Value | (0008,0100) | 1, Not sent | 1, Not sent |
| >>Coding Scheme Designator | (0008,0102) | 1, Not sent | 1, Not sent |
| >>Code Meaning | (0008,0104) | 1, Not sent | 1, Not sent |
| >Requested Procedure Description | (0032,1060) | 2, User input. If user does not provide input values, then empty | 2, Got from MWL. |
| >Scheduled Procedure Step ID | (0040,0009) | 2, User input. If user does not provide input values, then empty | 2, Got from MWL. |

OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X**GE Healthcare****DICOM CONFORMANCE STATEMENT**

2148033 REV 1

| | | | |
|---------------------------------------|-------------|----------------------------------------------------------------------------------------------------------------------------------------|------------------|
| >Scheduled Procedure Step Description | (0040,0007) | 2, User input. If user does not provide input values, then empty | 2, Got from MWL. |
| >Scheduled Protocol Code Sequence | (0040,0008) | 2, empty | 2, Got from MWL |
| >>Code Value | (0008,0100) | 1, Not sent | 1, Got from MWL |
| >>Coding Scheme Designator | (0008,0102) | 1, Not sent | 1, Got from MWL |
| >>Coding Scheme Version | (0008,0103) | 3, Not sent | 3, Got from MWL |
| >>Code Meaning | (0008,0104) | 3, Not sent | 3, Got from MWL |
| Patient's Name | (0010,0010) | 2, User input. If user does not provide input values, then UNNAMED xxxxx should be assigned to Last Name (xxxxx is a 5-digital number) | 2, Got from MWL. |
| Patient ID | (0010,0020) | 2, User input. If user does not provide input values, then empty | 2, Got from MWL. |
| Issuer of Patient ID | (0010,0021) | 3, empty | 3, Got from MWL. |
| Patient's Birth Date | (0010,0030) | 2, User input. If user does not provide input values, then empty | 2, Got from MWL. |
| Patient's Sex | (0010,0040) | 2, User input. If user does not provide input values, then empty | 2, Got from MWL. |
| Referenced Patient Sequence | (0008,1120) | 2, empty. | 2, empty. |
| >Referenced SOP Class UID | (0008,1150) | 1, Not sent | 1, Not sent |
| >Referenced SOP Instance UID | (0008,1155) | 1, Not sent | 1, Not sent |
| Admission ID | (0038,0010) | 3, Not sent | 3, Not sent |
| Issuer of Admission ID | (0038,0011) | 3, Not sent | 3, Not sent |
| Service Episode ID | (0038,0060) | 3, Not sent | 3, Not sent |
| Issuer of Service Episode ID | (0038,0061) | 3, Not sent | 3, Not sent |
| Service Episode Description | (0038,0062) | 3, Not sent | 3, Not sent |

7.4.3 Performed Procedure Step Information Module

**TABLE 7-4
PERFORMED PROCEDURE STEP INFORMATION MODULE ATTRIBUTES**

| Attribute Name | Tag | Type for SCU N-CREATE | Type for SCU N-SET | Use |
|---------------------------------------------------------------|-------------|------------------------------|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Performed Procedure Step ID | (0040,0253) | 1 | - | Set to the value of 'Scheduled Procedure Step ID' value from MWL. Or, Set to value entered by the user. Or, Set to value generated by the equipment. |
| Performed Station AE Title | (0040,0241) | 1 | - | Got from configure file. |
| Performed Station Name | (0040,0242) | 2 | - | Got from configure file. |
| Performed Location | (0040,0243) | 2 | - | Got from configure file. The value is Hospital Name. |
| Performed Procedure Step Start Date | (0040,0244) | 1 | - | Same as exam start date |
| Performed Procedure Step Start Time | (0040,0245) | 1 | - | Same as exam start time |
| Performed Procedure Step Status | (0040,0252) | 1 | 3 | When PPS start (N-CREATE) message is sent, this element will have the value "IN PROGRESS" . When PPS end (N-SET) message is sent, this element will have either "COMPLETED" or "DISCONTINUED" based on user selection. |
| Performed Procedure Step Description | (0040,0254) | 2 | 3 | This value is Scheduled Procedure Step Description |
| Performed Procedure Type Description | (0040,0255) | 2 | 3 | From "Requested Procedure Description" on screen which is come from MWL. |
| Procedure Code Sequence | (0008,1032) | 2 | 3 | Always set as Empty for local created exam. When using MWL, this value is Requested Procedure Code Sequence. |
| >Code Value | (0008,0100) | 1 | 1 | Only when using MWL sent. |
| >Coding Scheme Designator | (0008,0102) | 1 | 1 | Only when using MWL sent. |
| >Coding Scheme Version | (0008,0103) | 3 | 3 | Only when using MWL sent. |
| >Code Meaning | (0008,0104) | 3 | 3 | Only when using MWL sent. |
| Performed Procedure Step End Date | (0040,0250) | 2 | 3 | Date when PPS end got installed |
| Performed Procedure Step End Time | (0040,0251) | 2 | 3 | Time when PPS end got installed |
| Comments on the Performed Procedure Step | (0040,0280) | 3 | 3 | Not sent |
| Performed Procedure Step Discontinuation Reason Code Sequence | (0040,0281) | 3 | 3 | Used when exam is DISCONTINUED. Provide a way let user to select one |
| >Code Value | (0008,0100) | 1 | 1 | Follow CID 9301 |

| | | | | |
|---------------------------|-------------|---|---|-----------------|
| >Coding Scheme Designator | (0008,0102) | 1 | 1 | Follow CID 9301 |
| >Coding Scheme Version | (0008,0103) | 3 | 3 | Follow CID 9301 |
| >Code Meaning | (0008,0104) | 3 | 3 | Follow CID 9301 |

7.4.4 Image Acquisition Results Module

**TABLE 7-5
IMAGE ACQUISITION RESULTS MODULE ATTRIBUTES**

| Attribute Name | Tag | Type for SCU N-CREATE | Type for SCU N-SET | Use |
|-------------------------------------------------------|-------------|------------------------------|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Modality | (0008,0060) | 1 | - | When user create a local scheduled exam, default value will be set. Got from MWL server. |
| Study ID | (0020,0010) | 2 | - | Got from MWL server or generated by equipment. |
| Performed Protocol Code Sequence | (0040,0260) | 2 | 3 | Always set as Empty. |
| >Code Value | (0008,0100) | 1 | 1 | Not sent |
| >Coding Scheme Designator | (0008,0102) | 1 | 1 | Not sent |
| >Coding Scheme Version | (0008,0103) | 3 | 3 | Not sent |
| >Code Meaning | (0008,0104) | 3 | 3 | Not sent |
| Performed Series Sequence | (0040,0340) | 2 | 3 | N-Set: Will be sent with referenced instances. |
| >Performing Physician's Name | (0008,1050) | 2 | 2 | Got from MWL server. User input. |
| >Protocol Name | (0018,1030) | 1 | 1 | When using MWL, if Scheduled Protocol Sequence Code Meaning is not empty then this value is Scheduled Protocol Sequence Code Meaning, else this value is Scheduled Procedure Step Description. It will be set as "Unknown Protocol" for local created exam. |
| >Operator's Name | (0008,1070) | 2 | 2 | Always set as Empty. |
| >Series Instance UID | (0020,000E) | 1 | 1 | System creates. |
| >Series Description | (0008,103E) | 2 | 2 | This value is Protocol Name. |
| >Retrieve AE Title | (0008,0054) | 2 | 2 | Always set as Empty. |
| > Archive Requested | (0040,A494) | 3 | 3 | Not sent |
| >Referenced Image Sequence | (0008,1140) | 2 | 2 | Will be sent with referenced images. One item for each image created within the series |
| >>Referenced SOP Class UID | (0008,1150) | 1 | 1 | Referenced image SOP Class UID. |
| >>Referenced SOP Instance UID | (0008,1155) | 1 | 1 | Referenced image SOP Instance UID. |
| >Referenced Non-Image Composite SOP Instance Sequence | (0040,0220) | 2 | 2 | Will be sent with referenced instance |

| | | | | |
|-------------------------------|-------------|---|---|---------------------------------------|
| >>Referenced SOP Class UID | (0008,1150) | 1 | 1 | Referenced instance SOP Class UID. |
| >>Referenced SOP Instance UID | (0008,1155) | 1 | 1 | Referenced instance SOP instance UID. |

7.4.5 Radiation Dose Module

**Table 7-6
Radiation Dose Module Attributes**

| Attribute Name | Tag | Type for SCU N-CREATE | Type for SCU N-SET | Use |
|----------------------------------------------|-------------|------------------------------|---------------------------|------------------------------------------------------------------------------------------------------------|
| Anatomic Structure, Space or Region Sequence | (0008,2229) | 3 | 3 | Always set as empty. |
| Total Time of Fluoroscopy | (0040,0300) | 3 | 3 | Got value from DB. |
| Total Number of Exposures | (0040,0301) | 3 | 3 | Got value from DB. |
| Distance Source to Detector (SID) | (0018,1110) | 3 | 3 | Got value from configure file. |
| Distance Source to Entrance | (0040,0306) | 3 | 3 | Not sent |
| Entrance Dose | (0040,0302) | 3 | 3 | Got value from DB. |
| Entrance Dose in mGy | (0040,8302) | 3 | 3 | The same with (0040,0302), but unit is mGy. |
| Exposed Area | (0040,0303) | 3 | 3 | Not sent |
| Image Area Dose Product | (0018,115E) | 3 | 3 | Got value from DB. |
| Comments on Radiation Dose | (0040,0310) | 3 | 3 | Not Used |
| Exposure Dose Sequence | (0040,030E) | 3 | 3 | Zero or more items shall be included in this sequence. The items count will be the shot count. |
| >Radiation Mode | (0018,115A) | 3 | 3 | The value reported will be 'CONTINUOUS' for all X-Ray modes except pulsed mode which will report 'PULSED'. |
| >KVp | (0018,0060) | 3 | 3 | Got value from DB. |
| >X-ray Tube Current in μ A | (0018,8151) | 3 | 3 | Got value from DB. |
| >Exposure Time | (0018,1150) | 3 | 3 | Got value from DB. |
| >Filter Type | (0018,1160) | 3 | 3 | Not sent |
| >Filter Material | (0018,7050) | 3 | 3 | Not sent |

7.4.6 Billing and Material Management Codes Module

**TABLE 7-7
BILLING AND MATERIAL MANAGEMENT CODES MODULE ATTRIBUTES**

| Attribute Name | Tag | Type for SCU N-CREATE | Type for SCU N-SET | Use |
|---------------------------------|-------------|------------------------------|---------------------------|------------|
| Billing Procedure Step Sequence | (0040,0320) | 3 | 3 | Not sent |
| > 'Code Sequence Macro' | | | | |
| Film Consumption Sequence | (0040,0321) | 3 | 3 | Not sent |

| | | | | |
|---------------------------------------|-------------|---|---|----------|
| >Number of Films | (2100,0170) | 3 | 3 | Not sent |
| >Medium Type | (2000,0030) | 3 | 3 | Not sent |
| >Film Size ID | (2010,0050) | 3 | 3 | Not sent |
| Billing Supplies and Devices Sequence | (0040,0324) | 3 | 3 | Not sent |
| >Billing Item Sequence | (0040,0296) | 3 | 3 | Not sent |
| >> 'Code Sequence Macro' | | | | |
| >Quantity Sequence | (0040,0293) | 3 | 3 | Not sent |
| >>Quantity | (0040,0294) | 3 | 3 | Not sent |
| >>Measuring Units Sequence | (0040,0295) | 3 | 3 | Not sent |
| >>> 'Code Sequence Macro' | | | | |

7.5 STANDARD EXTENDED AND PRIVATE DATA ATTRIBUTES

The Workstation supports the Standard and Private Attributes defined in the following sections in Standard Extended MPPS Instances as Type 3 data elements.

7.5.1 Standard Attributes

The Workstation supports the following attributes, not specified in the MPPS IOD, in SOP Instances as Type 3 data elements.

**TABLE 7-8
STANDARD EXTENDED ATTRIBUTES**

| Attribute Name | Tag | Use |
|----------------|-----|-----|
| N/A | N/A | N/A |

7.6 STANDARD EXTENDED AND PRIVATE CONTEXT GROUPS

The Workstation supports coded terminology using Standard Extended, Private, and Configurable Context Groups defined in the following sections.

7.6.1 Standard Extended Context Groups

The Workstation supports the following extensions to standard Context Groups for SOP Instances created by this product. Extensions are indicated by **bold text**.

7.6.1.1 Context ID 9301 Modality PPS Discontinuation Reasons

**CONTEXT ID 9301
MODALITY PPS DISCONTINUATION REASONS
TYPE: EXTENSIBLE VERSION: 20090616**

| Coding Scheme Designator (0008,0102) | Code Value (0008,0100) | Code Meaning (0008,0104) |
|--------------------------------------|------------------------|------------------------------------|
| DCM | 110500 | Doctor cancelled procedure |
| DCM | 110501 | Equipment failure |
| DCM | 110502 | Incorrect procedure ordered |
| DCM | 110503 | Patient allergic to media/contrast |

OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X**GE Healthcare****DICOM CONFORMANCE STATEMENT**

2148033 REV 1

| | | |
|-----|--------|------------------------------------------|
| DCM | 110504 | Patient died |
| DCM | 110505 | Patient refused to continue procedure |
| DCM | 110506 | Patient taken for treatment or surgery |
| DCM | 110507 | Patient did not arrive |
| DCM | 110508 | Patient pregnant |
| DCM | 110509 | Change of procedure for correct charging |
| DCM | 110510 | Duplicate order |
| DCM | 110511 | Nursing unit cancel |
| DCM | 110512 | Incorrect side ordered |
| DCM | 110513 | Discontinued for unspecified reason |
| DCM | 110514 | Incorrect worklist entry selected |
| DCM | 110515 | Patient condition prevented continuing |
| DCM | 110516 | Equipment change |

8. STORAGE COMMITMENT PUSH MODEL IMPLEMENTATION

8.1 STORAGE COMMITMENT PUSH MODEL INFORMATION OBJECT DEFINITION

Please refer to DICOM Part 3 (Information Object Definitions) for a description of each of the attributes contained within the Storage Commitment Information Object.

The Storage Commitment Information Object is used both for N-ACTION Storage Commitment Requests by the SCU and N-EVENT-REPORT Storage Commitment Notifications by the SCP.

8.1.1 STORAGE COMMITMENT MODULE FOR N-ACTION

**TABLE 8-1
STORAGE COMMITMENT MODULE FOR N-ACTION**

| Attribute Name | Tag | SCU Use | SCP Use |
|------------------------------|-------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| Transaction UID | (0008,1195) | Created by equipment, using prefix: 1.2.840.113619.8.432.+MAC +transactionIndicator(70)+Ti meStamp_Current MAC: Number string of device MAC address. | Not used |
| Storage Media File-Set ID | (0088,0130) | Not used | Not used |
| Storage Media File-Set UID | (0088,0140) | Not used | Not used |
| Referenced SOP Sequence | (0008,1199) | Images/RDSR instances might be included here. | Not used |
| >Referenced SOP Class UID | (0008,1150) | | Not used |
| >Referenced SOP Instance UID | (0008,1155) | | Not used |
| >Storage Media File-Set ID | (0088,0130) | Not used | Not used |
| >Storage Media File-Set UID | (0088,0140) | Not used | Not used |

8.1.2 STORAGE COMMITMENT MODULE FOR N-EVENT-REPORT

**TABLE 8-2
STORAGE COMMITMENT MODULE FOR N-EVENT-REPORT**

| Attribute Name | Tag | SCU Use |
|----------------------------|-------------|---------------------------------------------------------------------|
| Transaction UID | (0008,1195) | Value received from SCP, if the value is empty, then return failed. |
| Retrieve AE Title | (0008,0054) | Not used |
| Storage Media File-Set ID | (0088,0130) | Not used |
| Storage Media File-Set UID | (0088,0140) | Not used |

| | | |
|------------------------------|-------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Referenced SOP Sequence | (0008,1199) | Success of storage commitment is reported to the user. This sequence will be used to mark storage status on screen. All the succeeded instance will marked as “Stored” in the system. |
| >Referenced SOP Class UID | (0008,1150) | |
| >Referenced SOP Instance UID | (0008,1155) | |
| >Retrieve AE Title | (0008,0054) | Not Used |
| >Storage Media File-Set ID | (0088,0130) | Not used |
| >Storage Media File-Set UID | (0088,0140) | Not used |
| Failed SOP Sequence | (0008,1198) | All the failed instance will not be marked as “Stored” in the system. |
| >Referenced SOP Class UID | (0008,1150) | |
| >Referenced SOP Instance UID | (0008,1155) | |
| >Failure Reason | (0008,1197) | See Section 8.1.2.1 for the list of processed values. |

8.1.2.1 Processing of Failure Reason when received in a N-Event-Report

When receiving an N-Event-Report request with an Event Type ID equal to 2, meaning that Storage Commitment is complete, but failure exists, following is the set of value that this Storage Commitment SCU AE is able to process:

| Failure Reason | Meaning | Application Behavior When Receiving Reason Code |
|-----------------------|------------------------------------|--------------------------------------------------------------------------|
| 0110H | Processing failure | Failure and reason is logged. A general failure is reported to the user. |
| 0112H | No such object instance | |
| 0213H | Resource limitation | |
| 0122H | Referenced SOP Class not supported | |
| 0119H | Class / Instance conflict | |
| 0131H | Duplicate transaction UID | |
| * | Other Failure Reason code values | |

9. BASIC DIRECTORY INFORMATION OBJECT IMPLEMENTATION

9.1 IOD MODULE TABLE

Table 9-1 identifies the defined modules within the entities which comprise the Basic Directory IOD. Modules are identified by Module Name.

See DICOM Part 3 for a complete definition of the entities, modules, and attributes.

TABLE 9-1
BASIC DIRECTORY IOD MODULES

| Entity Name | Module Name | Reference |
|-------------------------|-------------------------|-----------|
| File Set Identification | File Set Identification | 9.2.1 |
| Directory Information | Directory Information | 9.2.2 |

The FSC of this implementation creates a Directory Information Module for new media and creates XA IOD.

9.2 INFORMATION MODULE DEFINITIONS

Please refer to DICOM Standard Part 3 (Information Object Definitions) for a description of each of the entities and modules contained within the Basic Directory Information Object.

The following modules are included to convey Enumerated Values, Defined Terms, and Optional Attributes supported. Type 1 & Type 2 Attributes are also included for completeness and to define what values they may take and where these values are obtained from. It should be noted that they are the same ones as defined in the DICOM Standard Part 3 (Information Object Definitions). Also note that Attributes not present in tables are not supported.

9.2.1 File Set identification Module

**TABLE 9-2
FILE-SET IDENTIFICATION MODULE**

| Attribute Name | Tag | Type | Attribute Description |
|----------------------------------------------------|-------------|------|------------------------------------------------------------------------------------------------------------------------------------------------|
| File-set ID | (0004,1130) | 2 | Current DateTime stamp, format like:YYYYMMDDHHMMSS YYYY- year, MM- month, DD- day. HH – Hour 00-23, MM – minute 00-59, SS – second 00-59 |
| File-set Descriptor File ID | (0004,1141) | 3 | Not Used |
| Specific Character Set of File-set Descriptor File | (0004,1142) | 1C | Not Used |

9.2.2 Directory Information Module

**TABLE 9-3
DIRECTORY INFORMATION MODULE**

| Attribute Name | Tag | Type | Attribute Description |
|-------------------------------------------------------------------|-------------|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Offset of the First Directory Record of the Root Directory Entity | (0004,1200) | 1 | The number of bytes from the beginning of the file to the first byte of the first Directory Record. When the Root Directory entity contains no directory record, this offset is set to 00000000H |
| Offset of the Last Directory Record of the Root Directory Entity | (0004,1202) | 1 | The number of bytes from the beginning of the file to the first byte of the last Directory Record. When the Root Directory entity contains no directory record, this offset is set to 00000000H |
| File-set Consistency Flag | (0004,1212) | 1 | FSC writes 0000H: no known inconsistencies. |
| Directory Record Sequence | (0004,1220) | 2 | FSC creates items in this sequence (Directory Records). |
| >Offset of the Next Directory Record | (0004,1400) | 1 | Is set |
| >Record In-use Flag | (0004,1410) | 1 | Set the value to 0xFFFF for record is in use |
| >Offset of Referenced Lower-Level Directory Entity | (0004,1420) | 1 | Is set |
| >Directory Record Type | (0004,1430) | 1 | Enumerated Values are created by an FSC/FSU: PATIENT STUDY SERIES IMAGE SR DOCUMENT |
| >Private Record UID | (0004,1432) | 1C | Not Used |
| >Referenced File ID | (0004,1500) | 1C | Filename is generated to be unique on the media like: DCM<index>\IMG<number of image> Field is included if Directory Record Type is not PATIENT/STUDY/SERIES. |

| | | | |
|----------------------------------------------------|-------------|----|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| >Referenced SOP Class UID in File | (0004,1510) | 1C | Field is included if Directory Record Type is not PATIENT/STUDY/SERIES. Unique ID for the SOP Class of the Instance stored in the referenced File. |
| >Referenced SOP Instance UID in File | (0004,1511) | 1C | Field is included if Directory Record Type is not PATIENT/STUDY/SERIES. Unique Identifier for the SOP Instance stored in the referenced file. |
| >Referenced Transfer Syntax UID in File | (0004,1512) | 1C | Field is included if Directory Record Type is not PATIENT/STUDY/SERIES. Unique Identifier for the Transfer Syntax used to encode the Instance stored in the referenced file. |
| > Referenced Related General SOP Class UID in File | (0004,151A) | 1C | Not Used |
| >Record Selection Keys | | | See 9.2.3 |

9.2.3 Definition of Specific Directory Records

9.2.3.1 Patient Directory Record Definition

**TABLE 9-4
PATIENT KEYS**

| Key | Tag | Type | Attribute Description |
|------------------------|-------------|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Specific Character Set | (0008,0005) | 1C | ISO_IR 100 |
| Patient's Name | (0010,0010) | 2 | FSC/FSU will fill this attribute with the value in the referenced composite object instance. |
| Patient ID | (0010,0020) | 1 | FSC/FSU will fill this attribute with the value in the referenced composite object instance. FSC/FSU will fill this mandatory attribute with a value "DCMTKPAT + 6 digits" if it is empty in the referenced composite object instance. |

9.2.3.2 Study Directory Record Definition

**TABLE 9-5
STUDY KEYS**

| Key | Tag | Type | Attribute Description |
|------------------------|-------------|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Specific Character Set | (0008,0005) | 1C | ISO_IR 100 |
| Study Date | (0008,0020) | 1 | FSC/FSU will fill this mandatory attribute with the value in the referenced composite object instance. This value is mandatory for Workstation, can't be empty. |
| Study Time | (0008,0030) | 1 | FSC/FSU will fill this mandatory attribute with the value in the referenced composite object instance. This value is mandatory for Workstation, can't be empty. |

| | | | |
|--------------------|-------------|----|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Study Description | (0008,1030) | 2 | FSC/FSU will fill this attribute with the value in the referenced composite object instance. |
| Study Instance UID | (0020,000D) | 1C | FSC/FSU will fill this mandatory attribute with the value in the referenced composite object instance. This value is mandatory for Workstation, can't be empty. |
| Study ID | (0020,0010) | 1 | FSC/FSU will fill this mandatory attribute with the value in the referenced composite object instance. FSC/FSU will fill this mandatory attribute with a value "DCMTKSTUDY+ 6 digits" if it is empty in the referenced composite object instance. |
| Accession Number | (0008,0050) | 2 | FSC/FSU will fill this attribute using the value in the referenced composite object instance. |

9.2.3.3 Series Directory Record Definition

**TABLE 9-6
SERIES KEYS**

| Key | Tag | Type | Attribute Description |
|------------------------------|-------------|-------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Specific Character Set | (0008,0005) | 1C | ISO_IR 100 |
| Modality | (0008,0060) | 1 | This will be consistent with Image Content |
| Series Instance UID | (0020,000E) | 1 | FSC/FSU will fill this mandatory attribute with the value in the referenced composite object instance. This value is mandatory for Workstation, can't be empty. |
| Series Number | (0020,0011) | 1 | FSC/FSU will fill this mandatory attribute with the value in the referenced composite object instance. FSC/FSU will fill this mandatory attribute with a value start from 1 (+1 for each series) if it is empty in the referenced composite object instance. |
| Icon Image Sequence | (0088,0200) | 3 | Not Used |
| >Include 'Image Pixel Macro' | | | See 9.2.3.3.1 |

9.2.3.3.1 Icon Image Key Definition

Not applicable

9.2.3.4 Image Directory Record Definition

TABLE 9-7
IMAGE KEYS

| Key | Tag | Type | Attribute Description |
|------------------------------|-------------|------|--------------------------------------------------------------------------------------------------------|
| Specific Character Set | (0008,0005) | 1C | ISO_IR 100 |
| Image Type | (0008,0008) | 3 | FSC/FSU will fill this mandatory attribute with the value in the referenced composite object instance. |
| Instance Number | (0020,0013) | 1 | FSC/FSU will fill this mandatory attribute with the value in the referenced composite object instance. |
| Icon Image Sequence | (0088,0200) | 3 | FSC/FSU will create Icon Images at the image level. |
| >Include 'Image Pixel Macro' | | | See section 9.2.3.4.1 |

- Note: Please also refer to Section 3.2.1.2.1.1.1 for information about Additional Keys sent in DICOMDIR.

9.2.3.4.1 Icon Image Key Definition

- The FSC creates the icon image sequence from the referenced composite object instance. The FSR displays image icons to the user without reading the actual image pixel from the media.
 - Samples Per pixel (0028,0002) is set to 1.
 - Photometric Interpretations (0028,0004) is set to MONOCHROME2.
 - Row (0028,0010) and Column (0028,0011) are set to 128,128.
 - Bits Allocated (0028,0100) and Bits Stored (0028,0101) are created and set to 8.
 - High Bit (0028,0102) is created and set to 7.
 - Pixel Representation (0028,0103) is created and set to 0.
- Pixel Data (7EF0,0010) is set with the 128x128 icon image data.

9.2.3.5 SR Document Directory Record Definition

**TABLE 9-8
SR DOCUMENT KEYS**

| Key | Tag | Type | Attribute Description |
|---------------------------------------------|-------------|------|-----------------------------------------------------------------------------------------------------------------------------------|
| Specific Character Set | (0008,0005) | 1C | ISO_IR 100 |
| Instance Number | (0020,0013) | 1 | FSC/FSU will fill this mandatory attribute with the value in the referenced composite object instance. |
| Completion Flag | (0040,A491) | 1 | "COMPLETE" |
| Verification Flag | (0040,A493) | 1 | "Unverified". |
| Content Date | (0008,0023) | 1 | First Shot Date |
| Content Time | (0008,0033) | 1 | First Shot Time |
| Verification DateTime | (0040,A030) | 1C | Not Used |
| Concept Name Code Sequence | (0040,A043) | 1 | One item shall be included in this sequence. |
| >Code Value | (0008,0100) | 1 | "113701" |
| >Coding Scheme Designator | (0008,0102) | 1 | "DCM" |
| >Code Meaning | (0008,0104) | 1 | "X-Ray Radiation Dose Report" |
| Content Sequence | (0040,A730) | 1C | |
| >Relationship Type | (0040,A010) | 1 | All, and only, Content Items with the HAS CONCEPT MOD relationship from the root Content Item shall be included in this Sequence. |
| <i>>Include 'Document Content Macro'</i> | | | |
| | | | |

9.3 PRIVATE DATA DICTIONARY

Not applicable.

10. PRINT MANAGEMENT IMPLEMENTATION

10.1 INTRODUCTION

This section of the DICOM Conformance Statement specifies the implementation for the specific SOP Classes supported in the Basic Grayscale Print Management Meta SOP Class, the attributes supported for both IODs and services, and the valid range of values for mandatory and optional attributes.

10.2 BASIC FILM SESSION SOP CLASS

10.2.1 Basic Film Session N-Create Attributes

This table lists the attributes that are sent in the Basic Film Session N-Create Request:

| Attribute name | Tag | Use |
|------------------------|-------------|--------------------------------------------------------------------------------------------------|
| Specific Character Set | (0008,0005) | Not Used |
| Number of Copies | (2000,0010) | Range of this value is 1-10, default is 1. |
| Print Priority | (2000,0020) | Enumerated Terms can be sent: HIGH or MED or LOW |
| Medium Type | (2000,0030) | Enumerated Terms can be sent: PAPER, CLEAR FILM, BLUE FILM, MAMMO CLEAR FILM, MAMMO BLUE FILM |
| Film Destination | (2000,0040) | Enumerated Terms can be sent: MAGAZINE or PROCESSOR |
| Film Session Label | (2000,0050) | Set to "OEC One Print SCU" |
| Memory Allocation | (2000,0060) | Not Used |
| Owner ID | (2100,0160) | Set to "OEC One" |

10.2.2 Basic Film Session N-Delete

The N-DELETE is used to delete the complete Basic Film Session SOP Instance hierarchy. As a result, all references to Image SOP Instances within the film session are deleted.

10.3 BASIC FILM BOX SOP CLASS

10.3.1 Basic Film Box N-Create Attributes

This table lists the attributes that are sent to the SCP in the Basic Film Box N-Create Request, and that are received in the Basic Film Box N-Create Response from the SCP

| Attribute Name | Tag | Use |
|----------------------|-------------|--------------------|
| Image Display Format | (2010,0010) | User configurable. |

**OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | |
|------------------------------------------|-------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| | | User values =: STANDARD\1,1 STANDARD\1,2 STANDARD\2,1 STANDARD\2,2 STANDARD\2,3 STANDARD\3,3 STANDARD\3,4 STANDARD\4,4 STANDARD\4,5 STANDARD\5,4 Default value: First selection when declaring printer. |
| Referenced Film Session Sequence | (2010,0500) | Used |
| >Referenced SOP Class UID | (0008,1150) | Set to 1.2.840.10008.5.1.1.1 |
| >Referenced SOP Instance UID | (0008,1155) | This value should be returned from SCP in N-CREATE-RSP Basic Film Session SOP Class; Then this value is sent to SCP in N-CREATE-RQ Basic Film Box SOP Class requesting for Basic Grayscale Image Box SOP Class Instance UID which contained in Referenced Image Box Sequence from SCP. |
| Referenced Image Box Sequence | (2010,0510) | Used |
| >Referenced SOP Class UID | (0008,1150) | Set to 1.2.840.10008.5.1.1.4 |
| >Referenced SOP Instance UID | (0008,1155) | Received in the N-Create Response from SCP, this value should be sent to SCP in N-SET-RQ Basic Grayscale Image Box SOP Class |
| Referenced Basic Annotation Box Sequence | (2010,0520) | Not Used |
| >Referenced SOP Class UID | (0008,1150) | Not Used |
| >Referenced SOP Instance UID | (0008,1155) | Not Used |
| Film Orientation | (2010,0040) | Enumerated Terms can be sent: PORTRAIT or LANDSCAPE |
| Film Size ID | (2010,0050) | Enumerated Terms can be sent: 8INX10IN 10INX12IN 10INX14IN 11INX14IN 14INX14IN 14INX17IN 24CMX24CM 24CMX30CM |
| Magnification Type | (2010,0060) | Set to BILINEAR |
| Max Density | (2010,0130) | User configurable. The default is 300 |
| Configuration Information | (2010,0150) | User configurable. This value specifies the DICOM printer parameters. It is defined in the DICOM Printer's Conformance Statement file. |

| | | |
|--------------------------------------|-------------|-------------------------------------------------|
| Referenced Presentation LUT Sequence | (2050,0500) | Not Used |
| >Referenced SOP Class UID | (0008,1150) | Not Used |
| >Referenced SOP Instance UID | (0008,1155) | Not Used |
| Annotation Display Format ID | (2010,0030) | Not Used |
| Smoothing Type | (2010,0080) | Not Used |
| Border Density | (2010,0100) | Enumerated Terms can be sent: BLACK or WHITE |
| Empty Image Density | (2010,0110) | Enumerated Terms can be sent: BLACK or WHITE |
| Min Density | (2010,0120) | User configurable. The default is 20 |
| Trim | (2010,0140) | Set to NO |
| Illumination | (2010,015E) | Not Used |
| Reflected Ambient Light | (2010,0160) | Not Used |
| Requested Resolution ID | (2020,0050) | Not Used |
| ICC Profile | (0028,2000) | Not Used |

10.3.2 Basic Film Box N-Action Attributes

Following are the Action Reply arguments that are supported if present in the N-Action response of the Basic Film Box SOP Class

| Action Type Name | Action Type ID | Attribute | Tag | Usage SCU |
|------------------|----------------|-------------------------------|-------------|-----------|
| Print | 1 | Referenced Print Job Sequence | (2100,0500) | Not Used |
| | | >Referenced SOP Class UID | (0008,1150) | Not Used |
| | | >Referenced SOP Instance UID | (0008,1155) | Not Used |
| | | >Print Job ID | (2100,0010) | Not Used |

10.3.3 Basic Film Box N-Delete

The N-DELETE is used to delete the last created Basic Film Box SOP Instance hierarchy. As a result all the information describing the last film is deleted.

10.4 BASIC GRAYSCALE IMAGE BOX SOP CLASS

10.4.1 Basic Grayscale Image Box Pixel N-Set Attributes

This table lists the attributes that are sent in the Basic Grayscale Image Box N-Set Request:

| Attribute Name | Tag | Use |
|----------------|-------------|------------------------------------------------------------------------------|
| Image Position | (2020,0010) | Based on Image Display Format (2010,0010), range of values sent is [1 ~ 20]. |

**OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | |
|--------------------------------------|-------------|---------------------------------------------------------------------------------------|
| Basic Grayscale Image Sequence | (2020,0110) | A sequence which provides the content of the grayscale image pixel data to be printed |
| >Samples Per Pixel | (0028,0002) | 1 |
| >Photometric Interpretation | (0028,0004) | MONOCHROME2 |
| >Rows | (0028,0010) | Value is 1000 Note: For dose summary snapshot, this value will be 960. |
| >Columns | (0028,0011) | Value is 1000 Note: For dose summary snapshot, this value will be 960. |
| >Pixel Aspect Ratio | (0028,0034) | Not Used |
| >Bits Allocated | (0028,0100) | Value is 8 |
| >Bits Stored | (0028,0101) | Value is 8 |
| >High Bit | (0028,0102) | Value is 7 |
| >Pixel Representation | (0028,0103) | 0 |
| >Pixel Data | (7FE0,0010) | Image Pixels |
| Polarity | (2020,0020) | Not Used |
| Magnification Type | (2010,0060) | Not Used |
| Smoothing Type | (2010,0080) | Not Used |
| Min Density | (2010,0120) | Not Used |
| Max Density | (2010,0130) | Not Used |
| Configuration Information | (2010,0150) | Not Used |
| Requested Image Size | (2020,0030) | Not Used |
| Requested Decimate/Crop Behavior | (2020,0040) | Not Used |
| Referenced Presentation LUT Sequence | (2050,0500) | Not Used |
| > Referenced SOP Class UID | (0008,1150) | Not Used |
| > Referenced SOP Instance UID | (0008,1155) | Not Used |

10.5 PRINTER SOP CLASS

10.5.1 Printer N-Event-Report Attributes

The following table describes the Workstation behavior when receiving an **N-Event-Report** request from the Printer SCP depending on the Event Type ID value.

| Event Type Name | Event Type ID | Attribute | Tag | Use |
|-----------------|---------------|---------------------|-------------|---------|
| Normal | 1 | | | Ignored |
| Warning | 2 | Printer Status Info | (2110,0020) | Ignored |
| | | Film Destination | (2000,0040) | Ignored |
| | | Printer Name | (2110,0030) | Ignored |
| Failure | 3 | Printer Status Info | (2110,0020) | Ignored |
| | | Film Destination | (2000,0040) | Ignored |
| | | Printer Name | (2110,0030) | Ignored |

GE Healthcare

2148033 REV 1

10.5.2 Printer N-Get Attributes

| Attribute Name | Tag | Use |
|--------------------------|-------------|------------------------------------------------------------------------------------------------------------------------|
| Printer Status | (2110,0010) | When status is FAILURE, we'll not execute print any more and inform user print fail and the printer status is abnormal |
| Printer Status Info | (2110,0020) | Ignored |
| Printer Name | (2110,0030) | Ignored |
| Manufacturer | (0008,0070) | Ignored |
| Manufacturer Model Name | (0008,1090) | Ignored |
| Device Serial Number | (0018,1000) | Ignored |
| Software Versions | (0018,1020) | Ignored |
| Date Of Last Calibration | (0018,1200) | Ignored |
| Time Of Last Calibration | (0018,1201) | Ignored |

11. QUERY IMPLEMENTATION

11.1 WORKSTATION MAPPING OF DICOM ENTITIES

The Workstation maps DICOM Information Entities to local Information Entities in the product’s database and user interface.

TABLE 11-1
MAPPING OF DICOM ENTITIES TO WORKSTATION ENTITIES

| DICOM | Workstation Entity |
|---------|--------------------|
| Patient | Exam |
| Study | Exam |
| Series | Exam |
| Image | Image |

11.2 INFORMATION MODEL KEYS

Please refer to DICOM Standard PS 3.4 (Service Class Specifications) for a description of each of the levels contained within the Query/Retrieve Information Model.

11.2.1 Common Query Keys

The query key attributes specified in this section are used at all levels and in all classes of query.

TABLE 11-2
Q/R PATIENT LEVEL COMMON RETRIEVE ATTRIBUTES

| Attribute Name | Tag | Type | SCU Use |
|----------------------------|-------------|------|----------------------------------------------------------------------------------------|
| Specific Character Set | (0008,0005) | - | See 11.2.1.1.1 |
| Query Retrieve Level | (0008,0052) | - | Set to level of query: PATIENT STUDY SERIES IMAGE |
| Retrieve AE Title | (0008,0054) | - | Returned value will be used for Retrieve AE Title of retrieval of query response item. |
| Storage Media File-set ID | (0088,0130) | - | Attribute is not requested. |
| Storage Media File-set UID | (0088,0140) | - | Attribute is not requested. |

11.2.1.1 Q/R Common Attribute Descriptions

11.2.1.1.1 Specific Character Set

In Query requests the character set tag (0008,0005) is set to “ISO_IR 100” for all requests.

In processing of query responses the character set tag (0008,0005) is checked: absence of the tag is defined as ISO_IR 6 (7 bit ASCII) or the tag is present ISO_IR 100 (ISO-8859-1 Latin 1) will be accepted.

11.2.2 Patient Level

This section defines the keys at the Patient Level of the Patient Root Query/Retrieve Information Models that are supported by this implementation.

**TABLE 11-3
PATIENT LEVEL ATTRIBUTES FOR THE PATIENT ROOT
QUERY/RETRIEVE INFORMATION MODEL**

| Attribute Name | Tag | Type | SCU Use |
|-------------------------------------|-------------|-------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Patient's Name | (0010,0010) | R* | Single value and wildcard matching are supported. User could enter anything (including wildcard) except "^", "=", and "\". The value entered by the user will be used as matching key, but the response from server depends on server's capability. During the DICOM archive/sending, the "^" and "=" will be used as delimiter (added by system) in Patient Name components. |
| Patient ID | (0010,0020) | U* | Single value and wildcard matching are supported. User could enter anything (including wildcard) except "\". The value entered by the user will be used as matching key, but the response from server depends on server's capability. |
| Referenced Patient Sequence | (0008,1120) | O | Not Requested |
| >Referenced SOP Class UID | (0008,1150) | O | Not Requested |
| >Referenced SOP Instance UID | (0008,1155) | O | Not Requested |
| Patient's Birth Date | (0010,0030) | O | Universal Matching |
| Patient's Birth Time | (0010,0032) | O | Not Requested |
| Patient's Sex | (0010,0040) | O | Universal Matching |
| Other Patient IDs | (0010,1000) | O | Not Requested |
| Other Patient Names | (0010,1001) | O | Not Requested |
| Ethnic Group | (0010,2160) | O | Not Requested |
| Patient Comments | (0010,4000) | O | Not Requested |
| Number of Patient Related Studies | (0020,1200) | O | Universal Matching |
| Number of Patient Related Series | (0020,1202) | O | Universal Matching |
| Number of Patient Related Instances | (0020,1204) | O | Universal Matching |

Note: * in the *Type* column indicates that this information is displayed on screen, if available

11.2.3 Study Level – Patient Root

This section defines the keys at the Study Level of the Patient Root Query/Retrieve Information Models that are supported by this implementation.

TABLE 11-4
STUDY LEVEL ATTRIBUTES FOR THE PATIENT ROOT
QUERY/RETRIEVE INFORMATION MODEL

| Attribute Name | Tag | Type | SCU Use |
|------------------------------------|-------------|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Study Date | (0008,0020) | R* | Range value matching is supported. User could select one of "All", "Today", "Since Yesterday", "Last 3 days", "Last 7 days" and "Last 30 days" from the screen, and then date information will be sent as "yyyymmdd-yyyymmdd" format to the server for matching. |
| Study Time | (0008,0030) | R* | Universal Matching |
| Accession Number | (0008,0050) | R* | Single value and wildcard matching are supported. Leading and trailing wildcard matching for value entered. The value entered by the user will be used as matching key. |
| Study ID | (0020,0010) | R | Universal Matching |
| Study Instance UID | (0020,000D) | U | Universal Matching |
| Modalities in Study | (0008,0061) | O | Universal Matching |
| SOP Classes in Study | (0008,0062) | O | Not Requested |
| Referring Physician's Name | (0008,0090) | O | Not Requested |
| Study Description | (0008,1030) | O* | Universal Matching |
| Procedure Code Sequence | (0008,1032) | O | Not Requested |
| >Code Value | (0008,0100) | O* | Not Requested |
| >Coding Scheme Designator | (0008,0102) | O | Not Requested |
| >Coding Scheme Version | (0008,0103) | O | Not Requested |
| >Code Meaning | (0008,0104) | O | Not Requested |
| Name of Physician(s) Reading Study | (0008,1060) | O | Not Requested |
| Admitting Diagnoses Description | (0008,1080) | O | Not Requested |
| Referenced Study Sequence | (0008,1110) | O | Not Requested |
| >Referenced SOP Class UID | (0008,1150) | O | Not Requested |
| >Referenced SOP Instance UID | (0008,1155) | O | Not Requested |
| Issuer of Patient ID | (0010,0021) | O | Not Requested |

**OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | |
|-----------------------------------|-------------|---|--------------------|
| Patient's Age | (0010,1010) | O | Not Requested |
| Patient's Size | (0010,1020) | O | Not Requested |
| Patient's Weight | (0010,1030) | O | Not Requested |
| Occupation | (0010,2180) | O | Not Requested |
| Additional Patient History | (0010,21B0) | O | Not Requested |
| Other Study Numbers | (0020,1070) | O | Not Requested |
| Number of Study Related Series | (0020,1206) | O | Universal Matching |
| Number of Study related Instances | (0020,1208) | O | Universal Matching |
| Interpretation Author | (4008,010C) | O | Not Requested |

Note: * in the *Type* column indicates that this information is displayed on screen, if available

11.2.4 Study Level – Study Root

This section defines the keys at the Study Level of the Study Root Query/Retrieve Information Model that are supported by this implementation.

**TABLE 11-5
STUDY LEVEL ATTRIBUTES FOR THE STUDY ROOT
QUERY/RETRIEVE INFORMATION MODEL**

| Attribute Name | Tag | Type | SCU Use |
|-----------------------|-------------|-------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Study Date | (0008,0020) | R* | <p>Range value matching is supported.</p> <p>User could select one of "All", "Today", "Since Yesterday", "Last 3 days", "Last 7 days" and "Last 30 days" from the screen, and then date information will be sent as "yyyymmdd-yyyymmdd" format to the server for matching.</p> <p>All: All matches.</p> <p>Today: Date Range Matches "Today – Today"</p> <p>Since Yesterday: Date Range Matches "Yesterday - Today"</p> <p>Last 3 days: Date Range Matches "2 days before – Today"</p> <p>Last 7 days: Date Range Matches "6 days before – Today"</p> <p>Last 30 days: Date Range Matches "29 days before – Today"</p> <p>Default: Last 3 days.</p> |
| Study Time | (0008,0030) | R* | <p>Universal Matching.</p> <p>Always blank in the request.</p> |
| Accession Number | (0008,0050) | R* | <p>Single value and wildcard matching are supported.</p> <p>The value entered by the user will be used as matching key.</p> |

**OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | |
|------------------------------------|-------------|----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Study ID | (0020,0010) | R | Universal Matching |
| Patient's Name | (0010,0010) | R* | Single value and wildcard matching are supported. User could enter anything (including wildcard) except "^", "=", and "\". The value entered by the user will be used as matching key, but the response from server depends on server's capability. During the DICOM archive/sending, the "^" and "=" will be used as delimiter (added by system) in Patient Name components. |
| Patient ID | (0010,0020) | R* | Single value and wildcard matching are supported. User could enter anything (including wildcard) except "\". |
| Study Instance UID | (0020,000D) | U | Universal Matching |
| Modalities in Study | (0008,0061) | O | Universal Matching |
| SOP Classes in Study | (0008,0062) | O | Not Requested |
| Referring Physician's Name | (0008,0090) | O | Not Requested |
| Study Description | (0008,1030) | O* | Universal Matching |
| Procedure Code Sequence | (0008,1032) | O | Not Requested |
| >Code Value | (0008,0100) | O | Not Requested |
| >Coding Scheme Designator | (0008,0102) | O | Not Requested |
| >Coding Scheme Version | (0008,0103) | O | Not Requested |
| >Code Meaning | (0008,0104) | O | Not Requested |
| Name of Physician(s) Reading Study | (0008,1060) | O | Not Requested |
| Admitting Diagnoses Description | (0008,1080) | O | Not Requested |
| Referenced Study Sequence | (0008,1110) | O | Not Requested |
| >Referenced SOP Class UID | (0008,1150) | O | Not Requested |
| >Referenced SOP Instance UID | (0008,1155) | O | Not Requested |
| Issuer of Patient ID | (0010,0021) | O | Not Requested |
| Patient's Age | (0010,1010) | O | Not Requested |
| Patient's Size | (0010,1020) | O | Not Requested |
| Patient's Weight | (0010,1030) | O | Not Requested |
| Occupation | (0010,2180) | O | Not Requested |
| Additional Patient History | (0010,21B0) | O | Not Requested |
| Other Study Numbers | (0020,1070) | O | Not Requested |

**OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | |
|-------------------------------------|-------------|---|--------------------|
| Number of Study Related Series | (0020,1206) | O | Universal Matching |
| Number of Study related Instances | (0020,1208) | O | Universal Matching |
| Interpretation Author | (4008,010C) | O | Not Requested |
| Referenced Patient Sequence | (0008,1120) | O | Not Requested |
| >Referenced SOP Class UID | (0008,1150) | O | Not Requested |
| >Referenced SOP Instance UID | (0008,1155) | O | Not Requested |
| Patient's Birth Date | (0010,0030) | O | Universal Matching |
| Patient's Birth Time | (0010,0032) | O | Not Requested |
| Patient's Sex | (0010,0040) | O | Universal Matching |
| Other Patient IDs | (0010,1000) | O | Not Requested |
| Other Patient Names | (0010,1001) | O | Not Requested |
| Ethnic Group | (0010,2160) | O | Not Requested |
| Patient Comments | (0010,4000) | O | Not Requested |
| Number of Patient Related Studies | (0020,1200) | O | Not Requested |
| Number of Patient Related Series | (0020,1202) | O | Not Requested |
| Number of Patient Related Instances | (0020,1204) | O | Not Requested |

Note: * in the *Type* column indicates that this information is displayed on screen, if available

11.2.5 Series Level

This section defines the keys at the Series Level of the Patient Root and Study Root Query/Retrieve Information Models that are supported by this implementation.

**TABLE 11-6
SERIES LEVEL ATTRIBUTES FOR THE
QUERY/RETRIEVE INFORMATION MODEL**

| Attribute Name | Tag | Type | SCU Use |
|---------------------|-------------|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Modality | (0008,0060) | R* | Single value matching is supported. User could select one of "XA", "RF", "CR", "SC" and "ALL" for Modality Type from screen, then the corresponding string "XA", "RF", "CR", "OT" and "" will be sent to server for matching. Note: "ALL" means ignore the Modality Type, the Modality Type will leave empty during sending. |
| Series Number | (0020,0011) | R | Universal Matching |
| Series Instance UID | (0020,000E) | U | Universal Matching |

**OEC ONE WORKSTATION SOFTWARE RELEASE WKS-2.X
DICOM CONFORMANCE STATEMENT**

GE Healthcare
2148033 REV 1

| | | | |
|----------------------------------------------|--------------|---|--------------------|
| Number of Series Related Instances | (0020,1209) | O | Universal Matching |
| Series Description | (0008,103E) | O | Universal Matching |
| Performed Procedure Step ID | (0040, 0253) | O | Not Requested |
| Referenced Performed Procedure Step Sequence | (0008,1111) | O | Not Requested |
| >Referenced SOP Class UID | (0008,1150) | O | Not Requested |
| >Referenced SOP Instance UID | (0008,1155) | O | Not Requested |
| Request Attribute Sequence | (0040, 0275) | O | Not Requested |
| >Requested Procedure ID | (0040,1001) | O | Not Requested |
| >Scheduled Procedure Step ID | (0040,0009) | O | Not Requested |
| Performed Procedure Step Start Date | (0040,0244) | O | Not Requested |
| Performed Procedure Step Start Time | (0040,0245) | O | Not Requested |

Note: * in the *Type* column indicates that this information is displayed on screen, if available

11.2.6 Image Level

Currently no Image level query requests are made.