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Bevorzugter Name	Synonym
Alkylethylbenzoldimethylammoniumchlorid	NB
Alkoniumchlorid	Alkyldimethyl-3,4-dichlorbenzyl ammoniumchlorid
Benzalkoniumchlorid	Alkyldimethylbenzylammoniumchlorid
Cetalkoniumchlorid	Cetyldimethylbenzylammoniumchlorid
Didecyldimethylammoniumchlorid	Chloriddidecyldimethylammonium
Diocetyldimethylammoniumchlorid	Chloriddioctyldimethylammonium
Hexadecyldimethylbenzylammoniumchlorid	Chloridhexadecyldimethylbenzylammonium
Methyldodecylbenzyltrimethyl	Chloridmethyldodecylbenzyltrimethylammoniumchlorid Ammonium
Octadecyldimethylbenzylammoniumchlorid	Chloridoctadecyldimethylbenzyl ammonium
Octyldecyldimethylammoniumchlorid	Chloridoctyldecyldimethylammonium
Octyldimethylammoniumchlorid	Chloridoctyldimethylammonium

## Bibliographie

- 1) Angus, D. and R. Wax, „Epidemiology of sepsis: An update“, *Crit Care Med*, 2001. 29(Suppl): p. S109-S116
- 2) Martin, G., Mannino, D., Eaton, S. and Moss, M. „The epidemiology of sepsis in the United States from 1979 through 2000“, *N Engl J Med*. 348:1546-54, 2003
- 3) Parrillo, J.E., Parker, M.M., Natanson, C., Suffredini, A.F., Danner, R.L., Cunnion, R.E. and Ognibene, F.P. „Septic shock in humans“, *Annals of Internal Medicine*. 113:227-242, 1990
- 4) Burke, J. „Infection control - A problem for patient safety“, *N Engl J Med*. 348: 651-656, 2003
- 5) Weinstein, R. „Nosocomial infection update“, *Emerg Infect Dis*. 4:416-420, 1998
- 6) Jarvis, W. „Infection control and changing health-care delivery systems“, *Emerg Infect Dis*. 7:170-173, 2001
- 7) Maki, D. and Brookmeyer, P. „A survey of EKG telemetry harnesses as a reservoir of resistant noscomial pathogens“, *American Society for Microbiology* (Herausgeber). Zugriff: 9. Dezember 2003 unter [www.asmta.org/memonly/abstracts/AbstractView.asp?AbstractID=82845](http://www.asmta.org/memonly/abstracts/AbstractView.asp?AbstractID=82845)
- 8) Borer, A., Smolyakov, R., Eskira, S., Peled, N., Gilad, J., Man, I., Hyam, E., Riesenber, K. and Schlaeffer, F. „Cellular phones (CeP) of personnel as a potential source for nosocomial transmission of *Acinetobacter baumannii* (AB)“, *American Society for Microbiology* (Herausgeber). Zugriff: 9. Dezember 2003 unter [www.asmta.org/memonly/abstracts/AbstractView.asp?AbstractID=82844](http://www.asmta.org/memonly/abstracts/AbstractView.asp?AbstractID=82844)
- 9) Austin, D.J., Bonten, M.J.M., Weinstein, R.A., Slaughter, S. and Anderson, R.M. „Vancomycin-resistant enterococci in intensive-care hospital settings: Transmission dynamics, persistence, and the impact of infection control programs“, *Proc Natl Acad Sci*. 96: 6908-6913, 1999
- 10) Boyce, J.M. and Pittet, D. „Guideline for hand hygiene in health-care settings recommendations of the healthcare infection control practices advisory committee and the HICPAC/SHEA/ APIC/ IDSA hand hygiene task force“, *Centers for Disease Control* (Herausgeber). Zugriff: 9. Dezember 2003 unter [www.cdc.gov/mmwr/PDF/RR/RR5116.pdf](http://www.cdc.gov/mmwr/PDF/RR/RR5116.pdf)
- 11) „CDC releases new hand hygiene guidelines“, *Centers for Disease Control and Prevention* (Herausgeber). Zugriff: 10. Dezember 2003 unter [www.cdc.gov/handhygiene/pressrelease.htm](http://www.cdc.gov/handhygiene/pressrelease.htm)
- 12) Kampf, G. and Löffler, H. „Dermatological aspects of a successful introduction and continuation of alcohol-based hand rubs for hygienic hand disinfection“, *Journal of Hospital Infection*. 55:1-7, 2003
- 13) Rutala, W. „APIC guidelines for infection control practice“, *Association for Professionals in Infection Control and Epidemiology Inc.* (Herausgeber). Zugriff: 17. November 2003 unter [www.apic.org/pdf/gddisinf.pdf](http://www.apic.org/pdf/gddisinf.pdf)
- 14) „Category IV Monograph: Hard Surface Disinfectants“, *Health Canada Online* (Herausgeber). Zugriff: 18. November 2003 unter [www.hc-sc.gc.ca/hpfbdgpsa/tpd-dpt/hardsur\\_e.pdf](http://www.hc-sc.gc.ca/hpfbdgpsa/tpd-dpt/hardsur_e.pdf)

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© 2007 General Electric Company – Alle Rechte vorbehalten.
General Electric Company. GE Healthcare Oy., wird unter