

# Clinical excellence in computerized ECG analysis



“The fact is,  
the electrocardiogram  
is needed to deliver  
modern medical care  
now more than ever.”<sup>1</sup>

J. W. Hurst

## ECG and its role in modern medical care

“There is a misconception that newer diagnostic modalities have replaced the ECG. The fact is, the electrocardiogram (ECG) is needed now more than ever to deliver accurate and cost effective modern medical care.”<sup>2</sup> Furthermore, it is noted in the current American College of Cardiology (ACC)/American Heart Association (AHA) guidelines that: “Recording the resting 12-lead ECG continues to be the most commonly used laboratory procedure for the diagnosis of heart disease.” In addition, “The procedure is safe, simple, and reproducible; the ECG record lends itself to serial studies; and the relative cost is minimal.”<sup>3</sup>

The ECG continues to play a key role in the proper selection of the latest therapies. In the setting of an acute coronary syndrome, the ECG is used for treatment selection and expedited care.<sup>4</sup> The ECG continues to be central in the management of arrhythmias.<sup>5</sup> New advancements in electronic cardiac implants have solidified the role of the ECG in selection of appropriate target populations for such therapies as biventricular pacemakers<sup>6</sup> and implantable cardioverter defibrillators.<sup>7</sup> Even in the presence of superior diagnostic imaging methods, the ECG often provides valuable prognostic data showing that the ECG still has a significant role in determining patient care.<sup>8,9</sup>

## Quality digital waveforms driving clinical excellence

Many physicians can recall instances where a single ECG made all the difference in determining a critical diagnosis. To be the most effective and to maintain clinical excellence, it is critical the ECG enters the digital/paperless domain.

Advancements in computer technology brought digital writers and displays to the clinical environment. This allowed for the computerization (digitization) of the ECG. Digital acquisition of the ECG signal is important to reduce noise and artifact.<sup>10</sup> GE was the first to introduce a digital acquisition system and lead the way in developing sophisticated signal processing programs designed to improve ECG signal (waveform) quality. This processing is critical and helps to maximize the accuracy of computerized measurements, interpretation and waveform display.<sup>11,12</sup>

## Bringing continued innovation to the ECG

GE was the first company to provide digital electrocardiographs, which simultaneously acquired and analyzed all 10 seconds of all 12 leads of ECG waveforms.<sup>13</sup> Introduced in 1982, this innovation had a profound impact on the industry.<sup>14</sup> GE continues its history of new innovations and raising the standards of clinical excellence. (See our Breakthrough Innovation timeline brochure, identifying advances in ECG devices, analysis programs, and system connectivity, along with the impact they had on the establishment of new clinical guidelines and ECG standards.)

## ECG across the continuum of care

GE provides products and systems that incorporate digital ECG processing, in every patient care setting where the ECG is needed. Starting with out-of-hospital emergency care, GE's Marquette® 12SL™ program is implemented in ZOLL® and Medtronic LIFEPAK® pre-hospital defibrillators for recognition of acute coronary syndromes.<sup>15</sup> GE provides a complete set of diagnostic systems for resting ECG, stress testing, Holter, clinical patient monitoring, cardiac catheterization and electrophysiology (EP) testing, all of which include advanced analysis programs.<sup>16,17,18</sup>

## Information management for effective, broad support

GE uses industry-standard networking to connect ECG acquisition systems. All ECG types, whether it be Holter, stress, resting, signal averaged ECG, vectorcardiogram, etc., can be sent to the MUSE® Cardiology Information System. In addition, ECGs acquired from bedside monitors and telemetry units can also be transmitted via the Unity Network®. This network can include a secure wireless capability to provide automatic routing of the ECG to the point-of-decision. Furthermore, the network has HL7 functions for accurate administration of patient demographics, orders, test results and billing-related information. GE's network connectivity has been reported to reduce errors, cost, and report turnaround time resulting in expedited patient care.<sup>19,20</sup>

# Many physicians can recall instances where a single ECG made all the difference in determining a critical diagnosis.

Although GE offers a solution for an Electronic Medical Record (EMR),<sup>21</sup> it also can provide the ECG to other third-party information systems, including the Veterans Administration (VA) Hospital Vista System.<sup>22</sup> GE has a long history of providing quality ECG data to enterprise-wide hospital information, or office-based systems.

GE also supports multiple standards for graphical reports, as in the portable document format (PDF), for the export of ECG data (via SCP-ECG, HL7, DICOM and XML) and the import of ECG data (via GE's Information Interchange program called I<sup>2</sup> or XML).

## Streamlined workflow and real-time decision support

Throughout the continuum of care, GE provides real-time electrocardiograph-based applications that assist the user in obtaining high quality ECGs through tools, such as Marquette Hookup Advisor,<sup>TM23</sup> and clinical decision support tools, such as ACI-TIPI, for chest pain patient management.<sup>24</sup> Additionally, with systems such as MUSE and Unity Network, real-time decision support has been implemented to automatically capture, identify, prioritize, compare, and route critical ECGs.<sup>25</sup>

Real-time decision support is most valuable when a sophisticated analysis program is combined with an information system. This combination facilitates a practical, streamlined workflow, and provides the decision maker with timely patient information to make appropriate medical decisions.

## Arrhythmias and decision support

Major cardiac centers often store information about patient arrhythmias for the purpose of optimizing patient care decisions. These recordings of critical ECG events assist the electrophysiologist in making the most appropriate decisions regarding implant settings, ablation, or other forms of therapy.

Continuous ECG recordings from all bedside monitors or telemetry units are available via the Unity Network, and can be routed to an arrhythmia review workstation, known as MARS® (Multi-parameter Arrhythmia Review Station). At the MARS station, critical arrhythmias can be viewed, sorted and sent to the MUSE system for long-term storage.

Marquette EKPro<sup>TM</sup> is a program used for analyzing continuous ECG recordings from GE bedside monitors, telemetry units or Holter recorders. The Marquette EKPro program automatically labels the data and identifies significant arrhythmic episodes. In 2002, Marquette EKPro was enhanced to also analyze pulsatile data accompanying the ECG,<sup>26</sup> which "leads to a clinically significant reduction in the number of false positive heart rate alarms without reducing detection of serious arrhythmias."<sup>27</sup>

Atrial fibrillation is the most common arrhythmia that results in hospitalization in the United States. For this reason, the Marquette EKPro analysis program has been enhanced to include atrial fibrillation detection and trending in patient monitors. GE has developed advanced ECG processing techniques for the prediction,<sup>28,29</sup> detection<sup>30,31,32</sup> and quantification<sup>33</sup> of atrial fibrillation in many other care areas.

In addition to the detection, cataloging, and transmission of significant arrhythmias for decision support, GE is an industry-leading provider of automated measurements that have been correlated with Sudden Cardiac Death (SCD).

These measurements are generated from the following programs:

- Marquette 12SL ECG Analysis<sup>34</sup>
- Marquette Signal Averaged ECG (SAECG)<sup>35</sup>
- Marquette T-Wave Alternans (TWA) Analysis<sup>36</sup>
- Marquette Heart Rate Variability (HRV) spectral<sup>37</sup> and time domain<sup>38,39</sup>
- Marquette Heart Rate Turbulence (HRT) Analysis<sup>40,41,42</sup>

GE is also a provider of two analysis programs, TWA and HRT specifically used for the prediction of SCD.

GE is a leader in computerized ECG assessment of QT associated measurements, some of which have been correlated with the risk of SCD, such as QT dispersion<sup>43,44</sup> and principal component analysis (PCA).<sup>45,46</sup> To assist core ECG laboratories in the generation of more consistent, reproducible measurements, GE provided the first commercially available systems to offer high magnification, on-screen review and evaluation of QT/QTc measurements, using a digital editing tool which is supported by computer-assisted QT measurement analysis.<sup>47,48</sup>

### Acute coronary syndromes and decision support

The correlation between superior clinical outcomes and time-to-treatment is well established for acute myocardial infarction (AMI), whether the treatment is via primary percutaneous transluminal coronary angioplasty (PTCA)<sup>49</sup> or a thrombolytic.<sup>50</sup> As compared to thrombolytic therapy, primary PTCA has been shown to lower cost<sup>51</sup> and shorten length-of-hospital-stay.<sup>52</sup> However, expediting patients through the emergency medical system (EMS) to a prepared cath lab takes more planning and is more complex.<sup>53</sup> As more community hospitals adopt this new treatment path,<sup>54,55</sup> a need exists for decision support to ensure that balloon dilation in the cath lab occurs within the ACC benchmark of 90 ( $\pm$ 30) minutes once the patient enters the emergency medical system.<sup>56</sup>

The pre-cath ECG is the most common cited source of delay.<sup>57,58</sup> "There are numerous studies showing that the use of pre-hospital ECGs have been able to reduce door-to-balloon times."<sup>59</sup> By acquiring a pre-hospital ECG, admission and transport delays are reduced.<sup>60,61</sup> Studies also show that the ECG can be transmitted to the hospital before the patient arrives.<sup>62,63,64,65</sup>

Through the use of equipment with GE's Marquette 12SL program, it has been demonstrated that acquiring pre-hospital ECGs is easy to do,<sup>66</sup> significantly cuts total time-to-treatment,<sup>67,68</sup> and has "the potential to significantly increase the diagnostic accuracy in chest pain patients."<sup>69</sup> "Primary angioplasty programs employing pre-hospital 12-lead ECGs are reporting dramatic reductions in Emergency Department (ED) arrivals to reperfusion times. In particular, one hospital was able to cut its average time by more than 50%, from 87 minutes to only 40 minutes, with other case study programs either approaching or breaking through the 90-minute benchmark."<sup>70</sup>

Once the pre-hospital ECG is acquired by the MUSE system, GE's Marquette 12SL Serial Comparison program can identify significant changes in a patient's ECG history, including new left bundle branch block (LBBB),<sup>71</sup> an accepted indication for primary PTCA in the face of chest pain.<sup>72</sup>

The MUSE system can send a wireless electronic message, which includes the critical ECG, to the on-call cardiologist. Current and previous ECGs from the patient can be viewed on the physician's PDA.<sup>73,74</sup> This can be done automatically by the system, based on computerized recognition of an ST-elevated acute myocardial infarction (AMI),<sup>75,76,77,78</sup> or the emergency physician can manually request it to be sent to the cardiologist. Once confirmed by the cardiologist, the cath lab staff can be notified to prepare the cath lab for patient arrival.<sup>79</sup>

As the leading provider of cath lab monitoring equipment, GE offers the ability to continuously measure and trend the 12-lead ECG before, during, and after the cath lab intervention. These systems perform full 12-lead ST-segment analysis, and store all 12 leads for printing at a later time as specified in the latest AHA standard.<sup>80</sup> Continuous 12-lead trending can also be used to detect a failed thrombolysis and possibly prompt the decision as to whether or not rescue angioplasty is warranted.<sup>81</sup>

As GE Healthcare continues its leadership and expertise in providing assistance during difficult clinical care decisions, our understanding of disease at the molecular level is paving the path for breakthroughs that will transform healthcare. Working together, we'll help you to be better able to predict, diagnose, inform and treat.

## References

- <sup>1</sup>Hurst, J.W. "Current Status of Clinical Electrocardiography with Suggestions for the Improvement of the Interpretive Process." *Am J Cardiol.* 92 (Nov 2003): 1072-1079.
- <sup>2</sup>Hurst, J.W. "Current Status of Clinical Electrocardiography with Suggestions for the Improvement of the Interpretive Process." *Am J Cardiol.* 92 (Nov 2003): 1072-1079.
- <sup>3</sup>Kadish AH, et. al. ACC/AHA Clinical Competence Statement on Electrocardiography and Ambulatory Electrocardiography *Circulation.* 2001;104:3169.
- <sup>4</sup>Braunwald, E, et al. "ACC/AHA 2002 Guideline Update for the Management of Patients with Unstable Angina and Non-ST-Segment Elevation Myocardial Infarction – Summary Article: A Report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Committee on the Management of Patients With Unstable Anginal)." *J Am Coll Cardiol.* 40(7)(Oct 2002): 1366-74.
- <sup>5</sup>Tracy CM, et. al., American College of Cardiology/American Heart Association Clinical Competence Statement on Invasive Electrophysiology Studies, Catheter Ablation, and Cardioversion *Circulation.* 2000;102:2309.
- <sup>6</sup>Sandhu, R., and Bahler, R.C. "Prevalence of QRS Prolongation in a Community Hospital Cohort of Patients with Heart Failure and Its Relation to Left Ventricular Systolic Dysfunction." *Am J Cardiol.* 93(2)(Jan 2004): 244-6.
- <sup>7</sup>Gregory Engel, M.D., James G. Beckerman, M.D., Victor F. Froelicher, M.D., Takuya Yamazaki, M.D., Henry A. Chen, M.D., Electrocardiographic Arrhythmia Risk Testing *Curr Probl Cardiol.* July 2004.
- <sup>8</sup>Joseph P. Frolkis, M.D., Ph.D., Claire E. Pothier, M.S., Eugene H. Blackstone, M.D., and Michael S. Lauer, M.D. Frequent Ventricular Ectopy after Exercise as a Predictor of Death *New England Journal of Medicine* Volume 348:781-790 February 27, 2003.
- <sup>9</sup>Daniel B. Mark, Linda Shaw, F.E. Harrell, M.A. Hlatky, K.L. Lee, J.R. Bengton, C.B. McCants, R.M Califf, and D.B. Pryor. Prognostic value of a treadmill exercise score in outpatients with suspected coronary artery disease. *The New England Journal of Medicine*, 325(12):849-853, Sept. 1991.
- <sup>10</sup>ANSI/AAMI (Association for the Advancement of Medical Instrumentation) EC11:1991/(R)2001, Diagnostic electrocardiographic devices, 2ed.
- <sup>11</sup>Kaiser, W., and Findeis, M. "Artifact processing during exercise testing." *J Electrocardiol.* 32 (Suppl)(1999): 212-9.
- <sup>12</sup>Brodnick DE, inventor; GE Medical Systems Information Technologies, assignee. Method of removing signal interference from sampled data and apparatus for effecting the same. *US patent 6 351 664.* 2002 Feb 26.
- <sup>13</sup>Rowlandson, I. "The Marquette 12SL Program." In: *Common Standards for Quantitative Electrocardiography – 5th Progress Report.* J Williams, ed. Leuven, Belgium. 1985.
- <sup>14</sup>Drazen, E., et. al. "Survey of Computer-Assisted Electrocardiography in the United States." *J Electrocardiol.* 21(suppl)(1988): 598-104.
- <sup>15</sup>Aufferheide, T.P., et al. "Test of the Acute Cardiac Ischemia Time-Insensitive Predictive Instrument (ACI-TIPI) for Prehospital Use. *Ann Emerg Med.* 27(2)(Feb 1996): 193-198.
- <sup>16</sup>Hammill, S., et. al. "Establishment of Signal-Averaged Electrocardiographic Criteria with Frank XYZ Leads and Spectral Filter Used Alone and in Combination with Ejection Fraction to Predict Inducible Ventricular Tachycardia in Coronary Artery Disease." *Am J Cardiol.* 70(3)(Aug 1992): 316-320.
- <sup>17</sup>Nearing, B.D., and Verrier, R L "Modified Moving Average Analysis of T-Wave Alternans to Predict Ventricular Fibrillation with High Accuracy." *J Appl Physiol.* 92(2)(Feb 2002): 541-549.
- <sup>18</sup>Schmidt, Georg, Evaluation of electrocardiograms in the field of extrasystoles *U.S. Patent 6,496,722* December 17, 2002.
- <sup>19</sup>PDAs, Bluetooth And Wi-Fi: Beaming Healthcare Into The Future Edited Transcript – Forum for Decisionmakers The U.S. Medicine Institute For Health Studies Washington, D.C. December 2, 2003 U.S. Medicine Institute for Health Studies 2021 L St. NW, Suite 400, Washington DC 20036: with following quote "Electrocardiograms are now wireless in our hospital. We were one of the first to adopt this. Once you roll the cart up to the patient, you actually download from the hospital system the name and a lot of the demographics, and immediately, once you've taken the electrocardiogram, it is up on what's known as the MUSE System, which can be directly visualized by our hospital system, and that can be visualized from home. We can call in and see those tracings and find out whether the patient is having an impending myocardial infarction or not. We can actually send by e-mail that ECG to this Pocket PC."
- <sup>20</sup>Wake Forest University: Customer testimonial regarding HL7 and MUSE quote: Through an industry standard HIS interface from GE Medical Systems, Wake Forest University Baptist Medical Center has enhanced the flow of cardiology information between the Department of Cardiology and other points in the hospital. The results have been greatly improved efficiency, greater accuracy and easier billing."
- <sup>21</sup>Jerome S. Fischer, Lawrence Blonde. Impact of an Electronic Medical Record on Diabetes Practice Workflow *Clinical Diabetes*, April, 1999.
- <sup>22</sup>Fletcher, RD et. al. Computerized medical records in the Department of Veterans Affairs *Cancer* Volume 91, Issue S8 , Pages 1603 – 1606 Published Online: 17 Apr 2001.
- <sup>23</sup>Farrell RM, Young BJ, Effect of Lead Quality on Computerized ECG Interpretation *Computers in Cardiology* 2004;31:173-176.
- <sup>24</sup>Selker, H.P., et al. "Use of the Acute Cardiac Ischemia Time-Insensitive Predictive Instrument (ACI-TIPI) to Assist with Triage of Patients with Chest Pain or Other Symptoms Suggestive of Acute Cardiac Ischemia. A Multicenter, Controlled Clinical Trial." *Ann Intern Med.* 129(11)(Dec 1998): 845-55.
- <sup>25</sup>Leibrandt, P.N., et. al. "Validation of Cardiologist's Decisions to Initiate Reperfusion Therapy for Acute Myocardial Infarction with Electrocardiograms Viewed on Liquid Crystal Displays of Cellular Telephones." *Am Heart J.* 140(5)(Nov 2000): 747-52.
- <sup>26</sup>*US Patent 6,801,802 B2* System and method for selecting physiological data from a plurality of physiological data sources David A. Sitzman, Robert M. Farrell, Filed on Jun. 29, 2001, as Appl. No. 9/896,076. Prior Publication US 2003/0009106 A1, Jan. 09, 2003.
- <sup>27</sup>Shapiro, R.M, et al. "Reduction in Alarm Frequency with a Fusion Algorithm for Processing Monitor Signals." *Proceedings of the American Thoracic Society International Conference*, 2002.
- <sup>28</sup>GANG, YI et al. Preoperative Electrocardiographic Risk Assessment of Atrial Fibrillation After Coronary Artery Bypass Grafting. *Journal of Cardiovascular Electrophysiology* Volume 15 Issue 12 Page 1379 – December 2004.
- <sup>29</sup>Dhala, A., et al. "Signal Averaged P-Wave Analysis of Normal Controls and Patients with Paroxysmal Atrial Fibrillation: A Study in Gender Differences, Age Dependence, and Reproducibility." *Clin Cardiol.* 25(11)(Nov 2002): 525-31.
- <sup>30</sup>Farrell, R.M, Xue, J.Q., and Young, B.J. "Enhanced Rhythm Analysis for Resting ECG Using Spectral and Time-Domain Techniques." *Computers in Cardiology* 2003. Los Alamitos, CA: IEEE Computer Society Press. 30(Sept 2003): 733-736.
- <sup>31</sup>Reddy, B. R., Elko, P., and Swiryn, S. "A New Arrhythmia Analysis Program, MAC-RHYTHM, for Resting ECG. In: *Proceedings of the XXIII International Congress on Electrocardiology* 1996. Lieberman, J., ed. Cleveland, OH July-Aug 1996: 471-480.
- <sup>32</sup>Reddy, B.R., et al. "Prospective Evaluation of a Microprocessor-Assisted Cardiac Rhythm Algorithm: Results from One Clinical Center." *J Electrocardiol.* 30(Suppl)(1998): 28-33.
- <sup>33</sup>Young, B.J., Brodnick, D., and Rowlandson, I. "A Comparative Study of a Hidden Markov Model Detector for Atrial Fibrillation." *Proceedings of the 9th IEEE Workshop on Neural Networks for Signal Processing.* Madison, WI: IEEE Signal Processing Society. (19 Aug 1999).
- <sup>34</sup>Froelicher, V., et. al. "Prognostic Value of Computer Electrocardiography in Veteran Outpatients" *Federal Practitioner* (Mar 2004): 11-20.
- <sup>35</sup>Hammill, S., et. al. "Establishment of Signal-Averaged Electrocardiographic Criteria with Frank XYZ Leads and Spectral Filter Used Alone and in Combination with Ejection Fraction to Predict Inducible Ventricular Tachycardia in Coronary Artery Disease." *Am J Cardiol.* 70(3)(Aug 1992): 316-320.
- <sup>36</sup>Willem J. Kop, PhD; David S. Krantz, PhD; Bruce D. Nearing, PhD; John S. Gottdiener, M.D.; John F. Quigley, PhD; Mark O'Callahan, BS; Albert A. DelNegro, M.D.; Ted D. Friehling, M.D.; Pamela Karasik, M.D.; Sonia Suchday, PhD; Joseph Levine, M.D.; Richard L. Verrier, PhD. Effects of Acute Mental Stress and Exercise on T-Wave Alternans in Patients With Implantable Cardioverter Defibrillators and Controls *Circulation.* 2004 Apr 20;109(15):1864-9.
- <sup>37</sup>Sands, K.E.F., et al. "Power Spectrum Analysis of Heart Rate Variability in Human Cardiac Transplant Recipients." *Circulation.* 79(1)(Jan 1989): 76-82.
- <sup>38</sup>Kleiger RE, et al. Decreased heart rate variability and its association with increased mortality after acute myocardial infarction. *Am J Cardiol.* 1987 Feb 1;59(4):256-62.
- <sup>39</sup>Tsuji et al. Impact of Reduced Heart Rate Variability on Risk for Cardiac Events: The Framingham Heart Study *Circulation.* 1996;94:2850-2855.
- <sup>40</sup>Schmidt G, Malik M, Barthel P, Schneider R, Ulm K, Rolnitzky L, Camm AJ, Bigger JT Jr, Schomig A Heart-rate turbulence after ventricular premature beats as a predictor of mortality after acute myocardial infarction. *Lancet.* 1999 Apr 24;353(9162):1377-9.
- <sup>41</sup>Ghuran A. et al. Heart rate turbulence-based predictors of fatal and nonfatal cardiac arrest (The Autonomic Tone and Reflexes After Myocardial Infarction substudy). *Am J Cardiol.* 2002 Jan 15;89(2):184-90.
- <sup>42</sup>Barthel P. et al., Risk Stratification After Acute Myocardial Infarction by Heart Rate Turbulence *Circulation.* 2003;108:1221.
- <sup>43</sup>Okin, PM, et al. "Assessment of QT Interval and QT Dispersion for Prediction of All-Cause and Cardiovascular Mortality in American Indians: The Strong Heart Study." *Circulation.* 101(2000): 61-66.
- <sup>44</sup>Xue, Q., and Reddy, S. "New Algorithms for QT Dispersion Analysis." *Comput Cardiol.* (1996): 293-296.
- <sup>45</sup>Okin, PM, et al. "Principal Component Analysis of the T-Wave and Prediction of Cardiovascular Mortality in American Indians: The Strong Heart Study." *Circulation.* 105(2002): 714-719.
- <sup>46</sup>Milos K, et al. Principal Component Analysis of the T-Wave in Patients with Chest Pain and Conduction Disturbances. *Pacing and Clinical Electrophysiology* Vol. 27 Issue 10 Page 1378, October 2004.

- <sup>47</sup>Xue et. al., Clinical research workstation *US Patent 6,463,320* October 8, 2002.
- <sup>48</sup>Xue, Q., and Reddy, S. "Computerized QT Analysis Algorithms." *J Electrocardiol.* 30(Suppl)(1997): 181-186.
- <sup>49</sup>Patel, S., et. al. "Adverse Outcomes Accompanying Primary Angioplasty (PTCA) for Acute Myocardial Infarction (AMI) – Dangers of Delay." *J Am Coll Cardiol.* 27(suppl A)(1996): 62A.
- <sup>50</sup>Hermens, W.T., et.al. "Effect of Thrombolytic Treatment Delay on Myocardial Infarct Size. *Lancet.* 340(8830)(Nov 1992): 1297.
- <sup>51</sup>Mark, D.B., et. al. "Costs of Direct Angioplasty Versus Thrombolysis for Acute Myocardial Infarction: Results from the GUSTO-II Randomized Trial." *Circulation.* 94 (1996): 168.
- <sup>52</sup>Brodie, B.R. "Early Hospital Discharge After Acute Myocardial Infarction." *The Journal of Invasive Cardiology.* 7(suppl F)(1995): 22F-28F.
- <sup>53</sup>Kent DM, Lau J, Selker HP. Balancing the benefits of primary angioplasty against the benefits of thrombolytic therapy for acute myocardial infarction: the importance of timing. *Eff Clin Pract.* 2001 Sep-Oct;4(5):214-20.
- <sup>54</sup>Brush, J.E., et. al. "Retrospective Comparison of a Strategy of Primary Coronary Angioplasty Versus Intravenous Thrombolytic Therapy for Acute Myocardial Infarction in a Community Hospital Without Surgical Backup." *J Invasive Cardiol.* 8(2)(Mar 1996): 91-98.
- <sup>55</sup>McNamara, N.S., et. al. "Can Community Hospitals Provide Effective Primary PTCA Coverage at All Hours?" *J Am Coll Cardiol.* 29(suppl 1:91A)(1997): 1.91A.
- <sup>56</sup>Ryan TJ, Antman EM, Brooks NH, Califf RM, Hillis LD, Hiratzka LF, Rapaport E, Riegel B, Russell RO, Smith EE III, Weaver WD. ACC/AHA guidelines for the management of patients with acute myocardial infarction: 1999 update: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Committee on Management of Acute Myocardial Infarction).
- <sup>57</sup>Sharkey, S.W., et. al. "An Analysis of Time Delays Preceding Thrombolysis for Myocardial Infarction." *JAMA.* 262(22)(Dec 1989): 3171-4.
- <sup>58</sup>Kereiakes D., Weaver W., et. al., Time delays in the diagnosis and treatment of acute myocardial infarction: a tale of eight cities. *Amer Heart J* 1990;120(4):773-779.
- <sup>59</sup>Guerra, DR et. al. Door-to-Balloon Delays with PCI in Acute Myocardial Infarction *Current Treatment Options in Cardiovascular Medicine* Volume 6, Issue 1: February 2004:69-77.
- <sup>60</sup>Aufderheide, T.P., et. al. "Milwaukee Prehospital Chest Pain Project: Phase I: Feasibility and Accuracy of Prehospital Thrombolytic Candidate Selection." *Am J Cardiol.* 69(12)(Apr 1992): 991-996.
- <sup>61</sup>Karagounis L., et. al., Impact of field transmitted electrocardiography on time to in-hospital thrombolytic therapy in acute myocardial infarction. *Circulation* 1989;80:II-353.
- <sup>62</sup>Garza MA Cellular technology. Revolutionizing EMS communications. *JEMS.* 1990 May;15(5):46-7, 50-3, 55-6.
- <sup>63</sup>Sherrid M, Greenberg H, et. al. A pilot study of paramedic-administered, prehospital thrombolysis for acute myocardial infarction. *Clin Cardiol.* 1990 Jun;13(6):421-4.
- <sup>64</sup>Weaver, W.D., et. al. "Myocardial Infarction Triage and Intervention Project: Phase I: Patient Characteristics and Feasibility of Prehospital Initiation of Thrombolytic Therapy." *J Am Coll Cardiol.* 15(5)(Apr 1990): 925-931.
- <sup>65</sup>Gibler, W.B., et al. "Prehospital Diagnosis and Treatment of Acute Myocardial Infarction: A North-South Perspective. The Cincinnati Heart Project and the Nashville Prehospital TPA Trial." *Am Heart J.* 121(1 Pt 1)(Jan 1991): 1-11.
- <sup>66</sup>Aufderheide, T.P., et al. "A Prospective Evaluation of Prehospital 12-lead ECG Application in Chest Pain Patients." *J Electrocardiol.* 24 Suppl(1992): 8-13.
- <sup>67</sup>Kereiakes, D.J., et al. "Relative Importance of Emergency Medical System Transport and the Prehospital Electrocardiogram on Reducing Hospital Time Delay to Therapy for Acute Myocardial Infarction: A Preliminary Report from the Cincinnati Heart Project." *Am Heart J.* 23(4 Pt 1)(Apr 1992): 835-840.
- <sup>68</sup>Foster, D.B., et al. "Prehospital Recognition of AMI Using Independent Nurse/Paramedic 12-lead ECG Evaluation: Impact on In-Hospital Times to Thrombolysis in a Rural Community Hospital." *Am J Emerg Med.* 12(1)(Jan 1994): 25-31.
- <sup>69</sup>Aufderheide, T.P., et. al. "The Diagnostic Impact of Prehospital 12-lead Electrocardiography." *Annals of Emergency Medicine.* 19(11)(1990): 1280-1287.
- <sup>70</sup>The Advisory Board Company, Myocardial Infarction: Toward a Higher Standard of Care 1998:267.
- <sup>71</sup>Rowlandson, G. In Inventor *US Pat App: 20020087055* System and method for detecting new left branch bundle block for accelerating treatment of acute myocardial infarction July 4, 2002.
- <sup>72</sup>Edhouse J et. al. ABC of clinical electrocardiography Acute myocardial infarction-Part II *British Medical Journal* Vol. 324 20 Apr 2002.
- <sup>73</sup>Leibrandt, P.N., et. al. "Validation of Cardiologist's Decisions to Initiate Reperfusion Therapy for Acute Myocardial Infarction with Electrocardiograms Viewed on Liquid Crystal Displays of Cellular Telephones." *Am Heart J.* 140(5)(Nov 2000): 747-52.
- <sup>74</sup>Pettis, K.S., et. al. "Evaluation of the Efficacy of Hand-Held Computer Screens for Cardiologist's Interpretations of 12-Lead Electrocardiograms." *Am Heart J.* 138(4 pt 1)(Oct 1999): 765-70.
- <sup>75</sup>O'Rourke, M.F., Cook, A., Carroll, G., et al. Accuracy of a portable interpretive ECG machine in diagnosis of acute evolving myocardial infarction. *Aust N Z J Med* 1992; 22: 9-13.
- <sup>76</sup>Kudenchuk, P.J., et al. "Accuracy of Computer Interpreted Electrocardiography in Selecting Patients for Thrombolytic Therapy. MITI Project Investigators." *J Am Coll Cardiol.* 17(7)(Jun 1991): 1486-1491.
- <sup>77</sup>Xue, J., et. al. "A New Method to Incorporate Age and Gender Into the Criteria for the Detection of Acute Inferior Myocardial Infarction." *J Electrocardiol.* 34(4)(Part 2)(Oct 2001): 229-234.
- <sup>78</sup>Massel, D., Dawdy, J.A., Melendez, L.J. Strict reliance on a computer algorithm or measurable ST segment criteria may lead to errors in thrombolytic therapy eligibility. *Am Heart J.* 2000 Aug;140(2):221-6.
- <sup>79</sup>Rowlandson, Ian G., *US Patent 2002/0087355-A1*, "Automated Scheduling of Emergency Procedure Based On Identification of High-Risk Patient", Jul. 4, 2002.
- <sup>80</sup>Barbara J. Drew, RN, PhD, Chair AHA Scientific Statement Practice Standards for Electrocardiographic Monitoring in Hospital Settings An American Heart Association Scientific Statement From the Councils on Cardiovascular Nursing, Clinical Cardiology, and Cardiovascular Disease in the Young: Endorsed by the International Society of Computerized Electrocardiology and the American Association of Critical-Care Nurses *Circulation.* 2004;110:2721-2746.
- <sup>81</sup>Smith et. al., ACC/AHA Guidelines for Percutaneous Coronary Intervention (Revision of the 1993 PTCA Guidelines) *JACC VOL. 37, NO. 8, JUNE 2001:2339I-IXVI.*

©2005 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.

GE, GE Monogram, MUSE®, Marquette®, MARS® Unity Network®, EKPro™, Hookup Advisor™ and 12SL™ are trademarks of General Electric Company.

ZOLL® is a registered trademark of Zoll Medical Corporation.  
LIFEPAK® is a registered trademark of Medtronic.

GE Medical Systems Information Technologies,  
a General Electric company, going to market as GE Healthcare.

For more than 100 years, scientists and industry leaders have relied on General Electric for technology, services and productivity solutions. So no matter what challenges your healthcare system faces – you can always count on GE to help deliver the highest quality services and support.

For details, please contact your GE Healthcare representative today.

GE Healthcare  
3000 North Grandview  
Waukesha, WI 53188  
U.S.A.

[www.gehealthcare.com](http://www.gehealthcare.com)



imagination at work