			RSP6-16-RS	RIC5-9-RS RIC5-9W-RS	RAB2-5-RS	RAB4-8-RS	RNA5-9-RS	AB2-7-RS	SP10-16-RS	
	<b>GE Ultrasound Contact Gel</b>	Gel	GE Accessories	OK	OK	OK	OK	OK	OK	—
	<b>Aquasonics 100</b>	Gel	Parker Laboratory Inc.	OK	OK	OK	OK	OK	OK	—
	<b>Haiyin</b>	Gel	Wuxi Huali Industrial Company	⊗	⊗	⊗	⊗	⊗	OK	—
	<b>PROJELLY</b>	Gel	JEX Inc.	⊗	⊗	⊗	⊗	⊗	OK	—
	<b>Sonogel</b>	Gel	Sonogel Vertriebs GmbH	OK	OK	OK	OK	OK	⊗	—
	<b>Scan</b>	Gel	Parker Laboratory Inc.	OK	OK	OK	OK	OK	OK	—
	<b>Incidin Foam</b>	20%IPA, 10%Ethyl Alcohol	Ecolab	⊗	OK	⊗	⊗	⊗	⊗	< 5 minutes
	<b>T-Spray</b>	Quaternary Ammoniums	Pharmaceutical Innovations	OK	OK	OK	OK	OK	OK	< 20 minutes
	<b>T-Spray II</b>	Quaternary Ammoniums	Pharmaceutical Innovations	OK	OK	OK	OK	OK	OK	< 20 minutes
	<b>Transeptic Spray</b>	70% IPA	Parker Laboratory Inc.	⊗	OK	OK	OK	⊗	OK	< 15 minutes
	<b>Asepti-Wipes</b>	8%IPA, <1%Quat	Ecolab Inc.	⊗	OK	OK	OK	OK	⊗	< 15 minutes
	<b>Septiwipes</b>	Quaternary Ammoniums	Dr.Schumacher GmbH	OK	OK	OK	OK	OK	⊗	< 15 minutes
	<b>Cleanisept</b>	Quaternary Ammoniums	Dr.Schumacher GmbH	OK	OK	OK	OK	OK	⊗	< 15 minutes
	<b>MetriWipes</b>	8%IPA, <1%Quat	Metrex	⊗	OK	OK	OK	OK	⊗	< 15 minutes
	<b>NUCLEAN II</b>	Isopropanol / Ethanol	Alkapharm	⊗	OK	⊗	⊗	⊗	⊗	< 15 minutes
	<b>Sani-Cloth Active</b>	Quaternary Ammoniums	PDI	OK	⊗	⊗	⊗	⊗	⊗	< 15 minutes
	<b>Sani-Cloth Germicidal</b>	8%IPA, <1%Quat	PDI	⊗	OK	OK	OK	OK	⊗	< 15 minutes
	<b>Sani-Cloth HB</b>	Quaternary Ammoniums	PDI	OK	OK	OK	OK	OK	OK	< 30 minutes
	<b>PeraSafe</b>	Peracetic Acid	Antec International	OK	OK	OK	OK	OK	⊗	< 15 minutes
	<b>Alkazyme</b>	Proteolytic Enzymes	Alkapharm	⊗	⊗	⊗	⊗	⊗	⊗	< 10 minutes
	<b>Cidex</b>	Glutaraldehyde	Johnson & Johnson	OK	OK	OK	OK	OK	OK	< 50 minutes
	<b>Cidex OPA</b>	Ortho-phthalaldehyde	Johnson & Johnson	OK	OK	OK	OK	OK	OK	< 15 minutes
	<b>Cidex Plus</b>	Glutaraldehyde	Johnson & Johnson	⊗	OK	OK	OK	OK	⊗	< 25 minutes
	<b>Gigasept</b>	Succindialdehyde, formaldehyde	Schülke & Mayr	OK	OK	⊗	⊗	OK	OK	< 30 minutes
	<b>Gigasept AF</b>	Quaternary Ammoniums	Schülke & Mayr	⊗	OK	OK	OK	⊗	OK	< 30minutes
	<b>Gigasept FF</b>	Succindialdehyde	Schülke & Mayr	OK	OK	OK	OK	OK	⊗	< 30minutes
	<b>Klenzyme</b>	Proteolytic Enzymes	STERIS	⊗	OK	OK	OK	⊗	⊗	< 10 minutes
	<b>Metricide 14</b>	Glutaraldehyde	Metrex	OK	OK	OK	OK	OK	OK	< 50 minutes
	<b>Metricide 28</b>	Glutaraldehyde	Metrex	⊗	⊗	⊗	⊗	⊗	⊗	< 50 minutes
	<b>Metricide 30</b>	Glutaraldehyde	Metrex	⊗	⊗	⊗	⊗	⊗	⊗	< 25 minutes
	<b>Sekusept Plus</b>	Nonionic surfactant	Henkel-Ecolab	⊗	OK	⊗	⊗	⊗	⊗	< 30 minutes
	<b>Ster-Bac Blu</b>	Quaternary Ammoniums	Ecolab Inc.	OK	OK	OK	OK	OK	OK	< 50 minutes
	<b>Wavicide-01</b>	Glutaraldehyde	MCC	OK	OK	OK	OK	OK	OK	< 50 minutes

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