GE Healthcare

**Innova Vision, Innova EPVision, Innova
TrackVision
Conformance Statement for DICOM**

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WARNING**X-RAY EQUIPMENT IS DANGEROUS TO BOTH PATIENT AND OPERATOR UNLESS MEASURES OF PROTECTION ARE STRICTLY OBSERVED**

Though this equipment is built to the highest standards of electrical and mechanical safety, the useful x-ray beam becomes a source of danger in the hands of the unauthorized or unqualified operator.

Excessive exposure to x-radiation causes damage to human tissue.

Therefore, adequate precautions must be taken to prevent unauthorized or unqualified persons from operating this equipment or exposing themselves or others to its radiation.

Before operation, persons qualified and authorized to operate this equipment should be familiar with the Recommendations of the International Commission on Radiological Protection, contained in Annals Number 26 of the ICRP, and with applicable national standards.

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**WARNING
(EN)**

This service manual is available in English only.

- If a customer's service provider requires a language other than English, it is the customer's responsibility to provide translation services.
- Do not attempt to service the equipment unless this service manual has been consulted and is understood.
- Failure to heed this warning may result in injury to the service provider, operator or patient from electric shock, mechanical or other hazards.

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This page left intentionally blank.

Revision History

Part / Rev	Date	Reason for change	Pages
5392762-1-8EN rev 1	May 5, 2010	Initial release of direction 5392762-1-8EN	70
5392762-2-8EN rev 1	July 23, 2010	Initial release of direction 5392762-2-8EN	70
5392762-3-8EN rev 1	September 13, 2010	Initial release of direction 5392762-3-8EN	70

This page left intentionally blank.

Table of Contents

Chapter 1 CONFORMANCE STATEMENT OVERVIEW.....17

 1 Conformance Statement Overview.....17

Chapter 2 INTRODUCTION.....19

 1 Introduction.....19

 1.1 Overview.....19

 1.2 Overall Dicom Conformance Statement Document Structure.....20

 1.3 Intended Audience.....22

 1.4 Scope and Field Application.....22

 1.5 Important Remarks.....22

 1.6 References.....23

 1.7 Definitions.....23

 1.8 Symbols and Abbreviations.....26

Chapter 3 NETWORK CONFORMANCE STATEMENT.....29

 1 Network Conformance Statement.....29

 1.1 Introduction.....29

 1.2 Implementation Model.....30

 1.2.1 Application Data Flow Diagram.....30

 1.2.2 Presentation Context Table.....30

 1.2.3 Real-World Activities.....30

 1.2.4 SOP Instance UID.....30

 1.2.5 Support of Extended Character Sets.....31

**Chapter 4 INNOVA VISION APPLICATIONS DERIVED X-RAY ANGIOGRAPHY (XA)
 INFORMATION OBJECT IMPLEMENTATION.....33**

 1 Introduction.....33

 2 Innova Vision Applications Derived X-Ray Image IOD Implementation.....34

 3 Innova Vision Applications X-Ray Image IOD Entity-Relationship Model.....35

 4 Entities Description.....37

 5 Mapping of DICOM Entities.....38

6 IOD Module Table.....	39
7 Information Module Definitions.....	41
7.1 Common Information Module Definitions.....	41
7.2 Patient Entity.....	41
7.3 Study Entity.....	41
7.3.1 General Study Module.....	41
7.3.2 Patient Study Module.....	42
7.4 Series Entity.....	42
7.5 Equipment Entity.....	42
7.6 Image Entity.....	43
7.6.1 General Image Module.....	43
7.6.2 Image Pixel Module.....	44
7.6.3 Cine Module.....	44
7.6.4 Multi Frame Module.....	44
7.6.5 Frame Pointers Module.....	44
7.6.6 X-Ray Image Module.....	45
7.6.7 X-Ray Acquisition Module.....	45
7.6.8 X-Ray Table Module.....	45
7.6.9 XA Positioner Module.....	45
7.6.10 VOI LUT Module.....	46
7.6.11 SOP Common Module.....	46
7.7 Standard Extended and Private Data Attributes.....	46
7.7.1 Standard Attributes.....	46
7.7.2 Private Group DLX_SERIE_01.....	47
7.7.3 Private Group GEMS_DL_IMG_01.....	47
7.7.4 Private Group GEMS_XR3DCAL_01.....	51
7.7.5 Private Group GEMS_3D_INTVL_01.....	51
Chapter 5 INNOVA VISION APPLICATIONS SECONDARY INFORMATION OBJECT IMPLEMENTATION.....	55
1 Introduction.....	55
2 Innova Vision Applications Secondary Image IOD Implementation.....	56

3	Innova Vision Applications Secondary Image IOD Entity - Relationship Model.....	57
4	Entities Description.....	59
5	Mapping of DICOM Entities.....	60
6	IOD Module Table.....	61
7	Information Module Definitions.....	62
7.1	Common Information Module Definitions.....	62
7.2	Patient Entity.....	62
7.3	Study Entity.....	62
7.3.1	General Study Module.....	62
7.3.2	Patient Study Module.....	63
7.4	Series Entity.....	63
7.5	Equipment Entity.....	64
7.5.1	General Equipment Module.....	64
7.5.2	SC Equipment Module.....	64
7.6	Image Entity.....	64
7.6.1	General Image Module.....	64
7.6.2	Image Pixel Module.....	65
7.6.3	SC Image Module.....	65
7.6.4	SOP Common Module.....	65
7.7	Standard Extended and Private Data Attributes.....	66
7.7.1	Standard Attributes.....	66
7.7.2	Private Group DLX_SERIE_01.....	66
7.7.3	Private Group GEMS_DL_IMG_01.....	66

This page left intentionally blank.

Chapter 1 Conformance Statement Overview

1 Conformance Statement Overview

Innova Vision, Innova EPVision, Innova TrackVision are post processing applications that are installed on the same hardware platform as the base application, **Advantage Workstation**.

Innova Vision, Innova EPVision, Innova TrackVision do not have an intrinsic DICOM Network feature. They do not directly invoke the DICOM Server AE and they use the network services supported by the base application: Advantage Workstation.

Refer to the respective Conformance Statement - **Advantage Workstation** where Innova Vision, Innova EPVision, Innova TrackVision applications are running to understand the network services supported.

Innova Vision applications can be installed only on supported platform defined by Volume Viewer VS4 application (AW4.5 and above AW versions).

The application uses the following DICOM objects created in Innova system.

Table 1-1: (Input)

Sop Class Name	Sop Class UID
X-Ray Image Storage	1.2.840.10008.5.1.4.1.1.12.1

The application additionally uses the 3D model (3D volume can be saved in DICOM format) or 3D SaveState (a Volume Viewer object) that is created in AW Volume Viewer.

Table 1-2: (Input)

Sop Class Name	Sop Class UID
GEMS Private DICOM 3D Model	1.2.840.113619.4.26
Secondary Capture Image Storage (3D Save State)	1.2.840.10008.5.1.4.1.1.7

The application generates the derived instances of the following DICOM objects:

Table 1-3: (Output)

Sop Class Name	Sop Class UID
X-Ray Image Storage	1.2.840.10008.5.1.4.1.1.12.1
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7

This page left intentionally blank.

Chapter 2 Introduction

1 Introduction

1.1 Overview

This DICOM Conformance Statement is divided into Chapters as described below:

- **Chapter 1 (Conformance Statement Overview)**, which describes the purpose of this Conformance Statement.
- **Chapter 2 (Introduction)**, which describes the overall structure, intent, and references for this Conformance Statement.
- **Chapter 3 (Network Conformance Statement)**, which specifies the GEMS equipment compliance to the DICOM requirements for the implementation of Networking features.
- **Chapter 4 (Innova Vision Applications Derived X-Ray Angiography (XA) Information Object Implementation)**, which specifies the GEMS equipment compliance to DICOM requirements for the implementation of a Derived XA Image Information Object.
- **Chapter 5 (Innova Vision Applications Secondary Information Object Implementation)**, which specifies the GEMS equipment compliance to DICOM requirements for the implementation of a Secondary Capture Information Object.

The application uses the following DICOM objects created in Innova system. (For a complete description of the IOD, refer to the Innova DICOM conformance statement, direction 5394268-1-8EN.)

Table 2-1:

Sop Class Name	Sop Class UID
X-Ray Image Storage	1.2.840.10008.5.1.4.1.1.12.1

The application additionally uses the 3D model (3D volume can be saved in DICOM format) or 3D SaveState (a Volume Viewer object) that is created in AW Volume Viewer. (For a complete description of the IODs, refer to the AW Volume Viewer DICOM conformance statement, direction 5344139-100.)

Table 2-2:

Sop Class Name	Sop Class UID
GEMS Private DICOM 3D Model	1.2.840.113619.4.26
Secondary Capture Image Storage (3D Save State)	1.2.840.10008.5.1.4.1.1.7

The application generates the derived instances of the following DICOM objects:

Table 2-3:


Sop Class Name	Sop Class UID
X-Ray Image Storage	1.2.840.10008.5.1.4.1.1.12.1
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7

Innova Vision, Innova EPVision, Innova TrackVision is referred as Innova Vision applications in this document.

1.2 Overall Dicom Conformance Statement Document Structure

The Documentation Structure of the GEMS DICOM Conformance Statements is shown in the Illustration below.

Illustration 2-1: GEMS DICOM Conformance Statements

This DICOM Conformance Statement documents the DICOM Conformance Statement and Technical Specification required to interoperate with the GEMS network interface.

The GEMS Conformance Statement, contained in this document, also specifies the Lower Layer communications which it supports (e.g., TCP/IP). However, the Technical Specifications are defined in the DICOM Part 8 standard.

For more information regarding DICOM, copies of the Standard may be obtained on the Internet at <http://medical.nema.org>. Comments on the Standard may be addressed to:

DICOM Secretariat
 NEMA
 1300 N. 17th Street, Suite 1847

Rosslyn, VA 22209

USA

Phone: +1.703.841.3200

1.3 Intended Audience

The reader of this document is concerned with software design and/or system integration issues. It is assumed that the reader of this document is familiar with the DICOM Standard and with the terminology and concepts which are used in that Standard.

1.4 Scope and Field Application

It is the intent of this document to provide an unambiguous specification for GEMS implementations. This specification, called a Conformance Statement, includes a DICOM Conformance Statement and is necessary to ensure proper processing and interpretation of GEMS medical data exchanged using DICOM. The GEMS Conformance Statements are available to the public.

The reader of this DICOM Conformance Statement should be aware that different GEMS devices are capable of using different Information Object Definitions. For example, a GEMS CT Scanner may send images using the CT Information Object, MR Information Object, Secondary Capture Object, etc.

Included in this DICOM Conformance Statement are the Module Definitions which define all data elements used by this GEMS implementation. If the user encounters unspecified private data elements while parsing a GEMS Data Set, the user is well advised to ignore those data elements (per the DICOM standard). Unspecified private data element information is subject to change without notice. If, however, the device is acting as a "full fidelity storage device", it should retain and re-transmit all of the private data elements which are sent by GEMS devices.

1.5 Important Remarks

The use of these DICOM Conformance Statements, in conjunction with the DICOM Standards, is intended to facilitate communication with GE imaging equipment. However, **by itself, it is not sufficient to ensure that inter-operation will be successful**. The **user (or user's agent)** needs to proceed with caution and address at least four issues:

- **Integration** – The integration of any device into an overall system of interconnected devices goes beyond the scope of standards (DICOM v3.0), and of this introduction and associated DICOM Conformance Statements when interoperability with non-GE equipment is desired. The responsibility to analyze the applications requirements and to design a solution that integrates GE imaging equipment with non-GE systems is the **user's** responsibility and should not be underestimated. The **user** is strongly advised to ensure that such an integration analysis is correctly performed.
- **Validation** – Testing the complete range of possible interactions between any GE device and non-GE devices, before the connection is declared operational, should not be overlooked. Therefore, the **user** should ensure that any non-GE provider accepts full responsibility for all validation required for their connection with GE devices. This includes the accuracy of the

image data once it has crossed the interface between the GE imaging equipment and the non-GE device and the stability of the image data for the intended applications.

Such a validation is required before any clinical use (diagnosis and/or treatment) is performed. It applies when images acquired on GE imaging equipment are processed/displayed on a non-GE device, as well as when images acquired on non-GE equipment is processed/displayed on a GE console or workstation.

- **Future Evolution** – GE understands that the DICOM Standard will evolve to meet the user's growing requirements. GE is actively involved in the development of the DICOM Standard. DICOM will incorporate new features and technologies and GE may follow the evolution of the Standard. The GEMS protocol is based on DICOM as specified in each DICOM Conformance Statement. Evolution of the Standard may require changes to devices which have implemented DICOM. **In addition, GE reserves the right to discontinue or make changes to the support of communications features (on its products) described by these DICOM Conformance Statements.** The **user** should ensure that any non-GE provider, which connects with GE devices, also plans for the future evolution of the DICOM Standard. Failure to do so will likely result in the loss of function and/or connectivity as the DICOM Standard changes and GE Products are enhanced to support these changes.
- **Interaction** – It is the sole responsibility of the **non-GE provider** to ensure that communication with the interfaced equipment does not cause degradation of GE imaging equipment performance and/or function.

1.6 References

NEMA PS3:

Digital Imaging and Communications in Medicine (DICOM) Standard, available free at <http://medical.nema.org/>.

1.7 Definitions

Informal definitions are provided for the following terms used in this Conformance Statement. The DICOM Standard is the authoritative source for formal definitions of these terms.

Abstract Syntax

The information agreed to be exchanged between applications, generally equivalent to a Service/Object Pair (SOP) Class. Examples : Verification SOP Class, Modality Worklist Information Model Find SOP Class, Computed Radiography Image Storage SOP Class.

Application Entity (AE)

An end point of a DICOM information exchange, including the DICOM network or media interface software; i.e., the software that sends or receives DICOM information objects or messages. A single device may have multiple Application Entities.

Application Entity Title

The externally known name of an *Application Entity*, used to identify a DICOM application to other DICOM applications on the network.

Application Context

The specification of the type of communication used between *Application Entities*. Example: DICOM network protocol.

Association

A network communication channel set up between *Application Entities*.

Attribute

A unit of information in an object definition; a data element identified by a tag. The information may be a complex data structure (Sequence), itself composed of lower level data elements. Examples: Patient ID (0010,0020), Accession Number (0008,0050), Photometric Interpretation (0028,0004), Procedure Code Sequence (0008,1032).

Information Object Definition (IOD)

The specified set of *Attributes* that comprise a type of data object; does not represent a specific instance of the data object, but rather a class of similar data objects that have the same properties. The Attributes may be specified as Mandatory (Type 1), Required but possibly unknown (Type 2), or Optional (Type 3), and there may be conditions associated with the use of an Attribute (Types 1C and 2C). Examples: MR Image IOD, CT Image IOD, Print Job IOD.

Joint Photographic Experts Group (JPEG)

A set of standardized image compression techniques, available for use by DICOM applications.

Media Application Profile

The specification of DICOM information objects and encoding exchanged on removable media (e.g., CDs).

Module

A set of Attributes within an *Information Object Definition* that are logically related to each other. Example: Patient Module includes Patient Name, Patient ID, Patient Birth Date, and Patient Sex.

Negotiation

First phase of Association establishment that allows *Application Entities* to agree on the types of data to be exchanged and how that data will be encoded.

Presentation Context

The set of DICOM network services used over an Association, as negotiated between *Application Entities*; includes *Abstract Syntaxes* and *Transfer Syntaxes*.

Protocol Data Unit (PDU)

A packet (piece) of a DICOM message sent across the network. Devices must specify the maximum size packet they can receive for DICOM messages.

Security Profile

A set of mechanisms, such as encryption, user authentication, or digital signatures, used by an *Application Entity* to ensure confidentiality, integrity, and/or availability of exchanged DICOM data.

Service Class Provider (SCP)

Role of an *Application Entity* that provides a DICOM network service; typically, a server that performs operations requested by another *Application Entity* (*Service Class User*). Examples: Picture Archiving and Communication System (image storage SCP, and image query/retrieve SCP), Radiology Information System (modality worklist SCP).

Service Class User (SCU)

Role of an *Application Entity* that uses a DICOM network service; typically, a client. Examples: imaging modality (image storage SCU, and modality worklist SCU), imaging workstation (image query/retrieve SCU).

Service/Object Pair (SOP) Class

The specification of the network or media transfer (service) of a particular type of data (object); the fundamental unit of DICOM interoperability specification. Examples: Ultrasound Image Storage Service, Basic Grayscale Print Management.

Service/Object Pair (SOP) Instance

An information object; a specific occurrence of information exchanged in a *SOP Class*. Examples: a specific x-ray image.

Tag

A 32-bit identifier for a data element, represented as a pair of four digit hexadecimal numbers, the “group” and the “element”. If the “group” number is odd, the tag is for a private (manufacturer-specific) data element. Examples: (0010,0020) [Patient ID], (07FE,0010) [Pixel Data], (0019,0210) [private data element].

Transfer Syntax

The encoding used for exchange of DICOM information objects and messages. Examples: *JPEG* compressed (images), little endian explicit value representation.

Unique Identifier (UID)

A globally unique “dotted decimal” string that identifies a specific object or a class of objects; an ISO-8824 Object Identifier. Examples: Study Instance UID, SOP Class UID, SOP Instance UID.

Value Representation (VR)

The format type of an individual DICOM data element, such as text, an integer, a person’s name, or a code. DICOM information objects can be transmitted with either explicit identification of the type of each data element (Explicit VR), or without explicit identification (Implicit VR); with Implicit

VR, the receiving application must use a DICOM data dictionary to look up the format of each data element.

1.8 Symbols and Abbreviations

AE

Application Entity

AET

Application Entity Title

DICOM

Digital Imaging and Communications in Medicine

DPPS

Data Points Per Second

IOD

Information Object Definition

PACS

Picture Archiving and Communication System

SC

Secondary Capture

SCP

Service Class Provider

SCU

Service Class User

SOP

Service-Object Pair

VR

Value Representation

VM

Value Multiplicity

XA

X-ray Angiography

This page left intentionally blank.

Chapter 3 Network Conformance Statement

1 Network Conformance Statement

1.1 Introduction

This section of the DICOM Conformance Statement specifies the Innova Vision, Innova EPVision, Innova TrackVision compliance to DICOM requirements for Networking features.

Innova Vision, Innova EPVision, Innova TrackVision is a tool that is installed on the same hardware platform as the base application, **Advantage Workstation**. This base application is a Networked Medical Imaging Console dedicated to Examination Review and Diagnosis. The workstation uses DICOM services to import acquisition images for possible further analysis or processing, and to export images and radiotherapy data to other vendors.

Innova Vision, Innova EPVision, Innova TrackVision does not have an intrinsic DICOM Network feature. It does not directly invoke the DICOM Server AE. For some detailed information on DICOM features of Advantage Windows, refer to the respective Conformance Statement - Advantage Workstation where Innova Vision, Innova EPVision, Innova TrackVision application is running.

Innova Vision, Innova EPVision, Innova TrackVision is a post processing application and it creates DERIVED XA SOP Instances of X-Ray Image SOP Class.

The application uses the following DICOM objects created in Innova system. (For a complete description of the IOD, refer to the Innova DICOM conformance statement, direction 5394268-1-8EN.)

Table 3-1:

Sop Class Name	Sop Class UID
X-Ray Image Storage	1.2.840.10008.5.1.4.1.1.12.1

The application additionally uses the 3D model (3D volume can be saved in DICOM format) or 3D SaveState (a Volume Viewer object) that is created in AW Volume Viewer. (For a complete description of the IODs, refer to the AW Volume Viewer DICOM conformance statement, direction 5344139-100.)

Table 3-2:

Sop Class Name	Sop Class UID
GEMS Private DICOM 3D Model	1.2.840.113619.4.26
Secondary Capture Image Storage (3D Save State)	1.2.840.10008.5.1.4.1.1.7

The application generates the derived instances of the following DICOM objects:

Table 3-3:

Sop Class Name	Sop Class UID
X-Ray Image Storage	1.2.840.10008.5.1.4.1.1.12.1
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7

1.2 Implementation Model

1.2.1 Application Data Flow Diagram

Refer to the respective Conformance Statement - **Advantage Workstation** where Innova Vision, Innova EPVision, Innova TrackVision application is running.

1.2.2 Presentation Context Table

Refer to the respective Conformance Statement - **Advantage Workstation** where Innova Vision, Innova EPVision, Innova TrackVision application is running.

1.2.3 Real-World Activities

The user can load 3D model from AW Volume Viewer SaveState (with a DLO) object and start the Innova Vision application.

The user can request to store the reformatted series of X-Ray images (XA). Upon this users request, Innova Vision Applications XA Model SOP Instance is created and saved into Advantage Windows database.

The user can also request to store the single frame secondary capture of the reformatted X-Ray images. This creates the Innova Vision Applications SC Model SOP Instance and saved into Advantage Windows database.

The **goal of this document** is to give a detailed description of:

- Innova Vision Applications XA Model DICOM IOD
- Innova Vision Applications SC Model DICOM IOD

1.2.4 SOP Instance UID

Implementation UID assigned to Innova Vision (or) Innova EPVision (or) Innova TrackVision is:
1.2.840.113619.6.276

An UID generated by a product has 2 parts : <root>.<suffix>.

For a GE product root is 1.2.840.113619 where

- 1 identifies ISO
- 2 identifies the ISO member body branch
- 840 identifies the country code

- 113619 identifies GEMS as a specific organization.

For a Study , Series , Instances created in GE suffix is 2.Imp.id where

- Imp identifies a specific implementation and is registered by GCC
- id is an number or a substring (i.j or i.j.l...) defined by the implementation. In our implementation it means get UID from Advantage Windows (Conformance Statement - **Advantage Workstation** Conformance Statement for DICOM.)

InnovaVision (or) InnovaEPVision (or) InnovaTrackVision will generate UIDs for instances that looks like: **1.2.840.113619.2.276**

1.2.5 Support of Extended Character Sets

The Innova Vision applications can only process the Extended Char Set ISO_IR 100 (Latin alphabet Number 1 supplementary set).

The product user interface will allow the user to enter characters from the console keyboard that are within the Extended Char Set ISO_IR 100 character set.

The Innova Vision applications can only generate images with the Extended Char Set ISO_IR 100 (Latin alphabet Number 1 supplementary set).

This page left intentionally blank.

Chapter 4 Innova Vision Applications Derived X-Ray Angiography (XA) Information Object Implementation

1 Introduction


This section specifies the requirements for Derived XA IOD that is generated by the DICOM Innova Vision Applications.

2 Innova Vision Applications Derived X-Ray Image IOD Implementation

This section defines how X-Ray Image attributes are used within the implementation, and whether these attributes are mandatory or optional for the correct operation of the application.

3 Innova Vision Applications X-Ray Image IOD Entity-Relationship Model

Illustration 4-1: Innova Vision Applications Derived X-Ray Image Entity Relationship Diagram

The Entity-Relationship diagram is shown in Fig above. In this figure, the following convention is established to represent the information organization:

Each entity is represented by a rectangular box

Each relationship is represented by a diamond shaped box

The fact that a relationship exists between two entities is depicted by lines connecting the corresponding entity boxes to the relationship boxes.

The relationships are fully defined with the maximum number of possible entities in the relationship shown. See DICOM Part 3 Section for explanation of the entity-relationship notation.

4 Entities Description

Refer to DICOM Standard Part 3 (Information Object Definitions) for a description of each of the entities contained within the X-Ray Image information object.

5 Mapping of DICOM Entities

DICOM entities map to the *Innova Vision Applications Derived XA IOD* entities in the following manner:

Table 4-1: Mapping of DICOM Entities

DICOM	Innova Vision Applications Derived XA
Patient Entity	Patient Entity (Advantage Workstation)
Study Entity	Exam Entity (Advantage Workstation)
Series Entity	Series Entity (Advantage Workstation)
Equipment Entity	Equipment Entity (Advantage Workstation)
Image Entity	Image Entity (Advantage Workstation)

6 IOD Module Table

The Innova Vision Applications XA Information Object Definitions comprise the modules of the following tables and the Private Attributes.

Table 4-2: Innova Vision Applications Derived XA Image IOD Modules

Entity Name	Module Name	Usage	Reference
Patient	Patient	Mandatory	Section 7.2, Patient Entity
Study	General Study	Mandatory	Section 7.3.1, General Study Module
	Patient Study	User Option	Section 7.3.2, Patient Study Module
Series	General Series	Mandatory	Section 7.4, Series Entity
Equipment	General Equipment	Mandatory	Section 7.5, Equipment Entity
Frame of Reference	Synchronization	Not Used	N/A
Image	General Image	Mandatory	Section 7.6.1, General Image Module
	Image Pixel	Mandatory	Section 7.6.2, Image Pixel Module
	Contrast/Bolus	Not Used	N/A
	Cine	C - Required if pixel Frame Cine data	Section 7.6.3, Cine Module
	Multiframe	C - Required if pixel data is Multi-Frame Cine data	Section 7.6.4, Multi Frame Module
	Frame Pointers	User Option	Section 7.6.5, Frame Pointers Module
	Mask	Not Used	N/A
	Display Shutter	Not Used	N/A
	Device	Not Used	N/A
	Intervention	Not Used	N/A
	X-Ray Image	Mandatory	Section 7.6.6, X-Ray Image Module
	X-Ray Acquisition	Mandatory	Section 7.6.7, X-Ray Acquisition Module
	X-Ray Collimator	Not Used	N/A
	X-Ray Table	C - Required if Image is created with table motion, may be present otherwise	Section 7.6.8, X-Ray Table Module
	XA Positioner	Mandatory	Section 7.6.9, XA Positioner Module
	DX Detector	Not Used	N/A
Overlay Plane	Not Used	N/A	
Multi-frame Overlay	Not Used	N/A	
Modality LUT	Not Used	N/A	

Entity Name	Module Name	Usage	Reference
	VOI LUT	User Option	Section 7.6.10, VOI LUT Module
	SOP Common	Mandatory	Section 7.6.11, SOP Common Module

7 Information Module Definitions

Please refer to DICOM v3.0 Standard Part 3 (Information Object Definitions) for a description of each of the entities and modules contained within the XA Information Object.

The following modules are included to convey Enumerated Values, Defined Terms, and Optional Attributes supported. Type 1 & Type 2 Attributes are also included for completeness and to define what values they may take and where these values are obtained from. It should be noted that they are the same ones as defined in the DICOM v3.0 Standard Part 3 (Information Object Definitions).

7.1 Common Information Module Definitions

Please refer to DICOM Part 3 (Information Object Definitions) for a description of each of the entities, modules, and attributes contained within the XA Information Objects.

7.2 Patient Entity

Patient Module

Table 4-3: Patient Module Attributes

Attribute Name	Tag	Type	Attribute Description
Patient's Name	(0010,0010)	2	Copied from the 3D model without change
Patient ID	(0010,0020)	2	Copied from the 3D model without change
Patient's Birth Date	(0010,0030)	2	Copied from the 3D model without change
Patient's Sex	(0010,0040)	2	Copied from the 3D model without change
Other Patient ID	(0010,1000)	3	Copied from the 3D model without change

7.3 Study Entity

7.3.1 General Study Module

Table 4-4: General Study Module Attributes

Attribute Name	Tag	Type	Attribute Description
study date	(0008,0020)	2	Date the Innova Vision Applications exam is started in AW
study time	(0008,0030)	2	Time the Innova Vision Applications exam is started in AW
accession number	(0008,0050)	2	EMPTY
referring physician name	(0008,0090)	2	EMPTY
study description	(0008,1030)	3	"InnovaVision" OR "InnovaTrackVision" OR "InnovaEP-Vision"
study instance uid	(0020,000D)	1	Uniquely generated for each Innova Vision Applications exam started
study id	(0020,0010)	2	Hardcoded to "1"

7.3.2 Patient Study Module

Table 4-5: Patient Study Module Attributes

Attribute Name	Tag	Type	Attribute Description
Patient's Age	(0010,1010)	3	Copied from the 3D model without change
Patient's Size	(0010,1020)	3	Copied from the 3D model without change
Patient's Weight	(0010,1030)	3	Copied from the 3D model without change

7.4 Series Entity

General Series Module

Table 4-6: General Series Module Attributes

Attribute Name	Tag	Type	Attribute Description
series date	(0008,0021)	3	same as acquisition date of the original XA image that resulted in this series
series time	(0008,0031)	3	same as acquisition time of the original XA image that resulted in this series
series description	(0008,103E)	3	"InnovaVision" OR "InnovaTrackVision" OR "InnovaEP-Vision"
modality	(0008,0060)	1	"XA"
protocol name	(0018,1030)	3	"InnovaVision" OR "InnovaTrackVision" OR "InnovaEP-Vision"
patient position	(0018,5100)	2C	Copied from the original XA image without change
Series Instance UID	(0020,000E)	1	Generated uniquely in AW
series number	(0020,0011)	2	Generated in AW. Unique within the given Innova Vision Applications exam.
Related Series Sequence	(0008,1250)	3	One item only, to reference the Innova Series where X-Ray image (Fluoro or Record) was acquired
>Study Instance UID	(0020,000D)	1	Copied from the original XA image without change
>Series Instance UID	(0020,000E)	1	Copied from the original XA image without change
>Purpose of Reference Code Sequence	(0040,A170)	2	One item only, to describe the purpose of this reference - DCID 7210 that specifies the purpose "Simultaneously Acquired"
>>Code Value	(0008,0100)	1C	"122400"
>>Coding Scheme Designator	(0008,0102)	1C	"DCM"
>>Code Meaning	(0008,0104)	1C	"Simultaneously Acquired"

7.5 Equipment Entity

General Equipment Module

Table 4-7: General Equipment Module Attributes

Attribute Name	Tag	Type	Attribute Description
Manufacturer	(0008,0070)	2	Hardcoded to "GE MEDICAL SYSTEMS"
Institution Name	(0008,0080)	3	Institution where the equipment is located that produced the fused images (i.e, the AW system).
Station Name	(0008,1010)	3	User defined name identifying the machine that produced the fused images, typically the AETitle of the AW
Manufacturer's Model Name	(0008,1090)	3	Manufacturer's model name of the equipment, which is the Innova Vision Applications application name used to generate this image
Device Serial Number	(0018,1000)	3	Manufacturer's serial number of the equipment, which is the Serial Number of the AW machine
Software Version(s)	(0018,1020)	3	Has 2 values. Manufacturer's designation of the software version of the AW and the Innova Vision Applications application.

7.6 Image Entity

7.6.1 General Image Module

Table 4-8: General Image Module Attributes

Attribute Name	Tag	Type	Attribute Description
image type	(0008,0008)	3	"DERIVED\PRIMARYAUGFLUORO"
Acquisition Date	(0008,0022)	3	Same as acquisition date of the original XA image
Acquisition Time	(0008,0032)	3	Same as acquisition time of the original XA image
content date	(0008,0023)	2C	The date the F+ image pixel data creation started in AW.
content time	(0008,0033)	2C	The time the F+ image pixel data creation started in AW.
Referenced Instance Sequence	(0008,114A)	3	Sequence of one Item providing the referenced 3D model (3D image, planning information, markers) that is used for the fusion.
>Referenced sop class uid	(0008,1150)	1C	SOP Class UID of the referenced 3D model
>Referenced sop instance uid	(0008,1155)	1C	SOP Instance UID of the referenced 3D model
>Purpose of Reference Code Sequence	(0040,A170)	1	One item only, to describe the purpose of this reference
>>Code Value	(0008,0100)	1C	"INNOVA-001" - GE private Code value that specifies the purpose "Corresponding volumetric presentation"
>>Coding Scheme Designator	(0008,0102)	1C	"99GEMS" - GE private Code Scheme Designator that specifies the purpose "Corresponding volumetric presentation"
>>Code Meaning	(0008,0104)	1C	"Corresponding volumetric presentation" - GE private Code Meaning that specifies the purpose "Corresponding volumetric presentation"
instance number	(0020,0013)	2	starting at 1. Unique within a series
patient orientation	(0020,0020)	2C	sent as a zero-length empty field
image comments	(0020,4000)	3	User-defined comments about the image

7.6.2 Image Pixel Module

Table 4-9: Image Pixel Module Attributes

Attribute Name	Tag	Type	Attribute Description
Samples per Pixel	(0028,0002)	1	"1"
Photometric Interpretation	(0028,0004)	1	"MONOCHROME2"
Rows	(0028,0010)	1	"1024"
Columns	(0028,0011)	1	"1024"
Bits Allocated	(0028,0100)	1	"8"
Bits Stored	(0028,0101)	1	"8"
High Bit	(0028,0102)	1	"7"
Pixel Representation	(0028,0103)	1	"0"
Pixel Data	(7FE0,0010)	1	A data stream of the pixel samples which comprise the Image

7.6.3 Cine Module

Table 4-10: Cine Module Attributes

Attribute Name	Tag	Type	Attribute Description
start trim	(0008,2142)	3	"1"
stop trim	(0008,2143)	3	The frame number of the last frame of the Multiframe image to be displayed
recommended display frame rate	(0008,2144)	3	Copied from the original XA image without change
cine rate	(0018,0040)	3	Copied from the original XA image without change
frame time vector	(0018,1065)	1C	Values derived from the original XA image
frame delay	(0018,1066)	3	"0"

7.6.4 Multi Frame Module

Table 4-11: Multi Frame Module Attributes

Attribute Name	Tag	Type	Attribute Description
number of frames	(0028,0008)	1	Number of frames in a Multi-frame Image
frame increment pointer	(0028,0009)	1	"(0018,1065)"

7.6.5 Frame Pointers Module

Table 4-12: Frame Pointers Module Attributes

Attribute Name	Tag	Type	Attribute Description
representative frame number	(0028,6010)	3	The frame number selected for use as a pictorial representation (e.g. icon) of the Multi-frame Image

7.6.6 X-Ray Image Module

Table 4-13: X-Ray Image Module Attributes

Attribute Name	Tag	Type	Attribute Description
pixel intensity relationship	(0028,1040)	1	"DISP"
Image Type	(0008,0008)	1	Refer to Section 7.6.1
Samples per pixel	(0028,0002)	1	Refer to Section 7.6.2
Photometric Interpretation	(0028,0004)	1	Refer to Section 7.6.2
Bits Allocated	(0028,0100)	1	Refer to Section 7.6.2
Bits Stored	(0028, 0101)	1	Refer to Section 7.6.2
High Bit	(0028,0102)	1	Refer to Section 7.6.2
Pixel Representation	(0028, 0103)	1	Refer to Section 7.6.2

7.6.7 X-Ray Acquisition Module

Table 4-14: X-Ray Acquisition Module Attributes

Attribute Name	Tag	Type	Attribute Description
Kvp	(0018,0060)	2	EMPTY
exposure	(0018,1152)	2C	EMPTY
radiation setting	(0018,1155)	1	"GR" OR "SC"
intensifier size	(0018,1162)	3	Copied from the original XA image without change
field of view shape	(0018,1147)	3	"RECTANGLE"
field of view dimension(s)	(0018,1149)	3	Copied from the original XA image without change

7.6.8 X-Ray Table Module

Table 4-15: X-Ray Table Module Attributes

Attribute Name	Tag	Type	Attribute Description
table motion	(0018,1134)	2	"STATIC" OR "DYNAMIC"
table vertical increment	(0018,1135)	2C	Values derived from the original XA image
table lateral increment	(0018,1136)	2C	Values derived from the original XA image
table longitudinal increment	(0018,1137)	2C	Values derived from the original XA image
Table Angle	(0018,1138)	3	Values derived from the original XA image

7.6.9 XA Positioner Module

Table 4-16: XA Positioner Module Attributes

Attribute Name	Tag	Type	Attribute Description
distance source to detector	(0018,1110)	3	Copied from the original XA image without change
distance source to patient	(0018,1111)	3	Copied from the original XA image without change
positioner motion	(0018,1500)	2C	"STATIC" OR "DYNAMIC"

Attribute Name	Tag	Type	Attribute Description
positioner primary angle	(0018,1510)	2	Values derived from the original XA image
positioner secondary angle	(0018,1511)	2	Values derived from the original XA image
positioner primary angle increment	(0018,1520)	2C	Values derived from the original XA image
positioner secondary angle increment	(0018,1521)	2C	Values derived from the original XA image

7.6.10 VOI LUT Module

Table 4-17: VOI LUT Module Attributes

Attribute Name	Tag	Type	Attribute Description
window center	(0028,1050)	3	"128"
window width	(0028,1051)	3	"256"

7.6.11 SOP Common Module

Table 4-18: SOP Common Module Attributes

Attribute Name	Tag	Type	Attribute Description
specific character set	(0008,0005)	1C	"ISO_IR 100"
sop class uid	(0008,0016)	1	"1.2.840.10008.5.1.4.1.1.1 2.1"
sop instance uid	(0008,0018)	1	Uniquely identifies the SOP Instance.
Contributing Equipment Sequence	(0018,A001)	3	One item shall be present for Innova equipment
>Purpose of Reference Code Sequence	(0040,A170)	1	Describes the purpose for which the related equipment is being reference.
>>Code Value	(0008,0100)	1C	"109101"
>>Coding Scheme Designator	(0008,0102)	1C	"DCM"
>>Code Meaning	(0008,0104)	1C	"Acquisition Equipment"
>Manufacturer	(0008,0070)	1	"GE MEDICAL SYSTEMS"
>Manufacturer's Model Name	(0008,1090)	3	"DL"
>Contribution Description	(0018,A003)	3	"X-Ray Acquisition equipment"

7.7 Standard Extended and Private Data Attributes

The Product supports the Standard and Private Attributes defined in the following sections in Standard Extended X-Ray SOP Instances as Type 3 data elements.

7.7.1 Standard Attributes

The Product supports the following attributes in SOP Instances as Type 3 data elements.

Table 4-19: Standard Extended Attributes

Attribute Name	Tag	Information Entity Name - Image	Use
Curve Dimensions	(5000,0005)	3	"2"
Number of Points	(5000,0010)	3	Number of data points in this Curve
Type of Data	(5000,0020)	3	"ECG"
Data Value Representation	(5000,0103)	3	"0000H"
Curve Data Descriptor	(5000,0110)	3	"01"
Axis Units	(5000,0030)	3	"DPPSINONE"
Coordinate Start Value	(5000,0112)	3	"0"
Coordinate Step Value	(5000,0114)	3	"250"
Curve Data	(5000,3000)	3	Curve Data points in the curve, each dimension for the first point, followed by dimensions for second point

7.7.2 Private Group DLX_SERIE_01

Private Group DLX_SERIE_01 is modeled as part of the Image (X-Ray Table and X-Ray Positioner) Information Entity.

Table 4-20: Private Group DLX_SERIE_01

Attribute Name	Tag	VR	VM	Attribute Description and Use
Angle value 1	(0019,xx01)	DS	1	Copied from the original XA image without change
Angle value 2	(0019,xx02)	DS	1	Copied from the original XA image without change
Angle value 3	(0019,xx03)	DS	1	Copied from the original XA image without change
Table Vertical Position	(0019,xx21)	DS	1	Copied from the original XA image without change
Table Longitudinal Position	(0019,xx22)	DS	1	Copied from the original XA image without change
Table Lateral Position	(0019,xx23)	DS	1	Copied from the original XA image without change

7.7.3 Private Group GEMS_DL_IMG_01

Private Group GEMS_DL_IMG_01 is modeled as part of the Image (General Image, X-Ray Positioner, X-Ray Table and Curve) Information Entity.

Table 4-21: Private Group GEMS_DL_IMG_01

Attribute Name	Tag	VR	VM	Attribute Description and Use
patient position per image	(0019,xxC7)	CS	1	Patient position descriptor relative to the equipment. Defined terms are: HFP = Head First-Prone HFS = Head First-Supine HFDR = Head First-Decubitus Right HFDL = Head First-Decubitus Left FFDR = Feet First-Decubitus Right FFDL = Feet First-Decubitus Left FFP = Feet First-Prone FFS = Feet First-Supine Values derived from the original XA image
fov dimension double	(0019,xx0B)	DS	1-2	Value in floating point resolution, whose truncature is (0018,1149). Possible values are "400\400" OR "320\320" OR "300\300" OR "200\200" OR "172.8\172.8" OR "170\170" OR "160\160" OR "150\150" OR "147.2\147.2" OR "121.6\121.6" OR "120\120" Values derived from the original XA image.
image detector rotation angle	(0019,xx92)	DS	1	image rotation at the detector reading, before image flip. Values derived from the original XA image
image flip	(0019,xx95)	CS	2	Contain two defined terms, First value: <ul style="list-style-type: none"> • YES if horizontal flip • NO else Second value : <ul style="list-style-type: none"> • YES if vertical flip • NO else Values derived from the original XA image
can downscan 512	(0019,xxAA)	CS	1	Indicates the possibility to downscan the pixel data to 512x512 for exchange purposes: <ul style="list-style-type: none"> • YES • NO Values derived from the original XA image
Acquisition Mode Description	(0019,xxB1)	LO	1	The precise description of the "numerical code" (Adx acq mode). Values derived from the original XA image
Acquisition Mode Display Label	(0019,xxB2)	LO	1	Label that shall be displayed on the AW browser, for each sequence Values derived from the original XA image
Acquisition Region	(0019,xxBA)	CS	1	Coded String to determine whether the acquisition is Cardiac or Angio. Defined terms are CARDIAC, ANGIO and UNKNOWN Values derived from the original XA image
Acquisition SUB mode	(0019,xxBB)	CS	1	Coded String to determine whether the acquisition mode was designed for a subtracted or Non-subtracted review. Defined terms are SUB, NOSUB and UNKNOWN. Values derived from the original XA image

Attribute Name	Tag	VR	VM	Attribute Description and Use
table head tilt angle	(0019,xxEE)	FL	1	Angle of the head-feet axis of the table in degrees relative to the horizontal plane. Zero is defined when the head-feet axis is in the horizontal plane. Positive values are when the head of the table is upwards the horizontal plane. The valid range is from -45 to +45. Contains the value of the first frame. Values derived from the original XA image
table rotation angle	(0019,xxEA)	FL	1	Rotation angle in degrees of the table in the horizontal plane (clockwise when looking from above the table). The value of this attribute applies to the first frame of the multiframe image. Values derived from the original XA image
table cradle angle	(0019,xxBC)	FL	1	Angle of the left-right axis of the table in degrees relative to the horizontal plane. Positive values indicate that the left of the table is upwards. The value of this attribute applies to the first frame of the Multi-frame image. Values derived from the original XA image
table rotation status vector	(0019,xxBD)	CS	1-N	Status of the rotation of the table in the horizontal plane: <ul style="list-style-type: none"> • YES • NO This is a multi-valued attribute that contains the Table Rotation Status for each frame. Values derived from the original XA image
table rotation angle increment	(0019,xxC3)	FL	1-N	Incremental change in the rotation of the table in the horizontal plane (clockwise when looking from above the table) relative to the first frame of the Multi-frame image. Contains as many values as number of frames. Required if Table Motion is DYNAMIC. Values derived from the original XA image
Table X Position to Iso-center	(0019,xxEB)	FL	1	X position of the Table Reference Point with respect to the Isocenter (mm). Positive values are towards the LAO direction of the Isocenter. The value of this attribute applies to the first frame of the Multi-frame image. Values derived from the original XA image
Table Y Position to Iso-center	(0019,xxEC)	FL	1	Y position of the Table Reference Point with respect to the Isocenter (mm). positive values are downwards the horizontal plane in the vertical direction. The value of this attribute applies to the first frame of the Multi-frame image. Values derived from the original XA image
Table Z Position to Iso-center	(0019,xxED)	FL	1	Z position of the Table Reference Point with respect to the Isocenter (mm). Positive values are towards the CRA direction of the Isocenter. The value of this attribute applies to the first frame of the Multi-frame image. Values derived from the original XA image

Attribute Name	Tag	VR	VM	Attribute Description and Use
Table X Position to Iso-center increment	(0019,xxD7)	FL	1-N	Incremental change in X position of the Table Reference Point with respect to the Isocenter (in mm), relative to the first frame of the Multi-frame image. Positive values are towards the LAO direction of the Isocenter. Contains as many values as number of frames. Required if Table Motion is DYNAMIC. Values derived from the original XA image
Table Y Position to Iso-center increment	(0019,xxD8)	FL	1-N	Incremental change in Y position of the Table Reference Point with respect to the Isocenter (in mm), relative to the first frame of the Multi-frame image. Positive values are downwards the horizontal plane in the vertical direction. Contains as many values as number of frames. Required if Table Motion is DYNAMIC. Values derived from the original XA image
Table Z Position to Iso-center increment	(0019,xxD9)	FL	1-N	Incremental change in Z position of the Table Reference Point with respect to the Isocenter (in mm), relative to the first frame of the Multi-frame image. Positive values are towards the CRA direction of the Isocenter. Contains as many values as number of frames. Required if Table Motion is DYNAMIC. Values derived from the original XA image
Table Head Tilt Angle increment	(0019,xxDA)	FL	1-N	Incremental change in the angle of the headfeet axis of the table in degrees with respect to the horizontal plane, relative to the first frame of the Multi-frame image. Zero means when the headfeet axis is in the horizontal plane. Positive values are when the head of the table is upwards the horizontal plane. Contains as many values as number of frames. Required if Table Motion is DYNAMIC. Values derived from the original XA image
angle 1 increment	(0019,xx97)	DS	1-N	Incremental change in angle_value_1, sent if positioner motion is dynamic. Values derived from the original XA image
angle 2 increment	(0019,xx98)	DS	1-N	Incremental change in angle_value_2, sent if positioner motion is dynamic. Values derived from the original XA image
angle 3 increment	(0019,xx99)	DS	1-N	Incremental change in angle_value_3, sent if positioner motion is dynamic. Values derived from the original XA image
SID vector	(0019,xxBE)	FL	1-N	Distance in mm from source to detector center. This is a multi-valued attribute that contains the SID for each frame. Values derived from the original XA image
SOD vector	(0019,xxE9)	FL	1-N	Distance in mm from source to the system isocenter. This is a multi-valued attribute that contains the SOD for each frame. Values derived from the original XA image
Number of Points_before acq	(0019,xx60)	US	1	Number of data points in the ECG Curve before the acquisition.
Curve Data_before acq	(0019,xx61)	OW	1	Curve Data points in the curve before the XA acquisition

Attribute Name	Tag	VR	VM	Attribute Description and Use
Number of Points_trigger	(0019,xx62)	US	1	Number of data points of QRS trigger values in ECG data.
Curve Data_trigger	(0019,xx63)	OW	1	Curve Data points in the curve that describes the QRS trigger values.
ECG Synchronization	(0019,xx64)	SH	1	ECG synchronization ON / OFF
ECG Delay Mode	(0019,xx65)	SH	1	Prospective / Retrospective gating method used
ECG Delay vector	(0019,xx66)	IS	1-N	ECG % (integer value between 0 and 100) applied at frame acquisition time.

7.7.4 Private Group GEMS_XR3DCAL_01

Table 4-22: Private Group GEMS_XR3DCAL_01

Attribute Name	Tag	VR	VM	Attribute Description and Use
Generalized calibration	(0021,xx20)	LT	1	Augmented calibration string of Innova system

7.7.5 Private Group GEMS_3D_INTVL_01

Table 4-23: Private Group GEMS_3D_INTVL_01

Attribute Name	Tag	VR	VM	Attribute Description and Use
X-Ray Marker Sequence	(0023,xx01)	SQ	1	Sequence of markers
> Marker Id	(0023,xx02)	SH	1	Id of the Marker. Alphanumeric, limited to 16 characters
> Marker Type	(0023,xx03)	CS	1	Defined terms: POINTS, LINE (can be extended)
> Marker Size	(0023,xx04)	FL	1	Size of the marker as provided by the user. Float number.
> Marker Color CIELab Value	(0023,xx05)	US	3	A color that is chosen by the user. It is a triplet value , units in PCS-Values and values in CIE-Lab
> Marker Label	(0023,xx06)	LO	1	Displayed (translated) label. Alphanumeric, limited to 64 characters
> Marker Visible State	(0023,xx07)	CS	1	Visible state as chosen by the user during the store of F+ image. Enumerated terms (non extensible): • YES • NO
> Marker Description	(0023,xx08)	LO	1	Alphanumeric, limited to 64 characters
> Marker Points Sequence	(0023,xx10)	SQ	1	If present shall contain one or more items, each item is a point in the marker
>> Marker Point Id	(0023,xx11)	SH	1	Id of the point. Alphanumeric, limited to 16 characters

Attribute Name	Tag	VR	VM	Attribute Description and Use
>> Marker Point Position	(0023,xx12)	FL	3	Float coordinates x,y,z of the point
>> Marker Point Size	(0023,xx13)	FL	1	Size of marker. Float number.
>> Marker Point Color CIELab Value	(0023,xx14)	US	3	A color triplet value, units in PCSValues and values in CIELab
>> Marker Point Visible State	(0023,xx16)	CS	1	Visible state as chosen by the user during the store of F+ image. Enumerated terms (non extensible): <ul style="list-style-type: none"> • YES • NO
>> Marker Point Order	(0023,xx17)	IS	1	Order in the list of points. Integer.
Volume Manual Registration	(0023,xx18)	FL	3	Float coordinates x,y,z of the translation
Volumes Threshold	(0023,xx20)	IS	1-N	list of volume rendering thresholds in the augmented fluoro object (in mm) as set by user.
CutPlane activation flag	(0023,xx25)	CS	1	Enumerated values <ul style="list-style-type: none"> • YES • NO <p>Indicates if the clipping plane is activated by the user or not.</p>
CutPlane position Value	(0023,xx26)	IS	1	Clipping plane cursor position (selected by the user)
CutPlane normal Value	(0023,xx27)	FL	3	Clipping plane normal coordinates (Float coordinates x,y,z)
Volume scaling factor	(0023,xx28)	FL	1	Default value = 1 (no scaling). Magnification factor for the volume linked to the F+ sequence.
ROI to TableTop Distance	(0023,xx29)	FL	1	Estimated distance between the center of the anatomy / region of interest to the table top. Set by the user during storing Innova Vision Applications image.
DRR Threshold	(0023,xx30)	IS	1-N	list of DRR rendering thresholds in the augmented fluoro object
Volume Table position	(0023,xx31)	FL	3	Float coordinates x,y,z of the table position configured for the volume
Rendering mode	(0023,xx32)	IS	1	The rendering mode that was set by user. Values can be 0,1 or 2. Volume = 0, Surface = 1, DRR = 2.
3D Object Opacity	(0023,xx33)	IS	1	Set by the user. Value can range from 0-100.
Invert Image	(0023,xx34)	IS	1	Set by the user. 0 - Normal, 1 - Inverted
Enhance Full	(0023,xx35)	IS	1	Set by the user. 0 - Masked, 1 - Full
Zoom	(0023,xx36)	FL	1	Set by the user. Any positive value.
Roam	(0023,xx37)	IS	2	Set by the user. Horizontal and vertical direction Roam. Value depends on zoom factor and can be any integer value.

Attribute Name	Tag	VR	VM	Attribute Description and Use
window level	(0023,xx38)	IS	1	Set by the user. Value range from - 1024 to 1023
window width	(0023,xx39)	IS	1	Set by the user. Value range from 1 to 2047
BMC Setting	(0023,xx40)	CS	1	Set by User: <ul style="list-style-type: none"> • ON • OFF
Back View Setting	(0023,xx41)	CS	1	Set by User: <ul style="list-style-type: none"> • ON • OFF
Sub Volume Visibility	(0023,xx42)	CS	1-N	Set by User: <ul style="list-style-type: none"> • YES • NO
3D Landmarks Visibility	(0023,xx43)	CS	1	Set by User: <ul style="list-style-type: none"> • YES • NO
Ablation Point Visibility	(0023,xx44)	CS	1	Set by User: <ul style="list-style-type: none"> • YES • NO

This page left intentionally blank.

Chapter 5 Innova Vision Applications Secondary Information Object Implementation

1 Introduction


This section specifies the requirements for Secondary Capture IOD that is generated by the DICOM Innova Vision Applications.

2 Innova Vision Applications Secondary Image IOD Implementation

This section defines how X-Ray image attributes are used within the implementation, and whether these attributes are mandatory or optional for the correct operation of the application.

3 Innova Vision Applications Secondary Image IOD Entity - Relationship Model

Illustration 5-1: Innova Vision Applications Secondary Image Entity Relationship Diagram

The Entity-Relationship diagram is shown in Figure above. In this figure, the following convention is established to represent the information organization:

Each entity is represented by a rectangular box.

Each relationship is represented by a diamond shaped box.

The fact that a relationship exists between two entities is depicted by lines connecting the corresponding entity boxes to the relationship boxes.

The relationships are fully defined with the maximum number of possible entities in the relationship shown. See DICOM Part 3 Section for explanation of the entity-relationship notation.

4 Entities Description

Refer to DICOM Standard Part 3 (Information Object Definitions) for a description of each of the entities contained within the Secondary Image information object.

5 Mapping of DICOM Entities

DICOM entities map to the Innova Vision Applications Secondary IOD entities in the following manner:

Table 5-1: Mapping of DICOM Entities

DICOM	Innova Vision Applications Secondary
Patient Entity	Patient Entity (Advantage Workstation)
Study Entity	Exam Entity (Advantage Workstation)
Series Entity	Series Entity (Advantage Workstation)
Equipment Entity	Equipment Entity (Advantage Workstation)
Image Entity	Image Entity (Advantage Workstation)

6 IOD Module Table

The Innova Vision Applications Secondary Information Object Definitions comprise the modules of the following tables and the Private Attributes.

Table 5-2: Innova Vision Applications Secondary Image IOD Modules

Entity Name	Module Name	Usage	Reference
Patient	Patient	Mandatory	Section 7.2, Patient Entity
Study	General Study	Mandatory	Section 7.3.1, General Study Module
	Patient Study	User Option	Section 7.3.2, Patient Study Module
Series	General Series	Mandatory	Section 7.4, Series Entity
Equipment	General Equipment	Mandatory	Section 7.5.1, General Equipment Module
	SC Equipment	Mandatory	Section 7.5.2, SC Equipment Module
Image	General Image	Mandatory	Section 7.6.1, General Image Module
	Image Pixel	Mandatory	Section 7.6.2, Image Pixel Module
	SC Image	Mandatory	Section 7.6.3, SC Image Module
	Device	Not Used	N/A
	Overlay Plane	Not Used	N/A
	Modality LUT	Not Used	N/A
	VOI LUT	Not Used	N/A
	SOP Common	Mandatory	Section 7.6.4, SOP Common Module

7 Information Module Definitions

Please refer to DICOM V3.0 Standard Part 3 (Information Object Definitions) for a description of each of the entities and modules contained within the SC Information Object.

The following modules are included to convey Enumerated Values, Defined Terms, and Optional Attributes supported. Type 1 and Type 2 Attributes are also included for completeness and to define which values they may take, and where these values are obtained from. It should be noted that they are the same ones as defined in the DICOM V3.0 Standard Part 3 (Information Object Definitions).

Also note that Attributes not present in tables are not supported.

7.1 Common Information Module Definitions

Please refer to DICOM Part 3 (Information Object Definitions) for a description of each of the entities, modules, and attributes contained within the SC Information Objects.

7.2 Patient Entity

Patient Module

Table 5-3: Patient Module Attributes

Attribute Name	Tag	Type	Attribute Description
Patient's Name	(0010,0010)	2	Copied from the 3D model without change
Patient ID	(0010,0020)	2	Copied from the 3D model without change
Patient's Birth Date	(0010,0030)	2	Copied from the 3D model without change
Patient's Sex	(0010,0040)	2	Copied from the 3D model without change
Other Patient ID	(0010,1000)	3	Copied from the 3D model without change

7.3 Study Entity

7.3.1 General Study Module

Table 5-4: General Study Module Attributes

Attribute Name	Tag	Type	Attribute Description
study date	(0008,0020)	2	Date the Innova Vision Applications exam is started in AW
study time	(0008,0030)	2	Time the Innova Vision Applications exam is started in AW
accession number	(0008,0050)	2	EMPTY
referring physician name	(0008,0090)	2	EMPTY
study description	(0008,1030)	3	Innova Vision OR Innova TrackVision OR Innova EP-Vision
study instance UID	(0020,000D)	1	Uniquely generated for each Innova Vision Applications exam started

Attribute Name	Tag	Type	Attribute Description
study ID	(0020,0010)	2	Hard coded to "1"

7.3.2 Patient Study Module

Table 5-5: Patient Study Module Attributes

Attribute Name	Tag	Type	Attribute Description
Patient's age	(0010,1010)	3	Copied from the 3D model without change
Patient's size	(0010,1020)	3	Copied from the 3D model without change
Patient's weight	(0010,1030)	3	Copied from the 3D model without change

7.4 Series Entity

General Series Module

Table 5-6: General Series Module Attributes

Attribute Name	Tag	Type	Attribute Description
series date	(0008,0021)	3	same as acquisition date of the original XA image that resulted in this series
series time	(0008,0031)	3	same as acquisition time of the original XA image that resulted in this series
series description	(0008,103E)	3	Innova Vision OR Innova TrackVision OR Innova EP-Vision
modality	(0008,0060)	1	XA
protocol name	(0018,1030)	3	Innova Vision OR Innova TrackVision OR Innova EP-Vision
patient position	(0018,5100)	2C	Copied from the original XA image without change
series instance UID	(0020,000E)	1	Generated uniquely in AW
series number	(0020,0011)	2	Generated in AW. Unique within the given Innova Vision Applications exam
related series sequence	(0008,1250)	3	One item only, to reference the Innova Series where the X-Ray image (Fluoro or Record) was acquired
>Study Instance UID	(0020,000D)	1	Copied from the original XA image without change
>Series Instance UID	(0020,000E)	1	Copied from the original XA image without change
>Purpose of Reference Code Sequence	(0040,A170)	2	One item only, to describe the purpose of this reference
>>Code Value	(0008,0100)	1C	"122400"
>>Coding Scheme Designator	(0008,0102)	1C	"DCM"
>>Code Meaning	(0008,0104)	1C	"Simultaneously Acquired"

7.5 Equipment Entity

7.5.1 General Equipment Module

Table 5-7: General Equipment Module Attributes

Attribute Name	Tag	Type	Attribute Description
Manufacturer	(0008,0070)	2	Hard coded to "GE MEDICAL SYSTEMS"
Institution Name	(0008,0080)	3	Institution where the equipment is located that produced the fused images (i.e, the AW system)
Station Name	(0008,1010)	3	User defined name identifying the machine that produced the fused images, typically the AETitle of the AW
Manufacturer's Model Name	(0008,1090)	3	Manufacturer's model name of the equipment, which is the Innova Vision Applications application name used to generate this image
Device Serial Number	(0018,1000)	3	Manufacturer's serial number of the equipment, which is the Serial Number of the AW machine
Software Version(s)	(0018,1020)	3	Has 2 values: Manufacturer's designation of the software version of the AW and the Innova Vision Applications application

7.5.2 SC Equipment Module

Table 5-8: SC Equipment Module Attributes

Attribute Name	Tag	Type	Attribute Description
conversion type	(0008,0064)	1	WSD
SC manufacturer	(0018,1016)	3	GE MEDICAL SYSTEMS
SC manufacturer model name	(0018,1018)	3	Manufacturer's model number of the Secondary Capture Device (i.e, AW)

7.6 Image Entity

7.6.1 General Image Module

Table 5-9: General Image Module Attributes

Attribute Name	Tag	Type	Attribute Description
Image type	(0008,0008)	3	DERIVED\SECONDARY\AUGFLUORO
Acquisition Date	(0008,0022)	3	Same as content date
Acquisition Time	(0008,0032)	3	Same as content time
Content date	(0008,0023)	2C	Date the Innova Vision Applications SC image pixel data creation started in AW.
Content time	(0008,0033)	2C	Time the Innova Vision Applications image pixel data creation started in AW.
Instance number	(0020,0013)	2	starting at 1. Unique within a serie
Patient orientation	(0020,0020)	2C	Sent as a zero-length empty field
Image comments	(0020,4000)	3	User-defined comments about the image

Attribute Name	Tag	Type	Attribute Description
Burned in annotation	(0028,0301)	3	NO

7.6.2 Image Pixel Module

Table 5-10: Image Pixel Module Attributes

Attribute Name	Tag	Type	Attribute Description
Samples per Pixel	(0028,0002)	1	3
Photometric Interpretation	(0028,0004)	1	RGB
Rows	(0028,0010)	1	1000
Columns	(0028,0011)	1	1000
Bits Allocated	(0028,0100)	1	8
Bits Stored	(0028,0101)	1	8
High Bit	(0028,0102)	1	7
Pixel Representation	(0028,0103)	1	0
Planar Configuration	(0028,0006)	1C	0
Pixel Data	(7FE0,0010)	1	A data stream of the pixel samples which comprise the SC Image

7.6.3 SC Image Module

Table 5-11: SC Image Module Attributes

Attribute Name	Tag	Type	Attribute Description
SC date	(0018,1012)	3	Date of storing the SC image
SC time	(0018,1014)	3	Time of storing the SC image

7.6.4 SOP Common Module

Table 5-12: SOP Common Module Attributes

Attribute Name	Tag	Type	Attribute Description
Specific character set	(0008,0005)	1C	ISO_IR 100
SOP class UID	(0008,0016)	1	1.2.840.10008.5.1.4.1.1.7
SOP instance UID	(0008,0018)	1	Uniquely identifies the SOP Instance
Contributing Equipment Sequence	(0018,A001)	3	One item shall be present for Innova equipment
>Purpose of Reference Code Sequence	(0040,A170)	1	Describes the purpose for which the related equipment is being reference.
>>Code Value	(0008,0100)	1C	109101
>>Coding Scheme Designator	(0008,0102)	1C	DCM
>>Code Meaning	(0008,0104)	1C	Acquisition Equipment
>Manufacturer	(0008,0070)	1	GE MEDICAL SYSTEMS
>Manufacturer's Model Name	(0008,1090)	3	DL

Attribute Name	Tag	Type	Attribute Description
>Contribution Description	(0018,A003)	3	X-Ray Acquisition equipment

7.7 Standard Extended and Private Data Attributes

The Product supports the Standard and Private Attributes defined in the following sections in Standard Extended SC SOP Instances as Type 3 data elements.

7.7.1 Standard Attributes

The Product supports the following attributes in SOP Instances as Type 3 data elements.

Table 5-13: Standard Extended Attributes

Attribute Name	Tag	Information Entity Name - Image	Use
Kvp	(0018,0060)	3	EMPTY
Exposure	(0018,1152)	3	EMPTY
distance source to detector	(0018,1110)	3	Copied from the original XA image without change
distance source to patient	(0018,1111)	3	Copied from the original XA image without change
positioner motion	(0018,1500)	3	"STATIC" or "DYNAMIC"
positioner primary angle	(0018,1510)	3	Values derived from the original XA image
positioner secondary angle	(0018,1511)	3	Values derived from the original XA image

7.7.2 Private Group DLX_SERIE_01

Private Group DLX_SERIE_01 is modeled as part of the Image (X-Ray Positioner) Information Entity.

Table 5-14: Private Group DLX_SERIE_01

Attribute Name	Tag	VR	VM	Attribute Description and use
Angle value 1	(0019,xx01)	DS	1	Copied from the original XA image without change
Angle value 2	(0019,xx02)	DS	1	Copied from the original XA image without change
Angle value 3	(0019,xx03)	DS	1	Copied from the original XA image without change

7.7.3 Private Group GEMS_DL_IMG_01

Private Group GEMS_DL_IMG_01 is modeled as part of the Image (General Image, X-Ray Positioner, X-Ray Table and Curve) Information Entity.

Table 5-15: Private Group GEMS_DL_IMG_01

Attribute Name	Tag	VR	VM	Attribute Description and use
Patient position per image	(0019,xxC7)	CS	1	Patient position descriptor relative to the equipment. Defined terms are: HFP = Head First-Prone HFS = Head First-Supine HFDR = Head First-Decubitus Right HFDL = Head First-Decubitus Left FFDR = Feet First-Decubitus Right FFDL = Feet First-Decubitus Left FFP = Feet First-Prone FFS = Feet First-Supine Values derived from the original XA image
Image flip	(0019,xx95)	CS	2	Contain two defined terms: First value : YES if horizontal flip, else NO Second value : YES if vertical flip, else NO Values derived from the original XA image
SID vector	(0019,xxBE)	FL	1-N	Distance in mm from source to detector center. This is a multi-valued attribute that contains the SID for each frame. Values derived from the original XA image
SOD vector	(0019,xxE9)	FL	1-N	Distance in mm from source to the system iso-center. This is a multi-valued attribute that contains the SOD for each frame Values derived from the original XA image

This page left intentionally blank.

© 2010 General Electric Company.
GE Medical Systems, a General Electric Company, going to market as GE Healthcare.
283, rue de la Minière
78530, Buc
FRANCE

www.gehealthcare.com