

# Drawing Index

These sheets are a document set and should not be separated. Electrical information and references are contained on all sheets.

SITE READINESS	C1
EQUIPMENT LAYOUT	A1
(Equipment locations, heat loads, component weights, environmental specs)	
STRUCTURAL LAYOUT	S1
(Structural support/mounting locations for floor/wall/ceiling, wall support elevations)	
STRUCTURAL DETAILS	S2
(Floor and Ceiling loading information)	
ELECTRICAL LAYOUT	E1
(Contractor supplied wiring, interconnect methods, junction point locations and descriptions)	
ELECTRICAL SPECIFICATIONS	E2
(Maximum wiring run lengths, interconnect diagram, system power specifications)	
ELECTRICAL DETAILS	E3 THRU E4
EQUIPMENT DETAILS	D1 THRU D2

**\* REQUIRED REFERENCE \***

Innova 2000  
Preinstallation Manual  
5160944

A mandatory component of this drawing set is the GE Healthcare Preinstallation manual. Failure to reference the preinstallation manual will result in incomplete documentation required for site design and preparation.

Preinstallation documents for GE Healthcare products can be accessed on the web at:

<http://www.gehealthcare.com/company/docs/siteplanning.html>

# GE Healthcare



## Cardio-Vascular Site Planning



imagination at work

## Customer Site Readiness Requirements

- Any deviation from these drawings must be communicated in writing to and reviewed by your local GE Healthcare Installation Project Manager prior to making changes.
- Make arrangements for any rigging, special handling, or facility modifications that must be made to deliver the equipment to the installation site. If desired, your local GE Healthcare Installation Project Manager can supply a reference list of rigging contractors.
- New construction requires the following; 1. Secure area for equipment, 2. Power for drills and other test equipment, 3. Capability for image analysis, 4. Restrooms.
- Provide for refuse removal and disposal (e.g. crates, cartons, packing)
- Contact a radiation physicist or consultant to specify radiation containment requirements.

## GE Equipment Delivery Requirements

Items 1 through 8 on the GE Healthcare Site Readiness Checklist are REQUIRED to facilitate equipment delivery to the installation site. Equipment will not be delivered if these requirements are not satisfied.

GE Healthcare Site Readiness Checklist						
GEHC Global Order # : _____			Customer: _____			
GEHC On-site Representative : _____			MI Supplier: _____			
Name of customer reviewed with : _____			Lead Installer: _____			
GEHC PMI : _____			Phone Number: _____			
Target Site Prep Completion Date: _____			Helper: _____			
The customer is responsible for proper site preparation and site readiness regardless of any GEHC inspections/assessments.						
Item #	Inspection Date	Storage: Is item ready?	Predict (Pre-ship): Is this item ready? Will item be ready?	Verify (Delivery): Is item ready?	Validate (Mech Install): Is item ready?	Comments If "N", please enter in comments or action plan
1	Equipment installation drawings must match actual room size and must meet clearance requirements. Deviations that meet installation requirements may be red-lined, if red-lining is allowed by local code. Seismic requirements are identified on construction drawings.					
2	Delivery route to installation or storage area meets requirements and has been discussed and scheduled with the customer. Ensure floor protection is discussed, requirements identified, and will be available at time of delivery and installation.					
3	Rooms that will contain equipment, including storage areas, are dust free. Room security to prevent unauthorized access and theft has been discussed with customer. The customer is aware of these security issues, implications and responsibility.					
4	In room HVAC ductwork and units (in room) must be mechanically installed and dust free. Installation rooms appear to meet environmental conditions (see Further Definitions) and observed issues have been communicated to the customer. If being stored, storage area must meet PIM storage criteria.					
5	Ceiling grid is installed, Unistrut is located per the installation drawings, and permanent lighting is installed and operational.					
6	Floor is clean and prepared for final floor covering. Customer has verified floor leveling meets the equipment installation drawings and PIM specs and no visible defects are observed. Gantry and table baseplate are installed prior to delivery (if applicable)					
7	Access to a working phone at the facility for emergency use, including MR magnet delivery.					
8	All walls primed (final coat not needed on Day 1), and counter tops that will support equipment must be installed. No dust-producing cabinetry work in installation areas.					
9	Mechanical supplier has been provided with a set of equipment installation drawings for reference. For California, permitted construction drawings or PMI-specified installation drawings are required.					
10	Conduit/electrical cable ducting/dividers/ access flooring installed, with the exception of surface-mounted floor ducting. Wiring to the main disconnect panel is installed and compliant with equipment installation drawings or pre-installation manual.					

Issued Date: 7/9/07 Rev 11

GE Healthcare Technologies  
Installation Services Design Center  
Milwaukee, Wisconsin

SHEET TITLE: SITE READINESS  
MODALITY TYPE: INNOVA 2000  
THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, THE USER HAS BEEN MADE AWARE OF THE REQUIREMENTS FOR ACTUAL CONSTRUCTION PURPOSES, HOWEVER, AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:  
TYPICAL SPECIAL PROCEDURES 5-80F  
TYPICAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-80F	04
DATE:	10-08-07
DRAWN BY:	LLM
CHECKED BY:	TST

REVISION HISTORY:


SHEET  
C1



SCALE: 1/4" = 1'-0"

EQUIPMENT LAYOUT

RECOMMENDED CEILING HEIGHT = 9'-6"

GE EQUIPMENT LISTING


EQUIPMENT ON ORDER FROM GE HEALTHCARE, INSTALLED BY GE HEALTHCARE, PER : NEITHER A QUOTE OR GON WAS ISSUED AT THE DATE OF THESE DRAWINGS

NOTE: LOCAL CONDITIONS MAY DICTATE THAT ITEMS IDENTIFIED IN THIS CATEGORY BE INSTALLED BY OTHERS.

EQUIPMENT CROSS REFERENCE CHART  
 P = PRE-APPROVAL  
 C = CALCULATIONS/  
 PENDING APPROVAL  
 S = SPECIFICATIONS ONLY

ITEM NO.	QUANTITY ORDERED	REFER TO SHEET "D"	ITEM DESCRIPTION (* = EXISTING/REINSTALL)	WEIGHT	HEAT OUTPUT (PER HOUR)	DETAIL NO.	STRC PLAN	ELEC PLAN
1			OPERATORS CONSOLE	15 lbs	194 btu	A80066	---	WBC1 C
2			ADVANTX VMP 100 GENERATOR	1358 lbs	5119 btu	A0721	S11	MP2 C
3			REMOTE CONTROL FOR INJECTOR (OPTION)	4 lbs		B5028	---	IEC S
4			INJECTOR ELECTRONICS (OPTION)	37 lbs	320 btu	B5028	---	IE S
5			INJECTOR HEAD ON TABLE RAIL (OPTION)	15 lbs		B5030A	---	IH S
6			FOUR LCD MONITOR SUSPENSION ON 9 FT. 6 IN. XT INBOARD BRIDGE	485 lbs	1157 btu	B2004 B2010A	---	WBM1 C
7			OMEGA IV/V TABLE WITH ROTATING TOP	1300 lbs	614 btu	B5037	B5049M	LUS C
8			COUNTERBALANCED EYE AND THYROID SHIELD WITH LAMP (OPTION)	167 lbs	150 btu	B2064	B2064B	LMP S
9			POSITIONER CABINET	663 lbs	11126 btu	B5049E	S101	VPE1 C
10			INNOVA 2000 POSITIONER (REFERENCE TABLE BASE-PLATE DETAIL FOR FLOOR MOUNTING INFORMATION)	1543 lbs	1167 btu	B5049 B5049A B5049B	---	LC1 C
11			DLX DIGITAL CONSOLE AND KEYBOARD	17 lbs	870 btu	C7412J	---	DX2 C
12			DLX DIALOG DISPLAY	15 lbs		C7412G	---	---
13			DLX KEYPAD	4 lbs		C7412H	---	---
14			ATLAS CABINET	947 lbs	9556 btu	B0556	S100	ACAB C
15			DETECTOR CHILLER	33 lbs	706 btu	B5049F	---	DC S
16			WATER CHILLER	229 lbs	6143 btu	M0917A	---	CHLR C
17			LONGITUDINAL STATIONARY RAIL FOR XT SUSPENSION	68 lbs		B2004	---	---
18			JEDI GENERATOR CART	50 lbs		B2004	---	---

THE FOLLOWING ITEMS, WHICH HAVE BEEN ORDERED FROM GE HEALTHCARE, ARE TO BE INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.



ANCILLARY ITEMS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
60	MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 44 IN. W X 89 IN. H (1138mm X 2108mm). CONTINGENT ON A 96 IN. (2438mm) CORRIDOR WIDTH
61	X-RAY ON WARNING LIGHT - AVAILABLE FROM GE SUPPLY CALL 800-200-9760 GE CAT. NO. WX1A8VV-DF-XIU
62	CABLE DRAPE RAIL, CAT. NO. CP655 OR EQUIVALENT. CONTACT UNISTRUT WISCONSIN: 262-796-8710
63	CATHETER CABINETS
64	CONTROL WALL TO CEILING WITH LEAD GLASS VIEWING WINDOW.
65	COUNTER TOP WITH SINK, BASE AND WALL CABINETS
66	COUNTER TOP FOR EQUIPMENT-SHELVING MAY BE REQUIRED. PROVIDE DIMENSIONED OPENINGS AS REQUIRED TO ROUTE INTERCONNECT CABLES TO RACEWAY BELOW COUNTERTOP.
67	SLIDING EQUIPMENT ROOM DOORS
68	SHELF FOR EQUIPMENT

THE FOLLOWING ITEMS ARE AVAILABLE FROM GE HEALTHCARE TECHNOLOGIES. CONTACT YOUR LOCAL GE HEALTHCARE SERVICE REPRESENTATIVE FOR PRICING AND AVAILABILITY.

90	X-RAY ROOM WARNING LIGHT/ROOM LIGHTING CONTROL PANEL REFERENCE JUNCTION POINT 'XRLC' ON SHEET 'E1' FOR DETAILED DESCRIPTION -CAT. NO. E4500SS FOR WARNING LIGHT & ROOM LIGHT CONTROL
91	INNOVA MAIN DISCONNECT REFERENCE JUNCTION POINT 'A' ON SHEET 'E1' FOR DETAILED DESCRIPTION. CAT NO. E4502BK (24 W X 36 H X 11.2 IN. DEEP)

GENERAL SPECIFICATIONS

- o THE REQUIRED CEILING HEIGHT INDICATED ON THESE PLANS IS TO ENSURE EQUIPMENT FUNCTION IS NOT INHIBITED. CONSULT WITH YOUR LOCAL GEHC INSTALLATION SPECIALIST REGARDING ACCEPTABILITY OF OTHER CEILING HEIGHTS.
- o CHECK ALL DOOR OPENINGS AND HALLWAYS FROM DELIVERY LOCATION TO WHERE EQUIPMENT IS TO BE INSTALLED TO ENSURE THE ROUTE PHYSICALLY AND STRUCTURALLY WILL ACCOMMODATE THE EQUIPMENT AS SHIPPED.
- o RADIATION PROTECTION REQUIREMENTS ARE NOT INDICATED ON THIS PLAN. WHERE NEEDED PER NATIONAL OR LOCAL CODE THEY SHALL BE SPECIFIED BY A QUALIFIED RADIOLOGICAL PHYSICIST.
- o THE DEVELOPMENT OF THE EQUIPMENT LAYOUT, ROOM DIMENSIONS, MECHANICAL AND ELECTRICAL SUGGESTIONS IS PREDICATED UPON THE BEST INFORMATION OBTAINABLE FROM THE SITE, COUPLED WITH THE CUSTOMER'S KNOWN DESIRES. ARCHITECTURAL OR ELECTRICAL CHANGES INCLUDING RELOCATION OF EQUIPMENT ILLUSTRATED ON THIS DRAWING IS ALLOWED ONLY WITH NOTIFICATION, IN WRITING, AND REVIEW BY GEHC SERVICE DEPARTMENT. EQUIPMENT OPERATION, SERVICEABILITY, AND RESTRICTING CABLE LENGTHS, ETC., MAKE THIS ESSENTIAL FOR A PROPER INSTALLATION. GEHC RESERVES THE RIGHT TO MAKE ON THE JOB CHANGES BECAUSE OF CUSTOMER REQUIREMENTS AND/OR OBSTACLES IN CONSTRUCTION, ETC..
- o ALL WORK TO BE IN COMPLIANCE WITH NATIONAL AND LOCAL BUILDING SAFETY CODES.
- o DIMENSIONS ARE TO FINISHED SURFACES OF ROOM

SITE ENVIRONMENT SPECIFICATIONS

- o AMBIENT OPERATING TEMPERATURE - 58 TO 75 DEGREES (F), MAXIMUM ALLOWABLE TEMPERATURE CHANGE OF 15 DEGREES (F)/HOUR
- o HUMIDITY: REFER TO PREINSTALLATION MANUAL FOR THE EQUIPMENT ILLUSTRATED ON THIS DRAWING.
- o ALTITUDE: NOT TO EXCEED 8,000 FT. ABOVE SEA LEVEL.
- o THE ENVIRONMENT FOR THE ELECTRONICS CABINET MUST BE CONTROLLED SO THE ABOVE RESTRICTIONS ARE NOT EXCEEDED.
- o DO NOT RESTRICT THE AIR INTAKE AT THE LOWER FRONT OR AIR EXHAUST AT THE TOP OF THE ELECTRONICS CABINETS.

MAGNETIC INTERFERENCE SPECIFICATIONS

IMAGE INTENSIFIERS MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 1 GAUSS TO GUARANTEE SPECIFIED IMAGING PERFORMANCE.

X-RAY TUBES MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO GUARANTEE SPECIFIED PERFORMANCE.

SYSTEM ELECTRONICS MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO GUARANTEE DATA INTEGRITY.

OPERATORS CONSOLE EQUIPMENT MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO OBTAIN SPECIFIED GEOMETRIC LINEARITY.

**GE Healthcare Technologies**  
 Installation Services Design Center  
 Milwaukee, Wisconsin

SHEET TITLE: EQUIPMENT LAYOUT  
 MODALITY TYPE: INNOVA 2000

THIS PLAN IS LIMITED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CORRECTLY IDENTIFY ACTUAL CONSTRUCTION PURPOSES, HOWEVER, AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:  
**TYPICAL SPECIAL PROCEDURES 5-80F**  
 TYPICAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-80F	04
DATE:	10-08-07
DRAWN BY:	LLM
CHECKED BY:	TST

REVISION HISTORY:




SHEET  
**A1**

SCALE: 1/4" = 1'-0"

STRUCTURAL LAYOUT

RECOMMENDED CEILING HEIGHT = 9'-6"

TYPICAL WALL SUPPORT ELEVATIONS



STRUCTURAL SUPPORT METHODS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
1	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S11, FOR GENERATOR CABINET
2	AREA OCCUPIED BY GE SUPPLIED OMEGA TABLE BASE
3	UNISTRUT OR EQUIVALENT SUPPORTS FOR FASTENING THE OVERHEAD COUNTERPOISED SUSPENSION. SUPPORT TO BE LOCATED AS SHOWN. SUPPORT SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, BE PARALLEL SQUARE AND IN THE SAME HORIZONTAL PLANE. FLUSH WITH FINISHED CEILING. SUSPENSION REQUIRES 102 LBS/BOLT SUPPORT. METHODS OF SUPPORT THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION.
4	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S101, FOR POSITIONER CABINET.
5	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S108, FOR MAIN DISCONNECT CONTROL.
6	AREA OCCUPIED BY GE SUPPLIED POSITIONER BASEPLATE
7	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S100, FOR ATLAS CABINET.
8	UNISTRUT OR EQUIVALENT SUPPORT IN CEILING FOR FASTENING CEILING SUPPORTED EQUIPMENT. SUPPORTS TO RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL. RUN WALL TO WALL, BE PARALLEL SQUARE, AND IN THE SAME HORIZONTAL PLANE, FLUSH WITH THE FINISHED CEILING. RAILS ARE MOUNTED TO THESE SUPPORTS EVERY 2'-2" AND REQUIRE 50 LBS. (227 LBS. IN SEISMIC REGIONS) PER BOLT LOAD. METHODS OF SUPPORT THAT PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION.
9	>>COMPONENTS FLUSH WITH CEILING<< UNISTRUT OR EQUIVALENT SUPPORT IN CEILING FOR FASTENING CABLE DRAPE RAIL. SUPPORTS TO RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL. RUN WALL TO WALL, BE PARALLEL SQUARE, AND IN THE SAME HORIZONTAL PLANE, FLUSH WITH THE FINISHED CEILING. RAILS ARE MOUNTED TO THESE SUPPORTS EVERY 2'-2" AND REQUIRE 50 LBS. (227 LBS. IN SEISMIC REGIONS) PER BOLT LOAD. METHODS OF SUPPORT THAT PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION. >>COMPONENTS BELOW CEILING<< CABLE DRAPE RAIL, UNISTRUT CAT. NO. CPG55 or EQUIVALENT.

STRUCTURAL NOTES

- ALL STEEL WORK AND PARTS NECESSARY TO SUPPORT CEILING MOUNTED TUBE HANGER OR OTHER EQUIPMENT ARE TO BE SUPPLIED BY THE CUSTOMER OR HIS CONTRACTORS. THE UNISTRUT OR EQUIVALENT STRUCTURE SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, RUN WALL TO WALL, BE PARALLEL SQUARE AND IN THE SAME HORIZONTAL PLANE FLUSH WITH FINISHED CEILING. THE SYSTEM IS TO BE CROSS BRACED VERTICALLY, HORIZONTALLY AND DIAGONALLY TO ALLOW NO MOVEMENT AND A MAXIMUM OF 1.58mm(1/16") DEFLECTION. (10) 12.7mm (1/2") DIA. x 38.1mm (1 1/2") LONG BOLTS WITH UNISTRUT 12.7mm (1/2") NUTS WITH SPRINGS ARE TO BE PROVIDED BY CUSTOMER OR HIS CONTRACTORS FOR EACH STATIONARY AND AUXILIARY SUPPORT RAIL. CLOSURE STRIPS SHALL BE PROVIDED FOR AREAS OF UNISTRUT EXPOSED AND WITHOUT MOUNTING UNITS.
- METHODS OF SUPPORT FOR THE STEELWORK THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE CONCRETE OR MASONRY ANCHORS IN DIRECT TENSION.
- ALL UNITS THAT ARE WALL MOUNTED OR WALL SUPPORTED ARE TO BE PROVIDED WITH SUPPORTS WHERE NECESSARY. WALL SUPPORTS ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTORS. SEE PLAN AND DETAIL SHEETS FOR SUGGESTED LOCATIONS AND MOUNTING HOLE LOCATIONS.
- ALL CEILING MOUNTED FIXTURES, AIR VENTS, SPRINKLERS, ETC. TO BE FLUSH MOUNTED, OR SHALL NOT EXTEND MORE THAN 6.35mm (1/4") BELOW THE FINISHED CEILING.
- CONTROL WALLS WITH TUBE HANGER PASSAGE ABOVE SHALL BE CONSTRUCTED TO 2130mm (7'-0") HIGH.
- FLOOR SLABS ON WHICH EQUIPMENT IS TO BE INSTALLED MUST BE LEVEL TO 3.17mm (1/8") IN 3050mm (10'-0")
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.
- CUSTOMERS CONTRACTOR MUST PROVIDE ALL PENETRATIONS IN POST TENSION FLOORS.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL ANY NON-STANDARD ANCHORING. DOCUMENTS FOR STANDARD ANCHORING METHODS ARE INCLUDED WITH GE EQUIPMENT DRAWINGS FOR GEOGRAPHIC AREAS THAT REQUIRE SUCH DOCUMENTATION.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL HARDWARE FOR "THROUGH THE FLOOR" ANCHORING AND/OR ANY BRACING UNDER ACCESS FLOORS. THIS CONTRACTOR MUST ALSO PROVIDE FLOOR DRILLING THAT CANNOT BE COMPLETED BECAUSE OF AN OBSTRUCTION ENCOUNTERED WHILE DRILLING BY THE GE INSTALLER SUCH AS REBAR ETC.

**GE Healthcare Technologies**  
 Installation Services Design Center  
 Milwaukee, Wisconsin

SHEET TITLE: STRUCTURAL LAYOUT  
 MODALITY TYPE: INNOVA 2000

THIS PLAN IS LIMITED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, THE ENGINEER HAS BEEN MADE AWARE OF THE ACTUAL CONSTRUCTION PURPOSES, HOWEVER, AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:  
**TYPICAL SPECIAL PROCEDURES 5-80F**  
 TYPICAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-80F	04
DATE:	10-08-07
DRAWN BY:	LLM
CHECKED BY:	TST

REVISION HISTORY:


SHEET  
**S1**




SCALE: 1/4" = 1'-0"

ELECTRICAL PLAN

RECOMMENDED CEILING HEIGHT = 9'-6"

JUNCTION POINT DESCRIPTIONS



**ELECTRICAL OUTLET LEGEND**  
CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS. HEIGHT ABOVE FLOOR DETERMINED BY LOCAL CODES UNLESS OTHERWISE SPECIFIED.




- ⊕ DUPLEX HOSPITAL GRADE, DEDICATED WALL OUTLET 120-V, SINGLE PHASE POWER
- ⊕ DEDICATED TELEPHONE LINE(S) (SEE ELECTRICAL DETAIL ELEC-1 OR ELEC-67)
- ⊕ NETWORK OUTLET (SEE ELECTRICAL DETAILS ELEC-83 AND ELEC-84 OR ELEC-87)
- ⊕ 3-15R NEMA RECEPTACLE, DEDICATED OUTLET 120-V, SINGLE PHASE POWER
- ⊕ DUPLEX HOSPITAL GRADE, DEDICATED OUTLET 120-V EMERGENCY, SINGLE PHASE POWER, 15A

**DUCT HATCHING LEGEND**

- ▨ ABOVE CEILING DUCT
- ▨ UNDER FLOOR DUCT
- ▨ TRENCH DUCT (FLUSH FLOOR)
- ▨ SURFACE FLOOR DUCT
- ABOVE CEILING CONDUIT
- BELOW FLOOR CONDUIT

- JUNCTION POINT NOTES**
- ALL JUNCTION BOXES, CONDUIT, DUCT, DUCT DIVIDERS, SWITCHES, CIRCUIT BREAKERS, ETC., ARE TO BE SUPPLIED AND INSTALLED BY CUSTOMER'S ELECTRICAL CONTRACTOR.
  - CONDUIT AND DUCT RUNS SHALL HAVE SWEEP RADIUS BENDS.
  - CONDUITS AND DUCT ABOVE CEILING OR BELOW FINISHED FLOOR MUST BE INSTALLED AS NEAR TO CEILING OR FLOOR AS POSSIBLE TO REDUCE RUN LENGTH.
  - CEILING MOUNTED JUNCTION BOXES ILLUSTRATED ON THIS PLAN MUST BE INSTALLED FLUSH WITH FINISHED CEILING.
  - ALL DUCTWORK MUST MEET THE FOLLOWING REQUIREMENTS:
    - DUCTWORK SHALL BE METAL WITH DIVIDERS AND HAVE REMOVABLE, ACCESSIBLE COVERS.
    - DUCTWORK SHALL BE CERTIFIED/RATED FOR ELECTRICAL POWER PURPOSES.
    - DUCTWORK SHALL BE ELECTRICALLY AND MECHANICALLY BONDED TOGETHER IN AN APPROVED MANNER.
    - PVC AS A SUBSTITUTE MUST BE USED IN ACCORDANCE WITH ALL LOCAL AND NATIONAL CODES.
  - ALL OPENINGS IN ACCESS FLOORING ARE TO BE CUT OUT AND FINISHED OFF WITH GROMMET MATERIAL BY THE CUSTOMER'S CONTRACTOR.
  - GENERAL CONTRACTOR TO INSERT PULL CORDS FOR ALL CABLE RUN CONDUITS BETWEEN THE EQUIPMENT ROOM AND THE OPERATORS CONTROL ROOM.
  - 10 FOOT PITGALS AT ALL JUNCTION POINTS. NO ALUMINUM OR SOLID WIRES.
  - ALL WIRING MUST BE THHN OR TFFN STRANDED COPPER THERMOPLASTIC 600 VOLT OR EQUIVALENT UNLESS OTHERWISE STATED.
  - GROUNDING IS CRITICAL TO EQUIPMENT FUNCTION AND PATIENT SAFETY. SITE MUST CONFORM TO WIRING SPECIFICATIONS SHOWN ON THIS PLAN.

POINT	DESCRIPTION	QTY.	HARDWARE	DETAIL NO., SHT. E3
A	MAIN DISCONNECT AVAILABLE FROM GE CALL: 800-558-5102 OR LOCAL GE INSTALLATION PROJECT MGR.	1	INNOVA 150-AMP CIRCUIT BREAKER DISTRIBUTION PANEL GEMS CAT. NO. E4502K 150 EMERGENCY OFF PUSHBUTTON GE INSTALLATION PROJECT MGR.	ELEC-103
ACAB	ATLAS CABINET	1	32 IN. OF GROMMET MATERIAL FOR AN 8" X 8" OPENING IN DUCT COVER	ELEC-5 ELEC-6
CHLR	RECIRCULATING WATER CHILLER	1	12 IN. OF GROMMET MATERIAL FOR A 3" X 3" OPENING IN DUCT COVER	ELEC-5 ELEC-6
DC	DETECTOR CHILLER	1	12 IN. OF GROMMET MATERIAL FOR A 3" X 3" OPENING IN DUCT COVER	ELEC-5 ELEC-6
DX2	DLX CONTROL	1	12 IN. OF GROMMET MATERIAL FOR A 3" X 3" OPENING IN DUCT COVER	ELEC-5 ELEC-6
IE	INJECTOR ELECTRONICS	1	12 IN. OF GROMMET MATERIAL FOR A 3" X 3" OPENING IN DUCT COVER	ELEC-5 ELEC-6
IEC	INJECTOR CONTROL	1	12 IN. OF GROMMET MATERIAL FOR A 3" X 3" OPENING IN DUCT COVER	ELEC-5 ELEC-6
IH	INJECTOR HEAD	1	EXTERNALLY CONNECTED AT TABLE BASE	ELEC-9
LC1	INNOVA LC	1	GE SUPPLIED FITTING 12" X 12" X 6" IN. BOX 1 1/2" IN. DIA. LOCKNUT 4" IN. DIA. BUSHING 1/4" X 24" X 12" IN. BOX SUITABLE LENGTH OF 6" IN. DIA. THREADED CONDUIT OR PIPE 6" IN. DIA. LOCKNUTS 1/2" IN. DIA. BUSHING	ELEC-100
LMP	SURGICAL LAMP	1	COVERPLATE 1 1/2" X 4" X 4" IN. BOX 1/2" IN. DIA. CHASE NIPPLE	ELEC-8
LUS	OMEGA TABLE	1	12" X 12" X 6" IN. GROUND BAR WITH 1/4" IN. DIA. MACHINE SCREWS 2" X 8" IN. DIA. BUSHING & LOCKNUT	ELEC-9
MP1	ADVANTX GENERATOR	1	LOCATED WITHIN 'MP2'	ELEC-5 ELEC-6
MP2	POWER UNIT	1	32 IN. OF GROMMET MATERIAL FOR AN 8" X 8" OPENING IN DUCT COVER	ELEC-5 ELEC-6
RDS1	EMERGENCY OFF	1	PROVIDE A SINGLE GANG, 2 1/8" IN. DEEP, FLUSH MTD. WALL BOX.	ELEC-16
RDS2	EMERGENCY OFF	1	PROVIDE A SINGLE GANG, 2 1/8" IN. DEEP, FLUSH MTD. WALL BOX.	ELEC-16
RML1	ROOM LIGHTS AVAILABLE FROM GE. CALL: 800-558-5102	1	COVERPLATE SINGLE GANG BOX #E45055 24V X-RAY ROOM WARNING LIGHT AND ROOM LIGHT CONTROLLER OR EQUIVALENT.	ELEC-17
VPE1	POSITIONER CABINET	1	32 IN. OF GROMMET MATERIAL FOR AN 8" X 8" OPENING IN DUCT COVER	ELEC-5 ELEC-6
WBC1	ADVANTX CONSOLE	1	GE FURNISHED COVERPLATE	ELEC-4
WBM1	TV MONITOR	1	10" X 10" X 6" IN. FLUSH CEILING BOX 2 1/2" IN. DIA. CHASE NIPPLE 1 COVERPLATE	ELEC-8
WC	WATER CHILLER HOSE OUTLET	1	3" IN. CONDUIT STUBBED 2" IN. ABOVE FLOOR	ELEC-9
XRL1	WARNING LIGHT	1	COVERPLATE SINGLE GANG BOX #E45055 24V X-RAY ROOM WARNING LIGHT AND ROOM LIGHT CONTROLLER OR EQUIVALENT. DO NOT USE FLUORESCENT FIXTURES.	ELEC-17
XRLC	WARNING LIGHT CONTROLLER AVAILABLE FROM GE. CALL: 800-558-5102 OR LOCAL GE INSTALLATION PROJECT MGR.	1	#E45055 24V X-RAY ROOM WARNING LIGHT AND ROOM LIGHT CONTROLLER OR EQUIVALENT MAX 24V CONTROLLER	ELEC-17



NOTE: REFER TO CONDUIT LEGEND FOR ALL ADDITIONAL CONDUITS NOT SHOWN ON DRAWING.

**ADDITIONAL CONDUIT RUNS FOR INNOVA 2000 SYSTEMS**

CONDUITS REQUIRED FROM POINT "XRLC"

- XRLC TO XRL1 = ONE 1/2" CND.
- XRLC TO RML1 = ONE 1/2" CND.
- XRLC TO VPE1 = ONE 1/2" CND.
- XRLC TO 120-V 1 PHASE POWER CONDUIT AS REQUIRED

CONTACT YOUR LOCAL RADIO VASCULAR PROJECT MANAGER, INSTALLATIONS (CVPMI) FOR ANY MODIFICATIONS TO ROOM LAYOUT.

BEFORE PROCEEDING WITH INSTALLATION OF CEILING MOUNTED FIXTURES, PLEASE REFER TO STRUCTURAL SHEET S1 FOR LOCATIONS OF UNISTRUT AND OTHER STRUCTURAL SUPPORTED EQUIPMENT IN CEILING.

CONDUITS REQUIRED FROM POINT "WBM1"

- WBM1 TO ACAB = TWO 2 1/2" CND.

CONDUITS REQUIRED FROM POINT "LMP"

- LMP TO 120-V 1 PHASE POWER CONDUIT AS REQUIRED

**FEEDEER TABLE** REV. DATE: 02/26/06

\* CALCULATIONS BASED UPON NOMINAL VOLTAGE, WIRE SIZE IN AWG.  
\* RECOMMENDED FEEDER SIZES FROM DIST. TRNG. TO ROOM DISCONNECT. CALCULATIONS ARE AT NOMINAL VOLTAGE BASED UPON MAXIMUM 1/0 WIRE SIZE FROM ROOM DISCONNECT TO POWER CABINET WITH A MAXIMUM RUN OF 25 FT.  
\* NEUTRAL MUST BE TERMINATED INSIDE THE MAIN DISCONNECT PANEL AND NOT AT ANY GE CABINET.  
\* THE GROUNDING CONDUCTOR WILL BE OF SAME SIZE AS THE FEEDER WIRES. THIS GROUNDING WILL RUN FROM THE EQUIPMENT BACK TO THE FACILITY POWER SOURCE/MAIN GROUNDING POINT AND ALWAYS TRAVEL IN THE SAME CONDUIT WITH THE FEEDERS AND NEUTRAL.  
\* MINIMUM WIRE SIZE FOR CIRCUIT BREAKER, BASED ON RECOMMENDED OVERCURRENT PROTECTION.  
\* FOR A FULL SYSTEM UPS, REFER TO ELECTRICAL DETAILS FOR UPS FEEDER WIRES.

RUN LENGTH IN FEET	POWER SUPPLY VOLTAGE					
	324-396 360	342-418 360	360-440 400	378-462 400	398-484 440	414-506 480
	SIZE OF FEEDERS AND GROUND WIRES (AWG)					
50	1/0	1/0	1/0	1/0	1/0	1/0
100	1/0	1/0	1/0	1/0	1/0	1/0
150	3/0	2/0	2/0	1/0	1/0	1/0
200	4/0	4/0	3/0	3/0	2/0	1/0
250	30M	30M	25M	4/0	3/0	3/0
300	40M	35M	30M	25M	4/0	4/0
350	60M	50M	40M	35M	30M	25M
400	70M	60M	50M	40M	35M	30M

CONDUITS REQUIRED FROM POINT "A"

- A TO RDS1 = ONE 1/2" CND.
- A TO RDS2 = ONE 1/2" CND.
- A TO MP2 = ONE CND. AS REQ'D. FOR TWO CUSTOMER SUPPLIED POWER/GROUND RUNS.
- A TO ACAB = ONE 1" CND. FOR ONE GE SUPPLIED SIGNAL CABLES
- A TO 480-V 3 PHASE POWER CONDUIT AS REQUIRED

CONDUITS REQUIRED FOR MAC LAB

- PC/IVUS TO WBM1 = ONE 3" CND.
- PC/IVUS TO TRAM = ONE 3 1/2" CND.

CONDUITS REQUIRED FOR RADIO LAB

- PC/IVUS TO WBM1 = ONE 3" CND.
- PC TO TRAM = ONE 3" CND.
- IVUS TO TRAM = ONE 3" CND.

CONDUITS REQUIRED FOR COMBO LAB

- PC TO WBM1 = ONE 3" CND.
- PC TO TRAM/AMP = ONE 3" CND.
- PC TO RMOT = ONE 3" CND.
- IVUS TO TRAM/AMP = ONE 3" CND.


CONDUITS REQUIRED FOR GENERIC PHYSIO

- PC/IVUS TO WBM1 = ONE 3" CND.
- PC TO TRAM = ONE 3" CND.
- IVUS TO TRAM = ONE 3" CND.

CONTRACTOR SUPPLIED AND INSTALLED WIRING

ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS.

WIRE RUN, FROM - TO	QUANTITY, WIRE SIZE/COLOR
A > MP2	3-BLACK, 1-WHITE, 1-GREEN - REFER TO FEEDER TABLE
A > RDS2	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
XRLC > 1 PHASE	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
A > RDS1	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
3 PHASE > A	3-BLACK, 1-WHITE, 1-GREEN - REFER TO FEEDER TABLE
XRL1 > XRLC	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
120V > LMP	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
120-V > XRLC	1-BLACK, 1-WHITE, 1-GREEN - (SIZE AS REQUIRED)
XRLC > VPE1	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
RML1 > XRLC	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
A > ACAB	2-ND. 10 BLACK, 1-ND. 10 GREEN



**GE Healthcare Technologies**  
Installation Services Design Center  
Milwaukee, Wisconsin

SHEET TITLE: ELECTRICAL LAYOUT  
MODALITY TYPE: INNOVA 2000  
THIS PLAN IS LIMITED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS DRAWING, THE CONTRACTOR HAS BEEN MADE AWARE OF THE ACTUAL CONSTRUCTION PURPOSES, HOWEVER, AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:  
**TYPICAL SPECIAL PROCEDURES 5-80F**  
TYPICAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-80F	04
DATE:	10-08-07
DRAWN BY:	LLM
CHECKED BY:	TST

REVISION HISTORY:


SHEET  
**E1**





### ELECTRICAL DETAIL REMOVABLE COVERPLATE (TYPICAL)

ELEC-4  
REV. DATE: 10/20/94

6" [152 mm]  
6" [152 mm]  
6" [152 mm]  
6" [152 mm]

GE FURNISHED COVERPLATE

MOUNTED ON DUCT COVER  
MOUNTED ON WALL BOX

DETAIL NOT TO SCALE

### ELECTRICAL DETAIL HORIZONTAL WALL DUCT (TYPICAL)

ELEC-5  
REV. DATE: 03/19/04

TYPICAL WALL DUCT  
REMOVABLE DUCT COVER  
FINISHED FLOOR  
GROMMETED OPENING  
REMOVABLE SECTION OF WALL DUCT COVER  
REFER TO CHART FOR MINIMUM DIVIDER REQUIREMENT LOCAL CODES MAY REQUIRE ADDITIONAL DIVIDERS

ELECTRICAL DUCT  
COVER PLATE TO BE REMOVABLE  
RUBBER GROMMET  
EQUAL  
EQUAL  
DUCT WIDTH  
ELECTRICAL CONTRACTOR TO FURNISH AND INSTALL SCREWS AS SHOWN

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" [610mm]	2
18" [457mm]	2
10" [254mm]	2
6" [152mm]	1
4" [102mm]	1

DETAIL NOT TO SCALE

### ELECTRICAL DETAIL VERTICAL WALL DUCT (TYPICAL)

ELEC-6  
REV. DATE: 03/19/04

REFER TO CHART FOR MINIMUM DIVIDER REQUIREMENT LOCAL CODES MAY REQUIRE ADDITIONAL DIVIDERS  
REMOVABLE SECTION OF WALL DUCT  
REMOVABLE DUCT COVER  
ELECTRICAL DUCT  
DUCT WIDTH  
EQUAL  
EQUAL  
GROMMETED OPENING  
RUBBER GROMMET  
COVER PLATE TO BE REMOVABLE  
ELECTRICAL CONTRACTOR TO FURNISH AND INSTALL SCREWS AS SHOWN

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" [610mm]	2
18" [457mm]	2
10" [254mm]	2
6" [152mm]	1
4" [102mm]	1

DETAIL NOT TO SCALE

### ELECTRICAL DETAIL BOX WITH COVERPLATE (TYPICAL)

ELEC-8  
REV. DATE: 09/30/94

OUTLET BOX  
HARDWARE

DETAIL NOT TO SCALE

### ELECTRICAL DETAIL CONDUITS THRU-FLOOR (TYPICAL)

ELEC-9  
REV. DATE: 08/08/94

FINISHED FLOOR  
HARDWARE  
1.5" (38 mm) TYP.

DETAIL NOT TO SCALE

### ELECTRICAL DETAIL X-RAY WARNING LIGHT & ROOM LIGHT CONTROL PANEL

ELEC-17  
REV. DATE: 05/09/07

TO GE X-RAY ON SIGNAL  
MAXIMUM 24-VAC  
E4502SS  
X-RAY ROOM WARNING LIGHT / ROOM LIGHTING CONTROL PANEL  
X-RAY WARNING LIGHT  
0-VAC  
TO GE ROOM LIGHT SIGNAL  
MAXIMUM 24-VAC  
X-RAY WARNING LIGHT OR ROOM LIGHT ARE NOT PART OF THIS CAT. NO.  
ROOM LIGHTS  
0-VAC

THE R4500AL IS RECOMMENDED IF "X-RAY ON" WARNING LIGHT AND ROOM LIGHT CONTROL ARE UTILIZED  
THE R4500AM IS RECOMMENDED IF "X-RAY ON" WARNING LIGHT ONLY

E4502RL  
X-RAY ROOM WARNING LIGHT CONTROL PANEL  
X-RAY WARNING LIGHT IS NOT PART OF THIS CAT. NO.  
X-RAY WARNING LIGHT  
0-VAC  
TO GE X-RAY ON SIGNAL  
MAXIMUM 24-VAC  
120-VAC 10A MAXIMUM

CONTROL PANEL CAN BE LOCATED ABOVE THE CEILING NEAR THE WARNING LIGHT  
UNLESS SPECIFIED ON SHEET A1 AS BEING INCLUDED ON EQUIPMENT ORDER, ALL ITEMS ILLUSTRATED ARE TO BE FURNISHED AND INSTALLED BY CUSTOMER'S CONTRACTOR

### ELECTRICAL DETAIL EMERGENCY OFF BUTTON

ELEC-16  
REV. DATE: 08/22/05

PLAN VIEW  
FRONT VIEW  
SIDE VIEW  
SINGLE GANG BOX SUPPLIED BY CONTRACTOR  
PLATE & OFF BUTTON

DETAIL NOT TO SCALE

### ELECTRICAL DETAIL INNOVA MAIN DISCONNECT PANEL

ELEC-103  
REV. DATE: 9/24/03

CABINET DEPTH = 10" [254mm]  
System Main Disconnect  
24" [610mm]  
36" [914mm]  
18" [457mm]  
34.5" [876mm]

DETAIL NOT TO SCALE

### ELECTRICAL DETAIL POSITIONER INTERCONNECT DETAIL, UNDER FLOOR

ELEC-100  
REV. DATE: 03/30/04

THRU-FLOOR FITTING  
ISOCENTER  
4.5" [114mm]  
WATER LINES  
POSITIONER BASEPLATE  
TABLE SIDE  
WATER CABLE TROUGH OR CONDUIT  
0.9" [23mm]  
3.6" [91mm]  
5.0" [127mm]  
9" [228mm] DIA. OPENING THRU FLOOR.  
6" [152mm] I.D. PIPE OR CONDUIT  
ELECTRICAL CABLE TROUGH OR CONDUIT  
4.0" [102mm]  
24" x 24" x 12" BOX [610mm x 610mm x 305mm]

NOTE: PIPE, JUNCTION BOX AND DUCT OR CONDUIT ARE TO BE SUPPLIED AND INSTALLED BY CUSTOMER OR CUSTOMER'S CONTRACTOR.

DETAIL NOT TO SCALE

### ELECTRICAL DETAIL INSITE CONNECTION (TYPICAL)

ELEC-1  
REV. DATE: 04/24/02

ONE OF THE FOLLOWING TWO SELECTIONS MUST BE INSTALLED AT THE LOCATION SHOWN ON THE ELECTRICAL PLAN (SHEET E1) FOR GE INSITE CONNECTION BASED UPON SYSTEM CONFIGURATION.  
A) ONE INTERNET ACCESSIBLE VIRTUAL PRIVATE NETWORK (VPN) CONNECTION WITH A STATIC IP ADDRESS, AND ONE TELEPHONE LINE - DEDICATED-DIRECT-DIALING, VOICE GRADE.  
OR  
B) TWO TELEPHONE LINES - ONE DEDICATED DIRECT-DISTANCE-DIALING, VOICE GRADE AND ONE A DEDICATED DATA LINE.

FINISHED CEILING  
1" CONDUIT FROM J.B. TO ABOVE FINISHED CEILING.  
TO BE DETERMINED  
FINISHED FLOOR  
SINGLE GANG J.B.  
COVERPLATE WITH TWO TELEPHONE RECEPTACLES OR ONE TELEPHONE RECEPTACLE AND ONE NETWORK RECEPTACLE

ALL ITEMS ILLUSTRATED ARE TO BE FURNISHED AND INSTALLED BY CUSTOMER OR THEIR CONTRACTOR.

DETAIL NOT TO SCALE

### ELECTRICAL DETAIL BOX WITH COVERPLATE AND NETWORK JACK

ELEC-83  
REV. DATE: 10/06/98

BOX  
NETWORK JACK  
COVERPLATE

DETAIL NOT TO SCALE

### ELECTRICAL DETAIL NETWORK CONNECTION (TYPICAL)

ELEC-84  
REV. DATE: 03/06/04

LOCAL AREA NETWORK  
FINISHED CEILING  
1/2" CONDUIT FROM J.B. TO ABOVE FINISHED CEILING.  
TO BE DETERMINED  
FINISHED FLOOR  
SINGLE GANG J.B.  
COVERPLATE WITH NETWORK RECEPTACLE

FOR NUCLEAR SYSTEMS A DIRECT NETWORK CONNECTION IS TO BE MADE BETWEEN THE SYSTEM AND THE REVIEW WORKSTATION.

DETAIL NOT TO SCALE

**GE Healthcare Technologies**  
Installation Services Design Center  
Milwaukee, Wisconsin

SHEET TITLE: ELECTRICAL DETAILS  
MODALITY TYPE: INNOVA 2000

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, THE ENGINEER HAS BEEN MADE AWARE OF ALL LOCAL AND STATE REGULATIONS AND CODES. THE ENGINEER'S RESPONSIBILITY IS LIMITED TO THE ACTUAL CONSTRUCTION PURPOSES, HOWEVER, AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:  
**TYPICAL SPECIAL PROCEDURES 5-80F**  
TYPICAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-80F	04
DATE: 10-08-07	
DRAWN BY: LLM	
CHECKED BY: TST	


REVISION HISTORY:


SHEET  
**E3**

ELECTRICAL DETAIL  
J.B. / WALL DUCT DETAIL (TYPICAL)

ELEC-2

REV. DATE: 09/30/94



DETAIL NOT TO SCALE

SHEET TITLE: ELECTRICAL DETAILS

MODALITY TYPE: INNOVA 2000

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. PREPARING ENGINEER HAS BEEN MADE TO CORRECT DETAILS TO CONFORM WITH ACTUAL CONSTRUCTION PURPOSES, HOWEVER, AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:

TYPICAL SPECIAL  
PROCEDURES 5-80F  
TYPICAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-80F	04
DATE: 10-08-07	
DRAWN BY: LLM	
CHECKED BY: TST	

REVISION HISTORY:


SHEET

E4



**GE Healthcare Technologies**

Installation Services Design Center  
Milwaukee, Wisconsin



**EQUIPMENT DETAIL**  
CONTROL CONSOLE  
A80066  
REV. DATE: 03/10/99

TOP VIEW  
FRONT VIEW  
REAR VIEW

DRAWING NOT TO SCALE

**EQUIPMENT DETAIL**  
ADVANTX VMP GENERATOR CABINET  
A0721  
REV. DATE: 05/17/05

FRONT VIEW  
SIDE VIEW  
PLAN VIEW

DETAIL NOT TO SCALE

**EQUIPMENT DETAIL**  
INJECTOR REMOTE CONTROL AND ELECTRONICS  
B50-28

REMOTE CONTROL  
ELECTRONICS

DRAWING NOT TO SCALE

**EQUIPMENT DETAIL**  
INJECTOR ON TABLE RAIL  
B50-30A

PLAN VIEW  
SIDE VIEW

DRAWING NOT TO SCALE

**EQUIPMENT DETAIL**  
XT RADIOGRAPHIC SUSPENSION, INBOARD MOUNTING  
B2004  
REV. DATE: 12/07/94

XT STATIONARY RAIL  
CABLE DRAPE RAIL  
XT INBOARD BRIDGE

DETAIL NOT TO SCALE

**EQUIPMENT DETAIL**  
LCD MONITOR SUSPENSION, 4, 6 OR 8 MONITORS  
B2010A  
REV. DATE: 12/16/03

LCD MONITOR SUSPENSION  
XT INBOARD BRIDGE

DETAIL NOT TO SCALE

**EQUIPMENT DETAIL**  
OMEGA IV ANGIO/V LONG VASCULAR TABLE/STEPPING TOP  
B5037  
REV. DATE: 05/18/05

SIDE VIEW  
FRONT VIEW (FOOT END)  
PLAN VIEW

DETAIL NOT TO SCALE

**EQUIPMENT DETAIL**  
AADCO EYE & THYROID SHIELD WITH LAMP  
B20-64  
REV. DATE: 06/08/01

CEILING  
CARRIAGE TRACK  
LEAD GLASS SHIELD

DETAIL NOT TO SCALE

**EQUIPMENT DETAIL**  
INNOVA POSITIONER CABINET  
B5049E  
REV. DATE: 05/17/05

PLAN VIEW TOP / BOTTOM  
FRONT VIEW  
SIDE VIEW

DETAIL NOT TO SCALE

**EQUIPMENT DETAIL**  
INNOVA 2000 MAIN DISCONNECT  
E45-02K  
REV. DATE: 9/24/03

System Main Disconnect

DETAIL NOT TO SCALE

**EQUIPMENT DETAIL**  
INNOVA 2000 VASCULAR SYSTEM  
B5049  
REV. DATE: 04/26/04

PLAN VIEW

DRAWING NOT TO SCALE

**EQUIPMENT DETAIL**  
INNOVA 2000 VASCULAR SYSTEM  
B5049A  
REV. DATE: 04/26/04

FINISHED CEILING  
FINISHED FLOOR  
ISOCENTER  
FOCAL SPOT

DETAIL NOT TO SCALE

**GE Healthcare Technologies**  
Installation Services Design Center  
Milwaukee, Wisconsin

SHEET TITLE: EQUIPMENT DETAILS  
MODALITY TYPE: INNOVA 2000

THIS PLAN IS LIMITED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IT IS NOT TO BE USED FOR CONSTRUCTION OF THE FACILITY. THE USER SHALL BE RESPONSIBLE FOR OBTAINING ALL NECESSARY PERMITS AND APPROVALS FROM LOCAL, STATE AND FEDERAL AGENCIES. GE HEALTHCARE SHALL NOT BE RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:  
**TYPICAL SPECIAL PROCEDURES 5-80F**  
TYPICAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-80F	04
DATE:	10-08-07
DRAWN BY:	LLM
CHECKED BY:	TST

REVISION HISTORY:


SHEET  
**D1**

**EQUIPMENT DETAIL**  
INNOVA 2000 VASCULAR SYSTEM

B5049B  
REV. DATE: 04/26/04

FRONT VIEW

DRAWING NOT TO SCALE

**EQUIPMENT DETAIL**  
DLX CONSOLE AND KEYBOARD

C7412J  
REV. DATE: 12/22/03

DETAIL NOT TO SCALE

**EQUIPMENT DETAIL**  
DLX DIALOG DISPLAY

C7412G  
REV. DATE: 02/22/02

DETAIL NOT TO SCALE

**EQUIPMENT DETAIL**  
DLX or DL KEYPAD

C7412H  
REV. DATE: 09/03/03

DETAIL NOT TO SCALE

**EQUIPMENT DETAIL**  
INNOVA 2000 ATLAS CABINET

B0558  
REV. DATE: 05/17/05

SHIPPING DIMENSIONS:  
25.6"(650mm)D x 35.4"(900mm)W x 70.8"(1800mm)H

PLAN VIEW TOP / BOTTOM

FRONT VIEW

SIDE VIEW

DETAIL NOT TO SCALE

**EQUIPMENT DETAIL**  
INNOVA DETECTOR COOLER

B5049F  
REV. DATE: 8/01/00

PLAN VIEW

FRONT VIEW

SIDE VIEW

DETAIL NOT TO SCALE

**EQUIPMENT DETAIL**  
COOLIX 2200A RECIRCULATING CHILLER

M0917A  
REV. DATE: 05/17/05

SHIPPING DIMENSIONS:  
43.3"(1100mm)D x 32.2"(820mm)W x 49.6"(1260mm)H

PLAN VIEW

SIDE VIEW

FRONT VIEW

REAR VIEW

DETAIL NOT TO SCALE

**GE Healthcare Technologies**  
Installation Services Design Center  
Milwaukee, Wisconsin

SHEET TITLE: EQUIPMENT DETAILS  
MODALITY TYPE: INNOVA 2000

THIS PLAN IS LIMITED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. PREPARING THIS PLAN IS NOT INTENDED TO BE USED AS A CONSTRUCTION DOCUMENT. GE HEALTHCARE TECHNOLOGIES, INC. HAS BEEN MADE AWARE OF THE INTENT OF THIS PLAN AND ACCEPTS NO LIABILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:  
**TYPICAL SPECIAL PROCEDURES 5-80F**  
TYPICAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-80F	04
DATE:	10-08-07
DRAWN BY:	LLM
CHECKED BY:	TST

REVISION HISTORY:


SHEET  
**D2**