

Drawing Index

These sheets are a document set and should not be separated. Electrical information and references are contained on all sheets.

SITE READINESS	C1
EQUIPMENT LAYOUT (Equipment locations, heat loads, component weights, environmental specs)	A1
STRUCTURAL LAYOUT (Structural support/mounting locations for floor/wall/ceiling, wall support elevations)	S1
STRUCTURAL DETAILS (Floor and Ceiling loading information)	S2
ELECTRICAL LAYOUT (Contractor supplied wiring, interconnect methods, junction point locations and descriptions)	E1
ELECTRICAL SPECIFICATIONS (Maximum wiring run lengths, interconnect diagram, system power specifications)	E2
ELECTRICAL DETAILS	E3 THRU E4
EQUIPMENT DETAILS	D1 THRU D3

These drawings indicate the placement and interconnection of the listed equipment components. These drawings are not construction or site preparation drawings. Customer remains ultimately responsible for preparing the site to accommodate the operation of such equipment in compliance with GE Healthcare's written specifications and all applicable federal, state, and/or local requirements.

* REQUIRED REFERENCE *

Innova IGS
Pre Installation Manual
5421046-1-1EN

A mandatory component of this drawing set is the GE Healthcare Pre Installation manual. Failure to reference the preIS manual will result in incomplete documentation required for site design and preparation.

Pre Installation documents for GE Healthcare products can be accessed on the web at:

www.gehealthcare.com/siteplanning

GE Healthcare



Interventional Site Planning

CUSTOMER ACCEPTANCE



imagination at work

Customer Site Readiness Requirements

- Any deviation from these drawings must be communicated in writing to and reviewed by your local GE Healthcare Installation Project Manager prior to making changes.
- Make arrangements for any rigging, special handling, or facility modifications that must be made to deliver the equipment to the installation site. If desired, your local GE Healthcare Installation Project Manager can supply a reference list of rigging contractors.
- New construction requires the following; 1. Secure area for equipment, 2. Power for drills and other test equipment, 3. Capability for image analysis, 4. Restrooms.
- Provide for refuse removal and disposal (e.g. crates, cartons, packing)
- Contact a radiation physicist or consultant to specify radiation containment requirements.

GE Equipment Delivery Requirements

The items on the GE Healthcare Site Readiness Checklist are REQUIRED to facilitate equipment delivery to the IS site. Equipment will not be delivered if these requirements are not satisfied.

GE Healthcare Site Readiness Checklist Rev 19					
Before using this document ensure you have the latest Rev from MyWorkshop on D00422752					
GEHC Global Order #:		Customer:			
GEHC PMI:		FE / Installer:			
The customer is responsible for proper site preparation regardless of any GEHC measurements/inspections/assessments.					
	Inspection Date:	Storage is item ready?	PM is item ready?	FE is item ready?	Comments If 'N', enter comments or action plan
1					MR Magnet Delivery Requirements: Ensure oxygen venting system is available for magnet connection as defined by GEHC Pre-Installation Manual (PIM) requirements, exhaust fan system is installed and operational, 480V power, and chilled water supply is available 24x7 that meets system cooling requirements. External connectivity is available for magnet monitoring and phone service is available during delivery. Surface mount vitromat installed where required. Magnet room final flooring is in place.
2					MR RF Screen Room Requirements: RF Screen Room is tested with copy of Test Report, emailed to ISAdmin.COEMB@ge.com, that it is compliant with GEHC specifications. Dock Bolt and magnet anchors (if applicable) installed using 2 part anchor. For HDx systems, blower box mount bolts installed by RF vendor using 2 part anchors.
3					State Regulatory Requirements: Facility registration number provided for states of Ill, Ky, HI, RI, SC, TX, X-ray shielding plan and state acknowledgment letter provided to installer for AR, DC, NC, SC, CO & WA. Site Drawing Requirements: Final version of equipment network and antenna, installation drawings (including red lined versions) verified to match actual room and has been provided to installer.
4					Surface Penetration Requirements: Customer/Contractor scheduled to provide required drilling or cutting into floors, ceilings, and walls, OR surface penetration permit available and posted in the room when GEHC will perform the work.
5					Pre-Delivery Route Requirements: The equipment delivery route from the truck to the final destination within the facility has been reviewed with all key stakeholders to safely meet the minimum requirements for equipment access, and all communications/notifications have occurred. Arrangements have been made for special handling (elevator, rigging, floor protection, fork lift, rollback truck, etc).
6					Finished Room Requirements: Rooms that will contain equipment, including storage areas not in scan suite, are dust free. Provisions taken to maintain a dust free room. Precautions must be taken to prevent dust from entering rooms containing equipment when construction is incomplete in adjacent areas. All walls primed (final coat not needed on Day 1). Shielding, doors, and windows are to be installed. No contractor work being done during or after the installation that will cause dust in the installation areas or potential equipment damage. Room security to prevent unauthorized access and theft has been discussed with customer. The customer is aware of these security issues, implications and responsibility. For Storage: Room must meet PIM requirements for storage.
7					Electrical Requirements: Lockable (LOTO) Main Disconnect Panel (MDP) is installed per GE guidelines and system power is available. Conduits, electrical cable ducting/dividers/cable trays, and access flooring is installed in proper location and height. Surface floor duct and load-side wires can be installed at time of system installation. Validate outlet location and requirements meet specifications for device/equipment.
8					HVAC Requirements: The HVAC/Chilled Water systems designed to maintain the environment per spec/PMI is at running state and appears to provide the desired environmental conditions including location of vents, temperature and humidity for system operation.
9					Flooring Requirements: Floor is clean and prepared for final floor covering. Floor levelness/flatness is measured and within tolerance, and there are no visible defects per GEHC specifications. Confirm customer anchoring plan aligns with designed floor thickness. Final flooring installed where required for network racks.
10					Ceiling Requirements: Unistrut (or equivalent) location, levelness and spacing is measured (or vendor confirmed) and consistent with the requirement of the installation drawings. Ensure unistrut and rails are not used as mounting surfaces. Ceiling grid is installed. Remount lighting is installed and operational. HVAC diffusers are installed and connected to ductwork. Ceiling tiles installed per PMI discretion.

GE Healthcare
Healthcare Project Implementation - Design Center
Milwaukee, Wisconsin
Copyright 2009 General Electric Company - Proprietary to GE

SHEET TITLE: **SITE READINESS**
MODALITY TYPE: IGS 520, 530, 540
THIS PLAN IS SUBMITTED TO SURVEY LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE ACTUAL CONSTRUCTION PURPOSES AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
INTERVENTIONAL CARDIOLOGY (IC)
TYPICAL FINAL DRAWINGS

PROJECT	REVISION
5-106F	02
DATE:	18.Dec.13
DRAWN BY:	JPH
CHECKED BY:	TST
QT. NO.:	NONE
QT. DT.:	NONE

REVISION HISTORY:

SHEET
C1

PIM R2
RQ - 140198

GE EQUIPMENT LISTING

EQUIPMENT ON ORDER FROM GE HEALTHCARE, INSTALLED BY GE HEALTHCARE, PER QUOTE NONE DATED NONE

NOTE: LOCAL CONDITIONS MAY DICTATE THAT ITEMS IDENTIFIED IN THIS CATEGORY BE INSTALLED BY OTHERS.


ITEM NO.	QUANTITY ORDERED	REFER TO SHEET "D"	ITEM DESCRIPTION (* = EXISTING/REINSTALL)	WEIGHT	HEAT OUTPUT (PER HOUR)	DETAIL NO.	EQUIPMENT CROSS REFERENCE CHART	
							STRC PLAN	ELEC PLAN
1	3		KVA UPS CABINET (LARGE DISPLAY SUBSYSTEM OPTION)	77 lbs	546 btu	B2016	0	UPS3 C
2	1		LARGE DISPLAY MONITOR CABINET	253 lbs	3412 btu	B2014	0	LDC C
3	2		18 IN. MONITOR ON WALL SUPPORT	26 lbs	204 btu	C7617B	-	WBMB C
4	1		IVUS VOLCANO SS CONSOLE, INCLUDES FLAT PANEL MONITOR AND KEYBOARD (DESK MOUNTED)	68 lbs	1631 btu	B551	-	IVUS -
5	1		NURSING NOTES WORKSTATION	46 lbs	682 btu	-	-	S -
6	1		WORKSTATION CART	-	-	-	-	-
7	1		MAC-LAB CONSOLE, INCLUDES MONITORS AND KEYBOARD	566 lbs	2935 btu	-	-	PC S
8	1		TRAM NET RACK	8 lbs	-	B5047	-	TRAM S
9	1		COLOR PRINTER	-	-	-	-	S -
10	1		COUNTERBALANCED EYE AND THYROID SHIELD WITH R96 LAMP	143 lbs	-	B5031E	B5031F	LMP S
11	1		REMOTE CONTROL FOR INJECTOR	4 lbs	-	B5028	-	IEC S
12	1		INJECTOR ELECTRONICS	37 lbs	320 btu	B5028	-	IE S
13	1		INJECTOR HEAD ON TABLE RAIL	15 lbs	-	B5030A	-	IH S
14	2		LONGITUDINAL STATIONARY RAIL FOR XT SUSPENSION	68 lbs	-	B2007B	-	C
15	1		OMEGA IQ TABLE	1750 lbs	614 btu	B5049N	-	LUS C
16	1		XR BUZZER (LOCATED ABOVE CEILING)	2 lbs	-	B5150H	-	XR B
17	1		INNOVA POSITIONER (REFERENCE TABLE BASE-PLATE DETAIL FOR FLOOR MOUNTING INFORMATION)	1653 lbs	2416 btu	B5050A B5050B B5050C B5050D B5050E B5050F B5050G B5050H B5050J B5050P B5050R	-	LC1 C
18	1		UPS INTERFACE BOX	-	-	E4502IB	-	UI B
19	1		ATLAS CABINET (C2)	659 lbs	1825 btu	B0558C	\$100	C2 C
20	1		ATLAS CABINET (C1)	1115 lbs	3389 btu	B0558C	\$100	C1 C
21	1		UPS CABINET	1170 lbs	4061 btu	E4502SC	-	UPS -
22	1		DETECTOR CHILLER	33 lbs	706 btu	B5049F	-	DC S
23	1		WATER 4100 CHILLER	264 lbs	18730 btu	B-1G504 B-1G504	-	CHLR C
24	1		LARGE DISPLAY MONITOR ON SINGLE MONITOR SUSPENSION (9 FT 6 IN. INBOARD BRIDGE (MOUNT TWO GE MONITORS ON BACKSIDE OF LD MONITOR))	784 lbs	1706 btu	B2004 B2015	-	LDM C
25	1		CONTROL ROOM MONITOR WITH DL KEYPAD	22 lbs	204 btu	C7617H C7617	-	S
26	1		OPERATORS CONSOLE	22 lbs	546 btu	C7617 C7502 B5050C	-	WBC1 C
27	1		BOLUS CHASE HANDSWITCH	2 lbs	-	-	-	WBBC -
28	1		AW WORKSTATION	81 lbs	1201 btu	M1013AW A7617	-	C
29	1		COOLIX 4100 AUTOTRANSFORMER	66 lbs	238 btu	B-1G505	-	AT -

THE FOLLOWING ITEMS, WHICH HAVE BEEN ORDERED FROM GE HEALTHCARE, ARE TO BE INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.

30	6		VITALING SPEAKER	-	-	-	-	-
31	1		VITALING CONSOLE	-	-	B0566	-	-
32	1		VITALING MICROPHONE	-	-	-	-	-
33	1		VITALING MICROPHONE (ONE ON MONITOR BRIDGE IN EXAM ROOM)	-	-	-	-	-
34	1		INNOVA MAIN DISCONNECT, REFERENCE JUNCTION POINT PDB ON SHEET E1 FOR DETAILED DESCRIPTION.	326 lbs	1532 btu	E4502M	-	PDB -

EQUIPMENT LAYOUT RECOMMENDED CEILING HEIGHT = 9'-6"

This equipment layout indicates the placement and interconnection of the indicated equipment components. There may be federal, state, and/or local requirements that could impact the placement of these components. It remains the Customer's responsibility for ensuring the site and final equipment placement complies with all applicable federal, state, and/or local requirements.



ANCILLARY ITEMS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
60	BEARING BLOCK OUTLINE. SEE S1 FOR MORE INFORMATION.
61	COUNTER TOP WITH BASE AND WALL CABINETS
62	CATHETER CABINETS
63	COUNTER TOP FOR EQUIPMENT-SHELVING MAY BE REQUIRED PROVIDE GROMMETTED OPENINGS AS REQUIRED TO ROUTE INTERCONNECT CABLES TO RACEWAY BELOW COUNTERTOP.
64	CONTROL WALL TO CEILING WITH LEAD GLASS VIEWING WINDOW
65	SHELF - CUSTOMER TO PROVIDE ADEQUATE WALL SUPPORT
66	X-RAY ON WARNING LIGHT - AVAILABLE FROM GE SUPPLY CALL 800-200-9760 GE CAT. NO. WX1ABWW-DF-XIU
67	CABLE DRAPE RAIL.
68	MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 44 IN. W X 89 IN. H (1118mm X 2108mm), CONTINGENT ON A 96 IN. (2438mm) CORRIDOR WIDTH
69	150-AMP LOCAL SERVICE DISCONNECT FOR LOCK-OUT/TAG-OUT CAPABILITY. (MAY BE A FUSED DISCONNECT, CIRCUIT BREAKER OR SAFETY SWITCH.)
70	CUSTOMER SUPPLIED STORAGE CABINET
71	MED GASES IN CEILING
72	LEAD APRON RACK
73	SCRUB SINK

THE FOLLOWING ITEMS ARE AVAILABLE FROM GE HEALTHCARE TECHNOLOGIES. CONTACT YOUR LOCAL GE HEALTHCARE SERVICE REPRESENTATIVE FOR PRICING AND AVAILABILITY.

90	X-RAY ROOM WARNING LIGHT/ROOM LIGHTING CONTROL PANEL REFERENCE JUNCTION POINT 'XRLC' ON SHEET E1 FOR DETAILED DESCRIPTION -CAT. NO. E4502SS FOR WARNING LIGHT & ROOM LIGHT CONTROL.
----	---

GENERAL SPECIFICATIONS

- THE REQUIRED CEILING HEIGHT INDICATED ON THESE PLANS IS TO ENSURE EQUIPMENT FUNCTION IS NOT INHIBITED. CONSULT WITH YOUR LOCAL GEHC SPECIALIST REGARDING ACCEPTABILITY OF OTHER CEILING HEIGHTS.
- CHECK ALL DOOR OPENINGS AND HALLWAYS FROM DELIVERY LOCATION TO WHERE EQUIPMENT IS TO BE INSTALLED TO ENSURE THE ROUTE PHYSICALLY AND STRUCTURALLY WILL ACCOMMODATE THE EQUIPMENT AS SHIPPED.
- RADIATION PROTECTION REQUIREMENTS ARE NOT INDICATED ON THIS PLAN. WHERE NEEDED PER NATIONAL OR LOCAL CODE THEY SHALL BE SPECIFIED BY A QUALIFIED RADIOLOGICAL PHYSICIST.
- THE DEVELOPMENT OF THE EQUIPMENT LAYOUT, ROOM DIMENSIONS, MECHANICAL AND ELECTRICAL SUGGESTIONS IS PREDICATED UPON THE BEST INFORMATION OBTAINABLE FROM THE SITE, COUPLED WITH THE CUSTOMER'S KNOWN DESIRES. ARCHITECTURAL OR ELECTRICAL CHANGES INCLUDING RELOCATION OF EQUIPMENT ILLUSTRATED ON THIS DRAWING IS ALLOWED ONLY WITH NOTIFICATION, IN WRITING, AND REVIEW BY GEHC SERVICE DEPARTMENT. EQUIPMENT OPERATION, SERVICEABILITY, AND RESTRICTING CABLE LENGTHS, ETC., MAKE THIS ESSENTIAL FOR A PROPER IS. GEHC RESERVES THE RIGHT TO MAKE ON THE JOB CHANGES BECAUSE OF CUSTOMER REQUIREMENTS AND/OR OBSTACLES IN CONSTRUCTION, ETC..
- ALL WORK TO BE IN COMPLIANCE WITH NATIONAL AND LOCAL BUILDING SAFETY CODES.
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM

SITE ENVIRONMENT SPECIFICATIONS

- AMBIENT OPERATING TEMPERATURE: EQUIPMENT ROOM WITH FLUORO UPS OPTION 68° TO 77° F. (20° TO 25° C)
- AMBIENT OPERATING TEMPERATURE: CONTROL ROOM 68° TO 77° F. (20° TO 25° C)
- AMBIENT OPERATING TEMPERATURE: EXAM ROOM-DESIGN FOR PATIENT/OPERATOR COMFORT TARGET TEMPERATURE 64° F (18° C)
- HUMIDITY: 30° TO 75° FOR EQUIPMENT AND CONTROL ROOMS AND 30° TO 70° FOR EXAM ROOM
- ALTITUDE: NOT TO EXCEED 9,842 FT. (3000M) ABOVE SEA LEVEL.
- THE ENVIRONMENT FOR THE ELECTRONICS CABINET MUST BE CONTROLLED SO THE ABOVE RESTRICTIONS ARE NOT EXCEEDED.
- DO NOT RESTRICT THE AIR INTAKE OR AIR EXHAUST OF THE SYSTEM COMPONENTS.
- ENVIRONMENTAL CONDITIONS LISTED ABOVE MUST BE MAINTAINED AT ALL TIMES INCLUDING FOR EXAMPLE OVERNIGHT, WEEKENDS, AND HOLIDAYS.

MAGNETIC INTERFERENCE SPECIFICATIONS

IMAGE INTENSIFIERS MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 1 GAUSS TO GUARANTEE SPECIFIED IMAGING PERFORMANCE.

X-RAY TUBES MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO GUARANTEE SPECIFIED PERFORMANCE.

SYSTEM ELECTRONICS MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO GUARANTEE DATA INTEGRITY.

OPERATORS CONSOLE EQUIPMENT MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO OBTAIN SPECIFIED GEOMETRIC LINEARITY.

GE Healthcare
Healthcare Project Implementation - Design Center
Minneapolis

SHEET TITLE: **EQUIPMENT LAYOUT**
MODALITY TYPE: IGS 520, 530, 540

THIS PLAN IS SUBMITTED TO ASSIST IN THE LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO THE ACTUAL CONSTRUCTION PURPOSES AND THE USER'S INTENT. GE HEALTHCARE ACCEPTS NO LIABILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
INTERVENTIONAL CARDIOLOGY (IC) TYPICAL FINAL DRAWINGS

PROJECT	REVISION
5-106F	02


DATE: 18.Dec.13
DRAWN BY: JPH
CHECKED BY: TST
QT. NO.: NONE
QT. DT.: NONE

REVISION HISTORY:

SHEET
A1


TYPICAL WALL SUPPORT ELEVATIONS

S115




SUPPORT FOR UPS INTERFACE BOX
(NOT TO SCALE)

S107



SUPPORT FOR MAIN DISCONNECT CONTROL
(NOT TO SCALE)

S100




SUPPORT FOR ATLAS/SYSTEMS CABINET
(NOT TO SCALE)

SCALE: 1/4" = 1'-0"

STRUCTURAL LAYOUT

RECOMMENDED CEILING HEIGHT = 9'-6"



STRUCTURAL SUPPORT METHODS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

- | ITEM NO. | ITEM DESCRIPTION
(* INDICATES EXISTING) |
|----------|--|
| 1 | HATCHED AREA INDICATES MONITOR BRIDGE BEARING BLOCK PATH. |
| 2 | UNISTRUT OR EQUIVALENT SUPPORTS FOR FASTENING THE OVERHEAD COUNTERPOISED SUSPENSION. SUPPORTS TO BE LOCATED AS SHOWN. SUPPORT SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL. BE PARALLEL, SQUARE, AND IN THE SAME HORIZONTAL PLANE. FLUSH WITH FINISHED CEILING. SUSPENSION REQUIRES 102 LBS/BOLT SUPPORT. METHODS OF SUPPORT THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION. |
| 3 | UNISTRUT OR EQUIVALENT SUPPORT IN CEILING FOR FASTENING CEILING SUPPORTED EQUIPMENT. SUPPORTS TO RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL. RUN WALL TO WALL, BE PARALLEL, SQUARE, AND IN THE SAME HORIZONTAL PLANE. FLUSH WITH THE FINISHED CEILING. RAILS ARE MOUNTED TO THESE SUPPORTS EVERY 2'-2" AND REQUIRE 350 LBS. (<597 LBS. IN SEISMIC REGIONS) PER BOLT LOAD. METHODS OF SUPPORT THAT PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION. * TO ORDER, CALL UNISTRUT WISCONSIN AT 262-796-8710. |
| 4 | >>COMPONENTS FLUSH WITH CEILING<<
UNISTRUT OR EQUIVALENT SUPPORT IN CEILING FOR FASTENING CABLE DRAPE RAIL. SUPPORTS TO RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL. RUN WALL TO WALL, BE PARALLEL, SQUARE, AND IN THE SAME HORIZONTAL PLANE. FLUSH WITH THE FINISHED CEILING. RAILS ARE MOUNTED TO THESE SUPPORTS EVERY 2'-2" AND REQUIRE 30 LBS. PER BOLT LOAD. METHODS OF SUPPORT THAT PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION. * TO ORDER, CALL UNISTRUT WISCONSIN AT 262-796-8710. |
| 5 | AREA OCCUPIED BY GE SUPPLIED TABLE BASEPLATE ABOVE CEILING |
| 6 | MOUNT XR BUZZER BRACKET ON WALL. |
| 7 | AREA OCCUPIED BY GE SUPPLIED POSITIONER BASEPLATE |
| 8 | SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S115 FOR UPS INTERFACE BOX. |
| 9 | SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S107, FOR MAIN DISCONNECT CONTROL. |
| 10 | SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S100, FOR ATLAS CABINET. |

STRUCTURAL NOTES

- ALL STEEL WORK AND PARTS NECESSARY TO SUPPORT CEILING MOUNTED TUBE HANGER OR OTHER EQUIPMENT ARE TO BE SUPPLIED BY THE CUSTOMER OR HIS CONTRACTORS. THE UNISTRUT OR EQUIVALENT STRUCTURE SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, RUN WALL TO WALL, BE PARALLEL, SQUARE AND IN THE SAME HORIZONTAL PLANE FLUSH WITH FINISHED CEILING. THE SYSTEM IS TO BE CROSS BRACED VERTICALLY, HORIZONTALLY AND DIAGONALLY TO ALLOW NO MOVEMENT AND A MAXIMUM OF 1,58mm(1/16") DEFLECTION.
(10) 12,7mm (1/2") DIA. x 38,1mm (1 1/2") LONG BOLTS WITH UNISTRUT 12,7mm (1/2") NUTS WITH SPRINGS ARE TO BE PROVIDED BY CUSTOMER OR HIS CONTRACTORS FOR EACH STATIONARY AND AUXILIARY SUPPORT RAIL. CLOSURE STRIPS SHALL BE PROVIDED FOR AREAS OF UNISTRUT EXPOSED AND WITHOUT MOUNTING UNITS.
- METHODS OF SUPPORT FOR THE STEELWORK THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE CONCRETE OR MASONRY ANCHORS IN DIRECT TENSION.
- ALL UNITS THAT ARE WALL MOUNTED OR WALL SUPPORTED ARE TO BE PROVIDED WITH SUPPORTS WHERE NECESSARY. WALL SUPPORTS ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTORS. SEE PLAN AND DETAIL SHEETS FOR SUGGESTED LOCATIONS AND MOUNTING HOLE LOCATIONS.
- ALL CEILING MOUNTED FIXTURES, AIR VENTS, SPRINKLERS, ETC. TO BE FLUSH MOUNTED, OR SHALL NOT EXTEND MORE THAN 6,35mm (1/4") BELOW THE FINISHED CEILING.
- CONTROL WALLS WITH TUBE HANGER PASSAGE ABOVE SHALL BE CONSTRUCTED TO 2130mm (7'-0") HIGH.
- FLOOR SLABS ON WHICH EQUIPMENT IS TO BE INSTALLED MUST BE LEVEL TO 3,17mm (1/8") IN 3050mm (10'-0")
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.
- CUSTOMERS CONTRACTOR MUST PROVIDE ALL PENETRATIONS IN POST TENSION FLOORS.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL ANY NON-STANDARD ANCHORING. DOCUMENTS FOR STANDARD ANCHORING METHODS ARE INCLUDED WITH GE EQUIPMENT DRAWINGS FOR GEOGRAPHIC AREAS THAT REQUIRE SUCH DOCUMENTATION.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL HARDWARE FOR "THROUGH THE FLOOR" ANCHORING AND/OR ANY BRACING UNDER ACCESS FLOORS. THIS CONTRACTOR MUST ALSO PROVIDE FLOOR DRILLING THAT CANNOT BE COMPLETED BECAUSE OF AN OBSTRUCTION ENCOUNTERED WHILE DRILLING BY THE GE INSTALLER SUCH AS REBAR ETC.
- IT IS THE CUSTOMER'S RESPONSIBILITY TO PERFORM ANY FLOOR OR WALL PENETRATIONS THAT MAY BE REQUIRED. THE CUSTOMER IS ALSO RESPONSIBLE FOR ENSURING THAT NO SUBSURFACE UTILITIES (E.G., ELECTRICAL OR ANY OTHER FORM OF WIRING, CONDUITS, PIPING, DUCT WORK OR STRUCTURAL SUPPORTS (I.E. POST TENSION CABLES OR REBAR)) WILL INTERFERE OR COME IN CONTACT WITH SUBSURFACE PENETRATION OPERATIONS (E.G. DRILLING AND INSTALLATION OF ANCHORS/SCREWS) PERFORMED DURING THE INSTALLATION PROCESS. TO ENSURE WORKER SAFETY, GE INSTALLERS WILL PERFORM SURFACE PENETRATION OPERATIONS ONLY AFTER THE CUSTOMER'S VALIDATION AND COMPLETION OF THE "GE SURFACE PENETRATION PERMIT"

SHEET TITLE: STRUCTURAL LAYOUT

MODALITY TYPE: IGS 520, 530, 540

THIS PLAN IS SUBMITTED TO SUPPORT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE MOST RECENT GE DRAWINGS AND TO THE COMPANY'S ACTUAL CONSTRUCTION PRACTICES AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.




PROJECT TITLE:
**INTERVENTIONAL
CARDIOLOGY (IC)**
TYPICAL FINAL DRAWINGS

PROJECT TITLE:

PROJECT	REVISION
5-106F	02
DATE:	18.Dec.13
DRAWN BY:	JPH
CHECKED BY:	TST
QT. NO:	NONE
QT. DT:	NONE


REVISION HISTORY:

SHEET
S1



FLOOR MOUNTING : INNOVA 2100-3100-4100 (UNITY)/OMEGA V LONG TABLE (WITH IQ TILT TABLE BASEPLATE) INSTALLATION (TEMPLATE NO. 2360133) **B5049N**




REV. DATE: 06/04/09



WARNING!! THE RELATIONSHIP BETWEEN THE TABLE BASE AND THE POSITIONER BASEPLATE IS CRITICAL.

PRIOR TO DRILLING MOUNTING HOLES CONTACT LOCAL GE HEALTHCARE INSTALLATION PROJECT MANAGER OR LEAD FIELD ENGINEER TO VERIFY THAT THE PROPER FULL SIZE FLOOR MOUNTING TEMPLATE IS USED.

MEDICAL GAS FLOOR EXIT LOCATIONS



Customer/Contractor Alert: It is the responsibility of the Customer or their Contractor to drill all anchor/thru-bolting holes for anchoring the positioner and table to the floor. Refer to GEHC document no. *2290880-2-100 for installation preparation and procedures.

NOTE: THRU BOLTING IS HIGHLY PREFERRED FOR THE INSTALLATION OF THE POSITIONER BASEPLATE AND OMEGA TABLE. HARDENED BOLTS AND 4" x 4" [102mm x 102mm] STEEL PLATES TO BE USED ARE SUPPLIED BY GE HEALTHCARE AS INDICATED ON THE ACTUAL DETAIL DRAWING. BE ADVISED, HOWEVER, THAT ADDITIONAL SUPPORT STRUCTURES: STEEL BEAMS, PLATES, CORE BORING OF MOUNTING HOLES, ETC., ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.

NOTE: IF THRU BOLTING IS NOT POSSIBLE, FLOOR ANCHORS CAN BE USED IF APPROVED BY CUSTOMERS STRUCTURAL ENGINEER. FOR ON GRADE INSTALLATIONS, MOUNTING KIT CAT. NO. 2286398 SHOULD BE ORDERED. ANCHORS INCLUDED IN KIT SHOULD BE APPROVED BY CUSTOMERS STRUCTURAL ENGINEER.

NOTE: BASEPLATES MUST BE LEVEL WITHIN 1/32" [0.79mm]

NOTE: JOISTS MUST BE SPANNED WITH STEEL REINFORCING. SIZE AND THICKNESS OF STEEL REINFORCING ARE DETERMINED BY THE ACTUAL PAN CONSTRUCTION ON SITE. STEEL PLATES, CHANNELS OR BEAMS MAY BE USED.

NOTE: DETERMINE THE POSITION OF THE "REBARS" IN THE CONCRETE FLOOR SO ANCHOR HOLES WILL NOT RUN INTO THEM.

* DOCUMENT FURNISHED BY GEHC AS PART OF INSTALLATION KIT

POSITIONER BOLT FORCES FOR WORST CASE CONDITIONS

LOADS	BOLT TENSION (AT BOLT "A")
HORIZONTAL ACCELERATION = 625 lbs. [284 Kg]	MAXIMUM TENSION = 881 lbs. [400 Kg]
VERTICAL ACCELERATION = 209 lbs. [95 Kg]	BOLT SHEAR (U-ARM LOCKED)
	MAXIMUM SHEAR = 129 lbs. [54 Kg]/BOLT

OMEGA TABLE BOLT FORCES FOR WORST CASE CONDITIONS

LOADS	BOLT TENSION	BOLT SHEAR
	MAXIMUM TENSION = 1938 lbs. [880 Kg]/BOLT	MAXIMUM SHEAR = 407 lbs. [185 Kg]/BOLT

GE Healthcare
Healthcare Project Implementation - Design Center
Manufacture

SHEET TITLE: STRUCTURAL DETAILS
MODALITY TYPE: IGS 520, 530, 540

THIS PLAN IS SUBMITTED TO SUPPORT LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO DETAILS AND DIMENSIONS SHOWN ON THE DRAWING. GE HEALTHCARE SHALL NOT BE RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
INTERVENTIONAL CARDIOLOGY (IC)
TYPICAL FINAL DRAWINGS

PROJECT	REVISION
5-106F	02


DATE: 18.Dec.13
DRAWN BY: JPH
CHECKED BY: TST
QT. NO.: NONE
QT. DT.: NONE

REVISION HISTORY:

SHEET
S2

ELECTRICAL DETAIL
VERTICAL WALL DUCT (TYPICAL)


ELEC-6
REV. DATE: 03/19/04



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
J.B. / WALL DUCT DETAIL (TYPICAL)


ELEC-2
REV. DATE: 09/30/94



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
INNOVA PLUS WITH BOX AND CONDUIT RUNS

ELEC-177
REV. DATE: 11.Nov.11




DETAIL NOT TO SCALE

ELECTRICAL DETAIL
INSITE CONNECTION (TYPICAL)

ELEC-1
REV. DATE: 04/24/02

ONE OF THE FOLLOWING TWO SELECTIONS MUST BE INSTALLED AT THE LOCATION SHOWN ON THE ELECTRICAL PLAN (SHEET E1) FOR GE INSITE CONNECTION BASED UPON SYSTEM CONFIGURATION.
A) ONE INTERNET ACCESSIBLE VIRTUAL PRIVATE NETWORK (VPN) CONNECTION WITH A STATIC IP ADDRESS, AND ONE TELEPHONE LINE - DEDICATED-DIRECT-DIALING, VOICE GRADE.
OR
B) TWO TELEPHONE LINES - ONE DEDICATED DIRECT-DISTANCE-DIALING, VOICE GRADE AND ONE A DEDICATED DATA LINE.




ALL ITEMS ILLUSTRATED ARE TO BE FURNISHED AND INSTALLED BY CUSTOMER OR THEIR CONTRACTOR.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL
FLOOR BOX WITH NIPPLES (TYPICAL)


ELEC-13
REV. DATE: 09/30/94



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
BOX WITH COVERPLATE AND NETWORK JACK


ELEC-83
REV. DATE: 10/06/98



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
NETWORK CONNECTION (TYPICAL)


ELEC-84
REV. DATE: 03/06/04



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
BOX WITH COVERPLATE (TYPICAL)


ELEC-8
REV. DATE: 09/30/94



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
X-RAY WARNING LIGHT & ROOM LIGHT CONTROL PANEL


ELEC-157
REV. DATE: 04/23/09



CONTROL PANEL CAN BE LOCATED ABOVE THE CEILING NEAR THE WARNING LIGHT
UNLESS SPECIFIED ON SHEET A1 AS BEING INCLUDED ON EQUIPMENT ORDER,
ALL ITEMS ILLUSTRATED ARE TO BE FURNISHED AND INSTALLED BY CUSTOMER'S CONTRACTOR

ELECTRICAL DETAIL
EMERGENCY OFF BUTTON


ELEC-16
REV. DATE: 05/14/09



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
TABLE INTERCONNECTION - BOX BELOW FLOOR

ELEC-48
REV. DATE: 01/04/96



DETAIL NOT TO SCALE

SHEET TITLE: ELECTRICAL DETAILS
MODALITY TYPE: IGS 520, 530, 540

PROJECT TITLE:
INTERVENTIONAL
CARDIOLOGY (IC)
TYPICAL FINAL DRAWINGS

PROJECT	REVISION
5-106F	02
DATE:	18.Dec.13
DRAWN BY:	JPH
CHECKED BY:	TST
QT. NO:	NONE
QT. DT:	NONE

REVISION HISTORY:

SHEET
E3


PIM R2
RQ - 140198

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED

GE Healthcare
Healthcare Project Implementation - Design Center
Madison, WI


ELECTRICAL DETAIL
TABLE INTERCONNECT DETAIL, UNDER FLOOR

ELEC-134
REV. DATE: 05/10/04




ELECTRICAL DETAIL
POSITIONER INTERCONNECT DETAIL, UNDER FLOOR

ELEC-100
REV. DATE: 03/30/04




ELECTRICAL DETAIL
HORIZONTAL WALL DUCT (TYPICAL)

ELEC-5
REV. DATE: 03/19/04




ELECTRICAL DETAIL
HORIZONTAL WALL DUCT (TYPICAL)

ELEC-5A
REV. DATE: 06/16/08




ELECTRICAL DETAIL
CONDUITS THRU-FLOOR (TYPICAL)

ELEC-9
REV. DATE: 08/08/94



ELECTRICAL DETAIL
INNOVA PLUS MAIN DISCONNECT PANEL

ELEC-161
REV. DATE: 09/27/10



GE Healthcare
Healthcare Project Implementation - Design Center
Minneapolis, MN

SHEET TITLE: ELECTRICAL DETAILS
MODALITY TYPE: IGS 520, 530, 540

THIS PLAN IS SUBMITTED TO SUPPORT THE LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL LOCAL, STATE AND FEDERAL REGULATIONS AND STANDARDS. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.













PROJECT TITLE:
INTERVENTIONAL
CARDIOLOGY (IC)
TYPICAL FINAL DRAWINGS

PROJECT	REVISION
5-106F	02

DATE: 18.Dec.13
DRAWN BY: JPH
CHECKED BY: TST
QT. NO.: NONE
QT. DT.: NONE

REVISION HISTORY:

SHEET
E4



GE Healthcare

Healthcare Project Implementation - Design Center
Manufacture

SHEET TITLE: EQUIPMENT DETAILS

MODALITY TYPE: IGS 520, 530, 540

THIS PLAN IS SUBMITTED TO REQUEST FOR APPROVAL OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPLIANCE, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL APPLICABLE CODES AND STANDARDS. GE HEALTHCARE DOES NOT ACCEPT LIABILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:

INTERVENTIONAL CARDIOLOGY (IC)

TYPICAL FINAL DRAWINGS

PROJECT	REVISION
5-106F	02

DATE: 18.Dec.13
DRAWN BY: JPH
CHECKED BY: TST
QT. NO: NONE
QT. DT: NONE

REVISION HISTORY:

SHEET

D1

PIM R2
RQ - 140198

EQUIPMENT DETAIL
INNOVA VASCULAR SYSTEM

B5050B
REV. DATE: 06/07/05

FRONT VIEW

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
INNOVA VASCULAR SYSTEM

B5050
REV. DATE: 02/22/05

PLAN VIEW

SHIPPING DIMENSIONS:
110"L x 45.5"W x 77"H
(2790mmL x 1160mmW x 1950mmH)
(ON DOLLY)
WIDTH IS REDUCED TO 34" [865mm]
BY REMOVING SIDE RAILS

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
SHIPPING CREATE FOR OMEGA & INNOVA IQ TABLES

B5050E
REV. DATE: 12/07/09

OMEGA IV/V TABLE
OMEGA TABLE BASE ASSEMBLY ON PALLET
INNOVA IQ TILT TABLE
OMEGA 4/5 AND INNOVA IQ TILT TABLE
INNOVA IQ TILT TABLE BASE ASSEMBLY ON PALLET

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
INNOVA C2 CABINET ON SHIPPING DOLLY

B5050F
REV. DATE: 12/07/09

INNOVA C2 CABINET (ON DOLLY)

SHIPPING WEIGHT: 821 LBS. (372 kg)(C2-HARMONY SYSTEMS)
SHIPPING WEIGHT: 785 LBS. (356 kg)(C2-BIPLANE SYSTEMS)

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
SHIPPING DOLLY FOR INNOVA LC POSITIONER

B5050G
REV. DATE: 12/07/09

INNOVA LC POSITIONER (ON DOLLY)

NOTE:
BOTH ENDS OF THE DOLLY CAN BE REMOVED WHICH WILL SHORTEN LC GANTRY DOLLY DONE TO 86.22" (2190mm) RECOMMEND ONLY ONE SIDE BE REMOVED WHEN DELIVERY THROUGH HOSPITAL.

SHIPPING WEIGHT: 2340 lbs. (1060 kg)

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
INNOVA C1 & C3 CABINETS ON SHIPPING DOLLY

B5050H
REV. DATE: 12/07/09

INNOVA C1 CABINETS (ON DOLLY)

SHIPPING WEIGHT: 1277 lbs (579 kg) C1-HARMONY SYSTEMS)
SHIPPING WEIGHT: 1052 lbs (477 kg) C1 FRONTAL - BIPLANE SYSTEMS)
SHIPPING WEIGHT: 866 lbs (393 kg) C3 LATERAL - BIPLANE SYSTEMS)

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
LP AND LC GANTRY DELIVERY PATH

B5050J
REV. DATE: 12/07/09

LP GANTRY
LC GANTRY

CORRIDOR WIDTH (FOR GANTRY DELIVERY)	DOOR SIZE (FOR GANTRY DELIVERY)
8'-0" WIDE	4'-0" OPENING
7'-0" WIDE	5'-0" OPENING
6'-0" WIDE	6'-0" OPENING
5'-0" WIDE	7'-0" OPENING
4'-0" WIDE	8'-0" OPENING

NOTE:
WHEN DELIVERING GANTRY FROM CORRIDOR TO CORRIDOR THE SAME HOLDS TRUE AS CORRIDOR TO DOOR SIZE.

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
INNOVA RADIATION SCATTER PLOTS

B5050P
REV. DATE: 12/07/09

UNITS: RELATIVE AIR KERMA: $\mu\text{Gy}/\mu\text{Gym}^2$
DISTANCES: RADIUS AT 1, 2 AND 3 METERS

NOTE:
FOR REFERENCE ONLY. PLEASE REFER TO OPERATOR'S MANUAL FOR ADDITIONAL INFORMATION.

LATERAL 1 METER
LATERAL 1.5 METERS

4-4-4 GANTRY IN LATERAL POSITION - DOSE AT 1.5 METER FROM GROUND
4-4-3 GANTRY IN LATERAL POSITION - DOSE AT 1 METER FROM GROUND

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
INNOVA RADIATION SCATTER PLOTS

B5050R
REV. DATE: 12/07/09

UNITS: RELATIVE AIR KERMA: $\mu\text{Gy}/\mu\text{Gym}^2$
DISTANCES: RADIUS AT 1, 2 AND 3 METERS

NOTE:
FOR REFERENCE ONLY. PLEASE REFER TO OPERATOR'S MANUAL FOR ADDITIONAL INFORMATION.

VERTICAL 1 METER
VERTICAL 1.5 METERS

4-4-2 GANTRY IN VERTICAL POSITION - DOSE AT 1.5 METER FROM GROUND
4-4-3 GANTRY IN LATERAL POSITION - DOSE AT 1 METER FROM GROUND

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
18" FLAT PANEL MONITOR ON WALL SUPPORT

C76-17B
REV. DATE: 04/29/04

FRONT VIEW
SIDE VIEW

PLAN VIEW

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
IVUS VOLCANO S5i WORKSTATION

BS5I
REV. DATE: 04/04/08

IVUS SAFETY ISOLATION TRANSFORMER
IVUS PATIENT INTERFACE MODULE
IVUS CPU
IVUS CONTROL CONSOLE
IVUS JOYSTICK
IVUS VIDEO SWITCH
IVUS PRINTER
IVUS CONTROL ROOM MONITOR

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
TRAM-RAC 4A

B5047
REV. DATE: 05/26/04

TABLE RAIL MOUNT
FLOOR MOUNT

DRAWING NOT TO SCALE

GE Healthcare
Healthcare Project Implementation - Design Center
Manufacture, Wisconsin

SHEET TITLE: EQUIPMENT DETAILS
MODALITY TYPE: IGS 520, 530, 540

THIS PLAN IS SUBMITTED TO SUPPORT THE PURCHASE OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN. EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST REVISED SPECIFICATIONS AND TO THE LATEST GE HEALTHCARE DRAWINGS. GE HEALTHCARE ACCEPTS NO LIABILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
INTERVENTIONAL CARDIOLOGY (IC)
TYPICAL FINAL DRAWINGS

PROJECT	REVISION
5-106F	02

DATE: 18.Dec.13
DRAWN BY: JPH
CHECKED BY: TST
QT. NO: NONE
QT. DT: NONE

REVISION HISTORY:

SHEET
D2

EQUIPMENT DETAIL
VITALINQ COMMUNICATION AND MUSIC SYSTEM

B0566
REV. DATE: 06/14/05

PROCEDURE ROOM

SEE MANUFACTURER'S PREINSTALL MANUAL FOR PRODUCT PLACEMENT AND INSTALLATION.

COMMUNICATION MICROPHONE OVER TABLE END
LEFT MUSIC SPEAKER OVER TABLE END
RIGHT MUSIC SPEAKER OVER TABLE END
COMMUNICATION SPEAKER OVER DOCTOR
COMMUNICATION SPEAKER OVER OPERATOR
LEFT MUSIC SPEAKER ABOVE AND IN FRONT OF OPERATOR
CONSOLE

CONTROL ROOM

VIS-A-VIS, INC.
1-800-319-6014

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
MAVIG EYE & THYROID SHIELD WITH LAMP

B50-31E
REV. 00: 10/03/97

CEILING
CARRIAGE TRACK
65" MAX. ARC [1651mm]
65" MAX. ARC [1651mm]
LEAD GLASS SHIELD

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
INJECTOR REMOTE CONTROL AND ELECTRONICS

B50-28

REMOTE CONTROL
7.5" [191mm]
14.5" [368mm]
3.25" [83mm]
ELECTRONICS
10.5" [267mm]
17" [432mm]
12.75" [324mm]

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
INJECTOR ON TABLE RAIL

B50-30A

PLAN VIEW
SIDE VIEW

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
XR-BUZZER BRACKET

B5150H
REV. 00: 10/30/08

62" [165mm]
31" [85mm]
5 Ø5
3.12" [80mm]
2.5" [7mm]
3.7" [10mm]
3" [75mm]
19" [5mm]
2 Ø15
6.89" [170mm]
7.28" [180mm]
1.9" [5mm]
1.00mm [100mm]
75" [20mm]
1" [25mm]
5.50" [140mm]
5.6" [15mm]
12" [3mm]
6" [100mm]
12" [3mm]

NOTE: XR-BUZZER BRACKET IS MOUNTED ON WALL, ABOVE CEILING. PLACE SPEAKER ABOVE GRILLED CEILING TILE FOR SOUND PENETRATION.

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
UPS INTERFACE BOX

E45021B
REV. DATE: 07/11/05

BOTTOM VIEW
PLAN VIEW
FRONT VIEW
FRONT VIEW (COVER OFF)

(4) MOUNTING HOLES
4.3" [111mm]
3.1" [8mm]
11.3" [288mm]
9.0" [230mm]
9.0" [230mm]
11.0" [280mm]

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
DLX or DL KEYPAD

C7412H
REV. DATE: 09/03/03

3.23" [82mm]
9.37" [238mm]
11.82" [300mm]

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
18" FLAT PANEL MONITOR

C76-17
REV. DATE: 08/28/09

FRONT VIEW
SIDE VIEW
PLAN VIEW

16.5" [420mm]
18" [458mm]
13.2" [337mm]
3.7" [94mm]

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
RCIM WITH DL KEYBOARD CONSOLE

C75-02
REV. DATE: 10/25/10

6.5" [165mm]
17.35" [440.7mm]
17.32" [440mm]
8.3" [210mm]
17.9" [455mm]
20.7" [525mm]
16.5" [420mm]
18" [458mm]
13.2" [337mm]
3.7" [94mm]
18.1" [460mm]
6.5" [165mm]

PC TOWER XW6400
PC TOWER XW8200
FRONT VIEW
SIDE VIEW
PLAN VIEW
KEYBOARD

DETAIL NOT TO SCALE

TYPICAL CONTROL ROOM
INNOVA SINGLE PLANE

B5050C
REV. DATE: 08/26/08

LIVE MONITOR
REFERENCE MONITOR
ADVANTAGE WINDOWS WORKSTATION
INNOVA CONSOLE
DL CONTROL MONITOR
IVUS VOLCANO WORKSTATION
PHYSIO MONITORING WORKSTATION

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
WORKSTATION

M1013AW
REV. DATE: 04/25/01

6.93" [176mm]
15.75" [400mm]
17.1" [435mm]
18.7" [475mm]
19.5" [495mm]
18.5" [471mm]
18.0" [457mm]
7.25" [184mm]
2.25" [57mm]
17.6" [447mm]
18.1" [460mm]
8.6" [218mm]

ULTRASPAC STATION
KEYBOARD
21" COLOR MONITOR
FLAT PANEL MONITOR

DETAIL NOT TO SCALE

GE Healthcare
Healthcare Project Implementation - Design Center
Madison, Wisconsin

SHEET TITLE: EQUIPMENT DETAILS
MODALITY TYPE: ICS 520, 530, 540

THIS PLAN IS SUBMITTED TO REQUEST LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO THE DETAILS AND DIMENSIONS SHOWN. HOWEVER, THE USER SHALL VERIFY THE ACTUAL CONDITIONS OF THE PROJECT AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
INTERVENTIONAL CARDIOLOGY (IC)
TYPICAL FINAL DRAWINGS

PROJECT	REVISION
5-106F	02
DATE:	18.Dec.13
DRAWN BY:	JPH
CHECKED BY:	TST
QT. NO.:	NONE
QT. DT.:	NONE

REVISION HISTORY:

SHEET
D3

PIM R2
RQ - 140198