


LOGIQ P7 & P9

Physical Dimensions and Weight (mm) :

Please Note: The picture may vary according to order configuration.



External Inputs and Outputs :




Environmental Requirements

Operating Temperature	10 to 35 C (50 to 96F)
Operating Humidity	30 to 80 % rH, non-condensing , Pressure : 700 -1060 hPa
Heat Dissipation	680 BTU/hour +300 Btu for each person. Without lights, other equipment etc.
Storage Temperature	-10 to 50 C (-14 to 122 F),
Storage Humidity	10 % - 90 % rH, non-condensing , Pressure : 700 -1060 hPa
Lighting	Combination lighting system (dim/bright) is recommended. Keep in mind that lighting controls and dimmers can be a source of EMI which could degrade image quality.

Please Note: For more information please refer to the product's service & user manuals.

Recommended Floor Plan 4.3m x 5.2m (14ft x 17ft)



Electrical Specifications

Voltage	100 - 120 VAC ±10% 5.0 ... 4.1 A 50, 60 Hz (±2%) 220 - 240 VAC ±10% 2.2 ... 2.0 A 50, 60 Hz (±2%)
Power Consumption	500VA
Site Power Outlets	A dedicated AC power outlet must be within reach of the unit without extension cords
Inrush current	Inrush current is not a factor to consider due to the inrush current limiting properties of the power supplies.
Power outage may occur	Requires a dedicated single branch circuit. To avoid circuit overload and possible loss of critical care equipment, make sure you DO NOT have any other equipment operating on the same circuit.
Acoustic Noise output	Less than 50 dB(A)