

Drawing Index

These sheets are a document set and should not be separated. Electrical information and references are contained on all sheets.

SITE READINESS

C1

EQUIPMENT LAYOUT

A1

(Equipment locations, heat loads, component weights, environmental specs)

STRUCTURAL LAYOUT

S1

(Structural support/mounting locations for floor/wall/ceiling, wall support elevations)

STRUCTURAL DETAILS

S2

(Floor and Ceiling loading information)

ELECTRICAL LAYOUT

E1

(Contractor supplied wiring, interconnect methods, junction point locations and descriptions)

ELECTRICAL SPECIFICATIONS

E2

(Maximum wiring run lengths, interconnect diagram, system power specifications)

ELECTRICAL DETAILS

E3

EQUIPMENT DETAILS

D1 THRU D2

These equipment IS drawings indicate the placement and interconnection of the listed equipment components. These drawings are not construction or site preparation drawings. Customer remains ultimately responsible for preparing the site to accommodate the IS and operation of such equipment in compliance with GE Healthcare's written specifications and all applicable federal, state, and/or local requirements.

\* REQUIRED REFERENCE \*

INJECTORS CT & PET  
Pre Installation Manual  
N/A

A mandatory component of this drawing set is the GE Healthcare Pre Installation manual. Failure to reference the preIS manual will result in incomplete documentation required for site design and preparation.

Pre Installation documents for GE Healthcare products can be accessed on the web at:

www.gehealthcare.com/siteplanning

GE Healthcare



CT Site Planning




imagination at work

Customer Site Readiness Requirements

- Any deviation from these drawings must be communicated in writing to and reviewed by your local GE Healthcare Installation Project Manager prior to making changes.
- Make arrangements for any rigging, special handling, or facility modifications that must be made to deliver the equipment to the installation site. If desired, your local GE Healthcare Installation Project Manager can supply a reference list of rigging contractors.
- New construction requires the following; 1. Secure area for equipment, 2. Power for drills and other test equipment, 3. Capability for image analysis, 4. Restrooms.
- Provide for refuse removal and disposal (e.g. crates, cartons, packing)
- Contact a radiation physicist or consultant to specify radiation containment requirements.

GE Equipment Delivery Requirements

The items on the GE Healthcare Site Readiness Checklist are REQUIRED to facilitate equipment delivery to the IS site. Equipment will not be delivered if these requirements are not satisfied.



GE Healthcare Site Readiness Checklist Rev 19

Before using this document ensure you have the latest Rev from MyWorkshop on DOC0422752

GEHC Global Order #

Customer:

GEHC PMI

FE / Installer:

The customer is responsible for proper site preparation regardless of any GEHC measurements/inspections/assessments.

		Inspection Date:					Comments
			Storage is item ready?	PMI is item ready?	FE is item ready?		If "N", enter comments or action plan
	<div>GEHC Minimum Requirements</div>						
1	<div>MR Magnet Delivery Requirements:</div> <div>Ensure cryogen venting system is available for magnet connection as defined by GEHC Pre-Installation Manual (PIM) requirements, exhaust fan system is installed and operational, 480V power, and chilled water supply is available 24x7 that meets system cooling requirements. External connectivity is available for magnet monitoring and phone service is available during delivery. Surface mount vibromat installed where required. Magnet room final flooring is in place.</div>						
2	<div>MR RF Screen Room Requirements:</div> <div>RF Screen Room is tested with copy of Test Report, emailed to 3dskin@GE-Healthcare.com, that it is compliant with GEHC specifications. Dock Bolt and magnet anchors (if applicable) installed using 2 part anchor. For HDx systems, blower box mount bolts installed by RF vendor using 2 part anchors</div>						
3	<div>State Regulatory Requirements:</div> <div>Facility registration number provided for states of IL, KY, HI, RI, SC, TX, X-ray shielding plan and state acknowledgment letter provided to installer for AR, DC, NC, SC, CO &amp; WA.</div>						
4	<div>Site Drawing Requirements:</div> <div>Final version of equipment network and antenna, installation drawings (including red lined versions) verified to match actual room and has been provided to installer.</div>						
5	<div>Surface Penetration Requirements:</div> <div>Customer/Contractor scheduled to provide required drilling or cutting into floors, ceilings, and walls, OR surface penetration permit available and posted in the room when GEHC will perform the work.</div>						
6	<div>Pre-Delivery Route Requirements:</div> <div>The equipment delivery route from the truck to the final destination within the facility has been reviewed with all key stakeholders to safely meet the minimum requirements for equipment access, and all communications/notifications have occurred. Arrangements have been made for special handling (elevator, rigging, floor protection, fork lift, rollback truck, etc)</div>						
7	<div>Finished Room Requirements:</div> <div>Rooms that will contain equipment, including storage areas not in scan suite, are dust free. Provisions taken to maintain a dust free room. Precautions must be taken to prevent dust from entering rooms containing equipment when construction is incomplete in adjacent areas. All walls primed (final coat not needed on Day 1). Shielding, doors, and windows are to be installed. No contractor work being done during or after the installation that will cause dust in the installation areas or potential equipment damage. Room security to prevent unauthorized access and theft has been discussed with customer. The customer is aware of these security issues, implications and responsibility. For Storage: Room must meet PIM requirements for storage.</div>						
8	<div>Electrical Requirements:</div> <div>Lockable (LOTO) Main Disconnect Panel (MDPI) is installed per GE guidelines and system power is available. Conduits, electrical cable ducting/dividers/cable trays, and access flooring is installed in proper location and height. Surface floor duct and load-side wires can be installed at time of system installation. Validate outlet location and requirements meet specifications for device/equipment.</div>						
9	<div>HVAC Requirements:</div> <div>The HVAC/Chilled Water systems designed to maintain the environment per spec/PIM is at running state and appears to provide the desired environmental conditions including location of vents, temperature and humidity for system operation.</div>						
10	<div>Flooring Requirements:</div> <div>Floor is clean and prepared for final floor covering. Floor levelness/flatness is measured and within tolerance, and there are no visible defects per GEHC specifications. Confirm customer anchoring plan aligns with designed floor thickness. Final flooring installed where required for network racks.</div>						
11	<div>Ceiling Requirements:</div> <div>Unistrut (or equivalent) location, levelness and spacing is measured (or vendor confirmed) and consistent with the requirement of the installation drawings. Ensure unistrut and rails are not used as mounting surfaces. Ceiling grid is installed. Permanent lighting is installed and operational. HVAC diffusers are installed and connected to ductwork. Ceiling tiles installed per PMI discretion.</div>						
12	<div>Staging Requirements:</div> <div>Space has been identified to support the active installation process only. This area meets PIM/project book requirements. Storage space has been identified, if needed. This secured space would be used to store equipment indefinitely. If offsite, transportation plan has been developed at customer expense. This space must meet PIM requirements.</div>						
13	<div>Network Connectivity:</div> <div>Hardware for network connectivity/network drop is in place prior to delivery with specified network firewall configuration where required. Site Surveys for wireless mobile XR units have been completed.</div>						
14	<div>Medical Gases Requirements:</div> <div>Systems (hard piped or portable) in place to allow testing and calibration of equipment (anesthesia), including ventilation.</div>						

GE Healthcare

IS Services Design Center

Milwaukee, Wisconsin

Copyright 2009 General Electric Company – Proprietary to GE

SHEET TITLE: SITE READINESS

MODALITY TYPE: INJECTORS CT & PET

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS, IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST GEHC PIM. CUSTOMER TO VERIFY ALL DETAILS TO BE USED FOR ACTUAL CONSTRUCTION PURPOSES. HOWEVER, THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:

INJECTOR  
TYPICAL FINAL  
CEILING & PEDESTAL MOUNT

PROJECT	REVISION
TYPINJ	03
DATE: 01.Jul.15	
DRAWN BY: DJP	
CHECKED BY: TMS	

REVISION HISTORY:

SHEET

C1





GE EQUIPMENT LISTING									
EQUIPMENT ON ORDER FROM GE HEALTHCARE, INSTALLED BY GE HEALTHCARE, PER QUOTE CT & CT/PET DATED						EQUIPMENT CROSS REFERENCE CHART			
NOTE: LOCAL CONDITIONS MAY DICTATE THAT ITEMS IDENTIFIED IN THIS CATEGORY BE INSTALLED BY OTHERS.						SEISMIC STATUS	P = PREAPPROVAL C = CALCULATIONS/ PENDING APPROVAL S = SPECIFICATIONS ONLY		
ITEM NO.	QUANTITY ORDERED	REFER TO SHEET "D"				DETAIL NO.	STRC PLAN	ELEC PLAN	
⬢	ITEM DESCRIPTION (* = EXISTING/REINSTALL)	WEIGHT	HEAT OUTPUT (PER HOUR)						
①	NEMOTO								
1	MAIN UNIT	13 lbs			B8007C	-	IE	-	
②	INJECTOR CONTROL CONSOLE	8 lbs			B8007C	-	ICC	S	
1	INJECTOR HEAD ON PEDESTAL	26 lbs			B8007D	-	IH	-	
④	INJECTOR HEAD ON OVERHEAD COUNTERPOISED SUSPENSION.	22 lbs			B5031N	B5031B	IH	S	
⑤	MEDRAD STELLANT								
1	INJECTOR CONTROL AND ELECTRONICS	22 lbs	320 btu		E8007N	.	IE	S	
⑥	INJECTOR HEAD ON PEDESTAL	44 lbs			B8007A	.	IH	-	
⑦	INJECTOR HEAD ON OVERHEAD COUNTERPOISED SUSPENSION.	79 lbs			B5031	B5031I	IH	S	
⑧	MEDRAD PROVIS								
1	REMOTE CONTROL FOR INJECTOR	4 lbs			B5028	.	IEC	S	
⑨	INJECTOR ELECTRONICS	37 lbs	320 btu		B5028	.	IE	-	
⑩	INJECTOR HEAD ON OVERHEAD COUNTERPOISED SUSPENSION.	92 lbs			B5031W B5031X	B5031S	IH	S	
⑪	MEDRAD MARK IV/PLUS ECT								
1	REMOTE CONTROL FOR INJECTOR	4 lbs			B5028	.	IEC	S	
⑫	INJECTOR ELECTRONICS	37 lbs	320 btu		B5028	.	IE	S	
⑬	MEDRAD MARK V INJECTOR ON PEDESTAL	90 lbs	320 btu		B5030	.	IH	-	
⑭	INJECTOR HEAD ON OVERHEAD COUNTERPOISED SUSPENSION.	79 lbs			B5031	B5031S	IH	S	
⑮	MEDRAD ENVISION								
1	DISPLAY CONTROL UNIT	8 lbs			B8007	.	ICC	-	
⑯	INJECTOR ELECTRONICS	13 lbs	320 btu		B8007	.	IE	S	
⑰	INJECTOR HEAD ON PEDESTAL	44 lbs			B8007A	.	IH	-	
⑱	INJECTOR HEAD ON OVERHEAD COUNTERPOISED SUSPENSION.	79 lbs			B5031	B5031S	IH	S	
⑲	LF 9000								
1	REMOTE CONTROL FOR INJECTOR	6 lbs			B5029C	---	IEC	S	
⑳	INJECTOR ELECTRONICS	13 lbs	320 btu		B5029C	.	IE	S	
㉑	INJECTOR HEAD ON REMOTE STAND	44 lbs			---	---	IH	-	
㉒	INJECTOR HEAD ON OVERHEAD COUNTERPOISED SUSPENSION.	79 lbs			B5031	B5031S	IH	S	
㉓	EZ-EM PERCUPUMP								
1	REMOTE CONTROL FOR INJECTOR	4 lbs			.	.	ICC	S	
㉔	INJECTOR HEAD ON PEDESTAL	39 lbs			.	.	IH	S	
㉕	INJECTOR HEAD ON OVERHEAD COUNTERPOISED SUSPENSION.	79 lbs			B5031U	B5031S	IH	S	
㉖	EZ-EM EMPOWER								
1	EMPOWER INJECTOR POWER SUPPLY	6 lbs			B8009B	---	---	S	
㉗	EMPOWER INJECTOR REMOTE CONTROL	15 lbs			B8009C	---	ICC	S	
㉘	EMPOWER INJECTOR HEAD ON PEDESTAL	37 lbs			B8009A	---	IH	S	
㉙	EMPOWER INJECTOR HEAD ON OVERHEAD COUNTERPOISED SUSPENSION.	74 lbs			B8009	---	IH	S	

THE FOLLOWING ITEMS, WHICH HAVE BEEN ORDERED FROM GE HEALTHCARE, ARE TO BE INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.

--	--	--	--	--	--	--	--	--	--

SCALE: 1/4" = 1'-0"      EQUIPMENT LAYOUT      RECOMMENDED CEILING HEIGHT = 9'-0"

This equipment layout indicates the placement and interconnection of the indicated equipment components. There may be federal, state, and/or local requirements that could impact the placement of these components. It remains the Customer's responsibility for ensuring the site and final equipment placement complies with all applicable federal, state, and/or local requirements.



ANCILLARY ITEMS	
CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS	
ITEM NO. ⬢	ITEM DESCRIPTION (* INDICATES EXISTING)
THE FOLLOWING ITEMS ARE AVAILABLE FROM GE HEALTHCARE TECHNOLOGIES. CONTACT YOUR LOCAL GE HEALTHCARE SERVICE REPRESENTATIVE FOR PRICING AND AVAILABILITY.	

GENERAL SPECIFICATIONS	
<ul style="list-style-type: none"><li>THE REQUIRED CEILING HEIGHT INDICATED ON THESE PLANS IS TO ENSURE EQUIPMENT FUNCTION IS NOT INHIBITED. CONSULT WITH YOUR LOCAL GEHC SPECIALIST REGARDING ACCEPTABILITY OF OTHER CEILING HEIGHTS.</li><li>CHECK ALL DOOR OPENINGS AND HALLWAYS FROM DELIVERY LOCATION TO WHERE EQUIPMENT IS TO BE INSTALLED TO ENSURE THE ROUTE PHYSICALLY AND STRUCTURALLY WILL ACCOMMODATE THE EQUIPMENT AS SHIPPED.</li><li>RADIATION PROTECTION REQUIREMENTS ARE NOT INDICATED ON THIS PLAN. WHERE NEEDED PER NATIONAL OR LOCAL CODE THEY SHALL BE SPECIFIED BY A QUALIFIED RADIOLOGICAL PHYSICIST.</li><li>THE DEVELOPMENT OF THE EQUIPMENT LAYOUT, ROOM DIMENSIONS, MECHANICAL AND ELECTRICAL SUGGESTIONS IS PREDICATED UPON THE BEST INFORMATION OBTAINABLE FROM THE SITE, COUPLED WITH THE CUSTOMER'S KNOWN DESIRES. ARCHITECTURAL OR ELECTRICAL CHANGES INCLUDING RELOCATION OF EQUIPMENT ILLUSTRATED ON THIS DRAWING IS ALLOWED ONLY WITH NOTIFICATION, IN WRITING, AND REVIEW BY GEHC SERVICE DEPARTMENT. EQUIPMENT OPERATION, SERVICEABILITY, AND RESTRICTING CABLE LENGTHS, ETC., MAKE THIS ESSENTIAL FOR A PROPER IS. GEHC RESERVES THE RIGHT TO MAKE ON THE JOB CHANGES BECAUSE OF CUSTOMER REQUIREMENTS AND/OR OBSTACLES IN CONSTRUCTION, ETC..</li><li>ALL WORK TO BE IN COMPLIANCE WITH NATIONAL AND LOCAL BUILDING SAFETY CODES.</li><li>DIMENSIONS ARE TO FINISHED SURFACES OF ROOM</li></ul>	

SITE ENVIRONMENT SPECIFICATIONS	

MAGNETIC INTERFERENCE SPECIFICATIONS	

GE Healthcare

IS Services Design Center

Minwaukee, Wisconsin

EQUIPMENT LAYOUT

MODALITY TYPE: INJECTORS CT & PET

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS, IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO THE ACTUAL CONSTRUCTION. HOWEVER, THE USER SHALL BE RESPONSIBLE FOR ACCEPTING RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:

INJECTOR  
TYPICAL FINAL  
CEILING & PEDESTAL MOUNT

PROJECT	REVISION
TYPINJ	03

DATE:	01.Jul.15
DRAWN BY:	DJP
CHECKED BY:	TMS

REVISION HISTORY:

SHEET

A1

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED

TYPICAL WALL SUPPORT ELEVATIONS

SCALE: 1/4" = 1'-0"









STRUCTURAL LAYOUT

RECOMMENDED CEILING HEIGHT = 9'-0"

STRUCTURAL SUPPORT METHODS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
<div><div></div></div>	
1	THE PEDESTAL-CEILING MOUNT REQUIRES A FLUSH CEILING MOUNTING PLATE THAT IS STRUCTURALLY SUPPORTED TO HANDLE THE WEIGHT OF THE LOAD AS SHOWN IN DETAIL B50-318. A MAVID DESIGNED CEILING PLATE IS AVAILABLE AND CAN BE ORDERED - B5007P. THIS 14" x 14" PLATE HAS PRE DRILLED 1/2" HOLES IN EACH OF THE 4 MOUNTING CORNERS. IF AN EQUIVALENT PLATE IS USED THE SUPPLIED TEMPLATE SHOULD BE USED TO DRILL THE REQUIRED 4 - 1/2" HOLES IN A PATTERN AS SHOWN IN DETAIL B50-318. AN ADDITIONAL 1/32" HOLE IS REQUIRED FOR THE PEDESTAL-CEILING MOUNT SAFETY CHAIN.
2	UNISTRUT OR EQUIVALENT SUPPORTS FOR FASTENING THE OVERHEAD COUNTERPOISED SUSPENSION. SUPPORT TO BE LOCATED AS SHOWN. SUPPORT SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, BE PARALLEL SQUARE, AND IN THE SAME HORIZONTAL PLANE ABOVE FINISHED CEILING. ENSURE MOUNTING SURFACE IS INSTALLED LEVEL OR PLUMB WITHIN +/- 1 DEGREE, AND IS STRUCTURALLY SUFFICIENT TO MAINTAIN A LEVEL OR PLUMB CONDITION UNDER 110 LB (50KG) SYSTEM LOAD AND MAXIMUM SYSTEM MOMENT OF 4400 IN-LB (500N-M). METHODS OF SUPPORT THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION. 14" x 14" x 3/8" THICK STEEL PLATE PROVIDED BY MANUFACTURER. SEE DETAIL B50311 ON DETAIL SHEETS.
3	UNISTRUT OR EQUIVALENT SUPPORTS FOR FASTENING THE OVERHEAD COUNTERPOISED SUSPENSION. SUPPORT TO BE LOCATED AS SHOWN. SUPPORT SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, BE PARALLEL SQUARE, AND IN THE SAME HORIZONTAL PLANE ABOVE FINISHED CEILING. BODM REQUIRES 102 LBS/BOLT SUPPORT. METHODS OF SUPPORT THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION. 14" x 14" x 3/8" THICK STEEL PLATE PROVIDED BY THE CUSTOMER OR HIS CONTRACTOR. SEE DETAIL B50315 ON DETAIL SHEETS.
4	UNISTRUT OR EQUIVALENT SUPPORTS FOR FASTENING THE OVERHEAD COUNTERPOISED SUSPENSION. SUPPORT TO BE LOCATED AS SHOWN. SUPPORT SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, BE PARALLEL SQUARE, AND IN THE SAME HORIZONTAL PLANE ABOVE FINISHED CEILING. BODM REQUIRES 102 LBS/BOLT SUPPORT. METHODS OF SUPPORT THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION. 14" x 14" x 3/8" THICK STEEL PLATE PROVIDED BY THE CUSTOMER OR HIS CONTRACTOR. SEE DETAIL B50311 ON DETAIL SHEETS.

 <p>EZ-EM EMPOWER</p>	 <p>EZ-EM PERCUPUMP</p>	 <p>LF 9000</p>	 <p>MEDRAD ENVISION</p>
 <p>MEDRAD MARK IV/PLUS ECT</p>	 <p>MEDRAD PROVIS</p>	 <p>MEDRAD STELLANT</p>	 <p>NEMOTO</p>

STRUCTURAL NOTES

- ALL STEEL WORK AND PARTS NECESSARY TO SUPPORT CEILING MOUNTED EQUIPMENT IS TO BE SUPPLIED BY THE CUSTOMER OR HIS CONTRACTORS.
- METHODS OF SUPPORT FOR THE STEELWORK THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE CONCRETE OR MASONRY ANCHORS IN DIRECT TENSION.
- ALL UNITS THAT ARE WALL MOUNTED OR WALL SUPPORTED ARE TO BE PROVIDED WITH SUPPORTS WHERE NECESSARY. WALL SUPPORTS ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTORS. SEE PLAN AND DETAIL SHEETS FOR SUGGESTED LOCATIONS AND MOUNTING HOLE LOCATIONS.
- ALL CEILING MOUNTED FIXTURES, AIR VENTS, SPRINKLERS, ETC. TO BE FLUSH MOUNTED, OR SHALL NOT EXTEND MORE THAN 1/4" BELOW THE FINISHED CEILING.
- FLOOR SLABS ON WHICH EQUIPMENT IS TO BE INSTALLED MUST BE LEVEL TO 1/4" in 10'-0"
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.
- CUSTOMERS CONTRACTOR MUST PROVIDE ALL PENETRATIONS IN POST TENSION FLOORS.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL ANY NON-STANDARD ANCHORING. DOCUMENTS FOR STANDARD ANCHORING METHODS ARE INCLUDED WITH GE EQUIPMENT DRAWINGS FOR GEOGRAPHIC AREAS THAT REQUIRE SUCH DOCUMENTATION.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL HARDWARE FOR "THROUGH THE FLOOR" ANCHORING AND/OR ANY BRACING UNDER ACCESS FLOORS. THIS CONTRACTOR MUST ALSO PROVIDE FLOOR DRILLING THAT CANNOT BE COMPLETED BECAUSE OF AN OBSTRUCTION ENCOUNTERED WHILE DRILLING BY THE GE INSTALLER SUCH AS REBAR ETC.
- IT IS THE CUSTOMER'S RESPONSIBILITY TO PERFORM ANY FLOOR OR WALL PENETRATIONS THAT MAY BE REQUIRED. THE CUSTOMER IS ALSO RESPONSIBLE FOR ENSURING THAT NO SUBSURFACE UTILITIES (E.G., ELECTRICAL OR ANY OTHER FORM OF WIRING, CONDUITS, PIPING, DUCT WORK OR STRUCTURAL SUPPORTS (I.E. POST TENSION CABLES OR REBAR)) WILL INTERFERE OR COME IN CONTACT WITH SUBSURFACE PENETRATION OPERATIONS (E.G. DRILLING AND INSTALLATION OF ANCHORS/SCREWS) PERFORMED DURING THE INSTALLATION PROCESS. TO ENSURE WORKER SAFETY, GE INSTALLERS WILL PERFORM SURFACE PENETRATION OPERATIONS ONLY AFTER THE CUSTOMER'S VALIDATION AND COMPLETION OF THE "GE SURFACE PENETRATION PERMIT"

SHEET TITLE: STRUCTURAL LAYOUT

MODALITY TYPE: INJECTORS CT & PET

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO DETAILS IN THE GE EQUIPMENT DRAWINGS. HOWEVER, THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:

INJECTOR  
TYPICAL FINAL  
CEILING & PEDESTAL MOUNT

PROJECT	REVISION
TYPINJ	03
DATE:	01.Jul.15
DRAWN BY:	DJP
CHECKED BY:	TMS

REVISION HISTORY:

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SHEET

S1

GE Healthcare



IS Services Design Center


Minneapolis, Minnesota

Wisconsin

SUPPORT DETAIL  
OVERHEAD COUNTERPOISED SUSPENSION

B50-31B


REV. DATE: 09.Jul.13



SUPPORT DETAIL  
OVERHEAD COUNTERPOISED SUSPENSION

B50-31I


REV. DATE: 01.Nov.12



SUPPORT DETAIL  
INJECTOR & LEAD GLASS SHIELD  
ON OVERHEAD COUNTERPOISED SUSPENSION

B50-31S


REV. DATE: 03/05/07



SUPPORT DETAIL  
OVERHEAD COUNTERPOISED SUSPENSION

B50-31T

REV. DATE: 03/05/07



SHEET TITLE: STRUCTURAL DETAILS

MODALITY TYPE: INJECTORS CT & PET

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS, IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST EDITIONS OF THE NATIONAL ELECTRICAL CODE (NEC) AND THE NATIONAL FIRE PROTECTION ASSOCIATION (NFPA) CODES. HOWEVER, THE USER OR ACTING CONTRACTOR ASSUMES ALL RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:

INJECTOR  
TYPICAL FINAL  
CEILING & PEDESTAL MOUNT

PROJECT	REVISION
TYPINJ	03
DATE:	01.Jul.15
DRAWN BY:	DJP
CHECKED BY:	TMS

REVISION HISTORY:

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SHEET

S2








GE Healthcare


IS Services Design Center


Minwaukee,

Wisconsin

RECOMMENDED CEILING HEIGHT = 9'-0"



















 DUPLEX HOSPITAL GRADE, DEDICATED WALL OUTLET  
 120-V, SINGLE PHASE POWER

 DUPLEX HOSPITAL GRADE, DEDICATED CEILING OUTLET  
 120-V, SINGLE PHASE POWER

PLEASE SEE BELOW FOR ADDITIONAL REQUIRED  
CONDUIT RUNS AND SIZES.

- o ALL JUNCTION BOXES, CONDUIT, DUCT, DUCT DIVIDERS, SWITCHES, CIRCUIT BREAKERS, ETC., ARE TO BE SUPPLIED AND INSTALLED BY CUSTOMERS ELECTRICAL CONTRACTOR.
- o CONDUIT AND DUCT RUNS SHALL HAVE SWEEP RADIUS BENDS
- o CONDUITS AND DUCT ABOVE CEILING OR BELOW FINISHED FLOOR MUST BE INSTALLED AS NEAR TO CEILING OR FLOOR AS POSSIBLE TO REDUCE RUN LENGTH.
- o CEILING MOUNTED JUNCTION BOXES ILLUSTRATED ON THIS PLAN MUST BE INSTALLED FLUSH WITH FINISHED CEILING.
- o ALL DUCTWORK MUST MEET THE FOLLOWING REQUIREMENTS:
  - 1. DUCTWORK SHALL BE METAL WITH DIVIDERS AND HAVE REMOVABLE, ACCESSIBLE COVERS.
  - 2. DUCTWORK SHALL BE CERTIFIED/RATED FOR ELECTRICAL POWER PURPOSES.
  - 3. DUCTWORK SHALL BE ELECTRICALLY AND MECHANICALLY BONDED TOGETHER IN AN APPROVED MANNER.
  - 4. PVC AS A SUBSTITUTE MUST BE USED IN ACCORDANCE WITH ALL LOCAL AND NATIONAL CODES.
- o ALL OPENINGS IN ACCESS FLOORING ARE TO BE CUT OUT AND FINISHED OFF WITH GROMMET MATERIAL BY THE CUSTOMERS CONTRACTOR.
- o GENERAL CONTRACTOR TO INSERT PULL CORDS FOR ALL CABLE RUN CONDUITS BETWEEN THE EQUIPMENT ROOM AND THE OPERATORS CONTROL ROOM.
- o 10 FOOT PIGTAILS AT ALL JUNCTION POINTS.
- o ALL WIRING MUST BE THIN OR TFFN STRANDED COPPER THERMOPLASTIC 600 VOLT OR EQUIVALENT INSULATION. **ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.**
- o GROUNDING IS CRITICAL TO EQUIPMENT FUNCTION AND PATIENT SAFETY. SITE MUST CONFORM TO WIRING SPECIFICATIONS SHOWN ON THIS PLAN.



IE	TO	IH	ONE 2 1/2" CND.
----	----	----	-----------------

	DESCRIPTION	QTY.	HARDWARE	DELIV. NO., SHIP. C.
ICC	EZ-EM EMPOWER REMOTE CONTROL	1 1 1	COVERPLATE SINGLE GANG BOX 1 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
IH	INJECTOR HEAD	1 1 1	COVERPLATE SINGLE GANG BOX 1 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
ICC	EZ-EM PERCUPUMP INJECTOR CONTROL	1 1 1	4 X 4 X 4 IN. BOX COVERPLATE 2 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
IH	INJECTOR HEAD	1 1 1	4 X 4 X 4 IN. BOX 1 1/2 IN. DIA. CHASE NIPPLE COVERPLATE	ELEC-8
IE	LF 9000 INJECTOR ELECTRONICS	1 1 1	2 1/2 IN. DIA. CHASE NIPPLE 6 X 6 X 4 IN. BOX	ELEC-8
IEC	INJECTOR CONTROL	1	CONNECT EXTERNALLY	
IH	INJECTOR HEAD	1 1 1	6 X 6 X 4 IN. BOX 2 1/2 IN. DIA. CHASE NIPPLE COVERPLATE	ELEC-8
ICC	MEDRAD ENVISION INJECTOR CONTROL CONSOLE	1	EXTERNALLY CONNECTED	
IE	INJECTOR ELECTRONICS	1 1 1	6 X 6 X 4 IN. BOX COVERPLATE 2 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
IH	INJECTOR HEAD	1 1 1	2 1/2 IN. DIA. CHASE NIPPLE 6 X 6 X 4 IN. BOX COVERPLATE	ELEC-8
IE	MEDRAD MARK IV/PLUS INJECTOR ELECTRONICS	ECT 1 1 1	COVERPLATE 6 X 6 X 4 IN. BOX 2 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
IEC	INJECTOR ELECTRONICS CONTROL	1	CONNECT EXTERNALLY	
IH	INJECTOR HEAD	1 1 1	COVERPLATE 1/2 IN. DIA. CHASE NIPPLE 6 X 6 X 4 IN. BOX	ELEC-8
IED	INJECTOR OUTLET	1 1 1	SINGLE GANG BOX DUPLEX HOSPITAL GRADE RECEPTACLE AND COVERPLATE	
IE	MEDRAD PROVIS INJECTOR ELECTRONICS	1 1 1	COVERPLATE 6 X 6 X 4 IN. BOX 2 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
IEC	INJECTOR ELECTRONICS CONTROL	1	CONNECT EXTERNALLY	
IH	INJECTOR HEAD	1 1 1	6 X 6 X 4 IN. BOX COVERPLATE 2 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
IE	MEDRAD STELLANT INJECTOR ELECTRONICS	1 1 1	2 1/2 IN. DIA. CHASE NIPPLE COVERPLATE 6 X 6 X 4 IN. BOX	ELEC-8
IH	INJECTOR HEAD	1 1 1	COVERPLATE 2 1/2 IN. DIA. CHASE NIPPLE SINGLE GANG BOX	ELEC-8
IH	INJECTOR HEAD (PEDESTAL)	1 1 1	COVERPLATE 2 1/2 IN. DIA. CHASE NIPPLE 6 X 6 X 4 IN. BOX	ELEC-8
ICC	NEMOTO INJECTOR CONTROL CONSOLE	1	EXTERNALLY CONNECTED	
IE	INJECTOR ELECTRONICS	1 1 1	6 X 6 X 4 IN. BOX COVERPLATE 2 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
IH	INJECTOR HEAD	1 1 1	SINGLE GANG BOX COVERPLATE 2 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
IH	INJECTOR HEAD (PEDESTAL)	1 1 1	2 1/2 IN. DIA. CHASE NIPPLE 6 X 6 X 4 IN. BOX COVERPLATE	ELEC-8

WIRE RUN, FROM — TO	QUANTITY, WIRE SIZE/COLOR

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO DETAILS TO ACTUAL EQUIPMENT EXPECTED TO BE INSTALLED. IT IS NOT TO BE USED FOR ACTUAL CONSTRUCTION PURPOSES, HOWEVER, AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.


SHEET  
E 1

 **GE Healthcare**

---

S Services Design Center  
Milwaukee, Wisconsin

INTERCONNECT DIAGRAM



POWER SPECIFICATIONS

ELECTRICAL NOTES

- NOTE 1: ALL WIRES SPECIFIED SHALL BE COPPER STRANDED, FLEXIBLE, THERMO-PLASTIC, COLOR CODED, CUT 10 FOOT LONG AT OUTLET BOXES, DUCT TERMINATION POINTS OR STUBBED CONDUIT ENDS. ALL CONDUCTORS, POWER, SIGNAL AND GROUND, MUST BE RUN IN A CONDUIT OR DUCT SYSTEM. ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS. WIRE RUNS MUST BE CONTINUOUS COPPER STRANDED AND FREE FROM SPLICES. **ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.**
- NOTE 2: WIRE SIZES GIVEN ARE FOR USE OF EQUIPMENT. LARGER SIZES MAY BE REQUIRED BY LOCAL CODES.
- NOTE 3: IT IS RECOMMENDED THAT ALL WIRES BE COLOR CODED, AS REQUIRED IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 4: CONDUIT SIZES SHALL BE VERIFIED BY THE ARCHITECT, ELECTRICAL ENGINEER OR CONTRACTOR, IN ACCORDANCE WITH LOCAL OR NATIONAL CODES.
- NOTE 5: CONVENIENCE OUTLETS ARE NOT ILLUSTRATED. THEIR NUMBER AND LOCATION ARE TO BE SPECIFIED BY OTHERS. LOCATE AT LEAST ONE CONVENIENCE OUTLET CLOSE TO THE SYSTEM CONTROL, THE POWER DISTRIBUTION UNIT AND ONE ON EACH WALL OF THE PROCEDURE ROOM. USE HOSPITAL APPROVED OUTLET OR EQUIVALENT.
- NOTE 6: GENERAL ROOM ILLUMINATION IS NOT ILLUSTRATED. CAUTION SHOULD BE TAKEN TO AVOID EXCESSIVE HEAT FROM OVERHEAD SPOTLIGHTS. DAMAGE CAN OCCUR TO CEILING MOUNTING COMPONENTS AND WIRING IF HIGH WATTAGE BULBS ARE USED. RECOMMEND LOW WATTAGE BULBS NO HIGHER THAN 75 WATTS AND USE DIMMER CONTROLS (EXCEPT MR). DO NOT MOUNT LIGHTS DIRECTLY ABOVE AREAS WHERE CEILING MOUNTED ACCESSORIES WILL BE PARKED.
- NOTE 7: **ROUTING OF CABLE DUCTWORK, CONDUITS, ETC., MUST RUN DIRECT AS POSSIBLE OTHERWISE MAY RESULT IN THE NEED FOR GREATER THAN STANDARD CABLE LENGTHS (REFER TO THE INTERCONNECTION DIAGRAM FOR MAXIMUM USABLE LENGTHS POINT TO POINT).**
- NOTE 8: CONDUIT TURNS TO HAVE LARGE, SWEEPING BENDS WITH MINIMUM RADIUS IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 9: A SPECIAL GROUNDING SYSTEM IS REQUIRED IN ALL PROCEDURE ROOMS BY SOME NATIONAL AND LOCAL CODES. IT IS RECOMMENDED IN AREAS WHERE PATIENTS MIGHT BE EXAMINED OR TREATED UNDER PRESENT, FUTURE, OR EMERGENCY CONDITIONS. CONSULT THE GOVERNING ELECTRICAL CODE AND CONFER WITH APPROPRIATE CUSTOMER ADMINISTRATIVE PERSONNEL TO DETERMINE THE AREAS REQUIRING THIS TYPE OF GROUNDING SYSTEM.
- NOTE 10: THE MAXIMUM POINT TO POINT DISTANCES ILLUSTRATED ON THIS DRAWING MUST NOT BE EXCEEDED.
- NOTE 11: PHYSICAL CONNECTION OF PRIMARY POWER TO GE EQUIPMENT IS TO BE MADE BY CUSTOMERS ELECTRICAL CONTRACTOR WITH THE SUPERVISION OF A GE REPRESENTATIVE. THE GE REPRESENTATIVE WOULD BE REQUIRED TO IDENTIFY THE PHYSICAL CONNECTION LOCATION, AND INSURE PROPER HANDLING OF GE EQUIPMENT.

DIAGRAM KEY

- CUSTOMER/CONTRACTOR SUPPLIED WIRING. ROUTE IN ADEQUATE CONDUIT OR RACEWAY.
- \_\_\_\_\_ GE FURNISHED CABLE RUNS. ROUTE IN EMPTY CONDUIT OR RACEWAY.
- 59' [18M] MAXIMUM RUN LENGTH BETWEEN JUNCTION POINTS.  
Feet [Meters]

SHEET TITLE: ELECTRICAL SPECIFICATIONS

MODALITY TYPE: INJECTORS CT & PET

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS, IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST EDITIONS OF THE NATIONAL ELECTRICAL CODE, THE NATIONAL FIRE PROTECTION ASSOCIATION CONSTRUCTION PRACTICES, HOWEVER, THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:

INJECTOR  
TYPICAL FINAL  
CEILING & PEDESTAL MOUNT

PROJECT	REVISION
TYPINJ	03
DATE:	01.Jul.15
DRAWN BY:	DJP
CHECKED BY:	TMS

REVISION HISTORY:

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SHEET

E2



**GE Healthcare**


IS Services Design Center

Minwaukee, Wisconsin

ELECTRICAL DETAIL  
BOX WITH COVERPLATE (TYPICAL)

ELEC-8

REV. DATE: 09/30/94



DETAIL NOT TO SCALE

PROJECT TITLE:

INJECTOR  
TYPICAL FINAL  
CEILING & PEDESTAL MOUNT

PROJECT	REVISION
TYPINJ	03
DATE:	01.Jul.15
DRAWN BY:	DJP
CHECKED BY:	TMS

REVISION HISTORY:

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SHEET

E3

SHEET TITLE: ELECTRICAL DETAILS  
MODALITY TYPE: INJECTORS CT & PET

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS, IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST EDITIONS OF THE NATIONAL ELECTRICAL CODE, THE NATIONAL FIRE PROTECTION ASSOCIATION (NFPA) CODES, AND THE NATIONAL BUILDING CODES. HOWEVER, THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.



GE Healthcare

IS Services Design Center

Milwaukee, Wisconsin



EQUIPMENT DETAIL  
NEMOTO INJECTOR ELECTRONICS

B8007C  
REV. 00: 7/9/08

Diagram showing the Display Unit and Main Unit of the NEMOTO INJECTOR ELECTRONICS. The Display Unit has a height of 9.7" [248mm] and a width of 4.2" [107mm]. The Main Unit has a height of 7.8" [190mm] and a width of 9.4" [240mm]. The Display Unit is 10.6" [270mm] wide. The Main Unit is 3.9" [100mm] wide.

DISPLAY UNIT      DRAWING NOT TO SCALE      MAIN UNIT

EQUIPMENT DETAIL  
NEMOTO INJECTOR ON PEDESTAL

B80-07D  
REV. 00: 7/10/08

Diagram showing the NEMOTO INJECTOR ON PEDESTAL. The Plan View shows a height of 20.47" [520mm]. The Side View shows a height of 49.84" [1266mm].

PLAN VIEW      DRAWING NOT TO SCALE      SIDE VIEW

EQUIPMENT DETAIL  
NEMOTO INJECTOR ON OVERHEAD  
COUNTERPOISED SUSPENSION

B50-31N  
REV. DATE: 06/23/09

VERIFY MOUNTING ASSEMBLY DIMENSIONS  
WITH INJECTOR MANUFACTURER.

Diagram showing the NEMOTO INJECTOR ON OVERHEAD COUNTERPOISED SUSPENSION. The diagram includes a table for POST LENGTHS and GE ECAT #, and a table for MIN. PLATE HEIGHT and MAX. PLATE HEIGHT. The diagram also shows the TOP OF MOUNTING ASSEMBLY, 65" MAX. ARC [1660mm], HIGHEST HEAD POSITION, INJECTOR, and LOWEST HEAD POSITION.

POST	LENGTHS	GE ECAT #
SHORT	22.8" / 580mm	E8004NA
INTERMEDIATE	33.5" / 850mm	E8004NJ
LONG	N/A	N/A

POST	MIN. PLATE HEIGHT	MAX. PLATE HEIGHT
SHORT	107.75" [2737mm]	113.75" [2889mm]
INTERMEDIATE	118.5" [3010mm]	124.5" [3162mm]
LONG	-	-

RECOMMENDED HEIGHT TO FLOOR  
82"-88" [2083mm - 2235mm]

DRAWING NOT TO SCALE

EQUIPMENT DETAIL  
STELLANT INJECTOR ELECTRONICS

E80-07N  
REV. DATE: 19.SEP.13

Diagram showing the STELLANT INJECTOR ELECTRONICS. The Front View shows a height of 13.5" [343mm] and a width of 11.0" [279mm]. The Side View shows a height of 11.5" [292mm] and a width of 8.9" [225mm]. The Plan View shows a height of 8.8" [222mm] and a width of 12.5" [319mm].

FRONT VIEW      SIDE VIEW      PLAN VIEW

DRAWING NOT TO SCALE

EQUIPMENT DETAIL  
ENVISION INJECTOR ON PEDESTAL

B80-07A  
REV. 01: 7/10/08

Diagram showing the ENVISION INJECTOR ON PEDESTAL. The Plan View shows a height of 17.94" [456mm]. The Side View shows a height of 48.75" [1240mm].

PLAN VIEW      DRAWING NOT TO SCALE      SIDE VIEW

EQUIPMENT DETAIL  
MEDRAD INJECTOR ON OVERHEAD  
COUNTERPOISED SUSPENSION

B50-31  
REV. DATE: 01.Nov.12

VERIFY MOUNTING ASSEMBLY DIMENSIONS  
WITH INJECTOR MANUFACTURER.

Diagram showing the MEDRAD INJECTOR ON OVERHEAD COUNTERPOISED SUSPENSION. The diagram includes a table for POST LENGTHS and GE ECAT #, and a table for MIN. PLATE HEIGHT and MAX. PLATE HEIGHT. The diagram also shows the TOP OF MOUNTING ASSEMBLY, 65" MAX. ARC [1660mm], HIGHEST HEAD POSITION, INJECTOR, and LOWEST HEAD POSITION.

POST	LENGTHS	GE ECAT #
SHORT	22.8" / 580mm	E8007NB, ND, NG, NY, PB & PE
INTERMEDIATE	33.5" / 850mm	E8007PP, PN, PK, PR & PY
LONG	39.4" / 1000mm	E8007NE, NF, NJ, NN, PD & E8018AA

POST	MIN. PLATE HEIGHT	MAX. PLATE HEIGHT
SHORT	108" [2743mm]	114" [2896mm]
INTERMEDIATE	120" [3048mm]	
LONG	126" [3200mm]	132" [3353mm]

RECOMMENDED HEIGHT TO FLOOR  
82"-88" [2083mm - 2235mm]

DRAWING NOT TO SCALE

EQUIPMENT DETAIL  
INJECTOR REMOTE CONTROL AND ELECTRONICS

B50-28  
REV. DATE: 02.AUG.12

Diagram showing the INJECTOR REMOTE CONTROL AND ELECTRONICS. The Remote Control is 7.5" [191mm] high and 14.5" [368mm] wide. The Electronics are 10.5" [267mm] high and 17" [432mm] wide. The Remote Control is 3.25" [83mm] wide.

DRAWING NOT TO SCALE

EQUIPMENT DETAIL  
MEDRAD PROVIS INJECTOR

B5031W  
REV. DATE: 06/02/03

Diagram showing the MEDRAD PROVIS INJECTOR. The Standard Ceiling Mount shows a height of 32" min. [813mm] and a width of 7" [178mm]. The Short Ceiling Mount shows a height of 25" [635mm] and a width of 7" [178mm].

STANDARD CEILING MOUNT      SHORT CEILING MOUNT

DRAWING NOT TO SCALE

EQUIPMENT DETAIL  
MEDRAD PROVIS INJECTOR

B5031X  
REV. DATE: 06/02/03

Diagram showing the MEDRAD PROVIS INJECTOR. The diagram shows a height of 33.9" [860mm] and a width of 65.4" [1660mm].

DETAIL NOT TO SCALE

EQUIPMENT DETAIL  
INJECTOR ON MOBILE PEDESTAL

B50-30  
REV. DATE: 28.JUN.12

Diagram showing the INJECTOR ON MOBILE PEDESTAL. The Plan View shows a height of 26" [660mm] and a width of 52" [1321mm]. The Side View shows a height of 55" [1397mm] and a width of 26" [660mm].

PLAN VIEW      DRAWING NOT TO SCALE      SIDE VIEW

EQUIPMENT DETAIL  
ENVISION INJECTOR ELECTRONICS

B80-07  
REV. 00: 12/3/97

Diagram showing the ENVISION INJECTOR ELECTRONICS. The Display Control Console is 11.5" [292mm] high and 13.0" [330mm] wide. The Electronics Console is 12.31" [313mm] high and 13.75" [350mm] wide.

DISPLAY CONTROL CONSOLE      ELECTRONICS CONSOLE

DRAWING NOT TO SCALE

EQUIPMENT DETAIL  
LF 9000 INJECTOR

B50-29C  
REV. 00: 08/08/97

Diagram showing the LF 9000 INJECTOR. The Powerhead is 6" [150mm] high and 16" [410mm] wide. The Control is 12" [300mm] high and 8" [200mm] wide. The Electronics are 12" [300mm] high and 14" [360mm] wide.

POWERHEAD      CONTROL      ELECTRONICS

DRAWING NOT TO SCALE

GE Healthcare

IS Services Design Center

Minwaukee, Wisconsin

SHEET TITLE: EQUIPMENT DETAILS

MODALITY TYPE: INJECTORS CT & PET

THIS PLAN IS SUBMITTED TO REQUEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO THE ACTING CONSTRUCTION ENGINEER'S REQUIREMENTS. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

INJECTOR

TYPICAL FINAL

CEILING & PEDESTAL MOUNT

PROJECT TITLE:

PROJECT REVISION

TYPINJ 03

DATE: 01.Jul.15

DRAWN BY: DJP

CHECKED BY: TMS

REVISION HISTORY:

SHEET


D1



EQUIPMENT DETAIL  
INJECTOR ON OVERHEAD  
COUNTERPOISED SUSPENSION

B50-31U

REV. DATE: 08/03/98




DRAWING NOT TO SCALE

EQUIPMENT DETAIL  
EMPOWER INJECTOR POWER SUPPLY

B8009B

REV. DATE: 12/18/02




DETAIL NOT TO SCALE

EQUIPMENT DETAIL  
EMPOWER INJECTOR REMOTE CONTROL (DESKTOP)

B8009C

REV. DATE: 12/18/02




DETAIL NOT TO SCALE

EQUIPMENT DETAIL  
EMPOWER INJECTOR ON PEDESTAL

B8009A

REV. DATE: 12/18/02




DETAIL NOT TO SCALE

EQUIPMENT DETAIL  
EMPOWER CT INJECTOR ON CEILING MOUNT

B8009

REV. DATE: 12/18/02



DETAIL NOT TO SCALE

SHEET TITLE: EQUIPMENT DETAILS

MODALITY TYPE: INJECTORS CT & PET

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS, IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO THE ACTUAL CONSTRUCTION. HOWEVER, THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:

INJECTOR  
TYPICAL FINAL  
CEILING & PEDESTAL MOUNT

PROJECT	REVISION
TYPINJ	03
DATE:	01.Jul.15
DRAWN BY:	DJP
CHECKED BY:	TMS

REVISION HISTORY:

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SHEET

D2