




OPTIMA RF420 – NON-ELEVATING TABLE



ROOM DIMENSIONS	LENGTH & WIDTH	CEILING HEIGHT
MINIMUM:	20'-6" x 16'-0" [6.25m x 4.88m]	8'-6" [2.6m]

ROOM PLAN LOCATED ON PAGE 2-66

ITEM	QUANTITY	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
①	1	NON-ELEVATING TABLE	1322 lbs	1535 btu
②	1	WALL STAND	489 lbs	119 btu
③	1	GENERATOR <CONTROL CABINET>	595 lbs	3070 btu
④	1	STARTER	110 lbs	
⑤	1	HEAT EXCHANGER	22 lbs	
⑥	1	MONITOR CART	119 lbs	
⑦	1	OPERATOR CONSOLE	64 lbs	

ENVIRONMENT

Ambient Operation Temperature:


Allowable Temperature Change: 50° TO 95° (F), 10° TO 35° (C)

Humidity:

Allowable Humidity Change: 30% to 85% (no condensation)

Max - Altitude: 700 to 1060 hectopascals

ROOM PLAN LOCATED ON PAGE 2-66

ITEM 	QUANTITY	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
60		COUNTER TOP FOR EQUIPMENT - MINIMUM DEPTH 24" AND ADDITIONAL SHELVING MAY BE REQUIRED BELOW COUNTER TOP FOR PC TOWER. PROVIDE GROMMETED OPENINGS AS REQUIRED TO ROUTE CABLES.		
61		CONTROL WALL TO CEILING WITH LEAD GLASS VIEWING WINDOW.		
62		MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 47.2"W x 78.7"H [1200mm x 2000mm], CONTINGENT ON A 62.9" [1600mm] CORRIDOR WIDTH		
63		X-RAY ON WARNING LIGHT - AVAILABLE FROM GEXPRO CALL: 800-200-9760 GE CAT. NO. WXIABWW-0F-XIU		
64		DOOR LIMIT SWITCH <NEEDED ONLY IF REQUIRED BY STATE/LOCAL CODES>		
90		X-RAY ROOM WARNING LIGHT/ROOM LIGHTING CONTROL PANEL CAT. NO. E4502SS		
91		MAIN DISCONNECT,		