





INNOVA BIPLANE IGS 620/630 INTERVENTIONAL RADIOLOGY (IR)



NOTE: THIS IS A SUGGESTED LAYOUT BASED OFF OF GE HEALTHCARE PRE-INSTALL MANUAL. PLEASE CONTACT LOCAL PROJECT MANAGER OF INSTALLATION IF YOUR ROOM LAYOUT OR DIMENSIONS ARE DIFFERENT THAN THE LAYOUT SHOWN.

| ROOM DIMENSIONS | LENGTH & WIDTH | CEILING HEIGHT |
|-----------------|--------------------------------|----------------|
| RECOMMENDED: | 38'-0" x 26'-2" [11.5m x 8.0m] | 9'-4" [2.85m] |



SEE APPENDIX "A" FOR POWER SPECIFICATIONS

ROOM PLAN LOCATED ON PAGE 4-71

| ITEM | QUANTITY | ITEM DESCRIPTION | WEIGHT | HEAT OUTPUT |
|------|----------|---|----------|-------------|
| ① | 1 | 3 KVA UPS CABINET <LARGE DISPLAY SUBSYSTEM OPTION> | 77 lbs | 546 btu |
| ② | 1 | AP COOLIX 4100 AUTOTRANSFORMER | 66 lbs | 238 btu |
| ③ | 1 | AP COOLIX 4100 WATER CHILLER | 264 lbs | 18730 btu |
| ④ | 1 | AP DETECTOR CHILLER | 33 lbs | 706 btu |
| ⑤ | 1 | OPERATORS CONSOLE | 22 lbs | 546 btu |
| ⑥ | 1 | LATERAL CABINET <C3> | 703 lbs | 2945 btu |
| ⑦ | 1 | LATERAL POSITIONER BRIDGE MOUNT ASSEMBLY MOUNTED FROM CEILING SUPPORTS | 1421 lbs | 4126 btu |
| ⑧ | 2 | LONGITUDINAL STATIONARY RAIL FOR LATERAL GANTRY INNOVA POSITIONER | 68 lbs | |
| ⑨ | 1 | LARGE DISPLAY MONITOR ON SINGLE MONITOR SUSPENSION 7 ft. 9 in. INBOARD BRIDGE <MOUNT TWO GE MONITORS ON BACKSIDE OF LD MONITOR> | 784 lbs | 1706 btu |
| ⑩ | 1 | LARGE DISPLAY MONITOR CABINET | 253 lbs | 3412 btu |
| ⑪ | 2 | 18 in. MONITOR ON WALL SUPPORT | 26 lbs | 204 btu |
| ⑫ | 1 | IVUS VOLCANO S51 CONSOLE, INCLUDES FLAT PANEL MONITOR AND KEYBOARD <DESK MOUNTED> | 68 lbs | 1631 btu |
| ⑬ | 1 | NURSING NOTES WORKSTATION | 46 lbs | 682 btu |
| ⑭ | 1 | WORKSTATION CART | | |
| ⑮ | 1 | MAC-LAB CONSOLE, INCLUDES MONITORS AND KEYBOARD | 566 lbs | 2935 btu |
| ⑯ | 1 | TRAM NET RACK | 8 lbs | |
| ⑰ | 1 | COLOR PRINTER | | 1054 btu |
| ⑱ | 2 | COUNTERBALANCED EYE AND THYROID SHIELD WITH R96 LAMP | 143 lbs | |
| ⑲ | 1 | REMOTE CONTROL FOR INJECTOR | 4 lbs | |
| ⑳ | 1 | INJECTOR ELECTRONICS ON SHELF | 37 lbs | 320 btu |
| ㉑ | 1 | INJECTOR HEAD ON TABLE RAIL | 15 lbs | |
| ㉒ | 2 | LONGITUDINAL STATIONARY RAIL FOR XT SUSPENSION | 68 lbs | |
| ㉓ | 1 | INNOVA V TABLE | 1300 lbs | 614 btu |
| ㉔ | 1 | XR BUZZER <LOCATED ABOVE CEILING> | 2 lbs | |
| ㉕ | 1 | INNOVA POSITIONER <REFERENCE TABLE BASE-PLATE DETAIL FOR FLOOR MOUNTING INFORMATION> | 1653 lbs | 2416 btu |
| ㉖ | 1 | ATLAS CABINET <C2> | 621 lbs | 4570 btu |

ENVIRONMENT

Ambient Operation Temperature: 68°-77° F, (20°-25° C)
 Allowable Temperature Change: Not Available
 Humidity: 30%-75% non condensing
 Allowable Humidity Change: Not Available
 Max - Altitude: 9842 ft. (3000 M) above sea level

ROOM PLAN LOCATED ON PAGE 4-71

| ITEM | QUANTITY | ITEM DESCRIPTION | WEIGHT | HEAT OUTPUT |
|--|----------|---|----------|-------------|
| 27 | 1 | ATLAS CABINET <C1> | 888 lbs | 4413 btu |
| 28 | 1 | UPS CABINET | 1170 lbs | 4061 btu |
| 29 | 1 | BOLUS CHASE HANDSWITCH | 2 lbs | |
| 30 | 1 | AW WORKSTATION | 81 lbs | 1201 btu |
| 31 | 1 | 3 KVA UPS CABINET | 77 lbs | 546 btu |
| 32 | 1 | LATERAL COOLIX 4100 AUTOTRANSFORMER | 66 lbs | 238 btu |
| 33 | 1 | LATERAL COOLIX 4100 WATER CHILLER | 264 lbs | 18730 btu |
| 34 | 1 | LATERAL DETECTOR CHILLER | 33 lbs | 706 btu |
| 35 | 1 | LOGIQ I ULTRASOUND UNIT <ON CART> | | |
| 50 | 1 | INNOVA MAIN DISCONNECT, REFERENCE JUNCTION POINT 'A' ON SHEET E1 FOR DETAILED DESCRIPTION. | 899 lbs | 2215 btu |
| 51 | 2 | CABLE DRAPE RAIL FOR LP POSITIONER | | |
| 52 | 6 | VITALINQ SPEAKER | | |
| 53 | 1 | VITALINQ CONSOLE | | |
| 54 | 1 | VITALINQ MICROPHONE | | |
| 55 | 1 | VITALINQ MICROPHONE <ONE ON MONITOR BRIDGE IN EXAM ROOM> | | |
| ***** CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS ***** | | | | |
| 60 | | BEARING BLOCK OUTLINE, SEE S1 FOR MORE INFORMATION. | | |
| 61 | | CONTROL WALL TO CEILING WITH LEAD GLASS VIEWING | | |
| 62 | | SHELF - CUSTOMER TO PROVIDE ADEQUATE WALL SUPPORT | | |
| 63 | | X-RAY ON WARNING LIGHT - AVAILABLE FROM GE SUPPLY GE CAT. NO. WXIABW-OF-XIU CALL: 800-200-9760 | | |
| 64 | | CABLE DRAPE RAIL. | | |
| 65 | | MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 44 IN. W x 83 IN. H [1118mm x 2108mm], CONTINGENT ON A 96 IN. [2438mm] CORRIDOR WIDTH | | |
| 66 | | 150-AMP LOCAL SERVICE DISCONNECT FOR LOCK-OUT/TAG-OUT CAPABILITY. <MAY BE A FUSED DISCONNECT, CIRCUIT BREAKER OR SAFETY SWITCH> | | |
| 67 | | CUSTOMER SUPPLIED STORAGE CABINET | | |
| 68 | | CATHETER CABINETS | | |
| 69 | | MED GASES IN CEILING | | |
| 70 | | LEAD APRON RACK | | |
| 71 | | SCRUB SINK | | |
| 72 | | COUNTER TOP FOR EQUIPMENT- PROVIDE GROMMETED OPENINGS AS REQUIRED TO ROUTE INTERCONNECT CABLES TO RACEWAY BELOW COUNTERTOP. | | |
| 73 | | COUNTER TOP WITH BASE AND WALL CABINETS | | |
| 90 | | X-RAY ROOM WARNING LIGHT/ROOM LIGHTING CONTROL PANEL REFERENCE JUNCTION POINT 'XRLC' ON SHEET 'E1' FOR DETAILED DESCRIPTION -CAT. NO. E4502SS FOR WARNING LIGHT & ROOM LIGHT CONTROL. | | |