

ROOM PLAN LOCATED ON PAGE 1-138

ITEM	QUANTITY	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
①	1	SYSTEMS CABINET	679 lbs	1098 btu
②	1	XT RADIOGRAPHIC SUSPENSION WITH INBOARD MOUNTING.	628 lbs	501 btu
③	2	LONGITUDINAL STATIONARY RAIL FOR XT SUSPENSION	68 lbs	
④	1	FLEXI DT MOBILE TABLE	639 lbs	
⑤	1	DETECTOR SUPPORT ASSEMBLY	97 lbs	
⑥	1	OPERATORS CONSOLE	74 lbs	1235 btu
⑦	1	LONGITUDINAL DRIVE BELT 1 IN. WIDE		
⑧	2	ANCHOR RAILS		
⑨	1	GRID HOLDER (FIELD VERIFY IDEAL LOCATION)	35 lbs	
⑩	1	DIGITAL CHEST UNIT WITH EXTENDED RECEPTOR	595 lbs	822 btu
⑪	1	WEIGHT BEARING STAND	101 lbs	
⑫	1	IMAGE PASTING BARRIER	200 lbs	

ENVIRONMENT

Ambient Operation Temperature: 59-95° (F)
 Allowable Temperature Change: 18° (F) per hour
 Humidity: 20-80 percent non condensing
 Allowable Humidity Change: 10 percent per hour
 Max - Altitude: 8000 ft. above sea level

ROOM PLAN LOCATED ON PAGE

ITEM ◻	QUANTITY	ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
60		COUNTER TOP FOR EQUIPMENT- MINIMUM DEPTH 24 IN. AND ADDITIONAL SHELVING MAY BE REQUIRED BELOW COUNTER TOP FOR PC TOWER. PROVIDE GROMMETED OPENINGS AS REQUIRED TO ROUTE CABLES.		
61		COUNTER TOP WITH BASE CABINETS		
62		CONTROL WALL, 7 FT. HIGH WITH LEAD GLASS VIEWING WINDOW.		
63		CABLE DRAPE RAIL.		
64		X-RAY ON WARNING LIGHT - AVAILABLE FROM GE SUPPLY CALL: 800-200-9760 GE CAT. NO. WXIABWW-DF-XIU		
65		DOOR LIMIT SWITCH <NEEDED ONLY IF REQUIRED BY STATE/LOCAL CODES>		
66		MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 40 IN. W x 82 IN. H [1016mm x 2083mm], CONTINGENT ON A 72 IN. [1828mm] CORRIDOR WIDTH		
90		X-RAY ROOM WARNING LIGHT CONTROL PANEL REFERENCE JUNCTION POINT 'XRLC' ON SHEET 'E1' FOR DETAILED DESCRIPTION -E4502RL FOR WARNING LIGHT CONTROL ONLY.		
91		MAIN DISCONNECT, REFERENCE JUNCTION POINT 'A' ON SHEET E1 FOR DETAILED DESCRIPTION. CAT. NO. E4502RS OR WITH AUTO RESTART E4502SA. <20 W x 48 H x 6.68 in. D>		
		***** CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS *****		