





XR6000



ROOM DIMENSIONS	LENGTH & WIDTH	CEILING HEIGHT
MINIMUM:	VARIABLES: BASED ON EQUIPMENT AND APPLICATIONS	2500mm [8'-2.5"]

ROOM PLAN LOCATED ON PAGE 1-158

ITEM	QUANTITY			
ITEM		ITEM DESCRIPTION	WEIGHT	HEAT OUTPUT
①	1	STANDARD WALL STAND	107 kg 235 lbs	
②	1	XR6000 CONTROL CONSOLE PEDESTAL MOUNT	6 kg 13 lbs	6 W 20 btu
③	1	XR6000 TABLE WITH TUBESTAND	240 kg 529 lbs	725 W 2474 btu
<p>*****                      CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS                      *****</p>				
60		CONTROL WALL TO CEILING WITH LEAD GLASS VIEWING WINDOW.		
61		X-RAY ON WARNING LIGHT - AVAILABLE FROM GE SUPPLY CALL: 800-200-9760 GE CAT. NO. WXIABWW-OF-XIU		
62		DOOR LIMIT SWITCH <NEEDED ONLY IF REQUIRED BY STATE/LOCAL CODES>		
63		MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 1000MM W X 2250MM H CONTINGENT ON A 2500MM CORRIDOR WIDTH		
90		X-RAY ROOM WARNING LIGHT CONTROL PANEL REFERENCE JUNCTION POINT 'XRLC' ON SHEET 'E1' FOR DETAILED DESCRIPTION -E4502RL FOR WARNING LIGHT CONTROL ONLY.		
91		EMERGENCY OFF SWITCH <NEEDED ONLY IF REQUIRED BY STATE/LOCAL CODES>		
92		MAIN DISCONNECT, REFERENCE JUNCTION POINT 'A' ON SHEET E1. <16" W x 24" H x 6" D> [400MM W x 600MM H x 150MM D]		

ENVIRONMENT

Ambient Operation Temperature: 10-40° <C> / 50-104° <F>  
 Allowable Temperature Change: <10° <C> / 50° <F> per hour  
 Humidity: 30-75 percent non condensing  
 Allowable Humidity Change: -  
 Max - Altitude: 3000M <MIN: -100M> 9842' <MIN -328' >