

Vscan Extend

Protect patients. Staff. And your device.

Proper cleaning and disinfection between patients is necessary to help prevent disease transmission. By following recommended care and handling procedures, you can help maximize your system's performance and product life.

This document is intended to provide an overview of the necessary steps for probe and device reprocessing. Reference the Vscan Extend™ user manual for detailed instructions.

VSCAN EXTEND
USER MANUAL

Inspect the probe and device

Examine the cable connectors for mechanical defects and inspect the lens, probe housing and the cable for damage that could allow liquid into the probe. Be sure to check the entire length of the electrical cables for cracks, cuts and abrasions. If you notice any defects, do not use the Vscan Extend until it has been properly repaired.

Clean the device after each use; before disinfection

Cleaning is a separate, important step that should be performed prior to disinfecting the probe and device.

Visit the Device Cleaner/Disinfectant Compatibility Website for an overview of the various agents we have tested for compatibility with our device materials (compatibility testing only; not related to cleaning/disinfection efficacy).

CLEANER
COMPATIBILITY


Clean and Disinfect the probe after each use

Properly disinfecting the probe is a necessary step in protecting patients, staff and the probe. Disinfection procedures are based on the type of procedure the probe was used for.

Please reference the Vscan Extend User Manual for a complete list of steps and detailed instructions.

Visit the Transducer Cleaning and Disinfection Compatibility Page. You will need to search for "Vscan dual (sector/linear) or Vscan sector"

PROBE CLEANING
& DISINFECTION





GE Healthcare is a leading global medical technology and digital solutions innovator. GE Healthcare enables clinicians to make faster, more informed decisions through intelligent devices, data analytics, applications and services, supported by its Edison intelligence platform. With over 100 years of healthcare industry experience and around 50,000 employees globally, the company operates at the center of an ecosystem working toward precision health, digitizing healthcare, helping drive productivity and improve outcomes for patients, providers, health systems and researchers around the world.

Follow us on Facebook, LinkedIn, Twitter and Insights, or visit our website www.gehealthcare.com for more information.

For more information please visit:
vscan.gehealthcare.com

Imagination at work

© 2020 General Electric Company - All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, and Vscan Extend are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB74257XX(1)