



***Mac-Lab/CardioLab Anti-Virus Information
Software versions 6.5.3, 6.5.4, 6.5.6, 6.8, 6.8.1, 6.9, and 6.9.5***

Product Group: Interventional Invasive Products

Products: Mac-Lab IT/XT/XTi, CardioLab IT/XT/XTi,
SpecialsLab and ComboLab IT/XT/XTi
Recording Systems, Centricity Cardiology Data
Management Systems

Subject: Anti-Virus Information

Date: 6-July-2017



Contents

Mac-Lab/CardoLab v6.9/6.9.5 Anti-Virus Installation.....	3
Anti-Virus Requirements.....	3
Validated Anti-Virus Software.....	4
Customer-Provided Server Configuration.....	4
Anti-Virus Software Installation Settings.....	5
Symantec EndPoint Protection v12.1.2.....	7
Installation Overview.....	7
Pre-Installation Guidelines.....	7
Symantec EndPoint Protection Deployment Steps (Preferred Push Installation Method).....	7
Symantec EndPoint Protection Server Console Configurations.....	8
Symantec EndPoint Protection Post Installation Guidelines.....	10
McAfee VirusScan Enterprise v8.8 Patch 2.....	10
Installation Overview.....	10
McAfee VirusScan Enterprise Installation Procedure.....	11
McAfee VirusScan Enterprise Configuration.....	11
McAfee ePolicy Orchestrator v5.0.....	13
Installation Overview.....	13
Pre-Installation Guidelines.....	13
McAfee ePolicy Orchestrator Deployment Steps (Preferred Push Installation Method).....	13
McAfee ePolicy Orchestrator Server Console Configuration.....	15
McAfee ePolicy Orchestrator Post Installation.....	17
Trend Micro OfficeScan Client/Server Edition v10.6 SP2.....	17
Installation Overview.....	17
Pre-Installation Guidelines.....	17
Trend Micro OfficeScan Deployment Steps (Preferred Push Installation Method).....	18
Trend Micro OfficeScan Server Console Configuration.....	19
Trend Micro OfficeScan Post Installation Guidelines.....	21



Mac-Lab/CardoLab/INW Server v6.5.3/6.5.4/6.5.6/6.8/ 6.8.1/6.9/6.9.5 Anti-Virus Installation

Anti-virus software supports facilities in complying with privacy regulations, such as HIPAA.

Anti-Virus Requirements

WARNING:

ANTI-VIRUS SOFTWARE INSTALLATION

The System is delivered without anti-virus protection. It is recommended to have validated anti-virus software installed on the system before connecting to any network. Lack of validated virus protection could lead to system instability or failure.

Note the following requirements:

- The MLCL v6.5.3/6.5.4/6.5.6/6.8/6.8.1/6.9/6.9.5 Anti-Virus solutions listed below are in addition to the qualified Anti-Virus solutions listed in the applicable MLCL Security Guide.
- For MLCL v6.5.3/6.5.4/6.5.6/6.8/6.8.1/6.9/6.9.5 systems, **the listed Anti-Virus can only be installed after all applicable security patches are installed. The MLCL v6.5.3 and v6.5.4 Acquisition and Review systems also requires SP3 installed for Windows XP. The MLCL v6.9 INW Server also requires SP1 installed for Windows Server 2008 R2.**
- Anti-virus software is not provided with the Mac-Lab/CardioLab system and is the customer's responsibility to acquire, install, and maintain.
- The customer is responsible for updating anti-virus definition files.
- If a virus is found contact the facility System Administrator and GE Technical Support.
- Install only the anti-virus software packages listed in the listed in the Validated Anti-Virus Software section.
- Log in as an Administrator or member of that group to perform the activities in this document.
- Use a language version of the validated anti-virus software that matches the operating system language if possible. If there is no validated anti-virus software that matches the operating system language, install the English version of the anti-virus software.



Validated Anti-Virus Software

WARNING:

SYSTEM INSTABILITY

Do not install or use unvalidated anti-virus software (including unvalidated versions). Doing so may result in system instability or failure. Use only validated anti-virus software in the appropriate language version.

NOTE:

If the language specific anti-virus software is not available, install the English version of anti-virus software.

The Mac-Lab/CardioLab v6.9/6.95 systems have been validated to run with the software listed in the following table.

Supported Anti-Virus Software	Supported Languages	Supported Anti-Virus Software Version
McAfee VirusScan Enterprise	English	8.8 Patch 2
McAfee ePolicy Orchestrator (with McAfee VirusScan Enterprise 8.8 Patch 2)	English	v5.0
Symantec EndPoint Protection	English	12.1.2
Trend Micro OfficeScan Client/Server Edition	English	10.6 SP2, XG 12.0

The Mac-Lab/CardioLab v6.5.3/6.5.4/6.5.6/6.8/6.8.1 systems have been validated to run with the software listed in the following table.

Supported Anti-Virus Software	Supported Languages	Supported Anti-Virus Software Version
Trend Micro OfficeScan Client/Server Edition	English	XG 12.0

NOTE: Previously supported CA Total Defense Anti-Virus is no longer a commercially available product.

Customer-Provided Server Configuration

The anti-virus management console is required to be installed on the Customer-Provided Server.

The communication between Customer-Provided Server and Mac-Lab/CardioLab devices can be accomplished in different ways including:

1. Adding to INW domain.



2. Adding to Member Server domain.
3. Cross domain authentication.

NOTE: The customer-provided server should have two network ports. One network port to connect to the Centricity Cardiology INW network and the second network port to connect to the hospital network.

Anti-Virus Software Installation Settings

Disable Loopback Connection

Note: Perform the below steps only when Loopback Connection is present in Acquisition system.

On an Acquisition system connected to the Mac-Lab/CardioLab environment, disable the Loopback Connection to discover all client systems with the same subnet mask on the domain.

1. Log on as **Administrator** or a member of that group.
2. Right-click **My Network Places** on the desktop and select **Properties**.
3. Right-click **Loopback Connection** and select **Disable**.
4. Restart the Acquisition system.

NOTE: Disabling the Loopback connection on the Acquisition system is required to discover all client systems with same subnet mask on the domain.

Enable Loopback Connection

Note: Perform the below steps only when Loopback Connection is present in Acquisition system.

On an Acquisition systems connected to the Mac-Lab/CardioLab environment, enable the Loopback Connection using the steps below.

1. Log on as **Administrator** or member of that group.
2. Right-click **My Network Places** on the desktop and select **Properties**.
3. Right-click **Loopback Connection** and select **Enable**.
4. Restart the Acquisition system.

Configure Computer Browser Service Before Anti-Virus Installation

Check the Computer Browser service setting on networked Acquisition and Review systems to make sure it is configured correctly.

1. Log on as **Administrator** or member of that group.
2. Click **Start > Run**.
3. Type **services.msc** and press **Enter**.
4. Ensure the **Computer Browser** service is **Started** and **Automatic**.
5. If not then follow further instructions.
6. Double-click the **Computer Browser** service.



7. Change the **Startup type** to **Automatic**.
8. Click **Start**.
9. Click **OK**.
10. Close the **Services** window.

Configure Computer Browser Service After Anti-Virus Installation

Check the Computer Browser service setting on networked Acquisition and Review systems to make sure it is configured correctly.

1. Click **Start > Run**.
2. Type **services.msc** and press **Enter**.
3. Double-click the **Computer Browser** service.
4. Change the **Startup type** to **Manual**.
5. Click **OK**.
6. Close the **Services** window.



Symantec EndPoint Protection v12.1.2

Installation Overview

Install Symantec EndPoint Protection in a networked Mac-Lab/CardioLab environment only. In a networked environment, the Symantec EndPoint Protection must be installed on the customer-provided server and then deployed to the Centricity Cardiology INW server and Acquisition/Review workstation as clients. Use the following instructions to install and configure **Symantec EndPoint Protection** for English.

Virus updates are the responsibility of the facility. Update the definitions regularly to ensure that the latest virus protection is on the system.

Pre-Installation Guidelines

1. On the customer-provided server, make sure Symantec EndPoint Protection Manager is installed before continuing with these steps.
2. On the customer-provided server, open inbound port **8014** and name it **IN_8014_Symantec** and allow inbound connections for domain, public, and private.
3. Log on as Administrator or a member of that group on all client systems (Acquisition, Review, and INW Server) to install the anti-virus software.
4. Disable the Loopback Connection. Refer to Disable Loopback Connection for more information.
5. Configure the Computer Browser service. Refer to Configure Computer Browser Service Before Anti-Virus Installation for more information.

Symantec EndPoint Protection Deployment Steps (Preferred Push Installation Method)

1. Click **Start > All Programs > Symantec EndPoint Protection Manager > Symantec Endpoint Protection Manager**.
2. Enter the appropriate user name and password to log in to Symantec Endpoint Protection Manager. (Click **Yes** if a security prompt displays.)
3. Check **Do not show this Welcome Page again** and click **Close** to close the welcome screen.
4. Click **Home** in the **Symantec Endpoint Protection Manger** window.
5. Select **Install protection client to computers** from the **Common Tasks** drop-down list in the top-right of the **Home** window.
6. Select **New Package Deployment** and click **Next**.
7. Keep the default settings and click **Next**.
8. Select **Remote push** and click **Next**. Wait for the **Computer selection** screen to appear.
9. Expand **<Domain>** (example: INW). Systems connected to the domain are displayed in the **Computer selection** window.

NOTE: If all systems are not being recognized, click **Search Network** and click **Find Computers**. Use the **search by IP address** detection method to identify the client systems (Acquisition, Review, and INW Server).



10. Select all Mac-Lab/CardioLab client machines connected to the domain and click >>. The **Login Credentials** screen appears.
11. Enter the appropriate user name, password and domain name and click **OK**.
12. Make sure all selected machines appear under **Install Protection Client** and click **Next**.
13. Click **Send** and wait until the Symantec anti-virus software is deployed on all client systems (Acquisition, Review, and INW Server). When finished, the **Deployment Summary** screen appears.
14. Click **Next** and then click **Finish** to complete the Client Deployment Wizard.
15. Restart all the client machines (Acquisition, Review, and INW Server). Login with Administrator or as a member of that group on all client machines after the restart.

Symantec EndPoint Protection Server Console Configurations

1. Select **Start > All Programs > Symantec EndPoint Protection Manager > Symantec EndPoint Protection Manager**. The Symantec EndPoint Protection Manager log on window opens.
2. Enter the appropriate Symantec Endpoint Protection Manager Console password and click **Log On**.
3. Select the **Policies** tab and click **Virus and Spyware Protection** under **Policies**. The **Virus and Spyware Protection Policies** window opens.
4. Click **Add a Virus and Spyware Protection** policy under **Tasks**. The **Virus and Spyware Protection** window opens.
5. Under **Windows Settings > Scheduled Scans**, click **Administrator-Defined Scans**.
6. Select **Daily Scheduled Scan** and click **Edit**. The **Edit Scheduled Scan** window opens.
7. Change scan name and description to **Weekly Scheduled Scan** and **Weekly Scan at 00:00** respectively.
8. Select **Scan type** as **Full Scan**.
9. Select the **Schedule** tab.
10. Under **Scanning Schedule**, select **Weekly** and change the time to **00:00**.
11. Under **Scan Duration** uncheck **Randomize scan start time within this period (recommended in VMs)** and select **Scan until finished (recommended to optimize scan performance)**.
12. Under **Missed scheduled Scans** uncheck **Retry the scan within**.
13. Select the **Notifications** tab.
14. Uncheck **Display a notification message on the infected computer** and click **OK**.
15. Select the **Advanced** tab in the **Administrator-Defined Scans** window.
16. Under **Scheduled Scans** uncheck **Delay scheduled scans when running on batteries**, **Allow user-defined scheduled scans to run when scan author is not logged on** and **Display Notification about detections when the user logs on**.
17. Under **Startup and Triggered Scans** uncheck **Run an Active Scan when new definitions arrive**.
18. Under **Windows Settings > Protection Technology**, click **Auto-Protect**.



19. Select the **Scan Details** tab and select and lock **Enable Auto-Protect**.
20. Select the **Notifications** tab and uncheck and lock **Display a notification message on the infected computer** and **Display the Auto-Protect results dialog on the infected Computer**.
21. Select the **Advanced** tab and under **Auto-Protect Reloading and Enablement**, uncheck and lock the **When Auto-Protect is disabled, Enable after: option**.
22. Under **Additional Options** click **File Cache**. The **File Cache** window opens.
23. Uncheck **Rescan cache when new definitions load** and click **OK**.
24. Under **Windows Settings > Protection Technology**, click **Download Protection**.
25. Select the **Notifications** tab and uncheck and lock **Display a notification message on the infected computer**.
26. Under **Windows Settings > Email Scans**, click **Internet Email Auto-Protect**.
27. Select the **Notifications** tab and uncheck and lock **Display a notification message on the infected computer**, **Display a progress indicator when email is being sent**, and **Display a notification area icon**.
28. Under **Windows Settings > Email Scans**, click **Microsoft Outlook Auto-Protect**.
29. Select the **Notifications** tab and uncheck and lock **Display a notification message on the infected computer**.
30. Under **Windows Settings > Email Scans**, click **Lotus Notes Auto-Protect**.
31. Select the **Notifications** tab and uncheck and lock **Display a notification message on infected computer**.
32. Under **Windows Settings > Advanced Options**, click **Quarantine**.
33. Under **When New Virus Definitions Arrive**, select **Do nothing**.
34. Under **Windows Settings > Advanced Options**, click **Miscellaneous**.
35. Select the **Notifications** tab and uncheck **Display a notification message on the client computer** under the **When definitions are outdated** and **Display a notification message on the client computer** under **When Symantec Endpoint Protection is running without virus definitions** and **Display error messages with a URL to a solution**.
36. Click **OK** to close **Virus and Spyware Protection policies**.
37. Click **Yes** at the **Assign Policies** message box.
38. Select **My Company** and click **Assign**.
39. Click **Yes** at the message box.
40. Under **Policies** click **LiveUpdate**.
41. Select **LiveUpdate Settings policy** and under **Tasks**, click **Edit the policy**.
42. Under **Overview > Windows Settings**, click **Server Settings**.
43. Under **Internal or External LiveUpdate Server**, ensure **Use the default management server** is selected and uncheck **Use a LiveUpdate server**.



44. Click **OK**.
45. Click **Clients** from left pane and select the **Policies** tab.
46. Uncheck **Inherit policies and settings from parent group "My Company"** and click **Communications Settings** under **Location-Independent Policies and Settings**.
47. Under **Download**, make sure **Download policies and content from the management server** is checked and **Push mode** is selected.
48. Click **OK**.
49. **Click General Settings** under **Location-independent Policies and Settings**.
50. Select the **Tamper Protection** tab and uncheck and lock **Protect Symantec security software from being tampered with or shut down**.
51. Click **OK**.
52. Click **Admin** and select **Servers**.
53. Under **Servers**, select **Local Site (My Site)**.
54. Under **Tasks**, select **Edit Site Properties**. The **Site Properties for Locate Site (My Site)** window opens.
55. Select **LiveUpdate** tab and under **Download Schedule** ensure the schedule is set to **Every 4 hour(s)**.
56. Click **OK**.
57. Click **Policies** and select **Firewall**.
58. Select **Firewall policy** and under **Tasks** click **Edit the policy**. The **Firewall policy** window opens.
59. Click **Protection and Stealth** and under **Protection Settings** uncheck **Automatically block an attacker's IP address**.
- 60. Click OK.**
61. Click **Log Off** and close the Symantec EndPoint Protection Manager Console. Make sure Symantec Endpoint Protection Policies are pushed in client systems.

Symantec EndPoint Protection Post Installation Guidelines

1. Enable the Loopback Connection. Refer to Enable Loopback Connection for more information.
2. Configure the Computer Browser service. Refer to Configure Computer Browser Service After Anti-Virus Installation for more information.

McAfee VirusScan Enterprise v8.8 Patch 2

Installation Overview

McAfee VirusScan Enterprise should be installed on an individual Mac-Lab/CardioLab system and it should be managed individually. Use the following instructions to install and configure **McAfee VirusScan Enterprise** for English.

Virus updates are the responsibility of the facility. Update the definitions regularly to ensure that the latest virus protection is on the system.



McAfee VirusScan Enterprise Installation Procedure

1. Log on as Administrator or as a member of that group.
2. Insert the **McAfee VirusScan Enterprise 8.8 Patch 2 CD** into the CD drive.
3. Double-click **SetupVSE.Exe**. The Windows Defender dialog appears.
4. Click **Yes**. The **McAfee VirusScan Enterprise Setup** screen appears.
5. Click **Next**. The **McAfee End User License Agreement** screen appears.
6. Read the license agreement and complete any necessary fields, click **OK** when finished.
7. Select **Typical** and click **Next**.
8. Select **Standard Protection** and click **Next**.
9. Click **Install** and wait for the installation to complete. After successful installation of McAfee VirusScan Enterprise, the **McAfee Virus Scan Enterprise Setup has completed successfully** screen appears.
10. Uncheck the **Run On-Demand Scan** checkbox and click **Finish**.
11. If the **Update in Progress** window appears, click **Cancel**.
12. If a message box to restart the system appears, click **OK**.
13. Restart the system.
14. Log on as Administrator or as a member of that group.

McAfee VirusScan Enterprise Configuration

1. Right-click **McAfee** in the system tray and select **On-Access Scan Properties**.
2. Select the **ScriptScan** tab. The **ScriptScan** window opens.
3. Clear the **Enable scanning of scripts** check box.
4. Select the **Messages** tab. The **Messages** window opens.
5. Uncheck the **Show the messages dialog when a threat is detected and display the specified text in the message** check box.
6. Click **Apply**.
7. Click **OK** to close the **On-Access Scan Properties** window.
8. Select **Start > All Programs > McAfee > VirusScan Console**. The **VirusScan Console** window opens.
9. Select **Tools > Alerts**. The **Alert Properties** window opens.
10. Uncheck the **On-Access Scan, On-Demand Scan and scheduled scans, Email Scan and AutoUpdate** check boxes.
11. Click **Destination**. The **Alert Manager Client Configuration** window opens.
12. Select the **Disable alerting** check box.
13. Click **OK**. The **Alert Properties** window opens.



14. Select the **Additional Alerting Options** tab.
15. Select the **Suppress all alerts (severities 0 to 4)** option from the **Severity Filter** drop-down list.
16. Select the **Alert Manager Alerts** tab.
17. Clear the **Access Protection** check box.
18. Click **Apply**.
19. Click **OK** to close the **Alert Properties** window.
20. Right-click **AutoUpdate** on the VirusScan Console.
21. Click **Properties**. The McAfee **AutoUpdate Properties – AutoUpdate** window opens.
22. Click **Schedule**. The **Schedule Settings** window opens.
23. Clear the **Enable (scheduled task runs at specified time)** check box.
24. Click **Apply**.
25. Click **OK** to close the **Schedule Settings** window.
26. Click **OK** to close the **McAfee AutoUpdate Properties – AutoUpdate** window.
27. Right-click **Full Scan** on the VirusScan Console.
28. Click **Properties**. The **On Demand Scan Properties** window opens.
29. Click **Schedule**. The **Schedule Settings** window opens.
30. Check the **Enable (scheduled task runs at specified time)** check box.
31. Select the **Schedule** tab.
32. Select **Weekly** from **Run Task, 12:00 AM** from **Start Time**. check **Sunday** from **Schedule Task Weekly**.
33. Click **OK**. The **On Demand Scan Properties - Full Scan** window opens.
34. Select the **Exclusions** tab.
35. Click **Exclusions**. The Set **Exclusions** window opens.
36. Click **Add**.
37. Click **Browse** and navigate to **D:\GEMMS\Prucka** and **D:\GEData\Studies** folders one at a time and select the **Also exclude subfolders** checkbox.
38. Click **OK**.
39. In the **Set Exclusions** window, make sure the **D:\GEMMS\Prucka** and **D:\GEData\Studies** folders display.
40. Click **OK**.
41. Click **OK** to close the **On Demand Scan Properties – Full Scan** window.
42. Close the **VirusScan Console**.



McAfee ePolicy Orchestrator v5.0

Installation Overview

Install McAfee ePolicy Orchestrator on a networked Mac-Lab/CardioLab environment only. McAfee ePolicy Orchestrator must be installed on a customer-provided server and McAfee VirusScan Enterprise should be deployed to the Centricity Cardiology INW server and Acquisition/Review workstations as a client. Use the following instructions to install and configure **McAfee ePolicy Orchestrator** for English.

Virus updates are the responsibility of the facility. Update the definitions regularly to ensure that the latest virus protection is on the system.

Pre-Installation Guidelines

1. Make sure McAfee ePolicy Orchestrator Console is installed on the customer provided server before continuing with these steps.
2. On the customer-provided server, open inbound port **443** and name it **IN_443_McAfee** and allow inbound connections for domain, public, and private.
3. Log on as Administrator or a member of that group on all client systems (Acquisition, Review, and INW Server) to install the anti-virus software.
4. Disable the Loopback Connection. Refer to Disable Loopback Connection for more information.

McAfee ePolicy Orchestrator Deployment Steps (Preferred Push Installation Method)

1. Select **Start > All Programs > McAfee > ePolicy Orchestrator > Launch McAfee ePolicy Orchestrator 5.0.0 Console** to log on to the ePolicy Orchestrator console.

NOTE: Click **Continue with this website** if the **Security Alert** message box appears.

2. Enter the appropriate username and password and click **Log On**.
3. Select **Menu > Configuration > Server Settings > Port**.
4. Record the **Agent-to-server communication port** number.
5. On the customer-provided server, open the **Agent-to-server communication port** as an inbound port and name it **IN_<port number>_McAfee**.
6. Select **Start > All Programs > McAfee > ePolicy Orchestrator > Launch McAfee ePolicy Orchestrator 5.0.0 Console** to log on to the ePolicy Orchestrator console.

NOTE: Click **Continue with this website** if the **Security Alert** message box appears.

7. Enter the appropriate username and password and click **Log On**.
8. Select **Menu > Systems > Systems Tree**. The **System Tree** window opens.
9. Click **My Organization** and with the focus on **My Organization** click **System Tree Actions > New Systems** from the bottom left corner of the screen.
10. Select **Push agents and add systems to the current group (My Organization)** and click **Browse**.
11. Enter the **domain administrator** username and password and click **OK**.



12. Select the **INW** domain from the **Domain** drop-down list.
13. Select the client machines (Acquisition, Review, and INW Server) connected to the domain and click **OK**.
14. Select **Agent Version** as **McAfee Agent for Windows 4.8.0 (Current)**. Enter appropriate domain administrator username and password and click **OK**.
15. In client machines (Acquisition, Review, and INW Server) make sure the **C:\Program Files\McAfee\Common Framework** directory is present and **McAfee Agent** is installed in the same directory.

NOTE: For the INW Server make sure the **C:\Program Files (x86)\McAfee\Common Framework** directory is present and **McAfee Agent** is installed in the same directory.
16. Restart the client machines (Acquisition, Review, and INW Server).
17. Click **Start > All Programs > McAfee > ePolicy Orchestrator > Launch McAfee ePolicy Orchestrator 5.0.0 Console**.
18. Enter the appropriate username and password and click **Log On**.
19. Click **Menu > Systems > System Tree**.
20. Click **My Organization** and with the focus on **My Organization** click the **Assigned Client Tasks** tab.
21. Click **Actions > New Client Task Assignment** button at the bottom of the screen. The **Client Task Assignment Builder** screen appears.
22. Select the following:
 - a. **Product:** McAfee Agent
 - b. **Task Type:** Product Deployment
 - c. **Task name:** Create New Task
23. On the **Client Task Catalog: New Task- McAfee Agent: Product Deployment** screen, complete the fields as follows:
 - a. **Task Name:** Enter the appropriate task name
 - b. **Target platforms:** Windows
 - c. **Products and components:** VirusScan Enterprise
 - d. **Options:** Run at every policy enforcement (Windows only)
24. Click **Save**.
25. In the **1 select Task** screen, select the following:
 - a. **Product:** McAfee Agent
 - b. **Task Type:** Product Deployment
 - c. **Task Name:** Newly created task name



26. Click **Next**. The **2 Schedule** screen appears.
27. Select **Run immediately** from **Schedule type** drop-down list.
28. Click **Next**. The **3 Summary** screen appears.
29. Click **Save**. The **System Tree** screen appears.
30. Select the **Systems** tab and then select all the client machines (Acquisition, Review, and INW Server) which are connected to the domain.
31. Click **Wake up Agents** at bottom of the window.
32. Keep default settings and click **OK**.
33. Restart all the client machines (Acquisition, Review, and INW Server) and log in with Administrator or a member of that group on all client machines.
34. Click the **Log Off** link to close the **McAfee ePolicy Orchestrator Console**.

McAfee ePolicy Orchestrator Server Console Configuration

1. Select **Start > All Programs > McAfee > ePolicy Orchestrator > Launch McAfee ePolicy Orchestrator 5.0.0 Console** to log on to the ePolicy Orchestrator console.

NOTE: Click **Continue with this website** if the **Security Alert** message box appears.

2. Enter the appropriate Username and Password. The **ePO Summary** window opens.
3. Select **Menu > Systems > Systems Tree**. The **System Tree** window opens.
4. Click **My Organization**.
5. Select the **Assigned Policies** tab. The **Assigned Policies** screen opens.
6. From the **Product** drop-down list, select **VirusScan Enterprise 8.8.0**. The **Assigned Policies** window for VirusScan Enterprise 8.8.0 opens.
7. Click **My Default** for **On-Access General Policies**. The **General** window opens.
8. Select **Workstation** from the **Settings for** drop-down list. Click **ScriptScan** and uncheck **Enable scanning of scripts**.
9. Click **Messages**. The **Messages** window opens.
10. Uncheck **Show the messages dialog box when a threat is detected and display the specified text in the message**.
11. Select **Server** from the **Settings for** drop-down list.
12. Click **ScriptScan** and ensure **Enable scanning of scripts** is unchecked.
13. Click **Messages**. The **Messages** window opens.
14. Uncheck the **Show the messages dialog box when a threat is detected and display the specified text in the message**.
15. Click **Save**.
16. Select **My Default** for **Buffer Overflow Protection Policies**. The **Buffer Overflow Protection** window opens.



17. Select **Workstation** from the **Settings for** drop-down list and **uncheck Show the messages dialog box when a buffer overflow is detected.**
18. Select **Server** from the **Settings for** drop-down list and uncheck **Show the messages dialog box when a buffer overflow is detected.**
19. Click **Save.**
20. Click **My Default** for **Alert Policies.** The **Alert Manager Alerts** window opens.
21. Select **Workstation** from the **Settings for** drop-down list and uncheck **On-Access Scan, On- Demand Scan and scheduled scans, Email Scan** and **AutoUpdate.**
22. Check **Disable alerting.**
23. Click **Additional Alerting Options.** The **Additional Alerting Options** window opens.
24. From the **Severity Filters** drop-down menu, select **Suppress all alerts (severities 0 to 4).**
25. Select **Server** from the **Settings for** drop-down list and select the **Alert Manager Alerts** tab. The **Alert Manager Alerts** window opens.
26. Uncheck **On-Access Scan, On-Demand Scan and scheduled scans, Email Scan** and **AutoUpdate.**
27. Check **Disable alerting.**
28. Click **Additional Alerting Options.** The **Additional Alerting Options** window opens.
29. From the **Severity Filters** drop-down menu, select **Suppress all alerts (severities 0 to 4).**
30. Click **Save.**
31. Click **My Default** for **On-Access Default Processes Policies.** The **Processes** window opens.
32. Select **Workstation** from the **Settings for** drop-down list and ensure **Configure one scanning policy for all processes** is selected.
33. Click the **Exclusions** tab. The **Exclusions** window opens.
34. Click **Add** and select **By pattern.**
35. Enter the **D:\GEMMS\Prucka** and **D:\GEData\Studies** folder names and select **Also exclude subfolders.**
36. Click **OK.**
37. Select **Server** from **Settings for** drop-down list and select the **Processes** tab.
38. Ensure **Configure one scanning policy for all processes** is selected.
39. Select the **Exclusions** tab. The **Exclusions** window opens.
40. Click **Add** and select **By pattern.**
41. Enter **D:\GEMMS\Prucka** and select **Also exclude subfolders.**
42. Click **OK.**
43. Click **Save.**



44. From the **Product** drop-down menu, select **McAfee Agent**. The **Policies** window for McAfee Agent opens.
45. Click **My Default** for **Repository**. The **Repositories** window opens.
NOTE: Click **Close** for **Internet Explorer Security Message Box**.
46. Click **Proxy**. The **Proxy** window opens.
47. Select **Use Internet Explorer settings (For Windows)/System Preferences settings (For Mac OSX)**.
48. Click **Save**.
49. Click **Systems**.
50. Select all the client systems (Acquisition, Review and Centricity Cardiology INW server) into which the configured policies are to be deployed.
51. Select **Wake Up Agents**. The **Wake Up Agent** window opens.
52. Click **OK**.
53. Log off ePolicy Orchestrator.

McAfee ePolicy Orchestrator Post Installation

Enable the Loopback Connection. Refer to Enable Loopback Connection for more information.

Trend Micro OfficeScan Client/Server Edition v10.6 SP2

Installation Overview

Install Trend Micro OfficeScan Client/Server Edition on a networked Mac-Lab/CardioLab environment only. Trend Micro OfficeScan must be installed on the customer-provided server and then deployed to Centricity Cardiology INW server and Acquisition/Review workstation as clients. Use the following instructions to install **Trend Micro OfficeScan Client/Server Edition** for English.

Virus updates are the responsibility of the facility. Update the definitions regularly to ensure that the latest virus protection is on the system.

Pre-Installation Guidelines

1. On the customer-provided server, make sure Office Scan Web Console is installed before continuing with these steps.
2. On the customer-provided server, open inbound port **8080** and name it **IN_8080_OSCE** and allow inbound connections for domain, public, and private.
3. Log on as Administrator or member of that group on all client systems (Acquisition, Review, and INW Server) to install the anti-virus software.
4. Disable the Loopback Connection. Refer to Disable Loopback Connection for more information.
5. Configure the Computer Browser service. Refer to Configure Computer Browser Service Before Anti-Virus Installation for more information.



Trend Micro OfficeScan Deployment Steps (Preferred Push Installation Method)

1. Click **Start > All Programs > TrendMicro OfficeScan server - <server name> > Office Scan Web Console.**

NOTE: Continue by selecting **Continue to this website (not recommended)**. In the **Security Alert** window, check **In the future, do not show this warning** and click **OK**.

2. If you receive a certificate error indicating that the site is not trusted, manage your certificates to include Trend Micro OfficeScan.
3. If prompted, install the **AtxEnc** add-ons. The **Security Warning** screen will appear.
4. Click **Install**.
5. Enter appropriate username and password and click **Log On**.
6. If prompted, click **Update Now** to install new widgets. Wait until the new widgets are updated. The update is completed screen will appear.
7. Click **OK**.
8. From the left side menu bar, click **Networked Computers > Client Installation > Remote link**.
9. If prompted, install the **AtxConsole** add-ons. The **Security Warning** screen will appear.
10. Click **Install**.
11. Double-click **My Company** in the **Remote Installation** window. All domains will be listed under **My Company**.
12. Expand the appropriate domain (Example: INW) from the list. All systems connected to the domain appear.
13. If domains or systems are not listed in the **Domain and Computers** window, do the following on each of the client systems (Acquisition, Review, and INW Server):
 - a. Enter the client system IP address in "Search for Computers:"
 - b. Click the **Search** button.
 - c. Enter the appropriate <domain name>\username and password and click **Log on**.
 - d. Entered IP addresses will appear in the **Selected Computers** list
14. Select the client machines (Acquisition, Review, and INW Server) one at a time from the **Selected Computers** pane and click **Install**.
15. Click **Yes** at the confirmation box.
16. Click **OK** at the **Number of clients to which notifications were sent** message box.
17. Restart all the client machines (Acquisition, Review, and INW Server) and Log in as Administrator or a member of that group on all client machines and wait until the Trend Micro OfficeScan icon in system tray changes to blue with sin wave symbol.
18. Click the **Log Off** link to close the **OfficeScan Web Console**.



Trend Micro OfficeScan Server Console Configuration

1. Select **Start > All Programs > TrendMicro Office Scan server <servername> > Office Scan Web Console**. The **Trend Micro OfficeScan Login** window appears.
2. Enter the appropriate user name and password and click **Login**. The **Summary** window opens.
3. From the left side pane, select the **Networked Computers > Client Management** link.
4. On the right side, select **OfficeScan Server**.
5. From the **Settings** options, select **Privileges and Other Settings**.
6. Select only the following options in the **Privileges** tab and clear the remaining options:
 - **Scan Privileges > Configure Manual Scan Settings**.
 - **Scan Privileges > Configure Real-time Scan Settings**.
 - **Scan Privileges > Configure Scheduled Scan Settings**.
 - **Proxy Setting Privileges > Allow the client user to configure proxy settings**.
 - **Uninstallation > Require a password for the user to uninstall the OfficeScan Client**. Enter a suitable password.
 - **Unloading > Require a password for the user to unload the OfficeScan client**. Enter a suitable password.
7. Select the **Other Settings** tab.
8. Select **Client Security Settings > Normal** and clear the remaining options.
9. Click **Apply to All Clients**.
10. Click **Close** to close the **Privileges and Other Settings** window.
11. From the left side pane, select the **Client Management** link.
12. On the right side, select **OfficeScan Server**.
13. From the **Settings** options, select **Scan Settings > Scan Now Settings**.
14. Select only the following options in the **Target** tab and clear the remaining options:
 - **Files to Scan > File types scanned by IntelliScan**.
 - **Scan Settings > Scan Compressed files**.
 - **Scan Settings > Scan OLE objects**.
 - **Virus/Malware Scan Settings only > Scan boot area**.
 - **CPU Usage > Low**.
 - **Scan Exclusion > Enable scan exclusion**.
 - **Scan Exclusion > Apply scan exclusion settings to all scan types**.
 - **Scan Exclusion List (Directories) > Exclude directories where Trend Micro products are installed** and select "Add path to client Computers Exclusion list".
 - Enter the **D:\GEMMS\Prucka** and **D:\GEData\Studies** folders one at a time in the directory path for **Exclusion List** and click **Add**.
15. Click **Apply to All Clients**.
16. The exclusion list on this screen will replace the exclusion list on the clients or domains you selected in the client tree earlier. Do you want to proceed? Message will appear. Click **OK**.



17. Click **Close** to close the **Scan Now Settings** page.
18. From the left side pane, select the **Client Management** link.
19. On the right side, select **OfficeScan Server**.
20. From the **Settings** options, select **Settings->Scan Settings->Real-time Scan Settings**.
21. Select the **Target** tab. Select only the following options and clear the remaining options:
 - **Enable Virus/Malware scan.**
 - **User Activity on Files > created/modified and retrieved.**
 - **Files to Scan > File types scanned by IntelliScan.**
 - **Virus/Malware Scan Settings Only > Enable IntelliTrap.**
 - **Scan Exclusion > Enable scan exclusion.**
 - **Scan Exclusion > Apply scan exclusion settings to all scan types.**
 - **Scan Exclusion List (Directories) > Exclude directories where Trend Micro products are installed.**
 - Make sure the **D:\GEMMS\Prucka** and **D:\GEData\Studies** paths are present in the **Exclusion List**.
22. Click the **Action** tab.
23. Keep the default settings and clear **Virus/Malware > Display a notification message on the client computer when Virus/Malware is detected** option and **Spyware/Grayware->Display a notification message on the client computer when spyware/Grayware is detected**.
24. Click **Apply to All Clients**.
25. Click **Close** to close the **Real-time Scan Settings** page.
26. From the left side pane, select the **Client Management** link.
27. On the right side, select **OfficeScan Server**.
28. From the **Settings** options, select **Scan settings-> Scheduled Scan Settings**.
29. Select the **Target** tab. Select only the following options and clear the remaining options:
 - **Enable Virus/Malware scan.**
 - **Schedule > Weekly, every (Sunday).**
 - **Files to Scan > File types scanned by IntelliScan.**
 - **Virus/Malware Scan settings only > Scan boot area.**
 - **CPU Usage > Low.**
 - **Scan Exclusion > Enable scan exclusion.**
 - **Scan Exclusion > Apply scan exclusion settings to all scan types.**
 - **Scan Exclusion List (Directories) > Exclude directories where Trend Micro products are installed.**
 - Make sure the **D:\GEMMS\Prucka** and **D:\GEData\Studies** paths are present in the **Exclusion List**.
30. Click the **Action** tab.
31. Keep the default settings and uncheck the **Virus/Malware > Display a notification message on the client computer when Virus/Malware is detected** and **Spyware/Grayware-> Display a notification message on the client computer when spyware/Grayware is detected** options.
32. Click **Apply to All Clients**.



33. Click **Close** to close the **Scheduled Scan Settings** page.
34. From the left side pane, select the **Networked Computers > Global Client Settings** link.
35. Select only the following options and clear the remaining options:
 - **Scan Settings > Configure Scan settings for large compressed files.**
 - **Scan Settings > Do not scan files in the compressed file if the size exceeds 2 MB.**
 - **Scan Settings > In a compressed file scan only the first 100 files.**
 - **Scan Settings > Exclude the OfficeScan server database folder from Real-time Scan.**
 - **Scan Settings > Exclude Microsoft Exchange server folders and files from scans.**
 - **Reserved Disk Space > Reserve 60 MB of disk space for updates.**
 - **Proxy Configuration > Automatically detect settings.**

NOTE: It is important to clear the Alert Settings > Display a notification message if the client computer needs to restart to load a kernel driver.

36. Click **Save**.
37. Click **Log off** and close the **OfficeScan Web Console**.

Trend Micro OfficeScan Post Installation Guidelines

1. Enable the Loopback Connection. Refer to Enable Loopback Connection for more information.
2. Configure the Computer Browser service. Refer to Configure Computer Browser Service After Anti-Virus Installation on for more information.

Trend Micro OfficeScan Client/Server Edition XG 12.0

Installation Overview

Install Trend Micro OfficeScan Client/Server Edition on a networked Mac-Lab/CardioLab environment only. Trend Micro OfficeScan must be installed on the customer-provided server and then deployed to Centricity Cardiology INW server and Acquisition/Review workstation as clients. Use the following instructions to install **Trend Micro OfficeScan Client/Server Edition**.

Virus updates are the responsibility of the facility. Update the definitions regularly to ensure that the latest virus protection is on the system.

Pre-Installation Guidelines

Note: Internet Explorer 10 is the minimum IE browser required to run OfficeScan Web Console on customer-provided server.

1. On the customer-provided server, make sure Office Scan Web Console is installed before continuing with these steps.
2. On the customer-provided server, open inbound port **8080** and name it **IN_8080_OSCE** and allow inbound connections for domain, public, and private.
3. During installation of Trend Micro OfficeScan do the following on the customer-provided server:
 - a. Uncheck **Enable firewall** in the **Anti-virus Feature** window.



- b. Select **No, Please do not enable assessment mode** in the **Anti-spyware Feature** window.
- c. Uncheck **Enable web reputation policy** in the **Web Reputation Feature** window.
- 4. Log on as **Administrator** or a member of that group on all client systems (Acquisition, Review, and INW Server) to install the anti-virus software.
- 5. Disable the Loopback Connection, if applicable. See “Disable Loopback Connection” for more information.
- 6. Configure the Computer Browser service. See “Configure Computer Browser Service Before Anti-Virus Installation” for more information.
- 7. One at a time, install the required root and intermediate certificates.
 - a. Double-click the first certificate to install it on the MLCL systems.
 - b. Open the certificate and click **Install Certificate**.
 - c. Click **Next** when the **Certificate Import Wizard** appears.
 - d. On the **Certificate Store** window, select **Place all certificates in the following store** and click **Browse**.
 - e. Check **Show physical stores > Trusted Root Certification Authorities > Local Computer** and then click **OK**.
 - f. Click **Next** on **Certificate Import Wizard**.
 - g. Click **Finish**. The import was successful message should appear.
 - h. Repeat until all required certificates are installed.

Required Root and Intermediate Level Certificates
AddTrustExternalCARoot.crt
COMODOCodeSigningCA2.crt
UTNAddTrustObject_CA.crt
UTN-USERFirst-Object.crt
UTN-USERFirst-Object_kmod.crt

Note: Each certificate has an expiry date. Once the certificate has expired, they should be renewed and updated on the MLCL systems to ensure that the OfficeScan agent functions as expected.

Trend Micro OfficeScan Deployment Steps (Preferred Push Installation Method)

- 1. Click **Start > All Programs > TrendMicro OfficeScan server - <server name> > Office Scan Web Console**.

NOTE: Continue by selecting **Continue to this website (not recommended)**. In the **Security Alert** window, check **In the future, do not show this warning** and click **OK**.



2. If you receive a certificate error indicating that the site is not trusted, manage your certificates to include Trend Micro OfficeScan.
3. If prompted, install the **AtxEnc** add-ons. The **Security Warning** screen will appear.
4. Click **Install**.
5. Enter appropriate username and password and click **Log On**.
6. If prompted, click **Update Now** to install new widgets. Wait until the new widgets are updated. The update is completed screen will appear.
7. Click **OK**.
8. From the top menu bar, click **Agents > Agent Installation > Remote**.
9. If prompted, install the **AtxConsole** add-ons. The Security Warning screen displays.
10. Click **Install**.
11. Double-click **My Company** in the **Remote Installation** window. All domains will be listed under **OfficeScan Server**.
12. Double-click the domain (Example: INW) from the list. All systems connected to the domain appear.
13. If domains or systems are not listed in the **Domains and Endpoints** window, perform the steps listed in either a or b:
 - a. On the OfficeScan web console:
 - i. Enter the IP of each of the client systems in the **Search for endpoints** box one at a time and press Enter.
 - ii. Provide <domain name>\username and password and click **Log on**.
 - b. On each of the client systems (Acquisition, Review, and INW Server):
 - i. Log in as **Administrator** or a member of that group on all client machines.
 - ii. Click **Start > Run**.
 - iii. Type \\<Anti-Virus Management Console_server_IP_address> and press **Enter**.
 - iv. Enter the administrator username and password, when prompted.
 - v. Navigate to \\<Anti-Virus Management Console_server_IP_address>\ofsscan and double-click **AutoPcc.exe**. Enter the administrator username and password, when prompted.
 - vi. Restart the client systems when the installation is complete.
 - vii. Log in as **Administrator** or a member of that group on all client machines and wait until the Trend Micro OfficeScan icon in system tray changes to blue.
 - viii. Skip the remaining steps in this procedure and go to "Trend Micro OfficeScan Server Console Configuration."
14. Select the client machines (Acquisition, Review, and INW Server) and click **Add**.



15. Type the <domain name>\username and password and click **Log on**.
16. Select the client machines (Acquisition, Review, and INW Server) one at a time from the **Selected Computers** pane and click **Install**.
17. Click **Yes** at the confirmation box.
18. Click **OK** at the **Number of agents to which notifications were sent** message box.
19. Restart all the client machines (Acquisition, Review, and INW Server) and Log in as Administrator or a member of that group on all client machines and wait until the Trend Micro OfficeScan icon in system tray changes to blue with a green tick mark symbol.
20. Click the **Log Off** link to close the **OfficeScan Web Console**.

Trend Micro OfficeScan Server Console Configuration

1. Select **Start > All Programs > TrendMicro Office Scan server <servername> > Office Scan Web Console**. The **Trend Micro OfficeScan Login** window appears.
2. Enter the appropriate user name and password and click **Login**. The **Summary** window opens.
3. From the top pane, select **Agents > Agent Management**.
4. On the left side, select **OfficeScan Server**.
5. From the **Settings** options, select **Scan Settings > Manual Scan Settings**. The **Manual Scan Settings** screen appears.
6. Click the **Target** tab and select only the following options. Uncheck all other options.

Field	Selected Value
Files to Scan	File types scanned by IntelliScan
Scan Settings	Scan compressed files
Scan Settings	Scan OLE objects
Virus/Malware Scan Settings Only	Scan boot area
CPU Usage	Low

7. Click the **Scan Exclusion** tab and select only the following options. Uncheck all other options.

Field	Selected Value
Scan Exclusion	Enable scan exclusion
Scan Exclusion	Apply scan exclusion settings to all scan types



Scan Exclusion List (Directories)	Exclude directories where Trend Micro products are installed. Select Add path to agent Computers Exclusion list .
-----------------------------------	--

8. Select **Adds path** from the drop-down under **Saving the OfficeScan agent's exclusion list**.
 - a. Type **D:\GEMMS\Prucka** and click **Add**.
 - b. Type **D:\GEData\Studies** and click **Add**.
9. Click **Apply to All Agents**.
10. Click **OK** at the **The exclusion list on this screen will replace the exclusion list on the agents or domains you selected in the client tree earlier. Do you want to proceed?** message.
11. Click **Close** to close the **Manual Scan Settings** screen.
12. From the top pane, select **Agent > Agent Management**.
13. On the left side, select OfficeScan Server.
14. From the **Settings** options, select **Scan Settings > Real-time Scan Settings**. The **Realtime Scan Settings** screen appears.
15. Click the **Target** tab and select only the following options. Uncheck all other options.

Field	Selected Value
Real-Time Scan Settings	Enable virus/malware scan
Real-Time Scan Settings	Enable spyware/grayware scan
Files to Scan	File types scanned by IntelliScan
Scan Settings	Scan compressed files
Scan Settings	Scan OLE objects
Virus/Malware Scan Settings Only	Enable IntelliTrap

16. Click the **Scan Exclusion** tab and select only the following options. Uncheck all other options.

Field	Selected Value
Scan Exclusion	Enable scan exclusion
Scan Exclusion	Apply scan exclusion settings to all scan types
Scan Exclusion List (Directories)	Exclude directories where Trend Micro products are installed.

17. Confirm that **D:\GEMMS\Prucka** and **D:\GEData\Studies** are included in the Exclusion List.



18. Click the **Action** tab.

19. Keep the default settings and uncheck the following options:

Field	Unselected Value
Virus/Malware	Display a notification message on endpoints when virus/malware is detected
Spyware/Grayware	Display a notification message on endpoints when spyware/grayware is detected

20. Click **Apply to All Agents**.

21. Click **Close** to close the **Real-time Scan Settings** screen.

22. From the top pane, select **Agents > Agent Management**.

23. On the left side, select **OfficeScan Server**.

24. From **Settings**, select **Scan Settings > Scheduled Scan Settings**. The **Scheduled Scan Settings** screen appears.

25. Click the **Target** tab and select only the following options. Uncheck all other options.

Field	Selected Value
Scheduled Scan Settings	Enable virus/malware scan
Scheduled Scan Settings	Enable spyware/grayware scan
Schedule	Weekly, every Sunday, Start time: 00:00 hh:mm
Files to Scan	File types scanned by IntelliScan
Scan Settings	Scan compressed files
Scan Settings	Scan OLE objects
Virus/Malware Scan Settings Only	Scan boot area
CPU Usage	Low

26. Click the **Scan Exclusion** tab and select only the following options. Uncheck all other options.

Field	Selected Value
-------	----------------



Scan Exclusion	Enable scan exclusion
Scan Exclusion	Apply scan exclusion settings to all scan types
Scan Exclusion List (Directories)	Exclude directories where Trend Micro products are installed.

27. Confirm that **D:\GEMMS\Prucka** and **D:\GEData\Studies** are included in the Exclusion List.

28. Click the **Action** tab.

29. Keep the default settings and uncheck the following options:

Field	Unselected Value
Virus/Malware	Display a notification message on endpoints when virus/malware is detected
Spyware/Grayware	Display a notification message on endpoints when spyware/grayware is detected

30. Click **Apply to All Agents**.

31. Click **Close** to close the **Scheduled Scan Settings** screen.

32. From the top pane, select **Agents > Agent Management**.

33. On the left side, select **OfficeScan Server**.

34. From **Settings**, select **Scan Settings > Scan Now Settings**. The **Scan Now Settings** screen appears.

35. Click the **Target** tab and select only the following options. Uncheck all other options.

Field	Selected Value
Scan Now Settings	Enable virus/malware scan
Scan Now Settings	Enable spyware/grayware scan
Files to Scan	File types scanned by IntelliScan
Scan Settings	Scan compressed files
Scan Settings	Scan OLE objects
Virus/Malware Scan Settings Only	Scan boot area
CPU Usage	Low



36. Click the **Scan Exclusion** tab and select only the following options. Uncheck all other options.

Field	Selected Value
Scan Exclusion	Enable scan exclusion
Scan Exclusion	Apply scan exclusion settings to all scan types
Scan Exclusion List (Directories)	Exclude directories where Trend Micro products are installed.

37. Confirm that **D:\GEMMS\Prucka** and **D:\GEData\Studies** are included in the Exclusion List.

38. Click **Apply to All Agents**.

39. Click **Close** to close the **Scan Now Settings** screen.

40. From the top pane, select **Agents > Agent Management**.

41. On the left side, select **OfficeScan Server**.

42. From **Settings**, select **Web Reputation Settings**. The **Web Reputation Settings** screen appears.

43. Click the **External Agents** tab and uncheck **Enable Web reputation policy on the following operating systems**, if selected already during installation.

44. Click the **Internal Agents** tab and uncheck **Enable Web reputation policy on the following operating systems**, if selected already during installation.

45. Click **Apply to All Agents**.

46. Click **Close** to close the **Web Reputation** screen.

47. From the top pane, select **Agents > Agent Management**.

48. On the left side, select **OfficeScan Server**.

49. From **Settings**, select **Behavior Monitoring Settings**. The **Behavior Monitoring Settings** screen appears.

50. Uncheck **Enable Malware Behavior Blocking** and **Enable Event Monitoring**.

51. Click **Apply to All Agents**.

52. Click **Close** to close the **Behavior Monitoring** screen.

53. From the top pane, select **Agents > Agent Management**.

54. On the left side, select **OfficeScan Server**.



55. From **Settings**, select **Device Control Settings**. The **Device Control Settings** screen appears.

56. Click the **External Agents** tab and uncheck the following options:

Field	Unselected Value
Notification	Display a notification message on endpoints when OfficeScan detects unauthorized device access
NA	Block the AutoRun function on USB storage
NA	Enable Device Control

57. Click the **Internal Agents** tab. Uncheck the same options listed in step 56.

58. Click **Apply to All Agents**.

59. Click **Close** to close the **Device Control Settings** screen.

60. From **Settings**, select **Device Control Settings**. The **Device Control Settings** screen reopens.

61. Click the **External Agents** tab and uncheck **Enable Device Control**.

62. Click the **Internal Agents** tab and uncheck **Enable Device Control**.

63. Click **Apply to All Agents**.

64. Click **Close** to close the **Device Control Settings** screen.

65. From the top pane, select **Agents > Agent Management**.

66. On the left side, select **OfficeScan Server**.

67. From **Settings**, select **Privileges and Other Settings**.

68. Click the **Privileges** tab and select only the following options. Uncheck all other options.

Field	Selected Value/Operation
Scan Privileges	Configure Manual Scan Settings
Scan Privileges	Configure Real-time Scan Settings
Scan Privileges	Configure Scheduled Scan Settings
Proxy Setting Privileges	Allow the agent user to configure proxy settings



Uninstallation	Requires a password. Enter a suitable password and confirm password.
Unload and Unlock	Requires a password. Enter a suitable password and confirm password.

69. Click the **Other Settings** tab. Uncheck all options.

Note: It is important to clear the following options:

Field	Option to Clear
OfficeScan Agent Self-protection	Protect OfficeScan agent services
OfficeScan Agent Self-protection	Protect files in the OfficeScan agent installation folder
OfficeScan Agent Self-protection	Protect OfficeScan agent registry keys
OfficeScan Agent Self-protection	Protect OfficeScan agent processes

70. Click **Apply to All Agents**.

71. Click **Close** to close the **Privileges and Other Settings** screen.

72. From the top pane, select **Agents > Agent Management**.

73. On the left side, select **OfficeScan Server**.

74. From **Settings**, select **Additional Service Settings**.

75. Uncheck **Enable service on the following operating systems**.

76. Click **Apply to All Agents**.

77. Click **Close** to close the **Additional Service Settings** screen.

78. From the top pane, select **Agents > Global Agent Settings**.

79. Select only the following options. Uncheck all other options.

Field	Selected Value
Scan Settings for Large Compressed Files	Do not scan files in the compressed file if the size exceeds 2 MB. Follow this for Real-Time Scan and Manual Scan/Schedule Scan/Scan Now.
Scan Settings for Large Compressed Files	In a compressed file scan only the first 100 files. Follow this for Real-Time Scan and Manual Scan/Schedule Scan/Scan Now.



Scan Settings	Exclude the OfficeScan server database folder from Real-time Scan.
Scan Settings	Exclude Microsoft Exchange server folders and files from scans.

- 80. Click Save.
- 81. From the top pane, select **Updates > Agents > Manual Updates**.
- 82. Select **Manually select agents** and click **Select**.
- 83. Double-click the appropriate domain name under **OfficeScan Server**.
- 84. Select client system one at a time and click **Initiate Update**.
- 85. Click **OK** when prompted.
- 86. Click **Log off** and close the **OfficeScan Web Console**.

Trend Micro OfficeScan Post Installation Guidelines

- 1. Enable the Loopback Connection. See "Enable Loopback Connection" for more information.
- 2. Configure the Computer Browser service. See "Configure Computer Browser Service After Anti-Virus Installation" for more information.