

g

GE Medical Systems

Technical Publications

Document # **MM10001**
Revision **1.0**

GE PACS Broker v1.5.2 (Modality Worklist Interfaces) **CONFORMANCE STATEMENT** **for DICOM V3.0**

Copyright© 1998 By General Electric Co.

Do not duplicate

g

GE Medical Systems

Author	Date	Reason for Change	Rev
S. Douglas	1/18/99	Initial release (Mitra v 4.7)	1.0

Introduction

OVERVIEW

This DICOM Conformance Statement is divided into Sections as described below:

Section 1 (Introduction), which describes the overall structure, intent, and references for this Conformance Statement.

Section 2 (OEM Conformance Statement), is an attachment of the original OEM DICOM Conformance Statement.

This document specifies the DICOM v3.0 implementation. It is entitled:

GE PACS Broker v1.5.2 (Modality Worklist Interfaces)

Conformance Statement for DICOM v3.0

Document # MM10001

This DICOM Conformance Statement documents the DICOM v3.0 Conformance Statement and Technical Specification required to interoperate with the GEMS network interface. Introductory information, which is applicable to all GEMS Conformance Statements, is described in the document:

Introduction to the Integrated DICOM/Network v3.0 (ID/Net v3.0)

Conformance Statement

Direction: 2118780.

This Introduction familiarizes the reader with DICOM terminology and general concepts. It should be read prior to reading the individual products' GEMS Conformance Statements.

The OEM Conformance Statement, contained in this document, also specifies the Lower Layer communications which it supports (e.g., TCP/IP). However, the Technical Specifications are defined in the DICOM v3.0 Part 8 standard.

For more information including Network Architecture and basic DICOM concepts, please refer to the Introduction.

For more information regarding DICOM v3.0, copies of the Standard may be obtained by written request or phone by contacting:

NEMA Publication
1300 North 17th Street
Suite 1847
Rosslyn, VA 22209
USA
Phone: (703) 841-3200

INTENDED AUDIENCE

The reader of this document is concerned with software design and/or system integration issues. It is assumed that the reader of this document is familiar with the DICOM v3.0 Standards and with the terminology and concepts which are used in those Standards.

If readers are unfamiliar with DICOM v3.0 terminology they should first refer to the document listed below, then read the DICOM v3.0 Standard itself, prior to reading this DICOM Conformance Statement document.

Introduction to the Integrated DICOM/Network v3.0 (ID/Net v3.0)

Conformance Statement

Direction: 2118780

SCOPE AND FIELD OF APPLICATION

It is the intent of this document, in conjunction with the *Introduction to the Integrated DICOM/Network v3.0 (ID/Net v3.0) Conformance Statement, Direction: 2118780*, to provide an unambiguous specification for this OEM DICOM implementation. This specification, called a Conformance Statement, includes a

Document # MM10001, Rev 1.0

DICOM v3.0 Conformance Statement and is necessary to ensure proper processing and interpretation of GEMS medical data exchanged using DICOM v3.0. The GEMS Conformance Statements are available to the public.

The reader of this DICOM Conformance Statement should be aware that different GEMS devices are capable of using different Information Object Definitions. For example, a GEMS CT Scanner may send images using the CT Information Object, MR Information Object, Secondary Capture Object, etc.

If the user encounters unspecified private data elements while parsing a DICOM Data Set, the user is well advised to ignore those data elements (per the DICOM v3.0 standard). Unspecified private data element information is subject to change without notice. If, however, the device is acting as a "full fidelity storage device", it should retain and re-transmit all of the private data elements which are sent by these devices.

IMPORTANT REMARKS

The use of these DICOM Conformance Statements, in conjunction with the DICOM v3.0 Standards, is intended to facilitate communication with GE imaging equipment. However, **by itself, it is not sufficient to ensure that inter-operation will be successful.** The **user (or user's agent)** needs to proceed with caution and address at least four issues:

- **Integration** - The integration of any device into an overall system of interconnected devices goes beyond the scope of standards (DICOM v3.0), and of this introduction and associated DICOM Conformance Statements when interoperability with non-GE equipment is desired. The responsibility to analyze the applications requirements and to design a solution that integrates GE imaging equipment with non-GE systems is the **user's** responsibility and should not be underestimated. The **user** is strongly advised to ensure that such an integration analysis is correctly performed.
- **Validation** - Testing the complete range of possible interactions between any GE device and non-GE devices, before the connection is declared operational, should not be overlooked. Therefore, the **user** should ensure that any non-GE provider accepts full responsibility for all validation required for their connection with GE devices. This includes the accuracy of the image data once it has crossed the interface between the GE imaging equipment and the non-GE device and the stability of the image data for the intended applications.

Such a validation is required before any clinical use (diagnosis and/or treatment) is performed. It applies when images acquired on GE imaging equipment are processed/displayed on a non-GE device, as well as when images acquired on non-GE equipment is processed/displayed on a GE console or workstation.

- **Future Evolution** - GE understands that the DICOM Standard will evolve to meet the user's growing requirements. GE is actively involved in the development of the DICOM v3.0 Standard. DICOM v3.0 will incorporate new features and technologies and GE may follow the evolution of the Standard. The GEMS protocol is based on DICOM v3.0 as specified in each DICOM Conformance Statement. Evolution of the Standard may require changes to devices which have implemented DICOM v3.0. **In addition, GE reserves the right to discontinue or make changes to the support of communications features (on its products) reflected on by these DICOM Conformance Statements.** The **user** should ensure that any non-GE provider, which connects with GE devices, also plans for the future evolution of the DICOM Standard. Failure to do so will likely result in the loss of function and/or connectivity as the DICOM Standard changes and GE Products are enhanced to support these changes.
- **To be informed of the evolution of the implementation described in this document, the User is advised to regularly check the GE Internet Server, accessible via anonymous ftp (GE Internet Server Address: ftp.med.ge.com, 192.88.230.11).**
- **Interaction** - It is the sole responsibility of the **non-GE provider** to ensure that communication with the interfaced equipment does

not cause degradation of GE imaging equipment performance and/or function.

REFERENCES

A list of references which is applicable to all GEMS Conformance Statements is included in the *Introduction to the Integrated DICOM/Network v3.0 (ID/Net v3.0) Conformance Statement, Direction: 2118780.*

DEFINITIONS

A set of definitions which is applicable to all GEMS Conformance Statements is included in the *Introduction to the Integrated DICOM/Network v3.0 (ID/Net v3.0) Conformance Statement, Direction: 2118780.*

SYMBOLS AND ABBREVIATIONS

A list of symbols and abbreviations which is applicable to all GEMS Conformance Statements is included in the *Introduction to the Integrated DICOM/Network v3.0 (ID/Net v3.0) Conformance Statement, Direction: 2118780.*

Original Equipment Manufacturers DICOM Conformance Statement

The sections that follow are taken from the Original Equipment Manufacturers DICOM Conformance Statement for this product.

Note: These sections are provided as reference and are subject to change without notice.

LIMITATIONS

- The GE PACS Broker is a general purpose product that provides a wide range of DICOM services, some of which are not needed nor supported in the GE PACS environment.
- The GE PACS Broker is used to integrate GE PACS with a number of Radiology Information Systems, not documented here, as well as to integrate modalities supporting the DICOM Modality Worklist SOP Class as an SCU. Therefore, the only DICOM service classes supported by the GE PACS environment through the use of the GE PACS Broker is Modality Worklist SOP Class as an SCP.
- The GE PACS Broker may also be used in a standalone configuration, when a GE PACS system is not present, to integrate modalities supporting DICOM Modality Worklist SCU with Radiology Information Systems (RIS) that do not offer native support for DICOM Modality Worklist SCP.
- The Extensions/Specialization provided by the GE PACS Broker includes 6 Private SOP Classes that are not defined by the DICOM Standard but a proprietary definition of the OEM. GE Medical Systems does not support these SOP Classes and does not encourage their use.

Table of Contents

Introduction.....	4
OVERVIEW	4
INTENDED AUDIENCE	4
SCOPE AND FIELD OF APPLICATION	4
IMPORTANT REMARKS	6
REFERENCES.....	7
DEFINITIONS.....	7
SYMBOLS AND ABBREVIATIONS	7
Original Equipment Manufacturers DICOM Conformance Statement	8
LIMITATIONS	8
1. Introduction.....	15
1.1. Purpose of this Document	15
1.2. Sources for this Document	15
1.3. Acronyms and Abbreviations	15
1.4. Typographical Conventions.....	15
2. Implementation Model.....	16
2.1. Application Data Flow Diagram.....	16
2.1.1. N-EVENT-REPORTs	16
2.1.2. N-GET Requests.....	16
2.1.3. N-CREATE and N-SET Study Component.....	16
2.1.4. C-FIND Modality Worklist	16
2.1.5. C-FIND Report.....	17
2.1.6. C-FIND Study	17
2.1.7. C-MOVE Study.....	17
2.2. Functional Definitions of AEs.....	17
2.3. Sequencing of Real World Activities	17
3. AE Specifications	18
3.1. <i>Broker</i> Specifications	18
3.1.1. Association Establishment Policies	18
3.1.1.1. General	18
3.1.1.2. Number of Associations	18
3.1.1.3. Asynchronous Nature	18
3.1.1.4. Implementation Identifying Information.....	18
3.1.1.5. Called/Calling Titles.....	19
3.1.2. Association Initiation by Real World Activity	19
3.1.2.1. Real World Activity - Verification	19
3.1.2.2. Real World Activity - Detached Patient Management.....	19
3.1.2.3. Real World Activity - Detached Visit Management	21
3.1.2.4. Real World Activity - Detached Study Management.....	23
3.1.2.5. Real World Activity - Detached Results Management	25
3.1.2.6. Real World Activity - Detached Interpretation Management	26
3.1.2.7. Real World Activity - Study Component Management	27
3.1.2.8. Real World Activity - Modality Worklist Management	28
3.1.2.9. Real World Activity - Find.....	31
3.1.2.10. Real World Activity - Move	32
3.1.3. Association Acceptance Policy.....	32
3.1.3.1. Real World Activity - Verification	32
3.1.3.2. Real World Activity - Detached Patient Management.....	33
3.1.3.3. Real World Activity - Detached Visit Management	35
3.1.3.4. Real World Activity - Detached Study Management.....	37
3.1.3.5. Real World Activity - Detached Results Management	39

Document # MM10001, Rev 1.0

3.1.3.6. Real World Activity - Detached Interpretation Management	41
3.1.3.7. Real World Activity - Study Component Management	42
3.1.3.8. Real World Activity - Modality Worklist Management	43
4. Communications Profiles	46
4.1. TCP/IP Stack.....	46
4.2. Physical Medium Supported.....	46
5. Extensions / Specializations.....	47
5.1. <i>Broker</i> Extended Specifications	47
5.1.1. Association Acceptance Policy.....	47
5.1.1.1. Real World Activity – Mitra Report Management	47
5.1.1.2. Real World Activity – Mitra Detached Patient Management	49
5.1.1.3. Real World Activity – Mitra Detached Visit Management.....	51
5.1.1.4. Real World Activity – Mitra Detached Study Management	53
5.1.1.5. Real World Activity – Mitra Detached Results Management.....	55
5.1.1.6. Real World Activity – Mitra Detached Interpretation Management.....	56
6. Support for Extended Character Sets	57

Table of Tables

Table 1: Verification SOP Class.....	18
Table 2: Management SOP Classes	18
Table 3 Query/Retrieve SOP Classes	18
Table 4: Transfer Syntaxes	19
Table 5: Presentation Contexts	19
Table 6: Transfer Syntaxes	20
Table 7: Presentation Contexts	20
Table 8: Detached Patient Management Object N-Event-Report Attributes	20
Table 9: Transfer Syntaxes	21
Table 10: Presentation Contexts	21
Table 11: Detached Visit Management Object N-Event-Report Attributes.....	22
Table 12: Transfer Syntaxes	23
Table 13: Presentation Contexts	23
Table 14: Detached Study Management Object N-Event-Report Attributes	23
Table 15: Transfer Syntaxes	25
Table 16: Presentation Contexts	25
Table 17: Detached Results Management Object N-Event-Report Attributes.....	25
Table 18: Transfer Syntaxes	26
Table 19: Presentation Contexts	26
Table 20: Detached Interpretation Management Object N-Event-Report Attributes.....	27
Table 21: Transfer Syntaxes	28
Table 22: Presentation Contexts	28
Table 23: Study Component Management Object N-Create Attributes.....	28
Table 24: Transfer Syntaxes	29
Table 25: Presentation Contexts	29
Table 26: Modality Worklist Information Model Attributes	29

Document # MM10001, Rev 1.0

Table 27: Transfer Syntaxes31

Table 28: Presentation Contexts31

Table 29: C-Find Key Attributes31

Table 30: Transfer Syntaxes32

Table 31: Presentation Contexts32

Table 32: Move Extended Negotiation.....32

Table 33: Transfer Syntaxes33

Table 34: Presentation Contexts33

Table 35: Verification status codes.....33

Table 36: Transfer Syntaxes34

Table 37: Presentation Contexts34

Table 38: Detached Patient Management Object N-Get Attributes34

Table 39: Detached Patient Management status codes.35

Table 40: Transfer Syntaxes35

Table 41: Presentation Contexts36

Table 42: Detached Visit Management Object N-Get Attributes36

Table 43: Detached Visit Management status codes.....37

Table 44: Transfer Syntaxes37

Table 45: Presentation Contexts37

Table 46: Detached Study Management Object N-Get Attributes38

Table 47: Detached Study Management status codes.39

Table 48: Transfer Syntaxes40

Table 49: Presentation Contexts40

Table 50: Detached Results Management Object N-Get Attributes40

Table 51: Detached Results Management status codes.....40

Table 52: Transfer Syntaxes41

Table 53: Presentation Contexts41

Document # MM10001, Rev 1.0

Table 54: Detached Interpretation Management Object N-Get Attributes41

Table 55: Detached Interpretation Management status codes.....42

Table 56: Transfer Syntaxes42

Table 57: Presentation Contexts42

Table 58: Study Component Management Object N-Create Attributes43

Table 59: Transfer Syntaxes43

Table 60: Presentation Contexts43

Table 61: Modality Worklist Information Model Attributes44

Table 62: Transfer Syntaxes47

Table 63: Presentation Contexts47

Table 64: Mitra Report Information Model Attributes48

Table 65: Transfer Syntaxes50

Table 66: Presentation Contexts50

Table 67: Mitra Detached Patient Management Object N-Get Attributes50

Table 68: Mitra Detached Patient Management status codes.51

Table 69: Transfer Syntaxes51

Table 70: Presentation Contexts51

Table 71: Mitra Detached Visit Management Object N-Get Attributes.....52

Table 72: Mitra Detached Visit Management status codes.....53

Table 73: Transfer Syntaxes53

Table 74: Presentation Contexts53

Table 75: Mitra Detached Study Management Object N-Get Attributes53

Table 76: Mitra Detached Study Management status codes.55

Table 77: Transfer Syntaxes55

Table 78: Presentation Contexts55

Table 79: Mitra Detached Results Management Object N-Get Attributes.....55

Table 80: Mitra Detached Results Management status codes56

Document # MM10001, Rev 1.0

Table 81: Transfer Syntaxes	56
Table 82: Presentation Contexts	56
Table 83: Mitra Detached Interpretation Management Object N-Get Attributes.....	57
Table 84: Mitra Detached Interpretation Management status codes.....	57

1. Introduction

1.1. Purpose of this Document

This document is a provisional DICOM Conformance Statement for the software product *Broker*. *Broker* is a service class provider for DIMSE-N services relating to access to Hospital Information Systems (HIS) and Radiology Information Systems (RIS). *Broker* is intended for use with a wide range of HIS/RIS technologies. Actual combinations for *Broker* plus HIS/RIS may be subject to restrictions not noted in this conformance statement. (For example, some combinations may support certain functions or attributes, and some may not.) Unless otherwise stated, all features conform to the DICOM V3.0 specification; all mandatory elements are supported.

Where the HIS/RIS functions as a service class provider of DIMSE-N services as described above, *Broker* can be configured to function as a service class user and/or service class provider of these services.

1.2. Sources for this Document

- ACR-NEMA Digital Imaging and Communications in Medicine (DICOM) v3.0, Final Text, 1996.
- Health Level Seven Version 2.1 (HL7 V2.1), 1991.

1.3. Acronyms and Abbreviations

The following acronyms and abbreviations are used in this document.

- ACR American College of Radiology
- ANSI American National Standards Institute
- DICOM Digital Imaging and Communications in Medicine
- DIMSE DICOM Message Service Element
- DIMSE-C DICOM Message Service Element-Composite
- DIMSE-N DICOM Message Service Element-Normalized
- NEMA National Electrical Manufacturers Association
- PDU Protocol Data Unit
- SCP Service Class Provider
- SCU Service Class User
- SOP Service Object Pair
- TCP/IP Transmission Control Protocol/Internet Protocol
- UID Unique Identifier
- ISIS Information System - Imaging System interface

1.4. Typographical Conventions

This section is designed to assist the reader in understanding the terms and typographical conventions used in this document.


Formatting convention	Type of information
Bold type	DICOM SOP Class, or DIMSE Service
<i>Italic type</i>	Application Entity

2. Implementation Model

Broker includes the following components:

- There is a data repository that supports the ISIS model.
- There is one or more interfaces to devices in the hospital information system, including in most cases a RIS. Typically these interfaces will conform to a standard protocol such as HL7 or DICOM.
- *Broker* is a DICOM service class provider for the following service classes: Patient Management, Study Management, Results Management, and Basic Worklist Management. *Broker* both accepts DICOM requests and issues DICOM events.

2.1. Application Data Flow Diagram



2.1.1. N-EVENT-REPORTS

In most configurations, *Broker* will receive event notifications from the HIS/RIS. That information will be stored in the *Broker*'s data repository. Depending on how *Broker* is configured, *Broker* will also send matching DICOM events to specified devices in the PACS. The distribution pattern for each class of event, and the time at which events are issued, is all under the control of data tables that are configurable via the *Broker* GUI.

Where the HIS/RIS is a DICOM SCP, *Broker* can be configured to receive N-EVENT Reports as SCU.

2.1.2. N-GET Requests

As SCP, *Broker* responds to N-GET requests for patient, visit, study, study component, results and interpretation information. The information that is returned within the N-GET response comes from either *Broker*'s data repository or from secondary queries made to the HIS/RIS. Which method of information retrieval is dependent upon the capabilities of the HIS/RIS, and upon the configuration of *Broker*.

As SCU, *Broker* issues N-GET requests for patient, visit, study, results and interpretation information. *Broker* supplies the requested SOP instance UID on all N-Get requests, and will update its data repository accordingly with the response.

2.1.3. N-CREATE and N-SET Study Component

Broker treats N-CREATE and N-SET study component exactly the same.

Broker will accept requests to N-CREATE a study component. These study components will be stored within the *Broker* data repository as specified by the ISIS model.

Broker can be configured to issue HL7 Trigger Events to interested parties in the HIS when study components are created. (See *Broker* HL7 Conformance Statement.)

2.1.4. C-FIND Modality Worklist

Broker will accept requests to C-FIND a modality worklist. These modality worklists will be stored within the *Broker* data repository as specified by the ISIS model. *Broker* can be configured to generate requests to C-FIND a modality worklist when the RIS supports C-FIND Modality Worklist as an SCP.

Document # MM10001, Rev 1.0

2.1.5. C-FIND Report

Broker will accept requests to C-FIND a report. This functionality is offered as an alternative to fetching reports using the normalized services specified by the DICOM specification. *Broker* can be configured to generate requests to C-FIND a report when the RIS supports C-FIND Report as an SCP.

2.1.6. C-FIND Study

Broker can be configured to generate requests to a PACS to C-FIND a study when the RIS issues a request to prefetch a study from the PACS.

2.1.7. C-MOVE Study

Broker can be configured to generate requests to a PACS to C-MOVE a study when the RIS issues a request to prefetch a study from the PACS.

2.2. Functional Definitions of AEs

Broker is implemented as a single application entity for a service class provider. The same application entity may also be configured as a service class user.

2.3. Sequencing of Real World Activities

Broker must have an installed/working connection with a HIS/RIS system in order for it to function properly in a hospital environment.

Broker will use its internal database to store the information received from the HIS/RIS for later retrieval via the N-GET mechanism. This database is built up with patient, visit, study, results, and interpretation information from the triggers that the HIS/RIS sends to *Broker*. *Broker* can also send out unsolicited information via the N-EVENT-REPORT mechanism, based upon the triggers received from the HIS/RIS. *Broker* can also be configured to generate query requests and study move requests to a PACS system via the C-FIND, C-MOVE mechanisms, based upon a prefetch trigger received from the RIS. The C-FIND is issued to retrieve the PACS-based study UID which is then used in the C-MOVE request. *Broker*, based on the initial RIS trigger received, determines the prefetch destination. This is repeated for each study in the prefetch list provided by the RIS.

3. AE Specifications

3.1. Broker Specifications

Broker provides Standard Conformance to the following DICOM V3.0 **Verification** SOP Class as an SCU and an SCP.

Table 1: Verification SOP Class

SOP Class	SOP Class UID
Verification	1.2.840.10008.1.1

Broker provides Standard Conformance to the following DICOM V3.0 **Management** SOP Classes as an SCP.

Table 2: Management SOP Classes

SOP Class	SOP Class UID
Detached Patient Management	1.2.840.10008.3.1.2.1.1
Detached Visit Management	1.2.840.10008.3.1.2.2.1
Detached Study Management	1.2.840.10008.3.1.2.3.1
Detached Results Management	1.2.840.10008.3.1.2.5.1
Detached Interpretation Management	1.2.840.10008.3.1.2.6.1
Study Component Management	1.2.840.10008.3.1.2.3.2
Modality Worklist Management	1.2.840.10008.5.1.4.3.1

Broker provides Standard Conformance to the following DICOM V3.0 **Query/Retrieve** SOP Class as an SCU.

Table 3 Query/Retrieve SOP Classes

SOP Class	SOP Class UID
Study Root Query/Retrieve IM Find	1.2.840.10008.5.1.4.1.2.2.1
Study Root Query/Retrieve IM Move	1.2.840.10008.5.1.4.1.2.2.2

3.1.1. Association Establishment Policies

3.1.1.1. General

3.1.1.2. Number of Associations

The maximum number of simultaneous associations accepted by *Broker* is configurable at run time, based on the system resources available. By default, the maximum number of associations is set at 32. There is no inherent limit to the number of associations other than limits imposed by the computer operating system.

3.1.1.3. Asynchronous Nature

Broker allows a single outstanding operation on any association. Therefore, *Broker* does not support asynchronous operations window negotiation, other than the default as specified by the DICOM specification.

3.1.1.4. Implementation Identifying Information

Broker will respond with the following implementation identifying parameters:

Document # MM10001, Rev 1.0

- Implementation Class UID **1.2.124.113532.3320**
- Implementation Version Name **BROKER971.0**

The implementation version name policies are the following: product name “**BROKER**” followed by the year of the product “**97**”, and finally the version of the product, “**1.0**”.

3.1.1.5. Called/Calling Titles

The default calling title that *Broker* will use is “**BROKER**”. This parameter can be configured before application startup. *Broker* can be configured to validate the Called Title of the requesting SCU during association negotiation. This policy provides *Broker* with a rudimentary level of security, by preventing unknown SCU processes from accessing hospital information.

3.1.2. Association Initiation by Real World Activity

Broker attempts to issue a new association when an unsolicited event message needs to be sent to an SCU, or query/ retrieve issued to an SCP.

3.1.2.1. Real World Activity - Verification

3.1.2.1.1. Associated Real World Activity - Verification

Broker uses the verification service class to test communication with a remote entity.

3.1.2.1.2. Presentation Context Table - Verification

Broker supports the transfer syntaxes listed in Table 4.

Table 4: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 5: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Verification	1.2.840.10008.1.1	all from Table 4	SCU	None

3.1.2.1.3. SOP Specific Conformance - Verification

Broker provides standard conformance to the DICOM **Verification** Service Class.

3.1.2.2. Real World Activity - Detached Patient Management

3.1.2.2.1. Associated Real World Activity - Detached Patient Management

3.1.2.2.1.1. Activity as SCP

Broker will send DIMSE-N EVENT-REPORTs to indicate that a change in the status of patient information has occurred. The following events are supported:

- Patient Created - to signal that a new patient has been added to the HIS/RIS database
- Patient Deleted - to signal that a patient has been removed from the HIS/RIS database
- Patient Updated - to signal that the patient information has changed

3.1.2.2.1.2. Activity as SCU

Document # MM10001, Rev 1.0

Broker will query for patient information using the DIMSE-N GET Patient to complete its information where the patient information is not locally stored.

Only the requested SOP Instance UID is specified in the N-GET request, it is expected that the SCP will provide additional elements of interest.

3.1.2.2.2. Presentation Context Table - Detached Patient Management

Broker supports the transfer syntaxes listed in Table 6.

Table 6: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 7: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Detached Patient Management	1.2.840.10008.3.1.2.1.1	all from Table 6	SCU/SCP	None

3.1.2.2.3. SOP Specific Conformance - Detached Patient Management

Broker provides standard conformance to the DICOM Detached Patient Management Service Class.

Broker supports the following elements for this SOP class:

Table 8: Detached Patient Management Object N-Event-Report Attributes

Event Type Name	Attribute Name	Tag
Patient Created	Instance Creation Date	(0008,0012)
	Instance Creation Time	(0008,0013)
	Instance Creator UID	(0008,0014)
Patient Deleted/Updated	Specific Character Set	(0008,0005)
	Patient Name	(0010,0010)
	Patient ID	(0010,0020)
	Issuer of Patient ID	(0010,0021)
	Other Patient Ids	(0010,1000)
	Other Patient Names	(0010,1001)
	Patient Telephone Numbers	(0010,2154)
	Patient Address	(0010,1040)
	Patient Birth Date	(0010,0030)
	Patient Sex	(0010,0040)
	Patient Weight	(0010,1030)
	Ethnic Group	(0010,2160)
	Patient Religious Preference	(0010,21F0)
	Patient Data Confidentiality Constraint Desc.	(0040,3001)
	Patient State	(0038,0500)
Pregnancy Status	(0010,21C0)	

Medical Alerts	(0010,2000)
Contrast Allergies	(0010,2110)
Special Needs	(0038,0050)
Referenced Study Sequence	(0008,1110)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Visit Sequence	(0008,1125)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Reference Patient Alias Sequence	(0038,0004)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)

3.1.2.3. Real World Activity - Detached Visit Management

3.1.2.3.1. Associated Real World Activity - Detached Visit Management

3.1.2.3.1.1. Activity as SCP

Broker will send DIMSE-N EVENT-REPORTs to indicate that a change in the status of visit information has occurred. The following events are supported:

- Visit Created - to signal that this is a new visit for this patient
- Visit Scheduled - to signal that this patient will be arriving at some time in the future
- Patient Admitted - to signal that the patient has been admitted into the hospital
- Patient Transferred - to signal that a patient has moved to a new location
- Patient Discharged - to signal that the patient has been discharged from the hospital

3.1.2.3.1.2. Activity as SCU

Broker will query for patient information using the DIMSE-N GET Visit to complete its information where the patient information is not locally stored.

Only the requested SOP Instance UID is specified in the N-GET request, it is expected that the SCP will provide additional elements of interest.

3.1.2.3.2. Presentation Context Table - Detached Visit Management

Broker supports the transfer syntaxes listed in Table 9.

Table 9: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 10: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Detached Visit Management	1.2.840.10008.3.1.2.2.1	all from Table 9	SCU/SCP	None

Document # MM10001, Rev 1.0

3.1.2.3.3. SOP Specific Conformance - Detached Visit Management

Broker provides standard conformance to the DICOM **Detached Visit Management** Service Class.

Broker supports the following elements for this SOP class:

Table 11: Detached Visit Management Object N-Event-Report Attributes

Event Type Name	Attribute Name	Tag
Visit Created	Instance Creation Date	(0008,0012)
	Instance Creation Time	(0008,0013)
	Instance Creator UID	(0008,0014)
Visit Scheduled / Deleted / Updated, Patient Admitted / Transferred / Discharged	Specific Character Set	(0008,0005)
	Institution Name	(0008,0080)
	Institution Address	(0008,0081)
	Admission ID	(0038,0010)
	Issuer of Admission ID	(0038,0011)
	Visit Status ID	(0038,0008)
	Current Patient Location	(0038,0300)
	Patient's Institution Residence	(0038,0400)
	Admitting Date	(0038,0020)
	Admitting Time	(0038,0021)
	Admitting Diagnosis Description	(0008,1080)
	Discharge Date	(0038,0030)
	Discharge Time	(0038,0032)
	Discharge Diagnosis Description	(0038,0040)
	Route of Admissions	(0038,0016)
	Referring Physician's Name	(0008,0090)
	Referring Physician's Address	(0008,0092)
	Scheduled Admission Date	(0038,001A)
	Scheduled Admission Time	(0038,001B)
	Scheduled Patient Institution Residence	(0038,001E)
	Scheduled Discharge Date	(0038,001C)
	Scheduled Discharge Time	(0038,001D)
	Referenced Patient Sequence	(0008,1120)
	>Referenced SOP Class UID	(0008,1150)
	>Referenced SOP Instance UID	(0008,1155)
	Referenced Study Sequence	(0008,1110)
>Referenced SOP Class UID	(0008,1150)	
>Referenced SOP Instance UID	(0008,1155)	

3.1.2.4. Real World Activity - Detached Study Management

3.1.2.4.1. Associated Real World Activity - Detached Study Management

3.1.2.4.1.1. Activity as SCP

Broker will send DIMSE-N EVENT-REPORTs to indicate that a change in the status of study information has occurred. The following events are supported:

- Study Created - to signal that a new study has been created
- Study Scheduled - to signal that the study has been scheduled to occur
- Patient Arrived - to signal that the patient has arrived for this study
- Study Started - to signal that the study has begun
- Study Completed - to signal that a study has been completed
- Study Verified - to signal that this study was successful
- Study Read - to signal that this study has been reviewed
- Study Deleted - to signal that this study has been canceled
- Study Updated - to signal that the study information has changed

3.1.2.4.1.2. Activity as SCU

Broker will query for patient information using the DIMSE-N GET Study to complete its information where the patient information is not locally stored.

Only the requested SOP Instance UID is specified in the N-GET request, it is expected that the SCP will provide additional elements of interest.

3.1.2.4.2. Presentation Context Table - Detached Study Management

Broker supports the transfer syntaxes listed in Table 12.

Table 12: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 13: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Detached Study Management	1.2.840.10008.3.1.2.3.1	all from Table 12	SCU/SCP	None

3.1.2.4.3. SOP Specific Conformance - Detached Study Management

Broker provides standard conformance to the DICOM **Detached Study Management** Service Class.

Broker supports the following elements for this SOP class:

Table 14: Detached Study Management Object N-Event-Report Attributes

Event Type Name	Attribute Name	Tag
Study Created	Instance Creation Date	(0008,0012)
	Instance Creation Time	(0008,0013)
	Instance Creator UID	(0008,0014)
Study Updated / Deleted / Scheduled / Started / Completed / Verified / Read, Patient Arrived	Specific Character Set	(0008,0005)

Study ID	(0020,0010)
Study ID Issuer	(0032,0012)
Accession Number	(0008,0050)
Study Instance UID	(0020,000D)
Study Status ID	(0032,000A)
Study Priority ID	(0032,000C)
Scheduled Study Start Date	(0032,1000)
Scheduled Study Start Time	(0032,1001)
Scheduled Study Stop Date	(0032,1010)
Scheduled Study Stop Time	(0032,1011)
Scheduled Study Location	(0032,1020)
Scheduled Study Location Application Entity Title	(0032,1021)
Requesting Service	(0032,1033)
Requesting Physician	(0032,1032)
Requested Procedure Description	(0032,1060)
Requested Procedure Code Sequence	(0032,1064)
>Code Value	(0008,0100)
>Coding Scheme Designator	(0008,0102)
>Code Meaning	(0008,0104)
Study Arrival Date	(0032,1040)
Study Arrival Time	(0032,1041)
Study Date	(0008,0020)
Study Time	(0008,0030)
Study Completed Date	(0032,1050)
Study Completed Time	(0032,1051)
Study Verified Date	(0032,0032)
Study Verified Time	(0032,0033)
Study Read Date	(0032,0034)
Study Read Time	(0032,0035)
Name of Physician(s) Reading Study	(0008,1060)
Reason For Study	(0032,1030)
Referenced Patient Sequence	(0008,1120)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Visit Sequence	(0008,1125)
>Referenced SOP Class UID	(0008,1150)

>Referenced SOP Instance UID	(0008,1155)
Referenced Results Sequence	(0008,1100)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Study Component Sequence	(0008,1111)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)

3.1.2.5. Real World Activity - Detached Results Management

3.1.2.5.1. Associated Real World Activity - Detached Results Management

3.1.2.5.1.1. Activity as SCP

Broker will send DIMSE-N EVENT-REPORTs to indicate that a change in the status of Results information has occurred. The following events are supported:

- Results Created - to signal that the results have been created
- Results Updated - to signal that the results information has changed
- Results Deleted - to signal that the results information has been removed from the HIS/RIS database

3.1.2.5.1.2. Activity as SCU

Broker will query for patient information using the DIMSE-N GET Results to complete its information where the patient information is not locally stored.

Only the requested SOP Instance UID is specified in the N-GET request, it is expected that the SCP will provide additional elements of interest.

3.1.2.5.2. Presentation Context Table - Detached Results Management

Broker supports the transfer syntaxes listed in Table 15.

Table 15: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 16: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Detached Results Management	1.2.840.10008.3.1.2.5.1	all from Table 15	SCU/SCP	None

3.1.2.5.3. SOP Specific Conformance - Detached Results Management

Broker provides standard conformance to the DICOM **Detached Interpretation Management** Service Class.

Broker supports the following elements for this SOP class:

Table 17: Detached Results Management Object N-Event-Report Attributes

Event Type Name	Attribute Name	Tag
Results Created	Instance Creation Date	(0008,0012)
	Instance Creation Time	(0008,0013)
	Instance Creator UID	(0008,0014)
Results Updated / Deleted	Specific Character Set	(0008,0005)
	Results ID	(4008,0040)
	Results ID Issuer	(4008,0042)
	Impressions	(4008,0300)
	Referenced Study Sequence	(0008,1110)
	> Referenced SOP Class UID	(0008,1150)
	> Referenced SOP Instance UID	(0008,1155)
	Referenced Interpretation Sequence	(4008,0005)
	> Referenced SOP Class UID	(0008,1150)
	> Referenced SOP Instance UID	(0008,1155)

3.1.2.6. Real World Activity - Detached Interpretation Management

3.1.2.6.1. Associated Real World Activity - Detached Interpretation Management

3.1.2.6.1.1. Activity as SCP

Broker will send DIMSE-N EVENT-REPORTs to indicate that a change in the status of interpretation information has occurred. The following events are supported:

- Interpretation Created - to signal the creation of the interpretation
- Interpretation Recorded - to signal that the results have been dictated
- Interpretation Transcribed - to signal that the preliminary results are ready
- Interpretation Approved - to signal that the final results are ready
- Interpretation Updated - to signal that the interpretation information has changed

3.1.2.6.1.2. Activity as SCU

Broker will query for patient information using the DIMSE-N GET Interpretation to complete its information where the patient information is not locally stored.

Only the requested SOP Instance UID is specified in the N-GET request, it is expected that the SCP will provide additional elements of interest.

3.1.2.6.2. Presentation Context Table - Detached Interpretation Management

Broker supports the transfer syntaxes listed in Table 18.

Table 18: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 19: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			

Detached Interpretation Management	1.2.840.10008.3.1.2.6.1	all from Table 18	SCU/SCP	None
------------------------------------	-------------------------	-------------------	---------	------

3.1.2.6.3. SOP Specific Conformance - Detached Interpretation Management

Broker provides standard conformance to the DICOM Detached Interpretation Management Service Class.

Broker supports the following elements for this SOP class:

Table 20: Detached Interpretation Management Object N-Event-Report Attributes

Event Type Name	Attribute Name	Tag
Interpretation Created	Instance Creation Date	(0008,0012)
	Instance Creation Time	(0008,0013)
	Instance Creator UID	(0008,0014)
Interpretation Updated / Deleted / Recorded / Transcribed / Approved	Specific Character Set	(0008,0005)
	Interpretation ID	(4008,0200)
	Interpretation ID Issuer	(4008,0202)
	Interpretation Type ID	(4008,0210)
	Interpretation Status ID	(4008,0212)
	Interpretation Recorded Date	(4008,0100)
	Interpretation Recorded Time	(4008,0101)
	Interpretation Recorder	(4008,0102)
	Interpretation Transcription Date	(4008,0108)
	Interpretation Transcription Time	(4008,0109)
	Interpretation Transcriber	(4008,010A)
	Interpretation Author	(4008,010C)
	Interpretation Text	(4008,010B)
	Referenced Results Sequence	(0008,1100)
	> Referenced SOP Class UID	(0008,1150)
	> Referenced SOP Instance UID	(0008,1155)
	Interpretation Approver Sequence	(4008,0111)
	>Interpretation Approval Date	(4008,0112)
	>Interpretation Approval Time	(4008,0113)
	>Physicians Approving Interpretation	(4008,0114)

3.1.2.7. Real World Activity - Study Component Management

3.1.2.7.1. Associated Real World Activity - Study Component Management

3.1.2.7.1.1. Activity as SCP

Document # MM10001, Rev 1.0

Broker will receive DIMSE-N CREATE, and DIMSE-N SET of Study Component to indicate that a Study Component object has been created.

Broker may be queried by the SCU for existing study components using the DIMSE-N GET Study Component service class.

3.1.2.7.2. Presentation Context Table - Study Component Management

Broker supports the transfer syntaxes listed in Table 56.

Table 21: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 22: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Study Component Management	1.2.840.10008.3.1.2.3.2	all from Table 21	SCU/SCP	None

3.1.2.7.3. SOP Specific Conformance - Study Component Management

Broker provides standard conformance to the DICOM Study Component Management Service Class.

Broker supports the following elements for this SOP class:

Table 23: Study Component Management Object N-Create Attributes

Attribute Name	Tag
Specific Character Set	(0008,0005)
Referenced Study Sequence	(0008,1110)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Retrieve Application Entity Title	(0008,0054)
Modality	(0008,0060)
Study Description	(0008,1030)
Acquisition In Study	(0020,1004)
Study Status ID	(0032,000a)

3.1.2.8. Real World Activity - Modality Worklist Management

3.1.2.8.1. Associated Real World Activity - Modality Worklist Management

3.1.2.8.1.1. Extensions

Broker makes use of either the Study Status Id attribute from the Imaging Service Request module (or for some legacy systems, the private element Scheduled Procedure Step Status from the Scheduled Procedure Step module). This has been added so that a PACS may obtain the status of the study via the Modality Worklist response.

Removal of these elements is a configurable option.

3.1.2.8.1.2. Activity as SCU

Document # MM10001, Rev 1.0

Broker can be configured to send DICOM C-Find requests in response to an external device querying Broker for Worklist or as the result of Broker being configured to poll for Worklist automatically.

3.1.2.8.2. Presentation Context Table - Modality Worklist Management

Broker supports the transfer syntaxes listed in Table 24.

Table 24: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 25: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Modality Worklist Info Model –FIND	1.2.840.10008.5.1.4.31	all from Table 24	SCU	None

3.1.2.8.3. SOP Specific Conformance - Modality Worklist Management

Broker provides standard conformance to the DICOM Basic **Worklist Management** Service Class.

Broker supports all required matching key types:

Matching Key Types	
SV	single valued match
WC	wild card match
SQ	sequence match
DR	date range match

Subject to availability from the HIS, Broker supports all required return keys.

Broker supports a 'NOT' operator (!) for all attributes with a single valued match (SV) type.

Broker supports the following elements for this SOP class:

Table 26: Modality Worklist Information Model Attributes

Module	Attribute Name	Tag	Match	Return
SOP Common	Specific Character Set	(0008,0005)		1C
Scheduled Procedure Step	Scheduled Procedure Step Sequence	(0040,0100)	SQ	1
	>Scheduled Station AE Title	(0040,0001)	SV	1
	>Scheduled Procedure Step Start Date	(0040,0002)	DR	1
	>Scheduled Procedure Step Start Time	(0040,0003)	DR	1
	>Scheduled Procedure Step End Date	(0040,0004)		1
	>Scheduled Procedure Step End Time	(0040,0005)		1
	>Modality	(0008,0060)	SV	1
	>Scheduled Performing Physician	(0040,0006)	WC	2
	>Scheduled Procedure Step Desc.	(0040,0007)		1C
	>Scheduled Station Name	(0040,0010)	SV	2
	>Scheduled Procedure Step Location	(0040,0011)		2
	>Scheduled Action Item Code Seq.	(0040,0008)		1C
	>>Code Value	(0008,0100)		1C

Requested Procedure	>>Coding Scheme Designator	(0008,0102)		1C
	>Pre-Medication	(0040,0012)		2C
	>Scheduled Procedure Step ID	(0040,0009)		1
	>Scheduled Procedure Step Status	(0040,0020)	SV	3
	>Requested Contrast Agent	(0032,1070)		2C
	>Comments on the Scheduled Procedure Step	(0040,0400)		3
	Requested Procedure ID	(0040,1001)		1
Imaging Service Request	Requested Procedure Description	(0032,1060)		1C
	Reason for Requested Procedure	(0040,1002)		3
	Requested Procedure Code Sequence	(0032,1064)		1C
	>Code Value	(0008,0100)		1C
	>Coding Scheme Designator	(0008,0102)		1C
	Study Instance UID	(0020,000D)	SV	1
	Referenced Study Sequence	(0008,1110)		2
	>Referenced SOP Class UID	(0008,1150)		1C
	>Referenced SOP Instance UID	(0008,1155)	SV	1C
	Requested Procedure Priority	(0040,1003)		2
	Patient Transport Arrangements	(0040,1004)		2
	Procedure Location	(0040,1005)		3
	Procedure Placer Order	(0040,1006)		3
	Procedure Filler Order	(0040,1007)		3
	Requested Procedure Comments	(0040,1400)		3
	Accession Number	(0008,0050)	SV	2
	Requesting Physician	(0032,1032)		2
	Requesting Service	(0032,1032)		3
	Referring Physician Name	(0008,0090)		2
	Reason for Imaging Service Request	(0040,2001)		3
	Study Status Id	(0032,000a)	SV	3
	Study Priority Id	(0032,000c)		3
	Visit Identification	Admission ID	(0038,0010)	SV
Visit Status	Current Patient Location	(0038,0300)	SV	2
Visit Relationship	Referenced Patient Sequence	(0008,1120)		2
	>Referenced SOP Class UID	(0008,1150)		2
	>Referenced SOP Instance UID	(0008,1155)	SV	2
Patient Identification	Patient Name	(0010,0010)	WC	1
	Patient ID	(0010,0020)	SV	1

Patient Demographic	Other Patient Id	(0010,1000)		3
	Other Patient Name	(0010,1001)		3
	Patient Address	(0010,1040)		3
	Patient Telephone Numbers	(0010,2154)		3
	Patient Birth Date	(0010,0030)		2
	Patient Sex	(0010,0040)		2
	Patient Weight	(0010,1030)		2
	Confidentiality Constraint	(0040,3001)		2
Patient Medical	Patient State	(0038,0500)		2
	Pregnancy Status	(0010,21C0)		2
	Medical Alerts	(0010,2000)		2
	Contrast Allergies	(0010,2110)		2
	Special Needs	(0038,0050)		2

3.1.2.9. Real World Activity - Find

3.1.2.9.1. Associated Real World Activity - Find

3.1.2.9.1.1. Activity as SCU

Broker will query the SCP (typically a PACS) to get the study UID associated with the accession number provided by the RIS. This is done in preparation for a C-Move request as described in section Real World Activity - Move. Broker negotiates the Study Root query/ retrieve model only.

3.1.2.9.2. Presentation Context Table - Find

Broker supports the transfer syntaxes listed in Table 27. Broker will initiate any of the Presentation Contexts list in Table 28 for Query.

Table 27: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 28: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Study Root Query/Retrieve IM Find	1.2.840.10008.5.1.4.1.2.2.1	all from Table 27	SCU	None

3.1.2.9.3. SOP Specific Conformance - Find

SOP classes of the Query/Retrieve Service Class are implemented via the DIMSE C-FIND and C-MOVE services as defined in Part 7 of the DICOM standard.

Broker will include the following key attributes in its C-Find request issued to the PACS:

Table 29: C-Find Key Attributes

Attribute Name	Tag
----------------	-----

Document # MM10001, Rev 1.0

Specific Character Set	(0008,0005)
Accession Number	(0008,0050)
Study Instance UID	(0020,000D)

Broker will include the Study Instance UID attribute as an empty string indicating that Broker is requesting the study UID be returned.

3.1.2.10. Real World Activity - Move

3.1.2.10.1. Associated Real World Activity - Move

3.1.2.10.1.1. Activity as SCU

Broker will request the SCP (typically a PACS) move the study associated with the study UID provided by the C-Find (section Real World Activity - Find) to the destination AE which is determined based on the 'RIS-based prefetch event' received from the RIS. Broker negotiates the Study Root query/ retrieve model only.

3.1.2.10.2. Presentation Context Table - Move

Broker supports the transfer syntaxes listed in Table 30. Broker will initiate any of the Presentation Contexts list in Table 31 for Retrieval.

Table 30: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 31: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Study Root Query/Retrieve IM Move	1.2.840.10008.5.1.4.1.2.2.2	all from Table 30	SCU	see Note 1

Note 1: C-Move Extended Negotiation will be supported. Broker will include the following information in its association negotiation request:

Table 32: Move Extended Negotiation

Field Name	Value	Description of Field
Relational-retrieval	1	relational retrieval supported

3.1.2.10.3. SOP Specific Conformance - Move

SOP classes of the Query/Retrieve Service Class are implemented via the DIMSE C-FIND and C-MOVE services as defined in Part 7 of the DICOM standard.

3.1.3. Association Acceptance Policy

3.1.3.1. Real World Activity - Verification

3.1.3.1.1. Associated Real World Activity - Verification

Document # MM10001, Rev 1.0

Broker will respond to **Verification** requests to provide an SCU with the ability to determine if *Broker* is receiving DICOM requests.

3.1.3.1.2. Presentation Context Table - Verification

Broker supports the transfer syntaxes listed in Table 33. *Broker* will accept any of the Presentation Contexts listed in Table 34 for **Verification**.

Table 33: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 34: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Verification	1.2.840.10008.1.1	all from Table 33	SCP	None

3.1.3.1.3. SOP Specific Conformance - Verification

Broker provides standard conformance to the DICOM **Verification** Service Class. *Broker* returns one of the following status codes.

Table 35: Verification status codes.

Service Status	Further Meaning	Protocol Codes	Related Fields	Description
Success	Success	0000		Operation performed properly.

3.1.3.1.4. Presentation Context Acceptance Criterion - Verification

Broker will always accept a Presentation Context for the Verification SOP Class with the default DICOM transfer syntax listed in Table 33.

3.1.3.1.5. Transfer Syntax Selection Policies - Verification

Since no DICOM data object is associated with a **Verification** command, only the default DICOM transfer syntax is required/supported.

3.1.3.2. Real World Activity - Detached Patient Management

3.1.3.2.1. Associated Real World Activity - Detached Patient Management

3.1.3.2.1.1. Activity as SCP

Broker will respond to DIMSE N-GET requests to retrieve patient demographic information. *Broker* does not support any other DIMSE-N command as an SCP with this presentation context.

3.1.3.2.1.2. Activity as SCU

Broker will receive DIMSE-N EVENT-REPORTs to indicate that a change in the status of patient information has occurred. The following events are supported:

- Patient Created - to signal that a new patient has been added to the HIS/RIS database
- Patient Deleted - to signal that a patient has been removed from the HIS/RIS database

Document # MM10001, Rev 1.0

- Patient Updated - to signal that the patient information has changed

3.1.3.2.2. Presentation Context Table - Detached Patient Management

Broker supports the transfer syntaxes listed in Table 36.

Table 36: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 37: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Detached Patient Management	1.2.840.10008.3.1.2.1.1	all from Table 36	SCP	None

3.1.3.2.3. SOP Specific Conformance - Detached Patient Management

Broker provides standard conformance to the DICOM **Detached Patient Management** Service Class.

Broker supports the following elements for this SOP class:

Table 38: Detached Patient Management Object N-Get Attributes

Attribute Name	Tag
Specific Character Set	(0008,0005)
Patient Name	(0010,0010)
Patient ID	(0010,0020)
Issuer of Patient ID	(0010,0021)
Other Patient Ids	(0010,1000)
Other Patient Names	(0010,1001)
Patient Telephone Numbers	(0010,2154)
Patient Address	(0010,1040)
Patient Birth Date	(0010,0030)
Patient Sex	(0010,0040)
Patient Weight	(0010,1030)
Ethnic Group	(0010,2160)
Patient Religious Preference	(0010,21F0)
Patient Data Confidentiality Constraint Desc.	(0040,3001)
Patient State	(0038,0500)
Pregnancy Status	(0010,21C0)
Medical Alerts	(0010,2000)
Contrast Allergies	(0010,2110)
Special Needs	(0038,0050)
Referenced Study Sequence	(0008,1110)
>Referenced SOP Class UID	(0008,1150)

>Referenced SOP Instance UID	(0008,1155)
Referenced Visit Sequence	(0008,1125)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Reference Patient Alias Sequence	(0038,0004)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)

Broker returns one of the following status codes.

Table 39: Detached Patient Management status codes.

Service Status	Further Meaning	Protocol Codes	Related Fields	Description
Error	Failed	0110		The operation was not successful.
Success	Success	0000		Operation performed properly.

3.1.3.2.4. Presentation Context Acceptance Criterion - Detached Patient Management

Broker will always accept a Presentation Context for the Detached Patient Management SOP Class.

3.1.3.2.5. Transfer Syntax Selection Policies - Detached Patient Management

Broker supports only the Little Endian Implicit Transfer Syntax.

3.1.3.3. Real World Activity - Detached Visit Management

3.1.3.3.1. Associated Real World Activity - Detached Visit Management

3.1.3.3.1.1. Activity as SCP

Broker will respond to DIMSE N-GET requests to retrieve visit information. Broker does not support any other DIMSE-N command as an SCP with this presentation context.

3.1.3.3.1.2. Activity as SCU

Broker will receive DIMSE-N EVENT-REPORTs to indicate that a change in the status of visit information has occurred. The following events are supported:

- Visit Created - to signal that this is a new visit for this patient
- Visit Scheduled - to signal that this patient will be arriving at some time in the future
- Patient Admitted - to signal that the patient has been admitted into the hospital
- Patient Transferred - to signal that a patient has moved to a new location
- Patient Discharged - to signal that the patient has been discharged from the hospital

3.1.3.3.2. Presentation Context Table - Detached Visit Management

Broker supports the transfer syntaxes listed in Table 40.

Table 40: Transfer Syntaxes

Transfer Syntax	UID
-----------------	-----

DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2
---	-------------------

Table 41: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Detached Visit Management	1.2.840.10008.3.1.2.2.1	all from Table 40	SCP	None

3.1.3.3.3. SOP Specific Conformance - Detached Visit Management

Broker provides standard conformance to the DICOM **Detached Visit Management** Service Class.

Broker supports the following elements for this SOP class:

Table 42: Detached Visit Management Object N-Get Attributes

Attribute Name	Tag
Specific Character Set	(0008,0005)
Institution Name	(0008,0080)
Institution Address	(0008,0081)
Admission ID	(0038,0010)
Issuer of Admission ID	(0038,0011)
Visit Status ID	(0038,0008)
Current Patient Location	(0038,0300)
Patient's Institution Residence	(0038,0400)
Admitting Date	(0038,0020)
Admitting Time	(0038,0021)
Admitting Diagnosis Description	(0008,1080)
Discharge Date	(0038,0030)
Discharge Time	(0038,0032)
Discharge Diagnosis Description	(0038,0040)
Route of Admissions	(0038,0016)
Referring Physician's Name	(0008,0090)
Referring Physician's Address	(0008,0092)
Scheduled Admission Date	(0038,001A)
Scheduled Admission Time	(0038,001B)
Scheduled Patient Institution Residence	(0038,001E)
Scheduled Discharge Date	(0038,001C)
Scheduled Discharge Time	(0038,001D)
Referenced Patient Sequence	(0008,1120)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Study Sequence	(0008,1110)

>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)

Broker returns one of the following status codes.

Table 43: Detached Visit Management status codes.

Service Status	Further Meaning	Protocol Codes	Related Fields	Description
Error	Failed	0110		The operation was not successful.
Success	Success	0000		Operation performed properly.

3.1.3.3.4. Presentation Context Acceptance Criterion - Detached Visit Management

Broker will always accept a Presentation Context for the Detached Visit Management SOP Class.

3.1.3.3.5. Transfer Syntax Selection Policies - Detached Visit Management

Broker supports only the Little Endian Implicit Transfer Syntax.

3.1.3.4. Real World Activity - Detached Study Management

3.1.3.4.1. Associated Real World Activity - Detached Study Management

3.1.3.4.1.1. Activity as SCP

Broker will respond to DIMSE N-GET requests to retrieve study information. Broker does not support any other DIMSE-N command as an SCP with this presentation context.

3.1.3.4.1.2. Activity as SCU

Broker will receive DIMSE-N EVENT-REPORTs to indicate that a change in the status of study information has occurred. The following events are supported:

- Study Created - to signal that a new study has been created
- Study Scheduled - to signal that the study has been scheduled to occur
- Patient Arrived - to signal that the patient has arrived for this study
- Study Started - to signal that the study has begun
- Study Completed - to signal that a study has been completed
- Study Verified - to signal that this study was successful
- Study Read - to signal that this study has been reviewed
- Study Deleted - to signal that this study has been canceled
- Study Updated - to signal that the study information has changed

3.1.3.4.2. Presentation Context Table - Detached Study Management

Broker supports the transfer syntaxes listed in Table 44.

Table 44: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 45: Presentation Contexts

Abstract Syntax	Transfer Syntax	Role	Extended Negotiation

SOP Class	SOP Class UID			
Detached Study Management	1.2.840.10008.3.1.2.3.1	All from Table 44	SCP	None

3.1.3.4.3. SOP Specific Conformance - Detached Study Management

Broker provides standard conformance to the DICOM **Detached Study Management** Service Class.

Broker supports the following elements for this SOP class:

Table 46: Detached Study Management Object N-Get Attributes

Attribute Name	Tag
Specific Character Set	(0008,0005)
Study ID	(0020,0010)
Study ID Issuer	(0032,0012)
Accession Number	(0008,0050)
Study Instance UID	(0020,000D)
Study Status ID	(0032,000A)
Study Priority ID	(0032,000C)
Scheduled Study Start Date	(0032,1000)
Scheduled Study Start Time	(0032,1001)
Scheduled Study Stop Date	(0032,1010)
Scheduled Study Stop Time	(0032,1011)
Scheduled Study Location	(0032,1020)
Scheduled Study Location Application Entity Title	(0032,1021)
Requesting Service	(0032,1033)
Requesting Physician	(0032,1032)
Requested Procedure Description	(0032,1060)
Requested Procedure Code Sequence	(0032,1064)
>Code Value	(0008,0100)
>Coding Scheme Designator	(0008,0102)
>Code Meaning	(0008,0104)
Study Arrival Date	(0032,1040)
Study Arrival Time	(0032,1041)
Study Date	(0008,0020)
Study Time	(0008,0030)
Study Completed Date	(0032,1050)
Study Completed Time	(0032,1051)
Study Verified Date	(0032,0032)
Study Verified Time	(0032,0033)
Study Read Date	(0032,0034)
Study Read Time	(0032,0035)

Name of Physician(s) Reading Study	(0008,1060)
Reason For Study	(0032,1030)
Referenced Patient Sequence	(0008,1120)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Visit Sequence	(0008,1125)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Results Sequence	(0008,1100)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Study Component Sequence	(0008,1111)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)

Broker returns one of the following status codes.

Table 47: Detached Study Management status codes.

Service Status	Further Meaning	Protocol Codes	Related Fields	Description
Error	Failed	0110		The operation was not successful.
Success	Success	0000		Operation performed properly.

3.1.3.4.4. Presentation Context Acceptance Criterion - Detached Study Management

Broker will always accept a Presentation Context for the Detached Study Management SOP Class.

3.1.3.4.5. Transfer Syntax Selection Policies - Detached Study Management

Broker supports only the Little Endian Implicit Transfer Syntax.

3.1.3.5. Real World Activity - Detached Results Management

3.1.3.5.1. Associated Real World Activity - Detached Results Management

3.1.3.5.1.1. Activity as SCP

Broker will respond to DIMSE N-GET requests to retrieve results information. Broker does not support any other DIMSE-N command as an SCP with this presentation context.

3.1.3.5.1.2. Activity as SCU

Broker will receive DIMSE-N EVENT-REPORTs to indicate that a change in the status of results information has occurred. The following events are supported:

- Results Created - to signal that the results have been created
- Results Updated - to signal that the results information has changed
- Results Deleted - to signal that the results information has been removed from the HIS/RIS database

Document # MM10001, Rev 1.0

3.1.3.5.2. Presentation Context Table - Detached Results Management

Broker supports the transfer syntaxes listed in Table 48.

Table 48: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 49: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Detached Results Management	1.2.840.10008.3.1.2.5.1	all from Table 48	SCP	None

3.1.3.5.3. SOP Specific Conformance - Detached Results Management

Broker provides standard conformance to the DICOM **Detached Results Management** Service Class.

Broker supports the following elements for this SOP class:

Table 50: Detached Results Management Object N-Get Attributes

Attribute Name	Tag
Specific Character Set	(0008,0005)
Results ID	(4008,0040)
Results ID Issuer	(4008,0042)
Impressions	(4008,0300)
Referenced Study Sequence	(0008,1110)
> Referenced SOP Class UID	(0008,1150)
> Referenced SOP Instance UID	(0008,1155)
Referenced Interpretation Sequence	(4008,0005)
> Referenced SOP Class UID	(0008,1150)
> Referenced SOP Instance UID	(0008,1155)

Broker returns one of the following status codes.

Table 51: Detached Results Management status codes.

Service Status	Further Meaning	Protocol Codes	Related Fields	Description
Error	Failed	0110		The operation was not successful.
Success	Success	0000		Operation performed properly.

3.1.3.5.4. Presentation Context Acceptance Criterion - Detached Results Management

Broker will always accept a Presentation Context for the Detached Results Management SOP Class.

3.1.3.5.5. Transfer Syntax Selection Policies - Detached Results Management

Broker supports only the Little Endian Implicit Transfer Syntax.

3.1.3.6. Real World Activity - Detached Interpretation Management

3.1.3.6.1. Associated Real World Activity - Detached Interpretation Management

3.1.3.6.1.1. Activity as SCP

Broker will respond to DIMSE N-GET requests to retrieve interpretation information. Broker does not support any other DIMSE-N command as an SCP with this presentation context.

3.1.3.6.1.2. Activity as SCU

Broker will receive DIMSE-N EVENT-REPORTs to indicate that a change in the status of interpretation information has occurred. The following events are supported:

- Interpretation Created - to signal the creation of the interpretation
- Interpretation Recorded - to signal that the results have been dictated
- Interpretation Transcribed - to signal that the preliminary results are ready
- Interpretation Approved - to signal that the final results are ready
- Interpretation Updated - to signal that the interpretation information has changed

3.1.3.6.2. Presentation Context Table - Detached Interpretation Management

Broker supports the transfer syntaxes listed in Table 52.

Table 52: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 53: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Detached Patient Management	1.2.840.10008.3.1.2.6.1	all from Table 52	SCP	None

3.1.3.6.3. SOP Specific Conformance - Detached Interpretation Management

Broker provides standard conformance to the DICOM Detached Interpretation Management Service Class.

Broker supports the following elements for this SOP class:

Table 54: Detached Interpretation Management Object N-Get Attributes

Attribute Name	Tag
Specific Character Set	(0008,0005)
Interpretation ID	(4008,0200)
Interpretation ID Issuer	(4008,0202)
Interpretation Type ID	(4008,0210)
Interpretation Status ID	(4008,0212)
Interpretation Recorded Date	(4008,0100)
Interpretation Recorded Time	(4008,0101)
Interpretation Recorder	(4008,0102)
Interpretation Transcription Date	(4008,0108)

Interpretation Transcription Time	(4008,0109)
Interpretation Transcriber	(4008,010A)
Interpretation Author	(4008,010C)
Interpretation Text	(4008,010B)
Referenced Results Sequence	(0008,1100)
> Referenced SOP Class UID	(0008,1150)
> Referenced SOP Instance UID	(0008,1155)
Interpretation Approver Sequence	(4008,0111)
>Interpretation Approval Date	(4008,0112)
>Interpretation Approval Time	(4008,0113)
>Physicians Approving Interpretation	(4008,0114)

Broker returns one of the following status codes.

Table 55: Detached Interpretation Management status codes.

Service Status	Further Meaning	Protocol Codes	Related Fields	Description
Error	Failed	0110		The operation was not successful.
Success	Success	0000		Operation performed properly.

3.1.3.6.4. Presentation Context Acceptance Criterion - Detached Interpretation Management

Broker will always accept a Presentation Context for the Detached Interpretation Management SOP Class.

3.1.3.6.5. Transfer Syntax Selection Policies - Detached Interpretation Management

Broker supports only the Little Endian Implicit Transfer Syntax.

3.1.3.7. Real World Activity - Study Component Management

3.1.3.7.1. Associated Real World Activity - Study Component Management

3.1.3.7.1.1. Activity as SCP

Broker will receive DIMSE-N CREATE, and DIMSE-N SET of Study Component to indicate that a Study Component object has been created.

Broker may be queried by the SCU for existing Study Components using the DIMSE-N GET Study Component service class.

3.1.3.7.2. Presentation Context Table - Study Component Management

Broker supports the transfer syntaxes listed in Table 56.

Table 56: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 57: Presentation Contexts

Abstract Syntax	Transfer Syntax	Role	Extended
-----------------	-----------------	------	----------

SOP Class		SOP Class UID			Negotiation
Study Component Management		1.2.840.10008.3.1.2.3.2	all from Table 56	SCP	None

3.1.3.7.3. SOP Specific Conformance - Study Component Management

Broker provides standard conformance to the DICOM Study Component Management Service Class.

Broker supports the following elements for this SOP class:

Table 58: Study Component Management Object N-Create Attributes

Attribute Name	Tag
Specific Character Set	(0008,0005)
Referenced Study Sequence	(0008,1110)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Retrieve Application Entity Title	(0008,0054)
Modality	(0008,0060)
Study Description	(0008,1030)
Acquisition In Study	(0020,1004)
Study Status ID	(0032,000a)

3.1.3.8. Real World Activity - Modality Worklist Management

3.1.3.8.1. Associated Real World Activity - Modality Worklist Management

3.1.3.8.1.1. Extensions

Broker makes use of either the Study Status Id attribute from the Imaging Service Request module (or for some legacy systems, the private element Scheduled Procedure Step Status from the Scheduled Procedure Step module). This has been added so that the SCU may obtain the status of the study using the Modality Worklist response.

Removal of these elements is a configurable option.

3.1.3.8.1.2. Activity as SCP

Broker can be configured to respond to DICOM C-Find requests in response to an external device querying Broker for Worklist.

3.1.3.8.2. Presentation Context Table - Modality Worklist Management

Broker supports the transfer syntaxes listed in Table 59.

Table 59: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 60: Presentation Contexts

Abstract Syntax	Transfer Syntax	Role	Extended Negotiation

SOP Class	SOP Class UID			
Modality Worklist Info Model –FIND	1.2.840.10008.5.1.4.31	all from Table 59	SCP	None

3.1.3.8.3. SOP Specific Conformance - Modality Worklist Management

Broker provides standard conformance to the DICOM Basic **Worklist Management** Service Class.

Broker supports all required matching key types:

Matching Key Types	
SV	single valued match
WC	wild card match
SQ	sequence match
DR	date range match

Subject to availability from the HIS, Broker supports all required return keys.

Broker supports a 'NOT' operator (!) for all attributes with a single valued match (SV) type.

Broker supports the following elements for this SOP class:

Table 61: Modality Worklist Information Model Attributes

Module	Attribute Name	Tag	Match	Return
SOP Common	Specific Character Set	(0008,0005)		1C
Scheduled Procedure Step	Scheduled Procedure Step Sequence	(0040,0100)	SQ	1
	>Scheduled Station AE Title	(0040,0001)	SV	1
	>Scheduled Procedure Step Start Date	(0040,0002)	DR	1
	>Scheduled Procedure Step Start Time	(0040,0003)	DR	1
	>Scheduled Procedure Stop End Date	(0040,0004)		1
	>Scheduled Procedure Stop End Time	(0040,0005)		1
	>Modality	(0008,0060)	SV	1
	>Scheduled Performing Physician	(0040,0006)	WC	2
	>Scheduled Procedure Step Desc.	(0040,0007)		1C
	>Scheduled Station Name	(0040,0010)	SV	2
	>Scheduled Procedure Step Location	(0040,0011)		2
	>Scheduled Action Item Code Seq.	(0040,0008)		1C
	>>Code Value	(0008,0100)		1C
	>>Coding Scheme Designator	(0008,0102)		1C
	>Pre-Medication	(0040,0012)		2C
	>Scheduled Procedure Step ID	(0040,0009)		1
	>Scheduled Procedure Step Status	(0040,0020)	SV	3
	>Requested Contrast Agent	(0032,1070)		2C
>Comments on the Scheduled Procedure Step	(0040,0400)		3	
Requested Procedure	Requested Procedure ID	(0040,1001)		1
	Requested Procedure Description	(0032,1060)		1C
	Reason for Requested Procedure	(0040,1002)		3

Imaging Service Request	Requested Procedure Code Sequence	(0032,1064)		1C
	>Code Value	(0008,0100)		1C
	>Coding Scheme Designator	(0008,0102)		1C
	Study Instance UID	(0020,000D)	SV	1
	Referenced Study Sequence	(0008,1110)		2
	>Referenced SOP Class UID	(0008,1150)		1C
	>Referenced SOP Instance UID	(0008,1155)	SV	1C
	Requested Procedure Priority	(0040,1003)		2
	Patient Transport Arrangements	(0040,1004)		2
	Procedure Location	(0040,1005)		3
	Procedure Placer Order	(0040,1006)		3
	Procedure Filler Order	(0040,1007)		3
	Requested Procedure Comments	(0040,1400)		3
	Accession Number	(0008,0050)	SV	2
	Requesting Physician	(0032,1032)		2
	Requesting Service	(0032,1033)		3
	Referring Physician Name	(0008,0090)		2
	Reason for Imaging Service Request	(0040,2001)		3
Study Status ID	(0032,000a)	SV	3	
Study Priority ID	(0032,000c)		3	
Visit Identification	Admission ID	(0038,0010)	SV	2
Visit Status	Current Patient Location	(0038,0300)		2
Visit Relationship	Referenced Patient Sequence	(0008,1120)		2
Patient Identification	>Referenced SOP Class UID	(0008,1150)		2
	>Referenced SOP Instance UID	(0008,1155)		2
	Patient Name	(0010,0010)	WC	1
	Patient ID	(0010,0020)	SV	1
	Other Patient ID	(0010,1000)		3
	Other Patient Name	(0010,1001)		3
Patient Demographic	Patient Address	(0010,1040)		3
	Patient Telephone Numbers	(0010,2154)		3
	Patient Birth Date	(0010,0030)		2
	Patient Sex	(0010,0040)		2
	Patient Weight	(0010,1030)		2
	Confidentiality Constraint	(0040,3001)		2
Patient Medical	Patient State	(0038,0500)		2
	Pregnancy Status	(0010,21C0)		2

	Medical Alerts	(0010,2000)		2
	Contrast Allergies	(0010,2110)		2
	Special Needs	(0038,0050)		2

4. Communications Profiles

Broker provides DICOM V3.0 TCP/IP Network Communication Support as defined in Part 8 of the DICOM Standard.

4.1. TCP/IP Stack

Broker inherits its TCP/IP stack from the computer system upon which it executes.

4.2. Physical Medium Supported

Broker is indifferent to the physical medium over which TCP/IP executes; it inherits the medium from the computer system upon which it executes.

5. Extensions / Specializations

5.1. Broker Extended Specifications

Broker provides Conformance to the following private SOP Classes as an SCP.

Table 54: Extended SOP Classes

SOP Class	SOP Class UID
Mitra Report Management	1.2.840.113532.3500.8
Mitra Detached Patient Management	1.2.840.113532.3500.10
Mitra Detached Visit Management	1.2.840.113532.3500.11
Mitra Detached Study Management	1.2.840.113532.3500.13
Mitra Detached Results Management	1.2.840.113532.3500.15
Mitra Detached Interpretation Management	1.2.840.113532.3500.16
Mitra Detached Study Component Management	1.2.840.113532.3500.14

5.1.1. Association Acceptance Policy

5.1.1.1. Real World Activity – Mitra Report Management

5.1.1.1.1. Associated Real World Activity – Mitra Report Management

5.1.1.1.1.1. Activity as SCP

Broker can be configured to respond to DICOM C-FIND requests in response to an external device querying Broker for a Report.

5.1.1.1.1.2. Activity as SCU

Broker can be configured to request Reports via DICOM C-FIND in response to an external device querying Broker for a Report.

5.1.1.1.2. Presentation Context Table – Mitra Report Management

Broker supports the transfer syntaxes listed in Table 62.

Table 62: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 63: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Mitra Report Info Model –FIND	1.2.840.113532.3500.8	all from Table 62	SCP	None

5.1.1.1.3. SOP Specific Conformance – Mitra Report Management

Document # MM10001, Rev 1.0

Broker provides conformance to a Mitra Private SOP class used to fetch report information. Unlike the modality worklist query, all supported attributes are included in each response independent of which attributes are included in the query.

Subject to availability from the HIS, Broker supports the following elements for this SOP class:

Table 64: Mitra Report Information Model Attributes

Attribute Name	Tag	Match	Remark
Specific Character Set	(0008,0005)		
Patient Name	(0010,0010)		
Patient ID	(0010,0020)	SV	Required for some HISs
Issuer of Patient ID	(0010,0021)		
Other Patient Ids	(0010,1000)		
Other Patient Names	(0010,1001)		
Patient Telephone Numbers	(0010,2154)		
Patient Address	(0010,1040)		
Patient Birth Date	(0010,0030)		
Patient Sex	(0010,0040)		
Patient Weight	(0010,1030)		
Ethnic Group	(0010,2160)		
Patient Religious Preference	(0010,21F0)		
Patient Data Confidentiality Constraint Desc.	(0040,3001)		
Patient State	(0038,0500)		
Pregnancy Status	(0010,21C0)		
Medical Alerts	(0010,2000)		
Contrast Allergies	(0010,2110)		
Special Needs	(0038,0050)		
Study ID	(0020,0010)		
Study ID Issuer	(0032,0012)		
Accession Number	(0008,0050)	SV	One or more of the accession number , study instance uid, or patient id are required
Study Instance UID	(0020,000D)	SV	
Study Status ID	(0032,000A)		
Study Priority ID	(0032,000C)		
Requesting Service	(0032,1033)		
Requesting Physician	(0032,1032)		
Requested Procedure Description	(0032,1060)		
Requested Procedure Code Sequence	(0032,1064)		
>Code Value	(0008,0100)		
>Coding Scheme Designator	(0008,0102)		

>Code Meaning	(0008,0104)		
Study Arrival Date	(0032,1040)		
Study Arrival Time	(0032,1041)		
Study Date	(0008,0020)		
Study Time	(0008,0030)		
Study Completed Date	(0032,1050)		
Study Completed Time	(0032,1051)		
Study Verified Date	(0032,0032)		
Study Verified Time	(0032,0033)		
Study Read Date	(0032,0034)		
Study Read Time	(0032,0035)		
Name of Physician(s) Reading Study	(0008,1060)		
Reason For Study	(0032,1030)		
Results ID	(4008,0040)		
Results ID Issuer	(4008,0042)		
Impressions	(4008,0300)		
Interpretation ID	(4008,0200)		
Interpretation ID Issuer	(4008,0202)		
Interpretation Type ID	(4008,0210)		
Interpretation Status ID	(4008,0212)		
Interpretation Recorded Date	(4008,0100)		
Interpretation Recorded Time	(4008,0101)		
Interpretation Recorder	(4008,0102)		
Interpretation Transcription Date	(4008,0108)		
Interpretation Transcription Time	(4008,0109)		
Interpretation Transcriber	(4008,010A)		
Interpretation Author	(4008,010C)		
Interpretation Text	(4008,010B)		
Interpretation Approver Sequence	(4008,0111)		
>Interpretation Approval Date	(4008,0112)		
>Interpretation Approval Time	(4008,0113)		
>Physicians Approving Interpretation	(4008,0114)		

5.1.1.2. Real World Activity – Mitra Detached Patient Management

5.1.1.2.1. Associated Real World Activity – Mitra Detached Patient Management

5.1.1.2.1.1. Activity as SCP

Document # MM10001, Rev 1.0

Broker will respond to DIMSE N-GET requests to retrieve patient demographic information. *Broker* does not support any other DIMSE-N command as an SCP with this presentation context.

5.1.1.2.2. Presentation Context Table – Mitra Detached Patient Management

Broker supports the transfer syntaxes listed in Table 65.

Table 65: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 66: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Mitra Detached Patient Management	1.2.840.113532.3500.10	all from Table 65	SCP	None

5.1.1.2.3. SOP Specific Conformance – Mitra Detached Patient Management

Broker provides conformance to a Mitra Private SOP class used to fetch patient demographic information. The purpose for the private SOP class is to minimize the overhead that goes along with the normalized DICOM model (ie: you send something at the results level, and to obtain patient level info, the SCU has to perform 3 additional N-Gets).

If the SCU negotiates the Mitra Detached Patient Management SOP rather than, or in addition to, the standard DICOM Detached Patient Management, *Broker* will use the private SOP class which sends one response that also contains the info from the parent, all in the one DICOM object.

Broker supports the following elements for this SOP class:

Table 67: Mitra Detached Patient Management Object N-Get Attributes

Attribute Name	Tag
Specific Character Set	(0008,0005)
Patient Name	(0010,0010)
Patient ID	(0010,0020)
Issuer of Patient ID	(0010,0021)
Other Patient Ids	(0010,1000)
Other Patient Names	(0010,1001)
Patient Telephone Numbers	(0010,2154)
Patient Address	(0010,1040)
Patient Birth Date	(0010,0030)
Patient Sex	(0010,0040)
Patient Weight	(0010,1030)
Ethnic Group	(0010,2160)
Patient Religious Preference	(0010,21F0)
Patient Data Confidentiality Constraint Desc.	(0040,3001)
Patient State	(0038,0500)

Document # MM10001, Rev 1.0

Pregnancy Status	(0010,21C0)
Medical Alerts	(0010,2000)
Contrast Allergies	(0010,2110)
Special Needs	(0038,0050)
Referenced Study Sequence	(0008,1110)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Visit Sequence	(0008,1125)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Reference Patient Alias Sequence	(0038,0004)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)

Broker returns one of the following status codes.

Table 68: Mitra Detached Patient Management status codes.

Service Status	Further Meaning	Protocol Codes	Related Fields	Description
Error	Failed	0110		The operation was not successful.
Success	Success	0000		Operation performed properly.

5.1.1.3. Real World Activity – Mitra Detached Visit Management

5.1.1.3.1. Associated Real World Activity – Mitra Detached Visit Management

5.1.1.3.1.1. Activity as SCP

Broker will respond to DIMSE N-GET requests to retrieve visit information. Broker does not support any other DIMSE-N command as an SCP with this presentation context.

5.1.1.3.2. Presentation Context Table – Mitra Detached Visit Management

Broker supports the transfer syntaxes listed in Table 69.

Table 69: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 70: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Mitra Detached Visit Management	1.2.840.113532.3500.11	all from Table 69	SCP	None

Document # MM10001, Rev 1.0

5.1.1.3.3. SOP Specific Conformance – Mitra Detached Visit Management

Broker provides conformance to a Mitra Private SOP class used to fetch visit and patient demographic information. The purpose for the private SOP class is to minimize the overhead that goes along with the normalized DICOM model (ie: you send something at the results level, and to obtain patient level info, the SCU has to perform 3 additional N-Gets).

If the SCU negotiates the Mitra Detached Visit Management SOP rather than, or in addition to, the standard DICOM Detached Visit Management, *Broker* will use the private SOP class which sends one response that also contains the info from the parent, all in the one DICOM object.

Broker supports the following elements for this SOP class:

Table 71: Mitra Detached Visit Management Object N-Get Attributes

Attribute Name	Tag
Specific Character Set	(0008,0005)
ALL Attributes from Table 67	
Institution Name	(0008,0080)
Institution Address	(0008,0081)
Admission ID	(0038,0010)
Issuer of Admission ID	(0038,0011)
Visit Status ID	(0038,0008)
Current Patient Location	(0038,0300)
Patient's Institution Residence	(0038,0400)
Admitting Date	(0038,0020)
Admitting Time	(0038,0021)
Admitting Diagnosis Description	(0008,1080)
Discharge Date	(0038,0030)
Discharge Time	(0038,0032)
Discharge Diagnosis Description	(0038,0040)
Route of Admissions	(0038,0016)
Referring Physician's Name	(0008,0090)
Referring Physician's Address	(0008,0092)
Scheduled Admission Date	(0038,001A)
Scheduled Admission Time	(0038,001B)
Scheduled Patient Institution Residence	(0038,001E)
Scheduled Discharge Date	(0038,001C)
Scheduled Discharge Time	(0038,001D)
Referenced Patient Sequence	(0008,1120)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Study Sequence	(0008,1110)
>Referenced SOP Class UID	(0008,1150)

>Referenced SOP Instance UID	(0008,1155)
------------------------------	-------------

Broker returns one of the following status codes.

Table 72: Mitra Detached Visit Management status codes.

Service Status	Further Meaning	Protocol Codes	Related Fields	Description
Error	Failed	0110		The operation was not successful.
Success	Success	0000		Operation performed properly.

5.1.1.4. Real World Activity – Mitra Detached Study Management

5.1.1.4.1. Associated Real World Activity – Mitra Detached Study Management

5.1.1.4.1.1. Activity as SCP

Broker will respond to DIMSE N-GET requests to retrieve study information. Broker does not support any other DIMSE-N command as an SCP with this presentation context.

5.1.1.4.2. Presentation Context Table – Mitra Detached Study Management

Broker supports the transfer syntaxes listed in Table 73.

Table 73: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 74: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Mitra Detached Study Management	1.2.840.113532.3500.13	all from Table 73	SCP	None

5.1.1.4.3. SOP Specific Conformance – Mitra Detached Study Management

Broker provides conformance to a Mitra Private SOP class used to fetch study, visit and patient demographic information. The purpose for the private SOP class is to minimize the overhead that goes along with the normalized DICOM model (ie: you send something at the results level, and to obtain patient level info, the SCU has to perform 3 additional N-Gets).

If the SCU negotiates the Mitra Detached Study Management SOP rather than, or in addition to, the standard DICOM Detached Study Management, Broker will use the private SOP class which sends one response that also contains the info from the parent, all in the one DICOM object.

Broker supports the following elements for this SOP class:

Table 75: Mitra Detached Study Management Object N-Get Attributes

Attribute Name	Tag
Specific Character Set	(0008,0005)
ALL Attributes from Table 71	
Study ID	(0020,0010)
Study ID Issuer	(0032,0012)

Accession Number	(0008,0050)
Study Instance UID	(0020,000D)
Study Status ID	(0032,000A)
Study Priority ID	(0032,000C)
Scheduled Study Start Date	(0032,1000)
Scheduled Study Start Time	(0032,1001)
Scheduled Study Stop Date	(0032,1010)
Scheduled Study Stop Time	(0032,1011)
Scheduled Study Location	(0032,1020)
Scheduled Study Location Application Entity Title	(0032,1021)
Requesting Service	(0032,1033)
Requesting Physician	(0032,1032)
Requested Procedure Description	(0032,1060)
Requested Procedure Code Sequence	(0032,1064)
>Code Value	(0008,0100)
>Coding Scheme Designator	(0008,0102)
>Code Meaning	(0008,0104)
Study Arrival Date	(0032,1040)
Study Arrival Time	(0032,1041)
Study Date	(0008,0020)
Study Time	(0008,0030)
Study Completed Date	(0032,1050)
Study Completed Time	(0032,1051)
Study Verified Date	(0032,0032)
Study Verified Time	(0032,0033)
Study Read Date	(0032,0034)
Study Read Time	(0032,0035)
Name of Physician(s) Reading Study	(0008,1060)
Reason For Study	(0032,1030)
Referenced Patient Sequence	(0008,1120)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Visit Sequence	(0008,1125)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)
Referenced Results Sequence	(0008,1100)
>Referenced SOP Class UID	(0008,1150)

>Referenced SOP Instance UID	(0008,1155)
Referenced Study Component Sequence	(0008,1111)
>Referenced SOP Class UID	(0008,1150)
>Referenced SOP Instance UID	(0008,1155)

Broker returns one of the following status codes.

Table 76: Mitra Detached Study Management status codes.

Service Status	Further Meaning	Protocol Codes	Related Fields	Description
Error	Failed	0110		The operation was not successful.
Success	Success	0000		Operation performed properly.

5.1.1.5. Real World Activity – Mitra Detached Results Management

5.1.1.5.1. Associated Real World Activity – Mitra Detached Results Management

5.1.1.5.1.1. Activity as SCP

Broker will respond to DIMSE N-GET requests to retrieve results information. Broker does not support any other DIMSE-N command as an SCP with this presentation context.

5.1.1.5.2. Presentation Context Table – Mitra Detached Results Management

Broker supports the transfer syntaxes listed in Table 77.

Table 77: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 78: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Mitra Detached Results Management	1.2.840. 113532.3500.15	all from Table 77	SCP	None

5.1.1.5.3. SOP Specific Conformance – Mitra Detached Results Management

Broker provides conformance to a Mitra Private SOP class used to fetch results, study, visit and patient demographic information. The purpose for the private SOP class is to minimize the overhead that goes along with the normalized DICOM model (ie: you send something at the results level, and to obtain patient level info, the SCU has to perform 3 additional N-Gets). If the SCU negotiates the Mitra Detached Results Management SOP rather than, or in addition to, the standard DICOM Detached Results Management, Broker will use the private SOP class which sends one response that also contains the info from the parent, all in the one DICOM object.

Broker supports the following elements for this SOP class:

Table 79: Mitra Detached Results Management Object N-Get Attributes

Attribute Name	Tag
Specific Character Set	(0008,0005)

ALL Attributes from Table 75	
Results ID	(4008,0040)
Results ID Issuer	(4008,0042)
Impressions	(4008,0300)
Referenced Study Sequence	(0008,1110)
> Referenced SOP Class UID	(0008,1150)
> Referenced SOP Instance UID	(0008,1155)
Referenced Interpretation Sequence	(4008,0005)
> Referenced SOP Class UID	(0008,1150)
> Referenced SOP Instance UID	(0008,1155)

Broker returns one of the following status codes.

Table 80: Mitra Detached Results Management status codes.

Service Status	Further Meaning	Protocol Codes	Related Fields	Description
Error	Failed	0110		The operation was not successful.
Success	Success	0000		Operation performed properly.

5.1.1.6. Real World Activity – Mitra Detached Interpretation Management

5.1.1.6.1. Associated Real World Activity – Mitra Detached Interpretation Management

5.1.1.6.1.1. Activity as SCP

Broker will respond to DIMSE N-GET requests to retrieve interpretation information. Broker does not support any other DIMSE-N command as an SCP with this presentation context.

5.1.1.6.2. Presentation Context Table – Mitra Detached Interpretation Management

Broker supports the transfer syntaxes listed in Table 81.

Table 81: Transfer Syntaxes

Transfer Syntax	UID
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 82: Presentation Contexts

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
SOP Class	SOP Class UID			
Mitra Detached Interpretation Management	1.2.840.113532.3500.16	all from Table 81	SCP	None

5.1.1.6.3. SOP Specific Conformance – Mitra Detached Interpretation Management

Broker provides conformance to a Mitra Private SOP class used to fetch interpretation, results, study, visit and patient demographic information. The purpose for the private SOP class is to minimize the overhead that goes along with the normalized DICOM model (ie: you send

Document # MM10001, Rev 1.0

something at the results level, and to obtain patient level info, the SCU has to perform 3 additional N-Gets).

If the SCU negotiates the Mitra Detached Interpretation Management SOP rather than, or in addition to, the standard DICOM Detached Interpretation Management, *Broker* will use the private SOP class which sends one response that also contains the info from the parent, all in the one DICOM object.

Broker supports the following elements for this SOP class:

Table 83: Mitra Detached Interpretation Management Object N-Get Attributes

Attribute Name	Tag
Specific Character Set	(0008,0005)
ALL Attributes from Table 79	
Interpretation ID	(4008,0200)
Interpretation ID Issuer	(4008,0202)
Interpretation Type ID	(4008,0210)
Interpretation Status ID	(4008,0212)
Interpretation Recorded Date	(4008,0100)
Interpretation Recorded Time	(4008,0101)
Interpretation Recorder	(4008,0102)
Interpretation Transcription Date	(4008,0108)
Interpretation Transcription Time	(4008,0109)
Interpretation Transcriber	(4008,010A)
Interpretation Author	(4008,010C)
Interpretation Text	(4008,010B)
Referenced Results Sequence	(0008,1100)
> Referenced SOP Class UID	(0008,1150)
> Referenced SOP Instance UID	(0008,1155)
Interpretation Approver Sequence	(4008,0111)
>Interpretation Approval Date	(4008,0112)
>Interpretation Approval Time	(4008,0113)
>Physicians Approving Interpretation	(4008,0114)

Broker returns one of the following status codes.

Table 84: Mitra Detached Interpretation Management status codes.

Service Status	Further Meaning	Protocol Codes	Related Fields	Description
Error	Failed	0110		The operation was not successful.
Success	Success	0000		Operation performed properly.

6. Support for Extended Character Sets

Broker supports the following character sets:

- ISO-IR 6 (default) Basic G0 Set
- ISO-IR 100 Latin Alphabet No. 1

