LEUKOKIT

IT	Italiano	3
GB	English	8
DE	Deutsch	13
FR	Français	19
ES	Español	25
GR	Ελληνικά	31
PT	Português	37
NL	Nederlands	43

LEUKOKIT

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N.3 – Aghi a farfalla 19G – EMOTEC Srl – 
0051

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LEUKOKIT

Kit for autologous leukocytes purification and labelling

Instructions for use

Introduction

Leukokit is a medical device providing all reagents and materials needed to perform a complete procedure of autologous leukocytes purification from peripheral blood and their labelling with a radiopharmaceutical before re-injection into the patient.

The suggested procedure follows the recommendation of the Isorbe Consensus Protocol: Eur. J.Nucl.Med. (1998) 25:797-799 (Except for the use of the washing buffer)

All reagents and components are sterile and produced according to Good Manufacturing Practice; their use according to instructions allows aseptic manipulation during the cell labelling procedure.

The device has also been validated for its ability to preserve cell viability and chemotactic properties of manipulated leukocytes.

Radiopharmaceuticals should be prepared by the user in a manner which satisfies both radiation safety and pharmaceutical quality requirements. Appropriate aseptic precautions should be taken to comply with requirements of aseptic dispensing guidelines for hospitals.

Sterile Components Provided in the Kit

- 1 Sedimentation Syringe with sealing valve
- 1 Separation Vial
- 1 Vial (10 ml) of Anticoagulant Reagent (ACD)
- 1 Vial (10 ml) of Sedimentation Reagent (GSR)
- 1 Vial (14 ml) of Washing Buffer (NaCl)
- 2 Vented Spikes
- 1 Adaptor with Holder
- 3 Butterfly needles 19G
- 1 30ml Luer Lock Syringe
- 1 10 ml Luer Lock Syringe
- 2 5 ml Syringes
- 1 10 ml Syringe
- 1 2.5 ml Syringe
- 1 Sterile Cloth

Other Materials Provided in the Kit

- 1 Tray/Support
- 2 Luer Lock caps
- 3 Sets of colour-coding labels

Materials and Equipments required but not supplied

- Radiopharmaceutical for the labelling of leukocytes
- Centrifuge capable of spinning Falcon type tubes at 150g
- Radioprotective shields
- Sanitizing swabs pads soaked in 70% alcohol

Functional Characteristics of the Key Components

Separation Vial

The vial (Picture 3) provides a sterile closed environment with three different ports with specific functions.

- The injection port (Picture 3-A) can be repeatedly pierced with needles, and guarantees hermetic resealing following each withdrawal of the needle. This port is used to introduce the leukocyte rich plasma, the wash buffer and the radiopharmaceutical for leukocyte labelling into the vial.
- The air venting filter (green port Picture 3-B) allows equilibration of air pressure between the inside and the outside of the vial whilst preserving sterility of the entering air.
- The one way sealing valve (blue port Picture 3-C) connected to the plastic pipette, allows extraction of the solutions from the vial, preventing inadvertent introduction of air or liquids and re-sealing the vial following every syringe disconnection.
- The plastic pipette inside the separation vial allows accurate and progressive aspiration of the solutions from the vial. Its elastic supporting element allows to push down gradually the tip of the pipette inside the vial, thus allowing complete removal of supernatant without disturbing the cell pellet following radio-labelling (thus improving the radiochemical purity of the preparation) and to remove the labelled cells at the end of the procedure.

Sedimentation Syringe with sealing valve

- The syringe luer lock tip is sealed by its valve, thus providing a sterile closed environment.
- This special two-way valve can be opened only by connecting it to a female luer: i.e. to the vented spikes or to the butterfly needles provided (complete rotation is required).
- Disconnection from the female luer (by counter-clockwise rotation of the syringe), automatically reseals the valve.

Composition of reagents

Anticoagulation reagent (ACD)

Trisodium citrate x 2 H ₂ O	22 g/L
Citric acid x H ₂ O	8 g/L
Glucose	22.3 g/L
pH 4.7 – 5.3	

Sedimentation Reagent (GSR)

Succinylated Gelatine (modified fluid gelatine)	40 g/L
Sodium chloride	7 g/L
pH 7.1 – 7.7	

Washing buffer (NaCl)

Sodium chloride

pH 4.5 – 7.0

9 g/L

Procedure

Warnings and Precautions

Does not require any special storage conditions. Do not freeze.

“Leukokit” can be used only by adequately trained personnel.

Use the sterile components only if the packaging is intact.

Employ the kit components for a single separation and labelling procedure.

Use only the butterfly needles provided in the kit.

Do not use the “Leukokit” after the expiry date printed on the outer box.

Remember always to keep the plunger in the sedimentation syringe in a straight line during the various stages to avoid leaks, especially at either end of the stroke.

Single use device, do not re-use and discard according to current rules. Re-use of the device can be dangerous for the patient and/or harmful to the user, with risk of transmitting infections. Use of the device for a second time can compromise its function. Manufacturer waives any responsibility for damages to the patient deriving from misuse of the device or from improper use not complying with the instructions for use, or from use by personnel not properly qualified.

The use of the radiopharmaceutical product chosen for the leukocyte labelling, must be done following its specific instructions for use, enclosed with the product.

Use the following precautions during the whole procedure:

- Wear disposable plastic gloves
- Do not touch surfaces with the tip of the syringe sealing valve or of the needles.
- Always wipe with a sanitizing swab soaked in 70% alcohol the injection port of the separation vial and the stoppers of the reagent's vials and allow them to dry before piercing them.
- Appropriate aseptic techniques should be used during the labelling of the cells.

Blood collection

- Open and display the sterile cloth to obtain a clean area on the working surface.
- Open the sedimentation syringe package and identify it using a colour-coding label.
- Take note of the patient's name matched with the colour code.
- Pierce the stopper of the anticoagulation reagent vial using a vented spike.
- Remove the cap (light blue) of the vented spike and connect the syringe valve to the spike by rotating the syringe clockwise (Picture 1).
- Introduce 8 ml of anticoagulation reagent (ACD) in the syringe and disconnect it from the spike by rotating counter-clockwise (to avoid spillage, do not invert the vial whilst disconnecting the syringe).
- Connect a 19G butterfly needle to the syringe valve and perform the blood collection (approx. 40 ml). In case you would need to replace the butterfly needle due to difficulties during the blood drawing procedure, use the spare butterfly needle provided in the Kit (Picture 2).
- Pierce the stopper of the sedimentation reagent (GSR) vial using the second vented spike.
- Introduce in the syringe 8 ml of sedimentation reagent repeating the steps above, and mix thoroughly by repeatedly inverting the syringe.

- Position the syringe on the tray using the cardboard support (Picture 4).

Plasma transfer procedure

- Open the separation vial package (Picture 3), identify it with the same colour previously used, using a colour-coding label, and position it on the stand.
- Connect a new butterfly needle to the syringe valve and to the separation vial by piercing the injection port (Picture 4).
- Allow the red blood cells to settle (usually not less than 30 minutes are required).
- Open the cap of the venting filter (green) on the separation vial and keep it open throughout whole procedure.
- Transfer the leukocyte rich plasma to the separation vial (up to a maximum of 30 ml) by pushing down the syringe on the plunger. Do not disturb the red blood cells sediment.

Leukocytes separation and labelling

- Centrifuge the separation vial at 150 g for 5/10 minutes.
- Connect the 30 ml Luer lock syringe to the aspiration valve (blue) and aspirate the supernatant plasma (Picture 5).
- Remove the last part of the supernatant by pushing down the plastic needle (connected to the syringe) whilst aspirating (Picture 6). Leave in the vial approximately 0.5 ml of supernatant plasma. Allow the plastic pipette to return to its original position and continue aspiration until the plastic pipette is completely empty.
- Disconnect the syringe and close the syringe tip using one of the caps supplied.
- Gently shake the pellet until complete cell re-suspension.
- Prepare the radiopharmaceutical for leukocyte labelling according to the instructions of the manufacturer and introduce it in the separation vial via the injection port using the 2.5 ml syringe provided, or any other syringe supplied with the radiopharmaceutical compound (Picture 7).
- Incubate for the prescribed time and at the end, dilute the incubation mixture by adding 4-5 ml of washing buffer (NaCl) via the injection port using a 5 ml syringe
- Centrifuge at 150 g for 5/10 minutes and connect a 10 ml luer lock syringe to the aspiration valve (blue) and aspirate the supernatant containing the unbound radioactivity.
- Remove completely the supernatant by pushing down the plastic pipette (connected to the syringe) whilst aspirating. Do not disturb the cell pellet.
- Close the syringe tip with one of the caps provided. The syringe will contain the unbound radioactivity and can be used to calculate the labelling efficiency.

Patient's dose preparation

- Introduce 3-4 ml of washing buffer (NaCl) in the separation vial through the injection port using a 5 ml syringe
- Gently shake the pellet until complete cell re-suspension.
- Open the adaptor with holder package and the 10 ml syringe package; identify the syringe with same colour previously used, using a colour-coding label.
- Insert completely the needle into the adaptor with holder (Picture 8)
- Withdraw the syringe from the holder with a slight rotation to allow the adaptor to separate itself from the holder and remain positioned on the needle.
- Pierce the septum of the blue valve with the needle and insert it until the adaptor will be fully lodged into the neck of the blue valve (NOTE: pierce the septum only once at its centre).
- Aspirate completely the cell suspension (Picture 9). Continue aspiration until the plastic pipette will be completely empty.

- Eject any air from the syringe inside the separation vial (upside down) before extracting the needle from the blue valve.
- Extract the syringe with its needle to obtain the dose ready for injection.

Celltech srl
Via Alessandria 43/A
10098 Rivoli (TO)
Italy
Tel. (+39) 011 5368932
info@cell-tech.it
www.cell-tech.it

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N.3 - Vlindernaalden 19G – EMOTEC Srl –