Open Monitor Suspension

Service Instruction for Installation

5442831-4-1EN
Rev 1
© 2015, General Electric Company,
All Rights Reserved

Revision History

Part/Rev	Date	Change Description	Pages
5442831-1-1EN rev 1	November 7, 2012	Creation for Discovery IGS System	42
5442831-2-1EN rev 1	December 14, 2012	New release for FDR3A	44
5442831-3-1EN rev 1	December 11, 2014	Update for 58" LDM display	44
5442831-4-1EN rev 1	September 3, 2015	Update for Backout positions and Gantry Parking	48

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WARNING
(EN)

This service manual is available in English only.

- If a customer's service provider requires a language other than English, it is the customer's responsibility to provide translation services.
- Do not attempt to service the equipment unless this service manual has been consulted and is understood.
- Failure to heed this warning may result in injury to the service provider, operator or patient from electric shock, mechanical or other hazards.

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This page left intentionally blank.

Table of Contents

1	Introduction	11
1.1	Disclaimer	12
1.2	Pre-Installation	12
1.3	Spare Parts.....	12
1.4	Reference Documents	13
2	Safety & Regulatory Requirements	14
3	Interconnections Schematics	15
3.1	Power/Signal Schematic	15
3.2	Grounding Schematic	16
4	GE Discovery System Outputs.....	17
4.1	C1 Cabinet outputs to be used	17
4.2	Large Display cabinet outputs to be used	19
5	Installation Procedure.....	20
5.1	Pre-requisites.....	20
5.2	Suspension Installation.....	20
5.3	Monitors Installation	21
5.4	Cabling	24
5.5	Infrared Receiver Installation	29
5.6	Finalization.....	29
6	Grounding Tests.....	30
7	Leakage Current Measurement	31
7.1	Pre-requisites.....	31
7.2	Test Equipment.....	31
7.3	Test Procedure.....	32
7.4	Table to fill	33
8	Appendices.....	35
8.1	Appendix 1: Weight for the maximum load calculation.....	35
8.2	Appendix 2: 19" B&W Video Coax Cable Wiring Instructions	36

8.3	Appendix 3: 19" Power Cable Specification.....	37
8.4	Appendix 4: I/R Cable Specification	37
8.5	Appendix 5: Large Display Power & Video Cables	37
8.6	Appendix 5: Room layout overview	38
8.7	Appendix 7: Tables to fill for Logbook.....	45

1 Introduction

The Discovery™ IGS product with the Open Suspension option allows customers to power two B&W 19" monitors, up to 4 additional in-room 19" monitors (B&W or Color), and up to 2 Large Display (LD) monitors (8 Megapixels) on one or several Open Suspension of their choice, in addition or in replacement of the standard Mavig Suspension usually provided with the system. Third party monitors from external sources can also be installed on these suspensions, but are not powered by the system.

This document lists the GE mandatory requirements in order to provide safe system operation and provides the installation instructions. This Service manual shall be considered an addition to the Service Manual and Pre-Installation Manual of the system.

DISCLAIMER: The association of the Discovery™ IGS system and the purchaser's suspension(s) is not covered by the GEHC product certification.

NOTE: GE Healthcare is "GE Santé" in Province of Quebec - Canada.

General Information

1.1 Disclaimer

GE specifically disclaims any and all liability arising out of or relating to the use or performance of the monitor suspension (including cables) including, and without limitation, any liability or claims relating to patient injury, death, or the reliability of such monitors suspension.

The mechanical installation of the Open Suspension is fully under the customer and the installer responsibility.

The electrical installation of the third party monitors is fully under customer and installer responsibility.

The installer is responsible for ensuring that all requirements from this document are met.

1.2 Pre-Installation

It is the customer responsibility to ensure that the Open Suspension and its cables are installed prior to the GE equipment (gantry, table, cabinet) to ensure that standard GE Service Process can be followed during system installation. Monitors installation and connections to the GE equipment shall only be made in presence of a GE service representative.

It is recommended that the vendor contacts GE Service representative and reviews the site planning details before the suspension is installed. The position of each suspension in the exam room shall be planned following the recommendations from the chapter Room Layout of the Discovery™ IGS systems Pre-Installation Manual, in order to reduce the risk of collision between any fixed part of the suspension and the gantry or table (see illustrations in Appendix §9.6 for gantry and table range of motions).

In addition, the position of each suspension in the exam room shall be compatible with the maximum cable length that shall not exceed 36m (between C1 Cabinet and GE 19" monitors or between LD Cabinet and LD monitor).

Note: GE will not be responsible for any delay in installation if the suspension is not mounted and its cables routed before GE parts arrive on site.

1.3 Spare Parts

The customer is responsible for providing and replacing any parts of the Open Suspension and of the third party monitors.

1.4 Reference Documents

Document:
Discovery™ IGS systems - Service Manual
Discovery™ IGS systems - Pre-Installation Manual

2 Safety & Regulatory Requirements

It is the customer responsibility to ensure that the following requirements are met.

- Each suspension shall not be electrically motorized for up/down motion.
- Each suspension shall comply with the IEC 60601-1 standard and the applicable standards enforced in the country of installation (e.g. when installed in a European Community country, the suspension(s) shall be CE marked). In addition, for North America each suspension shall comply with UL/Canada deviations.
- Each suspension shall be manually adjustable in height and the force to be applied to lift the suspension when fully equipped shall not exceed 200 N in static in the vertical axis, in order to mitigate the risk of patient jammed between table and monitor suspension when the table is lifted up.
- Each suspension shall be installed in order to mitigate the risk of suspension fall on patients and the risk of collisions with the gantry, the table or any other suspension.
- Each suspension shall be attached to the ceiling in accordance to the manufacturer's instructions. It shall withstand the maximum suspension load with safety factors in accordance to applicable standards (at least 4x).
- Each suspension shall be compatible with the Environmental Requirements chapter of the Pre-Installation Manual of the system.
- When the system is installed in an operating room (OR configuration), each suspension shall be compatible with OR environmental constraints.

NOTICE

IT IS MANDATORY TO EXECUTE THE GROUNDING CONTINUITY MEASUREMENTS AS DEFINED IN SECTION 7.


NOTICE

IT IS MANDATORY TO EXECUTE THE LEAKAGE CURRENT MEASUREMENTS AS DEFINED IN SECTION 8.

3 Interconnections Schematics

3.1 Power/Signal Schematic

Note: The maximum length of the power and video cables between the C1 cabinet and each 19" monitor shall be 36 meters, including the cable length inside the suspension. For GE 19" monitors, 24m or 36m cables can be ordered.


Note: Power and video cables for Large Display Monitor are provided by GE and have only a 36 meters length between the Large Display cabinet and each 58" monitor.

The sets of cables provided with this option include the two cables required to connect the IR Receiver to the C1 Cabinet.

Refer to Chapter 6 for detailed cabling information.

3.2 Grounding Schematic


Here is an example of grounding schematic with 3 open suspensions and 2 optional Large Display monitors:


Note: Refer to Chapter 6 for detailed cabling information.

4 GE Discovery System Outputs


4.1 C1 Cabinet outputs to be used

C1 Cabinet Video Bulkhead

C1 Cabinet Power Distribution Unit

Video bulkhead outputs:

Connector	Description	Comments
J26	LIVE 1 (Monitor Suspension Live Fluoro Monitor 1)	Mandatory exam room monitor
J27	LIVE 2 (Control Booth Live Fluoro Monitor 2)	Mandatory control room monitor
J28	LIVE 3 (Optional)	
J24	Large Display monitor LIVE 1	
J30	ROADMAP 1 (Monitor Suspension Road Map Monitor)	Mandatory exam room monitor
J31	ROADMAP 2 (Control Booth Roadmap Slave Monitor)	
J32	ROADMAP 3 (Road Map Monitor Option)	
J54	ROADMAP 4 (for Large Display Option)	
J35	REVIEW 1 (Review monitor 1)	
J36	REVIEW 2 (Review monitor 2)	
J37	REVIEW 3 (Review monitor 3) (not used)	
J55	REVIEW 4 (for Large Display Option)	

4.2 Large Display cabinet outputs to be used

Refer to the following the GE Job Card for Large Display cabinet outputs details:

- IST0324 - Large Display Cabinet Installation

5 Installation Procedure

A MAXIMUM OF SIX 19" B&W AND COLOR MONITORS IN EXAM ROOM CAN BE POWERED BY THE SYSTEM (PDU J9 OUTPUT). DO NOT POWER MORE THAN SIX 19" B&W AND COLOR MONITORS IN EXAM ROOM.

NO OTHER EQUIPMENT (THIRD PARTY MONITORS, SURGICAL LIGHTS OR OTHER DEVICE) SHALL BE POWERED BY THE CONNECTOR J9 OR ANY OTHER POWER SOURCE FROM THE SYSTEM: OTHER CONNECTOR OF THE PDU, CONTROL ROOM POWER STRIP OR MULTIPLE SOCKET OUTLET OF THE SUSPENSION IF IT IS POWERED BY THE SYSTEM.

Up to two optional Large Display monitors can be connected to the system.

Note: Large Display monitors are connected to a Large Display cabinet, providing power and video signals to these monitors.

5.1 Pre-requisites

THE LOCKOUT TAGOUT (LOTO) PROCEDURES MUST BE APPLIED TO ENSURE FULL SAFETY DURING SERVICE INTERVENTION.

REFER TO THE LOCKOUT TAGOUT (LOTO) PROCEDURES OF EACH SUSPENSION TO DO THE LOTO FOR THE MECHANICAL AND ELECTRICAL ENERGIES.

REFER TO THE CORRESPONDING PROCEDURES OF THE DISCOVERY SERVICE MANUAL BEFORE INSTALLATION.

- REG0120 - Electrical Lockout Tagout procedures.
- REG0122 - Large Display Option Lockout Tagout Procedures.

5.2 Suspension Installation

The suspension installation shall be done under the customer's responsibility, following the suspension manufacturer's installation instructions (through the Service Manual or Operator Manual provided with the suspension).

The GE Field Engineer shall verify the following points:

- The load on each suspension must be compatible with the suspension acceptable maximum load. The weight of the GE monitors and IR receiver are given in Appendix 1.
- Ground studs shall not be put on an insulated surface such as anodized or painted one and shall be marked with the Protective Earth symbol (IEC 60417-5019 symbol).

5.3 Monitors Installation

The monitors shall be installed on the Open Suspension following the installation instructions of the suspension and fixed following the installation instructions of the monitors. After monitors' installation, the suspension balance shall be checked. If needed, it shall be adjusted following the installation instructions of the suspension.

NOTICE

FOR THE LARGE DISPLAY MONITOR INSTALLATION ON THE SUSPENSION, LIFTING THE MONITOR MANUALLY REQUIRES AT LEAST THREE PERSONS.


5.3.1 Type of fixation

19" monitors: standard VESA 100 x 100 mm. Use the supplied screws or four M4 x 10 screws.

4x M4 screws
Max insertion depth: 10 mm


Large Display (58") monitors: standard VESA 400 x 400 mm. Use the supplied screws or four M8 screws.

Top (14) and bottom (15) mounting holes for Large Display Monitor

Note: The installation of the LDM on the suspension frame can be made easier if you follow the recommendations below:

1. Have 3 persons lifting the LDM monitor while a fourth one inserts the screws and secure the LDM on the back.
2. Attach the suspension frame to the table without damaging the table top, UIF or accessories in order to avoid lifting the LDM monitor too high (see Illustration below).
3. Lift the LDM monitor using a tool such as the one recommended for the CMS installation (Genie Lift), the LDM being suspended to the forks using the parts provided with the install kit (p/n 5418782) delivered with the LDM.

5.4 Cabling

5.4.1 Interconnection Box

When an Interconnection box is present on the Open Suspension:

- The video cables shielding continuity shall be maintained.
- No live part shall be accessible.

FOR THE LARGE DISPLAY OPTION, NO INTERCONNECTION BOX SHALL BE USED FOR THE VIDEO CABLES.

5.4.2 19" Monitors

NOTICE

THE B&W 19" MONITORS PROVIDED WITH THE SYSTEM FOR THE DISPLAY OF THE REGULATORY INFORMATION (LIVE AND ROADMAP) SHALL BE ACCESSIBLE AND VISIBLE AT FLUOROSCOPIST WORKING POSITION.

NOTICE

WHEN THE LARGE DISPLAY OPTION IS PRESENT, THE TWO B&W 19" MONITORS PROVIDED WITH THE SYSTEM FOR THE DISPLAY OF THE REGULATORY INFORMATION (LIVE AND ROADMAP) SHALL BE ACCESSIBLE AND VISIBLE AT FLUOROSCOPIST WORKING POSITION IN LESS THAN ONE MINUTE. THOSE BACKUP MONITORS MAY BE MOUNTED EITHER AT THE BACK OF THE LARGE DISPLAY MONITOR SUSPENSION OR ON ANOTHER MONITOR SUSPENSION.

NOTICE

ONLY THE 19" B&W MONITORS (P/N 5492442) PROVIDED BY GE SHALL BE POWERED BY THE SYSTEM (C1 CABINET).

5.4.2.1 Power and grounding

Monitors power input:

Ground connection:

The power cables shall be provided by the installer and shall respect the following characteristics. Only this type of cable will guarantee the safety of the installation:

- The connector on the monitor side shall be IEC 320 Male Plug (over molded).
- The cable shall be 3 x 2.5 mm² minimum of any UL type 300V cable. The maximum length shall be 36 meters between the C1 cabinet and each monitor.
- The connector on cabinet side shall be: HUBELL Twist Lock HBL4570C. It shall be connected to J9 on the PDU of C1 cabinet.

See Appendix 3 of this document for the wiring of the cable on cabinet side.

**NOTICE**


EACH B&W OR COLOR 19" MONITOR SHALL BE GROUNDED WITH AN ADDITIONAL PROTECTIVE EARTH CABLE BETWEEN THE MONITOR AND THE GROUND BAR LOCATED AT THE BOTTOM OF THE C1 CABINET. THIS CABLE SHALL HAVE THE FOLLOWING CHARACTERISTICS:

- AWG10 minimum, 600V of any UL type.
- The maximum length shall be 36 meters between the C1 cabinet and each monitor.

5.4.2.2 Video (B&W Monitors)

It is recommended to use the video cables provided by GE from the cabinets to the monitors. If needed, additional cables shall be provided by the installer and shall respect the following characteristics to guarantee the image quality performances and EMC compliance.

- The cables shall be coax RG59U / KX-6 type, AWG24, double shielded, UL approved 300V. The maximum length shall be 36 meters between the C1 cabinet and each monitor.
- The connectors mounted at each extremity of these cables shall be BNC 75 ohms (ref 221185-9 AMP/TYCO or equivalent).
- These cables shall be preferably one-piece cables.
- The shield of the coax cable shall be grounded with clamps on the monitor side and on the cabinet side, as on the pictures below. Pay attention that if GE cables are used, clamps are not provided.


The video cable shall be connected to the input below:

See Appendix 2 of this document for the wiring instructions of the video cable.

5.4.2.3 Video (Color Monitors)

The video cables shall be provided by the installer, and shall be 15 pin VGA.

5.4.3 Large Display Monitors

NOTICE

ONLY THE CABLES PROVIDED BY GE SHALL BE USED TO CONNECT THE LARGE DISPLAY MONITOR. THE CONNECTION OF THE ADDITIONAL GROUND WIRE FROM THE LARGE DISPLAY CABINET TO THE MONITOR IS MANDATORY.

The GE Field Engineer shall connect the GE Large Display cabinet following the Job Card cabling instructions:

- IST0324 - Large Display Cabinet Installation

The Large Display Monitor shall be connected following the job card cabling instructions:

- IST0321 - Large Display Monitor installation on suspension (not all sections are applicable)

In case of a secondary Large Display Monitor, it shall be connected following the job card cabling instructions:

- IST0325 - Large Display Monitor Installation on articulating arm wall mount (not all sections are applicable)

To avoid damage to the video cables, the minimum bending radius of these cables shall be 35 mm.

The cables for Large Display are provided by GE. See Appendix in §9.5 for detailed characteristics.

5.4.4 Third Party Monitors

THIRD PARTY MONITORS INSTALLED ON THE SUSPENSION MUST NOT BE POWERED BY THE SYSTEM.

DANGEROUS VOLTAGE MAY BE PRESENT AT THE INPUT OF THESE MONITORS, EVEN WHEN PDB IS OFF.

THIRD PARTY MONITORS MUST BE GROUNDED TO C1 CABINET GROUND BAR BY AN ADDITIONAL GROUND WIRE, AWG10 MINIMUM, AND 600V OF ANY UL TYPE.

THIRD PARTY MONITORS SHALL BE CERTIFIED IEC60601-1 OR IEC60950.

Refer to the installation manual of the monitor or video source for power and video cabling requirements.

5.4.5 Suspension Grounding

ALL ACCESSIBLE CONDUCTIVE PARTS OF EACH SUSPENSION SHALL BE GROUNDED BY A DEDICATED GROUND WIRE CONNECTED TO THE C1 CABINET GROUND BAR.

5.5 Infrared Receiver Installation

If installed on the monitor suspension, the infrared remote control receiver should be installed on top of the suspension, with the front surface tilted down of 15 degrees as a default position. It is recommended to use self-adhesive Velcro tape (size: 50 mm x 200 mm) to fix the receiver on the suspension.

The cable for the infrared receiver is provided GE and only GE cable shall be used to guarantee the safety of the installation.

See Appendix in §9.4 for detailed characteristics.

Note: Always use the short IRK interface cable (P/N 2390547) to connect the IR Receiver. This cable shall be plugged on the suspension side and not on C1 cabinet side.

Note: If the short IRK cable cannot be routed inside the suspension, the IR Receiver can be mounted on another location independent of the suspension itself. In this case, it shall not be mounted in the patient environment.

5.6 Finalization

MAKE SURE THERE ARE NO POTENTIAL CABLE PINCHES AND SHARP EDGES WHEN ROUTING THE CABLES INSIDE THE SUSPENSIONS.

6 Grounding Tests

NOTICE

IT IS MANDATORY TO EXECUTE THESE GROUNDING CONTINUITY MEASUREMENTS AT THE END OF INSTALLATION IN ORDER TO GUARANTEE THE SAFETY OF THE PATIENT AND THE OPERATOR. THE RESULTS SHALL BE REPORTED INTO THE LOGBOOK.

For all the ground connections introduced by the open suspension:

- Check that the ground wires are screwed correctly on each monitor, on each suspension and on each ground bar (pay attention to the ground cables connected between cabinets and monitor suspensions).
- Check that the ground cabling is compliant with the ground schematic provided in this document.

Perform the ground measurement between the following points as defined in the paragraph 3.3 *Grounding impedance measurement of the CHK0103 JC - Grounding and Neutral Continuity Check for detailed instructions:*

- Between C1 Cabinet ground bar and each suspension metallic part.
- Between C1 Cabinet ground bar and each 19" monitor (B&W and color).
- Between C1 Cabinet ground bar and each Large Display monitor.
- Between C1 Cabinet ground bar and each third party monitor.

The measured value for each measurement shall be less than 0.1 ohms.

Report the results in the Discovery™ IGS systems Logbook.

7 Leakage Current Measurement

NOTICE

THIS LEAKAGE CURRENT MEASUREMENT IS A MANDATORY PROCEDURE THAT SHALL BE DONE AT THE END OF INSTALLATION IN ORDER TO GUARANTEE THE SAFETY OF THE PATIENT AND THE OPERATOR. THE TABLE SHALL BE FILLED WITH THE RESULTS AND BOUND INTO THE LOGBOOK.

7.1 Pre-requisites

Leakage current test results are dependent of cabling and ground connections, all system cabling must be checked with the System MISMAPS and Installation Job Cards.

The measure shall be executed in the worst conditions configuration:

- Make sure that the Table, Gantry and cabinets are not grounded by other cables than those provided by GE.
- In case of multiple exclusive options, the most critical option (regarding leakage current) must be used.
- No equipment shall be fixed to the table rail (TSUIs or other accessory): they shall be placed on the table. The table mount injectors shall be out of the table.
- When the system is equipped with the InnovaSense option, make sure that Capacitive Sensors are installed and 20°/s rotation speed is enabled.


The equipment shall be operated in Normal Conditions during the measurement. The complete system shall be energized and in full operating conditions with maximum functions activated.

7.2 Test Equipment

Tables to be filled:

Equipment	Model type and S/N	Next calibration date
True RMS Digital Voltmeter (DVM) ($\pm 1.5\%$)		
Standard IEC Measurement Device (MD) ($\pm 1\%$)		

Example of measurement device (MD cell):

Note: The Dale 601/601E (46-328406G1) from Dale Technology or the ESA612 Electrical Safety Analyzer from Fluke can be used for this measure (or any other equivalent device which measuring circuit conforms to above IEC60601 MD specifications).
 Be sure to read the tool user’s manual and understand all instructions, in particular those related to the selection of the measuring mode, current range and patient load which shall comply with IEC 60601-1 recommendations.

7.3 Test Procedure

The measurement shall be done when the system is in standby mode, and with the following movements:

- Move the table upwards to its highest position and simultaneously move the detector lift upwards (using the detector button on smart box).
- Move the table downwards to 0cm, then tilt table 15° down (Trendelenburg position). Move longitudinally towards head stop, then simultaneously move slowly the longitudinal axis to the table foot and tilt to the horizontal position.

It is not necessary to perform X-ray during the leakage current measurement.

In case a MD cell is used, the voltmeter shall be used in MAX mode following the current direction, successively in AC and DC mode.

For every measurement, the operator shall perform movements to reach the maximum absolute value. The leakage current value in μA is a direct reading of the voltmeter display in mV.

7.4 Table to fill

Measures between the following components and each suspension:

- **Stdby**= system in standby
- **Mvts**= system in movements as defined in the previous paragraph
- **Susp**= monitor suspension

Test criteria: the measured value shall be < 65µA.

Note: Do not fill the tables in this section. Refer to section 9.7 Appendix 7 to fill the tables, extract the corresponding pages from the binder and add them to the Discovery™ IGS systems Logbook.

From	To	Measured value (µA)								PASS /FAIL / NA
		Suspension 1		Suspension 2		Suspension 3		Suspension 4		
		Stdby	Mvts	Stdby	Mvts	Stdby	Mvts	Stdby	Mvts	
Right Table Rail	Each metallic part of the enclosure of each suspension.									
	Each monitor installed on each suspension.									
Gantry	Each metallic part of the enclosure of each suspension.									
	Each monitor installed on each suspension.									
TSSC	Each metallic part of the enclosure of each suspension.									
	Each monitor installed on each suspension.									

Measures between each suspension:

Test criteria: the measured value shall be < 65µA.

From	To	Measured value (µA)								PASS / FAIL / NA
		Suspension 1		Suspension 2		Suspension 3		Suspension 4		
		enclosure	monitors	enclosure	monitors	enclosure	monitors	enclosure	monitors	
Susp 1	Each metallic part of the enclosure of each suspension.									
	Each monitor installed on each suspension.									
Susp 2	Each metallic part of the enclosure of each suspension.									
	Each monitor installed on each suspension.									
Susp 3	Each metallic part of the enclosure of each suspension.									
	Each monitor installed on each suspension.									

Report the results in the Discovery™ IGS systems Logbook.

8 Appendices


8.1 Appendix 1: Weight for the maximum load calculation

Dimensions are given for information.


Part	Max Weight	Dimensions (W x H x D)
19" monitor (B&W)	6.3 Kg	422.5 mm x 348.5 mm x 95 mm
19" monitor (COLOR)	6.3 Kg	422.5 mm x 348.5 mm x 95 mm
Large Display monitor (58" monitor)	51 Kg	1317 mm x 775 mm x 145 mm
IR Receiver module	0.3 Kg	112 mm x 31 mm x 76 mm (+/-1 mm)

Note: The weight of the length of cables routed inside the monitor suspension has to be taken into account for the computation of overall equipment weight.


8.2 Appendix 2: 19" B&W Video Coax Cable Wiring Instructions

8.3 Appendix 3: 19" Power Cable Specification

8.4 Appendix 4: I/R Cable Specification


8.5 Appendix 5: Large Display Power & Video Cables

Large Display cables characteristics are:

- Video: 4x Ethernet cables with RJ45 connector
- Power: IEC 320 plug on monitor side


8.6 Appendix 5: Room layout overview

The following illustrations are given for informational purpose only. Please refer to the Pre-Installation Manual of Discovery™ IGS systems for more details.


All dimensions in mm (inches)


Potential collision between laser, detector lift and mast/chain


Backout positions (RIRP 1120 mm)


Backout positions (RIRP 1278 mm)


Backout positions (RIRP 1508 mm)

Gantry Parking (RIRP 1120 mm)

Gantry Parking (RIRP 1278 mm)


Gantry Parking (RIRP 1508 mm)

**FRONT VIEW
(HEAD END)**

All dimensions are in mm (inches)

NOTE

The table dimensions are correct for the position (longi=0, lat=0, lift=0, tilt=0, rot=0).

SIDE VIEW

TOP VIEW



Table dimensions

8.7 Appendix 7: Tables to fill for Logbook

Fill the following tables as required in section 8 Leakage Current Measurement, take out pages 41 to 42 from the binder and place these pages in the Discovery™ IGS systems Logbook, in Chapter 14 - Appendix.

Measures between the following components and each suspension:

- **Stdby**= system in standby
- **Mvts**= system in movements as defined in the previous paragraph
- **Susp**= monitor suspension

Test criteria: the measured value shall be < 65µA.

From	To	Measured value (µA)								PASS /FAIL / NA
		Suspension 1		Suspension 2		Suspension 3		Suspension 4		
		Stdby	Mvts	Stdby	Mvts	Stdby	Mvts	Stdby	Mvts	
Right Table Rail	Each metallic part of the enclosure of each suspension.									
	Each monitor installed on each suspension.									
Gantry	Each metallic part of the enclosure of each suspension.									
	Each monitor installed on each suspension.									
TSSC	Each metallic part of the enclosure of each suspension.									
	Each monitor installed on each suspension.									

Measures between each suspension:

Test criteria: the measured value shall be < 65µA.

From	To	Measured value (µA)								PASS / FAIL / NA
		Suspension 1		Suspension 2		Suspension 3		Suspension 4		
		enclosure	monitors	enclosure	monitors	enclosure	monitors	enclosure	monitors	
Susp 1	Each metallic part of the enclosure of each suspension.									
	Each monitor installed on each suspension.									
Susp 2	Each metallic part of the enclosure of each suspension.									
	Each monitor installed on each suspension.									
Susp 3	Each metallic part of the enclosure of each suspension.									
	Each monitor installed on each suspension.									

This page left intentionally blank.

© 2015 General Electric Company. All rights reserved
General Electric Company
283, rue de la Minière
78530, Buc
FRANCE

www.gehealthcare.com