



AHA compliant

SEER Holter Monitor ECG Cables



Quality ECG cables for your Holter devices play an important role in ensuring that they perform at the high level you have come to expect from GE Healthcare equipment.

SEER™ 12 ECG cables

10-leadwire cables are available for 12-channel Holter ECG recordings. Cables are tested and verified as compatible with the SEER 12 and CardioMem™ CM 3000-12 BT monitors.

SEER 1000 ECG cables

Cables are tested and verified as compatible with SEER 1000 24-hour, 48-hour and 7-day Holter monitors. Offerings are:

- 3-leadwire cables for 3-channel ECG recordings
- 5-leadwire cables for 2-channel ECG recordings
- 7-leadwire cables for 3-channel ECG recordings

Our equipment.
Our accessories.
Simply compatible.

You can trust GE Healthcare accessories and consumables to keep your clinical devices working at peak efficiency.

- Simplified ordering process
- Fast delivery
- One contact point for expert advice
- Quality support at every step

PRODUCT SPECIFICATIONS

Materials

Does not contain natural rubber latex, or DEHP
RoHS (Restriction of hazardous substances) compliant (Directive 2011/65/EU)
REACH conformity (EU Regulation 1907/2006)

Operating conditions

Recommended ambient humidity 5 to 95%, non-condensing
Recommended temperature 5° to 45°C (41° to 113° F)

Storage conditions

Recommended ambient humidity 0 to 90%, non-condensing
Recommended temperature -25° to 70°C (-13° to 158°F)

Cleaning solutions

Ethyl alcohol 70-80%
Isopropyl alcohol and n-propyl alcohol 70-80%
Sodium hypochlorite (5% household bleach): Minimum 1:500 dilution (minimum 100 ppm free chlorine) and maximum dilution 1:10
Glutaraldehyde 2% (pH 7.5-8.5)

Do not sterilize or autoclave

Compliance

American Heart Association (AHA) standards
CE Mark

Imagination at work

ORDERING INFORMATION

Part Number	Description
SEER 12	
2107228-001	12-Channel Holter ECG Cable, Round, AHA, 130 cm (51.2 in)
2107229-001	12-Channel Holter ECG Cable, Flat AHA, 105 cm (41.3 in)
SEER 1000	
2104821-001	SEER 1000 Leadwire set, 3 leadwire, 3 CH, 75 cm (29.5 in), AHA
2104822-001	SEER 1000 Leadwire set, 3 leadwire, 3 CH, 105 cm (41 in), AHA
2104814-001	SEER 1000 Leadwire set, 5 leadwire, 2 CH, 75 cm (29.5 in), AHA
2104816-001	SEER 1000 Leadwire set, 5 leadwire, 2 CH, 105 cm (41 in), AHA
2104815-001	SEER 1000 Leadwire set, 7 leadwire, 3 CH, 75 cm (29.5 in), AHA
2104817-001	SEER 1000 Leadwire set, 7 leadwire, 3 CH, 105 cm (41 in), AHA

Variations in length mean variations in the trunk cable. Lead length in all variants is 45 cm / 17.7 in (overall).

Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.com/promotional-locations

Distributed by:

GE Healthcare Finland Oy
Kuortaneenkatu 2
00510 Helsinki, Finland
Europe



Manufactured by:

GETEMED Medizin- und
Informationstechnik AG
Oderstraße 77
14513 Teltow, Germany