

# Imaging Protocol Manager

## The right exam. Every time. Effortlessly.



### Improve performance through standardization

#### The challenge

Standardize imaging device protocols to:



Reduce variation



Meet accreditation requirements



Help improve imaging consistency

Until now, there hasn't been a simple way to manage the process.

### The solution: Imaging Protocol Manager

#### How it works

Setting up and using Imaging Protocol Manager is simple:



- 1. Pull**—Import protocols to cloud library from CT and MR devices; allow basic edit
- 2. Publish**—View protocols and publish them as “Standard Protocol” for the fleet of same devices
- 3. Push**—Distribute protocols to scanners throughout the hospitals within the organization
- 4. Clinical instruction**—Verify positioning guidelines, contrast delivery, and technologist instructions
- 5. Deviation tracking**—Track and monitor protocol deviation between device and standard protocols on the cloud

This protocol management solution helps ensure all imaging devices are using the same protocols, so that providers can effortlessly deliver the right exam for each patient, every time.

- Cloud based
- Multi-modal
- Enterprise-wide<sup>1</sup>
- Simple to use

#### Outcomes

Imaging Protocol Manager can help:

- ✓ **Manage compliance**
  - Accreditation
  - Regulatory compliance
  - Payer reimbursements
- ✓ **Increase efficiency**
  - Fewer repeat exams
  - Reduced training time
  - Less manual protocol management
- ✓ **Transform culture**
  - Better morale
  - Increased focus on patient care

You can't be everywhere at once, but now your protocols can.

<sup>1</sup> Only available for CT and MR modalities only  
JB59811XXa