


CardioSoft Diagnostic System ECG Viewer

Efficient solution for ECG data viewing, editing and archiving.

CardioSoft™ Diagnostic System ECG Viewer transforms your GE MAC® resting electrocardiographs into an electronic storage system that efficiently organizes your ECGs and enables easy access for ECG viewing, editing and archiving.

- Easy-to-use yet comprehensive procedure review workstation that enhances ECG device workflow
- Access to a portfolio of analytical tools, including calipers, median display, arrhythmia review, vectors and lead measurements for detailed review of ECG test data
- Fast access to tens of thousands of patients and procedures using the searchable patient database
- Runs on most Windows®-based PCs and local area networks (LANs), maximizing your existing IT investment
- Scalable to meet the needs of your growing practice by adding clinical procedures and export capabilities





Workstation specifications

Microprocessor	Minimum: Pentium 4 \geq 1.6 GHz (Windows XP); 2 GHz (Windows Vista)
RAM	512 MB (Windows XP); 1 GB (Windows Vista)
Hard drive	4 GB (depending on number of tests to save); 50 MB of free memory
SW installation	CD-ROM drive
Pointer	Mouse
Graphics adapter	Minimum: SVGA 1024 x 768 Recommended: SXGA 1280 x 1024
Interfaces	Minimum: 2 USB ports, CD-RW, SD card, network interface card (recommended)
Operating system	Windows XP/Pro SP 2/3; Windows Vista Home Premium/Business SP 1; Windows Server 2000/2003
Printer	PLC 6 driver compatible; US letter or A4 paper; USB or parallel interface GE Healthcare recommends HP4250 network printers, HP4650 or HP LaserJet® P3005 dn for use with CardioSoft Diagnostic System and networked printing from CASE systems
Additional software for export functionality	Microsoft® Word and Excel® (optional)
Networking LAN	Wireless: 802.11G (optional) TCP/IP interface
Citrix®	Citrix application software not sold, installed or supported by GE Healthcare

Please note that some GE products may not be available in all regions due to local language or registration requirements. Please contact your GE representative in the region for the latest info on the availability of the products in your region.

GE Healthcare
Munzinger Straße 5
79111 Freiburg, Germany
Tel. +49 761 4543 0 • Fax +49 761 4543 233
www.gehealthcare.com



GE imagination at work

©2009 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.

GE, GE Monogram, CardioSoft and TONOPOINT are trademarks of General Electric Company.

GE Healthcare, a division of General Electric Company.

Citrix is a trademark of Citrix Systems, Inc.

HP LaserJet is a trademark of Hewlett-Packard Company.

Microsoft, Excel, Windows are trademarks of Microsoft Corporation.

Pentium is a trademark of Intel Corporation.

EMEA 2049702-002/1109
(Global 11926/MW)