



GE Medical Systems
Information Technologies

THE DINAMAP® DIFFERENCE™
– THINK 1, 2, 3

1.

**Accuracy monitoring
hypertensive patients**

With its unique algorithm, the DINAMAP PRO Series ensures exceptionally accurate BP readings in hypertensive situations.

2.

**Accuracy monitoring
in hypotensive patients**

Advanced dual hose technology enhances the sensitivity of the DINAMAP PRO Series providing precise readings in hypotensive cases.

3.

**Accuracy maintained even
if the patient moves**

The patented motion artefact rejection system built-in to every DINAMAP PRO Series monitor delivers complete accuracy even when the patient moves.

DINAMAP® PRO SERIES BLOOD PRESSURE, OXIMETRY & TEMPERATURE



Introduced to meet the ever more demanding requirements of today's healthcare professionals, the new DINAMAP® PRO Series incorporates the three market leading technologies for blood pressure, oximetry and temperature into a unique series of monitors.

Advanced technologies

The unmatched DINAMAP non-invasive blood pressure system, IVAC®* temperature monitoring and NELLCOR®** pulse oximetry are combined for the first time to produce a series of monitors that are as easy to use as they are powerful.

It's what's inside that counts

Earlier DINAMAP blood pressure monitors set the standard for accuracy, reliability and ease of use and the new PRO Series retains all of the features that have proven such a huge success:

- The patented motion artefact rejection system, which ensures accurate blood pressures determinations even when the patient moves;
- The advanced DINAMAP algorithm and systolic search facilities, which deliver accurate readings when measuring hypertensive patients;
- The dual hose technology that enhances monitor sensitivity in extreme hypotensive situations.

And of course, the intelligent inflation system maximises patient comfort at all times.

*IVAC is a trademark of Alaris Medical Systems
**NELLCOR is a trademark of Mallinckrodt, Inc.
Use only IVAC® probes and probe covers. The size, shape, and thermal characteristics of the probe covers can affect the performance of the instrument.
Inaccurate readings or retention problems may occur unless IVAC probes and probe covers are used.



DINAMAP® PRO SERIES BLOOD PRESSURE, OXIMETRY & TEMPERATURE

- One button press for a single blood pressure determination.
- One button press gives you a 'snap shot' reading of all parameters ('ALL OBS' function).
- One button press takes you to 'AUTO BP' cycle to set BP time intervals.
- NELLCOR pulse oximetry is activated automatically, as soon as the probe is attached to the patient.
- IVAC temperature is activated automatically, as soon as the probe is removed from its holder.
- User-adjustable target inflation pressure. You can now limit the cuff inflation pressure.
- Built-in printer on all models, so you can attach the print-out to the patient notes.
- Pole stand and basket included with all models.

Product selector

Model	Order code	Parameter set
DINAMAP® PRO 100	DP100UK	Blood pressure, heart rate
DINAMAP® PRO 200	DP200UK	Blood pressure, heart rate, temp
DINAMAP® PRO 300	DP300UK	Blood pressure, heart rate and oxygen saturation
DINAMAP® PRO 400	DP400UK	Blood pressure, heart rate, temp and oxygen saturation

BP Cuffs

DURA-CUF®	Order code	NHS Supplies code*	Approx. limb circumference
Thigh	2796	FFF056	38 - 50cm
Large adult	2791	FFF054	31 - 40 cm
Adult	2774	FFF055	23 - 33 cm
Small adult	2779	FFF051	17 - 25 cm
Child	2781	FFF050	12 - 19 cm
Infant	2783	FFF049	8.3 - 15 cm

(* = UK only)

Accessories

Item	Order code	Descriptions
SpO ₂ Sensor	DS100A	Oxygen saturation sensor (adult finger)
Printer paper	089100	Printer paper, 10 rolls
Oral temp probe	088012	IVAC Temp probe
Temp Covers	088015	IVAC Temp probe covers

Specifications

Dimensions: 25cm x 25cm x 18cm (H x W x D); Weight: 3.5 kg (including battery);
 NIBP accuracy: Meets or exceeds ANSI/AAMI standard SP-10;
 Default target: Adult: 160 ± 15mmHg; Cuff inflation, Neonatal: 110 ± 15mmHg;
 Target cuff inflation: Adult: 100 to 250mmHg, Neonatal: 100 to 140mmHg;
 Systolic range: Adult: 30 to 245mmHg, Neonatal: 40 to 140mmHg;
 Diastolic range: Adult: 10 to 195mmHg, Neonatal: 20 to 100mmHg;
 MAP range: Adult: 15 to 215mmHg, Neonatal: 30 to 115mmHg;
 Pulse rate accuracy: NIBP: ± 3.5%; SpO₂: ± 3beats/min;
 Pulse rate range: NIBP: Adult: 30-200 beats/min, NIBP: Neonatal: 30-200 beats/min, SpO₂: 20 to 250 beats/min;
 Temperature accuracy: Monitor mode: ±0.1°C, Predictive mode: ±0.6°C;
 Temperature range: Monitor mode: Max: 42.2°C, Min: 26.7°C;
 Predictive mode: Max: 42.2°C; Min: 31.6°C;
 Pulse oximeter accuracy Adult: 70% to 100% ± 3.5 digits, 0 to 69% unspecified;
 Neonatal: 70% to 100% ± 4.5 digits, 0 to 69% unspecified;
 Power requirements: AC input voltage: 100-240 VAC, 50/60 Hz (nominal); Internal battery: Sealed lead acid 12V 2.3Ah

GE Medical Systems *Information Technologies* reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation

For further information on DINAMAP® Blood Pressure monitors and accessories
 please call: Tel. +44 (0) 12 56 419 000



GE Medical Systems
Information Technologies

European Headquarters
 GE Medical Systems
 Information Technologies GmbH
 P.O. Box 60 02 65
 79032 Freiburg • Germany
 Tel. +49 761 45 43 - 0
 Fax +49 761 45 43 - 233

World Headquarters
 GE Medical Systems
 Information Technologies, Inc.
 8200 West Tower Avenue
 Milwaukee, WI 53223 • USA
 Tel. +1 414 355 5000
 Fax +1 414 355 3790

Asia Pacific
 GE Marquette Medical Systems
 11th Floor, The Lee Gardens, 33 Hysan Ave.
 Causeway Road • Hong Kong
 Tel. +852 2100 6300
 Fax +852 2100 6292

gemedicalsystems.com