

GE Entropy Sensor and GE Entropy Cable



GE Healthcare Entropy™ measurement is indicated for monitoring the state of the brain of adult and pediatric patients aged 2 years and older within a hospital.

The disposable GE Entropy Sensor and reusable GE Entropy Cable are accessories for Entropy measurement. The sensor acquires the electroencephalograph (EEG) and frontal electromyography (FEMG) signals which are processed into two Entropy parameters, State Entropy (SE) and Response Entropy (RE).

GE Entropy Sensor is backwards and forwards compatible with the GE Entropy measurement.

Features

- Single-point connection enables ergonomic and easy set-up
- Low and stable impedance provides good signal quality
- Wet-gel electrodes maximize patient comfort
- Leadwires offer flexible positioning of the sensor



GE Entropy Sensor and Cable

GE Entropy Sensor specifications

- Disposable, single-patient use
- 3 pre-attached leadwires and one connector
- For adults and children 2 years and older
- Soft wet-gel electrodes
- Skin preparation required
- GE Entropy Sensor must be used with GE Entropy Cable
- Can be safely used during electrosurgery and defibrillation
- Maximum usage time - 24 hours

Maintenance

GE Entropy Sensor for single use only – do not reuse.

GE Entropy Cable can be wiped with mild detergent or disinfectant. You may also use Ethylene Oxide (ETO) or Gas Plasma Technology (STERRAD).

Note: For detailed use information, always refer to the Instructions for Use accompanying the product.

Ordering Information

Description	Order code
GE Entropy Sensor (25/pkg)	M1038681
GE Entropy Cable, 3.5 m/11.5 ft (1/pkg)	M1050784
GE Entropy Start-up kit:	M1060609
GE Entropy Sensor (2 boxes)	
GE Entropy Cable (1pc)	



GE Healthcare
P.O. Box 900, FIN-00031 GE, Finland
Tel. +358 10 394 11
Fax +358 9 146 3310

www.gehealthcare.com

© 2008 General Electric Company – All rights reserved.
GE and GE Monogram are trademarks of
General Electric Company.

Entropy is a trademark of GE Healthcare Finland Oy.

GE Healthcare Finland Oy, a General Electric company,
doing business as GE Healthcare.

