

S/5™ Anesthesia and Critical Care Display Options

Multiple display options for maximum flexibility and optimal viewing



Features

- Each high-resolution display option shows up to 8 selectable waveforms in user-configurable colors
- Support for user-configurable second and third displays with different waveform and numerical information
- Secondary display user-configurable with or without main display's menus and alarms
- Five display options
 - Datex-Ohmeda 10 inch flat panel LCD display with VGA resolution with integrated Command Board and ComWheel™
 - Datex-Ohmeda 15 inch flat panel LCD display with XGA resolution. Removable ComBar™ and ComWheel
 - Datex-Ohmeda 15 inch video display with removable ComBar and ComWheel
 - 17 inch video display with integrated Command Board and ComWheel
 - Large 21 inch secondary display
- Unique command board and bars, specific to either Critical Care or Anesthesia, for intuitive access to monitoring functions
- Removable command bar brings monitor controls close to the user
- 10 inch or 15 inch flat panel LCD displays are a space-saving solution for main display or as a Record Keeping display

Technical specifications

Datex-Ohmeda S/5™ 15 inch color video display

15 inch color video display with removable ComBar with monitor controls, menu keys and ComWheel for menu selections and adjustments.

Display size:	15 inch diagonal
Display type:	Low radiation color video CRT
Display resolution:	High resolution, 1984 x 512
Number of traces:	Up to 8
Dimensions (W x D x H):	350 x 430 x 320 mm/ 14.6 x 17.7 x 13.0 in
Weight:	14.3 kg/31.6 lbs
Power supply to display:	Power supply from Central Unit

Datex-Ohmeda S/5™ 17 inch color video display

17 inch color video display with display base Command Board with monitor controls, menu keys and ComWheel for menu selections and adjustments.

Display size:	17 inch diagonal
Display type:	Low radiation color video CRT
Display resolution:	High resolution, 1984 x 512
Number of traces:	Up to 8
Dimensions (W x D x H):	412 x 445 x 415 mm/ 16.2 x 17.5 x 16.3 in
Weight:	19.0 kg/41.4 lbs
Mounting:	Only wall mounted
Power supply to display:	Power supply from Central Unit or external approved isolation transformer (local solutions). 120 VA
Rated voltages and frequencies:	100-120 V, 50/60 Hz or 220-240 V, 50/60 Hz

Datex-Ohmeda large 21 inch video display

Display size:	21 inch diagonal
Display type:	Low radiation, high contrast, color, flat face CRT with antiglare surface
Display resolution:	High resolution 1984 x 512
Number of traces:	Up to 8
Dimension (W x D x H):	515 x 544 x 490 mm/ 20.3 x 21.4 x 19.3 in
Weight:	32 kg/70 lbs
Power supply to secondary display:	Local solutions, approved isolation transformer 160 W
Operating voltage:	90 V-260 VAC
Operating frequency:	50/60 Hz

Datex-Ohmeda S/5™ 15 inch color LCD display

15 inch color LCD display with removable ComBar with monitor controls, menu keys and ComWheel for menu selections and adjustments.

Display size:	15 inch diagonal
Display type:	Active matrix color TFT LCD
Display resolution:	XGA, 1024 x 768
Number of traces:	Up to 8
Dimensions (W x D x H):	410 x 100 x 355 mm/ 16.1 x 3.9 x 14.0 in
Weight:	5 kg/11 lbs
Mounting:	Display mounting arm with adapters to poles or rails or desktop
Power supply to display:	Power supply from Central Unit

Datex-Ohmeda S/5™ 10 inch color LCD display

10 inch color LCD display with integrated Command Board with monitor controls, menu keys and ComWheel for menu selections and adjustments.

Display size:	10.4 inch diagonal
Display type:	Active matrix LCD
Display resolution:	VGA, 640x480
Number of traces:	Up to 8
Dimensions (W x D x H):	315 x 74 x 265 mm/ 12.4 x 2.9 x 10.5 in
Weight:	2.9 kg/6.4 lbs
Mounting:	Display mounting arm with adapters to poles or rails
Power supply to display:	Power supply from Central Unit

Ordering information

Description	Order code
S/5 15 inch video display	D-VMC15
S/5 15 inch LCD display	D-LCC15
S/5 Command Bar for Anesthesia	K-ANEB
S/5 Command Bar for Critical Care	K-ICUB
S/5 17 inch video display	D-VHC17-S5
S/5 Command Board for Anesthesia	K-ANE-S5
S/5 Command Board for Critical Care	K-ICU-S5
21 inch video display	D-445 PRO
S/5 10 inch LCD display (ANE)	D-LCC10A-S5



Datex-Ohmeda Division • Instrumentarium Corp.
P.O. Box 900, FIN-00031 Datex-Ohmeda, Finland
Tel. +358 10 39411 • Fax +358 9 146 3310

www.datex-ohmeda.com

Datex-Ohmeda, Inc.
P.O. Box 7550, Madison, WI 53707-7550, USA
Tel. 800 345 2700 • Fax 608 221 4384

In the USA: www.us.datex-ohmeda.com