



MUSE 7 Technical Service Class

Course Summary

The MUSE 7 Technical Service Class is an eight-day training course covering the function and interconnectivity of the MUSE 7 system, including MUSE client workstations, MUSE servers, MAC5000/5500 carts, MARS PC, and beside monitors, focusing on the Dash.

The course is designed to teach the BioMed, Field Engineer, System Administrator, and IT Professional practical skills needed for troubleshooting and System setup. The class is designed to be 40% instruction based and 60% hands-on labs.

Course Length

- 8 days

Required IT Skills

The following skills are required for successful completion of this class: IP Addressing, Subnetting, Hard Drives, Hardware Replacement, Hubs, Switch Port Configuration, Switch VLAN Configuration, Switch Troubleshooting, Router Configuration, Router Access Control Lists, Router Troubleshooting, POST, Access Point Configuration, Network Connectivity, Communication Paths, Network Troubleshooting and Back-up Systems.

Course Location

- Waukesha, WI

Course Objectives

Upon successful completion of the course, participants will be able to:

- Setup and configure the MUSE 7 Server
- Setup and configure the MUSE 7 Client workstation
- Configure the MAC 5000/5500 cart and send resting ECG's to the MUSE 7
- Setup Sites, Locations, Devices and Report Distribution
- Setup Muse Users
- Setup modem, fax, and printer devices on Server and Client
- Review system messages and event logs
- Perform system shutdown procedures
- Backup the MUSE file server
- Reload Muse Software on Server and Client Workstation
- Setup and configure the Patient Monitoring equipment to communicate with the MUSE file server.



- Setup and configure the Case Stress equipment to communicate with the MUSE file server.
- Setup and configure the MARS Holter equipment to communicate with the MUSE file server.
- Configure the MUSE CV Web option.
- Troubleshoot the configurations covered in the course
- Setup modem, fax, and printer devices on Server and Client
- Rebuild/upgrade the MUSE client workstation
- Review system messages and event logs
- Perform system shutdown procedures

Course Equipment

- MUSE Server
- MUSE Client Workstations
- MAC 5000/5500 with modem, wired, and wireless connections
- MARS PC server and clients
- Dash 3/4000
- Solar 8000
- Case and Case8000 Stress Systems

Related Courses

Case / Case 8000

MAC5000/5500/3500

MARS PC Server and Clients

DASH