

## Giraffe & Panda iRes Warmer Technical Education Class

### Course Summary

The Giraffe & Panda iRes Warmer technical education class is designed to teach service personnel the features of the warmers along with how to install, service, repair, and maintain them. Participants will be involved in basic system operations, identification of service tools and diagnostics, and disassembling the system. Participants will learn hands-on skills in labs designed to reinforce principles reviewed in class.

### FSE Requirements

None.

### Course Length

- 3 days

### Course Prerequisite

None

### Course Location

- Madison, WI

### Course Objectives

Upon successful completion of the course, participants will be able to:

- Provide an overview of the warmer, the various models and a description of the features.
- Describe the power distribution through the system.
- Locate, remove and replace the hardware components according to the procedures in the service manual.
- Perform the recommended maintenance on the warmer.
- Identify the variations of the Resuscitation Units and be able to service them.

### Course Equipment

- ❑ Giraffe & Panda Warmers
- ❑ Standard tool kit found in all classrooms
- ❑ 5.5 mm nut driver
- ❑ 13 mm socket/open end wrench
- ❑ Heater dish alignment tool (p/n M1118361)
- ❑ Laser holder (p/n M1118362)
- ❑ Null modem cable (p/n M1128795)



## GE Healthcare Education Services

- ❑ Test Fixture Kit (p/n M1132277)
- ❑ DISS O<sub>2</sub> Fitting (p/n 6700-0522-800)
- ❑ Certifier O<sub>2</sub> Cell (p/n 6051-0000-217)
- ❑ 18 mm wrench

### Related Courses

Giraffe & Panda Warmer Technical Overview

