



Giraffe Infant Care Technical Service Class

Course Summary

Students who participate in the Giraffe® technical class will receive instruction for both the Giraffe Omnibed® and the Giraffe Incubator. Upon completion, students will have an understanding of the differences between the two systems' functionality and hardware. Hands-on exercises will allow participants to perform expected annual preventive maintenance, adjustments, calibrations, and service checks. This course does not cover the Ohio® Infant Warmer Systems or the Care Plus® Incubators. A separate class is available for these Infant Care products.

Intended Audience

Bio-meds who will be servicing the Giraffe Infant Care Products.

Customer Requirements

Because of safety concerns, students are required to wear safety shoes with a steel toe.

Course Length

- 2 days

Course Prerequisites

- None

Course Location

- Because this class could be offered at different G.E. Training Centers, please check the schedule for location.



Course Objectives

Upon successful completion of the course, participants will be able to:

- Identify components which are common and unique to the Omni-bed and the Incubator.
- Identify and correct problems that result from operator error.
- Identify all circuit boards and list major functions of each board.
- Configure the software features of the units.
- Disassemble and assemble bed, base, and hood assemblies.
- Disassemble and assemble the Omni-bed rail assemblies.
- Perform required adjustments and calibrations to the Omni-bed.
- Calibrate Servo O2 and Humidifier.

Course Equipment

- Giraffe Omnibed
- Giraffe Incubator

Related Courses

Infant Care Class