

# S/5™ Collect

## S/5™ Collect4

S/5 Collect, a valuable research tool for professionals



S/5 Collect provides you with an easy and efficient way to collect data from Datex-Ohmeda monitors. Data can be collected from minutes to days at selectable time intervals, the only limit is your PC hard disk. With S/5 Collect you have ability to gather data straight from the monitor via an interface cable or connect directly from S/5 Central. You can add notes and markers during data collection. Analyzing the data is easy - all you need are known applications such as Excel®.

### Formats you are able to use


- .drc ( Datex Record Interface format) ; online
- ASCII; online/offline

### You can expand your applications with plug-ins such as

- LabVIEW
- C++
- DLL

### You can collect

- Waveforms
- Trends
- Alarms
- Notes
- Markers



Note: The S/5 Collect software is intended to be used as a research tool via Datex-Ohmeda monitors having no medical purpose and not affecting the intended use of the other medical devices.



**Datex-Ohmeda Division** ■ Instrumentarium Corp.  
 P.O. Box 900, FIN-00031 Datex-Ohmeda, Finland  
 Tél. +358 10 394 11 ■ Fax +358 9 146 3310

[www.datex-ohmeda.com](http://www.datex-ohmeda.com)

**Datex-Ohmeda, Inc.**  
 P.O. Box 7550, Madison, WI 53707-7550, USA  
 Tel. 800 345 2700 ■ Fax 608 221 4384

In the USA: [www.us.datex-ohmeda.com](http://www.us.datex-ohmeda.com)