

# Critikon Cuffs

## Total Solutions for Non-Invasive Blood Pressure

### Air Hose Tubing

GE Healthcare provides a comprehensive portfolio of air hoses for use with GE monitors connecting to Critikon® blood pressure cuffs. GE's air hoses are offered in a variety of sizes and styles – including Adult, Pediatric and Neonatal.

GE's air hoses are latex-free and prevent air leakage by providing a secure connection to the blood pressure cuff.



Adult-Pediatric, 12 ft, Gray Shroud to Screw Connector, Gray

Compatibility: DINAMAP® Pro, Compact, MPS Series Monitors  
**Part # 107365**



Neonatal, 12 ft, Gray Shroud to Female Slip Luer, Light Blue

Compatibility: DINAMAP® Pro, Compact, MPS Series Monitors  
**Part # 107368**



Adult-Pediatric, 12 ft, Gray Shroud to Submini Connector, Gray

Compatibility: DINAMAP® Pro, Compact, MPS Series Monitors  
**Part # 107363**



Adult-Pediatric, 24 ft, Gray Shroud to Screw Connector, Gray

Compatibility: DINAMAP® Pro, Compact, MPS Series Monitors  
**Part # 107366**



Air Hose Instrument Connector Assembly (Gray Shroud)

Compatibility: DINAMAP® Pro, Compact, MPS Series Monitors

**Part # 330073**



Adult-Pediatric, 12 ft, Rectangular Connector to Screw Connector, Gray

Compatibility: Eagle 3000, DASH 2000/3000/4000, Tram 400/450/451, Tram 800/850/851, Corometrics 120

**Part # 9461-217**



Neonatal, 8 ft, Rectangular Connector to Female Slip Luer, Light Blue

Compatibility: Eagle 3000, DASH 2000/3000/4000, Tram 400/450/451, Tram 800/850/851

**Part # 414874-001**





Adult-Pediatric, 12 ft, Screw Connector to Screw Connector, Blue

Compatibility: DINAMAP® 1846, 8100, 9300 and Plus Series Monitors

**Part # 008841**



Neonatal, 8 ft, Screw Connector to Female Slip Luer, Light Blue

Compatibility: DINAMAP® 1846, 8100, 9300 and Plus Series Monitors

**Part # 008840**



Adult-Pediatric, 12 ft, Rectangular Connector to Submin Connector, Gray

Compatibility: Eagle 3000, DASH 2000/3000/4000, Tram 400/450/451, Tram 800/850/851, Corometrics 120

**Part # 414873-001**



Adult-Pediatric, 12 ft, Round Connector to Submin Connector, Gray

Compatibility: Eagle 4000, Tram 100/200/250/300/350, Corometrics 120

**Part # 9461-203**



Adult-Pediatric, 12 ft, Screw Connector to Submin Connector, Blue

Compatibility: DINAMAP® 1846, 8100, 9300 and Plus Series Monitors

**Part # 088845**



Adult-Pediatric, 24 ft, Screw Connector to Screw Connector, Blue

Compatibility: DINAMAP® 1846, 8100, 9300 and Plus Series Monitors

**Part # 008842**



Adult-Pediatric, 12 ft, Round Connector to Screw Connector, Gray

Compatibility: Eagle 4000, Tram 100/200/250/300/350, Corometrics 118

**Part # 9461-214**



Neonatal, 8 ft, Round Connector to Female Slip Luer, Light Blue

Compatibility: Eagle 4000, Tram 100/200/250/300/350

**Part # 9461-205**



Neonatal, 12 ft, Rectangular Connector to Single Female Slip Luer, Gray

Compatibility: Eagle 3000, DASH 2000/3000/4000, Tram 400/450/451, Tram 800/850/851

**Part # 9461-216**



Neonatal, 12 ft, Rectangular Connector to Dual Female Slip Luer, Gray

Compatibility: Eagle 3000, DASH 2000/3000/4000, Tram 400/450/451, Tram 800/850/851

**Part # 9461-218**



Neonatal, 12 ft, Round Connector to Dual Female Slip Luer, Gray

Compatibility: Eagle 4000, Tram 100/200/250/300/350

**Part # 9461-215**



Adult-Pediatric, 10 ft, Round Connector to Round Connector, Black

Compatibility: Eagle 4000, Tram 100/200/250/300/350, Corometrics 118

**Part # 9461-201**

© 2005 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE and GE Monogram are trademarks of General Electric Company.

GE Medical Systems Information Technologies, Inc., a General Electric company, going to market as GE Healthcare.



imagination at work