

<b>100% Latex Free</b>	<b>Benefits of Latex Free</b>
<p>All Critikon® blood pressure cuffs, accessories, adaptors, connectors, bulbs, valves and air hoses are 100% Latex Free (as indicated below). Our 100% Latex Free product line allows us to assist your facility in limiting risk, controlling exposure and limiting latex sensitivity for both your patients and your staff.</p>	<p>Latex sensitivity affects over 6% of the population and can cause death. People with extended contact, like healthcare workers and patients who undergo multiple procedures are at an increased risk. Controlling exposure, ensure less risk, less paperwork and tracking and no need to maintain separate stock.</p>
<b>Latex Content By Product Grouping</b>	
<b>Dura-Cuf® Blood Pressure Cuffs</b>	<b>Soft-Cuf Blood® Pressure Cuffs</b>
<p>All Dura-Cuf® Blood Pressure Cuffs manufactured after march 1997 are latex free. The lot code is printed directly on the cuff and is read month over year (e.g. 3/7 is March 1997). To verify that the product in question does not contain latex, please look for the "latex free" symbol printed on the blood pressure cuff and "Contains No Latex" on the package. All packaging, bags, labels and boxes are latex free.</p>	<p>All Soft-Cuf® Blood Pressure Cuffs (which includes Soft-Cuf® Neonatal Blood Pressure Cuffs but excludes Complete Inflation Systems) have been latex free for over seven years. To verify that the product in question does not contain latex, please look for the "latex free" symbol printed on the blood pressure cuff and "Contains No Latex" on the package. All packaging, bags, labels and boxes are latex free.</p>
<b>Classic-Cuf® Blood Pressure Cuffs</b>	<b>Critikon® Accessories</b>
<p>All Classic-Cuf® Blood Pressure Cuffs (which includes Classic-Cuf® Neonatal Blood Pressure Cuffs but excludes Complete Inflation Systems) have been latex free for over seven years. To verify that the product in question does not contain latex, please look for the "latex free" symbol printed on the blood pressure cuff and "Contains No Latex" on the package. All packaging, bags, labels and boxes are latex free.</p>	<p>All Critikon® blood pressure cuff accessories (including connectors, adaptors, bulbs, valves and air hoses) manufactured after January 1998 are latex free. To verify that the product in question does not contain latex, please look for the "Contains No Latex" on the package. If in doubt, please discard the accessory in question and order a replacement. All packaging, bags, labels and boxes are latex free.</p>
<b>Complete Inflation Systems</b>	
<p>All Critikon® Complete Inflation Systems (including Dura-Cuf®, Soft-Cuf® and Classic-Cuf® Complete Inflation Systems) manufactured after January 1998 are latex free. To verify that the product in question does not contain latex, please look for the "latex free" symbol printed on the blood pressure cuff and "Contains No Latex" on the package. All packaging, bags, labels and boxes are latex free.</p>	<p>Products may come in contact with latex/natural rubber through incidental contact in the manufacturing process. Examples of this incidental contact include: protective clothing, gloved hands, gaskets, drive belts, adhesives, etc. For more information on Latex Sensitivity or Latex Content, please contact GE Healthcare Sales Representative.</p>