



ACCESSORIES CATALOG

A complete line of accessories specially designed for use with your Giraffe® OmniBed® and Giraffe® Incubator.






Giraffe Spot PT Lite™ Phototherapy System

- Focused light provides high intensity phototherapy (30 – 40 $\mu\text{W}/\text{cm}_2/\text{nm}$)
- Cooler spot light minimizes effect on incubator or baby's temperature
- Long-lasting metal halide bulb provides up to 2500 hours life
- Mounts to any Ohmeda Medical unit equipped with dovetail rail system, including Giraffe, Care Plus, Panda Baby Warmer or Ohio Infant Warmer System



Giraffe Monitor Shelf

- Mounts to dovetail rail behind control panel to provide space for monitors and other accessories
- Holds up to 50 lbs (22 kg)
- Size: 12 in x 26 in (30.5 cm x 66 cm)



Giraffe Exam Light

- Two brightness switch for high or low light intensity
- May be plugged directly into Giraffe or wall outlet
- Adjustable gooseneck positions light anywhere you need it
- Mounts to any Ohmeda Medical unit equipped with dovetail rail system, including Giraffe, Care Plus, Panda Baby Warmer or Ohio Infant Warmer System



Rotating Dovetail Rail Extension

- Increases the amount of available mounting space
- Ideal for suction and oxygen therapy products and other lighter equipment
- Holds up to 11 lbs (5 kg)



Giraffe Instrument Shelf

- Swivels for ease of positioning and mounts to dovetail rail for smaller accessories
- Holds up to 20 lbs (9 kg)
- Size: 12 in x 12 in (30.5 cm x 30.5 cm)



Rotating IV Pole (Adjustable Height)

- May be loosened and rotated behind the bed allowing the bed to fit through the door
- Mounts to any Ohmeda Medical unit equipped with dovetail rail system, including Giraffe, Care Plus, Panda Baby Warmer or Ohio Infant Warmer System
- Holds up to 20 lbs (9 kg)
- Min. length: 37 in (94 cm); Max. length: 64 in (163 cm)



**Easy-Load
E-Cylinder
Holder**

- Allows cylinders to be replaced easily without lifting the length of the tank
- Mounts to any Ohmeda Medical unit equipped with dovetail rail system, including Giraffe, Care Plus, Panda Baby Warmer or Ohio Infant Warmer System



Corner Trays

- Fit securely into mattress translation deck for storing small supplies
- Two trays included per pack




**Tubing
Management
Arm**

- Secures tubes and lines to reduce stress on infant
- Accepts multiple size vent tubes and lines
- Can be placed in any corner of mattress deck



Hood Cover

- Custom designed cover protects infant from light and sound
- Side flaps allow easy access to doors without removing cover
- Top flap provides opening for phototherapy
- Flame retardant and machine washable



**Silo Support
Arm**

- Designed to suspend gastroschisis silos or other care items over a baby



**High
Frequency
Oscillatory
Vent (HFOV)
Porthole
Cover**

- Helps maintain temperature, humidity, and O2 levels during HFOV use



**Fitted
Mattress
Sheet**

- Fitted sheet designed to conform to rounded mattress
- Flame retardant and machine washable

ORDERING INFORMATION

Details for **highlighted** products are shown on previous pages.

Giraffe® Family Accessories

(accessories can be used on any Giraffe® family product)

Giraffe Monitor Shelf	6600-0824-800
Giraffe Instrument Shelf	6600-0865-700
Tubing Management Arm	6600-0837-800
High Frequency Vent Porthole Cover	6600-0838-800
Giraffe Easy-Load E-cylinder Holder	6600-0836-800
Foam Mattress	6600-1116-500
Porthole wristlets (Box of 8)	6600-0211-850
Filter (package of 10)	6600-0207-850
Humidification Reservoir	6600-0216-850
In-Bed Scale	6600-0501-900
Silo Support	6600-0853-800
Fitted Sheet	6600-0850-800
Giraffe Disposable Patient Probe (Box of 10)	6600-0873-700
Giraffe Disposable Patient Probe (Box of 50)	6600-0874-700
Giraffe Reusable Patient Probe	6600-0875-700
Heat Reflecting Patch (Box of 50)	0203-1980-300
IRIS Cuff (Box of 10)	6600-0499-800
IRIS Cuff (Box of 100)	6600-0500-800
Corner Tray - SE and NW	6600-1793-500
Corner Tray - SW and NE	6600-1794-500
DIN Rail	6600-0855-800

Giraffe® OmniBed®-Specific Accessories

Servo Oxygen Upgrade Kit with DISS Inlet	6600-0678-800
Servo Oxygen Upgrade Kit with NIST Inlet	6600-0678-801
Servo Oxygen Upgrade Kit with DIN Inlet	6600-0678-802
Servo Oxygen Upgrade Kit w/ Air Liquide Inlet	6600-0678-803
OmniBed Hood Cover	6600-0825-800
Service Manual, North America	6600-0834-801
Service Manual, International A4	6600-0834-802

Giraffe® Incubator-Specific Accessories

Servo Oxygen Upgrade Kit with DISS Inlet	6600-0849-800
Servo Oxygen Upgrade Kit with NIST Inlet	6600-0849-801
Servo Oxygen Upgrade Kit with DIN Inlet	6600-0849-802
Servo Oxygen Upgrade Kit w/ Air Liquide Inlet	6600-0849-803
Incubator Hood Cover	6600-0846-800
Service Manual, North America	6600-0847-801
Service Manual, International A4	6600-0847-802

Giraffe® Spot PT Lite™

English - Giraffe plug (IEC plug)	6600-0833-800
French - Giraffe plug (IEC plug)	6600-0833-801
Spanish - Giraffe plug (IEC plug)	6600-0833-802
English - US plug	6600-0845-800
French - US plug	6600-0845-801
Spanish - US plug	6600-0845-802
German - CE plug	6600-0845-803
Italian - Italian plug	6600-0845-804
Dutch - CE plug	6600-0845-809
English - UK plug	6600-0845-811
Replacement Metal Halide Bulb (1 qty)	6600-0235-850
Portable Roll Stand	6600-0894-216

Giraffe® Exam Light

120V - North American Plug	6600-0832-800
120V - IEC Plug	6600-0832-801
230V - UK Plug	6600-0832-802
230V - IEC Plug	6600-0832-803
230V - CE Plug	6600-0832-804
100V - Japanese Plug	6600-0832-805
100V - IEC Plug	6600-0832-806
Replacement Bulb	6600-0912-200

General Accessories

(accessories can be used on any Ohmeda Medical product with dovetail rail system)

Dovetail Rail Extension	6600-0852-800
Rotating IV Pole	6600-0851-800
IV Pole Dual Hook - 12"	0217-5378-800
IV Pole Dual Hook - 24"	6600-0491-801
IV Pump Mounting Post (20" x 1")	0217-5376-800
Ventilator Mounting Pole	0217-5357-800
Utility Post (3.5" x 1")	0217-5374-800
Universal Manifold Bracket	0217-5359-800
Vacuum Slide Bracket	0217-5367-800
Retaining Clips (6)	6600-0055-851
Universal Rail Adapter Plate	0217-5363-800
Standard Gas Cylinder Holder	6600-0422-800



Making life a little easier

Ohmeda Medical
8880 Gorman Road
Laurel MD 20723 USA
www.ohmedamedical.com

Ohmeda Medical Customer Service
Tel 800 345 2700 in the US
Tel 410 888 5200 outside the US

Latin America, Caribbean
Ohmeda Medical
8880 Gorman Road
Laurel MD 20723 USA
Tel 410 888 5220
Fax 301 483 8340
icss@ohmedamedical.com

Europe, Africa, Middle East
Ohmeda Medical
8880 Gorman Road
Laurel MD 20723 USA
Tel 410 888 5200
Fax 301 483 8340
icss@ohmedamedical.com

Asia Pacific
Datex-Ohmeda Pte Ltd
152 Beach Road
No 12-05/07 Gateway East
Singapore 189721
Tel 65 6391 8618
Fax 65 6291 6618

Canada
Dynamed Health Care Systems
235 Shields Court
Markham, Ontario L3R 8V2
Tel 905 752 3300
Fax 905 752 3304