

BIOPSY GUIDE



gemedicalsystems.com

LOGIQ 700 EXPERT/PRO Series



Fixed angle biopsy device compatible with the 618E endocavity probe.

LOGIQ 500 PRO Series



Fixed angle biopsy device compatible with the E721 endocavitary probe.

LOGIQ 400 PRO Series



Mechanical biopsy device and disposable needle receptacles compatible with convex, sector, and linear array probes.

	Probe	Biopsy Starter Kit Catalog Number	Biopsy Guide Type
LOGIQ 700 EXPERT/PRO Series	Convex 348C	E8385MG*	Single Angle, Reusable
	Convex 548C	E8385MG*	Single Angle, Reusable
	Microconvex 618C	E8385NJ*	Single Angle, Reusable
	Microconvex 618E (TV/TR)	E8386CF	Single Angle, Disposable
	Linear 546L	E8385ML*	Single Angle, Reusable
	Linear 739L	E8385MC*	Single Angle, Reusable
	Linear LA39	E8385MM*	Single Angle, Reusable
LOGIQ 400/500 PRO Series	Convex C358	E8386RK	Multi Angle, Reusable
	Convex C551	E8385R** / E8385MF*	Multi/Single Angle, Reusable
	Convex C364	E8385PA** / E8385MD*	Multi/Single Angle, Reusable
	Microconvex E721 (TV/TR)	E8385MJ / H46222AD	0°/Single Angle, Disposable/Reusable
	Microconvex E721 (TV/TR)	H4550BG	5° Single Angle, Disposable
	Microconvex C721	E8385RB**	Multi Angle, Reusable
	Sector S317	E8385PD** / E8385SA*	Multi/Single Angle, Reusable
	Sector S611	E8386RH**	Multi Angle, Reusable
	Linear 546L	E8385PC** / E8385ML*	Multi/Single Angle, Reusable
	Linear 739L	E8386RJ** / E8385MC*	Multi/Single Angle, Reusable
	Linear LA39	E8385MM*	Single Angle, Reusable
	Linear L764	E8385RA** / E8385LA*	Multi/Single Angle, Reusable
	T-Type T739	E8386RG**	Multi Angle, Reusable

*E8385LC – Additional 24 Disposable Procedure Kits for Single Angle Biopsy

**E8385RC – Additional 24 Disposable Procedure Kits for Multi Angle Biopsy

E8386CB – Additional 25 needle guides for 618E Probe

E8385RD – Additional 24 Disposable Procedure Kits for C358 Probe

Starter Kit includes one reusable bracket (except for disposable type) and five Disposable Procedure Kits.

Disposable Procedure Kit includes needle guides, probe covers, sterile gel, and bands.