

TEE Ultrasound Probe Accessories

Protect your investment

You spend a lot of time and money investing in your GE Healthcare ultrasound equipment. Take the extra step and protect your investment – and reduce your repair cost – by using our cardiac accessories and supplies. Most TEE failures can be prevented if recommended accessories are used.

Here are a few things you can do to maximize your uptime:

- Conduct weekly probe inspections.
- Inspect probe cables and scan heads for nicks and cuts. If you find damage, report it quickly.
- Use OEM-recommended prepping and cleaning materials.
- Use OEM-recommended mouth protectors and storage devices to increase probe reliability.

Less downtime. More cost effective.



Adult Reusable Mouth Protector
Standard TEE probe for adults.

P/N: KX307217



Pediatric Mouth Protector
Perfect for protecting TEE probes when used on children or on small patients.

P/N: KX307216



Clip-on Mouth Protector

Clip on/off feature for rapid mounting and removal during the intubation and extubation process.

P/N: KX307214



**Operating Room
Clip-on Mouth Protector**

Designed for operating room. Use with clip on/off feature for rapid mounting and removal during the intubation and extubation process.

P/N: KX307215



Foam Scan Head Protector

Slips on over scan head. Disposable, inexpensive protection for delicate scan heads.

P/N: KX307218 – Adult

P/N: KX307219 – Pediatric



Leak Path Indicator Tester

Quick and easy-to-use tester identifies penetrating holes, tears or fractures while promoting patient safety.

P/N: KZ200800



Simplify Your Ordering

Now our ordering process has become even easier. **Simply make copies of the order form below** and use them as needed. Fax in your order or call the 800 number to get all of your GE ultrasound accessories and supplies.

Order form

Description	Part Number	Number of Packages	Price per Pkg.**	Ext. Price	Date
Adult <i>Reusable</i> Mouth Protector Package*	KX307217		\$94.95		
Mouth Protector – Pediatric Package*	KX307216		\$57.95		
Clip-on Mouth Protector Package*	KX307214		\$75.95		
OR Clip-on Mouth Protector Package*	KX307215		\$75.95		
Foam Scan Head Protection Covers – Adult Package*	KX307218		\$39.99		
Foam Scan Head Protection Covers – Pediatric Package*	KX307219		\$39.99		
Leak Path Indicator Tester	KZ200800		\$766.15		

* 1 package equals 25 parts

**Prices subject to change

Shipping Information

Facility Name: _____

Customer Name: _____

Address: _____

City, State, Zip: _____

Phone number: _____

E-mail: _____

To place an order please call 800-558-2040, press option 2, or fax 262-548-2011, 7a.m. – 6 p.m., CST.

GE Healthcare
9900 Innovation Drive
Wauwatosa, WI 53226
U.S.A.

www.gehealthcare.com

©2008 General Electric Company – All rights reserved.
GE and GE Monogram are trademarks of General Electric Company.

General Electric Company, doing business as GE Healthcare.



imagination at work