

GE Healthcare



Silhouette VR

Analog X-ray system





Value. Simplicity.

Silhouette VR.

The Silhouette VR radiographic system from GE Healthcare delivers maximum performance and productivity in a compact, easy-to-use package.

This fully integrated system offers intuitive controls with advanced applications for maximum performance – all in a space-saving design backed by the best service in the industry.

Designed for ease. Designed for you.

GE designed the Silhouette VR based on technologists' input to make it better than ever. This sleek system is easy to use, which helps boost productivity and enhance throughput in your facility. With an automatic bucky-tray coupling between on- and off-table applications, it helps minimize operator intervention to smooth workflow.

Designed with a user-friendly console, the Silhouette VR improves throughput over a full range of on- and off-table applications, while the pedestal-style table with a four-way floating top speeds positioning.

Functional and flexible in any space.

The freestanding Silhouette VR is compact to save space in your facility. With a three-phase generator neatly tucked under the table, it offers a small footprint to fit in any exam room. A fully counterbalanced tube stand rotates 180 degrees for superior flexibility and head-to-toe coverage for most patients. It also allows anatomical programming (APR) by morphology, view and body region.

Easy front-cover access to the generator speeds service calls, and the state-of-the-art generator is both microprocessor-controlled and self-calibrating to reduce maintenance.

Standard configuration:

- Three-phase, high-frequency 32 kW generator (optional 50 kW)
- Anatomical Programming (APR)
- Dual, 0.6 mm and 1.2 mm focal spots
- Tube stand integrated with table
- Manual collimator
- Four-way floating table top

Options and accessories:

- Automatic Exposure Control (AEC)
- Wall stand
- 100 cm grid for wall stand
- Hand grips
- Lateral cassette holder
- Patient support band

Get the details today

For details and complete specifications, please contact your GE Representative today.



Space-saving generator integrated into table.



Silhouette VR table and tube stand.



Rotation of tube stand.

©2007 General Electric Company.
GE and GE Monogram are trademarks of General Electric Company.

General Electric Company, doing business as GE Healthcare.

Healthcare Re-imagined

GE is dedicated to helping you transform healthcare delivery by driving critical breakthroughs in biology and technology. Our expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, and biopharmaceutical manufacturing technologies is enabling healthcare professionals around the world to discover new ways to predict, diagnose and treat disease earlier. We call this model of care "Early Health." The goal: to help clinicians detect disease earlier, access more information and intervene earlier with more targeted treatments, so they can help their patients live their lives to the fullest. Re-think, Re-discover, Re-invent, Re-imagine.

GE Healthcare
3000 North Grandview
Waukesha, WI 53188
U.S.A.

www.gehealthcare.com



imagination at work