




Improving your image.

ILUMA™ Ultra Cone Beam CT scanner



Authorized Distributor
GE Healthcare

Striving to help you better understand disease, diagnose and offer treatment options for your patients.


Opportunity knocks.

That's why we are proud to offer the ILUMA™. An Ultra Cone Beam CT scanner manufactured by IMTEC Imaging, Inc. Giving you the opportunity to improve your practice, and enhance your patient care experience by utilizing high-resolution diagnostic CT imaging from the convenience of your own office.


- Complete in-office patient care
- Advanced visualization software
- Highly detailed head and neck imaging
- Low dose

ILUMA







Head and neck




Inner ear



Multiplanar imaging




Virtual fly-through of sinuses



Temporal mandibular joint

Sharp edges.
Superb image quality.
Large field of view.
The results?



Surgical navigation

Clarity.

ILUMA uses proprietary Flash CT technology, the latest generation of advancement in ultra cone beam volumetric tomography to produce state-of-the-art sinus, skull base, head and neck, and temporal bone images in a simple 20 second scan.

- Images as thin as 0.09 mm isotropic reconstructed voxels. For typical ENT and surgical navigation applications, 0.4 mm image reconstructions are standard.
- Reconstruction time is two-and-a-half minutes at 0.4 mm voxels.
- Reconstructed 14-bit (16,384 levels of grayscale) provides optimal resolution.
- Focal spot of 0.3 mm enhances image quality.



Continuity. Convenience. Success.

Intuitive software design allows for easy set-up and completion of exams. With ILUMA you can scan a patient, instantly review the results and immediately discuss treatment options in a single visit —saving time and anxiety.

Capitalize on the additional revenue ILUMA can bring to your practice. When your patient needs a head or neck scan, ILUMA can allow your office the revenue opportunity of doing the images rather than sending your patient elsewhere.

An ILUMA CT scanner allows you a highly marketable advantage within your practice. You are not

only offering your patients a time savings by allowing for in-office imaging, but also superb image quality resulting in the diagnosis they need in one convenient office visit.

The results? Better utilization of time for both you and your patient. Increased patient throughput and plenty of opportunities for return on your CT scanner investment.



©2007 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, and GE Monogram are trademarks of General Electric Company.

ILUMA is a trademark of IMTEC Corporation.

GE OEC Medical Systems, Inc., a General Electric company, going to market as GE Healthcare.

Distributed by:

GE Healthcare
384 Wright Brothers Drive
Salt Lake City, Utah 84116
USA

**GE Healthcare,
Surgery – Americas:**
Phone 801-328-9300
Fax 801-328-4300

GE Healthcare – Europe:
Paris, France
Fax 33-1-30-70-94-35

GE Healthcare – Asia:
Tokyo, Japan –
Fax: +81-452-85-5490
Hong Kong –
Fax: +852-2559-3588



Authorized Distributor
GE Healthcare