



SENSA-CUF

Blood Pressure Cuffs

The CRITIKON™ SENSA-CUF™ is a highly versatile solution for a wide range of extended-care settings. Made using strong, breathable, easy-to-clean material, the SENSA-CUF is designed to maximize patient comfort in low-acuity areas in a hospital where risk of contamination are diminished, and durable enough for multi-patient use.



A trusted choice for clinical excellence

- Cuffs are held to the same quality standard as all GE patient monitors, with repeated inflations to ensure accurate readings
- Large, bold print helps you quickly choose the right cuff, enabling a more accurate measurement

Constructed for patient safety and durability

- Soft, breathable, easy-to-clean material is not made with PVC, making them safer for the environment, or natural rubber latex, making them safer for your patients
- DINACLICK™ connector meets international standards and prevents misconnections to luer fittings

Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the monitor faster than other dual-tube connectors
- Cuffs compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with automated monitors and a variety of manual devices, including wall manometers, enabling a universal connection system throughout the facility

Assortment packs available

- For low-volume cuff utilizers, several types of assortment packs are available for easier ordering and more optimal inventory management



SENSA-CUF			ISO 80369-5 compliant	ISO 80369-1 compliant	Does not conform to ISO 80369 standard	
						
Size	Color	Limb Circumference	1-Tube Bayonet Connector	2-Tube DINACLICK Connector	2-Tube Screw Connector	2-Tube Submin Connector
Infant	Orange	8-13 cm	2470	SEN-P1-2A	2458	2482
Child	Green	12-19 cm	2472	SEN-P2-2A	2460	2484
Small Adult	Light Blue	17-25 cm	2474	SEN-A1-2A	2462	2486
Small Adult Long	Light Blue	17-25 cm	2475	SEN-A1-2A-L	2463	2487
Adult	Navy	23-33 cm	2476	SEN-A2-2A	2464	2488
Adult Long	Navy	23-33 cm	2477	SEN-A2-2A-L	2465	2489
Large Adult	Rose	31-40 cm	2478	SEN-A3-2A	2466	2490
Large Adult Long	Rose	31-40 cm	2479	SEN-A3-2A-L	2467	2491
Thigh	Brown	38-50 cm	2480	SEN-T1-2A	2468	2492

Note: 5 cuffs/box

DINACLICK assortment packs

Connector Type	Part Number	Description
2-Tube DINACLICK Connectors	2059300-001	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/box)

Other connector assortment packs

Connector Type	Part Number	Description
2-Tube Screw Connectors	2494	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/box)

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com



© 2020 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, CRITIKON, SENSА-CUF and DINACLICK are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB00417XX 8/20