

Vscan Extend™ with Sector Probe

Indications Reference Guide

Deep scanning transducer

Vscan Extend with sector probe is indicated for ultrasound imaging, measurement, and analysis of the human body in clinical applications that include...

CLEARED APPLICATION	ANATOMY	EVALUATION
Abdominal ultrasound (adult & pediatric)	<ul style="list-style-type: none"> • Aorta • Gall bladder • Kidney • Liver • Spleen 	<ul style="list-style-type: none"> • Aortic aneurysm • Gall stone • Kidney stones or fluid • Hepatomegaly • Splenomegaly • Peritoneal fluid
Urology ultrasound (adult & pediatric)	<ul style="list-style-type: none"> • Urinary bladder • Kidney/ureter 	<ul style="list-style-type: none"> • Bladder dysfunction • Post-renal obstruction • Pre-post bladder volume
Cardiac ultrasound (adult & pediatric)	<ul style="list-style-type: none"> • Left ventricular function • Myocardium septum wall • Inferior vena cava • Mitral valve function • Aortic valve function • Pericardium 	<ul style="list-style-type: none"> • Left ventricular dysfunction • Left ventricular hypertrophy • Dilation • Valve open & closure • Pericardial fluid
Fetal/obstetrics & gynecology	<ul style="list-style-type: none"> • Fetal position • Fetal heart beat • Amniotic Fluid • Placenta • Uterus 	<ul style="list-style-type: none"> • Breech • Visualize heart motion • Fluid levels • Position • Fibroids
Thoracic/lung ultrasound (adult & pediatric)	<ul style="list-style-type: none"> • A, B, E lines • Lung sliding • Lung point 	<ul style="list-style-type: none"> • Pleural effusion • Acute lung injury including pneumo- or hemothorax • Alveolar interstitial • Lung consolidation <ul style="list-style-type: none"> – Pneumonia or pneumonitis – Pulmonary fibrosis • Acute respiratory distress syndrome
Cephalic	<ul style="list-style-type: none"> • Cerebral vessels 	<ul style="list-style-type: none"> • Flow
Procedure guidance	<ul style="list-style-type: none"> • Heart • Lung • Abdomen • Uterus • Thorax 	<ul style="list-style-type: none"> • Fluid detection including: pericardial, pleural, peritoneal, amniotic • Procedures: thoracentesis, pericardiocentesis-, paracentesis, and amniocentesis

Products mentioned in the material may be subject to government regulation and may not be available in all countries. Shipment and the effective sale can only occur after approval from the regulator. Please check with your local GE Healthcare representative for details.

Data subject to change.

© 2021 General Electric Company.

GE, the GE Monogram, and Vscan Extend are trademarks of General Electric Company.

