

CARESCAPE Monitor B850

Engineered to help provide better care

The CARESCAPE™ Monitor B850 is a high-acuity monitor that can help you manage your patient flow by providing a dependable level of data continuity and integration across care areas. Its intuitive design allows for easy training for all staff experience levels, and upgradability protects your long-term system investments.

Providing clinical excellence, dependable data continuity and integration

- iPanel™ integration enables easy access to labs, X-rays, charting and other data at the bedside
- CARESCAPE Patient Data Module provides consistent hemodynamic measurement during intra-hospital transport and also transfers trend data
- Innovative algorithms aid in accurate diagnoses, including GE EK-Pro, GE DINAMAP® non-invasive blood pressure, and 12 SLTM diagnostic ECG with direct 2-way MUSE® communication
- Specifically designed airway modules provide for complete respiratory monitoring with gas exchange measurement
- Various parameters can help you assess the adequacy of anesthesia in the OR, and give insight to patients' readiness for ventilator weaning in the ICU
- Auto View on Alarm (AVOA) automatically shares high priority alarms within the care unit
- High-definition visualization of real-time and historic patient data, for beat-to-beat and breath-to-breath trending
- Integration with the Aware® Gateway enables communications to EMR systems through standard HL7 protocol

Intuitive design

- Pages & Profiles functionality can increase flexibility and workflow efficiency by configuring monitors to unit standards and patient populations
- Integrated alarm light on top of the display provides a visual indication of alarm status, enabling viewing from a greater distance
- User-controlled views, from eight to 14 waveforms and up to 22 parameter boxes, allow for flexibility and customization, based on caregivers' and patients' needs
- Integrated USB ports allow for a keyboard, mouse, barcode readers and other data entry device attachment
- Trim Knob® control and optional hard keys provide one-touch activation of key features
- Dedicated software and parameters for the OR, PACU, CCU, NICU and ED care areas make it easy to configure to your work patterns

Protecting long-term investments

- Combines the best of Marquette® and Datex-Ohmeda® legacies into one platform
- Backwards compatible with many existing Marquette and Datex-Ohmeda components
- Remote service through InSite ExC provides updates and predictive maintenance
- Structured upgrade programs can lower costs of ownership and simplify long-term capital equipment planning



Technical specifications

CPU

Dimensions (H x W x D) 9.1 x 40.1 x 34 cm
(3.6 x 15.8 x 13.4 in)

Weight 7.5 kg (16.5 lb)

Power consumption 300 Watts (max)

Power requirements 100 to 240 VAC, 50-60 Hz

Heat dissipation 1024 BTU/hour (max)

Intel 423 embedded

Degree of protection against harmful ingress of water

IPX1

Alarms

Processor

Priority 3 levels - High, Medium and Low

High-priority alarm volume in accordance with IEC 60601-1-8 D15K* 81 dB(A)

D15K* 81 dB(A) D19KT* 82 dB(A)

Visual alarm notification Red, yellow, cyan

Audio silence indicator General alarm indicator

Audible alarm notification Internal speakers

Audio pause Configurable for 2 or 5 minutes

I/O connectors

Ethernet 2 RJ45 for CARESCAPE Network MC and CARESCAPE Network IX

and Carescape Network IX

E-Port/Tram-Net® 2 DB9F for CARESCAPE Patient Data

Module, Frame, or Tram-Rac®

Serial port 2 DB9M

Video 1 DVI-I, 1 DVI-D

Optional video 1 DVI-I with optional third video card

USB port 4 USB 2.0 standard

M-port 4 RJ45

Display characteristics

| Colors | 16 million |
|--------------------------------|---|
| Type | Active matrix color TFT LCD |
| Brightness D15K D19KT | 210 cd/m ² 180 cd/m ² |
| Resolution | 1024 x 768 @ 60 Hz (XGA) |
| Viewing angle D15K D19KT | (horizontal/vertical) 170°/170° 178°/178° |

^{*} D15K = 15-inch display with keypad;

D19KT = 19-inch display with keypad and touchscreen

| - : | 1.0 | |
|------------|------|-------|
| PIXE | dime | nsion |

| D15K | $0.297 \times 0.297 \text{ mm}$ |
|-------|---------------------------------|
| D19KT | 0.294 x 0.294 mm |

Contrast ratio

D15K 700:1 D19KT 1300:1

Response time

D15K 25 ms D19KT 20 ms

Waveforms

Based on configuration. Standard view includes up to 8 individual waveforms, optional views include overlays or insets for up to 14 waveforms per display. Additional clinical display allows for 16-28 waveforms.

Numerics

16 to 22 parameter numerics based on configuration, per display. Additional clinical display allowed for 32 to 44 parameter numerics.

See separate technical specifications document for displays. Displays are ordered separately.

Frames

Target use and support modules

| F7 | Carestation® Aisys®, Avance®, Aestiva®, Aespire®, ADU anesthesia machines, Patient Side Module, E-modules |
|--------------------|--|
| F5 | Stand-alone use; E-modules, Patient Side Module, CARESCAPE Patient Data Module |
| Power requirements | 16.7 V @ 4.6 A supplied by CPU-C1 Host through ePort interface cable |
| Power consumption | < 3 Watts idle (no modules) |
| F7 | 65 Watts (max) |
| F5 | 77 Watts (max) |
| | |

Frames are ordered separately.

Parameters and modules

| Parameters | Patient Side Module | CARESCAPE Patient Data Module /Tram® |
|-------------------|------------------------|--|
| ECG | 3, 5, 6, 12-lead ECG | 3, 5, 6**, 12-lead ECG |
| SpO ₂ | GE SpO ₂ | Masimo SET®, Nellcor OxiMax® |
| NIBP | GE | GE DINAMAP® SuperSTAT™ algorithm** |
| InvBP | 0 or 2 | 0 or 4 |
| Temp | 2 | 2, Optional with C.O. |
| Cardiac output | - | Optional with temperature |

^{**}CARESCAPE Patient Data Module only

| Parameters | E-Modules‡ | Tram and Tram-Rac® Modules |
|---|---|--|
| Multi parameter | modules | |
| | E-PRESTN E-RESTN E-PRETN | Tram models: 451, 451M & 451N5 851, 851M & 851N5 |
| InvBP & Temp | E-P [§] E-PP [§] E-PT | BP/ Dual Temp Module [†] BP Module [§] Dual BP Module [§] |
| SvO ₂ & C.O. | E-COP E-COPSv | - - |
| Single paramete | r modules | |
| SpO ₂ | E-NSATX E-MASIMO | SpO ₂ Module Masimo SpO ₂ |
| NMT | E-NMT | - |
| EEG | E-EEG | - |
| BIS | E-BIS | - |
| Entropy™ | E-ENTROPY | - |
| Gas | | |
| ms CO ₂ | - | Capnostat Mainstream CO ₂ Module |
| ss CO ₂ | E-miniC | - |
| ms/ss CO ₂ | - | CapnoFlex LF Module |
| ss CO ₂ , O ₂ & N ₂ O | E-CO E-COVX E-CAIO E-CAIOV E-CAIOVX | - |
| Agents | E-CAiO E-CAiOV E-CAiOVX | - |
| Patient Spirometry | E-COV E-COVX E-CAiOV E-CAiOVX | - |
| Gas exchange/ metabolics | E-COVX E-CAiOVX | - |

Parameter modules are ordered separately.

NOTE: For a complete list of compatible devices, including M-series modules, please refer to the CARESCAPE Monitor B850 Addendum for Device Compatibility

- § Module measures invasive blood pressure only
- † Module compatible for blood pressure only, not temperature measurement
- ‡ Acronyms for the parameters measured are as follows: P=invasive pressure; R=respiration rate, E=ECG; S=SpO₂; T=temperature; N=NIBP; C=CO₂ and N₂O; Ai=anesthetic agents and nitrous oxide with agent identification; O=O₂; V=Spirometry; X=gas exchange

Software options

| Software options | |
|-------------------------|---|
| Main software | Care-area specific software versions to optimize workflows: |
| | OR, PACU, CCU, ED, NICU |
| Other software options | Extended software options specific to each main software package. |
| Updgradeable with the 0 | CARESCAPE Life Upgrade Program |
| Networking | |
| Compatibility | CARESCAPE Network |
| Features | Peer-to-peer communication with AVOA functionality, centralized viewing and remote alarm management |
| Network type | LAN |
| Interfaces | For a complete list of compatible networked devices, please refer to the CARESCAPE Monitor B850 Addendum for Device Compatitiblity |
| Mounting | |
| GCX compatible | |
| Paper recorder | |
| 0 0001 50 14 1: :: | |

Optional, PRN 50-M digital writer or network laser printer. See separate technical specifications for PRN 50-M digital writer. Paper recorders are ordered separately.

Environmental specifications

| Operating temperature | 10 to 35°C (50 to 95°F) |
|---|--|
| Storage/transport temperature | -20 to 60°C (-4 to 140°F) |
| Operating humidity | 15 to 90% (non-condensing) |
| Storage humidity | 10 to 90% (non-condensing) |
| Operating altitude | 700 to 1060 hPa; 3011 to -382 meters |
| Storage/transport altitude | 500 hPa to 1060 hPa; 5,572 to -382 meters |
| Warranty | |
| One uear. Extended service contracts available. | |

Certifications

IEC 60601-1

CAN/CSA C22.2 No. 601.1-M90

UL 60601-1

IEC 60601-1-2

CE marking according to the Medical Devices Directive 93/42/EEC

For accessory information, please refer to CARESCAPE Monitor B850 Supplies and Accessories document.

For upgrade information, please refer to the CARESCAPE Life Upgrade Program.

© 2009 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, GE Monogram and CARESCAPE are trademarks of General Electric Company.

12SL, Aware, DINAMAP, Marquette, MUSE, Tram, Tram-Net, Tram-Rac and Trim Knob are trademarks of GE Medical Sustems Information Technologies, Inc.

Aespire, Aestiva, Aisys, Avance and Carestation are trademarks of Datex-Ohmeda.

Entropy and iPanel are trademarks of GE Healthcare Finland Ou.

BIS is a trademark of Aspect Medical Systems, Inc. Capnostat is a trademark of Respironics Novametrix, Inc. Masimo SET is a registered trademark of Masimo Corporation Nellcor and OxiMax are trademarks of Tyco Healthcare Group LP or an affiliate.

GE Medical Systems Information Technologies, Inc. a General Electric Company, doing business as GE Healthcare.

GE Healthcare Finland Oy, a General Electric Company, doing business as GE Healthcare.

About GE Healthcare

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement and performance solutions services help our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems.

Our "healthymagination" vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access and improving quality and efficiency around the world.

GE Healthcare 8200 W. Tower Ave. Milwaukee, WI 53223 USA GE Healthcare Finland Oy Kuortaneenkatu 2 00510 Helsinki Finland

GE Healthcare 3/F Building # 1, GE Technology Park 1 Hua Tuo Road Shanghai 201203 China

www.gehealthcare.com

