



# CRITIKON SOFT-CUF

## Blood Pressure Cuffs

Our most popular cuff, the CRITIKON SOFT-CUF™ is designed with soft, absorbent material with rounded corners for maximizing patient comfort during extended stays at the hospital. These cuffs are ideal for patients with longer term stays or for protocols where the cuff travels with the patient through their total length of time in the hospital. It is economically priced to help support your infection control efforts.



### Designed for patient safety, comfort and accuracy

- Priced economically for pairing with a single patient throughout their hospital stay to help with infection control efforts
- Soft, absorbent material is not made with DEHP or natural rubber latex to help enhance patient safety
- Designed with rounded corners for patient comfort
- Designed and clinically validated for performance accuracy to help you make well-informed decisions
- DINACLICK™ connector meets ISO 80369-5 international standards, lowering the risk of improper connections to other types of medical equipment
- A full range of cuff sizes are available, from neonate, pediatric, adult and thigh, ensuring that the measurement you receive is from proper fit

### Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to connect the cuff to the monitor faster than other dual-tube connectors
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with GE and non-GE automated BP devices, as well as on a variety of manual devices, including wall manometers and hand-held gauges, enabling a universal connection system throughout the facility and optimization of inventory space

### Cuff assortment packs available

- For low-volume cuff utilizers, three types of assortment packs are available for easier ordering and more optimal inventory management



SOFT-CUF			ISO 80369-5 compliant		ISO 80369-1 compliant	Does not conform to ISO 80369 standard					
			2-Tube DINACLICK Connector	1-Tube Bayonet Connector	2-Tube DINACLICK Connector	2-Tube Screw Connector	2-Tube Submin Connector	2-Tube Mated Submin Connector	1-Tube Screw Connector	1-Tube Submin Connector	
Size	Color	Limb Circumference									
Infant	Orange	8-13 cm	SFT-P1-2A-5S	2360	SFT-P1-2A	2500	2401	2450 <sup>1</sup>	002331	2320	
Child	Green	12-19 cm	SFT-P2-2A-5S	2361	SFT-P2-2A	2501	2402	2451	002332	2321	
Child Long	Green	12-19 cm	SFT-P2-2A-L-5S		SFT-P2-2A-L	2506	2400				
Small Adult	Light Blue	17-25 cm	SFT-A1-2A-5S	2362	SFT-A1-2A	2502	2403	2452	002333	2322	
Small Adult Long	Light Blue	17-25 cm	SFT-A1-2A-L-5S		SFT-A1-2A-L	2507	2407		002334		
Adult	Navy	23-33 cm	SFT-A2-2A-5S	2363	SFT-A2-2A	2503	2404	2453	002335	2323	
Adult Long	Navy	23-33 cm	SFT-A2-2A-L-5S	2140	SFT-A2-2A-L	2604	2116	2454	002336		
Large Adult	Rose	31-40 cm	SFT-A3-2A-5S	2364	SFT-A3-2A	2504	2405	2455	002337	2324	
Large Adult Long	Rose	31-40 cm	SFT-A3-2A-L-5S	2141	SFT-A3-2A-L	2644	2117	2456	002338		
Thigh	Brown	38-50 cm	SFT-T1-2A-5S	2365	SFT-T1-2A	2505	2406	2457	002339	2325	

Note: 20 cuffs/box

### DINACLICK assortment packs

Connector Type	Part Number (ISO 80369-5 compliant)	Part Number (ISO 80369-1 compliant)	Description
2-Tube DINACLICK	2059303-501	2059303-001	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)
	2059303-502	2059303-002	6 Infant, 10 Child, 4 Small Adult (20 cuffs/box)
	2059303-503	2059303-003	3 Small Adult, 7 Adult, 3 Adult Long, 3 Large Adult, 4 RADIAL-CUF (forearm) (20 cuffs/box)
	2059303-504	2059303-004	2 Infant, 2 Child, 2 Small Adult, 6 Adult, 2 Adult Long, 2 Large Adult, 2 RADIAL-CUF (forearm), 2 Thigh (20 cuffs/box)
	2059303-506	2059303-006	1 Small Adult, 1 Adult, 1 Large Adult, 1 RADIAL-CUF (4 cuffs/pack)

### Other connector assortment packs

Connector Type	Part Number	Description
2-Tube Screw	2695	1 Infant, 1 Child, 2 Small Adult, 2 Adult, 1 Adult Long, 2 Large Adult, 1 Thigh (10 cuffs/box)
	2059303-005	1 Small Adult, 1 Adult, 1 Large Adult, 1 RADIAL-CUF (4 cuffs/pack)
2-Tube Submin	2056304-002	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)
	2298	1 Infant, 1 Child, 2 Small Adult, 2 Adult, 1 Adult Long, 2 Large Adult, 1 Thigh (10 cuffs/box)
	2059303-007	Small Adult, 1 Adult, 1 Large Adult, 1 RADIAL-CUF (4 cuffs/pack)

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at [cuffs.gehealthcare.com](https://cuffs.gehealthcare.com)



1. Infant BP Cuff for Datex-Ohmeda™ monitors has Male Slip Luer Connector.

© 2020 General Electric Company - All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, CRITIKON, DURA-CUF, DINACLICK and Datex-Ohmeda are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB00418XX 8/20