



EMR Gateway Pro for MAC 2000

Simple, flexible ECG-to-EMR connectivity



Complex healthcare environments depend on simple, reliable connectivity that helps clinicians be productive while providing excellent patient care. EMR Gateway Pro is an easy-to-use ECG-to-EMR connectivity solution that provides seamless clinical workflow and simple, secure, and efficient data access. This software utility allows physician offices to effectively move clinical procedure data into the EMR using industry standard formats and protocols.

The EMR Gateway Pro solution seamlessly connects the MAC™ 2000 to your existing electronic medical records system, immediately making it easier to acquire, transmit, store, and access ECG data.

- Bi-directional connectivity with EMR. Orders and patient demographics are sent seamlessly to the MAC 2000, and results and test data are sent back to your EMR.
- Simplified workflow. After installation and configuration, the software operates in the background with no need for user intervention, saving time and minimizing errors.
- Industry-standard technology. EMR Gateway Pro uses well-established secure communication and data transfer protocols, such as HL7®, XML, and PDF, so you can seamlessly connect your MAC 2000 to EMRs.
- Easy installation. Uses familiar Microsoft® Windows®-based install technology, including flexible templates and straightforward configuration tools to simplify the process.
- Dependable service and support. Our experienced support team is always just a call away. Fast, remote diagnostics and support are available through InSite ExC, GE Healthcare's encrypted data connection for service needs.



EMR Gateway Pro: Operating Requirements
Minimum Specifications for Physical or Virtual Machine Server

Feature	Specification
Hardware	
CPU	2.0 GHz or faster 32-bit (x86) or 64-bit (x64) processor
Memory	2 GB minimum reserved
Hard Disk Drive	30 GB
Optical Drive	1 DVD-ROM
Network Port	1 100 Mbps or faster
Video Port	1
Keyboard Port	1
Mouse Port	1
Internet access	Required
Software	
Operating System	Microsoft® Windows® 7 SP1 32-bit or 64-bit, or Microsoft Windows Server® 2008 R2 SP1
Components installed by EMR Gateway Pro	C Runtime Libraries (multiple versions) .NET Framework 4.0 Full Java 7 Microsoft SQL Server® 2012 Express Edition with Tools IIS 7 with ASP.NET
Internet Browsers*	Microsoft Internet Explorer® IE 9, IE 10, IE 11 Mozilla Firefox® v26

EMR Gateway Pro HL7 Options

HL7 interface:

- Inbound ADT
- Inbound ORM (orders)
- Outbound results ORU with embedded pdf
- Outbound results MDM with UNC path to the pdf file

*Note: The Internet Browsers that are listed are those that GE Healthcare has tested the EMR Gateway Pro Application with. When using later versions of these Browsers, or when using other Internet Browsers, reference the EMR Gateway Pro Administrator's Manual to insure that all functionality is accessible.

Imagination at work

GE Healthcare
 P.O. Box 900,
 FIN-00031 GE, Finland
 GE Direct United Kingdom: +44 (0) 800 0329201
 www.gehealthcare.com

©2015 General Electric Company - All rights reserved.
 GE, GE Monogram and SEER are trademarks of General Electric Company.
 GE Healthcare, a division of General Electric Company.
 GE, the GE Monogram, Centricity, MAC, and imagination at work are trademarks of General Electric Company.
 HL7 is a registered trademark of Health Level Seven, Inc.
 Microsoft, Internet Explorer, SQL Server, Windows, and Windows Server are trademarks of Microsoft Corporation.
 Mozilla and Firefox are trademarks of Mozilla Corporation.
 JB31884GB 06/15