



CRITIKON DURA-CUF

Blood Pressure Cuffs

CRITIKON™ DURA-CUF™ is gentle and designed for durability. Made with tough, flexible nylon, they are convenient and easy-to-clean with curved corners to enhance your patient's comfort. They are constructed for multiple patient use in low-acuity care settings and ideal for physicians offices where the risk of contamination is diminished.










Designed for patient safety, comfort and accuracy

- Made with durable materials that can withstand extensive use, and multiple disinfections and cleanings
- Designed with rounded corners to enhance patient comfort
- Constructed with materials not made with PVC or natural rubber latex, making them safer for your patient
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- DINACLICK™ connector meets ISO 80369-5 international standards, lowering the risk of improper connections to other types of medical equipment

Suitable for virtually every physician's office

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the NIBP hose quicker than a pair of traditional tube connectors
- A full range of adapters and connectors are compatible with GE and non-GE automated BP devices, as well as on a variety of manual devices, including wall manometers and hand-held gauges, enabling a universal connection system throughout the facility
- Three types of cuff assortment packs enable easier ordering and optimization of inventory space



| DURA-CUF | | | ISO 80369-5 compliant | | ISO 80369-1 compliant | Does not conform to ISO 80369 standard | | | |
|------------------|------------|--------------------|---|---|---|---|---|---|---|
| | | |  |  |  |  |  |  |  |
| Size | Color | Limb Circumference | 2-Tube DINACLICK Connector | 1-Tube Bayonet Connector | 2-Tube DINACLICK Connector | 2-Tube Screw Connector | 2-Tube Submin Connector | 2-Tube Mated Submin Connector | 1-Tube Screw Connector |
| Infant | Orange | 8-13 cm | DUR-P1-2A-5S | 002274 | DUR-P1-2A | 002783 | 002200 | 2750 ¹ | 002760 |
| Child | Green | 12-19 cm | DUR-P2-2A-5S | 002275 | DUR-P2-2A | 002781 | 002201 | 2751 | 002761 |
| Child Long | Green | 12-19 cm | | | DUR-P2-2A-L | 2785 | | | 002769 |
| Small Adult | Royal Blue | 17-25 cm | DUR-A1-2A-5S | 002276 | DUR-A1-2A | 002779 | 002202 | 2752 | 002762 |
| Small Adult Long | Royal Blue | 17-25 cm | DUR A1-2A-L-5S | | DUR A1-2A-L | | | | 002763 |
| Adult w/ Hanger | Navy | 23-33 cm | | | | 002771 | | | |
| Adult | Navy | 23-33 cm | DUR-A2-2A-5S | 002277 | DUR-A2-2A | 002774 | 002203 | 2753 | 002764 |
| Adult Long | Navy | 23-33 cm | DUR-A2-2A-L-5S | 002130 | DUR-A2-2A-L | 002772 | 002206 | 002756 | 002765 |
| Large Adult | Wine | 31-40 cm | DUR-A3-2A-5S | 002278 | DUR-A3-2A | 002791 | 002204 | 2754 | 002766 |
| Large Adult Long | Wine | 31-40 cm | DUR-A3-2A-L-5S | 002131 | DUR-A3-2A-L | 002784 | 002207 | 002757 | 002767 |
| Thigh | Brown | 38-50 cm | DUR-T1-2A-5S | 002279 | DUR-T1-2A | 002796 | 002205 | 2755 | 002768 |

Note: 5 cuffs/box

DINACLICK assortment packs

| Connector Type | Part Number (ISO 80369-5 compliant) | Part Number (ISO 80369-1 compliant) | Description |
|-----------------------------|-------------------------------------|-------------------------------------|---|
| 2-Tube DINACLICK Connectors | 2059301-501 | 2059301-001 | 1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/box) |
| | 2059301-502 | 2059301-002 | 2 Infant, 3 Child, 1 Small Adult (6 cuffs/box) |
| | 2059301-503 | 2059301-003 | 1 Small Adult, 3 Adult, 1 Adult Long, 1 Large Adult (6 cuffs/box) |
| | 2059301-504 | 2059301-004 | 1 Infant, 1 Child, 1 Small Adult, 1 Adult, 1 Large Adult, 1 Thigh (6 cuffs/box) |

Other connector assortment packs

| Connector Type | Part Number | Description |
|--------------------------|-------------|---|
| 2-Tube Screw Connectors | 002699 | 1 Infant, 1 Child, 1 Small Adult, 1 Adult, 1 Large Adult, 1 Thigh (6 cuffs/box) |
| | 002697 | 2 Infant, 3 Child, 1 Small Adult (6 cuffs/box) |
| | 002698 | 1 Small Adult, 2 Adult, 1 Adult Long, 1 Large Adult, 1 Large Adult Long (6 cuffs/box) |
| | 2059301-005 | 1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/pack) |
| 2-Tube Submin Connectors | 2299 | 1 Infant, 1 Child, 1 Small Adult, 1 Adult, 1 Large Adult, and 1 Thigh (6 cuffs/box) |
| | 2296 | 2 Infant, 3 Child, 1 Small Adult (6 cuffs/box) |
| | 2297 | 1 Small Adult, 2 Adult, 1 Adult Long, 1 Large Adult, 1 Large Adult Long (6 cuffs/box) |
| | 2059301-006 | 1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/pack) |

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com

1. Infant BP Cuff for Datex-Ohmeda™ monitors has Male Slip Luer Connector.

© 2021 General Electric Company - All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, CRITIKON, DURA-CUF, DINACLICK and Datex-Ohmeda are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB00414XX Rev 2 1/21

