

# g

*GE Healthcare*

---

## **Technical Publications**

Document Number: DOC0403010

Revision: 1

***Centricity HIE Services***

**Version 1.00**

## **Integrating the Healthcare Enterprise (IHE) Integration Statement**

Copyrightã 2007 by General Electric Company



# INTRODUCTION

## OVERVIEW

This IHE Integration Statement describes the intended conformance of this GE product with the IHE Technical Framework. The IHE Technical Framework defines a subset of the functional components of the healthcare enterprise, called IHE Actors, and specifies their interactions in terms of a set of coordinated, standards-based transactions. These transactions are organized into functional units called Integration Profiles that highlight their capacity to address specific clinical needs.

Users familiar with these concepts can use Integration Statements as an aid to determine what level of integration this product supports with other systems and what clinical and operational benefits such integration might provide.

Other GE Integration Statements as well as further information on IHE may be found on Internet at: [ge.com/ihe](http://ge.com/ihe)

## IMPORTANT REMARKS

The use of these IHE Integration Statements, **by itself, it is not sufficient to ensure that inter-operation will be successful**. The **user** needs to proceed with caution and ensure that a qualified systems' integration (such as that provided by GE's IT Professional Services) address at least four issues:

- **Integration** - The integration of any system into an integrated department or enterprise requires an analysis of the applications requirements. The design of a solution that integrates GE systems with non-GE systems is the **user's** responsibility and should not be underestimated. Special care should be exercised in defining the expected workflow, the special conditions, and the clinical practice to be supported to ensure a safe and effective operation.
- **Validation** - Testing the complete range of possible interactions between any GE system and non-GE system, before the integration is declared operational, should not be overlooked. Therefore, the **user** should ensure that any non-GE provider accepts full responsibility for all validation required for their connection with GE systems. This includes the accuracy of the data once it has crossed the interface and the stability of the data for the intended applications. Such a validation is required before any clinical use (diagnosis and/or treatment) is performed.
- **Future Evolution** - GE understands that the IHE Technical Framework and the underlying communication standards will evolve to meet the user's growing requirements. GE is actively involved in the development of the IHE Technical Framework and the underlying standards. Such evolutions may require changes to already installed systems. **In addition, GE reserves the right to discontinue or make changes to product features including those reflected by this IHE Integration Statement.** The **user** should ensure that any non-GE system that connects with GE systems also plans for the future evolution of IHE. Failure to do so will likely result in the loss of function and/or connectivity as the IHE Integration profiles evolve and GE Products are enhanced to support these changes.

<b>General Electric - IHE Integration Statement</b>		<b>Date</b>
		11 Dec 2007
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>
GE Healthcare	Centricity HIE Services	1.00
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Cross Enterprise Document Sharing (XDS)	Document Repository	
	Document Registry	
	Document Consumer	Basic Patient Privacy Proof Basic Patient Privacy Enforcement
	Document Source	Basic Patient Privacy Enforcement Document Life Cycle Management
	Patient Identity Source	
Basic Patient Privacy Consents (BPPC)	Content Creator	Basic Patient Privacy Acknowledgement
	Content Consumer	Basic Patient Privacy Acknowledgement View
Patient Identity Cross Referencing (PIX)	Patient Identifier Cross Reference Manager	
	Patient Identity Source	
	Patient Identifier Cross Reference Consumer	PIX Update Notification
Patient Demographic Query (PDQ)	Patient Demographics Supplier	
	Patient Demographics Consumer	
Cross Enterprise Document Sharing of Medical Summaries (XDS-MS)	Content Creator	Discharge Summary Referral Summary
	Content Consumer	View Option Discrete Data Import Option
Exchange of Personal Health Record Information (XPHR)	Content Creator	
	Content Consumer	View Option Discrete Data Import Option
Emergency Department Referral (EDR)	Content Consumer	View Option Discrete Data Import Option
Emergency Department Encounter Summary (EDES)	Content Consumer	View Option Discrete Data Import Option
Antepartum Summary (APS)	Content Consumer	View Option
Sharing Laboratory Reports (XD-LAB)	Content Creator	
	Content Consumer	View Option Discrete Data Import Option
Audit Trail and Node Authentication (ATNA)	Secure Node	Encryption Option
	Audit Repository	
Consistent Time (CT)	Time Client	SNTP
<b>Internet address for GE Healthcare IHE information:</b> <a href="http://www.ge.com/ihe">http://www.ge.com/ihe</a>		
<b>Links to general information on IHE</b>		
In North America: <a href="http://www.ihe.net">www.ihe.net</a> <a href="http://www.himss.org/ihe">www.himss.org/ihe</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>

This page is left blank

g

## **GE Healthcare**

*GE Healthcare – Americas : Fax 262-544-3384  
P.O.Box 414, Milwaukee, Wisconsin, 53201 U.S.A  
Internet – <http://www.gehealthcare.com>*

*GE Healthcare – Europe : Fax 33-1-30-70-95-95  
Paris, France  
Internet - <http://www.gehealthcareeurope.com>*

*GE Healthcare – Asia  
Tokyo, Japan : Fax 81-3-3223-8560  
Singapore : Fax 65-291-7006*