



GE Healthcare

BiliSoft™ Disposable Pad Covers and Pad Nests

Intensive therapy, as easy as wrapping a baby in a blanket

The BiliSoft Pad Covers and Pad Nests are designed for use with the BiliSoft Phototherapy Systems, which are used to treat both premature and full-term infants for hyperbilirubinemia.

Features

- Ultra-soft, cushioned, disposable BiliSoft covers are made of flame-retardant, skin-friendly fabric
- Available in small and large sizes to match the size of the light pad
- Materials used are engineered to allow therapy light to pass through the fabric
- Covers are x-ray compatible

**Our equipment.
Our accessories.
Simply compatible.**

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- Simplified ordering process
- Fast delivery
- One contact point for expert advice
- Quality support at every step



[gehealthcare.com/products/accessories-and-supplies](https://www.gehealthcare.com/products/accessories-and-supplies)

PRODUCT SPECIFICATIONS

Product use

BiliSoft Pad covers are gentle on the baby's skin and provide a soft cushioned covering for all BiliSoft products.

Materials

Top Layer:	Non-woven polypropylene, no dye
Inner Padding:	Polyester fiber batting, no dye
Back Panel:	Polyester Mesh, no dye
Stitching:	Polyester thread

Materials composition

Not made with any BPA, DEHP, natural rubber latex, PVC or other phthalates.

Operating environment

Temperature:	5° to 40°C (41 to 104°F)
Humidity:	up to 95% RH (non-condensing)

Storage

Temperature:	-25° to 70°C (-13° to 158°F)
Humidity:	up to 90% RH (non-condensing)

Disposal considerations

The BiliSoft disposable covers are for single-patient-use only. Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.com.

Packaging materials

Primary: Polybag – Polyethylene (PE)

Regulatory standards

Pad Covers Meet: ISO 10993-5 and ISO 10993-10 for Biocompatibility Requirements and 16 CFR Part 1632.6 for Flammability Requirements.

Compliance

CE Mark, regulation (EU) 2017/745, REACH regulation 1907/2006 compliant, including amendments.

PART LIST

Part Number	Description
2105422-001	BiliSoft 2.0 Disposable Pad Covers, Flat, Small, 20pk
2105418-001	BiliSoft 2.0 Disposable Pad Covers, Flat, Small, 50pk
2107201-001	BiliSoft Disposable Pad Covers, Flat, Small, 20pk
2107200-001	BiliSoft Disposable Pad Covers, Flat, Small, 50pk
2105423-001	BiliSoft 2.0* Disposable Pad Covers, Flat, Large, 20pk
2105419-001	BiliSoft 2.0* Disposable Pad Covers, Flat, Large, 50pk
2105420-001	BiliSoft 2.0 Disposable Nest Pad Covers, Flat, Large, 15pk
2107202-001	BiliSoft Disposable Nest Pad Covers, Small, 15pk
2105421-001	BiliSoft 2.0* Disposable Nest Pad Covers, Large, 15pk

*Note: the BiliSoft 2.0 Disposable Pad (flat and nest) in size large also fits the BiliSoft Phototherapy System (legacy) large pad.

Distributed by:

GE Healthcare

Manufactured by:

Lumitex Medical Devices Inc.

8443 Dow Circle

Strongsville, OH 44136

© 2020 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram and BiliSoft are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2392468 07/20