



CRITIKON CLASSIC-CUF

Blood Pressure Cuffs

CRITIKON™ CLASSIC-CUF™ BP cuffs are designed and clinically validated for limited reuse on multiple patients. Made with durable vinyl construction, this time-tested CRITIKON blood pressure cuff is both easy-to-clean and strong to withstand repeated inflations – yet economical enough for single patient use.



A trusted choice for clinical excellence

- Cuffs are held to the same quality standard as all GE patient monitors, with repeated inflations to ensure accurate readings
- Large, bold print helps you quickly choose the right cuff, enabling more accurate measurement

Constructed for patient safety and durability

- Easy-to-clean, vinyl material does not contain natural rubber latex or DEHP, making them safer for your patients
- A full range of cuff sizes are available, from neonate, pediatric, adult and thigh, ensuring that the measurement you receive is from proper fit
- DINACLICK™ connector meets ISO 80369-5 international standards, lowering the risk of improper connections to other types of medical equipment







Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the NIBP hose quicker than a pair of traditional tube connectors
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with automated monitors and a variety of manual devices, including wall manometers and hand-held gauges, enabling a universal connection system throughout the facility

Cuff assortment packs available

- For low-volume cuff utilizers, several types of assortment packs are available for easier ordering and more optimal inventory management



CLASSIC-CUF			ISO 80369-5 compliant		ISO 80369-1 compliant	Does not conform to ISO 80369 standard		
								
Size	Color	Limb Circumference	2-Tube DINACLICK Connector	1-Tube Bayonet Connector	2-Tube DINACLICK Connector	2-Tube Screw Connector	2-Tube Submin Connector	2-Tube Mated Submin Connector
Infant	Orange	8-13 cm	CLA-P1-2A-5S	2300	CLA-P1-2A	2618	2354	2170 ¹
Child	Green	12-19 cm	CLA-P2-2A-5S	2301	CLA-P2-2A	2613	2355	2171
Small Adult	Light Blue	17-25 cm	CLA-A1-2A-5S	2302	CLA-A1-2A	2608	2356	2172
Adult	Navy	23-33 cm	CLA-A2-2A-5S	2303	CLA-A2-2A	2603	2357	2173
Adult Long	Navy	23-33 cm	CLA-A2-2A-L-5S	2306	CLA-A2-2A-L	2647	2352	2174
Large Adult	Rose	31-40 cm	CLA-A3-2A-5S	2304	CLA-A3-2A	2643	2358	2175
Large Adult Long	Rose	31-40 cm	CLA-A3-2A-L-5S	2307	CLA-A3-2A-L	2649	2353	2176
Thigh	Brown	38-50 cm	CLA-T1-2A-5S	2305	CLA-T1-2A	2648	2359	2177

Note: 20 cuffs/box

DINACLICK assortment packs

Connector Type	Part Number (ISO 80369-5 compliant)	Part Number (ISO 80369-1 compliant)	Description
2-Tube DINACLICK Connectors	2059304-501	2059304-001	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)
	2059304-502	2059304-002	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/pack)

Other connector assortment packs

Connector Type	Part Number	Description
2-Tube Screw Connectors	2692	1 Infant, 1 Child, 2 Small Adult, 2 Adult, 1 Adult Long, 2 Large Adult, 1 Thigh (10 cuff/box)
2-Tube Submin Connectors	2292	1 Infant, 1 Child, 2 Small Adult, 2 Adult, 1 Adult Long, 2 Large Adult, 1 Thigh (10 cuffs/box)

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com



1. Infant BP Cuff for Datex-Ohmeda™ monitors has Male Slip Luer Connector.

© 2020 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, CLASSIC-CUF, CRITIKON, DINACLICK and Datex-Ohmeda are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.