



INVASIVE CARDIOLOGY SECURITY WEBSITE – SECURITY PATCH BASE CONFIGURATIONS FOR ALTiX / ALTiX BT21

Interventional - Invasive Cardiology

| | |
|-----------------------|--|
| Product Group: | Interventional Invasive Products |
| Products: | Mac-Lab , CardioLab and ComboLab Recording Systems, Centricity Cardiology Data Management Systems |
| Subject: | Security Information – Security Patch Base Configurations for ALtiX / ALtiX BT21 |
| Date: | April 14, 2021 |

Summary

The following information is provided to GE Healthcare Technologies customers in regard to the base security patch configurations for Mac-Lab Hemodynamic, CardioLab Electrophysiology and ComboLab IT Recording Systems for Cath Lab, EP Lab and other interventional labs.

Base Configuration

When a new Mac-Lab, CardioLab or ComboLab IT Recording System version is released the product is validated with all applicable security patches up to the development cut-off date. This document contains the security patch base configuration at release.

For security issues identified after the product release, the Mac-Lab, CardioLab and ComboLab IT Recording Systems go through a monthly security analysis process described within the Invasive Cardiology Security Website. This website lists the applicable security patches qualified in addition to the base configuration presented in this document.



| | | |
|-----------|---|----------|
| 1. | VALIDATED MICROSOFT PATCHES INCLUDED ON 7.0 / 7.1 INSTALLATION | 3 |
| 1.1. | DELL 5820 7.0 / 7.1 WITH MICROSOFT WINDOWS 10 ENTERPRISE 2016 LTSB..... | 3 |
| 1.2. | DELL 5060 7.0 / 7.1 WITH MICROSOFT WINDOWS 10 ENTERPRISE 2016 LTSB..... | 5 |
| 1.3. | DELL 5070 7.0 / 7.1 WITH MICROSOFT WINDOWS 10 ENTERPRISE 2016 LTSB..... | 7 |
| 1.4. | DELL 7080 7.0 / 7.1 WITH MICROSOFT WINDOWS 10 ENTERPRISE LTSC 2018..... | 9 |



1. VALIDATED MICROSOFT PATCHES INCLUDED ON 7.0 / 7.1 INSTALLATION

1.1. DELL 5820 7.0 / 7.1 with Microsoft Windows 10 Enterprise 2016 LTSC

| # | Windows Update Description |
|----|---|
| 1 | Hotfix 3335 for Microsoft SQL Server Browser (KB4557397) |
| 2 | Definition Update for Microsoft Office 2016 (KB3115407) 64-Bit Edition |
| 3 | Security Update for Microsoft Access 2016 (KB4484167) 64-Bit Edition |
| 4 | Security Update for Microsoft Office 2016 (KB3085538) 64-Bit Edition |
| 5 | Security Update for Microsoft Office 2016 (KB3213551) 64-Bit Edition |
| 6 | Security Update for Microsoft Office 2016 (KB4011574) 64-Bit Edition |
| 7 | Security Update for Microsoft Office 2016 (KB4022162) 64-Bit Edition |
| 8 | Security Update for Microsoft Office 2016 (KB4022176) 64-Bit Edition |
| 9 | Security Update for Microsoft Office 2016 (KB4484258) 64-Bit Edition |
| 10 | Security Update for Microsoft Outlook 2016 (KB4484433) 64-Bit Edition |
| 11 | Security Update for Microsoft Project 2016 (KB4484441) 64-Bit Edition |
| 12 | Security Update for Microsoft Publisher 2016 (KB4011097) 64-Bit Edition |
| 13 | Security Update for Microsoft Visio 2016 (KB4484244) 64-Bit Edition |
| 14 | Security Update for Microsoft Word 2016 (KB4484438) 64-Bit Edition |
| 15 | Update for Microsoft Excel 2016 (KB4484437) 64-Bit Edition |
| 16 | Update for Microsoft Office 2016 (KB2920678) 64-Bit Edition |
| 17 | Update for Microsoft Office 2016 (KB2920712) 64-Bit Edition |
| 18 | Update for Microsoft Office 2016 (KB2920717) 64-Bit Edition |
| 19 | Update for Microsoft Office 2016 (KB2920720) 64-Bit Edition |
| 20 | Update for Microsoft Office 2016 (KB2920724) 64-Bit Edition |
| 21 | Update for Microsoft Office 2016 (KB3114528) 64-Bit Edition |
| 22 | Update for Microsoft Office 2016 (KB3114852) 64-Bit Edition |
| 23 | Update for Microsoft Office 2016 (KB3114903) 64-Bit Edition |
| 24 | Update for Microsoft Office 2016 (KB3115081) 64-Bit Edition |
| 25 | Update for Microsoft Office 2016 (KB3115276) 64-Bit Edition |
| 26 | Update for Microsoft Office 2016 (KB3118262) 64-Bit Edition |
| 27 | Update for Microsoft Office 2016 (KB3118263) 64-Bit Edition |
| 38 | Update for Microsoft Office 2016 (KB3118264) 64-Bit Edition |
| 39 | Update for Microsoft Office 2016 (KB3141456) 64-Bit Edition |
| 40 | Update for Microsoft Office 2016 (KB3178666) 64-Bit Edition |
| 41 | Update for Microsoft Office 2016 (KB3191929) 64-Bit Edition |
| 42 | Update for Microsoft Office 2016 (KB3213650) 64-Bit Edition |
| 42 | Update for Microsoft Office 2016 (KB4011035) 64-Bit Edition |
| 43 | Update for Microsoft Office 2016 (KB4011225) 64-Bit Edition |
| 44 | Update for Microsoft Office 2016 (KB4011259) 64-Bit Edition |



| | |
|----|---|
| 45 | Update for Microsoft Office 2016 (KB4011629) 64-Bit Edition |
| 46 | Update for Microsoft Office 2016 (KB4011634) 64-Bit Edition |
| 47 | Update for Microsoft Office 2016 (KB4011669) 64-Bit Edition |
| 48 | Update for Microsoft Office 2016 (KB4022161) 64-Bit Edition |
| 49 | Update for Microsoft Office 2016 (KB4022193) 64-Bit Edition |
| 48 | Update for Microsoft Office 2016 (KB4032254) 64-Bit Edition |
| 49 | Update for Microsoft Office 2016 (KB4461435) 64-Bit Edition |
| 50 | Update for Microsoft Office 2016 (KB4462119) 64-Bit Edition |
| 51 | Update for Microsoft Office 2016 (KB4464535) 64-Bit Edition |
| 52 | Update for Microsoft Office 2016 (KB4464538) 64-Bit Edition |
| 53 | Update for Microsoft Office 2016 (KB4464586) 64-Bit Edition |
| 54 | Update for Microsoft Office 2016 (KB4475580) 64-Bit Edition |
| 55 | Update for Microsoft Office 2016 (KB4484145) 64-Bit Edition |
| 56 | Update for Microsoft Office 2016 (KB4484171) 64-Bit Edition |
| 57 | Update for Microsoft Office 2016 (KB4484174) 64-Bit Edition |
| 58 | Update for Microsoft Office 2016 (KB4484327) 64-Bit Edition |
| 59 | Update for Microsoft Office 2016 (KB4484328) 64-Bit Edition |
| 60 | Update for Microsoft Office 2016 (KB4484335) 64-Bit Edition |
| 61 | Update for Microsoft Office 2016 (KB4484392) 64-Bit Edition |
| 62 | Update for Microsoft Office 2016 (KB4484394) 64-Bit Edition |
| 63 | Update for Microsoft Office 2016 (KB4484439) 64-Bit Edition |
| 64 | Update for Microsoft OneDrive for Business (KB4022219) 64-Bit Edition |
| 65 | Update for Microsoft OneNote 2016 (KB4484329) 64-Bit Edition |
| 66 | Update for Microsoft PowerPoint 2016 (KB4484337) 64-Bit Edition |
| 67 | Update for Microsoft Visio Viewer 2016 (KB2920709) 64-Bit Edition |
| 68 | Update for Skype for Business 2016 (KB4484326) 64-Bit Edition |
| 69 | Hotfix 3006 for SQL Server 2017 (KB4038634) (64-bit) |
| 70 | Hotfix 3035 for SQL Server 2017 (KB4293805) (64-bit) |
| 71 | Hotfix 3076 for SQL Server 2017 (KB4484710) (64-bit) |
| 72 | Hotfix 3103 for SQL Server 2017 (KB4494352) (64-bit) |
| 73 | Hotfix 3192 for SQL Server 2017 (KB4505225) (64-bit) |
| 74 | Hotfix 3335 for SQL Server 2017 (KB4557397) (64-bit) |
| 75 | Microsoft Visual C++ 2010 x64 Redistributable (KB2565063) |
| 76 | Update for Microsoft Visual Studio 2015 (KB3095681) |
| 77 | Hotfix 3335 for Microsoft SQL Server VSS Writer (64-bit) (KB4557397) |
| 78 | Security Update for Adobe Flash Player (KB4561600) |
| 79 | Security Update for Microsoft Windows (KB4485447) |
| 80 | Security Update for Microsoft Windows (KB4498947) |
| 81 | Security Update for Microsoft Windows (KB4565511) |



| | |
|----|---|
| 82 | Security Update for Microsoft Windows (KB4565912) |
| 83 | Update for Microsoft Windows (KB2693643) |
| 84 | Update for Microsoft Windows (KB4054590) |
| 85 | Update for Microsoft Windows (KB4494175) |
| 86 | Update for (KB2504637) |

1.2. DELL 5060 7.0 / 7.1 with Microsoft Windows 10 Enterprise 2016 LTSC

| # | Windows Update Description |
|----|---|
| 1 | Hotfix 3335 for Microsoft SQL Server Browser (KB4557397) |
| 2 | Definition Update for Microsoft Office 2016 (KB3115407) 64-Bit Edition |
| 3 | Security Update for Microsoft Access 2016 (KB4484167) 64-Bit Edition |
| 4 | Security Update for Microsoft Office 2016 (KB3085538) 64-Bit Edition |
| 5 | Security Update for Microsoft Office 2016 (KB3213551) 64-Bit Edition |
| 6 | Security Update for Microsoft Office 2016 (KB4011574) 64-Bit Edition |
| 7 | Security Update for Microsoft Office 2016 (KB4022162) 64-Bit Edition |
| 8 | Security Update for Microsoft Office 2016 (KB4022176) 64-Bit Edition |
| 9 | Security Update for Microsoft Office 2016 (KB4484258) 64-Bit Edition |
| 10 | Security Update for Microsoft Outlook 2016 (KB4484433) 64-Bit Edition |
| 11 | Security Update for Microsoft Project 2016 (KB4484441) 64-Bit Edition |
| 12 | Security Update for Microsoft Publisher 2016 (KB4011097) 64-Bit Edition |
| 13 | Security Update for Microsoft Visio 2016 (KB4484244) 64-Bit Edition |
| 14 | Security Update for Microsoft Word 2016 (KB4484438) 64-Bit Edition |
| 15 | Update for Microsoft Excel 2016 (KB4484437) 64-Bit Edition |
| 16 | Update for Microsoft Office 2016 (KB2920678) 64-Bit Edition |
| 17 | Update for Microsoft Office 2016 (KB2920712) 64-Bit Edition |
| 18 | Update for Microsoft Office 2016 (KB2920717) 64-Bit Edition |
| 19 | Update for Microsoft Office 2016 (KB2920720) 64-Bit Edition |
| 20 | Update for Microsoft Office 2016 (KB2920724) 64-Bit Edition |
| 21 | Update for Microsoft Office 2016 (KB3114528) 64-Bit Edition |
| 22 | Update for Microsoft Office 2016 (KB3114852) 64-Bit Edition |
| 23 | Update for Microsoft Office 2016 (KB3114903) 64-Bit Edition |
| 24 | Update for Microsoft Office 2016 (KB3115081) 64-Bit Edition |
| 25 | Update for Microsoft Office 2016 (KB3115276) 64-Bit Edition |
| 26 | Update for Microsoft Office 2016 (KB3118262) 64-Bit Edition |
| 27 | Update for Microsoft Office 2016 (KB3118263) 64-Bit Edition |
| 38 | Update for Microsoft Office 2016 (KB3118264) 64-Bit Edition |
| 39 | Update for Microsoft Office 2016 (KB3141456) 64-Bit Edition |
| 40 | Update for Microsoft Office 2016 (KB3178666) 64-Bit Edition |
| 41 | Update for Microsoft Office 2016 (KB3191929) 64-Bit Edition |



| | |
|----|---|
| 42 | Update for Microsoft Office 2016 (KB3213650) 64-Bit Edition |
| 42 | Update for Microsoft Office 2016 (KB4011035) 64-Bit Edition |
| 43 | Update for Microsoft Office 2016 (KB4011225) 64-Bit Edition |
| 44 | Update for Microsoft Office 2016 (KB4011259) 64-Bit Edition |
| 45 | Update for Microsoft Office 2016 (KB4011629) 64-Bit Edition |
| 46 | Update for Microsoft Office 2016 (KB4011634) 64-Bit Edition |
| 47 | Update for Microsoft Office 2016 (KB4011669) 64-Bit Edition |
| 48 | Update for Microsoft Office 2016 (KB4022161) 64-Bit Edition |
| 49 | Update for Microsoft Office 2016 (KB4022193) 64-Bit Edition |
| 48 | Update for Microsoft Office 2016 (KB4032254) 64-Bit Edition |
| 49 | Update for Microsoft Office 2016 (KB4461435) 64-Bit Edition |
| 50 | Update for Microsoft Office 2016 (KB4462119) 64-Bit Edition |
| 51 | Update for Microsoft Office 2016 (KB4464535) 64-Bit Edition |
| 52 | Update for Microsoft Office 2016 (KB4464538) 64-Bit Edition |
| 53 | Update for Microsoft Office 2016 (KB4464586) 64-Bit Edition |
| 54 | Update for Microsoft Office 2016 (KB4475580) 64-Bit Edition |
| 55 | Update for Microsoft Office 2016 (KB4484145) 64-Bit Edition |
| 56 | Update for Microsoft Office 2016 (KB4484171) 64-Bit Edition |
| 57 | Update for Microsoft Office 2016 (KB4484174) 64-Bit Edition |
| 58 | Update for Microsoft Office 2016 (KB4484327) 64-Bit Edition |
| 59 | Update for Microsoft Office 2016 (KB4484328) 64-Bit Edition |
| 60 | Update for Microsoft Office 2016 (KB4484335) 64-Bit Edition |
| 61 | Update for Microsoft Office 2016 (KB4484392) 64-Bit Edition |
| 62 | Update for Microsoft Office 2016 (KB4484394) 64-Bit Edition |
| 63 | Update for Microsoft Office 2016 (KB4484439) 64-Bit Edition |
| 64 | Update for Microsoft OneDrive for Business (KB4022219) 64-Bit Edition |
| 65 | Update for Microsoft OneNote 2016 (KB4484329) 64-Bit Edition |
| 66 | Update for Microsoft PowerPoint 2016 (KB4484337) 64-Bit Edition |
| 67 | Update for Microsoft Visio Viewer 2016 (KB2920709) 64-Bit Edition |
| 68 | Update for Skype for Business 2016 (KB4484326) 64-Bit Edition |
| 69 | Hotfix 3006 for SQL Server 2017 (KB4038634) (64-bit) |
| 70 | Hotfix 3035 for SQL Server 2017 (KB4293805) (64-bit) |
| 71 | Hotfix 3076 for SQL Server 2017 (KB4484710) (64-bit) |
| 72 | Hotfix 3103 for SQL Server 2017 (KB4494352) (64-bit) |
| 73 | Hotfix 3192 for SQL Server 2017 (KB4505225) (64-bit) |
| 74 | Hotfix 3335 for SQL Server 2017 (KB4557397) (64-bit) |
| 75 | Update for Microsoft Visual Studio 2015 (KB3095681) |
| 76 | Hotfix 3335 for Microsoft SQL Server VSS Writer (64-bit) (KB4557397) |
| 77 | Security Update for Adobe Flash Player (KB4561600) |



| | |
|----|---|
| 78 | Security Update for Microsoft Windows (KB4485447) |
| 79 | Security Update for Microsoft Windows (KB4498947) |
| 80 | Security Update for Microsoft Windows (KB4565511) |
| 81 | Security Update for Microsoft Windows (KB4565912) |
| 82 | Update for Microsoft Windows (KB2693643) |
| 83 | Update for Microsoft Windows (KB4054590) |
| 84 | Update for Microsoft Windows (KB4494175) |
| 85 | Update for (KB2504637) |

1.3. DELL 5070 7.0 / 7.1 with Microsoft Windows 10 Enterprise 2016 LTSB

| # | Windows Update Description |
|----|---|
| 1 | Hotfix 3335 for Microsoft SQL Server Browser (KB4557397) |
| 2 | Definition Update for Microsoft Office 2016 (KB3115407) 64-Bit Edition |
| 3 | Security Update for Microsoft Access 2016 (KB4484167) 64-Bit Edition |
| 4 | Security Update for Microsoft Office 2016 (KB3085538) 64-Bit Edition |
| 5 | Security Update for Microsoft Office 2016 (KB3213551) 64-Bit Edition |
| 6 | Security Update for Microsoft Office 2016 (KB4011574) 64-Bit Edition |
| 7 | Security Update for Microsoft Office 2016 (KB4022162) 64-Bit Edition |
| 8 | Security Update for Microsoft Office 2016 (KB4022176) 64-Bit Edition |
| 9 | Security Update for Microsoft Office 2016 (KB4484258) 64-Bit Edition |
| 10 | Security Update for Microsoft Outlook 2016 (KB4484433) 64-Bit Edition |
| 11 | Security Update for Microsoft Project 2016 (KB4484441) 64-Bit Edition |
| 12 | Security Update for Microsoft Publisher 2016 (KB4011097) 64-Bit Edition |
| 13 | Security Update for Microsoft Visio 2016 (KB4484244) 64-Bit Edition |
| 14 | Security Update for Microsoft Word 2016 (KB4484438) 64-Bit Edition |
| 15 | Update for Microsoft Excel 2016 (KB4484437) 64-Bit Edition |
| 16 | Update for Microsoft Office 2016 (KB2920678) 64-Bit Edition |
| 17 | Update for Microsoft Office 2016 (KB2920712) 64-Bit Edition |
| 18 | Update for Microsoft Office 2016 (KB2920717) 64-Bit Edition |
| 19 | Update for Microsoft Office 2016 (KB2920720) 64-Bit Edition |
| 20 | Update for Microsoft Office 2016 (KB2920724) 64-Bit Edition |
| 21 | Update for Microsoft Office 2016 (KB3114528) 64-Bit Edition |
| 22 | Update for Microsoft Office 2016 (KB3114852) 64-Bit Edition |
| 23 | Update for Microsoft Office 2016 (KB3114903) 64-Bit Edition |
| 24 | Update for Microsoft Office 2016 (KB3115081) 64-Bit Edition |
| 25 | Update for Microsoft Office 2016 (KB3115276) 64-Bit Edition |
| 26 | Update for Microsoft Office 2016 (KB3118262) 64-Bit Edition |
| 27 | Update for Microsoft Office 2016 (KB3118263) 64-Bit Edition |
| 38 | Update for Microsoft Office 2016 (KB3118264) 64-Bit Edition |



| | |
|----|---|
| 39 | Update for Microsoft Office 2016 (KB3141456) 64-Bit Edition |
| 40 | Update for Microsoft Office 2016 (KB3178666) 64-Bit Edition |
| 41 | Update for Microsoft Office 2016 (KB3191929) 64-Bit Edition |
| 42 | Update for Microsoft Office 2016 (KB3213650) 64-Bit Edition |
| 42 | Update for Microsoft Office 2016 (KB4011035) 64-Bit Edition |
| 43 | Update for Microsoft Office 2016 (KB4011225) 64-Bit Edition |
| 44 | Update for Microsoft Office 2016 (KB4011259) 64-Bit Edition |
| 45 | Update for Microsoft Office 2016 (KB4011629) 64-Bit Edition |
| 46 | Update for Microsoft Office 2016 (KB4011634) 64-Bit Edition |
| 47 | Update for Microsoft Office 2016 (KB4011669) 64-Bit Edition |
| 48 | Update for Microsoft Office 2016 (KB4022161) 64-Bit Edition |
| 49 | Update for Microsoft Office 2016 (KB4022193) 64-Bit Edition |
| 48 | Update for Microsoft Office 2016 (KB4032254) 64-Bit Edition |
| 49 | Update for Microsoft Office 2016 (KB4461435) 64-Bit Edition |
| 50 | Update for Microsoft Office 2016 (KB4462119) 64-Bit Edition |
| 51 | Update for Microsoft Office 2016 (KB4464535) 64-Bit Edition |
| 52 | Update for Microsoft Office 2016 (KB4464538) 64-Bit Edition |
| 53 | Update for Microsoft Office 2016 (KB4464586) 64-Bit Edition |
| 54 | Update for Microsoft Office 2016 (KB4475580) 64-Bit Edition |
| 55 | Update for Microsoft Office 2016 (KB4484145) 64-Bit Edition |
| 56 | Update for Microsoft Office 2016 (KB4484171) 64-Bit Edition |
| 57 | Update for Microsoft Office 2016 (KB4484174) 64-Bit Edition |
| 58 | Update for Microsoft Office 2016 (KB4484327) 64-Bit Edition |
| 59 | Update for Microsoft Office 2016 (KB4484328) 64-Bit Edition |
| 60 | Update for Microsoft Office 2016 (KB4484335) 64-Bit Edition |
| 61 | Update for Microsoft Office 2016 (KB4484392) 64-Bit Edition |
| 62 | Update for Microsoft Office 2016 (KB4484394) 64-Bit Edition |
| 63 | Update for Microsoft Office 2016 (KB4484439) 64-Bit Edition |
| 64 | Update for Microsoft OneDrive for Business (KB4022219) 64-Bit Edition |
| 65 | Update for Microsoft OneNote 2016 (KB4484329) 64-Bit Edition |
| 66 | Update for Microsoft PowerPoint 2016 (KB4484337) 64-Bit Edition |
| 67 | Update for Microsoft Visio Viewer 2016 (KB2920709) 64-Bit Edition |
| 68 | Update for Skype for Business 2016 (KB4484326) 64-Bit Edition |
| 69 | Hotfix 3006 for SQL Server 2017 (KB4038634) (64-bit) |
| 70 | Hotfix 3035 for SQL Server 2017 (KB4293805) (64-bit) |
| 71 | Hotfix 3076 for SQL Server 2017 (KB4484710) (64-bit) |
| 72 | Hotfix 3103 for SQL Server 2017 (KB4494352) (64-bit) |
| 73 | Hotfix 3192 for SQL Server 2017 (KB4505225) (64-bit) |
| 74 | Hotfix 3335 for SQL Server 2017 (KB4557397) (64-bit) |



| | |
|----|--|
| 75 | Update for Microsoft Visual Studio 2015 (KB3095681) |
| 76 | Hotfix 3335 for Microsoft SQL Server VSS Writer (64-bit) (KB4557397) |
| 77 | Security Update for Adobe Flash Player (KB4561600) |
| 78 | Security Update for Microsoft Windows (KB4485447) |
| 79 | Security Update for Microsoft Windows (KB4498947) |
| 80 | Security Update for Microsoft Windows (KB4565511) |
| 81 | Security Update for Microsoft Windows (KB4565912) |
| 82 | Update for Microsoft Windows (KB2693643) |
| 83 | Update for Microsoft Windows (KB4054590) |
| 84 | Update for Microsoft Windows (KB4494175) |
| 85 | Update for (KB2504637) |

1.4. DELL 7080 7.0 / 7.1 with Microsoft Windows 10 Enterprise LTSC 2018

| # | Windows Update Description |
|----|--|
| 1 | Definition Update for Microsoft Office 2016 (KB3115407) 64-Bit Edition |
| 2 | Hotfix 3370 for SQL Server 2017 (KB4583457) (64-bit) |
| 3 | Hotfix 3381 for Microsoft SQL Server Browser (KB5000685) |
| 4 | Hotfix 3381 for Microsoft SQL Server VSS Writer (64-bit) (KB5000685) |
| 5 | Hotfix 3381 for SQL Server 2017 (KB5000685) (64-bit) |
| 6 | Microsoft .NET Framework 4 KB2468871 |
| 7 | Microsoft .NET Framework 4 KB2468871v2 |
| 8 | Microsoft .NET Framework 4 KB2478063 |
| 9 | Microsoft .NET Framework 4 KB2504637 |
| 10 | Microsoft .NET Framework 4 KB2533523 |
| 11 | Microsoft .NET Framework 4 KB2544514 |
| 12 | Microsoft .NET Framework 4 KB2600211 |
| 13 | Microsoft .NET Framework 4 KB2600217 |
| 14 | Security Update for Microsoft Access 2016 (KB4484509) 64-Bit Edition |
| 15 | Security Update for Microsoft Office 2016 (KB3085538) 64-Bit Edition |
| 16 | Security Update for Microsoft Office 2016 (KB3213551) 64-Bit Edition |
| 17 | Security Update for Microsoft Office 2016 (KB4022162) 64-Bit Edition |
| 18 | Security Update for Microsoft Office 2016 (KB4022176) 64-Bit Edition |
| 19 | Security Update for Microsoft Office 2016 (KB4486755) 64-Bit Edition |
| 20 | Security Update for Microsoft Office 2016 (KB4493168) 64-Bit Edition |
| 21 | Security Update for Microsoft Office 2016 (KB4011574) 64-Bit Edition |
| 22 | Security Update for Microsoft Windows (KB4601345) |
| 23 | Security Update for Microsoft Windows (KB4601393) |
| 24 | Security Update for Microsoft Word 2016 (KB4493156) 64-Bit Edition |
| 25 | Update for (KB2504637) |



| | |
|----|---|
| 26 | Update for Microsoft Excel 2016 (KB4493196) 64-Bit Edition |
| 27 | Update for Microsoft Office 2016 (KB2920678) 64-Bit Edition |
| 28 | Update for Microsoft Office 2016 (KB2920712) 64-Bit Edition |
| 29 | Update for Microsoft Office 2016 (KB2920717) 64-Bit Edition |
| 30 | Update for Microsoft Office 2016 (KB2920720) 64-Bit Edition |
| 31 | Update for Microsoft Office 2016 (KB2920724) 64-Bit Edition |
| 32 | Update for Microsoft Office 2016 (KB3114852) 64-Bit Edition |
| 33 | Update for Microsoft Office 2016 (KB3114903) 64-Bit Edition |
| 34 | Update for Microsoft Office 2016 (KB3115081) 64-Bit Edition |
| 35 | Update for Microsoft Office 2016 (KB3115276) 64-Bit Edition |
| 36 | Update for Microsoft Office 2016 (KB3118262) 64-Bit Edition |
| 37 | Update for Microsoft Office 2016 (KB3118263) 64-Bit Edition |
| 38 | Update for Microsoft Office 2016 (KB3118264) 64-Bit Edition |
| 39 | Update for Microsoft Office 2016 (KB3141456) 64-Bit Edition |
| 40 | Update for Microsoft Office 2016 (KB3178666) 64-Bit Edition |
| 41 | Update for Microsoft Office 2016 (KB3191929) 64-Bit Edition |
| 42 | Update for Microsoft Office 2016 (KB3213650) 64-Bit Edition |
| 43 | Update for Microsoft Office 2016 (KB4011035) 64-Bit Edition |
| 44 | Update for Microsoft Office 2016 (KB4011225) 64-Bit Edition |
| 45 | Update for Microsoft Office 2016 (KB4011259) 64-Bit Edition |
| 46 | Update for Microsoft Office 2016 (KB4011629) 64-Bit Edition |
| 47 | Update for Microsoft Office 2016 (KB4011634) 64-Bit Edition |
| 48 | Update for Microsoft Office 2016 (KB4011669) 64-Bit Edition |
| 49 | Update for Microsoft Office 2016 (KB4022161) 64-Bit Edition |
| 50 | Update for Microsoft Office 2016 (KB4022193) 64-Bit Edition |
| 51 | Update for Microsoft Office 2016 (KB4032254) 64-Bit Edition |
| 52 | Update for Microsoft Office 2016 (KB4461435) 64-Bit Edition |
| 53 | Update for Microsoft Office 2016 (KB4462119) 64-Bit Edition |
| 54 | Update for Microsoft Office 2016 (KB4464535) 64-Bit Edition |
| 55 | Update for Microsoft Office 2016 (KB4464538) 64-Bit Edition |
| 56 | Update for Microsoft Office 2016 (KB4475580) 64-Bit Edition |
| 57 | Update for Microsoft Office 2016 (KB4475584) 64-Bit Edition |
| 58 | Update for Microsoft Office 2016 (KB4484145) 64-Bit Edition |
| 59 | Update for Microsoft Office 2016 (KB4484171) 64-Bit Edition |
| 60 | Update for Microsoft Office 2016 (KB4484174) 64-Bit Edition |
| 61 | Update for Microsoft Office 2016 (KB4484251) 64-Bit Edition |
| 62 | Update for Microsoft Office 2016 (KB4484395) 64-Bit Edition |
| 63 | Update for Microsoft Office 2016 (KB4486668) 64-Bit Edition |
| 64 | Update for Microsoft Office 2016 (KB4486716) 64-Bit Edition |
| 65 | Update for Microsoft Office 2016 (KB4486747) 64-Bit Edition |



| | |
|----|---|
| 66 | Update for Microsoft Office 2016 (KB4493189) 64-Bit Edition |
| 67 | Update for Microsoft OneDrive for Business (KB4022219) 64-Bit Edition |
| 68 | Update for Microsoft OneNote 2016 (KB4484329) 64-Bit Edition |
| 69 | Update for Microsoft Outlook 2016 (KB4493190) 64-Bit Edition |
| 70 | Update for Microsoft PowerPoint 2016 (KB4493164) 64-Bit Edition |
| 71 | Update for Microsoft Project 2016 (KB4493159) 64-Bit Edition |
| 72 | Update for Microsoft Publisher 2016 (KB4484334) 64-Bit Edition |
| 73 | Update for Microsoft Visio 2016 (KB4486709) 64-Bit Edition |
| 74 | Update for Microsoft Visio Viewer 2016 (KB2920709) 64-Bit Edition |
| 75 | Update for Microsoft Visual Studio 2015 (KB3095681) |
| 76 | Update for Microsoft Windows (KB2693643) |
| 77 | Update for Microsoft Windows (KB4601060) |
| 78 | Update for Removal of Adobe Flash Player (KB4577586) |
| 79 | Update for Skype for Business 2016 (KB4486745) 64-Bit Edition |

Contact Information

If you have any additional questions, please contact our Technical Support Department.

*Trademark of General Electric Company.

**All third party trademarks are the property of their respective owners.