

## 1.1 APPLICATION – ADVANTAGE WORKSTATION 4.4

Any application claiming to conform to the Standard shall provide the following certificate and conform to the requirements that are checked.

### Health Level Seven Context Management Standard

#### Declaration of Conformance

#### Version 1.4

#### Application

**It is required that all the boxes below be checked.**

- √ Implements all of the policies, protocols, rules, and constraints defined for applications as specified in the Architecture Specification; Implements the ContextParticipant interface, and uses the ContextManager interface as well as the ContextData, ContextFilter, SecureContextData, SecureBinding, and ContextAction interfaces as necessary as defined in the Architecture Specification.
- √ Implements the user interface policies, protocols, rules, and constraints defined for applications as specified in the User Interface Specification.
- √ Can function as a fully CCOW compliant application without any dependencies on the implementation of any other CCOW applications, including session management applications, or CCOW components.

**At least one of the boxes below must be checked.**

- √ Implements links for one or more standard identity subjects.  
Names of links implemented: \_\_\_\_\_

**The following boxes below may be checked.**

- ∅ Implements links for one or more custom subjects.  
Names of links implemented: \_\_\_\_\_
- ∅ Requests (i.e., is a client of) at least one of the following standard or custom action subjects.  
Names of actions requested: \_\_\_\_\_
- ∅ Uses an authentication repository, including the AuthenticationRepository interface.

**At least one of the boxes below must be checked.**

- √ Employs the technology, subject to all of the policies, protocols, rules, and constraints, defined for applications as specified in the ActiveX Technology Mapping.
- ∅ Employs the technology, subject to all of the policies, protocols, rules, and constraints, defined for applications as specified in the Web Technology Mapping.

