



YELLOW CRITIKON CLASSIC-CUF

Blood Pressure Cuffs

The yellow CRITIKON™ CLASSIC-CUF™ is designed for single patient use within an isolation environment. These vinyl cuffs are priced as disposable-use cuffs and support your infection prevention protocols. They provide the same quality and assurance that you expect from the entire line of CRITIKON cuffs, while tailored for isolation areas.



A trusted choice for clinical excellence

- Cuffs are held to the same quality standard as all GE patient monitors, with repeated inflations to ensure accurate readings
- Large, bold print helps you quickly choose the right cuff, enabling a more accurate measurement
- Yellow color helps make it easy to identify in isolation rooms




Constructed for patient safety and durability

- Easy-to-clean, vinyl material is not made with natural rubber latex or DEHP, making them safer for your patients
- A full range of cuff sizes are available, from infant to adult and thigh, ensuring that the measurement you receive is from a proper fit
- DINACLICK™ connector meets ISO 80369-5 international standards, lowering the risk of improper connections to other types of medical equipment

Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the monitor quicker than other dual-tube connectors
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with automated monitors and a variety of manual devices, enabling a universal connection system throughout the facility inventory management



YELLOW CLASSIC-CUF			ISO 80369-5 compliant	ISO 80369-1 compliant	Does not conform to ISO 80369 standard	
						
Size	Color	Limb Circumference	2-Tube DINACLICK Connector	2-Tube DINACLICK Connector	2-Tube Screw Connector	2-Tube Submin Connector
Infant	Yellow/Blue	8-13 cm	CLA-P1-2A-X-5S	CLA-P1-2A-X	2650	2670
Child	Yellow/Blue	12-19 cm	CLA-P2-2A-X-5S	CLA-P2-2A-X	2651	2671
Small Adult	Yellow/Blue	17-25 cm	CLA-A1-2A-X-5S	CLA-A1-2A-X	2607	2672
Adult	Yellow/Blue	23-33 cm	CLA-A2-2A-X-5S	CLA-A2-2A-X	2602	2673
Large Adult	Yellow/Blue	31-40 cm	CLA-A3-2A-X-5S	CLA-A3-2A-X	2642	2674
Thigh	Yellow/Blue	38-50 cm	CLA-T1-2A-X-5S	CLA-T1-2A-X	2652	2675

Note: 20 cuffs/box

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com



© 2020 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, CRITIKON, CLASSIC-CUF and DINACLICK are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB00415XX 8/20