



AHA compliant

# ECG Cables and Leadwires

## Easy connection. Accurate measurements.

GE ECG cables and leadwires help you reliably make the correct connection to GE devices. These high-quality reusable cables complement your equipment, enabling accurate measurements. GE cables and leadwires undergo strict testing to help ensure quality, reliability and safety. They have been tested and verified as compatible with GE devices to deliver optimum performance and a high standard of care.

### Durable

Resistant jacketing extends cable life, facilitates cleaning and remains flexible to withstand bending and twisting and enhance patient comfort.

### Easy to use

The leadwires and adaptors include color coding and labels on both ends to reduce the errors in setting up the leadwires and save time to the clinician to focus on patient care.

Integrated cable and leadwire set has unique grouping of leadwires that will help in untangling the cables easily and quick application on patient

### Product compatibility

The ECG trunk cables and leadwires are compatible and verified with the following GE equipment:

- MAC™ VU360, MAC 5500, MAC 3500, MAC 2000, MAC 800, MAC 600
- CASE® Cardiac Testing System
- CardioSoft® Cardiac Testing System

### Our equipment. Our accessories. Simply compatible.

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- One contact point for expert advice
- Simplified ordering process
- Fast delivery
- Quality support at every step



## FEATURES

GE ECG leadwires are available in 10- and 14-lead configurations.

Leadwires are electrically shielded to protect ECG signal integrity.

## PRODUCT SPECIFICATIONS

### Materials

---

Does not contain natural rubber latex, DEHP, BPA or soft PVC

ROHS (Restriction of hazardous substances) compliant

Bio-compatibility compliant

### Storage conditions

---

Recommended temperature      -40° to 70°C (41°F to 158°F)

Recommended relative humidity      10% to 95%

### Cleaning solutions

---

Ethanol      96%

Hydrogen peroxide      3%

Sodium Hypochlorite      0.5% NaOCl

Isopropyl Alcohol      70%

Mikrozid®      sensitive wipes

Kleenaseptic®

Cidex® Plus

Sporicidin®

Never autoclave or steam-clean

### Compliance

---

CE Mark

American Heart Association (AHA) standards

## PART LIST

Part Number	Description	MAC VU360	MAC 5500	MAC 3500	MAC 2000	MAC 800	MAC 600	CASE	Cardio-Soft
<b>ECG Trunk Cables - AHA</b>									
2104723-001	Trunk Cable 10-lead AHA			x	x	x	x		
<b>Integrated cable and leadwire set - AHA</b>									
2104727-001	ECG Value Cable 10-lead Cable / Leadwire AHA				x	x	x		
<b>ECG Leadwires - banana - AHA</b>									
2104724-001	Leadwire Set, Base 10, banana, AHA	x	x	x	x	x	x	x	x
2104751-001	Leadwire Set, 4 Add-on, banana, AHA-HEIM	x	x					x	x
2104752-001	Leadwire Set, 3 Add-on, banana, AHA-PED (V3R, V4R, V7)	x	x					x	
<b>ECG Leadwires - grabber - AHA</b>									
2104749-001	Leadwire Set, Base 10, grabber, AHA	x	x	x	x			x	x
2104725-001	Leadwire Set, 4 Add-on, grabber, AHA-HEIM	x	x					x	x
<b>Accessories</b>									
2104729-001	Adapter Clip, Wide Mouth Clip, 14 / pack	x	x	x	x	x	x	x	x
2104767-001	ECG Leadwire port plug 4 / pack	x	x					x	x

**Product may not be available in all countries and regions.**  
 Contact a GE Healthcare Representative for more information.  
 Please visit [www.gehealthcare.com/promotional-locations](http://www.gehealthcare.com/promotional-locations)

*Distributed by:*

GE Medical Systems Information Technologies, Inc.



*Manufactured by:*

Shenzhen Launch  
 Electrical Co., Ltd.

No.F Building, Zhonggangxing  
 Industrial Estate, zhangge Community, Guanlan street,  
 Longhua New district,  
 Shenzhen, China

Imagination at work

© 2018 General Electric Company - All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE and the GE Monogram, are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE, the GE monogram, Imagination at work, MAC, Case and CardioSoft are trademarks of General Electric Company. MikroZid is a trademark of Schülke & Mayr GmbH. Kleenaseptic is a trademark and brand of Metrex Research, LLC. Cidex is a trademark of Johnson & Johnson. Sporidicin is a trademark of Contec, Inc. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2180253