



Thermally sensitive GE ECG Recording Papers

Thermally sensitive chart papers from GE enable accurate tracings from your GE devices. Our comprehensive product line covers a wide range of formats to fit different printer configurations.

- Papers available in rolls or z-fold boxes, blank or with gridlines
- Some papers feature multiple chart versions offered with BPM scales for various patient types
- Large space for entering patient data on certain papers
- Indicators for end of pack and end of roll signal time to change paper

Product compatibility

GE ECG papers are tested and compatible with GE Healthcare equipment.

| Part number | Archiving Years | Diagnostic Cardiology | | | | | | | | | Patient Monitoring | | | | Fetal Monitoring | |
|-------------|-----------------|-----------------------|--------|--------|----------|----------|-----------|-----------|-----------|------------|--------------------|--------|------|------------|------------------|-------|
| | | CASE | MAC™ 5 | MAC™ 7 | MAC™ 600 | MAC™ 800 | MAC™ 2000 | MAC™ 3500 | MAC™ 5500 | MAC™ VU360 | Bx50 | Bx0(i) | B1x5 | V100 VC150 | Corometrics™ 170 | 250cx |
| 2104769-001 | 10 | x | | | | | | x | x | | | | | | | |
| 2104772-001 | 10 | | x | x | | | x | | | x | | | | | | |
| 2104768-001 | 10 | x | | | | | | x | x | | | | | | | |
| 2104773-001 | 8 | | | | x | | | | | | | | | | | |
| 2104775-001 | 8 | | | | | x | | | | | | | | | | |
| 2104770-001 | 10 | | | | | | | x | x | x | | | | | | |
| 2104771-001 | 10 | | x | x | | | x | x | | x | | | | | | |
| 5684683 | 10 | | x | | | | | | | | | | | | | |
| 2106825-001 | 25 | | | | | | | | | | | | x | | | |
| 2106824-001 | 10 | | | | | | | | | | x | | | | | |
| 2106823-001 | 10 | | | | | | | | | | x | x | x | | | |
| 2104907-001 | 10 | | | | | | | | | | | | | | x | x |
| 2104908-001 | 10 | | | | | | | | | | | | | | x | x |

PRODUCT SPECIFICATIONS

Materials

Thermally sensitive paper for accurate tracings and durability for archiving

Medium weight for easy passage through printing devices

Recycled materials used for carton and paper stiffeners

BPA Free

Storage conditions

Before usage

Recommended temperature 10° to 32°C (50° to 90°F)

Paper should be kept dry.

Archive conditions

Recommended temperature < 25°C (77°F)

Relative humidity < 70%

Store in a dark environment, avoiding natural or artificial light.

Avoid contact with chemicals such as plasticizers, oils, solvents, water and adhesives.

Product may not be available in all countries and regions. Full product technical specification is available upon request. Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.com/promotional-locations

Distributed by:

GE Healthcare Finland Oy
Kuortaneenkatu 2
00510 Helsinki, Finland
Europe



Manufactured by:

Graphic Controls
400 Exchange Street
Buffalo, NY 14204
USA

ORDERING INFORMATION

Diagnostic Cardiology and Patient Fetal Monitoring

| Part Number | Description |
|---|--|
| Size A4 Paper ¹ | |
| 2104769-001 | Thermal paper A4, red grid 195mm wide, z-fold, hole queue, 300 sheets, 8 packs |
| 2104772-001 | Thermal paper A4, red grid 205mm wide, z-fold, block queue, 150 sheets, 10 packs |
| Size A5 Paper ¹ | |
| 5684683 | Thermal paper A5, red grid 205mm wide, z-fold, block queue, 150 sheets, 20 packs |
| Size 8.5 x 11 inch Red Grid Paper | |
| 2104768-001 | Thermal paper 8.5" x 11", red grid 200mm wide, z-fold, hole queue, 300 sheets, 8 packs |
| Size 8.5 x 11 inch Red Grid/White Data Area | |
| 2104770-001 | Thermal paper 8.5" x 11", white patient data area, red grid 155mm wide, z-fold, hole queue, 300 sheets, 8 packs |
| 2104771-001 | Thermal paper 8.5" x 11", white patient data area, red grid 155mm wide, z-fold, hole queue, 150 sheets, 16 packs |
| Blank Paper Rolls | |
| 2106825-001 | Thermal paper 57.2 mm x 24.4 m (2.25 in x 80 ft), blank, 10 rolls |
| 2106824-001 | Thermal paper 50 mm x 25.6 m (2 in x 84 ft), blank, 20 rolls |
| 2106823-001 | Thermal paper, 50 mm x 27.4 m (2 in x 90 ft), blank, 48 rolls |
| Fetal Monitor Paper | |
| 2104907-001 | Thermal GE paper fetal monitors 30-240 bpm scaling, 170 series and 250cx series 40 packs/case |
| 2104908-001 | Thermal GE paper fetal monitors 50-210 bpm scaling, 170 series and 250cx series 40 packs/case |
| Size Custom Red Grid Paper | |
| 2104773-001 | Thermal paper 80mm wide, red grid 75mm wide, z-fold, block queue, 280 sheets, 10 packs |
| 2104775-001 | Thermal paper 110mm wide, red grid 100mm wide, z-fold, block queue, 200 sheets, 10 packs, 40 packs /case |

¹ Not available for sale in the US.

© 2022 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE and the GE Monogram, are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2146753 rev5 2/22