



CRITIKON SOFT-CUF

Blood Pressure Cuffs

Our most popular cuff, the CRITIKON SOFT-CUF™ is designed with soft, absorbent material with rounded corners for maximizing patient comfort during extended stays at the hospital. These cuffs are ideal for patients with longer term stays or for protocols where the cuff travels with the patient through their total length of time in the hospital. A limited reuse cuff, providing flexibility for single patient use supporting infection control or cleaned and reused for multi-patient protocols.



Designed for patient safety, comfort and accuracy

- Limited reuse, yet priced economically for pairing with a single patient throughout their hospital stay to help with infection control efforts
- Soft, absorbent material is not made with DEHP or natural rubber latex to help enhance patient safety
- Designed with rounded corners for patient comfort
- Designed and clinically validated for performance accuracy to help you make well-informed decisions
- DINACLICK™ connector meets ISO 80369-5 international standards, lowering the risk of improper connections to other types of medical equipment
- A full range of cuff sizes are available, from neonate, pediatric, adult and thigh, ensuring that the measurement you receive is from proper fit

Designed to improve workflow



- Easy-to-use, fast DINACLICK connector enables nurses to connect the cuff to the monitor faster than other dual-tube connectors
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with GE and non-GE automated BP devices, as well as on a variety of manual devices, including wall manometers and hand-held gauges, enabling a universal connection system throughout the facility and optimization of inventory space

Cuff assortment packs available

- For low-volume cuff utilizers, three types of assortment packs are available for easier ordering and more optimal inventory management

Scan the QR code below to view our
DINACLICK Connection System video



SOFT-CUF			ISO 80369-5 compliant		ISO 80369-1 compliant	Does not conform to ISO 80369 standard				
										
Size	Color	Limb Circumference	2-Tube DINACLICK Connector	1-Tube Bayonet Connector	2-Tube DINACLICK Connector	2-Tube Screw Connector	2-Tube Submin Connector	2-Tube Mated Submin Connector	1-Tube Screw Connector	1-Tube Submin Connector
Infant	Orange	8-13 cm	SFT-P1-2A-5S	2360	SFT-P1-2A	2500	2401	2450 ¹	002331	2320
Child	Green	12-19 cm	SFT-P2-2A-5S	2361	SFT-P2-2A	2501	2402	2451	002332	2321
Child Long	Green	12-19 cm	SFT-P2-2A-L-5S		SFT-P2-2A-L	2506	2400			
Small Adult	Light Blue	17-25 cm	SFT-A1-2A-5S	2362	SFT-A1-2A	2502	2403	2452	002333	2322
Small Adult Long	Light Blue	17-25 cm	SFT-A1-2A-L-5S		SFT-A1-2A-L	2507	2407		002334	
Adult	Navy	23-33 cm	SFT-A2-2A-5S	2363	SFT-A2-2A	2503	2404	2453	002335	2323
Adult Long	Navy	23-33 cm	SFT-A2-2A-L-5S	2140	SFT-A2-2A-L	2604	2116	2454	002336	
Large Adult	Rose	31-40 cm	SFT-A3-2A-5S	2364	SFT-A3-2A	2504	2405	2455	002337	2324
Large Adult Long	Rose	31-40 cm	SFT-A3-2A-L-5S	2141	SFT-A3-2A-L	2644	2117	2456	002338	
Thigh	Brown	38-50 cm	SFT-T1-2A-5S	2365	SFT-T1-2A	2505	2406	2457	002339	2325

Note: 20 cuffs/box

DINACLICK assortment packs

Connector Type	Part Number (ISO 80369-5 compliant)	Part Number (ISO 80369-1 compliant)	Description
2-Tube DINACLICK	2059303-501	2059303-001	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)
	2059303-502	2059303-002	6 Infant, 10 Child, 4 Small Adult (20 cuffs/box)
	2059303-503	2059303-003	3 Small Adult, 7 Adult, 3 Adult Long, 3 Large Adult, 4 RADIAL-CUF (forearm) (20 cuffs/box)
	2059303-504	2059303-004	2 Infant, 2 Child, 2 Small Adult, 6 Adult, 2 Adult Long, 2 Large Adult, 2 RADIAL-CUF (forearm), 2 Thigh (20 cuffs/box)
	2059303-505	2059303-005	1 Small Adult, 1 Adult, 1 Large Adult, 1 RADIAL-CUF (4 cuffs/pack)

Other connector assortment packs

Connector Type	Part Number	Description
2-Tube Screw	2695	1 Infant, 1 Child, 2 Small Adult, 2 Adult, 1 Adult Long, 2 Large Adult, 1 Thigh (10 cuffs/box)
	2059303-005	1 Small Adult, 1 Adult, 1 Large Adult, 1 RADIAL-CUF (4 cuffs/pack)
2-Tube Submin	2056304-002	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)
	2298	1 Infant, 1 Child, 2 Small Adult, 2 Adult, 1 Adult Long, 2 Large Adult, 1 Thigh (10 cuffs/box)
	2059303-007	Small Adult, 1 Adult, 1 Large Adult, 1 RADIAL-CUF (4 cuffs/pack)

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com

1. Infant BP Cuff for Datex-Ohmeda™ monitors has Male Slip Luer Connector.

© 2022 General Electric Company.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, CRITIKON, DURA-CUF, DINACLICK, and Datex-Ohmeda are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB00418XX 5/22

